

#### **Installation Guide**

VHM Arm with 8" Extension & Swivel Stop for Mounting Philips IntelliVue MP60/70, MX400/450/500/550/600/700/800 on Dräger Apollo Left Side



The purpose of this guide is to describe mounting of Arm on left side of machine.



**WARNING**: USE OF MOUNTING HARDWARE AND MONITORING COMPONENTS OTHER THAN THOSE DESCRIBED IN THIS DOCUMENT MAY RESULT IN SERIOUS INJURY DUE TO TIPPING OF THE ANESTHESIA MACHINE.

**WARNING**: ALL DEVICES AND MOUNTING EQUIPMENT SHOULD BE POSITIONED INWARD, CLOSE TO THE APOLLO, DURING TRANSPORT.

Installer: When installation is completed, provide these instructions and all related documentation to the end user for future reference.

#### **Parts Reference**

The following parts and hardware are included with this installation kit (hardware not shown):

Item #	Description	Qty
1	VHM Arm with 8" Extension	1
2	Nut Plate	2
3	M6 x 25mm Flat Head Socket Cap Screw (FHSCS)	2
4	Locking Lever, Height Adjustment	1
5	Cable Guide	1
6	Down Post, 6"	1
7	#10-32 x 3/8" Socket Head Cap Screw (SHCS)	3
8	M6 x 8mm FHMS	2
9	M6 x 12mm FHMS	1
10	5/32" Hex Wrench	1
11	1/2" Socket Wrench	1

VHM Arm with 8" Extension

Height Adjustment Locking Lever

With 8" Extension

Table Mount (Supplied by Philips)


\*\*Cable Guide


Tools Required: 5/32" hex wrench (provided), Phillips Screwdriver 1/2" socket wrench (provided).

<sup>\*\*</sup> See DU-WS-0001 page 6 for Cable Management

## Attaching Down Post to VHM Arm

- 1. Remove plastic bolt cap from Swivel Cup (not shown).
- 2. Using the 5/32" hex wrench provided, fasten Post to Swivel Cup with three (3) #10-32 x 3/8" SHCS as shown below.






# **Mounting the Arm in the Accessory Track**

- 1. Thread one (1) M6 x 25mm FHSCS a few turns into each Nut Plate (right).
- M6 x 25mm FHSCS (2)

  2. Insert Nut Plate/screw assemblies in top and bottom slots of Channel Slide (rear of Arm), holding bottom screw in place to prevent it from falling from the slot.

  Nut Plate (2) -
  - IL


- 3. Hold Nut Plates in slots (below left) and guide entire assembly to top of track.
- 4. Using the 5/32" hex wrench provided, tighten two (2) flat head socket cap screws (below right).



WARNING: Do not attempt to remove the Arm from the track while an instrument is mounted.


While \*supporting the weight of the mounted instrument, loosen the two (2) flat head socket cap screws (FHSCS) in the Slide. A 5/32" hex wrench is provided for this task. Push up or pull down at the end of the arm nearest the Slide. Tighten screws to secure Arm in track.

\* It may be necessary to dismount instrument before adjusting Arm in track.



### Attaching the Philips Table Mount to the Mounting Adapter


1. Fasten the Philips-supplied Table Mount to the Mounting Adapter on the VHM Arm with one (1) M6 x 12mm FHMS and two (2) M6 x 8mm FHMS as shown below.



### **Mounting the Monitor on the Arm (Table Mount)**

- 1. Thread the Height Adjustment Locking Lever into the mounting hole in the right side of the Arm (right).
- 2. Mount monitor according to Philips table mount instructions.
- Perform counterbalance adjustment in accordance with Section 3.0 of "Installation & Operation Manual for VHM™ Series Arms" (DU-WS-0001).

  Height Adjustment



#### Mounting the MX400/450/500/550 Monitor (Fixed Mount Adapter)

- 1. Thread the Height Adjustment Locking Lever into the mounting hole in the right side of the Arm (right).
- **2.** Mount the monitor in accordance with separate Installation Guide DU-PH-0073-15 provided with the Mounting Adapter, ordered separately.
- **3.** Perform counterbalance adjustment in accordance with Section 3.0 of "Installation & Operation Manual for VHM™ Series Arms" (DU-WS-0001).

#### VHM Arm Operation and Adjustment

Refer to "Installation & Operation Manual for VHM™ Series Arms" (DU-WS-0001).

### **Routine Maintenance**

Periodically check all mounting hardware. Tighten as necessary for optimal operation.

### **Cleaning the Mounting Assembly**

**CAUTION:** GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted instruments or accessory equipment, refer to the specific instructions delivered with those products.

- 1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
- 2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.
- 3. Do not use steel wool or other abrasive material to clean the mounting assembly.
- 4. Damage caused by the use of unapproved substances or processes will not be covered by warranty. We recommended that you test any cleaning solution on a small area of the arm that is not visible to verify compatibility.
- **5.** Never submerge or allow liquids to enter the arm. Wipe any cleaning agents off of the arm immediately, using a water-dampened cloth. Dry the arm thoroughly after cleaning.