PHILIPS

Foley Catheter Temperature Probe With 400 Series Thermistor temperature probe

Application Note

Description

The Philips Medical Systems Foley Catheter temperature probe utilizes a standard Foley Catheter with a secondary lumen in which a temperature sensor is permanently placed. The sensing end of the temperature sensor is near the tip of the catheter. The sensor's lead wires are terminated by an electrical connector which facilitates interconnection with the temperature monitoring system. An integral part of the catheter is the retention balloon which is inflated through a self sealing valve housed inside the luer lock adapter. The catheter's tip is a smooth bullet form configuration to minimize trauma during insertion. The assembly is sterile and available in sizes 14 French, 16 French, and 18 French.

Leakage current of the temperature probe is of a magnitude well within clinically acceptable non-hazardous levels. The accuracy of the thermistor is ± 0.2 °C in the range of 0 - 60 °C.

Indications

The Foley Catheter temperature probe is intended for temperature monitoring in patients which are catheterized because of fluid management problems, evaluation of urinary output, urine clearance following surgery or trauma involving pelvic organs, or obstruction and paralysis.

The Philips Foley Catheter temperature probe is designed for insertion through the urethra. It is compatible with all Philips Temperature Monitoring Modules via extension cables having either 1/4 inch phone jack or two pin connectors.

Directions for Use and Precautions

1. Use currently accepted aseptic procedures for the catheter care and patient preparation.
2. Lubricate the catheter with recommended lubricant and insert it using the currently accepted medical techniques. Do not use petroleum based lubricants.
3. Do not inflate catheter before verifying proper drainage.
4. To inflate balloon proceed as follows: Fill syringe with 5cc of sterile water. Insert luer end of syringe into the luer housing of the self sealing valve. Push the plunger slowly at a steady rate. Remove the syringe and attach the catheter's open end to the collection bag line. Do not over inflate the balloon. The balloon volume is 5cc with maximum 10cc over inflation.
5. Secure the catheter and drainage tubing using currently accepted procedures. Do not kink or squeeze the tubing.

6. Align the sensor's thermistor connector with the monitor cable connector and push firmly to assure a full contact. Forced mating of the connectors can interrupt electrical continuity.
7. To prevent accidental dislodgement of the probe, attach the cable clip to the bedding.
8. Three minutes of equilibration time is necessary prior to the initiation of monitoring.
9. Prior to removing the catheter, deflate the balloon by using an empty syringe to withdraw the water completely through the luer of the self releasing valve.
10. Remove the catheter using currently accepted medical procedure, disconnect it from the cable and discard it.

Electromagnetic Compatibility

The electromagnetic compatibility (EMC) validation of the 21096A, 21097A and M2255A included testing performed according to the international standard for EMC with medical devices. See the Manufacturer's Declaration for details.

EMC Testing

During the test program the 21096A, 21097A and M2255A used with Philips Patient Monitors were subjected to many EMC tests, both international standard and Philips proprietary tests. During most of the testing no anomalies were observed. For the tests, IEC 801-3 Radiated Immunity some reduced performance was observed. IEC 801-3 specifies that the product be subjected to a field of 3V/m over a frequency range of 26 to 1000 MHz with no degradation of performance. At many of the test frequencies over the specified range no anomalies were observed. However, at a number of test points, a change in the temperature display on the Philips Patient Monitor was observed. For the test points the radiated field was reduced to the level at which the trace returned to normal. These reduced levels are 0.96 - 3V/m from 26-1000 MHz.

Avoiding Electromagnetic Interference

When electromagnetic interference (EMI) is encountered there are a number of things that can be done to mitigate the problem.

- Eliminate the source. Possible sources of EMI can be turned off or moved away to reduce their strength.
- Attenuate the coupling. If coupling is through accessory leads, the interference may be reduced by moving and/or rearranging the leads.
- Add external attenuators. If EMI becomes an unusually difficult problem, external devices may be of help. A Philips Customer Engineer can be of help in determining the need for external devices.

Adverse Reactions

Reported adverse reactions associated with the use usage of the Foley Catheter with or without temperature sensor are: Septicemia, Urethritis, Urethrocutaneous Fistulae, Urinary tract infection and encrustations, bladder spasms, tissue burns due to aberrant electro-cautery RF current pathways.

Warnings

Do not over inflate the balloon. Balloon volume is 5cc, with 10cc maximum over inflation of water. Deflate balloon completely prior to catheter removal. Do not use petroleum based lubricants. During surgical dissection which employs electro-cautery, use currently accepted procedures to minimize conditions of the thermistor and lead wires functioning as an alternate path for radio-frequency current to return to ground, causing localized tissue burns. Procedures which may minimize risk of electro-cautery are:

Keep both active and ground electrodes of the electro-cautery system in close proximity so that the temperature sensor is outside of the radio-frequency current field.

Keep temperature monitor with its associated cables separated from the electro-cautery system.

Unusual, fast artificial variations in temperature readings may occur with concomitant applications of the electro-cautery system.

Contraindications

There are no known contraindications.

Sterile-Disposable Single Patient Use Only

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<http://shop.medical.philips.com>

Asia Pacific

Philips Medical Systems
24F Cityplaza One
1111 King's Road
Taikoo Shing, Hong Kong
(+852) 3197 7777

Canada

Philips Medical Systems
2660 Matheson Blvd. E
Mississauga, ON L4W 5M2
(800) 291-6743

United States

Philips Medical Systems
3000 Minuteman Road
Andover, MA 01810
(800) 225-0230

Latin America

Philips Medical Systems
5200 Blue Lagoon Drive
9th Floor
Miami, FL 33126
(305) 267-4220

Europe, Middle East and Africa

Philips Medizinsysteme Boeblingen GmbH
Medical Consumables
Hewlett-Packard Str. 2
71034 Boeblingen
Germany
Fax: (+49) 7031 463 1552

Austria:

(+43) 1 25125-600
office_medical@philips.com

Belgium/Luxembourg

(+32) 2 404-9105

Finland

(+358) 10 855 2450

France

(+33) 0825 01 07 16

Italy

(+39) 02 9260 8265

Medical Distribution

Europe, Middle East, Africa
(+41) 22 780-6162

Netherlands

(+31) 20 547-2510

Poland

(+48) 22 549 1400

Russia

(+7) 095 797 3973
F (+7) 095 797 3903

Spain

(+34) 91-631.31.00
F (+34) 90 230 40 50

Sweden

(+46) 8 506 488 12

Switzerland

(+41) 735-9111 (German Swiss)
(+41) 22 780-6111 (Suisse Romande)

United Kingdom

(+44) 870 607 76 77

© Copyright 2000 Philips Medical Systems
Printed in U.S.A. 11/01
AW-5980-1275-2 Rev. B
5980-1275E

21096 application note: Foley Catheter Temperature Probe w/400 series thermistor

PHILIPS PROPRIETARY:

THIS DOCUMENT AND/OR ELECTRONIC FILE IS THE PROPERTY OF
PHILIPS MEDICAL SYSTEMS AND CONTAINS CONFIDENTIAL AND PROPRIETARY
INFORMATION OF PHILIPS. NEITHER IT NOR THE INFORMATION CONTAINED
THEREIN SHALL BE DISCLOSED TO OTHERS OR DUPLICATED OR USED BY
OTHERS EXCEPT AS AUTHORIZED BY PHILIPS MEDICAL SYSTEMS IN WRITING.

PART NO.: **5980-1275E** REV: **B**
DRAWING NO.: **AW-5980-1275-2 Rev. B**
REV.HISTORY: **Rebranding from Agilent**
REF.SPEC.DRG.: **to Philips.**
NAME OF ORIG.: **Lapham/Miller Associates**
C.O.: