PHILIPS

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Skin Temperature Probe 2109IA

With 400 Series Thermistor temperature probe

Application Note

Description

Monitoring of skin surface temperature is a well accepted procedure for the detection of hypothermic or hyperthermic conditions. The Philips Medical Systems skin temperature probe is designed for placement on the surface of the skin and consists of a thermistor, plastic film, and foam covered with a metallic layer. The surface of the plastic is coated with medical grade hypoallergenic adhesive suitable for skin application. The hypoallergenic adhesive holds the sensor in situ. The foam thermally insulates the sensor and the metallic layer reflects external infrared heat energy, thus providing accurate measurement of body temperature. Leakage currents of the device is of magnitude well within clinically acceptable non-hazardous levels. The device is disposable, single patient use only and it is packaged sterile.

Range and Accuracy

The accuracy of the thermister is $\pm 0.1^{\circ}\text{C}$ from 25°C to 45°C .

Indications

Philips skin temperature probe is intended to be used with Philips Temperature Monitoring Modules via extension cables. Skin probe is specially designed for monitoring neonates and infants.

Directions for Use and Precautions

1. Dry skin completely in area of intended sensor placement.
2. Remove protective paper backing and place the sensor over previously dried skin area.
3. Align the sensor's connector with the monitor cable's connector and push firmly to assure a full contact. Forced mating of the connectors without proper alignment may cause damage to the connectors and interruption in electrical continuity.
4. During surgical procedures which employ electro-cautery, use currently acceptable procedures to minimize conditions of the thermistor and a lead wire functioning as an alternate path for radio-frequency current to return to ground, causing localized tissue burns.

Here are several procedures which may minimize risk of electro-surgical burns:

Keep both active and ground electrodes of the electro-cautery system in close proximity so that the skin temperature sensor is outside of the radio-frequency current field.

Keep the temperature monitor with its associated cables separated from electrocautery systems.

Unusual, fast variations in temperature reading may occur with concomitant applications of electro-cautery systems.

5. Reuse may cause degradation of adhesive or electrical properties. Do not reuse on another patient due to risk of cross-infection.

Disposal

Follow local regulations for disposal.

Sterile-Disposable

0123

	Do not reuse
	Do not use if seal is broken
	Caution see documentation
	STERILE R Sterilized via irradiation

Electromagnetic Compatibility

The electromagnetic compatibility (EMC) validation of the 2109IA included testing performed according to the international standard for EMC with medical devices. See the Manufacturer's Declaration for details.

Avoiding Electromagnetic Interference

If electromagnetic interference (EMI) is encountered there are a number of things that can be done to mitigate the problem.

- Eliminate the source. Possible sources of EMI can be turned off or moved away to reduce their strength.
- Attenuate the coupling. If coupling is through accessory leads, the interference may be reduced by moving and/or rearranging the leads.
- Add external attenuators. If EMI becomes an unusually difficult problem, external devices may be of help. A Philips Customer Engineer can be of help in determining the need for external devices.

Adverse Reactions

Adverse reactions reported during applications of skin temperature probes include skin abrasion reaction and tissue burns due to aberrant electro-cautery RF current pathways.

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Copyright © 2012
Böblingen GmbH