

Installation Guide


Dual Flat Panel M-Series Arm Wall Mount Kit for Philips MP90

The purpose of this guide is to describe installation of the Dual Flat Panel M-Series Arm Wall Mount Kit.

Parts Reference

The following parts and hardware are included with this installation kit (see photos for parts; hardware not shown)

Item #	Description	Qty
1	M-Series Arm (12" Arm with 6" Riser)	1
2	Dual Display Bar (includes hardware & installation guide)	1
3	Keyboard Angle Mount (includes hardware & installation guide)	1
4	Keyboard Plate (includes hardware & installation guide)	1
5	Down Post Kit, 12" (includes hardware & installation guide)	1
6	Display/Remote Speedpoint Mounting Kit (includes hardware & installation guide)	1
7	Handle (includes hardware & installation guide)	1




12" Down Post




12" M-Series Arm with 6" Riser

Keyboard Plate



Keyboard Angle Mount




Dual Display Bar

Tools Required: Phillips screwdriver (not provided), 5/32" and 1/8" hex wrenches (provided).



Installing the Dual Display Wall Mount Kit

Installation Note: Each component of this installation kit (see page 1) is packaged separately with its own mounting hardware and installation guide. Please refer to the separate installation guides as you install the mounting assembly in the numerical order shown below.



Mounting FMS on Down Post

1. Mount FMS on Down Post by placing the Clamp around the Post and tightening the Knob.



Mounting Display, Remote Speedpoint, External Alert Device, and Power Supply

Please refer to separate installation guide: *Philips MP90 Remote Speedpoint Mounting Bracket and Adapter, and External Alert Device Mounting Bracket (DU-AG-0019-80)*. **Installation Note:** Flat Panel Display mounting is also covered in referenced installation guide. Mount Display according to Section 2.0 for a 15" Display, or Section 3.0 for a 17" Display.

M-Series Arm Operation and Adjustment

See Installation Guide provided with the M-Series Arm (*DU-WMM-0001*).

Routine Maintenance

Periodically check all mounting hardware. Tighten as necessary for optimal operation.

Cleaning the Mounting System

CAUTION: GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted instruments or accessory equipment, refer to the specific instructions delivered with those products.

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.
3. Do not use steel wool or other abrasive material to clean the mounting assembly.
4. Damage caused by the use of unapproved substances or processes will not be covered by warranty. We recommended that you test any cleaning solution on a small area of the arm that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the arm. Wipe any cleaning agents off of the arm immediately, using a water-dampened cloth. Dry the arm thoroughly after cleaning.