Using the Twist-Push Handle Fetal Double-Spiral Electrode (FDSE)

Insert

FDSE guide tube and gently hold it flat against fetal presenting part while performing steps below



Twist

handle 45-degrees CLOCKWISE to unlock Protection Tab



Push

handle forward, sliding it over the outer guide tube until you feel spiral tips contact the presenting part




Turn

handle CLOCKWISE until mild
resistance indicates full attachment
WARNING: NEVER TURN MORE
THAN 1/2-TURN AFTER CONTACTING
PRESENTING PART



Pull

outer guide tube back to slide both applicator tubes off wires



Connect

FDSE connector to Adapter Cable - refer to *Using the* DECG Legplate Adapter Cable on this card



FDSE Connector

Copyright © 2006 Koninklijke Philips Electronics N.V. All rights reserved Printed in USA December 2006 453564053781



Objectives Fetal Double-Spiral Electrode (FDSE)

Instructor Notes

Use the following objectives to introduce clinicians to the Fetal Double-Spiral Electrode (FDSE):

- Explain how to set up the FDSE, attachment electrode, and legplate adapter cable.
- Refer to the *Instructions for Use* shipped with the FDSE and legplate adapter cable for set-up details.
- Explain that yellow color-coded components indicate prouct compatibility.
- Explain how to insert the FDSE.
- Demonstrate how to attach the FDSE using the Philips simulator (or a similar tool that represents the fetal presenting part).
- Demonstrate safe removal of FDSE from the fetal simulator: grasp the wires as close as possible to the fetal presenting part, turning counter-clockwise until spiral tips come free. WARNING: Do not pull the spiral tips out of the fetal skin, and do not pull the wires apart.

Participant Learning Objectives

After completing this training, you should be able to:

- Identify FDSE product safety features.
- Demonstrate how to use the FDSE twist-push handle.
- Demonstrate how to *attach* and *detach* the FDSE using the fetal simulator.
- Explain how to connect the FDSE, legplate adapter cable, and attachment electrode.


2. 6, 4, 2, 1, 5, 3 3. A

Answers to Check Your Knowledge:


1. a. guide tube b. spiral tips c. protection tab d. twist-push

Check Your Knowledge Fetal Double-Spiral Electrode (FDSE)

1. Fill in the blanks using the terms <i>spiral tips, twist-push, protection tab, guide tube.</i>
a. Removed by sliding it over the lead wires:
b. Inserted into the presenting part:
c. Secures the spiral tips in their protected
position until ready for attachment:
d. Describes how to use the FDSE handle:
2. Number the following steps in the order in which yo should apply the FDSE:
STEP DESCRIPTION Connect the electrode to the legplate adapter cable
Turn handle clockwise until mild resistance indicates full attachment
Twist handle 45-degrees clockwise to unlock protection tab
Insert guide tube and hold it flat against the presenting part
Pull guide tube back to slide both application tubes off the electrical wires
Push handle forward until you feel spiral tips contact fetus
3. Which illustration is correct?
A B



DECG Reusable Legplate Adapter Cable





9898 031 37651 Philips, Agilent, HP



9898 031 37791 GE Corometrics 115, 116, 118, 119, 173, 174 Adapters 1442AAO, 1442BAO



9898 031 37801 GE Corometrics 120, 126, 128, 129



9898 031 37811 Spacelabs IM77

Using the DECG Reusable Legplate Adapter Cable

FDSE Connector

Check

that the cleaning port cover is closed

Connect

legplate adapter cable to fetal module/monitor

Snap

legplate adapter cable onto attachment electrode

Press

attachment electrode onto the mother's leg

Verify

that yellow FDSE connector is dry, wipe if needed

Plug

FDSE connector into legplate adapter cable

After Use

If needed, use the cleaning port to flush and clean the legplate adapter cable (refer to DECG Reusable Legplate Adapter Cable Instructions for Use)