Philips TRANSPAC® 4 Disposable **Transducer Pressure Monitoring Kits**

Introduction to Steps

Transducer Set Up



 Slide the transducer into the mounting bracket.





- Activate the flush device by compressing the tags of the Squeeze flush device.
- Flush until fluid exits the side port of the zeroing stopcock.



Connect the cable extension of the TRANSPAC[®] 4 transducer to the reusable cable.

· Connect the other end of the reusable cable to the patient monitor.

6



• Turn the zeroing stopcock handle "off" to the zeroing port. Cover with a yellow non-vented cap.



· Ensure that all luer connections are tightened.

7



Activate the flush device to flush the remaining patient line. If additional stopcocks are in line, ensure that fluid exits each side port.



· Connect the administration set to the fluid source.

· Gravity prime recommended.

8



- · Verify that all stopcock handles are turned "off" to the side ports and that side ports are covered by yellow non-vented caps.
- Pressurize the IV fluid source to 300 mm Hg. Close the pressure cuff clamp.

Leveling and Zeroing the Transducer



· Position the pole mount so that the zeroing port is level with the patient's phlebostatic axis.



- Turn the zeroing stopcock handle "off" to the patient.
- Remove the yellow non-vented cap from the zeroing port.
- · Zero balance and calibrate per the monitor manufacturer's instructions.
- When completed, turn the zeroing stopcock handle "off" to the zeroing port.
- Cover the zeroing port with a yellow non-vented cap.



· Activate the flush device to flush the remaining patient line. If additional stopcocks are in line, ensure that fluid exits each side port.

