PHILIPS

Esophageal Stethoscope

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Esophageal Stethoscope Temperature Probe 21093A, 21094A, 21095A With 400 Series Thermistor

Application Note

Description

The Philips Medical Systems esophageal stethoscope with temperature probe is a disposable device which provides for accurate measurement of core body temperature as well as transmission of heart and lung sounds. Heart and lung sounds are transmitted across a thin cuff wall, through the side and distal openings of the tube. Made of a special material which optimizes sound transmission and remains both durable and flexible, the cuff maintains its integrity and prevents secretion from entering the tube. For connection to any standard ear piece, a male luer adapter is attached to the proximal end of the esophageal stethoscope. The stethoscope is 19 in. (48.3 cm) long. The tube, the cuff and luer adapter are made of light weight, nontoxic, implant tested, pyrogen free material.

The temperature transducer is in the form of an electrically insulated thermistor; which is permanently secured within the lumen of the esophageal stethoscope. The thermistor itself is placed toward the distal end of the tube but not outside the distal end. Leakage current of the Philips temperature probe (thermistor and stethoscope) is of a magnitude well within accepted nonhazardous levels.

Range and Accuracy

The accuracy of the thermister is $\pm 0.1^{\circ}\text{C}$ from 25°C to 45°C .

Indications

The Philips esophageal stethoscope temperature probe is suitable for usage with Philips Temperature Monitoring Modules via extension cables.

Directions for Use and Precautions

1. If a patient has to be intubated with an E.T. tube, perform intubation prior to placing the esophageal stethoscope.
2. Lubricate the stethoscope prior to insertion, and place the probe in accordance with currently acceptable medical procedures.
3. Verify position of the stethoscope by direct laryngoscopy or other medical techniques.

4. Align the stethoscope's connector with cable connector and push firmly to assure a full contact. Forced mating of the connectors without proper alignment may cause damage to the connectors and interruption in electrical continuity.

5. During surgical procedures which employ electro-cautery, use currently acceptable procedures to minimize the situation in which the thermistor and a lead wire may function as an alternate path for radiofrequency current to return to ground, which may cause localized tissue burns.

Here are several procedures which may minimize risk of electro-surgical burns:

Keep both active and ground electrodes of electro-cautery system in close proximity so that the esophageal stethoscope is outside of the radiofrequency current field. Keep the temperature monitor with associated cables separated from electro-cautery systems.

Unusual, fast artificial variations in temperature reading may occur with concomitant applications of electrocautery systems.

6. Reuse may cause degradation of physical or electrical properties. Do not reuse on another patient due to risk of cross-infection.

Electromagnetic Compatibility

The electromagnetic compatibility (EMC) validation of the 21093A, 21094A and 21095A included testing performed according to the international standard for EMC with medical devices. See the Manufacturer's Declaration for details.

Avoiding Electromagnetic Interference

If electromagnetic interference (EMI) is encountered there are a number of things that can be done to mitigate the problem.

- Eliminate the source. Possible sources of EMI can be turned off or moved away to reduce their strength.
- Attenuate the coupling. If coupling is through accessory leads, the interference may be reduced by moving and/or rearranging the leads.
- Add external attenuators. If EMI becomes an unusually difficult problem, external devices may be of help. A Philips Customer Engineer can be of help in determining the need for external devices.

Adverse Reactions

Adverse reactions reported during applications of esophageal stethoscope with or without temperature probes include accidental tracheal; or bronchial intubation accompanied with airway obstruction, esophageal abrasion, and tissue burns due to aberrant electro-cautery RF current pathways.

Contraindications

The use of the esophageal stethoscope may be contraindicated in neonates and small infants undergoing tracheostomy or internal jugular vein catheterization.

Disposal

Follow local regulations for disposal.

Sterile-Disposable

CE 0123

	Do not reuse		Do not use if seal is broken
	Caution see documentation		STERILE Sterilized via irradiation

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Copyright © 2012
Böblingen GmbH