

3875 Cypress Drive Petaluma, CA 94954 800.228.2555 707.773.1100 Fax 707.773.1180 www.gcx.com INSTRUMENT MOUNTING SYSTEM

# Installation Guide

VHM Arm with 8" Extension & Swivel Stop for Mounting Philips MP60/70 on Dräger Apollo Left Side



The purpose of this guide is to describe mounting of Arm on left side of machine.



**WARNING**: USE OF MOUNTING HARDWARE AND MONITORING COMPONENTS OTHER THAN THOSE DESCRIBED IN THIS DOCUMENT MAY RESULT IN SERIOUS INJURY DUE TO TIPPING OF THE ANESTHESIA MACHINE.


**WARNING**: ALL DEVICES AND MOUNTING EQUIPMENT SHOULD BE POSITIONED INWARD, CLOSE TO THE APOLLO, DURING TRANSPORT.

### Installer: <u>When installation is completed, provide these instructions and all related documentation to the end</u> <u>user for future reference.</u>

### **Parts Reference**

The parts list below includes parts and hardware that will be used in this installation procedure (see photos for parts; hardware not shown).

Item #	Description	Qty
1	VHM Arm with 8" Extension	1
2	Nut Plate	2
3	M6 x 25mm Flat Head Socket Cap Screw (FHSCS)	2
4	Locking Lever, Height Adjustment	1
5	Cable Clip	2
6	M6 x 8mm FHMS	2
7	M6 x 12mm FHMS	1
8	5/32" Hex Wrench	1
9	1/2" Socket Wrench	1




Tools Required: 5/32" hex wrench (provided), 1/2" socket wrench (provided).

Table Top Mount (Supplied by Philips)

# Mounting the Arm in the Accessory Track

- 1. Thread one (1) M6 x 25mm FHSCS a few turns into each Nut Plate (right).
- 2. Insert Nut Plate/screw assemblies in top and bottom slots of Channel Slide (rear of Arm), holding the bottom screw in place to prevent it from falling from the slot.
- 3. Hold Nut Plates in slots (below left) and guide entire assembly to top of track.
- 4. Using the 5/32" hex wrench provided, tighten two (2) flat head socket cap screws (below right).




## **Positioning Arm in Track**

**WARNING**: Do not attempt to remove the Arm from the track while an instrument is mounted.


While \*supporting the weight of the mounted instrument, loosen the two (2) flat head socket cap screws (FHSCS) in the Slide. A 5/32" hex wrench is provided for this task. Push up or pull down at the end of the arm nearest the Slide. Tighten screws to secure Arm in track.

\* It may be necessary to dismount instrument before adjusting Arm in track.



# Attaching the Philips Table Top Mount to the Mounting Adapter

1. Fasten the Philips-supplied table top mount to the Mounting Adapter on the VHM Arm with one (1) M6 x 12mm FHMS and two (2) M6 x 8mm FHMS as shown below.



### Mounting the Monitor on the Arm

- 1. Thread the Height Adjustment Locking Lever into the mounting hole in the right side of the Arm (right).
- 2. Mount monitor according to Philips table top mount instructions.
- 3. Perform counterbalance adjustment in accordance with Section 3.0 of "Installation & Operation Manual for VHM™ Series Arms" (DU-WS-0001).



### VHM Arm Operation and Adjustment

Refer to "Installation & Operation Manual for VHM™ Series Arms" (DU-WS-0001).

### **Routine Maintenance**

Periodically check all mounting hardware. Tighten as necessary for optimal operation.

### **Cleaning the Mounting Assembly**

**CAUTION:** GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted instruments or accessory equipment, refer to the specific instructions delivered with those products.

- 1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
- 2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.
- 3. Do not use steel wool or other abrasive material to clean the mounting assembly.
- 4. Damage caused by the use of unapproved substances or processes will not be covered by warranty. We recommended that you test any cleaning solution on a small area of the arm that is not visible to verify compatibility.
- 5. Never submerge or allow liquids to enter the arm. Wipe any cleaning agents off of the arm immediately, using a waterdampened cloth. Dry the arm thoroughly after cleaning.