## OmniLab Advanced +

## Sleep lab competency checklist

	Assessed		ssessor:		Date	:
--	----------	--	----------	--	------	---

Goals		Date goals were met	Assessor
Principles of operation			
Describe the intended use of the intended	the device.		
☐ Designated patient population	ons		
☐ Describe appropriate delive			
2. Explain the purpose of each ite	em:		
☐ CPAP mode	☐ Rate		
$\square$ Spontaneous mode	$\square$ Inspiratory time		
$\square$ Spontaneous/timed mode	$\square$ Breath rate setting		
$\square$ Timed therapy mode	$\square$ Volume setting		
☐ autoSV mode	☐ IPAP <sub>max</sub> , IPAP <sub>min</sub>		
☐ PC mode	☐ Maximum pressure		
☐ Auto CPAP	☐ EPAP <sub>max</sub> , EPAP <sub>min</sub>		
☐ Auto Bi-level	☐ Min pressure support		
☐ AVAPS	☐ Max pressure support		
☐ AVAPS AE	☐ Ramp		
☐ Flex options	☐ Unintentional leak/mask type		
☐ Rise time	☐ Total leak		
	☐ System One		
	Resistance Control		
Setup and operation			
1. Identify the circuit components	of the OmniLab Advanced + system:		
☐ Heated tube			
$\square$ Non-heated tube			
$\square$ Power cord			
$\square$ White and gray filters			
☐ Humidifier			
2. Demonstrate assembly of the			
and patient circuit, including th	e indications, setup, and		
humidifier adjustment.			
3. Demonstrate the steps to access and connect the device			
to the PSG system.			
,			
4. Demonstrate how to access and	d modify OmniLab Direct.		



Goals	Date goals were met	Assessor
The clinician will be able to:		
1. Demonstrate how to utilize OmniLab Direct to modify settings		
specific to each mode of therapy.		
2. Describe how to navigate and interpret the following settings		
from OmniLab Direct:		
☐ Total leak ☐ Respiratory rate		
☐ Unintentional leak ☐ Exhaled tidal volume		
☐ Minute ventiliation		
3. Describe the system alerts and actions to resolve.		
5. Describe the system are is and actions to resolve.		
4. Describe the method and issues involved with adding		
supplemental oxygen.		
Alert conditions and troubleshooting		
1. Demonstrate understanding of device alarms/alerts, including		
possible causes, and corrective action.		
Specific alarms/alerts:		
☐ Apnea ☐ Patient disconnect		
☐ Low minute ventilation ☐ Low tidal volume		
☐ Reset therapy meter		
Cleaning and care		
1. Demonstrate the method of cleaning the device and		
replacement of the filters.		
Assessed signature: As	ssessor signature:	
Assessed signature. As	osessor signature	

autoSV,AVAPS, OmniLab and Respironics are trademarks of Philips Respironics and its affiliates. All rights reserved.

## Please visit www.philips.com/omnilab



©2013 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

 $\label{eq:caution:cauchy} \textbf{CAUTION: US federal law restricts these devices to sale by or on the order of a physician.}$ 

edoc SB 9/30/13 MCI 4105864

Philips Healthcare is part of Royal Philips

Respironics Asia Pacific +65 6298 1088 Respironics Australia +61 (2) 9666 4444 Respironics Europe, Middle East, Africa +33 1 47 52 30 00 Respironics United Kingdom +44 800 1300 845 Philips Respironics 1010 Murry Ridge Lane Murrysville, PA 15668

Customer Service +1 724 387 4000 +1 800 345 6443 (toll free, US only) www.philips.com/respironics