Instructions for Use

BiPAP Disposable Circuit

- ⚠ Caution:** U. S. federal law restricts this device to sale by or on the order of a physician.
- 🚫 Note:** This circuit does not contain natural rubber latex or DEHP.

Symbols

Warning or Caution

Note

Do Not Reuse

⚠ Warnings:

- The exhaust vent on the Exhalation Port is designed to exhaust CO₂ from the patient circuit. Continuous flow is required for safe operation. Do not block or otherwise try to seal the exhaust vent on the exhalation port.
- Before instituting mechanical ventilation, confirm the proper connection of the circuit, and verify ventilator operation. Patient, ventilator and circuit must be monitored on a regular basis per established hospital standards of care.
- Do not pull or stretch the circuit tubing. This could result in circuit leakage.
- This circuit is not recommended for use with a heated humidifier since it does not allow for the use of water traps.
- This circuit is for single use only and is not intended to be disinfected or cleaned after patient use. Dispose of the circuit after patient use. Multi-patient use without disinfection between patients increases the risk of infection.
- This circuit is not designed for reuse. Respironics cannot guarantee the performance specifications once this circuit has been disinfected or cleaned after patient use.

Assembly Instructions


The BiPAP Disposable Circuit is made up of Circuit Tubing and a Proximal Pressure Line Assembly.

⚠ Caution: The length of time a circuit is used should be determined by the institution's accepted infection control procedures.

1. Connect the circuit tubing (Figure 1).
 - Connect a low-resistance, mainflow bacteria filter on the outlet port of the BiPAP unit.
 - Connect the cuffed end (Figure 1-B) of the 6' Tubing to the outlet port of the filter.
2. Connect the Proximal Pressure Line Assembly.

🔍 Figure 1

- A = Tee Connector
 B = Tubing Cuff
 C = Proximal Pressure Line Assembly
 D = Hose Clip
 E = 6' Circuit Tubing
 F = Hose Hanger
 G = Hose Clip
 H = Exhalation Port Exhaust Vent


- BiPAP S-D or S/T-D Units with Airway Pressure Monitor (APM) (Figure 2)
 - a. Connect the 24" length of line to the Pressure Port on the front of the BiPAP unit.
 - b. Connect the 2" length of line to the Airway Connection Port on the rear of the APM.

- BiPAP S-D or S/T-D Units without Airway Pressure Monitor (APM) (Figure 3)
 - a. When the Airway Pressure Monitor is not used, only the 84" length of proximal pressure line tubing is needed. Disconnect the 2" and 24" length from the assembly by cutting the 84" tubing near the Tee Connector. Discard the Tee Connector, 2" length and 24" length of tubing.
 - b. Connect the 84" length to the Pressure Port on the front of the BiPAP unit.


- BiPAP S or S/T Units with Airway Pressure Monitor (APM) (Figure 4)
 - a. When a non-diagnostic unit is used (BiPAP S unit or BiPAP S/T unit), only the 84" length of proximal pressure line tubing is needed. Disconnect the 2" and 24" length from the assembly by cutting the 84" tubing near the Tee Connector. Discard the Tee Connector, 2" length and 24" length of tubing.
 - b. Connect the 84" length of line to the Airway Connection Port on the rear of the APM.

3. Connect the Exhalation Port End of the Tubing to the Mask Rotating Connector.
 - The vent on the Exhalation Port should be positioned away from the patient's face.


▲ Figure 2

- A = Pressure Port
- B = BiPAP Outlet Port
- C = Mainflow Bacteria Filter
- D = Hose Hanger
- E = Exhalation Port (Do not block)
- F = Rotating Connector
- G = Proximal Pressure Line Assembly
- H = 6' Circuit Tubing
- I = Hose Clips

▲ Figure 3

- A = Pressure Port
- B = BiPAP Outlet Port
- C = Mainflow Bacteria Filter
- D = Hose Hanger
- E = Exhalation Port (Do not block)
- F = Rotating Connector
- G = Proximal Pressure Line Assembly
- H = 6' Circuit Tubing
- I = Hose Clips

▲ Figure 4

- A = BiPAP Outlet Port
- B = Mainflow Bacteria Filter
- C = Hose Hanger
- D = Exhalation Port (Do not block)
- E = Rotating Connector
- F = Proximal Pressure Line Assembly
- G = 6' Circuit Tubing
- H = Hose Clips

Reorder Numbers

312107	BiPAP Disposable Circuit, Box of 10
312111	BiPAP Disposable Humidifier Tubing, 36", Box of 10
312112	Disposable Proximal Pressure Line Assembly, Box of 10
312121	Disposable Proximal Pressure Line, 84", Box of 10
312106	Disposable Headstrap, Box of 10
312151	Disposable Tubing, 6', Box of 10
312153	Disposable Tubing, 18", Box of 10

Specifications

Disposal

Dispose of in accordance with local regulations.

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Respironics Inc.