


**REF 9960STD-00**  
REORDER

**REF 9960LGE-00**  
REORDER

**REF 9960PED-00**  
REORDER



## CAPNO<sub>2</sub>mask

### O<sub>2</sub> Delivery/CO<sub>2</sub> Mainstream Monitoring Mask



Single  
Patient  
Use Only



NOTE: These directions are guidelines intended for use by qualified medical personnel.

- 1 Select appropriate size mask for patient.
- 2 Attach oxygen supply tubing to gas source.
- 3 Set oxygen to desired flow and verify flow through mask.
- 4 Attach CAPNOSTAT CO<sub>2</sub> sensor to airway adapter.
- 5 Place mask on patient's face - pull straps to secure mask and minimize leaks, mold metal strap to fit nose.

O <sub>2</sub> LPM	FiO <sub>2</sub> Range
6	25-45%
7, 8, 9	30-60%
10, 11, 12	40-65%
13, 14, 15	45-75%

NOTE: Delivered FiO<sub>2</sub> is dependent on patients tidal volume and respiratory rate.

**WARNING:** Always verify that the opening of the airway adapter assembly is unobstructed. There is no provision for entrainment of room air should the opening become occluded.

#### CAUTION

- Minimum oxygen flow is 6 LPM.
- Federal law (U.S.) restricts this device to sale by or on the order of a medical practitioner.
- Monitor patient ventilation during use.
- Inspect mask, ensure all components are secure and undamaged.
- Verify proper mask, tubing connection and placement during use.
- Refer to CO<sub>2</sub> monitor User's Manual prior to use.
- Verify flow.
- Clean or replace mask if occluded from patient secretions or fluids.
- Do not use mask for O<sub>2</sub> delivery without active CO<sub>2</sub> monitoring.
- Mask is not to be used for any respiratory application other than described.

NOTE: Turn on the CO<sub>2</sub> monitor prior to using the mask to minimize delay in CO<sub>2</sub> measurement.

Respironics Novamatrix, LLC  
5 Technology Drive  
Wallingford, CT USA 06492



Respironics Deutschland  
Gewerbestr. 17



82211 Herrsching  
Germany

9966-02 ver B

Signature/Date: EDMS