## Instructions for Use

### AF811 Gel Full Face Mask

#### Intended Use

The AF811 Gel Full Face Mask is intended to provide an interface for application of CPAP or bi-level therapy to patients. The mask is for single use in the hospital/institutional environment only. The mask is to be used on patients (>66lbs/30kg) for whom CPAP or bi-level therapy has been prescribed.

**Note:** An exhalation port is not built into the AF811 Gel Full Face Mask. A separate exhalation device must be used with this mask.

**Note:** This mask does not contain natural rubber latex or DEHP.

**Caution:** U. S. federal law restricts this device to sale by or on the order of a physician.

---

#### Symbols

Warning or  
Caution

Note

Consult instructions  
for use

Does not contain  
natural rubber latex

Do Not  
Reuse

Leak symbol  
and value

#### Warnings:

- This mask is not suitable for providing life support ventilation.
- This mask requires a separate exhalation device.
- This mask is designed for use with CPAP or bi-level systems recommended by your health care professional or respiratory therapist. Do not wear this mask unless the CPAP or bi-level system is turned on and operating properly. Do not block or try to seal the exhalation port. Explanation of the Warning: CPAP systems are intended to be used with special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. This warning applies to most models of CPAP systems. Rebreathing of exhaled air for longer than several minutes can, in some circumstances, lead to suffocation.
- If oxygen is used with the device, the oxygen flow must be turned off when the device is not operating. Explanation of the Warning: When the device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the device enclosure. Oxygen accumulated in the device enclosure will create a risk of fire.
- Oxygen supports combustion. Oxygen should not be used while smoking or in the presence of an open flame.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate. This warning applies to most types of CPAP and bi-level machines.
- Some users may experience skin redness, irritation, or discomfort. If this happens, discontinue use and contact your healthcare professional.
- The patient's physician should be contacted if the patient experiences the following symptoms while using the masks or after removing it: Unusual chest discomfort, shortness of breath, stomach distension, belching, or severe headache; drying of the eyes, eye pain, or eye infections; blurred vision. (Consult an ophthalmologist if symptoms persist.)
- At low CPAP or EPAP pressures the flow through the exhalation

- port may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- A minimum of 3cm H<sub>2</sub>O (hPa) must be maintained when using this mask.
  - This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask.
  - This mask is not recommended if the patient is taking a prescription drug that may cause vomiting.
  - If an additional exhalation device is added to the patient circuit, you may need to adjust the pressure level to compensate for the additional leak of the exhalation device.
  - The gel mask is not customizable. Do not attempt to heat.
  - Hand wash prior to first use. Inspect the mask for damage or wear (cracking, crazing, tears, cushion damage resulting in gel exposure, etc). Discard and replace any components as necessary.
  - This mask is for single use only and is not intended to be disinfected or cleaned after patient use. Dispose of the mask after patient use. Multi-patient use without disinfection between patients increases the risk of infection.
  - This mask is not designed for reuse. Respiromics cannot guarantee the performance specifications once this mask has been disinfected or cleaned after patient use.
  - Use of a nasal or full face mask may cause tooth, gum, or jaw soreness or aggravate an existing dental condition. Consult your physician or dentist if symptoms occur.
  - Do not block or try to seal the entrainment valve.

### Contraindications

This mask may not be suitable for persons with the following conditions: impaired cardiac sphincter function, excessive reflux, impaired cough reflex, and hiatal hernia. It should not be used if the patient is uncooperative, obtunded, unresponsive, or unable to remove the mask by themselves.

### Before Use

- Read and understand the instructions completely.
- Hand wash the mask.
- Clean the patient's face.
- Verify that the mask and headgear are the correct size.
- If a nasogastric (NG) tube or similar device is in place, use the optional NG tube sealing pad. Position the pad so that its flat surface is against the patient's face and the C-shaped opening surrounds the tube.
- Verify that the therapy device, i.e., ventilator, including the alarms and safety systems, has been validated prior to use.
- Verify that the air entrainment valve functions correctly, as noted below.
- Inspect the mask and replace it if the cushion has hardened or is torn, or if any parts are broken.
- Verify therapy device pressure(s).

### Cleaning Instructions

1. Hand wash the mask in warm water with a mild dishwashing detergent.
- ⚠ Caution:** Use a mild liquid dishwashing detergent only. Do not use bleach, alcohol, or cleaning solutions containing bleach or alcohol, or cleaning solutions containing conditioners or moisturizers.
2. Rinse thoroughly and air dry completely before use.


Figure 1

- A = Capstrap  
B = Forehead Bracket and Forehead Pad  
C = Swivel Clip (ball and socket)  
D = Mask Leak Symbol and Value  
E = Entrainment Valve with Fresh Air Inlet  
F = Faceplate with Gel Cushion  
G = Four Position Forehead Arm  
H = Headgear Strap

### Air Entrainment Valve

1. Locate the air entrainment valve flapper inside the mask elbow.
2. With the airflow turned off, the entrainment valve flapper will lie flat, so that room air can flow through the fresh air inlet in the valve.
3. With the airflow on, the flapper will now cover the fresh air inlet and air from the CPAP or bi-level device should flow into the mask.
4. Do not block or try to seal the fresh air inlet in the elbow. Ensure that the entrainment valve is not blocked with secretions and that the flapper is dry. If the flapper does not function properly, replace the mask.

### Leak Symbol and Port Settings

Some ventilators may incorporate the use of a leak symbol and value in the mask selection setup procedures. The leak characteristics of this mask is leak symbol (K1). The leak symbol and value represents the intentional leak characteristics of the interface. On ventilators equipped with a Mask Selection control, enter the leak symbol value (K1) that corresponds with the leak symbol value on the mask.

### Achieving the Right Fit

1. Hold the mask lightly against the patient's face and slide the Capstrap headgear over their head. Connect the ball-and-socket swivel clips after the headgear is in position with the headgear tabs facing outward.
2. Adjust the position of the forehead arm to ensure minimal mask contact pressure on the nasal bridge while achieving minimal mask cushion leak. The lowest position provides the highest contact pressure at the bridge of the nose and provides the best mask cushion seal. The highest position provides the lowest contact pressure at the bridge of the nose and provides lower mask cushion seal.
3. Connect the patient circuit (exhalation device and flexible tubing) to the mask elbow with air entrainment valve.
4. Unfasten the headgear tabs on the top strap of the Capstrap and gradually adjust to ensure the mask is tight enough to achieve minimal mask cushion leak and a comfortable fit. **DO NOT OVERTIGHTEN.** Unfasten the headgear tabs at the swivel clips and gradually adjust to ensure the mask is tight enough to achieve minimal mask cushion leak and a comfortable fit. **DO NOT OVERTIGHTEN.**
5. Turn on the noninvasive ventilator. Instruct the patient to breathe normally.
6. Make final adjustments to all adjustment points to manage leaks and to provide a comfortable fit.


Figure 2

### Notes:

- Re-adjust the strap tension if leaks occur as the patient's position changes.
- Do not overtighten the straps. Overtightening can cause or worsen leaks.


### Connecting and Disconnecting the Headgear

The Capstrap headgear design provides the opportunity for greater headgear positioning and comfort. This mask can also be used with the Four Point Headstrap with swivel clips (REF 1015788).

## Removing the Capstrap Headgear

First, undo the bottom swivel clips by pulling the ball out and away from the socket (Figure 3). Next, gently twist the Capstrap hook (I) so that it disconnects from the forehead support bracket.


Figure 3

## Replacing the Capstrap Headgear

Press the Capstrap hook (I) into the back of the slot on the forehead support bracket and press until they snap together (Figure 4). On the lower straps, press the balls of the swivel clips into the sockets on the faceplate.


Figure 4

## Removing the Four Point Headstrap

First, undo the bottom swivel clips by pulling the ball out and away from the socket. Next, open the headgear tabs and pull them out of the slots on the forehead support arm.


Figure 4

## Replacing the Four Point Headstrap

On the top strap, thread the headgear tabs through the slots in the forehead support arm (Figure 5). Fold the tabs back so that they catch onto the headgear straps. On the lower straps, press the balls of the swivel clips into the sockets on the faceplate.


Figure 5


## Patient Access and Removing the Mask

Disconnect the bottom strap swivel clips and swing the mask upward (Figure 6). This allows easy access to the patient. To remove completely, disconnect the bottom strap swivel clips and pull the entire assembly over the patient's head.

Figure 6


If using the Four Point Headstrap, disconnect the bottom strap swivel clips and pull the entire assembly off over the patient's head.

## Cushion and Flap Assembly

This mask is intended for single patient use and is fully assembled and ready to use. Place the gel cushion onto the faceplate (Figure 7). Then place the flap over the gel cushion and push the retaining ring onto the faceplate until the connection tabs click into place. Tug on the cushion to make sure all of the parts are securely in place.

Figure 7


## Reorder Numbers

- 1063015 Respiromics AF811 Mask w/ Capstrap Headgear, L, 1 EA
- 1063016 Respiromics AF811 Mask w/ Capstrap Headgear, M, 1 EA
- 1063047 Respiromics AF811 Mask w/ Capstrap Headgear, S, 1 EA
- 1063051 Respiromics AF811 Mask w/ Capstrap Headgear, L, 5 PK
- 1063052 Respiromics AF811 Mask w/ Capstrap Headgear, M, 5 PK
- 1063053 Respiromics AF811 Mask w/ Capstrap Headgear, S, 5 PK
- 1019547 Respiromics Capstrap Headgear Only, 5 PK
- 1015788 Respiromics Four Point Headgear w/ Swivel Clips, 1 EA
- 452040 NG Tube Sealing Pad, 10 PK

## Specifications

**Warning:** The technical specifications of the mask are provided for your healthcare professional to determine if it is compatible with your CPAP or bi-level therapy device. If used outside these specifications, or if used with incompatible devices, the mask may be uncomfortable, the seal of the mask may not be effective, optimum therapy may not be achieved, and leak, or variation in the rate of leak, may affect device function.

### Intentional Leak

### Pressure Drop cm H<sub>2</sub>O (hPa)

50 SLPM 100 SLPM

All mask sizes: 0.4 1.0

### Deadspace

Small 260 mL

Medium 280 mL

Large 300 mL

### Disposal

Dispose of in accordance with local regulations.

### Storage Conditions

Temperature: -4° F to +140° F (-20° C to +60° C)

Relative Humidity: 15% to 95%, non-condensing

### AUSTRALIAN SPONSOR DETAILS:

Philips Electronics Australia Ltd.  
65 Epping Road, North Ryde, NSW 2113  
Australia

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Respironics Inc.