



# TORQMAX<sup>®</sup> SHEATH GRIP ACCESORY

Ergonomic Grip • Easy Side-Loading • Flexible Sizing

# TORQMAX® II SHEATH GRIP ACCESSORY

## Ergonomic Grip


- Contoured polymer shell with soft over-mold construction designed for comfortable user interface
- Provides mechanical advantage to more easily rotate the body of an associated sheath
- Distributes forces along the sheath body to aid in sheath advancement

## Easy Side-Loading

- Loads quickly and easily from the side – permitting easy removal or repositioning during procedures
- Stays on the sheath where it is positioned, until the user chooses to move it

## Flexible Sizing

- One size spans a wide range of sheath outer diameters – from 11.9F to 22.5F
- Compatible with all Spectranetics Laser Sheaths and associated outer sheaths



## TORQMAX SPECIFICATIONS

Model Number	Minimum Sheath Outer Diameter (F / in / mm)	Maximum Sheath Outer Diameter (F / in / mm)	Sheath Grip Length (in / mm)
502-001	11.9 / 0.155 / 4.0	22.5 / 0.296 / 7.5	2.5 / 64

### Indications for Use

#### TorqMax Sheath Grip Accessory

The Spectranetics The TorqMax Sheath Grip Accessory is intended for use in providing ergonomic grip on outer support sheaths, dilator sheaths, and Spectranetics laser sheaths.

### Contraindications

None known.

### Warnings

Observe all warnings for the sheath to be used as indicated in the associated "Instructions for Use."

### Precautions

For Single use only. The TorqMax Sheath Grip Accessory must not be resterilized and/or reused.

Do not alter the TorqMax Sheath Grip Accessory from its original state prior to use.

Observe all precautions for the sheath to be used as indicated in the associated "Instructions for Use."



#### Corporate Headquarters

The Spectranetics Corporation  
9965 Federal Dr., Colorado Springs,  
CO 80921  
Tel: 719-447-2000 • Fax: 719-447-2022  
Customer Service: 800-231-0978

#### Spectranetics International B.V.

Plesmanstraat 6, 3833 LA Leusden  
The Netherlands  
Tel: +31 33 4347 050  
Fax: +31 33 4347 051

#### German Office

Spectranetics Deutschland GmbH  
Schweinfurter Str. 7  
97080 Würzburg  
Germany  
Tel: +49 931 4520080 • Fax: +49 931 45200811