

Fast CMR methods

for anatomy, function and more

Cardiac Expert extension is an add-on to the comprehensive Cardiac Expert option. It provides additional techniques for cardiac MR studies like fast black blood imaging and functional imaging.

Cardiac Zoom is a small FOV imaging technique that accelerates black blood TSE of the heart and great vessels. It decreases the required breath hold duration by up to 30% without changing spatial resolution by enabling single beat (shot)imaging, which is challenging for conventional (multi-beat) imaging approaches.¹


3D Non-selective delivers 3D bFFE with reduced banding artifacts compared to Philips 3D Selective 3D FFE imaging.

Retrospective EPI combines retrospective triggering with EPI sampling. kt-SENSE is a spatio-temporal acceleration technique that offers all the benefits of k-t BLAST in addition to enhanced image uniformity. 2


Cardiac Expert extension

Field strength	1.5T and 3.0T Ingenia MR systems, excluding Ingenia Prodiva ¹
Prerequisite	k-t BLAST SW option
Main applications	Cardiac, Vascular (2D/4D flow)
Sequence	Black Blood TSE, 2D and 3D FFE/TFE, 3D balanced FFE multi-dynamic or multi-phase acquisitions, Phase Contrast Angiography
Image contrast	Available for a range of contrasts, including T1w, T2w, PDw and Quantitative Flow
Speed	 Cardiac Zoom reduces scan time up to 30% with equal spatial resolution³ Retrospective EPI delivers faster scanning compared to non-EPI Q-flow
Image quality	- Cardiac Zoom reduced sensitivity to motion due to the small FOV and shorter scan time³ - 3D Non-selective reduces black banding artefacts in cardiac bFFE, enables bFFE on 3.0T⁴ - k-t SENSE enhances image uniformity due to dS SENSE reconstruction with CLEAR


1 Ingenia Prodiva 1.5T is not for sale in the USA



Non-selective 3D bFFE Resolution: 1.5 x 1.5 x 1.5 mm Scan time: 2:48 min Ingenia 1.5T



Non-selective 3D bFFE Resolution: 1.5 x 1.5 x 1.5 mm Scan time: 2:48 min Ingenia 1.5T



Cardiac Zoom T1w TSE Black Blood Resolution: 0.9 x 1.3 x 8.0 mm Scan time: 9.8 sec Ingenia 1.5T

- 3 Compared to conventional Philips black blood imaging.
- 4 Compared to Philips 3D Selective 3D bFFE imaging.

© 2019 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.



How to reach us healthcare@philips.com www.philips.com/healthcare

More information www.philips.com/mrclinicalapplications