
Register your product and get support at / 请登录以下网址并得到相应的帮助




[www.philips.com/welcome](http://www.philips.com/welcome)

**PHILIPS**

**ENGLISH 6**

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## Introduction

This booklet contains important information on the PulseRelief device. The complete user manual for this device is included in the app and can also be downloaded from **[www.philips.com/support](http://www.philips.com/support)**.

Read chapter 'Contra-indications, reactions and side effects' and chapter 'Important safety information' before you start using the device.

## Intended use

---

The PulseRelief is intended to be used by adult consumers who experience mild to moderate chronic pain or acute post-surgical pain. The purpose of this device is to provide pain relief through transcutaneous electrical nerve stimulation (TENS). You can also use it for electrical muscle stimulation (EMS). The device is suitable for home use. If you recently had surgery or are in the care of a doctor, do not use this device without consulting your doctor first. Do not place the electrodes over open wounds or rashes.

## Important safety information

### Contra-indications

---

Do not use the device if:

- you have a cardiac pacemaker, implanted defibrillator or other implanted metallic or electronic device (e.g. a drug delivery system), as use in these cases could cause electric shock, burns, electrical interference, or even death.
- you have a heart disease.
- you have epilepsy.
- you are in the first trimester of your pregnancy, as the effects of TENS on the development of the foetus are as yet unknown. Avoid use over the uterus during pregnancy to reduce the risk of inducing

labour. If you are pregnant, always consult your doctor or midwife before use.

- you have a cognitive impairment.
- you have venous or arterial thrombosis or thrombophlebitis. If you use the device for EMS, there is a risk of mobilising the blood clot.
- you have undiagnosed pain, unless your doctor recommends use of this device.

## **Reactions**

---

Stop using the device and consult your doctor if the pain worsens or if you experience adverse reactions from the device.

Possible adverse reactions may include the following:

- skin irritation beneath the electrodes, although the gel used on the electrodes is not known to cause allergic reactions
- burns beneath the electrodes
- headaches or other painful sensations
- nausea or feeling faint

## **Side effects**

---

- After EMS, some muscle soreness may occur some time after the treatment.
- Some reddening of the skin under or around the electrodes may occur during and shortly after EMS treatment. Skin reddening usually disappears within two hours after treatment. If the skin continues to be red for more than a day, please consult your doctor.

## Warnings

---

### General warnings

- Read this user manual carefully and always adhere to the treatment instructions.
- Stop using the device and consult a doctor if the pain worsens.
- This is a medical device. Keep the device out of the reach of children. Children may not be able to use it according the instructions in this user manual. This device should not be used on or by children.
- This appliance is not intended for use by adults with reduced sensory or mental capabilities. They may not be able to use it according the instructions in this manual and they may become agitated by the treatment.
- Do not use the device on children, as the device is not designed to use on children.
- Do not modify the device or the electrodes. This could cause improper functioning.
- Only use the adapter supplied to charge the device. Use of another adapter may cause a hazardous situation.
- This device is designed for use by and on a single adult person.

### Do not apply to or treat the following parts or areas

- Stimulation on the sides of the neck over the carotid sinus may have severe adverse effects on your heart rhythm or blood pressure.

- Stimulation on the front of the neck may cause severe muscle spasms that could close your airway and result in breathing difficulties.
- Do not apply stimulation across your chest because the introduction of electrical current into the chest may cause rhythm disturbances to your heart, which could be lethal.
- Stimulation across the head or on opposite sides of the head must be avoided, as the effects on the brain are unknown.
- Do not place the electrodes on skin that is not intact, clean or healthy. Open, irritated or otherwise damaged skin may lead to the delivery of too much current to the area, which may potentially cause burns.
- If you recently had surgery, do not place one or both electrodes on the cut. Place the electrodes at a distance of 3 to 5 cm from the cut or along the proximal nerves or near the spine where the skin is intact.
- Do not place the electrodes near cancerous lesions, as this may have a negative impact on the cancerous lesion.
- Do not place the electrodes on areas of skin that lack normal sensation. This may cause you to be unaware of the high intensity current administered, which may result in electrical skin burns.
- Do not place the electrodes over open wounds or rashes. Open wounds may lead to the delivery of too much current to the area, potentially leading to

burns. It might also cause substances in the electrode to enter the skin.

- Do not place the electrodes over swollen, red, infected or inflamed areas or skin eruptions (e.g. phlebitis, thrombophlebitis and varicose veins). Stimulation must not be applied to areas with known thrombosis or thrombophlebitis, as stimulation may increase circulation, resulting in a greater risk of emboli.
- Do not place the electrodes inside body cavities, e.g. the mouth. This device is not designed for internal application.
- For EMS, only place the electrodes on healthy, uninjured muscles.

### **Do not use the device in the following conditions**

- Do not use the device if you have recently had surgery or if you are in the care of a doctor and have not consulted your doctor first.
- Do not use the device if you have a tendency to bleed internally due to any impact or injury.
- Do not use EMS stimulation to contract a muscle if contraction of the muscle may disrupt healing. For instance, if the muscle or tendon is torn, muscle contraction may exacerbate the tear, just like a voluntary contraction. This may also be the case after a recent surgical procedure or following acute trauma or fracture. Muscle contraction in case of tendinitis may also worsen symptoms.

- Do not use the device while driving, operating machines or performing any other activity in which electrical stimulation can put you at risk of injury.
- Do not use if you are likely to fall asleep during the treatment, as this may cause you to sense injuries too late. If you use the device at bedtime, set the timer to make it switch off automatically
- Do not use while in the bath, shower or swimming pool. This increases the risk of electric shocks and burns to the skin.
- Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Do not use the device when you are connected to high-frequency surgical equipment. This may result in burns on the skin under the electrodes and may damage the device.
- Do not use the device within less than 1 metre from shortwave or microwave medical equipment. Close proximity to this equipment may cause unstable device output.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

## 12 ENGLISH

- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

### **Precautions**

---

- If you had surgery recently, consult your doctor before you use the device. Using the device may disrupt the healing process.
- Use Burst TENS programs while resting, as these programs can result in muscle contraction in the treatment area.
- Make sure that you end the treatment in the Philips Treatment App or press the on/off button on the device before you remove the device or the electrodes. If you do not end the treatment, you may experience an unpleasant sensation in your fingers when you touch the magnets. This sensation is not harmful, but can be unpleasant.
- Only use this device with the adapter, cord and accessories recommended by the manufacturer.
- Always check the device and the electrodes for damage before use. Do not use the device or an electrode if it is damaged.
- The device is not waterproof. Do not use the device in wet surroundings and prevent it from getting wet.

- Although you can use the device indoors and outdoors, it does not withstand all weather conditions.
- The device does not withstand high and low temperatures. Please check the operating conditions in chapter 'Specifications'.
- The electrodes have a limited shelf life. Please check the packaging for the use-by date prior to use. Do not use electrodes whose use-by date has expired. For ordering replacement electrodes, see chapter 'Ordering accessories'.
- Do not use plaster or tape to attach the electrodes to the skin.
- Always use and store the electrodes according to the instructions in chapters 'Using the device' and 'After use'.
- If the device does not function as described here, stop using it and contact the Consumer Care Centre. For contact details, see chapter 'Guarantee and support'.

## **Compliance with standards**

---

- The device meets the relevant standards for this type of Class IIa electrical medical appliance and appliances using electrical stimulation for home use.
- This Philips appliance complies with all applicable standards and regulations regarding exposure to electromagnetic fields.

### Product overview (Fig. 1)

- 1 On/off button
- 2 Status indicator
- 3 Detachable connector
- 4 Micro-USB socket of device
- 5 Battery indicator
- 6 Adapter
- 7 Socket for standard USB plug
- 8 Standard USB plug
- 9 Micro-USB plug
- 10 Self-adhesive hydrogel electrodes

### Replacement

#### Electrodes

Replace the electrodes if:

- they are damaged or torn.
- they are past the use-by date indicated on the resealable bag.
- they have lost their adhesive power. Never use plaster or tape to attach them to your skin.
- when stimulation is uncomfortable, i.e. when you experience an unpleasant stinging or biting sensation.

*Note: Always replace the electrodes with electrodes recommended for this device by the manufacturer. You can order new electrodes at [www.philips.com](http://www.philips.com).*

### Guarantee and support

Your device has been designed and developed with the greatest possible care to guarantee an expected service life of 5 years.

If you need service or information or if you have a problem, please visit our website at [\*\*www.philips.com/support\*\*](http://www.philips.com/support). You can also contact the Philips Consumer Care Centre in your country. Its telephone number and contact URL is:

- United Kingdom: 0844 338 04 89 (5p/min. from a BT landline, other landline and mobile providers may charge more)
- Hong Kong: 852 2619 9663 (free of charge for landline calls).

## Guarantee restrictions

The terms of the guarantee do not cover the electrodes, as the electrodes are subject to wear and have to be replaced regularly.

## Manufacturer's legal address

Philips Consumer Lifestyle B.V.

Tussendiepen 4, 9206 AD Drachten, Netherlands

## Recycling

- This symbol on a product means that the product is covered by European Directive 2012/19/EU (Fig. 2).
- This symbol means that the product contains a built-in rechargeable battery covered by European Directive 2006/66/EC which shall not be disposed of with normal household waste. Follow the instructions in section 'Removing the rechargeable battery' to remove the battery (Fig. 3).
- Follow your country's rules for the separate collection of electrical and electronic products and rechargeable batteries. Correct disposal helps prevent negative consequences for the environment and human health.
- All plastic parts are marked with recycle symbols (Fig. 4).

## Removing the rechargeable batteries

*Note: Always remove the rechargeable batteries before you hand in the device at an official collection point. If you have trouble removing the batteries, you can also take the device to a Philips service centre, which will remove the battery pack for you and will dispose of it in an environmentally safe way.*

- 1** Make sure the rechargeable batteries are empty. Switch on the device as often as is necessary until you are no longer able to switch it on because the batteries have run out completely.
- 2** Insert the tip of a flat-head screwdriver between upper and lower shell at the micro-USB socket.

*Note: This action destroys the device and must only be executed if you intend to discard the device.*

- 3** Turn the screwdriver to separate the two shells.

*Note: You may have to insert and turn the screwdriver a few times to break the two shells loose from each other.*

- 4** When the shells have come apart, use the screwdriver to remove the batteries. Dispose of the plastic parts and the batteries separately according to local rules

### Specifications

Model	PR3840
Rated voltage(V) - adapter	100-240 ~ 150mA
Rated frequency (Hz) - adapter	50-60
Rated output - adapter	5V  500 mA
Class	Medical device Class IIa
Ingress of water - device	IP22
Battery type - device	Li-Ion
Adapter	The adapter is specified as part of the medical system.

### Pulse output parameters

Frequency range	1-100 Hz
Pulse duration	60-350 $\mu$ s for programs 1 to 20
Maximum current output	60mA at 500-1000 Ohm
Maximum output voltage	60V
Current pulse shape	Biphase symmetrical (net current 0 DC)

### Operating conditions

Temperature	from +5°C to +40°C
Relative humidity	from 15% to 93% (non-condensing)
Atmospheric pressure	from 700hPa to 1060hPa

## Storage and transport conditions

Temperature - electrodes	from 0°C to +40°C
Temperature - device	from -10°C to +50°C
Relative humidity	less than 93% (non-condensing)
Atmospheric pressure	from 700hPa to 1060hPa

*Note: If you want to store the electrodes for more than a month, keep them at temperatures between +5°C and +27°C.*

*Note: Do not store the electrodes in the freezer or refrigerator. Do not expose them to high temperatures, nor immerse them in water, or leave them outside of the plastic bag.*

## Explanation of symbols

### Symbols on the adapter and the device

- This symbol indicates that the adapter is double insulated (Class II), according to IEC 60601-1 (Fig. 5).
- This symbol means: Do not throw away with the normal household waste. For further instructions, see chapter 'Recycling' (Fig. 2).

### Symbols on the device

- This is the symbol for serial number. This symbol is followed by the manufacturer's serial number (Fig. 6).
- This is the symbol for reference number. This symbol is followed by the manufacturer's reference number, which is the type number of the device (Fig. 7).
- This symbol on the device means: protected against access to hazardous parts with a finger and against vertically falling water drops when tilted up to 15 degrees (Fig. 8).
- This symbol means that the part of the device that comes into physical contact with the user (also known as the applied part) is of type BF (Body Floating) according to IEC 60601-1. The applied parts are the electrodes (Fig. 9).

## 18 ENGLISH

- This symbol means that this device emits non-ionising radiation. All devices and systems that include RF transmitters or that intentionally apply RF electromagnetic energy must have a label with this symbol (Fig. 10).
- This symbol means: caution TENS output. You find this symbol near all electrode connections (Fig. 11).
- This symbol means 'Manufactured by' and appears next to the address of the legal manufacturer (Fig. 12).
- This symbol means: Conforms to EC Directives. CE stands for 'Conformité Européenne'. 0344 is the number of the notified body (Fig. 13).
- This symbol means: read the user manual before you start using the device (Fig. 14).
- This symbol is the mark that identifies the standby button (Fig. 15).
- This is the pulse symbol. It appears next to the status indicator (Fig. 16).
- This is the battery symbol. It appears next to the battery indicator (Fig. 17).

### **Symbols on the adapter**

- This symbol means: alternating current (Fig. 18).
- This symbol means: direct current (Fig. 19).
- This symbol means: Conforms to EC Directives. CE stands for 'Conformité Européenne' (Fig. 20).
- This symbol means: may only be used indoors (Fig. 21).
- This symbol means: USB connector (Fig. 22).
- This symbol means: micro-USB connector (Fig. 23).
- This symbol is the certification logo that includes identification of the accredited compliance test house. (Fig. 24)
- This symbol is the certification mark for technical equipment. The GS mark is based on the German Equipment and Product Safety Act (Fig. 25).
- This symbol shows that the no-load power CEC (California Energy Commission) Efficiency Level is V (five) in order to meet EU requirements (Fig. 26).

### **Symbols on electrode bag**

- This symbol means: Do not use while driving, operating machines or performing other activities that may present a risk of injury (Fig. 27).
- This symbol means: Do not use damaged or worn electrodes or electrodes that have lost their adhesive power (Fig. 28).
- This symbol means: The electrodes should be used by one person and not shared with others (Fig. 29).
- This symbol means: Do not place the electrodes on top of each other or so close to each other that they touch each other (Fig. 30).
- This symbol means: Do not place electrodes on your head or on your neck (Fig. 31).
- This symbol means: Do not place the electrodes on your chest (Fig. 32).
- This symbol means: Do not put the electrodes on open wounds, cancerous lesions, rashes and infected or inflamed skin. If you recently had surgery, do not place electrodes on the cut (Fig. 33).
- This symbol indicates the minimum and maximum storage temperature for electrodes if they are stored for more than a month (Fig. 34).
- This is the symbol for 'use by'. It indicates after which date the electrodes must not be used anymore (Fig. 35).

### **Symbols on the package**

- This symbol means that the package and its contents can be stored at temperatures between -10°C and +50°C (Fig. 36).
- This symbol means that the 2-year worldwide guarantee applies to this device (Fig. 37).
- This symbol means that the device contains a Li-Ion battery (Fig. 38).
- This symbol refers to the Green Dot scheme, which is covered under the European Packaging and Packaging Waste Directive (94/62/EC). Under this scheme, industry has to pay for recycling the packaging materials of consumer goods (Fig. 39)

## **Electromagnetic emissions and immunity**

### **Electromagnetic Compatibility (EMC)**

The PulseRelief device is approved according to EMC safety standard EN 60601-1-2. It is designed to be used in typical domestic or clinical environments.

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