DANSK

Nu introduserer vi PulseRelief TENS

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Introducing PulseRelief

TENS

Transcutaneous electrical nerve stimulation (TENS) is a method to relieve pain. It works via electrically stimulating nerves in the painful area.

EMS

Electrical muscle stimulation (EMS), also called neuro-muscular electrical stimulation (NMES), is a method to strengthen muscles. In healthy people, it works via electrically stimulating the motor nerves that during movement transport the signals to the muscles to contract.

Electrode placement

The electrode placement depends on the type of stimulation you would like to use - TENS or EMS - and the body part you would like to treat.

For suggestions regarding placement of the electrodes, see the ‘Electrode Placement Guide’.

For both TENS and EMS, the optimal electrode positions may vary slightly from person to person. So try moving the electrodes around until you get the most effective results. In any case, the electrodes should be placed at least 3 cm apart from each other.

For TENS (programs 1-15)

As a rule, place the electrodes on or near the area where you have pain.

Try to find an electrode position where you experience a strong but comfortable sensation from the stimulation. Stimulation may result in visible muscle contraction, especially in burst TENS programs (programs 7 to 11). This is not harmful, but if you find these contractions uncomfortable, try moving the electrodes to a different place, reduce the stimulation intensity or choose a different program.

For EMS (programs 16-20)

As a rule, place the electrodes on the muscle that you would like to strengthen so that the two electrodes are alongside of the muscle. One electrode should always be placed on the middle of the muscle belly.

When using EMS programs, muscle contraction is a desired effect and should always occur. If no muscle contraction occurs in an EMS program, you have not placed the electrodes correctly or you have chosen a too low intensity. If the muscle contractions feel uncomfortable, try varying the placement of the electrodes.

Attention

- Do not place electrodes over the carotid sinus on the sides of your neck or on the front of your neck, across your chest or on your head.
- Do not place electrodes on red or inflamed skin and skin with open wounds, on cancerous lesions or rashes.

Pre-set programs

The device comes standard with 15 pre-set programs for TENS and 5 pre-set programs for EMS.

Please note that the order of the programs in the list does not signify a recommendation of specific programs for specific conditions. As the effect of a treatment program may strongly differ from person to person, we advise you to try out different treatment programs to find out which one works best for you.

The treatment app will keep track of all your treatments to make it easy for you to filter out which one works for you.

Program specifications

The table below gives information on the frequency, pulse width and pre-set time of each program. The current pulse in all programs has a symmetrical bi-phase waveform.

Program name	Program type	Frequency	Pulse width	Pre-set time
Conventional	Conventional	100 Hz	200 µs	continuous
Thorough stimulus	Conventional	40 Hz	200 µs	continuous
Pre-activity	Conventional	100 Hz	150 µs	continuous
Soft stimulus	Conventional	80 Hz	60 µs	continuous
Deep stimulus	Conventional	80 Hz	150 µs	continuous
Radiating pulse	Conventional	60 Hz	150 µs	continuous
Continuous burst	Burst	100 Hz, 2 Hz	250 µs	30 min.
Slow kneading	Burst	80 Hz, 1 Hz	200 µs	30 min.
Mild kneading	Burst	80 Hz, 3 Hz	200 µs	30 min.
Deep kneading	Burst	80 Hz, 2 Hz	150 µs	30 min.

Diffuse burst	Burst	80 Hz, 2 Hz	180 µs	30 min.
Mixed frequency	Freq. Mod.	8 - 100 Hz	250 ms	30 min.
Pulsing massage	Freq. Mod.	2 - 8 Hz	250 µs	30 min.
Deep massage	Freq. Mod.	2 - 80 Hz	180 µs	30 min.
Gentle massage	Freq. Mod.	8 - 80 Hz	1150 µs	30 min.
Muscle stimulation pulse	Ampl. Mod.	40 Hz	200 µs	30 min.
Muscle stimulation flow	Ampl. Mod.	50 Hz	200 ms	20 min.
Muscle stimulation wave	Ampl. Mod.	65 Hz	200 µs	30 min.
Muscle stimulation flux	Ampl. Mod.	40 Hz	350 µs	25 min.
Muscle stimulation bounce	Ampl. Mod.	50 Hz	300 ms	20 min.

Program overview

Conventional TENS programs

Conventional TENS is a common TENS program with high-frequency stimulation. It is based on the Gate Control Theory that states that this type of stimulation will block transmission of pain signals in the nervous system. It is comparable to ‘rubbing the pain away’. Most people use Conventional TENS for 30 minutes per session, but you can use it as often and as long as it feels comfortable. The effect is usually most perceptible during stimulation and generally wears off rather quickly after stimulation. You have to adjust the intensity of the program so that you clearly feel the stimulation. Stimulations must not feel unpleasant or painful.

Name of program	Description	Suggested treatment time
Conventional	May give a continuous tingling stimulation	No limit
Thorough stimulus	May give a more thorough stimulating sensation	No limit
Pre-activity	May give a more local and soft stimulating sensation	No limit
Soft stimulus	May give a softer stimulating sensation, e.g. for sensitive body parts	No limit
Deep stimulus	May give a deeper stimulating sensation	No limit
Radiating pulse	May give you a more radiating stimulation sensation	No limit

Burst TENS programs

Burst TENS offers a high-frequency stimulation that comes in trains of pulses. This type of stimulation is known to trigger endorphin release. Endorphins are painkillers produced by the body. You can use Burst TENS around 30 minutes multiple times a day, for example 3 times. We advise you to adjust the intensity of the Burst program so that you feel the stimulation strongly. Burst programs generally give a sharper sensation, which can feel somewhat unpleasant. It is also advised to adjust the intensity up to a level that produces visible muscle contractions in the treatment area. These contractions are not harmful. If you find these contractions too uncomfortable, you can decrease the intensity, choose a different program or move the electrodes to a different place. Apply Burst TENS programs while resting, as they can result in muscle contraction in the treatment area.

Name of program	Description	Suggested treatment time
Continuous burst	May give a pulsing burst stimulation to release endorphins	30 min. *
Slow kneading	May give a slower and massaging burst stimulation to release endorphins	30 min. *
Mild kneading	May give a milder burst stimulation to release endorphins	30 min. *
Deep kneading	May give a more local and deeper burst stimulation to release endorphins	30 min. *
Diffuse burst	May give a more diffuse burst stimulation to release endorphins	30 min. *

*optionally multiple times a day, can be changed in the app.

Frequency-modulated TENS programs

Frequency-modulated TENS is a program in which the frequency of the pulse varies from high to low. This can be more comfortable or effective than the other programs in some cases. You can use frequency-modulated TENS for around 30 minutes multiple times a day, for example 3 times.

Name of program	Description	Suggested treatment time
Mixed frequency	May give a flowing stimulation sensation	30 min. *
Pulsing massage	May give a more pulsing stimulation sensation	30 min. *
Deep massage	May give a slower and deeper flowing sensation	30 min. *
Gentle massage	May give a mild flowing stimulation sensation	30 min. *

* optionally multiple times a day, can be changed in the app.


Electrical muscle stimulation (EMS) programs

Electrical muscle stimulation (EMS), also known as neuromuscular electrical stimulation (NMES) or electromyostimulation, uses pulses to evoke muscle contraction. This can for example help to strengthen muscles. Use EMS only on healthy, uninjured muscles. Apply EMS programs while resting (sitting, lying), as they result in muscle contraction. First-time users can experience electrical muscle stimulation as an unfamiliar and unusual sensation. We therefore advise you to start at a low intensity setting to get used to the sensation before you choose a higher intensity setting. In general, stimulation must always feel comfortable. The muscle contractions must not feel exhausting, strenuous or tiring. After strenuous exercise or exertion, we advise you to use a lower setting to avoid muscle fatigue.

Name of program	Description	Suggested treatment time
Muscle stimulation pulse	May give nice soothing sensation	30 min. *
Muscle stimulation flow	May give comfortable tingling sensation	20 min.
Muscle stimulation wave	May give mild contracting and relaxing sensation	30 min. *
Muscle stimulation flux	May give energizing muscle sensation	25 mon.
Muscle stimulation bounce	May give rhythmically relaxing sensation	20 min. *

* optionally multiple times a day, can be changed in the app

SV: Elektrodpliceringsguide för TENS
UK: Electrode placement guide for TENS

SV: Elektrodpliceringsguide för EMS
UK: Electrode placement guide for EMS