



PHILIPS

Magnetic Resonance

SmartPath upgrade

Ingenia 1.5T Evolution

Step into the future,
with confidence

Step into the future, with confidence

The healthcare industry is under tremendous pressure to hold costs down while delivering exceptional patient care. As a result, there is a clear trend towards keeping existing medical equipment.

With Philips MR SmartPath you can convert your existing Ingenia 1.5T system to an Ingenia 1.5T Evolution. It enables you to boost your performance with innovative SmartWorkflow solutions that include touchless patient sensing technology, in-room guidance on patient set-up and initiation of the exam at the patient's side. Adding SmartSpeed allows you to perform MRI exams up to 3 times faster¹ in both 2D- and 3D scanning and for all anatomies. MR SmartPath also includes the ComfortPlus mattress, which delivers a more comfortable table experience². A positive patient experience can be further supported through installation of an immersive audio-visual experience that calms patients and guides them through MR exams. Your MR SmartPath conversion also gives you access to the latest scanning techniques for confident diagnosis of the most challenging indications.

**Discover how Philips SmartPath program can help you step into the future with confidence.
It's your path to longer system lifetime**

¹ Compared to SENSE imaging.
² Compared to the standard mattress.

Extend the lifetime of your equipment	4
Boost your performance with latest workflow solutions	6
Scan up to 3 times faster, with equal or better image quality¹	17
Up to 50% faster MRI exams with virtually equal image quality¹	25
A more comfortable experience for your patients	33
Access to the latest scanning techniques	39
Enhance the value of your MR investment	46

Extend the lifetime of your equipment

SmartPath allows you to enhance your investment, extend the lifetime of your equipment and easily upgrade to the latest technology for long term success. With this program, you completely renew your trusted Ingenia MR system, just as if you had bought a new MR, extending the lifetime of your equipment and improving your total cost of ownership. Furthermore, it alleviates the delays and expense of installing a new MR magnet, that can involve breaking down walls and ceilings and using heavy cranes to transport the magnet in and out of the hospital.



Your current Ingenia 1.5T...



... will be dismantled and upgraded to become ...



... your new Ingenia 1.5T Evolution

A sustainable choice for substantially less cost

Re-using your existing magnet and converting to the next generation of MR costs substantially less than purchasing a new system. When you convert your MR instead of buying a new system, you also make the sustainable choice . Saving the CO₂ output and energy usage that would be required to manufacture a new magnet , and saving the costs of transporting, lifting and installing a new magnet weighing up to 3000 kilograms.



Make the sustainable choice



For less cost than purchase of a new system

Boost your performance with **latest workflow solutions**

With a growth in the elderly population and constant demands to do more with less, the pressure on healthcare providers is immense. This pressure is also evident in radiology departments and imaging centers. The increasing use of MR to diagnose a variety of conditions and illnesses has led to demands for greater efficiency, even as departments try to manage a shortage of MR operators and variability in staff expertise. Too often, it seems that productivity is at odds with giving patients the time and attention they desire.

With SmartWorkflow, you can achieve high productivity while enabling your staff to focus on patients. It reduces and simplifies the number of steps needed in a conventional MR exam workflow, using technology to guide and coach where required, and automate where possible. An end-to-end workflow solution that directly boosts efficiency through reduced variability and task automation, while supporting a better patient and staff experience, resulting in patient-centered productivity.



Smart Workflow in the exam room

“The entire workflow is smooth: Patient positioning and setup; launching the scan as soon as we leave the exam room; the intuitive touchscreen on the gantry; Touchless patient sensing... All of these things are much better than on our old system.”

Laura Barlow, RTMR MRI Technologist
Supervisor at the University of British Columbia





Guided exam set-up

Coaching and visual guidance are provided at the front of the magnet façade




Auto patient centering

Region of interest is automatically placed in the iso-center of the magnet



Touchless respiratory-triggering

Patient's breathing is detected without any operator interaction



Auto coil element selection

Optimal elements are selected automatically based on the anatomy planning



In-room exam start

Exam start can be initiated with a single touch of at the patient's side

* According to the definition of AI from the EU High-Level Expert Group.

Smart Workflow in the control room

“ We don’t have to manually direct the patient to breathing and not breathing. We can go ahead and let the machine do the work of the breathing instructions while we continue our planning of the exam.”

Carlos Avila, RT Technologist at Miami Cardiac & Vascular Institute




Confidence for MR Conditional implants
Step-by-step guidance to enter the condition values as specified by the implant manufacturer



Automated planning, scanning and processing
Fully automated geometry planning, coil element selection and execution of complete MR exams



Up to 3 times faster imaging¹
Breakthrough acceleration technique delivering image quality and speed without compromise



Automated patient coaching
Patients are guided via announcements of scan duration, table movements and breath hold instructions



Plan your day in advance
Dashboard to plan examinations before patient arrival, allowing you to stay on schedule

¹ Compared to SENSE imaging.
* According to the definition of AI from the EU High-Level Expert Group.



A virtual coach **guiding exam set-up**

Increase staff confidence and speed up patient set-up through automated real-time guidance and insights on the details of the current patient study. Achieve high quality results, independent from staff's expertise level. VitalScreen provides guidance at your staff's fingertips. Two 12-inch interactive touchscreens on the scanner provide coaching and visual guidance on recommended patient position, study laterality, coil and accessory placement. Moreover, feedback is provided on important exam details, including physiology signals (both VCG and respiratory) and – if applicable – contrast usage and breath-hold guidance.

Put your patients at ease, while manual steps in the workflow are automated

Free up your staff from monotonous, manual steps and enable them to focus on the patient through automatic placement of the region of interest in the scanner iso-center. The manual use of a laser light visor for iso-center positioning has become obsolete. VitalScreen automatically detects landmarks for selected anatomies and places the region of interest in the iso-center of the magnet. Once the patient is positioned on the table, only the push of a button is required to position the patient in the center of the bore.

Increase productivity and free up time for other tasks

Start exams as soon as possible, eliminating extra steps for your staff and decreasing the time the patient has to spent in the magnet, resulting in a more positive patient experience. VitalScreen allows staff to initiate the exam with a single touch of a button at the patient side. The exam starts immediately after the operator has closed the exam room door, so no time is wasted.

Continuous and robust respiratory signal providing superior image quality

Relieve your staff from the burden of positioning – and re-positioning – a respiratory belt. Positioning a belt shifts the operator's focus from the patient to the technology at a moment when it is critical that the patient is comfortable and reassured. Enjoy optical sensing and AI¹ to automatically detect patient respiratory patterns. VitalEye touchless patient sensing provides a fast detection of patient's breathing without any operator interaction. With VitalEye, the technologist no longer needs to set up an old-fashioned respiratory belt but receives a continuous and robust respiratory signal without any interaction. This revolution in touchless patient sensing helps your staff to keep a caring eye on your patient. The quality of the physiology signal detected by VitalEye is better than a belt-based approach providing superior image quality, for a broad range of patient sizes.

“ It always works, and it's always there.”


MR operator, University of Bonn, Germany



¹ AI stands for Artificial Intelligence, according to the definition of AI from the EU High-Level Expert Group.

Touchless patient sensing

Superior image quality with VitalEye¹, consistently



¹ Compared to Philips belt-based signal. Requires an unobstructed line of sight. Results from case studies are not predictive of results in other cases. Results in other cases may vary.



Up to 3 times faster MRI exams with no loss in image quality¹


Time is one of the most precious commodities you have in your MR department. What if we told you there was a way to recover time you have been losing during your MR examinations? And use the time you do have more wisely? Imagine how that could help you make better use of your scarce resources and better meet the demands of referring physicians. That's exactly what our acceleration technologies such as Philips SmartSpeed can do for your MR department.

SmartSpeed AI can speed up scan time nearly 3 times with no loss in image quality¹, free up time to improve your patient experience. You can use the time gained to scan more patients and reduce the cost per scan, to add unplanned patients to the schedule or free up time to improve your patient experience. It can also provide higher image quality¹ to enhance diagnostic confidence.


¹ Compared to SENSE imaging.

Fast scan times with SmartSpeed

SmartSpeed




Sagittal T2w TSE
SmartSpeed factor 2.2
0.6 x 0.7 x 3.0 mm
1:44 min




Axial PDw TSE
SmartSpeed factor 2.1
0.5x0.6x 3mm
2:10 min



Coronal PDw TSE
SmartSpeed factor 3
1.0 x 1.3 x .0 mm
2:15 min




Axial PDw TSE SPAIR
SmartSpeed factor 2.5 1.3
x 1.6 x 3.0 mm
2:37 min




Coronal T1w TSE
SmartSpeed factor 2.5
0.7 x 0.9 x 3.0 mm
1:42 min


SmartSpeed




Sagittal PDw TSE SPAIR
SmartSpeed factor 2
0.5 x 0.6 x 3.0 mm
1:50 min



Coronal T1w TSE
SmartSpeed factor 2.2
0.4 x 0.5 x 3.0 mm
1:14 min



Axial T2w TSE
SmartSpeed factor 2.2
0.6 x 0.7 x 3.0 mm
1:24 min




Coronal T2w TSE
SmartSpeed factor 2.4
0.6 x 0.7 x 3.0 mm
2:12 min




Sagittal PDw TSE
SmartSpeed factor 2.5 0.3
x 0.4 x 3.0 mm
2:10 min


C-SENSE




Sagittal T2w TSE
C-Sense factor 1.6
0.6 x 0.7 x 3.0 mm
2:19 min



Axial PDw TSE
C-Sense factor 1.6
0.5 x 0.6 x 3.0 mm
2:50 min



Coronal PDw TSE
C-Sense factor 2.2
1.0 x 1.3 x 3.0 mm
2:55 min




Axial PDw TSE SPAIR
C-Sense factor 1.8
1.3 x 1.7 x 3.0 mm
3:24 min




Coronal T1w TSE
C-Sense factor 2.1
0.7 x 0.9 x 3.0 mm
2:07 min


C-SENSE




Sagittal PDw TSE SPAIR
Sense 1.2
0.5 x 0.6 x 3.0 mm
3:00 min




Coronal T1w TSE
Sense 1.2
0.4 x 0.5 x 3.0 mm
2:13 min



Axial T2w TSE
Sense 1.2
0.6 x 0.7 x 3.0 mm
2:30 min



Coronal T2w TSE
Sense 1.4
0.3 x 0.4 x 3.0 mm
3:36 min



Sagittal PDw TSE
Sense 1.5
0.3 x 0.4 x 3.0 mm
3:40min

Motion-Free scans with SmartSpeed

SmartSpeed



Axial T2w FLAIR Motion free
SmartSpeed factor 1.4
0.4 x 0.4 x 5.0 mm,
4:00 min



Axial T1w MotionFree
SmartSpeed factor 1.8
0.4 x 0.4 x 5mm
2:58 min



Axial T2w MotionFree
SmartSpeed factor 1.8
0.4 x 0.4 x 5.0 mm,
2:08 min

C-SENSE



Axial T2w FLAIR MultiVane XD
SENSE factor 1.4
0.6 x 0.6 x 5.0 mm,
4:00 min




Axial T1w MultiVane XD
SENSE factor 1.8
0.4 x 0.4 x 5.0 mm,
2:58 min




Axial T2w MultiVane XD
SENSE factor 1.5
0.4 x 0.4 x 5.0 mm,
2:32 min

Prostate imaging in faster scan times with SmartSpeed


SmartSpeed



Axial T2w TSE PI-RADS
SENSE factor 1.3
0.4 x 0.7 x 3.0 mm
5:55 min




Axial T2w TSE PI-RADS
SmartSpeed factor 1.4
0.4 x 0.7 x 3.0 mm
3:37 min




Axial T2w TSE High Resolution
SmartSpeed factor 2.5
0.4 x 0.4 x 3.0 mm
3:37min


C-SENSE



Coronal T2w TSE PI-RADS
SENSE factor 1.3
0.4 x 0.7 x 3.0 mm
5:55 min



Coronal T2w PI-RADS
SmartSpeed factor 1.5
0.4 x 0.7 x 3.0 mm
3:27 min



Coronal T2w High Resolution
SmartSpeed factor 2.5
0.4 x 0.4 x 3.0 mm
4:10 min



Up to 50% faster MRI exams with virtually equal image quality¹

Time is one of the most precious commodities you have in your MR department. What if we told you there was a way to recover time you have been losing during your MR examinations? And use the time you do have more wisely? Imagine how that could help you make better use of your scarce resources and better meet the demands of referring physicians.


That's exactly what Compressed SENSE can do for your MR department. It accelerates your existing MR scans by up to 50% with virtually equal image quality, frees up time to improve your patient experience and can provide up to 60% higher resolution to enhance diagnostic confidence.¹



Compressed SENSE is suitable for all anatomies and can be used for all anatomical contrasts, in both 2D- and 3D scanning.


¹ Compared to Philips scans without Compressed SENSE

Up to 50% faster MRI exams with virtually equal IQ¹

Brain ExamCard with 2D and 3D protocols




	DWI (b1000)	1.4 x 2.0 x 3.0 mm	2:20 min
	T2w MultiVane XD	0.6 x 0.6 x 4.0 mm	2:32 min
	3D T2w FLAIR	1.1 x 1.1 x 1.1 mm	5:02 min
	3D T1w TFE	1.0 x 1.0 x 1.0 mm	5:36 min
	3D PCA	0.9 x 0.9 x 1.6 mm	6:49 min
	3D Inflow	0.6 x 0.9 x 1.4 mm	3:58 min
	 17:24 min		

	DWI (b1000)	1.4 x 2.0 x 3.0 mm	2:20 min	0%
	T2w MultiVane XD	0.6 x 0.6 x 4.0 mm	2:32 min	0%
	3D T2w FLAIR	1.1 x 1.1 x 1.1 mm	3:41 min	32%
	3D T1w TFE	1.0 x 1.0 x 1.0 mm	2:55 min	52%
	3D PCA	0.9 x 0.9 x 1.6 mm	3:31 min	49%
	3D Inflow	0.6 x 0.9 x 1.4 mm	2:25 min	37%

Spine ExamCard with 2D protocols




	Sag T2w TSE	0.8 x 1.0 x 4.0 mm	3:10 min
	Sag T1w TSE	0.8 x 1.1 x 4.0 mm	2:01 min
	Cor T2w SPAIR	0.9 x 1.2 x 4.5 mm	1:53 min
	Ax T2w TSE	0.6 x 0.8 x 4.0 mm	3:23 min
	 10:27 min		


	Sag T2w SE	0.8 x 1.0 x 4.0 mm	2:09 min	33%
	Sag T1w TSE	0.8 x 1.1 x 4.0 mm	1:18 min	41%
	Cor T2w SPAIR	0.9 x 1.2 x 4.5 mm	1:19 min	22%
	Ax T2w TSE	0.6 x 0.8 x 4.0 mm	2:08 min	35%
	 6:54 min			

¹ Compared to Philips scans without Compressed SENSE. Results from case studies are not predictive of results in other cases. Results in other cases may vary.
Courtesy: Kantonsspital Winterthur, Switzerland, Ingenia 1.5T

Results from case studies are not predictive of results in other cases. Results in other cases may vary.
Courtesy: Radiologie am St.Josef Stift, Bremen, Germany, Ingenia 1.5T


Fast push button exams

 Brain 6:17 min



Sag 3D T1w TFE	1.1 x 1.1 x 1.1 mm	1:26 min
Ax DWI	2.0 x 2.0 x 5.0 mm	0:36 min
Ax T2w TSE	0.7 x 0.9 x 6.0 mm	0:31 min
Cor T2w FLAIR	0.9 x 1.1 x 6.0 mm	1:28 min
Cor 3D TOF	0.8 x 0.8 x 1.2 mm	2:16 min


 L-Spine 4:21 min



Sag T2w TSE	0.9 x 1.2 x 4.0 mm	1:04 min
Sag T1w TSE	0.9 x 1.2 x 4.0 mm	0:56 min
Ax T2w TSE	0.7 x 1.0 x 4.0 mm	1:04 min
Cor STIR	1.0 x 1.5 x 4.0 mm	1:17 min


Enabled by Compressed SENSE and SmartWorkflow

 Knee 4:17 min



Sag T1w TSE	0.4 x 0.6 x 3.0 mm	0:57 min
Sag PDw TSE	0.4 x 0.5 x 3.0 mm	0:53 min
Cor PDw FatSat	0.5 x 0.9 x 3.0 mm	1:16 min
Ax PDw FatSat	0.5 x 0.8 x 3.0 mm	1:11 min

 Ankle² 7:46 min




Sagittal PDw TSE mDIXON (W+IP)	0.5 x 0.7 x 3.0 mm	2:39 min
Axial PDw TSE FatSat	0.4 x 0.6 x 3.0 mm	2:36 min
Coronal PDw TSE FatSat	0.5 x 0.7 x 3.0 mm	1:24 min
Coronal T1w TSE	0.4 x 0.5 x 3.0 mm	1:07 min

Add extra patient slots to your daily MRI schedule


Many radiology departments and imaging centers are looking for ways to increase the utilization of their MR equipment to meet the rising demand for MRI services. A full MRI exam performed with Compressed SENSE, for example, can save minutes compared to a conventional MRI exam. This could free up one or two extra exam slots in your daily schedule, which can result in much higher productivity and shorter waitlists without adding more operator hours.

Scan >5 patients/day extra within same scanning hours



Radiologie Dr Wagner in Gottingen , Germany can accommodate > 5 more patients per day, within the same scanning hours, after the introduction of Compressed SENSE.

Allowing you to reduce your backlog faster




Easily fit in unplanned patients

Do unscheduled patients disrupt your daily schedule and put extra stress on your staff? With Compressed SENSE you can create a buffer to easily handle emergency cases or urgent patients that are referred on the same day. This extra capacity can help you serve patients and referring physicians faster and make daily workflow go smoother.

“We can now provide a more flexible and faster MRI service to our patients and referring physicians. For instance, when a referring physician is requesting it, we can now quite smoothly insert an additional MRI examination without previous appointment on the same day.”


Hideki Koyasu, MD, Neurosurgical Clinic in Kanagawa, Japan

Easily fit in unplanned patients




Reduction in overtime, while maintaining same patient throughput per day

Having to work overtime is a recurring issue for many radiology departments and imaging centers that can impact staff satisfaction and run-up operational costs. The stress caused by heavy workloads and overtime hours greatly contribute to burnout among radiology technologists, not to mention long-term mental and physical health issues.¹ By reducing MRI scan times and improving scheduling flexibility, Compressed SENSE helps patients and staff to get home on time. This can improve the experience for all involved.



ComputerTomography Institut in Innsbruck, Austria has been able to reduce overtime by more than one hour, keeping the same patient throughput per day, after the introduction of Compressed SENSE.

¹ Vinu, Raj. Occupational stress and Radiography. NCBI. Nov-Dec 2006. <https://www.ncbi.nlm.nih.gov/pubmed/17119177>
Results from case studies are not predictive of results in other cases. Results in other cases may vary.





A more comfortable experience for your patients

Your patients are at the heart of Ingenia 1.5T Evolution – which includes an MR experience that enhances comfort and compliance. With up to 80% acoustic noise reduction¹, voice guidance, immersive in-bore visuals and a comfortable table mattress, Ingenia 1.5T Evolution helps your patients feel at ease, resulting in smooth, fast exams.

¹ Compared to scanning without ComforTone.

Provide an **immersive visual experience**

Your patients' scanning experience is significantly enhanced with Ingenia 1.5T Evolution. Designed to offer a relaxing sensory experience, Ambient Experience provides positive distractions for patients by incorporating dynamic lighting, projection and sound, contributing to a positive, engaging environment to benefit quality of care. From the moment a patient is moved into the scanner (the point at which people report the most stress), through completion of the scan, the In-Bore Connect solution can help patients to relax, follow directions and minimize motion. In a study, conducted using our in-bore solution, Herlev Gentofte University Hospital in Denmark managed to reduce the number of rescans by up to 70%¹. A case study at Radiologisches Zentrum am Kaufhof, Lübeck, Germany showed that the number of patients needing sedation was reduced by 80%².

***We've had a lot of patients provide compliments on the environment.
We have the Ambient solution in there that creates a soothing environment.”***

“

Carol Melvin, COO, Miami Cardiac and Vascular Institute



¹ Compared to the average of the other 5 Philips Ingenia MR scanners without Ambient Experience and In-Bore Connect. Results from case studies are not predictive of results in other cases. Results in other cases may vary.

² Results from case studies are not predictive of results in other cases. Results in other cases may vary. *The tranquilizer referred to is a valium-based derivative called "Diazepam".





Comfort in every detail

Because no detail is too small when it comes to helping your patients feel comfortable, Ingenia 1.5T Evolution includes the ComfortPlus mattress. On average, 90% of patients in severe discomfort find it easy to lie still on the ComfortPlus mattress. Overall comfort for this group of patients can increase by up to 36%.¹

“The most frequent comment we are getting from our technologists, is that for patients who have had scans on other Philips scanners, this new mattress is really significantly more comfortable.”

Dr. Oswood, Hennepin County Medical Center

¹ Compared to using a standard mattress.

90%

of patients in severe discomfort found it easy to lie still on the Comfort Plus Mattress during the exam

87%

of the patients indicated it was easy to lie still on the Comfort Plus Mattress



Access to the **latest scanning techniques**

Your SmartPath conversion also gives you access to the latest scanning techniques for confident diagnosis of the most challenging indications.


Robust suppression of motion artifacts with DWI TSE and MultiVane XD. Reduced geometrical distortion¹ with Zoom Diffusion and reduced susceptibility artifacts² caused by metal implants³ with O-MAR XD. Address patients who have difficulty holding their breath with 4D FreeBreathing. Enhance your clinical workflow by generating synthetic high b-value images with Computed DWI. Non-invasive quantification of liver fat fraction with mDIXON Quant and assessment of differences in liver tissue stiffness with MR Elastography.

¹ Due to reduced EPI echo train length in DWI-EPI compared to conventional Philips full FOV DWI-EPI.


² Compared to standard high bandwidth spin-echo based techniques.

³ Only for use with MR Safe or MR Conditional implants by strictly following the Instructions for Use.


Enhance your diagnostic confidence for Brain imaging



3D BrainVIEW – T2w FLAIR with Compressed SENSE
Reformats in any plane without loss of resolution
1.1 x 1.1 x 1.1 mm, 3:55 min



Zoom Diffusion, b1000 & ADC
Small FOV imaging with reduced geometrical distortion¹
1.2 x 1.2 x 3.0 mm, 5:40 min



DWI XD TSE
Robust supression of motion artifacts²
1.5 x 1.5 x 4.0 mm, 5:30 min



bFFE
0.6 x 0.6 x 0.6 mm, 4:38 min



bFFE XD
Better visualization of fine structures
0.6 x 0.6 x 0.6 mm, 3:35 min

Uniform and complete fat-free Spine imaging




T2w TSE mDIXON XD Compressed SENSE
Water only + In Phase
0.8 x 1.1 x 4.0 mm, 5:20 min
Courtesy : Radiologie am St Josef Stift, Bremen, Germany



T2w TSE
0.8 x 1.0 x 3.0 mm
4:10 min
Courtesy: Kurashiki Central Hospital, Japan



Zoom Diffusion, b1000
2.5 x 2.5 x 3.0 mm
3:45 min



3D NerveVIEW (MIP)
1.2 x 1.2 x 1.0 mm, 5:50 min


40 ¹ Due to reduced EPI echo train length in DWI-EPI compared to conventional Philips full FOV DWI-EPI | ² Compared to Philips multi shot DWI TSE
Results from case studies are not predictive of results in other cases. Results in other cases may vary.

Results from case studies are not predictive of results in other cases. Results in other cases may vary.

Fast and robust MSK imaging





PDw TSE SPIR
Compressed SENSE
0.3 x 0.4 x 2.0 mm
3:32 min




3D VIEW – PDw SPAIR
Compressed SENSE
0.8 x 0.8 x 0.8 mm
4:47 min



T1w TSE
Compressed SENSE
0.3 x 0.4 x 2.0 mm
2:01 min




PDw TSE
1.1 x 1.4 x 3.0 mm
2:38 min




PDw TSE – O-MAR XD
1.3 x 1.6 x 3.0 mm
7:55 min


New possibilities for Body and Breast imaging




T2w SPIR
MultiVane XD
1.1 x 1.1 x 3.0 mm
3:36 min
Courtesy: Toyonaka Municipal Hospital, Ingenia 1.5T




T2w TSE
MultiVane XD
1.0 x 1.0 x 3.0 mm
3:18 min




mDIXON Quant
Fat Fraction
2.0 x 2.0 x 2.3 mm
15 sec
Courtesy: Kurashiki Central Hospital, Japan



MR Elastography
Stiffness map
5.2 x 5.2 x 10.0 mm
14 sec (8 slices)




3D VIEW - T2w TSE
Compressed SENSE
0.8 x 0.8 x 1.0 mm
3:22 min
Courtesy: Radiologie am St.Josef Stift, Bremen, Germany



3D mDIXON XD
T1w FFE
0.8 x 0.8 x 1.0 mm
2:16 min


Diagnose liver lesions without breath holds



4D FreeBreathing
1.7 x 1.7 x 1.8 mm

Multi-phase contrast-enhanced MRI
Temporal resolution of 3 seconds per phase
Real-time reconstruction = images available for viewing during the acquisition

Generate synthetic high b-value images



Acquired DWI, b1000 + ADC
3.1 x 2.6 x 3.0 mm, 2:15 min

Computed DWI, b2000

Computed DWI, b3000

Computed DWI, b5000

Enhance the value of your MR investment

Imaging is both a clinical and an economic challenge. You need to manage a host of financial obligations and opportunities, all while keeping your focus on your patients. We can help, by putting together a package of offerings that keep total cost of ownership in check while providing you with tailored solutions for maintenance, fleet management, cybersecurity, education and financing.

Prevent issues before they occur

Scanner downtime can disrupt your schedule and delay patient care. We offer maintenance agreements that are suited to your needs, enabled by the latest service innovations and including an uptime guarantee. We prevent issues before they occur through proactive remote monitoring, remote diagnostics and remote and field service support. With e-Alerts and other remote data, we monitor more than 500 parameters of your MR system from a distance, detecting and resolving issues without impacting your department's operations. In fact, more than 50% of MR service cases are resolved remotely.¹ Our Philips-qualified service experts can also proactively resolve issues on-site, fix your system before it causes any disruption, and provide reliable and knowledgeable support.²

500
parameters are
monitored on an MR




+90,000
remote connections across
25,000
healthcare facilities in

139
countries³



25%
connected Philips
MR service cases² are
resolved before they
cause downtime, due
to proactive monitoring

>50%
of MR service cases
are resolved remotely¹



¹ For the Philips diagnostic imaging installed base

² Based on data collected between July '18 and July '19 on all service events registered on remotely connected Philips MR systems [globally]. Downtime does not include time due to planned maintenance

³ Based on global Philips only data.


Protecting your MR equipment from patient data breaches and cyber-attacks

Protecting patient health information requires constant vigilance. To keep health information and medical devices secure, we employ best practices in medical device security. Our multi-layered defense barriers include security policies, procedures, access controls, technical measures, training, and risk assessments. The Technology Maximizer Plus subscription program conveniently keeps your MR systems up-to-date through access to the latest cybersecurity patches and mandatory safety fixes via regular and ongoing software upgrades and hardware refreshes.

Standardize your MR fleet at a fixed annual cost


If you own more than one Philips scanner, standardizing under the same software release can enhance efficiency through one user interface for operators to learn and use the same ExamCards across multiple scanners. The Ingenia 1.5T Evolution is delivered with the latest available software release, providing a perfect opportunity to upgrade your fleet to this release and enter into a Technology Maximizer Plus subscription program.³ Under the program, your Ingenia 1.5T Evolution and the rest of your fleet will receive software updates whenever available, giving you the benefits of software improvements and cyber-security advances while maintaining all your MR systems on the same level.

Cost of maintenance and upgrades in previous years




- Large variation in annual cost
- Upgrading once in 4 years

Cost of maintenance and upgrades with Technology Maximizer program



- ✓ Fixed annual cost
- ✓ Yearly updating

Reduction in accumulated cost of maintenance and upgrades over 4 years Before vs with Technology Maximizer



- ✓ Technology Maximizer saves cost and provides more frequent updates

“It was consistently a challenge to plan for annual upgrades and predict their costs. Thanks to Technology Maximizer, we can now continue to have the latest versions of software for all of our MRI systems.”

Eliseo Vañó Galván, MD, Cardiovascular radiologist, Chairman of the CT & MR Department at Hospital Nuestra Senora del Rosario, Madrid, Spain

1 Based on global Philips-only data.
2 Requires minimum maintenance contract. Conditions apply. Offerings are available in selected countries and for selected products only.
3 Check for compatibility with your Philips representative.

Achieve excellence through **ongoing education**

Delivering consistent healthcare day-in and day-out is a challenge, particularly when faced with staff shortages and the need to cross-train department personnel. Our Philips MR Healthcare Education can help unlock the full potential of your staff, technology, and organization through innovative and meaningful healthcare education, delivered on-site or as e-Learning. For example, the Philips MR Technologist Development Program at Burjeel Hospital for Advanced Surgery (BHAS), a leading orthopedic and joint care center in Dubai, UAE, resulted in an average of 30% improvement in image quality across all procedures.¹ Team knowledge increased 30–40% in the key areas of patient care¹, imaging procedures, data acquisition and physics of image formation. The comprehensive, clinically-relevant courses, programs, and learning paths are designed to support clinical excellence, enhance operational efficiency and provide high-quality patient care.

Tailored financing solutions in line with your cash flow needs, budgets, and business strategy

Providing access to best-in-class healthcare is a leading priority for facilities like yours around the globe. At the same time, financial security and protecting your assets over time are also high on the agenda. To manage your financial challenges, you need to know whether your healthcare investments are sustainable – and how to get the most from your equipment. Financing your Ingenia 1.5T Evolution helps you exchange variability and unpredictability for visibility and certainty. This helps you avoid the burden and risk of upfront expenditures and benefit from transparent, predictable cost structures. As a result, you can manage and plan budgets more efficiently and free up capital that would otherwise be tied up in fixed assets.

¹ Results from case studies are not predictive of results in other cases. Results in other cases may vary





How to reach us

Please visit www.philips.com/mri
healthcare@philips.com