# PHILIPS

*SmartSleep*

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
## Introduction

Congratulations on your purchase and welcome to Philips! To fully benefit from the support that Philips offers, register your product at **[www.philips.com/welcome](http://www.philips.com/welcome)**.

The first Philips solution designed to alert you before your snoring starts. The positional monitor detects when you're moving onto your back and delivers adaptive vibrations that prompt you to sleep on your side – leading to quieter nights for you and your partner<sup>1</sup>.

<sup>1</sup>The product is for position-dependent back snorers for whom snoring ceases when sleeping on their side.

Tip: Watch the how-to videos of the Snoring Relief Band at

**<https://www.usa.philips.com/c-e/smartsleep/snoring-relief-band>**

## General description

- |   |                            |    |                              |
|---|----------------------------|----|------------------------------|
| 1 | Sensor                     | 10 | Clip                         |
| 2 | Activation button          | 11 | Window for display           |
| 3 | Display                    | 12 | Window for activation button |
| 4 | Adaptation phase indicator | 13 | Pouch for sensor             |
| 5 | Battery status indicator   | 14 | Insertion slit for sensor    |
| 6 | Response rate              | 15 | Strap                        |
| 7 | Length of sleep session    | 16 | Strap adjuster               |
| 8 | Hours on your back         | 17 | Adapter                      |
| 9 | Socket for USB             | 18 | USB cable                    |

## Intended use

The Philips SmartSleep Snoring Relief Band is intended to reduce snoring disruptions, leading to an improved sleep experience. This product is not a medical device and is not intended for use in the mitigation and treatment of disease or conditions including POSA (Positional Obstructive Sleep Apnea).

## When not to use

- The Snoring Relief Band must not be used if your breathing repeatedly stops and starts during sleep while snoring. You might have obstructive sleep apnea which is potentially a serious sleep disorder. Consult a doctor for diagnosis.
- Do not use the Snoring Relief Band if you are younger than 18 year.

- The Snoring Relief Band should not be used by a person who has a medical implant (e.g. pacemaker, insulin pump) that can be affected by mild vibrational stimuli on the chest, without consulting a doctor first.
- During pregnancy, consult a doctor before use.
- The Snoring Relief Band is not recommended for persons who are supposed to sleep on their back due to a physical condition (e.g. because of shoulder or back problems).
- The Snoring Relief Band cannot be used by persons who sleep in an upright position or require more than 2 pillows during sleep.
- Do not use the Snoring Relief Band if you are not able to switch sleeping positions independently due to physical limitations.
- Do not use the Snoring Relief Band in case of existing allergies to Polyester, Nylon or Spandex.

## **Important safety information**

Read this important information carefully before you use the product and save it for future reference.

# DANGER

## **To reduce the risk of electric shock:**

- The Snoring Relief Band is for indoor use only.
- Keep the adapter and the sensor dry and away from liquids and moisture (running tap, bathroom etc.).
- Never immerse the adapter and sensor in water nor rinse it under the tap.

# WARNING

## **To reduce the risk of burns, fire, electric shock, or injury to persons:**

- Only use this product for its intended purpose as shown in this user manual.
- Do not place the strap on damaged skin.
- This device is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the device by a person responsible for their safety.
- Do not let children play with the Snoring Relief Band, the strap, the adapter or the USB cable as these may present a strangulation hazard.
- Always check the Snoring Relief Band before you use it. Do not use the Snoring Relief Band or adapter if it is damaged, became wet or hot as this may cause injury.
- Be aware of the kinetic energy stored in elastic strap when stretched and released suddenly, it may cause bruises.
- Only use the adapter WAA1001 which is provided with the Snoring Relief Band for charging the battery. Using

another adapter may jeopardize system safety.

- If the adapter and/or USB cord is damaged or defective, always seek a replacement from Philips in order to avoid a hazard (see 'Accessories').
- The sensor contains a non-replaceable Lithium battery.
- Do not open the sensor, there are no serviceable parts inside. It will void warranty and may induce a battery related hazard.
- Do not dispose of batteries in fire or extreme heat. Batteries may explode or leak.

## CAUTION

- The sensor is not washable. Never immerse the sensor in water or rinse it under the tap.
- Never use compressed air, scouring pads, abrasive cleaning agents or aggressive liquids such as petrol or acetone to clean the Snoring Relief Band.
- Check if the voltage indicated on the adapter corresponds to the local mains voltage before you connect the adapter.
- Do not use the adapter in or near wall sockets that contain or have contained an electric air freshener to prevent irreparable damage to the adapter.
- Use and store the product under the right environmental conditions as indicated in this user manual (see 'Specifications '). If not, this could affect the performance and lifetime of the product.
- The strap should feel snug but comfortable below your chest.
- Reset the sensor if you intend to change the user as the sensor has a self-learning algorithm which is optimized to your personal sleep behavior.

## General


- For hygienic reasons, the strap can be washed and air-dried. Tumble dry and/or excessive washing will deteriorate the aesthetic appearance of the strap.
- For optimal battery performance of the product life time, charge the battery when the battery low indicator blinks orange.
- Do not store the sensor when the battery is empty. Charge the sensor for two hours and then store it at room temperature to ensure optimal optimal battery performance over the product life time.
- After charging, remove the adapter from the wall socket to preserve energy.
- After the product has reached end-of-life, dispose the sensor battery as described in chapter 'recycling' to preserve our environment.

# SAVE THESE INSTRUCTIONS


## Electromagnetic fields (EMF)

This Philips product complies with all applicable standards and regulations regarding exposure to electromagnetic fields.

## Preparing for use

The Snoring Relief Band consists of a sensor and a strap. Wear the Snoring Relief Band when you sleep to reduce snoring caused by sleeping on your back. The sensor will gently vibrate to encourage you to turn from your back to your side.


## Charging

Charge the sensor before first time use and when the battery low indicator blinks orange.

Charging takes up to 8 hours. A fully charged sensor has an operating time of up to 8 sleep sessions.

### Charging with the adapter

Tip: Charge the sensor when it is in the pouch of the strap.

- 1 Insert the small plug of the USB cable into the socket on the back of the sensor.
- 2 Insert the large plug of the USB cable into the socket of the adapter.
- 3 Insert the adapter into the wall socket.  
Note: You cannot use the sensor while charging.

Note: Put the sensor on a flat surface while charging.

Note: The USB socket on the sensor can only be used for charging.

Note: During charging, the battery indicator blinks slowly white.

### **Battery fully charged**

When the battery is fully charged, the battery indicator lights up white continuously and switches off automatically after 30 minutes.

### **Quick charge**

When the battery is low or empty, you can quickly charge the sensor for 1 hour. This allows you to use the sensor during one sleeping session.

### **Battery low**


The sensor can only be used for one more sleep session.

### **Battery empty**


When the battery is almost empty, the battery indicator blinks orange quickly. Charge the sensor when the battery low indication lights up before using it.

## Inserting the sensor into the pouch

- 1 Push open the slit on the back of the pouch.


- 2 Insert the sensor into the pouch.

Note: Make sure that the sensor is completely inserted into the pouch.


- 3 The display and the activation button of the sensor should be visible through the pouch windows.


## Wearing the Snoring Relief Band during the night

- 1 Center the pouch with the sensor in front of your torso, just below the chest.


- 2 Do not wear the strap too high or too low.

- 3 Fully push the clip into the loop to securely fasten the strap.

Note: Make sure that the activation button is pointing towards your chin and the display is visible through the window of the pouch.

Note: For best results, wear the Snoring Relief Band on bare skin. In case you use body lotions make sure that these are absorbed by the skin before assembling the strap.

- 4 Adjust the length of the strap with the strap adjuster. The strap should feel snug but comfortable below your chest. The Snoring Relief Band should stay in the same position during the night. The strap is tight enough if you are able to feel the vibrations.

## Using the product


### Starting a sleep session

When you are ready to go to sleep, put the strap with the sensor around your chest (see 'Wearing the Snoring Relief Band during the night'). For men, the strap should be worn below the upper chest. For women, it should be worn below the breasts.

#### **You need to manually start a sleep session.**

Start a sleep session when you are in one of the data screens or when the sensor is off. Press and hold the activation button for two seconds.

- The sensor will vibrate briefly and show alternating dashes '- -' to indicate that the sleep session has started. These dashes are visible on the display during the sleep session.

#### **Notes:**


- Only start a sleep session when you have the intention to fall asleep (don't start a session when you first want to watch television).
- The sensor allows you to fall asleep in any position. This means that during

the first 30 minutes of a sleep session the sensor will not vibrate. Once the 30 minutes have passed, the sensor will vibrate when you are lying on your back.

- A sleep session can last up to 12 hours. After 12 hours the sensor automatically stops the sleep session and switches off.
- If a sleep session has been stopped and started again within a two hour timespan, both sessions are counted as one full sleep session.

## Automatic pause mode

A sleep session is automatically paused, when the sensor detects you sat up or got out of bed (vertical position). After you lie down again, the sensor will not vibrate for the first 30 minutes to allow you to fall asleep. Once the 30 minutes have passed, the sensor encourages you to turn from your back to your side.

## Temporarily interrupt vibrations


If the sleep session is active and you are awake and lying on your back, but you do not want the sensor to vibrate at that moment, tilt the sensor to a vertical position and back to start the automatic pause mode. The pause mode is not shown in the sensor data.

## Stopping a sleep session

**You need to manually stop the sleep session when you wake up and get out of bed.**

To stop a sleep session, press and hold the activation button for two seconds. If you do not stop the sleep session, the sensor keeps vibrating when the Snoring Relief Band is taken off and the sensor is placed on a horizontal surface. Not stopping a session also influences the sensor data.

After the sleep session has ended, the sensor will vibrate briefly and show the data of your last sleep session.

## Sensor data

The display shows three types of sensor data.

### Response rate:

How many times (in percentage) you turned from your back to your side after the sensor vibrates.

### Length of sleep session:


How many hours you spent in bed.

### Hours on your back:

How many hours you spent on your back.


## Display

- 1 Press the activation button briefly to switch on the sensor.


- 2 The sensor switches on and quickly shows an animation of the sensor data and battery indicator.

- 3 Sensor data point 1:  Response rate. The response rate indicates how many times (in percentage) you turned from your back to your side after the sensor vibrates.

Press the activation button briefly again to see the length of your sleep session.

- 4** Sensor data point 2:  Length of sleep session. The sensor shows how many hours you spent in bed.

Press the activation button briefly again to see the hours spent on your back.

- 5** Sensor data point 3:  Hours on your back.

The display shows the number of hours you have spent on your back.

**Note:** The sensor will only save the data of your last sleep session. If you use the sensor for the first time the display shows '-.-' to indicate that there is no data available yet.

**Note:** When the sensor shows one of the data screens, the sensor switches off automatically after eight seconds of inactivity.

## Error code

The sensor may show A1 or A2 on the display. You do not need to do anything. The code informs you why the sensor does not show data.

A1 indicates a session of less than one hour. Sessions shorter than one hour are not saved as a sleep session. The sensor will not show any data for response rate and the hours on your back.

A2 indicates a session of more than twelve hours. After twelve hours, the sensor automatically stops the sleep session and switches off. The sensor will not show any data for response rate and the hours on your back.



# Learning more about your Snoring Relief Band

## Adaptation phase

You may need some time to get used to the Snoring Relief Band. To help you get used to the vibrations, the sensor starts with an adaptation phase that consists of 7 sleep sessions.

During this phase, the Snoring Relief Band does not vibrate every time you lie on your back. With each next sleep session it vibrates more often. From the eighth session onwards, the sensor always encourages you to turn from your back to your side.

The sensor shows a small dot on the left side of the display to indicate that the adaptation phase is active.

## Vibration level

The sensor has a self-learning algorithm. The sensor optimizes the vibration level based on your reaction time.

## Vibration patterns

Every five days, the vibration pattern changes to make sure that you do not get used to the vibrations and that you will keep responding to them.

## Cleaning

### Cleaning the sensor

**Warning: Do not immerse the sensor in water and do not rinse it under the tap.**

**Warning: Never clean the sensor with abrasive cleansing products or (aggressive) liquids.**

Clean the sensor with a dry cloth after use, if needed.

### Cleaning the strap

**Warning: Always remove the sensor from the strap before you put the strap in the washing machine.**

- 1 You can wash the strap in the washing machine according to the washing instructions.

**Caution:** To prevent damage to the strap and other clothes in the washing machine, put the strap in a fine-mesh washing bag or zippered/tied pillowcase.

Note: Select a gentle wash cycle and a maximum temperature of 30 °C (86 °F).

**2** Let the strap air-dry.

- Do not tumble dry the strap
- Do not iron the strap.
- Do not stretch the strap excessively.

## Storage

Store the sensor and strap in a clean, dry place at room temperature.

## Assistance

For assistance, visit our website:  
[www.philips.com/support](http://www.philips.com/support).

## Troubleshooting

This chapter summarizes the most common problems you could encounter with the product. If you are unable to solve the problem with the information below, visit [www.philips.com/support](http://www.philips.com/support).

Problem	Possible cause / solution
The Snoring Relief Band doesn't work for me.	You do not only snore when you are sleeping on your back, but also when you sleep on your side and/or stomach. The sensor is most effective for people who ONLY snore when sleeping on their back. Use the Snoring Relief Band for at least 4 weeks to determine if it is effective for you. If you sleep on your stomach, the strap might feel uncomfortable because the sensor needs to be on your front to function properly.

Problem	Possible cause / solution
I'm lying on my back, but the sensor does not vibrate.	<p>When the adaptation phase is active (session 1-7), it may take longer than 30 minutes before the sensor vibrates. The sensor slowly builds up the number of vibrations each night. During this phase the sensor will not vibrate every time you sleep on your back (see 'Adaptation phase').</p> <p>Check if a dot is shown on the display, indicating the active adaptation phase active (see 'Adaptation phase'). From the eighth session onwards, the sensor always encourages you to turn from your back to your side.</p>
	<p>The sensor is switched off or the sleep session is not active.</p> <p>Start a session when when the sensor is off or if you are in one of the data screens. Press and hold the activation button for two seconds. The sensor vibrates briefly and shows '- -' on the display to indicate that the sleep session has started (see 'Starting a sleep session').</p>

Problem	Possible cause / solution
The sensor vibrates when I'm lying on my back awake, but I am not snoring.	The sensor is designed for people who snore sleeping on their back. It doesn't record your snoring, but it measures your position. Try to fall asleep on your side, so the sensor doesn't prevent you to fall asleep due to the vibrations.
Problem	Possible cause / solution
I don't feel the vibrations	<p>The Snoring Relief band is wrongly positioned (e.g. on your back or too low/high on your chest) (see 'Wearing the Snoring Relief Band during the night').</p> <p>Center the pouch with the sensor in front of your torso, just <b>below the chest</b>. Close the strap with the clip and adjust the strap so it's snug, but comfortable around your torso.</p> <hr/> <p>The sensor is not positioned correctly in the pouch (backwards or upside down) (see 'Inserting the sensor into the pouch').</p> <p>Make sure that the sensor is fully inserted into the pouch, with the activation button pointing towards your chin. The display should be visible through the window of the pouch.</p>

<b>Problem</b>	<b>Possible cause / solution</b>
The sensor data shows incorrect information.	The battery was almost empty when you started the sleep session. If the battery runs out during a sleep session, the sensor switches off and stops collecting data.
	The sleep session lasted more than 12 hours and automatically stopped (to save battery). The display shows error code A2 (see 'Error code'). The sleep session has exceeded the maximum of 12 hours or was not stopped properly. Make sure you manually stop a session (see 'Stopping a sleep session').
<b>Problem</b>	<b>Possible cause / solution</b>
The sensor data shows incorrect information.	Within two hours of stopping a sleep session, you started another sleep session. If you start another sleep session within two hours after switching off the previous session, the sensor will register these two sessions as one single session.
	The sleep session was shorter than one hour. The display shows no data for response rate and hours on your back. It shows 'A1' for sessions shorter than 1 hour (see 'Error code').

Problem	Possible cause / solution
The reset does not work.	<p>The long press was too short (&lt;5 seconds). Press and hold the activation button for five seconds, until you see the countdown (3-2-1) on the display and feel a vibration. Keep pressing even when you see '- -' on the display (after 2 seconds of long pressing). The final vibration indicates that the reset was successful. The adaptation phase is now active again.</p> <hr/> <p>The sensor is off or the sensor is in a sleep session '- -'. The long press (five seconds) only works when the display is active and shows one of the sensor data screens (see 'Sensor data').</p>

## Full Two-Year Warranty

Philips Japan warrant this new product against defects in materials or workmanship for a period of two years from the date of purchase, and agree to repair or replace any defective product without charge.

**IMPORTANT:** This warranty does not cover damage resulting from accident, misuse or abuse, lack of reasonable care, or the affixing of any attachment not provided with the product or loss of parts or subjecting the product to any but the specified voltage or batteries.\*

**NO RESPONSIBILITY IS ASSUMED FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES.**

In order to obtain warranty service, go to **[www.philips.com/support](http://www.philips.com/support)** for assistance. It is suggested that for your protection you return shipments of product by insured mail, insurance prepaid. Damage occurring during shipment is not covered by this warranty.

**NOTE:** No other warranty, written or oral, is authorized by Philips Japan, Ltd.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state, province to province or country to country. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above exclusion and limitations may not apply to you.

\* Read enclosed instructions carefully.

PHILIPS and Philips shield are registered trademarks of Koninklijke Philips N.V.

©2020 Philips Japan, Ltd. All rights reserved.

## Disposal

- This product contains a rechargeable lithium-ion battery which must be disposed of properly.
- Contact your local town or city officials for battery disposal information.
- For assistance, visit our website **[www.philips.com/support](http://www.philips.com/support)**.
- Your product is designed and manufactured with high quality materials and components, which can be recycled and reused. For recycling information, please contact your local waste management facilities or visit **[www.recycle.philips.com](http://www.recycle.philips.com)**.

## Safety and compliance

The Snoring Relief Band meets the safety requirements in the USA and Canada:

- UL/CSA 62368: Audio/video, information and communication technology equipment
- UL1642: Batteries
- UL 1310, CAN/CSA C22.2 No 223-M91: Adapters

## FCC Information

- This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
  - Reorient or relocate the receiving antenna.
  - Increase the separation between the equipment and the receiver.
  - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
  - Consult the dealer or an experienced radio/TV technician for help.

- This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
- This Class B digital apparatus complies with Canadian ICES-003.

## Specifications

### Adapter Philips WAA1001

Safety classification	Class II power supply (double insulation)
Mode of operation	Continuous
Ingress Protection marking	IPX4 (protection against splashing water)
Input characteristics	100-240 Vac, 50/60 Hz, 3.5 W
No-load power consumption	0.03 W
Output characteristics	5.0 Vdc, 2.5 W, 0.5 A
Weight (EU / US)	39 gr (1.38 oz.) / 41 gr (1.44 oz.)
Length of detachable USB cord	60 cm (23.6 in.)

### Sensor

Safety classification	Class III device (sensor battery is charged by a Class II adapter)
Dimensions Sensor (height x width x Depth)	Approx. 79 x 43 x 13 mm (3.1 x 1.7 x 0.5 in.)
Weight Sensor	Approx. 29 gr (1 oz.)
Ingress Protection marking	IP 51
Internal non-replaceable battery	Lithium Polymer 3.7 V, 300 mAh, UL1642

### Strap

Material strap	See textile label on the strap for materials used.
Material hook / slider	Coated steel
Length (min - max)	Approx. (72 cm - 130 cm) / (28 in. - 51 in.)
Weight	Approx. 33 gr (1.16 oz.)

### Operating conditions

Temperature	+5 °C to +40 °C (41 °F to 104 °F)
Relative humidity	15 % to 95 % (no condensation)
Atmospheric pressure	86 - 106 kPa

## Storage conditions

Short term (<1 month ) storage temperature	-20 °C to +60 °C (-4 °F to +140 °F)
Long term (>1 month) storage temperature	+15 °C to +25 °C (+59 °F to +77 °F)
Relative humidity	≤85 % (no condensation)
Atmospheric pressure	86 - 106 kPa

## Explanation of symbols

Symbol	Explanation
	Manufacturer logo.
	UL listing mark for Canada and the United States.
	For indoor use only.
	Do not dispose of the product with normal household waste.
<b>IP51</b>	Protected against dust and falling drops of water.
	Refers to the user manual. Read the user manual before you use the product.
	DC power.
<b>SN3710</b>	Model number.

	Serial number.
	Detachable supply unit.
	Stand-by.
$t_a=40^{\circ}\text{C}$	Maximum ambient temperature.
	Efficiency class.
	SMPS (Switch Mode Power Supply).
	SMPS incorporating a short-circuit-proof safety isolating transformer.
<b>IPX4</b>	Ingress protection against water splashes from any direction.
 	Machine wash in water not exceeding 30 °C / 86 °F, at delicate/gentle setting.
 	Do not bleach.
 	Do not tumble dry.
 	Do not iron.
 	Do not dry clean.

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