**PHILIPS**

*SmartSleep*

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## Welcome

**Congratulations** on the decision to commit yourself to better sleep for you and your partner. Many people get concerned because of the effects of snoring on those around them. Philips is here to help.

**The device has proven quite successful** but some people report **it takes time to adjust** to the vibrations and new sleeping position. **This adjustment period is normal.**

**Be patient** and don't give up. Remind yourself you are working toward **quieter nights** for you and your partner.

Sincerely,

A handwritten signature in black ink, appearing to read 'Mark Aloia', with a stylized, cursive script.

Mark Aloia, Ph.D.  
Sleep Expert at Philips

## How it works


The Philips Snoring Relief Band is for **people who snore sleeping on their back**. Research shows that **sleeping on your side reduces snoring**.

Worn around the chest, the sensor uses **clinically proven** technology to detect when you're laying on your back and **vibrates to prompt you to roll to your side**.


**86% of users reported** reduced snoring after four weeks of use.

## Set up strap and sensor

- **Charge the sensor before you use it for your first week.** A full charge takes eight hours.
- **Insert the sensor into the gray pouch** that is part of the Snoring Relief Band.
- **Put on the strap** just **below your chest** and securely fasten the strap. The strap should be tight enough so you feel the vibrations and the strap does not move during the night. The strap should not be so tight that it hurts.


- **Press and hold the activation button for two seconds to start the sleep session.** The sensor will vibrate briefly and show '--'.

- Remember to **turn off the sensor after every sleep session.** Press and hold the button again for **two seconds.** The sensor will **vibrate** briefly and show the data.

## Your first nights

During the first week, the number of **vibrations starts low and increases** over time. **From the eighth session onwards, the sensor always encourages you to turn from your back to your side.**

A built-in **30-minute timer** allows **you to fall asleep** in your favorite position – back, stomach, or side – before the sensor starts vibrating.

30 Minutes after you start a sleep session, **the sensor will be activated.**

**Vibrations will stop** when you roll onto your side.


During the night, you are in control with the **Pause Function**: Flip the sensor 90 degrees to **pause the vibrations for 30 minutes.**

If you sleep on your stomach, the strap might feel uncomfortable because the sensor needs to be on your front to function properly.

# Keep track of **your progress**

Each night, the sensor collects data on **your sleep**. Those values give you **insight into your progress**, and how much your body is adjusting to the new sleep routine.

Press the activation button briefly to show **three types of data** collected that night.

## Response rate in %

The number increases as you roll more often onto your side after the sensor vibrates.

## Length of sleep session

How many hours from start to finish of your sleep session.

## Hours on your back

The more you roll over, the less time you lay on your back.

Chart your data each morning to see your progress.

	Response rate	Length of sleep session	Hours on your back
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
11			
12			
13			
14			
15			
16			
17			
18			
19			
20			
21			
22			
23			
24			
25			
26			
27			
28			

# Review the process

## Week 1

### The **initial adaptation phase**

The number of vibrations starts low and increases each night over the course of the week. If the vibrations feel too strong, try moving the strap a little up or down to increase the wearing comfort, or try wearing the strap over a t-shirt.

Ask yourself, do the vibrations nudge you to turn on your side?

## Week 2

### **Rolling to my side** when prompted

From this week onwards, the sensor starts the normal vibration phase and encourages you to turn from your back to your side. Remember to use the pause function if you need a break in vibrations.

Has your partner noticed a change?

## Week 3

### You are learning your **new routine**

Rolling over when the sensor vibrates feels easy and goes mostly unnoticed.

Ask yourself how you incorporated the vibrations and rolling over this week?

## Week 4

### **Put your snoring to rest, for good**

You should feel comfortable rolling onto your side as part of your sleep routine.

Are you and your partner enjoying a better night's sleep?

The notes section at the end of this welcome guide is to write you and your partner's additional motivations and observations.

## Notes

[illegible]

## Notes

This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

© 2020 Koninklijke Philips N.V.  
All rights reserved /  
3000.056.2731.2 (28/10/2020)