PHILIPS


AVENT

Single/Double
Electric Breast
Pump

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



English 6
Français 32
العربية 60

Contents

Introduction	6
General description	7
Intended use and indications for use	7
Important safety information	7
Before first use	9
Cleaning and disinfecting	10
Using the breast pump	13
Using the breast pump when you have symptoms of engorgement	19
After use	19
Feeding breast milk with the bottle	20
Compatibility	22
Ordering accessories and spare parts	22
Recycling	22
Battery removal instructions	23
Warranty and support	23
Troubleshooting	23
Undesirable side effects	25
Storage conditions	26
Electromagnetic fields (EMF)	26
Technical information	26
Explanation of symbols	29

Introduction

Congratulations on your purchase and welcome to Philips Avent! To fully benefit from the support that Philips Avent offers, register your breast pump at **www.philips.com/mybreastpump**.

The Philips Avent single/double electric breast pump Premium is inspired by babies and their natural suction movement. The soft silicone cushion stimulates the nipple with gentle massaging movements, helping you to express milk comfortably and effectively. The combination of suction and nipple stimulation triggers fast milk flow. The no-lean-forward design enables you to sit in a relaxed, upright position while expressing.

Healthcare professionals and the World Health Organization recommend to exclusively feed children with breast milk during the first six months of life and to continue giving them breast milk for two years next to offering additional nutritious foods. This is because your breast milk is specially adapted to your baby's needs and helps to protect your baby against infection and allergies. To breastfeed longer, you can express and store your milk so that your baby can still enjoy the benefits of it, even if you cannot be there to provide it yourself. As the pump is compact, quiet and discrete to use, you can take it with you anywhere, allowing you to express milk at your own convenience and to maintain your milk supply.

If breastfeeding goes well, it is advisable (unless recommended otherwise by your healthcare professional) to wait until your milk supply and breastfeeding schedule are established (normally at least 2 to 4 weeks after giving birth) before you start expressing.

General description

Note: The numbers mentioned below refer to the numbers in the figure on the front foldout page of this user manual.

Note: The double breast pump comes with two pieces of the items shown in Figures B and C.

Product description (Fig. A)

A1 Motor unit

A2 Breast pump kit

A3 Bottle

A4 Adapter*

A5 Double breast pump

* The adapter supplied may vary for different countries.

Breast pump kit (Fig. B)

B1 Silicone tube and cap

B2 Silicone diaphragm

B3 Pump body

B4 Cover

B5 Cushion

B6 White valve

Philips Avent Natural bottle (Fig. C)

C1 Bottle cap

C2 Screw ring

C3 Teat

C4 Sealing disc

C5 Bottle

Accessories

There are multiple versions of this breast pump package, which all come with a different combination of the following Philips Avent accessories:

- Breast pads
- Travel bag
- Insulation pouch
- Cleaning brush
- Storage cups with adapter*
- Breast milk storage bags*
- Natural bottles

* For these accessories, a separate user manual is provided.

Intended use and indications for use

The Philips Avent single/double electric breast pump Premium is intended to express and collect milk from the breast of a lactating woman, and indicated to alleviate the symptoms of engorgement. The device is intended for a single user.

Important safety information

Read this user manual carefully before you use the breast pump and save it for future reference.

This user manual can also be found online at: www.philips.com/support.

This product is operated by the lactating woman herself.

You can use this breast pump in an everyday environment.

Contraindication

- Never use the breast pump while you are pregnant, as pumping can induce labor.

Warning

Warnings to avoid choking, strangulation and injury:

- Do not allow children or pets to play with the motor unit, the adapter, packaging materials or accessories.

- Always unplug the breast pump after use and when fully charged. Store the breast pump when not in use.

- This breast pump is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge. Such persons can only use this breast pump if they are supervised by or have received proper instruction concerning use of the appliance by a person responsible for their safety.

Warning to prevent transmitting infectious diseases:

- If you are a mother who is infected with Hepatitis B, Hepatitis C or Human Immunodeficiency Virus (HIV), expressing breast milk will not reduce or remove the risk of transmitting the virus to your baby through your breast milk.

Warnings to avoid electric shock:

- Inspect the breast pump, including the adapter, for signs of damage before each use. Do not use the breast pump if the adapter or plug is damaged, if it does not work properly or if it has been dropped or submerged in water.
- Only use the adapter (S009AHZ050yyyy) provided with this breast pump.
- Always ensure that it is easy to disconnect the device from the mains in case of malfunction. The adapter is considered the disconnection device.

Warnings to avoid contamination and to ensure hygiene:

- For hygienic reasons, the breast pump is only intended for repeated use by a single user.
- Clean and disinfect all parts before first use and after each use.
- Do not use antibacterial or abrasive cleaning agents to clean the breast pump parts as this may cause damage.
- Wash your hands thoroughly with soap and water before you touch breast pump parts and breasts to prevent contamination. Avoid touching the inside of containers or lids.

Warnings to avoid breast and nipple problems including pain:

- Do not use the breast pump if the silicone diaphragm appears to be damaged or broken. See 'Ordering accessories and spare parts' for information on how to obtain replacement parts.
- Only use accessories and parts recommended by Philips Avent.
- There are no user serviceable parts inside the motor unit of the breast pump.
- No modification of the breast pump is allowed. If you do so, your warranty becomes invalid.
- Never use the breast pump while you are sleepy or drowsy to avoid lack of attention during use.
- Always switch off the breast pump before you remove the pump body from your breast, in order to release the vacuum.
- If the vacuum is uncomfortable or causes pain, switch off and remove the breast pump from your breast.

- Do not continue pumping for more than five minutes at a time if you do not succeed in expressing any milk. Try to express at another time during the day.
- If the process becomes very uncomfortable or painful, stop using the breast pump and consult your healthcare professional.

Warnings to prevent damage and malfunction of the breast pump:

- Portable radio frequency (RF) communications equipment (including mobile phones and peripherals such as antenna cables and external antennas) should not be used closer than 30 cm (12 in) to any part of the breast pump, including adapter. This could have a negative effect on the performance of the breast pump.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation.

Warning to avoid increased electromagnetic emissions or decreased electromagnetic immunity of this device, resulting in improper operation:

- Do not use any other adapter or charging cable than supplied with the device.

Battery safety instructions

- Only use the product for its intended purpose. Read this important information carefully before you use the product and its batteries and accessories, and save it for future reference. Misuse can lead to hazards or serious injuries. Accessories supplied may vary for different products.
- To charge the battery, only use the detachable adapter (S009AHZ050yyy) provided with the product. Do not use an extension cord.
- Charge, use and store the product at the temperature ranges indicated in the technical specifications table in this user manual.
- Keep product and batteries away from fire and do not expose them to direct sunlight or high temperatures.
- If the product becomes abnormally hot or smelly, changes color or if charging takes longer than usual, stop using and charging the product and contact Philips.
- Do not place products and their batteries in microwave ovens or on induction cookers.
- Do not open, modify, pierce, damage or dismantle the product or battery to prevent batteries from heating up or releasing toxic or hazardous substances. Do not short-circuit, overcharge or reverse charge batteries.
- If batteries are damaged or leaking, avoid contact with the skin or eyes. If this occurs, immediately rinse well with water and seek medical care.

Caution

Cautions to prevent damage to and malfunctioning of the breast pump:

- Prevent the adapter and motor unit from coming into contact with water.
- Never put the motor unit or the adapter in water or in a dishwasher, as this causes permanent damage to these parts.
- Keep the adapter and the silicone tubes away from heated surfaces to avoid overheating and deformation of these parts.

Before first use


Disassemble, clean and disinfect all parts that come into contact with breast milk before you use the breast pump for the first time and after every use.

Fully charge the device before you use it for the first time or after a long period of non-use.

Cleaning and disinfecting

Overview


Clean and disinfect the parts that come into contact with the breast and breast milk as described below:

	When	How
	Clean and disinfect before first use and after every use.	Disassemble all parts and then clean them as described in chapter 'Cleaning' and disinfect them as described in chapter 'Disinfecting'.
		

The following parts do not come into contact with the breast and breast milk; clean them as described below:

	When	How
	Clean when needed.	Wipe with a clean, damp cloth with water and mild detergent.

Step 1: Disassembling

Disassemble the breast pump, bottle and storage cup completely.
Also remove the white valve from the pump body (Fig. B3).

Step 2: Cleaning

The parts that come in contact with milk can be cleaned manually or in the dishwasher.

⚠ Warning: Do not use antibacterial or abrasive cleaning agents to clean the breast pump parts as this may cause damage.


⚠ Caution: Never put the motor unit or the adapter in water or in a dishwasher, as this causes permanent damage to these parts.

⚠ Caution: Be careful when you remove and clean the white valve. If it gets damaged, your breast pump does not function properly. To remove the white valve, gently pull at the ribbed tab on the side of the valve. To clean the white valve, rub it gently between your fingers in warm water with some dishwashing liquid. Do not insert objects into the white valve, as this may cause damage.

Step 2A: Manual cleaning

Supplies needed:


- Mild dishwashing liquid
- Drinking-quality water
- Soft, clean brush
- Clean tea towel or drying rack
- Clean sink or bowl


1. Rinse all parts under a running tap with lukewarm water.

2. Soak all parts for 5 minutes in warm water with some mild dishwashing liquid.

3. Clean all parts with a cleaning brush in warm water with some mild dishwashing liquid. Remove any food residues. If you use a brush to clean the tip of the feeding teat, clean it as carefully as possible to avoid damage.

4. Thoroughly rinse all parts under a running tap with cold clear water.

5. Leave all parts to air-dry on a clean tea towel or drying rack.

Step 2B: Cleaning in the dishwasher

Supplies needed:

- Mild dishwashing liquid or a dishwashing tablet
- Drinking-quality water

Note: Food colorings may discolor parts.

1. Place all parts on the top rack of the dishwasher.

2. Put dishwashing liquid or a tablet in the machine and run a standard program.

3. Leave all parts to air-dry on a clean tea towel or drying rack.

Step 3: Disinfecting


Supplies needed:

- A household pot
- Drinking-quality water

⚠ Caution: During disinfection with boiling water, prevent the bottle or other parts from touching the sides of the pot. This can cause irreversible product deformation or damage that Philips cannot be held liable for.

1. Fill a household pot with enough water to cover all parts and place all parts in the pot. Let the water boil for 5 minutes. Make sure the parts do not touch the sides of the pot.

2. Allow the water to cool down.

3. Gently remove the parts from the water. Place the parts on a clean surface/tea towel or drying rack to air-dry.

4. Store the dry items safely in a clean, protected area until needed. Ensure that the clean parts have air-dried completely before storing them to help prevent germs and mold from growing.


Using the breast pump

Cushion size

The Philips Avent single/double electric breast pump Premium has a soft, active cushion. There is only one cushion size. It gently stimulates your nipple to trigger milk flow. The cushion is made of flexible silicone, fitting nipple sizes up to 30 mm.

Charging the breast pump


Fully charge the battery before you use the breast pump for the first time and when the battery is low. Fully charging the battery takes up to 4 hours.


1. Switch off the device.


2. Insert the adapter into the wall socket and insert the small plug at the other end into the motor unit.

3. The battery lights light up one by one and then all go out: the battery is charging.


4. All lights are steady on: the battery is fully charged.

5. Remove the adapter from the wall socket and pull the small plug out of the motor unit.

Battery indicator light description

All lights are steady on. The battery is fully charged.
Charged battery operates approximately 3 pumping sessions.*


Two lights are steady on.
The battery is charged sufficiently for maximum two sessions.*

First light flashes slowly.
The battery is charged sufficiently for maximum one session.*

First light flashes quickly.
The battery is very low, the device will soon stop operating, charging is needed.

Battery lights light up one by one and go out all together.
The battery is charging.

* One session is the equivalent of 15 minutes (double electric breast pump) or 30 minutes (single electric breast pump).

Assembling the breast pump

Make sure you have cleaned and disinfected the parts of the breast pump that come into contact with milk.

⚠ Warning: Be careful, when you have disinfected the breast pump parts by boiling them, they can be very hot. To prevent burns, only start assembling the breast pump when disinfected parts have cooled down.

⚠ Warning: Wash your hands thoroughly with soap and water before you touch breast pump parts and breasts to prevent contamination. Avoid touching the inside of containers or lids.

Check the breast pump kit parts for wear or damage before use and replace them if necessary.

📌 Note: Correct placement of the cushion, silicone diaphragm and silicone tube is essential for the breast pump to form a proper vacuum.

1. Wash your hands thoroughly with soap and water.

2. Push the white valve in the pump body as far as possible.

3. Screw the pump body onto the bottle.


4. Place the cushion in the pump body and make sure that the rim covers the pump body.

4A. Push the inner part of the cushion into the funnel against the line (indicated with an arrow).


5. Place the silicone diaphragm in the pump body.


5A. Make sure that the silicone diaphragm fits securely around the rim by pressing it down with your thumbs.


6. Attach the silicone tube and cap onto the silicone diaphragm. Push down the cap until it is securely in place.

7. Connect the tube(s) to the motor unit.

8A. When the breast pump is charged, you can use it cordless. Press and hold the on/off button until the device switches on.

8B. You can also use the breast pump connected to the wall socket. If the battery is low, insert the adapter into the wall socket and insert the small plug at the other end into the motor unit.


9. The breast pump is now ready for use.


① Note: You can place the cover over the cushion to keep the breast pump clean when you take it with you and while you are preparing for expressing.

Motor unit part description

- 1 USB power inlet
- 2 On/off button with pause/play function
- 3 Mode selection button
- 4 'Level down' button
- 5 'Level up' button
- 6 LED display
- 7 Tubing port for single breast pump
- 8 Tubing ports for double breast pump

LED display description

- 1 Battery indicator lights
- 2 Timer (minutes) and vacuum level indication
- 3 Expression icon
- 4 Stimulation icon

Explanation of modes

The breast pump has two modes. Below you can find an explanation of these modes.

Mode icons	Mode	Explanation	No. of suction levels
	Stimulation mode 	Mode to stimulate the breast to start up the milk flow.	8 suction levels
	Expression mode 	Mode for efficient milk removal after milk has started flowing. After 90 sec. of stimulation, the device automatically switches to the expression mode.	16 suction levels
		The number on the display shows the running time in minutes.	
		When you change the suction level, the display shows the selected suction level for a few seconds and then shows the running time again.	

Pumping instructions

⚠ Warning: Always switch off the breast pump before you remove the pump body from your breast, in order to release the vacuum.

⚠ Warning: Do not continue pumping for more than five minutes at a time if you do not succeed in expressing any milk. Try to express at another time during the day.

⚠ Warning: If the process becomes very uncomfortable or painful, stop using the breast pump and consult your healthcare professional.

📌 Note: If you regularly express more than 125 ml/4 fl oz per session, you can purchase and use a 260 ml/9 fl oz Philips Avent Natural bottle to prevent overfilling and spillage.


📌 Note: The breast pump can be used cordless if it is sufficiently charged, but you can also use it connected to the wall socket.

1. Wash your hands thoroughly with soap and water and make sure that your breast is clean.

2. Place the breast pump kit onto your breast. Make sure that you center your nipple.


3. Press and hold the on/off button until the device switches on. The indicator lights go on to indicate that the device is in stimulation mode.

4. The breast pump starts in the stimulation mode and slowly increases its vacuum to the last used stimulation suction level (or level 5 when used for the first time).


5. The display shows the running time.


6. Use the 'level up' and 'level down' buttons to change the suction level to your preference.


7. After 90 seconds, the breast pump will switch smoothly to expression mode and slowly increase its vacuum to the last used expression suction level (or level 11 when used for the first time).


8. If your milk starts flowing earlier or if you want to switch from expression to stimulation modes with the mode selection button.

9. Briefly press the on/off button if you want to pause. To continue, briefly press the on/off button again.

10. Press and hold the on/off button to switch off.

11. Unplug the adapter from the wall socket and pull the small plug out of the motor unit.

① Tip: For optimal expression, choose the highest setting that still feels comfortable. This setting may differ across different sessions.

① Tip: You can use the double breast pump as a single breast pump by connecting only the tube of one pump kit to the motor unit.

Note: The breast pump will switch off automatically after 30 minutes in order to prevent overheating.

Using the breast pump when you have symptoms of engorgement


Breast engorgement means your breasts are overfilled with milk. This is a painful and unpleasant condition. Frequent extraction of milk and completely emptying the breasts of milk can reduce engorgement. When you have symptoms of engorgement, make sure you use a lower vacuum level to minimize stretching of breast tissue.

Note: When symptoms of engorgement are present, make sure you empty your breasts completely during each expression session.


After use

1. Unscrew the bottle.


2A. To store the breast milk: close the bottle with the sealing disc and screw ring.

2B. To prepare the bottle for feeding your baby: place a teat in a screw ring and screw the assembled screw ring onto the bottle. Seal the teat with the bottle cap (see 'Assembling the bottle').

3. Remove the silicone tube and cap from the breast pump kit.

4. Disassemble and clean all parts that have been in contact with breast milk (see 'Cleaning and disinfecting').

5. For easy storage, wrap the silicone tube around the motor unit and clip the cap onto the tube.

Warning: Refrigerate or freeze expressed breast milk immediately or keep it at a temperature of 16-29 °C for a maximum of 4 hours before you feed your baby.

 **Note:** For more information about the breast pump and tips for expressing milk, visit www.philips.com/support.

Storing breast milk

Below you find guidelines for storing breast milk:

Storage location	Temperature	Max. storage duration
Room	16-29 °C	4 hours
Refrigerator	4 °C	4 days
Freezer	< -4 °C	6 months

 **Warning:** Never refreeze thawed breast milk to avoid a decrease in milk quality.

 **Warning:** Never add fresh breast milk to frozen breast milk to avoid a decrease in milk quality and unintended thawing of the frozen milk.

Feeding breast milk with the bottle

For your child’s safety and health

WARNING!

- Continuous and prolonged sucking of fluids will cause tooth decay.
- Always check food temperature before feeding.
- Throw away at the first signs of damage or weakness.
- Keep components not in use out of the reach of children.
- Never attach to cords, ribbons, laces or loose parts of clothing. The child can be strangled.
- Never use feeding teats as a soother.
- Always use this product with adult supervision.
- Keep the bottle cap away from children to avoid suffocation.
- Do not allow children to play with small parts or walk/run while using bottles.
- Inspect all parts before use and pull the feeding teat in all directions. Throw away any part at the first signs of damage or weakness.
- Always discard any breast milk that is left over at the end of a feed.

Caution: To avoid damage to the bottle


- Do not place in a heated oven, plastic can melt.
- Plastic material properties may be affected by disinfection and high temperatures. This can affect the fit of the bottle cap.
- Do not leave a feeding teat in direct sunlight or heat, or leave in disinfectant for longer than recommended because this may damage the product.

Before using the bottle

- Inspect the bottle and feeding teat before each use and pull the feeding teat in all directions to prevent a choking hazard. Throw away at the first signs of damage or weakness.
- Clean and disinfect the bottle before first use and after every use (see 'Cleaning and disinfecting').


Assembling the bottle

1. The teat is easier to assemble if you move it back-and-forth while pulling it upwards instead of pulling it up in a straight line.

2. Make sure you pull the teat through until its lower part is aligned with the screw ring.

3. Make sure you place the cap vertically onto the bottle so that the teat sits upright.

4. Screw the screw ring with teat and cap on the bottle.

5. To remove the cap, place your hand over the cap and your thumb in the dent of the cap.

Heating breast milk

If you use frozen breast milk, let it defrost completely before you heat it.

Note: If you need to feed your baby quickly, you can defrost the milk in a bowl of hot water.

Heat the bottle with defrosted or refrigerated breast milk in a bowl of hot water or in a bottle warmer. Remove the screw ring and sealing disc from the bottle.

Warning: It is not recommended to use a microwave to warm up breast milk. Microwaves might alter the quality of breast milk and destroy some valuable nutrients and may produce localized high temperatures. If you heat up breast milk in the microwave, only place the container without screw ring, teat and cap in microwave. Always stir heated breast milk to ensure even heat distribution and check the temperature before serving.

Choosing the right teat for your baby

Philips Avent Natural Response teats are available with different flow rates. The breast pump comes with teat number 2. This is a 0m+ teat, which is recommended to start with for a newborn.

Use a higher flow rate if your baby falls asleep during feeding, gets frustrated, or plays with the teat instead of drinking. Use a lower flow rate if your baby is gulping the milk down or if milk drips out of its mouth during feeding.

Philips Avent Natural Response teats are clearly numbered on the side, to indicate flow rate. For more information, see www.philips.com/avent.

Storing the bottles

- Store all parts in a dry, clean and covered container.
- Do not leave a feeding teat in direct sunlight or heat, or leave in disinfectant for longer than recommended because this may damage the product.
- For hygiene reasons, we recommend replacing the teats after 3 months.

Compatibility

The Philips Avent single/double electric breast pump Premium is compatible with all baby bottles in the Philips Avent Natural range and Philips Avent storage cups with adapter. We do not recommend using the Philips Avent Natural glass bottles with the breast pump.

Ordering accessories and spare parts

To buy accessories or spare parts, visit www.philips.com/support or go to your Philips dealer. You can also contact the Philips Consumer Care Center in your country.

When you order spare teats, make sure that you use a teat with the correct flow rate for your baby (see 'Choosing the right teat for your baby') and do not mix Philips Avent Anti-colic bottle parts and teats with parts of the Philips Avent Natural bottles. They do not fit and could cause leakage or other issues.

Recycling

- Do not throw away the product with the normal household waste at the end of its life, but hand it in at an official collection point for recycling. By doing this, you help to preserve the environment.
- This product contains a built-in rechargeable battery which shall not be disposed of with normal household waste. Please take your product to an official collection point or a Philips service center to have a professional remove the rechargeable battery.
- Follow your country's rules for the separate collection of electrical and electronic products and rechargeable batteries. Correct disposal helps prevent negative consequences for the environment and human health.

Battery removal instructions


- When you handle batteries, make sure that your hands, the product and the batteries are dry.
- To avoid accidental short-circuiting of batteries after removal, do not let battery terminals come into contact with metal objects (e.g. coins, hairpins, rings). Do not wrap batteries in aluminum foil. Tape battery terminals or put batteries in a plastic bag before you discard them.

ⓘ Note: Only remove the rechargeable battery when you discard the product. Before you remove the battery, make sure that the product is disconnected from the wall socket and that the battery is completely empty.


ⓘ Note: Take any necessary safety precautions when you handle tools to open the product and when you dispose of the rechargeable battery.


1. Remove the rubber feet from the bottom of the motor unit.


2. Remove the screw with a torx screwdriver.


3. Open the motor unit by pressing at the sides of the device.

4. Disconnect the connectors from the printed circuit board.

5. Remove the frame from the motor unit.

6. Cut or fold the edges of the plastic battery frame. Remove the battery from the plastic battery frame.

ⓘ Note: Dispose of the rechargeable battery at a separate local collection point. Other parts can be disposed at a specific place for disposal of electric equipment.

Warranty and support

If you need information or support, please visit www.philips.com/support or read the international warranty leaflet.

Troubleshooting

This chapter summarizes the most common problems you could encounter with the breast pump. If you are unable to solve the problem with the information below, visit www.philips.com/support for a list of frequently asked questions or contact the Consumer Care Center in your country.

Problem	Solution
I experience pain when I use the breast pump.	- Choose a suction level that is comfortable for you. - Check if the silicone diaphragm is not damaged (does not have small cracks, holes etc.). - Make sure that your nipple is placed in the center of the cushion. - If the pain persists, stop using the breast pump and consult your healthcare professional.
The breast pump does not switch on while it is connected to the wall socket.	- Only use the adapter supplied. - Make sure that the adapter has been inserted correctly into the wall socket and that the small plug has been inserted correctly into the motor unit. - Connect the adapter to another wall socket. Then press and hold the on/off button to check if the breast pump can be switched on. - If the problem persists, contact the Consumer Care Center via www.philips.com/support.
The breast pump does not switch on (indicator lights do not go on).	- Make sure that the battery is charged (see 'Charging the breast pump') or connect the breast pump to the wall socket.
The breast pump is not charging.	- Only use the adapter supplied. - Make sure that the adapter has been inserted correctly into the wall socket and that the small plug has been inserted correctly into the motor unit. - The breast pump is too cold or too hot for charging. Let the breast pump acclimatize to room temperature for 30 minutes. - If the problem persists, contact the Consumer Care Center via www.philips.com/support.
Charging the breast pump takes a very long time.	- Only use the adapter supplied. - Make sure that the adapter has been inserted correctly into the wall socket and that the small plug has been inserted correctly into the motor unit. - The breast pump is too cold or too hot. It is recommended to charge the breast pump at room temperature. - If the problem persists, contact the Consumer Care Center via www.philips.com/support.
I do not feel any suction/the suction level is too low.	- If possible, try a higher suction level. - Check if the silicone diaphragm is assembled correctly (see 'Assembling the breast pump'). - Check if the white valve is not damaged and is assembled correctly (see 'Assembling the breast pump'). - Check if the other parts of breast pump are intact and assembled correctly (see 'Assembling the breast pump'). - Make sure that the breast pump is positioned correctly on the breast to allow a seal between breast pump and breast. - If you still do not feel any suction, please contact the Consumer Care Center via www.philips.com/support.

Problem	Solution
The suction of the breast pump is too strong.	- If possible, try a lower suction level. When you express for the first time, start with the default suction level (level 5 in stimulation mode, level 11 in expression mode) and increase/decrease the level when needed. During different expression sessions, different levels may feel most comfortable for you. - Make sure that you only use the breast pump parts provided with this breast pump. - Check if the silicone diaphragm is not damaged (does not have small cracks, holes etc.). - If the problem persists, contact the Consumer Care Center via www.philips.com/support.
I express little/no milk when I use my breast pump.	- If you do not feel enough suction, please check Troubleshooting item 'I do not feel any suction/the suction level is too low'. - If you are using the Philips Avent breast pump for the first time, you may need some practice before you are able to express any milk. More expressing tips can be found on the Philips website www.philips.com/mybreastpump.
When I switch on the breast pump after storage, it does not start working immediately.	- Make sure that the battery is fully charged (see 'Charging the breast pump') or connect the breast pump to the wall socket. - Perhaps you stored the breast pump at a temperature below the minimum storage temperature or above the maximum storage temperature. Let the breast pump acclimatize in a room of 20 °C for 30 minutes to adjust its temperature to the operating range (5 °C to 40 °C).
The display shows 'Er' and the device switches off automatically.	- If you are operating the breast pump on battery power, let it cool down for 30 minutes and try to switch it on again. - If you are operating the breast pump with the adapter connected to the wall socket, check if you are using the adapter that was supplied with the device. - If the problem persists, contact the Consumer Care Center via www.philips.com/support.

Undesirable side effects

Any serious incident that has occurred in relation to this device should be reported to Philips via **www.philips.com/support** and to the competent authority of the country in which the user and/or patient is established.

When using the breast pump, the below undesirable side effects may occur. If you experience any of these symptoms, contact a healthcare professional or breastfeeding specialist.

Undesirable side effect	Description
Pain sensation	Pain felt on the breast or nipple, or felt as a response to the suction being above a comfortable level.
Sore nipples	Persistent pain in the nipples at the beginning of the pumping session, or that lasts throughout the entire pumping session, or pain between sessions.
Erythema (redness)	Redness of the skin of the breast.

Undesirable side effect	Description
Bruise	A reddish-purple discoloration that does not fade when pressed. When a bruise fades, it may become green and brown.
Thrombus	Thrombus can result into tenderness, pain, redness, swelling and/or heat.
Injured tissue on nipple (nipple trauma)	- Fissures or cracked nipples. - Skin tissue peeling off the nipple. Normally occurs in combination with cracked nipples and/or blisters. - Nipple tear. - Blister. Looks like small bubbles on the surface of the skin. - Bleeding. Cracked or torn nipples can lead to bleeding of the affected area.

Storage conditions

Keep the breast pump out of direct sunlight as prolonged exposure may cause discoloration. Store the breast pump and its accessories in a safe, clean and dry place.

If the device has been stored in a hot or cold environment, place it in an environment with a temperature of 20 °C for 30 minutes to let it reach a temperature within the usage conditions (5 °C to 40 °C) before you use it.

Fully charge the battery before storing your breast pump for an extended period of time.

Electromagnetic fields (EMF)

This Philips device complies with all applicable standards and regulations regarding exposure to electromagnetic fields.

Technical information

	Single electric breast pump	Double electric breast pump
Vacuum level	Stimulation: -60 to -200 mbar (-45 to -150 mmHg) Expression: -60 to -360 mbar (-45 to -270 mmHg) Adjustable in steps of 20 mbar Intermittent vacuum	
Cycle speed	Stimulation: 105-120 cycles/min, depending on vacuum level Expression: 53-85 cycles/min, depending on vacuum level	
Motor unit rated input	5 V direct current / 1.1 A	5 V direct current / 1.8 A
Mains adapter rated input	100-240 V  50/60 Hz 400 mA	
Mains adapter rated output	5 V direct current / 1.1 A	5 V direct current / 1.8 A

	Single electric breast pump	Double electric breast pump
Adapter type number	S009AHZ050yyyy The letters "yyyy" represent the output current from 0100 (1000 mA) to 0180 (1800 mA), increments in steps of 100 mA. The letter "z" denotes type of plug, it can be: U, V, B, S, C, A, D, K, T, I, E, F.	
Battery	Rechargeable battery: 3.6 V 2600 mAh (non-replaceable)	
Power port	Micro-USB	
Operating conditions	Temperature: 5 °C to 40 °C Relative humidity: 15 to 90 % (non-condensing)	
Storage and transportation conditions	Temperature: -25 °C to 70 °C Relative humidity: up to 90 % (non-condensing) For optimal performance of the built-in rechargeable battery, it is recommended to store the device at room temperature.	
Operating and storage pressure	700-1060 hPa (<3000 m altitude)	
Net weight	Approximately 275 g	Approximately 355 g
External dimensions	Motor unit: 145 mm x 95 mm x 45 mm (LxWxH)	
Device classification	Applicable classifications: Degree of protection against electric shock: - Class II ME equipment - Type BF applied part The breast pump and the adapter are specified as medical electrical system according to IEC 60601-1.	
Materials	ABS, silicone (motor unit) Polypropylene, silicone (other parts)	
Service life	500 hours	
Mode of operation	Continuous operation	Non-continuous operation 55 min on, 180 min off at 22 °C (duty cycle 23 %) 20 min on, 90 min off at 40 °C (duty cycle 18 %)
Automatic power-off	After 30 minutes (after last user interaction)	
Ingress protection	IP22 (motor unit), IP20 (adapter)	

Electromagnetic compatibility - compliance information

This device has been tested and found to comply with IEC 60601-1-2:2014 Ed 4.0 for electromagnetic compatibility (EMC) as Class B according to CISPR 11:2009. Tests levels are listed in the accompanying tables. Only use this Philips Avent single/double electric breast pump Premium in combination with the adapter supplied.

Emissions	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The device does not use RF energy for its intended function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in locations in residential environments and in establishments directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Complies Class A Category	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The Philips Avent single/double electric breast pump Premium is intended for use in the electromagnetic environment listed below. The user of the device should ensure that it is used in such an environment to prevent improper operation.

Immunity test	Immunity test level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	Floors can be a source of static electricity. When the air is very dry (<20 % RH) the user is advised to prefer to use the device in a room with concrete or wooden floor (instead of using it in a room with a synthetic floor).
Radiated RF EM fields IEC 61000-4-3	10 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	385 (27), 450 (28), 710 (9), 745 (9), 780 (9), 810 (28), 870 (28), 930 (28), 1720 (28), 1845 (28), 1970 (28), 2450 (28), 5240 (9), 5500 (9), 5785 (9) MHz (V/m)	The device is immune against any RF wireless communication device such as smart phone, Wi-Fi or Bluetooth device at a distance of 30 cm or more. Avoid a closer distance as this may result in improper operation.
Rated power frequency magnetic fields IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	
Electrical fast transient/burst IEC 61000-4-4	±2 kV 100 kHz repetition frequency	Mains power quality expected of a typical commercial or residential environment.
Surges IEC 61000-4-5	±1 kV Line-to-line ±2 kV Line-to-ground	Mains power quality expected of a typical commercial or residential environment.

Immunity test	Immunity test level	Electromagnetic environment - guidance
Conducted disturbances induced by RF fields IEC 61000-4-6	3 Vrms, before modulation is applied 0.15 MHz – 80 MHz 6 Vrms, before modulation is applied, in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles single phase at 0° 0 % UT; 250/300 cycle	Mains power quality expected of a typical commercial or residential environment. At power interruptions the breast pump will automatically switch-over to battery mode without interruption of the function.

Note: UT is the AC mains voltage prior to application of the test level. If in the unlikely event the Philips Avent single/double electric breast pump Premium is locked inoperative in a function by disturbances from other devices beyond the levels stated in the table above, the device will shut-off completely and restart. In that case make sure the distance to the other device causing this disturbance is enlarged.

Explanation of symbols

The warning signs and symbols are essential to ensure that you use this device safely and correctly and to protect you and others from injury. Below you find the meaning of the warning signs and symbols on the label and in the user manual.

Symbol	Explanation
	Indicates to follow the instructions for use.
	Indicates important information such as warnings and cautions.
	Indicates usage tips, additional information or a note.
	Indicates the manufacturer.
	Indicates the manufacturing date.
	Indicates the compliance with the essential requirements of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

Symbol	Explanation
	Indicates manufacturer's catalog number of the breast pump.
	Indicates the batch number of the medical device.
	Indicates that the device is a medical device.
	Indicates that the part of the appliance that comes into physical contact with the user (also known as the applied part) is of type BF (Body Floating) according to IEC 60601-1. The applied part is the breast pump kit.
IP22	Indicates the protection against ingress of solid foreign objects greater than 12.5 mm and against harmful effects due to the ingress of dripping water when tilted at 15 degrees.
IP20	Indicates the protection against ingress of solid foreign objects greater than 12.5 mm and no protection against ingress of water.
	Indicates 'alternating current'.
	Indicates 'direct current'.
	Indicates 'Class II Equipment'. The adapter is double-insulated (Class II).
	Indicates USB.
	Indicates the storage and transportation temperature limits to which the device can be safely exposed: -25 °C to 70 °C.
	Indicates the relative humidity upper limits to which the device can be safely exposed: up to 90 %.
	Indicates the operating pressure range: 700 to 1060 hPa.
	Indicates that the breast pump must be kept dry.
	Indicates that the breast pump must be kept out of direct sunlight.
	Indicates 2-year Philips worldwide guarantee.

Symbol	Explanation
	Indicates Forest Stewardship Council. The FSC trademarks enable consumers to choose products that support forest conservation, offer social benefits, and enable the market to provide an incentive for better forest management.
	Indicates that the material used in the product is safe for food contact.
	Identifies the UDI carrier, including the AIDC and human readable information.
	Indicates that the material used for the packaging is corrugated board.
	Indicates that the material used for the packaging is non-corrugated cardboard/fiberboard.
	Indicates that the material used for the plastic bag is low density polyethylene.
	Indicates indoor use only.
	Indicates that the part is certified with six level energy efficiency rating.
	Indicates the Moroccan Conformity Mark.

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Trademarks are the property of Koninklijke Philips N.V. / Les marques commerciales sont la propriété de Koninklijke Philips N.V.
All rights reserved. / Tous droits réservés.
Philips Consumer Lifestyle B.V.