
Register your product and get support at
www.philips.com/welcome

HF3330

PHILIPS

ENGLISH 6

DANSK 25

SUOMI 44

NORSK 62

SVENSKA 81

Introduction

Congratulations on your purchase and welcome to Philips! To fully benefit from the support that Philips offers, register your product at www.philips.com/welcome.

Intended use

The goLITE BLU energy light is intended to make people feel more energetic, to adjust circadian rhythms and to provide mood relief from the winter blues. The medical intended use of the goLITE BLUE energy light is to treat winter depression (seasonal affective disorder). This treatment should always take place under supervision of a specialist.

Light, health and circadian rhythms

The human body uses light cues, such as those provided by the sun, to time certain internal functions. Properly timed rhythms regulate mood, sleep, energy, appetite and digestion. These daily internal cycles, called circadian rhythms, sometimes fall out of sync, resulting in a dysregulated body clock.

Unfortunately, modern living has dramatically altered nature's cues. A modern day no longer starts at the crack of dawn and ends at 8:00 p.m. Workdays are getting longer and many people face shift-work schedules. Additionally, the advent of electric lighting allows social gatherings and personal activities to extend well into the night. These factors have diminished the body's natural ability to regulate the body clock.

Rapid travelling across time zones is another source of disharmony between the internal clock and the external time, which is referred to as jet lag. Proper timing of bright light exposure allows quick adjustment to the new time zone.

Besides regulating the body clock, light has direct alerting and energising effects and improves our performance. Bright light therapy is also successfully used for mood relief, especially for seasonal mood swings. In countries with long winters, many people experience loss of energy and low spirits due to lack of light. These symptoms are referred to as the winter blues. They typically start in autumn and may last for several weeks. Exposure to bright light sources can help alleviate these symptoms and restore the summer mood and summer energy level.

Research shows that a special photo receptor in our eyes is responsible for regulating our energy and sleep/wake cycles. This receptor responds to the blue light of the summer sky. We do not get as much of this light in the autumn and winter, and indoor light does not emit enough of this essential colour. The goLITE BLU energy light helps to shift circadian rhythms back to their normal pattern by delivering the wavelength, colour and intensity of light that our body responds most efficiently to. The goLITE BLU energy light naturally improves your mood and energy levels by a relatively low light intensity compared to traditional white light therapy devices.

General description (Fig. 1)

- A LED panel
 - B Socket for appliance plug
 - C Interactive touchscreen display
 - 1 Alarm touchscreen button
 - 2 Alarm chime indication
 - 3 Alarm light indication
 - 4 Timer-on indication
 - 5 Clock time
 - 6 Demo mode indication
 - 7 Charge indication
 - 8 Plus (+) touchscreen button
 - 9 Timer touchscreen button
 - 10 Timer countdown time/Alarm time
 - 11 Intensity indication
 - 12 Minus (-) touchscreen button
 - D Small plug
 - E Adapter
 - F Stand
 - G Reset button
 - H Stand storage cavity
 - I Cover of built-in rechargeable battery
 - J On/off button
- Not shown: Pouch

Important

- Read this user manual carefully before you use the appliance and save it for future reference.
- The adapter is double insulated according to Class II (Fig. 2).

Danger

- Keep the adapter dry.

Warning

- Never use the adapter if it is damaged in any way. Always have the adapter replaced with one of the original type in order to avoid a hazard.
- The adapter contains a transformer. Do not cut off the adapter to replace it with another plug, as this causes a hazardous situation.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance.
- Keep the appliance away from water. Do not place or store it over or near water contained in a bathtub, washbasin, sink etc. Do not immerse it in water or any other liquid. This might impair the safety of the appliance.
- Use the goLITE BLU energy light in vertical position with the stand in place. Do not use it in horizontal position without the stand, as this might cause the appliance to overheat.
- Do not use the appliance in a room that is not illuminated by any other lamp, as this may cause eyestrain. Use the light of this appliance as an addition to the normal light (sunlight, electric light) in your home.
- Do not look straight into the light too long, as this may be uncomfortable for your eyes.
- Do not use the goLITE BLU energy light where its light may compromise your or other people's ability to perform essential tasks such as driving or operating heavy equipment.

- Do not sit closer to the appliance than recommended to obtain a higher light intensity. Increased light intensity does not improve the result and may increase the risk of headache or eyestrain.
- Do not use the appliance adjacent to or on top of other equipment, nor put any equipment on top of the appliance. However, if this is unavoidable, observe the appliance to make sure that it operates normally.
- When the goLITE BLU energy light is connected to the mains, do not position it in a place where it is difficult to remove the adapter from the wall socket.
- No modification of the goLITE BLU energy light is allowed.
- The appliance is not suitable for use in the presence of flammable anaesthetic mixtures with air, oxygen or nitrous oxide.
- Always consult your doctor before you start using the goLITE BLU energy light if:
 - 1 you suffer or have suffered from depression;
 - 2 you have bipolar disorder;
 - 3 you are hypersensitive to light (e.g. due to diabetes or epilepsy);
 - 4 you have an eye complaint or have had eye surgery;
 - 5 you are on photo-sensitising medications (e.g. certain antidepressants, psychotropic drugs or malaria tablets).

Caution

- Use, charge and store the appliance at a temperature between 0°C and 35°C.
- Only use the adapter supplied (AC/DC adapter KSAD0950120W1UV-1) to charge the appliance. If you use a different adapter than the adapter specified and provided by Philips, this may result in increased electromagnetic emission or decreased immunity of the goLITE BLU energy light.
- The appliance has an integrated rechargeable lithium-ion battery. Do not incinerate, disassemble or expose the rechargeable battery to temperatures above 60°C. The rechargeable battery may explode if overheated.
- Do not use the appliance outdoors or near heated surfaces or heat sources (Fig. 3).

- Stop using the appliance if it is damaged in any way. If repair is needed, consult chapter 'Guarantee and service' for information.
- Do not place heavy objects on the appliance, the adapter or the cord.
- Do not pull at the cord when you remove the adapter from the wall socket. Always grasp the adapter.
- Use the appliance according to the guidelines for electromagnetic compatibility (EMC) provided in this user manual.

Safety and compliance

- This Philips appliance complies with all standards regarding electromagnetic fields (EMF). If handled properly and according to the instructions in this user manual, the appliance is safe to use based on scientific evidence available today.
- This appliance meets the requirements of Medical Device Directive MDD 93/42/EEC and 2007/47/EC for medical appliances.


General

The adapter is equipped with an automatic voltage selector and is suitable for mains voltages ranging from 100 to 240 volts.

Symbol key

The following symbols may appear on the appliance:

Symbol	Description
	Consult the user manual supplied
	DC Power
	DC Jack Polarity
	Serial Number
HF3330	Model Number

Symbol	Description
 0344	European Declaration of Conformity
	Compliant with the Waste Electrical and Electronic Equipment/Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment (WEEE) recycling directives.
	Philips Shield
	Manufactured for: Philips Consumer Lifestyle B.V, Tussendiepen 4, 9206 AD Drachten, The Netherlands. Fax: +31 (0)512594316
 YYYYYWW 	Date Code

Preparing for use

When to use the goLITE BLU energy light

Knowing what time to use the goLITE BLU energy light can mean the difference between feeling better in just a few days instead of weeks, or perhaps not responding at all. Follow these guidelines to understand how light affects your body clock and when to use the goLITE BLU energy light for the best results.

Duration and intensity

Most people get beneficial results from 15 to 30 minutes of usage each day. You can start using the goLITE BLU energy light at 50% intensity and increase the light intensity if you feel no response after a week.

Sleep and energy

Using the goLITE BLU energy light at any time of day directly affects your alertness. You can use the goLITE BLU energy light as light therapy for the winter blues throughout the day. The best time to use it is in the morning. Do not use the goLITE BLU energy light late in the evening. If you use it within 2 hours of going to bed, it could disrupt your sleep.

- If you find it difficult to get up each day without an alarm clock or if you often sleep in, you need morning light. Use the goLITE BLU energy light within the first hour or two after waking up.
- If you find it difficult to stay awake in the late evening and wake up too early in the morning, you need evening light. Use the goLITE BLU energy light before you start to feel drowsy.
- If mornings or evenings do not seem to be the problem, but you feel generally tired during the day, then the time at which you use the goLITE BLU energy light is not as important. You can use the goLITE BLU energy light in the morning or around midday. For example, when you wake up in the morning, use it for fifteen to thirty minutes and then again at midday or in the evening for approximately the same amount of time. Place the goLITE BLU energy light further away from you during the evening treatment. This keeps your body clock from confusing the stronger morning signal with evening time. Do not use the goLITE BLU energy light within two hours of your desired bedtime. Blue light can suppress the nighttime hormones for a while, and you may not be able to fall asleep.

Winter depression

If your seasonal mood swings are more severe and you are suffering from depression, you should consult a therapist and only use light therapy under supervision. Your therapist may prescribe 45 minutes of treatment per day at level 1 or 2 over several weeks. The timing of light exposure can be determined as described above, but should be aligned with your therapist.

Jet lag

Rapidly travelling across several time zones creates disharmony between your internal body clock and the new external time. This not only causes sleeping problems, but influences your mood and energy level as well. If

you stay for more than a few days, you may also want to shift your body clock to the new time zone. This will help you sleep better and be more alert and energetic when you need to. Because special blue light is effective for shifting circadian rhythms, your goLITE BLU energy light can quickly adjust your rhythm to the new time zone. One of the best strategies to combat jet lag is to shift your body clock a day or two before you leave.

Travelling east: When you travel east, start using the goLITE BLU energy light in the morning two to three days before you leave. The number of time zones you travel corresponds to the hours your body clock has to shift. If you are only travelling a few time zones, one day's preparation is sufficient. Add a day of light therapy for each additional 2 time zones. For example, for a trip from the US to Europe, use the following schedule: On the first day, wake up an hour early and use the goLITE BLU energy light for 30–45 minutes. On the next day or days, you should wake up a couple of hours earlier and repeat this process. You should advance your bedtime accordingly. Once you arrive at your destination, avoid light before 10 a.m. local time (wear sunglasses) and use the goLITE BLU energy light in the midmorning to fully adjust to the local time. When returning, use your goLITE BLU energy light in the late evening for a couple of nights before the return flight. On the last night of the trip, try to stay up and use your goLITE BLU energy light as late as reasonably possible past midnight. On the day of the flight, avoid any bright light until after 2:00 p.m. (see below under 'Travelling west')

Travelling west: Since most people's body clocks adjust more easily to later time zones, you should not have to use the goLITE BLU energy light when you travel 3 time zones or less. You should only need to use your goLITE BLU energy light for one day before travelling 4 time zones, and add an additional preparation day for each 3 time zones. You should not have to use the goLITE BLU energy light for more than 3 days, even on the longer westbound trips. For example, when travelling from London to Los Angeles, use your goLITE BLU energy light for 30–40 minutes in the late evening (and again a couple of hours later, if you can) each night, starting 2–3 days before you leave. Use the goLITE BLU energy light 2–3 hours later each successive night. When you arrive at your destination, try to stay awake until bedtime and avoid bright evening light (wear sunglasses) for the first few days. Before you travel home, avoid any bright

afternoon and evening light and use your goLITE BLU energy light in the morning for a few days (see above under 'Travelling east').

For more instructions on how to avoid jet lag, visit www.philips.com/golite.

Charging

Note: You can also use the appliance while it is charging.

The goLITE BLU energy light comes with a 9.5V AC/DC adapter that operates on voltages ranging from 100 to 240 volts to allow worldwide use. The adapter comes with a set of exchangeable adapter plugs that fit most common worldwide wall sockets.

- 1** Put the adapter in the wall socket.
- 2** Insert the appliance plug into socket on the appliance (Fig. 4).
 - Only the outline of the charge indication is visible.
 - As the battery charges, the four bars of the charge indication flash one by one and fill as the charge increases.
 - When the battery is fully charged, the charge indication stops flashing.

Note: Charge the goLITE BLU energy light at least 8 hours before you use it for the first time. Charging normally takes 3 to 4 hours.

Note: When the battery is fully charged, it contains enough energy for 3 to 5 treatments.

Installation

- 1** Remove the stand from the stand storage cavity (Fig. 5).
- 2** Insert the bevelled end of the stand into the hole on the back of the appliance. Magnets hold the stand in place (Fig. 6).
- 3** Place the goLITE BLU energy light approximately 50-75cm from your eyes in such a way that the light bathes your face from the side (Fig. 7).
 - You can read, eat, work at the computer, watch TV or exercise while you use your goLITE BLU energy light.

- For the best results, place the goLITE BLU energy light at the same level as your midriff (stomach area), e.g. place it on the table at which you are sitting.

Note: *Do not stare into the light of the LED panel. It is sufficient that the blue light reaches your eyes indirectly from the side.*

- Use your goLITE BLU energy light in a well-lit room to minimise eye strain.

Setting the clock

- 1** Press the on/off button briefly to switch on the appliance (Fig. 8).
 - ▶ The touchscreen display goes on.
- 2** Press the alarm button and the timer button simultaneously for 2 seconds (Fig. 9).
 - ▶ The clock time starts to flash on the display (Fig. 10).
- 3** Set the clock time with the plus or minus button. If you keep the button pressed, the time moves forward or backward quickly (Fig. 11).
- 4** To switch between the 12-hour (a.m./p.m.) and 24-hour time indication, press the timer button and minus button simultaneously for 2 seconds.
- 5** Press either the alarm button or the timer button or wait 10 seconds to exit the time-setting mode.

Setting the timer

You can use the timer function to set the treatment time.

- 1** Press the on/off button briefly to switch on the appliance (Fig. 8).
 - ▶ The touchscreen display goes on.
- 2** Press the timer button (Fig. 12).
 - ▶ The timer on-indication and the timer countdown time start to flash (Fig. 13).

- 3** Press the plus or minus button to set the treatment time. If you keep the button pressed, the time moves forward or backward quickly (Fig. 11).

You can set a time from 1 to 60 minutes.

- ▶ The set session time is stored in the memory.
- ▶ The timer-on indication and the timer countdown time stop flashing.

Note: When you switch on the appliance, the timer counts down the set time. When the countdown is finished, the LED panel switches off automatically.

Demo mode

The demo mode is used to demonstrate the goLITE BLU energy light in shops. When the appliance is switched on in demo mode, customers of the shop can:

- adjust the light intensity
- switch the appliance on and off
- set the alarm
- set the treatment time

If the appliance is not operated for 3 minutes, it reverts to its demo mode default settings.

To activate the demo mode:

- 1** Make sure the appliance is switched off.

Note: If you want to switch off the appliance, press and hold the on/off button for 2 seconds.

- 2** Put the adapter in the wall socket and insert the appliance plug into socket on the appliance.

- 3** Press and hold the on/off button for 10 seconds (Fig. 8).

- ▶ The appliance switches on and the word 'DEMO' appears on the display.
- ▶ The LED panel goes on at a light intensity of 25%.
- ▶ The display shows a timer countdown time of 15 minutes, but the timer does not count down. It only starts to count down when the timer time is adjusted.

- ▶ The display shows an alarm time of 12:00 a.m. and shows the normal clock time, but the alarm mode is off.

Note: The appliance remains switched on in demo mode until it is switched off.

To deactivate the demo mode, press and hold the on/off button for 10 seconds.

Using the appliance

- 1 Press the on/off button briefly to switch on the appliance (Fig. 8).
 - ▶ The touchscreen display goes on.

Note: If you do not touch the display or press any button for 5 minutes when the display is switched on and the appliance operates on battery power, the appliance switches off automatically to save battery power.

Note: If you want to switch off the appliance, press and hold the on/off button for 2 seconds.

- 2 Briefly press the on/off button again to switch on the LED panel.
 - ▶ The timer counts down the set time.

Note: The default treatment time is 15 minutes. Use the timer to change the treatment time (see section 'Setting the timer' in chapter 'Preparing for use').

- ▶ When the countdown is finished, the LED panel switches off automatically.

Note: If you want to switch off the LED panel before the countdown is finished, press the on/off button briefly.

- 3 Press the plus or the minus button to change the light intensity (Fig. 11).

You can increase or decrease the intensity in increments of 25%.

Pausing a session

- 1 Press the on/off button or press the minus button to reduce the intensity to 0%.
 - ▶ The remaining minutes of the session are still shown on the display.

- 2 Press the on/off button to continue the session or press any of the touchscreen buttons to end the session.

Note: If you pause a session while the appliance operates on battery power and do not continue the session within 5 minutes, the appliance switches off automatically.

Setting the alarm

The alarm has three modes: chime only, light only or both chime and light.

- 1 Press the on/off button briefly to switch on the appliance (Fig. 8).
 - ▶ The touchscreen display goes on.
- 2 Press the alarm button (Fig. 14).
 - ▶ The alarm indications start to flash (Fig. 15).
- 3 Press the plus or minus button to switch between the three alarm modes: chime only, light only or both chime and light (Fig. 11).
- 4 Press the alarm button twice or wait 10 seconds to store the selected mode.
- 5 To set the alarm time, press the alarm button twice.
 - ▶ The alarm time starts to flash.
- 6 Press either the plus or the minus button to set the desired time. If you keep the button pressed, the time moves forward or backward quickly.
- 7 Press the alarm button or wait 10 seconds to store the set alarm time.

Note: The set alarm mode and time remain stored in the memory as long as there is energy in the battery.

Switching off the alarm sound

- 1 Press any of the touchscreen buttons or the on/off button to switch off the alarm sound when the alarm goes off.

Disabling the alarm

- 1 To disable the alarm, press and hold the alarm button until the alarm chime and the alarm light indications as well as the alarm time disappear from the display.

Cleaning

Never immerse the appliance in water or any other liquid, nor rinse it under the tap.

Never use scouring pads, abrasive cleaning agents or aggressive liquids such as petrol or acetone to clean the appliance.

- 1 Clean the appliance with a dry cloth.

Environment

- Do not throw away the appliance with the normal household waste at the end of its life, but hand it in at an official collection point for recycling. By doing this, you help to preserve the environment (Fig. 16).
- The built-in rechargeable battery contains substances that may pollute the environment. Always remove the rechargeable battery before you discard and hand in the appliance at an official collection point. Dispose of the rechargeable battery at an official collection point for batteries. If you have trouble removing the rechargeable battery, you can also take the appliance to a Philips service centre. The staff of this centre will remove the rechargeable battery for you and will dispose of it in an environmentally safe way. (Fig. 17)

Removing the rechargeable battery

Only remove the rechargeable battery when you discard the appliance. Make sure the battery is completely empty when you remove it.

- 1 Disconnect the appliance from the mains and remove the appliance plug.

- 2** Switch on the LED panel and leave it on until it goes out automatically.
- 3** Undo the 2 screws with a screwdriver (Fig. 18).
- 4** Remove the cover of the built-in rechargeable battery. (Fig. 19)
- 5** Lift the blue rectangular rechargeable battery out of the appliance. The battery is connected to the appliance through a red wire and a black wire to a socket inside the appliance (Fig. 20).
- 6** Grab the wires and pull the socket out of the appliance. (Fig. 21)
D The rechargeable battery has now been removed.

Guarantee and service

If you need service or information or if you have a problem, please visit the Philips website at www.philips.com or contact the Philips Consumer Care Centre in your country (you find its phone number in the worldwide guarantee leaflet). If there is no Consumer Care Centre in your country, go to your local Philips dealer.

Frequently asked questions

This chapter lists the questions most frequently asked about the appliance. If you cannot find the answer to your question, contact the Consumer Care Centre in your country.

Question	Answer
Do I have to take off my glasses when I use the goLITE BLU energy light?	No, you can use the goLITE BLU energy light with your glasses on. Do not wear tinted glasses or sunglasses when you use the goLITE BLU energy light, as these glasses reduce its effectiveness.

Question	Answer
Does the appliance have any side effects?	The goLITE BLU energy light is safe if used according to the instructions. Temporary headaches or nausea have been reported. These effects are usually solved by sitting further away from the appliance.
Can I get sunburn from using the goLITE BLU energy light?	No, the technology used in the goLITE BLU energy light eliminates all UV light.
I have heard that light can be dangerous. Can the goLITE BLU energy light damage my eyes?	The goLITE BLU energy light emits safe levels of blue light, according to current standards. However, if you suffer from an eye disease or have a history of eye disease in your family, consult your ophthalmologist before you use the goLITE BLU energy light. If you have a temporary eye problem, wait until it is resolved before you use the goLITE BLU energy light.
	The goLITE BLU energy light is safe if used according to the instructions. The appliance is free from any ultraviolet light risk to the eye and it produces a safe, diffused field of light. The appliance has been tested for ocular safety.
	As with any bright light source, do not look into the light too long. The beneficial reaction happens when the photoreceptors in the eye are stimulated. Since these photoreceptors lie in the periphery of the eye, you get a better result if you place the goLITE BLU energy light at a slight angle to your face.

Question	Answer
Should I use the goLITE BLU energy light every day?	You can use the goLITE BLU energy light daily, especially during the winter months. Many people find that they do not need the light in summer. Others use it year-round because they spend most of their time indoors and do not get the right kind of light signals to regulate their circadian rhythm properly.
Can the light cause harm to children or pets?	No. The light is safe to be used around children and pets. However, do not let children play with the appliance.
Can I use the goLITE BLU energy light more than once a day?	For most people, it is sufficient to use the goLITE BLU energy light once a day. Keep in mind that you may have trouble falling asleep if you use the goLITE BLU energy light too late in the evening.
How long do the LEDs last?	If you use the appliance 30 minutes a day throughout the year, the LEDs have a lifetime of over 50 years.
What if one or more LEDs stop working prematurely?	In the unlikely event that one or two LEDs fail prematurely, you can continue to use your goLITE BLU energy light normally because its effectiveness is not affected. However, if several LEDs stop functioning, have the goLITE BLU energy light repaired to ensure that you get the right amount of light. See chapter 'Guarantee and service'.

Question	Answer
What does 'DEMO' mean on the display?	When 'DEMO' is on the display, the appliance is in demo mode. This mode is used to demonstrate the goLITE BLU energy light in shops. To switch off this mode, press and hold the on/off button for 10 seconds when the goLITE BLU energy light is switched off. The goLITE BLU energy light switches on and the word 'DEMO' disappears from the display.
How do I use the reset button?	The reset button is located in the stand storage cavity on the back of the appliance. If the goLITE BLU energy light locks up or does not work correctly, reset the appliance by inserting the end of a paperclip into the reset button hole. This resets the goLITE BLU energy light and clears its memory.

Specifications

Model	HF3330
-------	--------

Electrical

Rated voltage appliance (VAC)	100-240
Rated voltage adapter (V DC)	9.5
Rated frequency (Hz)	50/60
Rated input power adapter (A)	1.2
Power consumption appliance (A)	1.1
Power consumption adapter (A)	0.4

IEC 60601-1 classification

- Insulation class	II
- Ingress of water	IPX0

Model	HF3330
--------------	---------------

- Mode of operation	Continuous
---------------------	------------

Rechargeable battery type

2-cell lithium-ion pack	7.4V, 850mAh
-------------------------	--------------

IEC 62133 and UL 2054 compliant

Operating conditions

Temperature	from 0°C to +35°C
-------------	-------------------

Relative humidity	from 15% to 90% (no condensation)
-------------------	-----------------------------------

Storage conditions

Temperature	from -20°C to +50°C
-------------	---------------------

Relative humidity	from 15% to 90% (no condensation)
-------------------	-----------------------------------

Physical characteristics

Dimensions	14 x 14 x 2.5cm
------------	-----------------

Weight	0.4kg (less than 14 ounces)
--------	-----------------------------

Peak wavelength	468nm +/- 8nm
-----------------	---------------

Half-peak bandwidth	20nm
---------------------	------

Typical irradiance (100% intensity)	300μW/cm ² +/- 10% at 50cm
-------------------------------------	---------------------------------------

Circuit diagrams, component part lists and other technical descriptions are available on request.

As user of this appliance, you should take measures to prevent the occurrence of electromagnetic interference or other types of interference between this appliance and other appliances.

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