

Register your product and get support at
www.philips.com/welcome

HF3319

PHILIPS

1 — O/I

2 — ☀

3 — ●●●●●

ENGLISH	6
DANSK	12
DEUTSCH	18
SUOMI	25
FRANÇAIS	31
ITALIANO	38
NEDERLANDS	44
NORSK	50
SVENSKA	56

Introduction

Congratulations on your purchase and welcome to Philips! To fully benefit from the support that Philips offers, register your product at www.philips.com/welcome.

Intended use

This appliance helps to counteract the tiredness, low spirits and lack of energy that many people experience in autumn and winter, when the days are darker and shorter. These symptoms are often referred to as the winter blues. A shortage of light can also have a negative effect on our sleep-wake rhythm and our biological clock.

The medical intended use of the EnergyLight is to treat winter depression (seasonal affective disorder). This treatment should always take place under supervision of a specialist.

The EnergyLight is suitable for both light baths and normal lighting.

EnergyLight and the winter blues

The winter blues are characterised by a combination of the following symptoms: tiredness, lack of energy, a need for more sleep, unexpected mood changes, trouble to concentrate and trouble to get out of bed in the morning. Sometimes the winter blues can also cause you to gain weight, as your body craves food that is rich in sugar and carbohydrates, such as pasta and chocolate. The winter blues occur during the darker period of the year, i.e. from September until March. Many people are troubled by one or more symptoms of the winter blues. Many health experts agree that the winter blues can be counteracted by exposure to bright light. 'Bright light' means daylight with a certain intensity. In winter, the light intensity is often much lower than in summer: on a gloomy winter day the light intensity varies between 2,500 and 10,000 lux, whereas this may be as much as 100,000 lux on a bright summer day. Moreover, in winter most people spend a lot of time indoors, either at home or at work, where artificial light often has an intensity of as little as 500 lux and sometimes even less. Together with scientists and medical specialists in the field of winter blues treatment, Philips has developed this special lighting appliance. This appliance allows you to safely get the amount of light you need to make the winter blues disappear. Taking light baths is generally accepted as one of the safest and most effective ways of fighting the winter blues.

EnergyLight and the sleep-wake rhythm

A shortage of light plays a major role in the disruption of our sleep-wake rhythm. Hormones are no longer produced at exactly the right time of day, it becomes harder to sleep soundly throughout the night and the need for short naps during the day increases. An extra dose of light restores our sleep-wake rhythm and improves the quality of our sleep as well as the way we function during the day. The EnergyLight provides this extra dose of light.

For more information about the health benefits, see www.philips.com, search word "EnergyLight".

General description (Fig. 1)

- A Energy lamp
- B Lamp housing with angle adjustment
- C Base
- D On/off knob and treatment time indicator
 - ① On/off symbol
 - ② Light symbol
 - ③ Light intensity settings from low to high
- E Treatment time symbol

Important

Read this user manual carefully before you use the appliance and save it for future reference.

Danger

- Water and electricity are a dangerous combination. Do not use this appliance in wet surroundings (e.g. in the bathroom or near a shower or swimming pool).

Warning

- Check if the voltage indicated on the appliance corresponds to the local mains voltage before you connect the appliance.
- If the mains cord is damaged, you must have it replaced by Philips, a service centre authorised by Philips or similarly qualified persons in order to avoid a hazard.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance.
- No modification of this appliance is allowed.

Caution

- Always unplug the appliance after use.
- Do not use the appliance in a room that is not illuminated by any other lamp, as this may cause eyestrain. Use the light of this appliance as an addition to the normal light (sunlight, electric light) in your home.
- Do not look straight into the light too long, as this may be uncomfortable for your eyes.
- Consult your doctor before you start using the Philips EnergyLight if:
 - 1 you suffer or have suffered from severe depression;
 - 2 you are hypersensitive to light (e.g. due to diabetes or epilepsy);
 - 3 you have an eye complaint;
 - 4 you are on specific medication (i.e. certain antidepressants, psychotropic drugs or malaria tablets);
 - 5 you have high blood pressure.
- If use of the appliance causes you to become too active, we advise you to reduce the exposure to bright light. If there is no improvement, consult your doctor.
- If you feel exhausted, distressed or restless and/or suffer from sleeping problems (e.g. insomnia) after you have used the appliance, we advise you to consult your doctor or psychotherapist.
- This appliance does not have an on/off switch. If you want to make sure that the appliance is fully isolated from the mains supply, remove the plug of the appliance from the wall socket.


General

- The duration of the light bath depends on the distance at which the appliance is used. See chapter 'Using the appliance' for the appropriate distance and duration.
- After the first few light baths, you may experience some eye strain and/or a light headache. These side effects tend to disappear quickly.
- This appliance meets the requirements of Medical Device Directive MDD 93/42/EEC and 2007/47/EC for medical appliances.

Symbol key

The following symbols may appear on the appliance:

Symbol	Description
	Consult the user manual supplied
HF3319	Model number
 0344	European Declaration of Conformity
O/I	Off/On

Symbol	Description
	Class II equipment
	Compliant with the Waste Electrical and Electronic Equipment/Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment (WEEE) recycling directives.
	General warning symbol
	Philips Shield
	Manufactured for: Philips Consumer Lifestyle B.V, Tussendiepen 4, 9206 AD Drachten, The Netherlands. Fax: +31 (0)512594316
	Date Code

Using the appliance

- It is advisable to start using the EnergyLight as soon as you feel the winter blues coming on or when you feel like you need extra light.
- Preferably take your light baths between 6 o'clock in the morning and 8 o'clock in the evening. (Fig. 2)

Tip: Visit www.philips.com/EnergyLight for more information on when to use the EnergyLight in your case.

- 1 Place the appliance on a tabletop.
- 2 Put the plug in the wall socket.
- 3 To switch on the appliance, turn the on/off knob clockwise to the light intensity you prefer.

Duration

The advisable duration of the light bath depends on the distance at which the appliance is used:

- approx. 2 hours at a distance of 50cm
- approx. 1 hour at a distance of 30cm
- approx. 30 minutes at a distance of 15-20cm

Treatment time indicator

The treatment time indicator shows how long you have been using the appliance.

- When you switch on the appliance, the LEDs go on in quick succession to show how the LED sequence is built up. Then the first LED lights up at half intensity during the first 15 minutes (Fig. 3).
- At the end of the first 15 minutes, the first LED is on at full intensity and the second LED goes on at half intensity.
- After 2 hours, all 8 LEDs light up at full intensity until you switch off the appliance.

Note: The LEDs only show the treatment time that has passed. You do not have to take any action when the LEDs go on.

- 4 Position yourself near the appliance and turn the screen in such a way that you can comfortably look in the direction of the light (Fig. 4).

Make sure there is a distance between your face and the screen of 15 to 50cm.

- You do not need to sit right in front of the appliance
- You do not have to take the whole light bath without interruption. If it suits your schedule better, you can interrupt your light bath and continue later.

- 5 Look into the light every now and then.

Do not look into the light continuously. You can simply engage in other activities such as reading, writing, embroidering or handicraft while taking a light bath and look into the light every now and then.

6 Repeat the light baths until you feel better.

Generally, improvements occur within 1 week.

You may repeat the light bath as often as you like during the dark period of the year. Take a light bath on at least 5 successive days.

7 After use, turn the on/off knob anticlockwise to switch off the appliance.

Note: Do not turn the on/off knob clockwise to the OFF setting.

8 Remove the mains plug from the wall socket.

Cleaning

Never use scouring pads, abrasive cleaning agents or aggressive liquids such as petrol or acetone to clean the appliance.

Do not let water run into the appliance or spill water onto the appliance (Fig. 5).

1 Unplug the appliance and let it cool down.

2 Clean the outside of the appliance with a damp cloth (Fig. 6).

Replacement

Except for the tube lamps, the appliance does not contain any serviceable parts. The tube lamps have to be replaced when they stop burning or after a period of 1,700 burning hours (i.e. after many years of normal use). The tube lamps must be replaced with original Philips EnergyLight PL-L 36W tube lamps.

Have the tube lamps replaced by a service centre authorised by Philips. Only they have the knowledge and skills required for this job and the original spare parts for the appliance.

Environment

- Tube lamps contain substances that may pollute the environment. Always remove the tube lamps before you discard and hand in the appliance at an official collection point. Dispose of the tube lamps at an official collection point for chemical waste.
- Do not throw away the appliance with the normal household waste at the end of its life, but hand it in at an official collection point for recycling. By doing this, you help to preserve the environment (Fig. 7).

Guarantee & service

If you need service or information or if you have a problem, please visit the Philips website at www.philips.com or contact the Philips Consumer Care Centre in your country (you find its phone number in the worldwide guarantee leaflet). If there is no Consumer Care Centre in your country, go to your local Philips dealer.

Technical specification

Model	HF3319
Rated voltage	230V
Rated frequency	50Hz

Model	HF3319
Rated input	72W
Class	II
Classification	Not classified
Ingress of water	IPX0

Operating conditions	
Temperature	+10°C to +30°C
Relative humidity	from 30% to 90%

Storage conditions	
Temperature	from -20°C to +50°C
Relative humidity	from 30% to 90% (no condensation)

Transport conditions	
Temperature	from -20°C to +50°C
Relative humidity	from 30% to 90% (no condensation)

Circuit diagrams, component part lists and other technical descriptions are available on request. As user of this appliance, you should take measures to prevent the occurrence of electromagnetic interference or other types of interference between this appliance and other appliances.

Frequently asked questions

This chapter lists the questions most frequently asked about the appliance. If you cannot find the answer to your question, contact the Consumer Care Centre in your country.

Question	Answer
Can I look directly at the appliance when it is switched on?	Yes, you can look into the light every now and then, but we advise you not to look into the light continuously. It is sufficient to sit close to the lamp.
Are there any side effects to using the EnergyLight?	During the first few days of use, some people may experience minor headaches or strained or watery eyes. These symptoms generally disappear within a few days. If not, consult your doctor before you proceed with your light treatment plan.
How often should I use the EnergyLight?	We recommend that you start taking daily sessions for 5 to 7 days in a row, preferably in the morning. Most people start to feel more energetic within 5 days. If you do not notice any improvement within that time, continue using the EnergyLight for another 7 to 14 days and make the sessions longer (e.g. 3 hours at 50cm). Do not give up too soon, since it may take up to 2 weeks before you notice any significant improvement. However, if using the EnergyLight causes you to become too active, we advise you to reduce the light dose.

Question	Answer
Can I use the EnergyLight if I am on medication?	If you are on prescribed medication, we advise you to consult your doctor first before you use the EnergyLight or undergo any other form of light therapy. Certain medicines cannot be used in combination with light therapy. Anti-depressants combined with light therapy may result in eye damage.
Can I use the EnergyLight if I have eye problems?	If you have eye problems, you should consult your doctor or ophthalmologist before you use the EnergyLight or undergo any other form of light therapy. There are certain eye and eye-related conditions that may respond unfavourably to light therapy.
Can I use the EnergyLight if I have a skin allergy?	If you are allergic to light or have a diagnosed chronic skin disease, it is safest to consult your doctor before you use the EnergyLight or undergo any other form of light therapy.
Can I use the EnergyLight during pregnancy?	Yes, there are no known side effects from using the EnergyLight during pregnancy.
Can I use the EnergyLight if I am allergic to light?	If you are allergic to light, we recommend that you consult your doctor before you use the EnergyLight or undergo any other form of light therapy.
Is it harmful to sit in front of a switched-on EnergyLight for a prolonged period of time?	No, using the EnergyLight longer than the recommended time is not harmful. It may, however, disturb your sleep pattern, particularly if you have your light baths after 8 p.m.
Do I have to wear any sunscreen lotion?	No, the EnergyLight is not a tanning device. It has a special screen that filters out all UV rays; it is therefore safe for your skin and does not give you a tan. If you use the Energy Light in accordance with the instructions in the user manual, you do not have to worry about tanning or other skin reactions.

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