
PHILIPS

HF3418

English 6
Español 17
Français (Canada) 29

IMPORTANT SAFETY INSTRUCTIONS

Read this user manual carefully before you use the appliance and save it for future reference.

Warning (a warning indicates the possibility of injury to the user or operator)

- Water and electricity are a dangerous combination. Do not use this appliance in wet surroundings (e.g. in the bathroom or near a shower or swimming pool). Do not immerse appliance in water or any other liquid.
- Check if the voltage indicated on the appliance corresponds to the local voltage before you connect the appliance.
- The adapter contains a transformer. Do not cut off the adapter to replace it with another adapter or plug, as this may cause a hazardous situation.
- Never use the adapter if it is damaged in any way. Always have the adapter replaced with one of the original type in order to avoid a hazard.
- Do not use the adapter in or near a power outlet that contains an electric air freshener to prevent damage to the adapter.
- Never open the appliance to prevent electric shock, finger burns or temporary black spots in the eye.
- If the power cord is damaged, do not use the appliance. Call 1-866-832-4361 for assistance.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- There are no user-serviceable parts inside the appliance. No modification of this appliance is allowed.
- Never use the appliance if the lamp housing is damaged, broken or missing.
- Do not use the appliance where the light may compromise your or other people's ability to perform essential tasks, such as driving or operating heavy equipment.
- Consult your doctor before you start using the Philips EnergyLight if:
 - you suffer or have suffered from severe depression;
 - you have an eye condition for which your physician advised you to avoid bright light;
 - you are on specific medication (i.e. certain antidepressants, psychotropic drugs or malaria tablets) as it might cause side effects (e.g. skin irritation).

Caution (a caution indicates the possibility of damage to the appliance)

- Always unplug the adapter if not used for longer time.
- Do not use the appliance outdoors or near heated surfaces or heat sources.
- Place the appliance on a stable, level and non-slippery surface, with the stand in place

- Do not place heavy objects on the appliance, the adapter or the cord.
- When you remove the adapter from the wall outlet, do not pull at the cord but grasp the adapter directly to avoid cord damage.
- When you remove the small plug from the appliance socket, do not pull at the cord but grasp the small plug directly to avoid cord damage.
- The appliance has no on/off switch. To disconnect the appliance, remove the plug from the wall socket. The adapter must remain easily accessible at all times.
- Do not impair the cooling of the appliance by covering the appliance with items such as blankets, curtains, clothes, papers, etc.
- Use this appliance for its intended household use as described in this manual. Do not use attachments not recommended by Philips North America LLC.
- If this appliance is damaged in any way or it is not working, it should not be operated. Please visit philips.com/support or call 1-866-832-4361 for assistance.

Contraindications

The following conditions may be contraindications for use of this appliance:

- bipolar disorder
- recent eye surgery or a diagnosed eye condition for which your physician advised you to avoid bright light
- use of photo-sensitizing medications (e.g. certain antidepressants, psychotropic drugs or malaria tablets) as it might cause side effects (e.g. skin irritation)

If any of the above conditions apply to you, always consult your physician before you start using the appliance.

Possible reactions

Possible adverse reactions to light from this appliance are relatively mild and of a transient nature. They may include:

- Headache
- Eyestrain
- Nausea

Measures to avoid, minimize or alleviate these adverse reactions:

- Always use the appliance in a well-lit room.
- Reduce the light intensity if these reactions occur.
- Stop using the appliance for a few days to make the reactions disappear and then try again.

SAVE THESE INSTRUCTIONS

Safety and compliance

- This appliance meets the requirements in the USA and Canada for portable electric luminaires: UL153, UL8750 and CSA 22.2 no.12.

8 English

- This device complies with the EMC requirements of the USA and Canada: FCC part 15, ICES-005.
- This is not a general illumination product.

FCC Information

- This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
 - Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and the receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - Consult the dealer or an experienced radio/TV technician for help.
- This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
- This Class B digital apparatus complies with Canadian ICES-005.

Introduction

Congratulations on your purchase and welcome to Philips! To fully benefit from the support that Philips offers, register your product at www.philips.com/welcome.

Intended use

The Philips EnergyLight is intended to make people feel more energetic, to increase alertness, and to provide mood relief from the winter blues. If you suffer from severe mood swings or depression, make sure to consult a professional, and do not use this device for self-treatment.

Light - a natural energizer

Today's busy modern lifestyles demand peak performance all day long. Even the fittest person suffers energy dips, particularly when deprived of natural daylight, such as in an office in winter. Daylight is essential to well-being, keeping fit, and feeling energized throughout the day – just like healthy nutrition and regular exercise. Philips EnergyLight mimics the natural energizing power of daylight on a bright sunny day. It stimulates special receptors in the eye to trigger the body's natural response to sunlight. This helps you to fight energy dips, to feel more energetic, active and alert, and to improve your mood.

Sleep patterns

The human body uses light cues, such as those provided by the sun, to time certain internal functions. Properly timed rhythms regulate energy, mood, sleep, appetite and digestion. These daily internal cycles, called circadian rhythms, sometimes fall out of sync with nature, which can have impact on our well-being. Nature's cues no longer govern our modern lifestyle. Our day no longer starts at sunrise and ends at sunset, working days are getting longer, and many people work in shifts. Moreover, our activities often extend well into the night. Many people tend to get too little sleep during the week, which diminishes their energy levels during the day. Using Philips EnergyLight at the right time of day can help you establish and reinforce regular sleep patterns.

Light via eyes

Research shows that bright light via the eyes influences the regulation of our energy level, mood and sleep/wake cycle. Indoor spaces generally do not provide enough light intensity. This appliance helps you to be at your best by delivering high quality intense light whenever and wherever you need it.

Winter blues

Lack of light during the darker winter months causes many people to experience a noticeable drop in energy level and low spirits as autumn sets in. These symptoms, known as 'winter blues', may cause you to feel low, listless and fatigued for weeks or even months. This affects both personal performance and mood. Bright light therapy is a successful method to fight these winter blues symptoms. It provides your body with summer vitality throughout the dark winter months.

Features (Fig. 1)

- 1 Light panel
- 2 Socket for small plug
- 3 Time indicator bar
- 4 Standby button
- 5 Stand
- 6 Intensity increase and decrease buttons
- 7 Adapter
- 8 Small plug

Note: Product illustrations are located on the inside of the front and back covers of this instruction booklet.

Using the appliance

When to use

Knowing at what time to use the EnergyLight is essential. It may mean the difference between feeling better in a few days instead of weeks and between experiencing benefits or experiencing none.

Follow these guidelines to understand how light affects us and when to use the EnergyLight for the best results.

Duration

Most people experience benefits if they use the appliance 20 to 30 minutes each day. It is safe to use the EnergyLight longer, provided you do so at the correct time of day. You may notice some energizing effects of light immediately. Other effects become noticeable after several days of regular use.

Indoor life and energy dips

If you spend a lot of time indoors with little daylight and artificial light, you can experience the energizing effects of the Philips EnergyLight very quickly after you switch it on. Many people experience a natural dip in the afternoon. The EnergyLight can help you counteract this dip.

Winter blues

If you experience lower energy levels, an increased need for sleep, and a lower mood in the dark season compared to spring and summer, use the energy light for 20-30 minutes a day. Preferably use it in the morning to help you start the day. You may also use it later during the day, but do not use it in the last few hours before bedtime. If you do, the light could delay your sleep and make you even less energetic the next day.

Sleep patterns

- Morning persons tend to wake up early in the morning and usually find it difficult to stay awake in the evening. If you want to get more out of your evening, avoid bright light directly after waking up and use the EnergyLight in the evening.
- Evening persons find it difficult to get up in the morning and are full of energy in the evening. If this applies to you, you need morning light. Use the EnergyLight within the first 1 to 2 hours after you get out of bed. Avoid bright light in the last 4 hours before bedtime. Daily use of the EnergyLight in the morning may help you establish a regular sleep pattern that allows you to function better in the morning.

Note: If your desired wake-up time is less than six hours after your normal bedtime, do not use the EnergyLight immediately after you wake up. This could give your body a signal that is opposite to the one intended.

Recommended times to use the appliance

When to use

Energy need because of:	Morning	Lunchtime	Afternoon	Evening
Lack of daylight indoors	+++	+++	+++	---
Being a morning person	---	+	++	+++
Being an evening person	+++	++	+	---

When to use

Energy need because of:	Morning	Lunchtime	Afternoon	Evening
An afternoon energy dip	+	+++	+++	---
Fighting the winter blues	+++	+	+	---

- The more pluses (+++) there are in a column, the more suitable the time.
- --- means: do not use the appliance at this time.

Light intensity

Use the EnergyLight in addition to the normal room lighting. The EnergyLight is more comfortable to use in a well-lit room. Adjust the brightness to a comfortable level. Allow a few minutes for your eyes to get used to the bright light. If the highest setting is too bright for you, start using the EnergyLight at a lower setting. Increase the light intensity if you feel no effect after a week.

How to use the appliance

- 1 Insert the small plug in the socket on the back of the appliance (Fig. 2).
- 2 Fold out the stand. The stand can be set to a variable angle.
 - If you are sitting straight up or standing and the top edge of the appliance is below eye level, place the appliance in a tilted position so that its light points further upwards and reaches your eyes (Fig. 3).
 - If you are sitting down at a table or desk and the top edge of the appliance is above eye level, place the appliance in a more upright position (Fig. 4).
- 3 Insert the adapter in a wall outlet (Fig. 5).
- 4 Place the appliance at arm's length (at approximately 20-30 in./50-75cm from your eyes) within your field of vision. The light must bathe your face (Fig. 6).
 - You can read, eat, work at the computer, watch TV or exercise while you use the appliance.
 - For the best results, place the appliance at the same level as your midriff (stomach area), e.g. place it on the table at which you are sitting.
 - Make sure that the surroundings of the appliance are well-lit for more visual comfort.

Note: Do not stare into the light of the appliance. It is sufficient that the bright light reaches your eyes indirectly.

- 5 Press the standby button briefly to switch on the appliance (Fig. 7).
 - The first time you switch on the appliance, the light comes on at the default middle intensity setting (setting 3) of 5 intensity settings.
 - The next time you switch on the appliance, the light comes on at the last used intensity setting.
- 6 You can choose five intensity settings, respectively 20%, 30%, 45%, 67% and 100% intensity.
 - To choose a lower setting, touch the - button briefly (Fig. 8).
 - To choose a higher setting, touch the + button briefly (Fig. 9).

Note: For the best results, use the highest light setting that is comfortable to your eyes.

- 7 The timer indicates the time elapsed from the moment you switch on the appliance.

12 English

- The indicator has 6 segments. Each segment represents 10 minutes. All segments together represent one hour.
- After you switch on the light, the first segment of the time indicator bar flashes (Fig. 10).
- After 10 minutes, the first segment lights up solid and the second segment starts flashing (Fig. 11).
- After 10 more minutes, the second segment lights up solid and the third segment starts flashing.
- This process continues, until all segments of the time indicator bar light up solid after 1 hour (Fig. 12).

Note: If you want to turn off the appliance, press the standby button briefly.

Cleaning

Do not immerse the appliance in water or any other liquid and do not rinse it under the tap.

Never use scouring pads, abrasive cleaning agents or aggressive liquids such as gasoline or acetone to clean the product.

- 1 Clean the appliance with a dry or moist cloth. Do not use a wet cloth (Fig. 13).

This product has no other user-serviceable parts. For assistance call 1-866-832-4361.

Storage

- If the appliance will not be used for an extended period of time, remove the power cord from the outlet and store in a safe, dry location where it will not be crushed, banded, or subject to damage.
- Do not wrap cord around the appliance when storing.

Disposal

Your product is designed and manufactured with high quality materials and components, which can be recycled and reused. For recycling information, please contact your local waste management facilities or visit www.recycle.philips.com.

Assistance

For assistance visit our website: www.philips.com/support or call toll free: 1-866-832-4361. Online information is available 24 hours a day, 7 days a week.

90-Day Money-Back Guarantee

If you are not fully satisfied with your product send the product back and we'll refund you the full purchase price.

The product must be shipped prepaid by insured mail, insurance prepaid, have the original sales receipt, indicating purchase price and date of purchase, and the money-back guarantee return authorization form enclosed.

We cannot be responsible for lost mail.

The product must be postmarked no later than 90 days after the date of purchase. Philips reserves the right to verify the purchase price of the product and limit refunds not to exceed suggested retail price.

To obtain money-back guarantee return authorization form, call 1-866-832-4361 for assistance.

Delivery of refund check will occur 6–8 weeks after receipt of returned product.

Full Two-Year Warranty

Philips North America LLC (USA) and Philips Electronics Ltd (Canada) warrant this new product against defects in materials or workmanship for a period of two years from the date of purchase, and agree to repair or replace any defective product without charge.

IMPORTANT: This warranty does not cover damage resulting from accident, misuse or abuse, lack of reasonable care, or the affixing of any attachment not provided with the product or loss of parts or subjecting the product to any but the specified voltage or batteries.*

NO RESPONSIBILITY IS ASSUMED FOR ANY SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES.

In order to obtain warranty service, go to www.philips.com/support for assistance. It is suggested that for your protection you return shipments of product by insured mail, insurance prepaid. Damage occurring during shipment is not covered by this warranty.

NOTE: No other warranty, written or oral, is authorized by Philips North America LLC.

This warranty gives you specific legal rights, and you may also have other rights which vary from state to state, province to province or country to country. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above exclusion and limitations may not apply to you.

* Read enclosed instructions carefully.

In US manufactured for: Philips Personal Health, a division of Philips North America LLC, P.O. Box 10313, Stamford, CT 06904.

In Canada imported by: Philips Electronics LTD, 281 Hillmount Road, Markham, Ontario L6C 2S3.

PHILIPS and Philips shield are registered trademarks of Koninklijke Philips N.V.

©2018 Philips North America LLC. All rights reserved.

Frequently asked questions

This chapter lists the most frequently asked questions about the appliance. If you cannot find the answer to your question, call 1-866-832-4361 for assistance.

Question	Answer
I have heard that LED light can be dangerous. Can this appliance damage my eyes?	Philips EnergyLight complies with the photobiological safety standard IEC62471. According to this standard, the appliance represents no risk to healthy eyes and can be safely used.
Can I get sunburn from using this appliance?	No, sunburn is caused by UV light. The light produced by this appliance does not contain UV.
Is the Philips EnergyLight safe to use around children and pets?	The light is safe to be used around children and pets. However, do not let children play with the appliance.
Can I use the Philips EnergyLight when I am pregnant?	Yes, there are no known side effect from using the Philips EnergyLight during pregnancy.
Do I have to take off my glasses when I use the appliance?	No, you can use the appliance with your glasses on. Do not wear tinted glasses or sunglasses when you use the appliance, as these glasses reduce its effectiveness.
How often do I have to use the EnergyLight?	Use the Philips EnergyLight every day at the right time to achieve the result you want. For more information on the duration of a session, read chapter 'Using the appliance'. Except when you try to adjust your sleep time, we advise you to use the EnergyLight at the same time every day. It is usually sufficient to use it 20 to 30 minutes a day. It is absolutely safe to use it longer, if you want to.
When do I start to notice improvement from using the Philips EnergyLight?	You may start to feel the energizing effect any time between first use and 1 to 2 weeks of daily use.
Do I have to sit close to the appliance to notice its effect?	No, we advise you to place the appliance at arm's length. If you use the appliance at a distance of 20-30 in. (50-75 cm), you can expect the full benefit of its bright light.
Is it harmful to sit in front of a switched-on EnergyLight for a long time?	No, the Philips EnergyLight is absolutely safe, also when you use it throughout the day. However, in general avoid using the light in the last 4 hours before you go to bed, as the energizing effect of the light may interfere with your sleep. Only 'morning persons' who want to get more out of their evening can use it in the evening, see chapter 'Using the appliance', section 'Sleep patterns'.
Can the lamp of the EnergyLight be replaced?	No, the Philips EnergyLight has a durable LED lamp that lasts many years.

Specifications

Model HF3418

Electrical

Rated input voltage adapter 100–240 Vac

Rated input frequency adapter 50/60 Hz

No-load power adapter < 0.1 W

Rated output voltage adapter 24 Vdc

Rated output power adapter 24 W

Standby power appliance < 0.25 W

Classifications

Portable Electric Luminaires UL153, UL8750 and CSA 22.2 no.12

Adapter Class 2 power supply

Ingress of an object and water IP31 (main unit only)

Mode of operation Continuous

Operating conditions

Temperature from 5°C to +35°C (41°F to 95°F)

Relative humidity from 15% to 90% (no condensation)

Storage and transport conditions

Temperature from -25°C to +70°C (-13°F to +158°F)

Relative humidity from 15% to 90% (no condensation)

Physical characteristics

Dimensions 29.4 x 37.6 x 4.1 cm (11.6 x 14.8 x 1.6 in.)

Weight 2.0kg (70 oz.)

Color temperature 4000–5000K

Color rendering index > 70

Nominal output (100% intensity) 10.000 lux at 17.5 cm (10.000 lux at 6.8 in.)

Nominal output (100% intensity) 2750 lux at 50 cm (2750 lux at 20 in.)

Symbol key

The following symbols may appear on the appliance:

Symbol	Description
	Consult the user manual supplied
	DC power
	DC jack polarity
IP31	Degree of protection against liquid and particles.
HF3418	Model number
	UL listing mark for Canada and the United States.
	Philips shield
	In China manufactured for Philips Personal Health, a division of Philips North America LLC, P.O. Box 10313, Stamford, CT 06904. In Canada imported for: Philips Electronics LTD, 281 Hillmount Road, Markham, Ontario L6C 2S3.
	Date code and serial number of this appliance
	Standby button
	Decreasing light intensity button
	Increasing light intensity button
	Compliant with the Waste Electrical and Electronic Equipment/Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment (WEEE) recycling directives.
CAN ICES-005 (B)/NMB-5 (B)	Appliance complies to Canadian EMC directive ICES-005
$t_a=35^{\circ}\text{C}$	Maximum ambient temperature while operating.

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Importé au Canada pour : Philips Electronics LTD, 281 Hillmount Road, Markham, Ontario L6C 2S3.