
PHILIPS

HF 3301


PHILIPS

- English** **Page 4**
• Keep page 3 open when reading these operating instructions.
- Français** **Page 7**
• Pour le mode d'emploi: dépliez la page 3.
- Deutsch** **Seite 10**
• Schlagen Sie beim Lesen der Gebrauchsanweisung Seite 3 auf.
- Nederlands** **Pagina 13**
• Sla bij het lezen van de gebruiksaanwijzing pagina 3 op.
- Italiano** **Pagina 16**
• Prima di usare l'apparecchio, aprite la pagina 3e leggete attentamente le istruzioni per l'uso.
- Dansk** **Side 19**
• Hold side 3 opslået mens De læser brugsanvisningen.
- Norsk** **Side 22**
• Slå opp på side 3 før De leser videre.
- Svenska** **Sid 25**
• Ha sidan 3 utvikt när ni läser bruksanvisningen.
- Suomi** **Sivu 28**
• Pidä sivu 3 auki samalla kun luet käyttöohjetta.

1

Important

- The appliance is suitable for both normal lighting and a light bath. The light bath is used to counteract the tiredness and lowered mood which can be caused by a derangement of the “biologic clock”. This derangement can be due to lack of bright light and the shortening of the days during winter and autumn.
- Before connecting the appliance to the mains socket, check if the voltage indicated on the rating plate corresponds with the mains voltage in your home.
- Water and electricity are a dangerous combination! Therefore, do not use the appliance in wet surroundings (e.g. in the bathroom or shower, or near a swimming pool).
- Remove the plug from the wall socket after use.
- Take care that children cannot play with the appliance.
- The advisable duration of the lightbath depends on the selected operational distance. (Please refer to section “Light bath”.)

Preferably have your lightbaths during the day or early in the evening, between 6 o'clock in the morning and 8 o'clock in the evening.

Repeat the lightbath on at least five successive days.



Besides, you can take lightbaths now and then whenever you feel like you need them, or start new courses (of five successive days or more).

- Use the light of this appliance as additional to the normal light conditions (sunlight, electric light) in your home. Do not use the appliance in a furthermore dark room.
- It is not advisable to keep staring into the lamps.
- After using the appliance for the first few times, some eye-stress and/or a light headache may occur. Usually these side-effects disappear soon.
- If using specific medicines (i.e. certain antidepressants, other psychotropic drugs) more serious complaints may happen. First consult your doctor if you take such drugs.
- With eye complaints or hypersensitivity to light it is advisable to contact your doctor first.
- If using the appliance should raise your activity too much, we advise you to reduce the light dose.
- If there is no improvement we advise you to consult your doctor.
- If you feel exhausted, distressed or restless and/or suffer from sleeping problems (e.g. insomnia), it is always advisable to consult your doctor or psychotherapist. Your problem could also be caused by something else.


General description (fig. 1)

- A Tube lamps (type Philips Original Bright Light PL-L 55 W)
- B Reflector
- C Opalescent/satin screen (acrylic)
- D On/Off switch
- E Suspension hole
- F Stand (fold-out)
- G Handle (for transport)
- H Screw and plug (for hanging the appliance)



Making the appliance ready for use


- Hang the appliance on the wall (fig. 2) or unfold the stand (fig. 3) and place it on a table top.
- Insert the plug into a wall socket (220-240V) (fig. 4).

Light bath

- Switch on (fig. 5).
- Position yourself in front of the appliance (fig. 6). For a light bath, the distance between your face and the opalescent/satin screen should be between 25 and 75 cm.


- The advisable duration of the lightbath depends on the selected operational distance:
 - approx. **2 hours** for a distance of approx. **75 cm** (fig. 7);
 - approx. **1 hour** for a distance of approx. **50 cm** (fig. 8);
 - approx. **½ hour** for a distance of approx. **25 cm** (fig. 9).


- You need not take the whole lightbath without a break. If it suits your schedule better, you may interrupt your lightbath and continue it later. Even the 30 minute session may be split into two sessions if you wish.
- You may look into the lamps every now and then, not continuously. So you are free to engage in such activities as reading, writing, fancy-work or handiwork, and look into the lamps every now and then (figs. 10, 11).

- Preferably have your lightbaths between 8 o'clock in the morning (a.m.) and 8 o'clock in the evening (p.m.) (fig. 12).



- Repeat the lightbaths until you feel better. Generally, improvements occur within one week. So as a rule, we advise you to repeat the lightbath on at least 5 successive days (fig. 13).

- You may repeat the light bath as much as you like, during the 'dark period' of the year. Do repeat the lightbath on at least five successive days.


After use



- Switch the appliance off (fig. 14) and remove the plug from the wall socket.
- You can lift the appliance by its handle (fig. 15).

Cleaning

Always unplug the appliance and let it cool down before cleaning.

- You can clean the outside of the appliance with a damp cloth (fig. 16). Make sure that no water runs into the appliance (fig. 17). Do not use abrasives or scourers, or liquids such as alcohol, petrol or acetone.


- From time to time the tube lamps, reflectors and the inner side of the opalescent screen need cleaning (figs. 18-21).
 - Remove the screen by pressing the clips gently with a screwdriver (fig. 20).

Replacing the tube lamps

After approx. 10.000 hours (many years in normal use) we advise you to replace the tube lamps.

Have the tube lamps replaced only by Philips Original Bright Light PL-L 55W.

You can have the tube lamps replaced by authorized service personnel or replace them yourself. If you wish to replace the lamps yourself:

- Remove the screen (figs. 18 and 19).
- Bend the bracket, which keeps the lamp to its place, a little sideways and remove it (fig. 22).
- Press the red button and pull the tube lamp forward, out of the lampholder (fig. 23).
- When *placing* a tube lamp, follow these instructions in reverse order (figs. 24 and 25).

Note: While sliding the tube light into the lampholder you need not press the red button. As soon as you hear a 'Click!', the lamp is placed correctly.

The tube lamps contain substances which may pollute the environment. When discarding lamps, please ensure that they are separated from the normal household refuse and that they are disposed of at an officially assigned place.

Replacing the mains cord

If the mains cord of this appliance is damaged, it must only be replaced by Philips or their service representative, as special tools and/or parts are required.

This appliance meets the requirements of MDD 93/42/EEC for medical appliances.

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