

References

1. Pokorney et al. Outcomes Associated with Extraction versus Capping and Abandoning Pacing and Defibrillator Leads Circulation 2017 Oct 10;136(15):1387-1395. doi: 10.1161/CIRCULATIONAHA.117.027636. Epub 2017 Aug 22.
2. Wilkoff, Bruce L., et al. How to treat and identify device infections. Heart Rhythm, Vol 4, No 11, 2007, 1467-1470.
3. Hussein, Ayman A., et al. "Cardiac Implantable Electronic Device Infections: Added Complexity and Suboptimal Outcomes With Previously Abandoned Leads." JACC: Clinical Electrophysiology (2016).
4. Byrd, CL, et al. Intravascular extraction of problematic or infected permanent pacemaker leads: 1994-1996. U.S. Extraction Database, MED Institute. PACE 1999; 22:1348-1357.
5. Maytin, M. MD, et al. Multicenter Experience with Extraction of the Riata™/Riata™ ST ICD Lead, Heart Rhythm, <http://dx.doi.org/10.1016/j.hrthm.2014.05.014>.
6. de Cock CC, et al. Long-term outcome of patients with multiple (> or = 3) noninfected leads: a clinical and echocardiographic study. PACE, Vol 23, No 4, 2000, 423-6.
7. Mattei, E., Gentili, G., Censi, F., Triventi, M. and Calcagnini, G. (2015), Impact of capped and uncapped abandoned leads on the heating of an MR-conditional pacemaker implant. Magn Reson Med, 73: 390-400. doi: 10.1002/mrm.25106.
8. Kalin R, Stanton MS. Current clinical issues for MRI scanning of pacemaker and defibrillator patients. PACE. April 2005;28(4):326-328.
9. Barakat, Amr F., et al. "Transvenous lead extraction at the time of cardiac implantable electronic device upgrade: Complexity, safety, and outcomes." Heart rhythm 14.12 (2017): 1807-1811.
10. Wazni, O et. al. Lead Extraction in the Contemporary Setting: The LExiCon Study: A Multicenter Observational Retrospective Study of Consecutive Laser Lead Extractions, J Am Coll Cardiol, 55:579-586.
11. Document on file D027561. When deployed, the Bridge Occlusion Balloon reduces blood loss by up to 90%, on average, in an animal model of an SVC tear. Testing was conducted in a heparinized porcine model which has shorter SVC length than is typical in humans. A balloon design scaled for use specifically in the porcine model was used in generating this data.
12. Document on file, D026197. In an animal model with SVC tears up to 3.5 cm, with 2 pacing leads and 1 ICD lead.
13. Roger G. Carrillo, MD; Darren C. Tsang, BS; Ryan Azarrafy, BA; Thomas A. Boyle, BS. Multi-Year Evaluation of Compliant Endovascular Balloon in Treating Superior Vena Cava Tears During Transvenous Lead Extraction. EHRA late-breaking trial, March 19, 2018.
14. Kusumoto et al. 2017 HRS Expert Consensus Statement on Cardiovascular Implantable Electronic Device Lead Management and Extraction. Heart Rhythm, 2017.

Important safety information

Bridge occlusion balloon

The Bridge occlusion balloon is indicated for use for temporary vessel occlusion of the superior vena cava in applications including perioperative occlusion and emergency control of hemorrhage. Use of the Bridge occlusion balloon in procedures other than those indicated is not recommended.

The adverse events associated with an occlusion balloon procedure include, but are not limited to allergic reactions, death, embolization, hematoma, hemorrhage, sepsis/infection, short-term hemodynamic deterioration, thromboembolic episodes, vascular thrombosis, vessel dissection, vessel perforation, vessel spasm.

In order to facilitate rapid delivery, it is recommended that a guidewire is in place in the superior vena cava prior to beginning the lead extraction procedure. Attempting to place the guidewire after a tear has occurred may:

- Result in an inability to traverse the superior vena cava with the guidewire
- Result in the guidewire exiting the vasculature at the tear site
- Result in an inability to place the Bridge occlusion balloon catheter
- Delay or prevent the ability to achieve occlusion

This information is not intended to replace a discussion with your healthcare provider on the benefits and risks of this procedure to you.

GlideLight laser sheath

The GlideLight laser sheath is intended for use with other lead extraction tools in patients who are suitable candidates for removal of implanted pacemaker and defibrillator leads. The use of the GlideLight laser sheath may be unsafe in some patients, or with certain leads, or when the leads cannot be extracted through the superior veins (that is, when groin or surgical extraction is required). Rarely a patient undergoing lead extraction may require urgent surgical treatment for a complication; therefore, patients should not undergo lead extraction with a laser sheath in centers where emergency surgical procedures cannot be performed. Leads not intended for extraction may be damaged during the procedure and may require replacement. Ask your doctor if you are a candidate for lead extraction with the GlideLight laser sheath.

Potential minor adverse events associated with lead extraction procedures that may or may not require medical or surgical treatment include: a tear or damage to the blood vessels, the heart or its structures; bleeding at the surgical site; or collapsed lung.

Rare but serious adverse events that require emergency medical or surgical procedures may include: a tear or damage to the blood vessels, the heart, lungs or their structures; blood clot or obstruction of the blood vessels or lungs by debris or lead fragments. Other serious complications may include: irregular heartbeat, weakened heart muscle, infection, respiratory failure or complications associated with anesthesia, stroke or death.

This information is not intended to replace a discussion with your healthcare provider on the benefits and risks of this procedure to you.

Caution: Federal law restricts this device to sale by or on the order of a physician.



Philips

3721 Valley Centre Drive, Suite 500, San Diego, CA 92130 USA
www.philips.com/IGTdevices

©2018 Koninklijke Philips N.V. All rights reserved. Some or all products manufactured by Spectranetics, a Philips company. Approved for external distribution. D016277-02 102018

PHILIPS

Lead Management



Risks and considerations

of cardiac lead extraction

Risks of capping a lead

Increased risk of infection^{2,3,4}

- Increased risk of infection five years post-procedure¹
- Risk of infection increases 2-7% at each device change²

Worse outcomes and increased difficulty of future extraction

- 2.6 times more likely to have a major adverse event (MAE)³
- 4 times more likely to have infected lead material retention, leading to a 2x increase in one-month mortality³
- Risk of failed lead removal doubles every three years⁴

Increased risk of venous occlusion and tricuspid regurgitation^{5,6}

- Can prevent access for new leads
- Risk of increased tricuspid regurgitation and resulting atrial fibrillation or right-sided heart failure

MRI contraindication⁷

- 75% of device patients will need an MRI in their lifetime⁸
- Abandoned leads are a risk for tissue damage or inappropriate cardiac stimulation⁷

Considerations of lead extraction

Extraction safety at device upgrade⁹

- Less complex and more likely to be successful, with a 97.2% clinical success rate
- Lower complication rates, despite significantly worse clinical profiles, with a 0.4% mortality rate and a 1% MAE rate

Overall procedural safety of lead extraction²

- Laser-assisted lead extraction is clinically proven safe, with a 97.7% clinical success rate, 1.4% procedural MAE rate, and 0.28% procedural mortality rate¹⁰

New innovations designed for safety

- In the rare event of an superior vena cava (SVC) tear, the Philips Bridge occlusion balloon can reduce blood loss by 90%¹¹ on average and provide 30 minutes of acceptable hemostasis¹²
- With proper use of the device, SVC tear survival has gone from 56.4% to 91.7%¹³

Overall, lead extraction has a

99.72%
procedural
survival
rate¹⁰

Shared decision-making

- It is an HRS Class I indication¹⁴ for physicians and patients to discuss the risks of lead abandonment and the risks of lead extraction.
- Operator-specific information about success rates, case volume, and complication rates should be discussed with patients prior to deciding to proceed with an extraction procedure (Class I).
- Extraction may be considered after shared decision-making process with patients (Class IIb).¹⁴

No better time to extract than now

Capping and abandoning leads poses significant risks that can be mitigated proactively with safe lead extraction.

