

**PHILIPS**

Stellarex

Drug-coated 0.035"  
angioplasty balloon

# The clear DCB choice



Featuring EnduraCoat  
technology



# Treatment that endures even in complex patients


Stellarex exhibits consistent three-year patency across trials—even in the most complex cases. Now, all patients can experience the power of proven results:

- 64.2% patency in ILLUMENATE Pivotal proves only durable three-year RCT results in complex patients
- 43.9% severe calcium is four or five times the rate of severe calcium studied in competitive trials<sup>4</sup>

Top-tier primary patency in complex patients with up to

43%


less drug load<sup>1,3</sup>



Competitor studies are independent clinical trials with different protocols and definitions. Therefore, they are not head-to-head comparisons, and data presented cannot be directly compared. Calcium definitions may vary from study to study, and the rates presented here are based on those used and reported in each respective study. Complex patients refers to high rates of severe calcium, diabetes and renal insufficiency. Primary patency based on Kaplan-Meier estimates.

## Show complex patients a better way

Stellarex continues to perform in patients with the highest rates of complex comorbidities—even through three years.




Data overview for informational purposes only and not for head-to-head comparison. Calcium definitions may vary from study to study, and the rates presented here are based on those used and reported in each respective study.

## Calcium can be beat


Stellarex is the only DCB reported to have durable patency at 2 years in severely calcified lesions—showing patients with severe calcium can have a meaningful treatment effect with Stellarex.

- 67.8% two-year patency in severely calcified lesions
- Demonstrates a robust treatment effect over PTA in severely calcified lesions<sup>1</sup>
- 82% of severely calcified segments are  $\geq 5$  cm long<sup>4</sup>



## The smart treatment choice for today's PAD patients

The prevalence of calcified lesions increases with age and diabetes.<sup>5</sup> Shown to work in calcium, Stellarex is the clear treatment choice for PAD.



## Pay for the DCB that works in calcium

Only Stellarex gives you confidence that your DCB use and spend in severely calcified lesions will be worthwhile and beneficial to patients.

	Incidence of severe calcium in RCT	Two-year patency rate in severe Ca++	DCB shown effective in RCT with severe Ca++	DCB provides treatment confidence in severe Ca++
Stellarex <sup>1</sup>	43.9%	67.8%	Yes	Yes
In.Pact <sup>3</sup>	8.1%	N/A	?	?
Lutonix <sup>6</sup>	10.4%	N/A	?	?
Ranger <sup>14</sup>	11.5%	N/A	?	?

Overview is provided for informational purposes only and not for head-to-head comparisons. Protocols and definitions may vary from study to study.




# No mortality difference for long quality of life

Overlapping survival curves demonstrate Stellarex has no statistical mortality difference compared to PTA at any follow up time point through five years.<sup>2</sup>

Only DCB with  
no mortality  
difference through


5  
Years



Stellarex treated patients with higher rates of comorbidities and likely sicker patients and still had no difference in mortality compared to PTA.

## Efficiency as easy as 1-2-3

**1. Limit wire exchanges**—Stellarex tracks easily over your choice of guide wire, including 0.035, 0.018 and 0.014, potentially reducing the need for guide wire exchanges.<sup>10</sup>



**Stellarex**  
Over guide wire

**2. Fast tracking**—Stellarex distorts less in tight bends, and the strong low-profile tip inhibits flare, minimizing catch on lesions and maximizing pushability and trackability.<sup>10</sup>



**Stellarex**  
Lateral view



**Stellarex**  
Axial view



**In.Pact**  
Axial view



**In.Pact**  
Lateral view

**3. Less waiting**—At just 60 seconds, Stellarex recommends the shortest inflation hold time of any DCB

**Stellarex**  
1-minute  
inflation hold<sup>8</sup>

**Lutonix**  
2-minute  
inflation hold<sup>6</sup>

**In.Pact and Ranger**  
3-minute  
inflation hold<sup>14,17</sup>



# Differentiated technology— next-generation EnduraCoat


Stellarex EnduraCoat was designed for performance in complex and severely calcified lesions and patients with multiple comorbidities.

- Hybrid paclitaxel offers prompt drug transfer and sustained tissue residency through 28 day restenotic window<sup>12</sup>
- Excipient polyethylene glycol (PEG) offers excellent adhesion and durability to protect low dose paclitaxel<sup>11</sup>
- Reduces drug loss during transit, relieving clinicians of transit time requirements<sup>8,10</sup>

## Designed for performance in calcium



## High transfer efficiency and effective residency<sup>13,15</sup>



PEG forms strong ionic bonds with hydroxyl apatite (HAp), the primary component of calcified atherosclerotic lesions.<sup>9</sup>

PEG's affinity for HAp may result in limited PTX washout in the presence of calcium.

PEG may protect PTX, giving it time to be absorbed into vessel when calcium is present.

Hybrid  
paclitaxel



PEG  
excipient



Top-tier clinical  
outcomes

## Why an effective low drug dose matters

Dose excess and particulate downstream possibly results in a delay of wound healing, loss of microcirculation and creation of aneurysms. Stellarex is the only low dose DCB with a statistically significant treatment effect at two years.



- In.Pact has a 75% higher drug dose than Stellarex<sup>3,8</sup>
- Compared to Stellarex, In.Pact loses 2.7 times more drug ( $\mu\text{g}$ ) during tracking to the deployment site<sup>10</sup>
- In.Pact coating visually flakes off during device prep<sup>10</sup>
- Lutonix low dose is mostly amorphous paclitaxel, which may lead to short-term tissue residency<sup>11</sup>
- Ranger mostly crystalline paclitaxel<sup>16</sup> requires loading tool to limit drug loss<sup>14</sup>




Image on file.

## Track

PEG offers exceptional durability during handling, tracking and inflation, helping prevent premature drug loss<sup>10,11</sup>

## Deliver

EnduraCoat achieves uniform and efficient drug transfer<sup>13</sup>

## Sustain

Hybrid paclitaxel provides prompt drug availability and sustained tissue residency<sup>12</sup>

# Stellarex 0.035" OTW drug-coated angioplasty balloon

Product catalog number	Sheath size (Fr)	Balloon diameter (mm)	Balloon length (mm)	Shaft length (cm)	Nominal pressure (atm)	Rated burst pressure (atm)
AB35SX040040135	6	4	40	135	10	20
AB35SX040060135	6	4	60	135	10	20
AB35SX040080135	6	4	80	135	10	20
AB35SX040100135	6	4	100	135	10	20
AB35SX040120135	6	4	120	135	10	20
AB35SX040150135	6	4	150	135	10	20
AB35SX040200135	6	4	200	135	10	20
AB35SX050040135	6	5	40	135	10	18
AB35SX050060135	6	5	60	135	10	18
AB35SX050080135	6	5	80	135	10	18
AB35SX050100135	6	5	100	135	10	18
AB35SX050120135	6	5	120	135	10	16
AB35SX050150135	6	5	150	135	10	16
AB35SX050200135	6	5	200	135	10	16
AB35SX060040135	6	6	40	135	8	14
AB35SX060060135	6	6	60	135	8	14
AB35SX060080135	6	6	80	135	8	14
AB35SX060100135	6	6	100	135	8	14
AB35SX060120135	6	6	120	135	8	12
AB35SX060150135	6	6	150	135	8	12
AB35SX060200135	6	6	200	135	8	11

## Important safety information

The Stellarex 0.035" OTW drug-coated angioplasty balloon is indicated for percutaneous transluminal angioplasty (PTA), after appropriate vessel preparation of de novo or restenotic lesions up to 180 mm in length in native superficial femoral or popliteal arteries with reference vessel diameters of 4-6 mm.

The Stellarex 0.035" OTW drug-coated angioplasty balloon is contraindicated for use in:

- Patients with known hypersensitivity to paclitaxel or structurally related compounds
- Patients who cannot receive recommended antiplatelet and/or anticoagulation therapy
- Women who are breastfeeding, pregnant or are intending to become pregnant, or men intending to father children
- Coronary arteries, renal arteries and supra-aortic/cerebrovascular arteries
- Patients judged to have a lesion that prevents complete inflation of an angioplasty balloon or proper placement of the delivery system

Possible adverse effects associated with the balloon dilation procedure include, but are not limited to: Abrupt vessel closure; Allergic reaction to contrast medium, antiplatelet therapy or catheter system components (drug, excipients and materials); Amputation/Loss of limb; Arrhythmias; Arterial aneurysm; Thrombosis; Arterio-venous fistula (AVF); Bleeding; Death; Embolism/Device embolism; Fever; Hematoma; Hemorrhage; Hypertension/Hypotension; Infection or pain at insertion site; Inflammation; Ischemia or infarction of tissue/organ; Occlusion; Pain or tenderness; Peripheral edema; Pseudoaneurysm; Renal insufficiency or failure; Restenosis; Sepsis or systemic infection; Shock; Stroke/Cerebrovascular accident; Vessel dissection, perforation, rupture, spasm or recoil; Vessel trauma that requires surgical repair; Balloon rupture; Detachment of a component of the balloon and/or catheter system; Failure of the balloon to perform as intended; Failure to cross the lesion.

Additional complications that may be associated with the addition of paclitaxel to the balloon include, but may not be limited to the following: Allergic/Immunologic reaction to paclitaxel; Alopecia; Anemia; Gastrointestinal symptoms (diarrhea, nausea, pain, vomiting); Hematologic dyscrasia (including neutropenia, leukopenia, thrombocytopenia); Hepatic enzyme changes; Histologic changes in vessel wall including inflammation, cellular damage or necrosis; Myalgia/Arthralgia; Myelosuppression; Peripheral neuropathy. Caution: Federal law restricts this device to sales by or on the order of a physician.

## References

1. Mathews J. Stellarex in the Treatment of the SFA and Popliteal: Late-Breaking 3-Year Data. NCVH, 2019. May 29, 2019. New Orleans, LA.
2. Gray, W. MD. Long-Term Safety of Stellarex DCB Across ILLUMENATE RCTs: 5-year results from a patient-level meta-analysis of the ILLUMENATE RCTs. Charing Cross 2021. April 20, 2021.
3. Schneider, P. Treatment Effect of Drug-Coated Balloons Is Durable to 3 Years in the Femoropopliteal Arteries: Long-Term Results of the IN.PACT SFA Randomized Trial. Circ. Card Int. 2018;11:1: e005891
4. Holden. Comparing Trials Data in the Management of Calcified Arteries. Charing Cross 2018. April 24-26, 2018; London, UK.
5. Rocha-Singh, et al. Peripheral Arterial Calcification: Prevalence, Mechanism, Detection and Clinical Implications. Cath and Cardiovasc Intervent. 2014;83:E212-E220.
6. Bard Lutonix Instructions for Use, BAW1387400r3.
7. Lyden, S. ILLUMENATE Trial Update: 3-year ILLUMENATE Global Stellarex DCB Data. AMP 2019. Chicago, IL.
8. Stellarex IFU. P011966-D.
9. Venkatasubbu GD, et al. Surface modification and paclitaxel drug delivery of folic acid modified polyethylene glycol functionalized hydroxyapatite nanoparticles. Powder Technology. 2013;235:437-442.
10. Data on file. D044595-00.
11. Mark J, et al. Physical properties of polymers. Cambridge University Press. 3rd ed. 2004.
12. Granada JF. Future directions, clinical applications and local drug delivery technologies. Presented at the Transcatheter Cardiovascular Therapeutics (TCT) 25th Annual Scientific Symposium; October 27-November 1, 2013; San Francisco, California. Slide 7, 15.
13. Superimposed PK curves from different datasets: Philips internal data on file D029433. Yazdani SK, Pacheco E, Nakano M, et al. Vascular, downstream, and pharmacokinetic responses to treatment with a low dose drug-coated balloon in a swine femoral artery model. Catheter Cardiovasc Interv. 2014;83:132-140. Melder R. Presented at: EuroPCR; May 15-18, 2012; Paris: France.
14. Ranger IFU 50589598-01.
15. Gongora, CA. Impact of Paclitaxel Dose on Tissue Pharmacokinetics and Vascular Healing: A Comparative Drug-Coated Balloon Study in the Familial Hypercholesterolemic Swine Model of Superficial Femoral In-Stent Restenosis. JACC Cardiovasc Interv. JACC: Cardiovascular Interventions, Volume 8, Issue 8, 2015, Pages 1115-1123.
16. <https://www.bostonscientific.com/en-EU/products/balloons--drug-coated/agent/trans-pax-efficient-coating-technology.html> Retrieved 3/11/2021.
17. Medtronic In.Pact Instructions for Use, M052624T001 Rev 1F.

**StellarexDCB.com**

