

Magnetic Resonance

Integrated clinical solutions

MR Workspace and clinical application portfolio

Extending the **power of MR**

For more than 130 years, we have been creating meaningful innovations to improve people's lives and make the world more sustainable. We are inspired to continually advance the state of precision diagnosis with customer- and patientcentric solutions that deliver clear care pathways and predicable outcomes. In MR, our mission is to achieve a fast, fully automated, and personalized exam for every patient, while acting responsibly towards our planet and society. With AI*-driven smart connected imaging, optimized workflows, and integrated clinical solutions, we strive to improve MR department's productivity, enhance patient and staff experience, and deliver high quality diagnostic imaging.

Our portfolio of MR clinical solutions can help to break down diagnostic boundaries by delivering speed, comfort and confidence, all with the aim to help improve patient care.

Disclaimer: All products, features and applications described in this brochure are subject to availability depending on country regulatory approvals and system compatibility

^{*}According to the definition of AI from the EU High-Level Expert Group

MR Workspace 4

- dS Scan Suite Essential 3
- AV MR Essential
- SmartExams
 8
- Smart WorkFlow
 9
- SmartSpeed Plus 10

11

15

- SmartSpeed Pro Body
- SmartSpeed Pro Cardiac 12
- SmartSpeed Premium 13
- dS IQ Suite 14
- dS Vascular Suite Pro

• dS Diffusion Suite F	Pro 16
• AV Onco Diffusion	Pro 17
• dS Spine Suite Pro	18
• dS Neuro Suite	19
• AV Neuro	20
• dS Spectro Suite Pi	ro 21
• dS Cardiac Suite	22
• dS Breast Suite Pro	23
• dS Liver Suite Prem	nium 24
• dS Pediatric Suite	25
• A la Carte	



With MR Workspace we aim to support you to:



Increase schedule efficiency



Deliver consistent image quality



Improve staff experience



Reduce training time



Provide faster time to results

MR Workspace

MR Workspace supports efficiency and staff satisfaction in the control room through intelligence, guidance and ease of use. Technologists can prepare exams before patients arrive and aim to achieve consistent quality regardless of experience, by using Protocol Assistant, an Al¹-driven solution that learns your protocol preferences and suggests the most appropriate ones based on clinical indication. Advanced visualization includes step-by-step guidance so technologists can begin post-processing, delivering more² diagnostic information to the radiologist, thus saving reporting time. The intuitive interface, large display of clinical images and essential parameter reveal contribute to outstanding ease-of-use.

MR Workspace also supports a large collection of clinical application suites. It is the Philips platform of choice for existing and future clinical innovations.

1 According to the definition of AI from the EU High-Level Expert Group.

2 The addition of step-by-step guidance and automation of routine and complex post-processing applications can now be performed

by the technologist on the console, saved via bookmark functionality, and handed off to the radiologist, which reduces time to results.

dS Scan Suite Essential



ScanTools Pro Clinical imaging performance



AutoVoice Guiding your patients through the MR examination



ComforTone Reduce acoustic noise and enhance MR experience



ScanWise Implant A key to confidence with MR Conditional implants



3D MSK VIEW View your 3D TSE imaging data in any plane



3D BrainVIEW View your 3D TSE imaging data in any plane



Whole Body Get comfortable body imaging with head-to-toe coverage



bFFE XD Expanding clinical applications of balanced FFE



2k imaging High resolution imaging



O-MAR Efficient near-metal soft tissu and bone imaging

AV MR Essential



MR MobiView Easy viewing of multistation studies



MR T1 Perfusion Review blood supply characteristics



MR Subtraction Perform basic calculations between two volumes



MR Diffusion Analyze diffusion and anisotropic properties of tissue



MR Echo Accumulation Calculation of new images based on series with multiple echoes

SmartExams

Plus

Pro



SmartExam Brain Standardized exams forconsistent MRI results



SmartExam Shoulder Standardized exams for consistent MRI results



SmartExam Spine Standardized exams for consistent MRI results



SmartExam Knee Standardized exams for consistent MRI results



Smart Workflow





VitalScreen Guidance at your fingertips

SmartExam Brain Standardized exams forconsistent MRI result



VitalScreen Guidance at your fingertips



SmartExam Brain Standardized exams forconsistent MRI result



SmartExam Spine Standardized exams for consistent MRI results





VitalScreen Guidance at your fingertips



VitalEye Touchless respiratorytriggering



SmartExam Brain Standardized exams forconsistent MRI results



SmartExam Spine Standardized exams for consistent MRI results



SmartExam Shoulder Standardized exams for consistent MRI results



SmartExam Knee Standardized exams for consistent MRI results





SmartSpeed Plus



SmartSpeed Neuro/Spine



SmartSpeed MSK



SmartSpeed MotionFree



SmartSpeed Implant



SmartSpeed DWI

SmartSpeed Pro Body



SmartSpeed Neuro/Spine



SmartSpeed MSK



SmartSpeed Body



SmartSpeed MotionFree



SmartSpeed MotionFree Body



SmartSpeed 3D FreeBreathing



SmartSpeed Implant



SmartSpeed DWI



SmartSpeed DWI Body

SmartSpeed Pro Cardiac



SmartSpeed Neuro/Spine



SmartSpeed MSK



SmartSpeed Cardiac



SmartSpeed MotionFree



SmartSpeed Implant



SmartSpeed DWI

SmartSpeed Premium



SmartSpeed Neuro/Spine



SmartSpeed MSK



SmartSpeed Body



SmartSpeed Cardiac



SmartSpeed MotionFree



SmartSpeed MotionFree Body



SmartSpeed 3D FreeBreathing



SmartSpeed Implant



SmartSpeed DWI



SmartSpeed DWI Body



Plus



MultiVane XD Motion-free imaging in short scan time



mDIXON XD TSE Replace all your FatSat by one single fat-free imaging solutior





mDIXON XD FFE Improve your fat-free imaging performance



3D VANE XD Free breathing abdominal imaging





MultiVane XD Motion-free imaging in short scan time



3D VANE XD Free breathing abdomina imaging



mDIXON XD TSE Replace all your FatSat by one single fat-free imaging solution



O-MAR XD Efficient near-metal soft tissue and bone imaging



mDIXON XD FFE Improve your fat-free imaging performance



dS Vascular Suite Pro



4D-TRAK XD Flexibility in your MR Angiography studies



mDIXON XD MultiStation Non-subtraction peripheral MR Angiography

dS Diffusion Suite Pro



Diffusion Excellence Pack Step up your diffusion performance



DWI XD TSE Speed up and improve the quality of your diffusion TSE



Zoom Diffusion Small FOV diffusion imaging for improved image quality

AV Onco Diffusion Pro



Advanced Diffusion Analyses Calculate and generate diffusion maps

॑



3D SpineVIEW View your 3D TSE imaging data in any plane



IRIS Zoom Improved small FOV spine diffusion imaging

T

dS Neuro Suite



Pro



SWIp Exquisite susceptibility contrast



Bold Motion-free imaging in short scan time





Spectroscopy XD More precise and more robust MR brain spectroscopy



Black Blood imaging Enhance your diagnostic confidence for brain imaging



DTI FiberTrak Fast, easy clinical fiber tracking



MultiBand SENSE High acceleration for your fMRI and DTI sequences



FiberTrak Extension High definition fiber tracking in the brain or spine



NeuroScience Explore brain connectivit



3D ASL Reproducible contrast-free brain perfusion



NeuroScience extension Extend your diffusion MRI studies



Plus

Π



MR Neuro Perfusion MR parametric brain analysis





MR FiberTrak Visualization of white matter tracts



MR iViewBold Real-time functional activation maps

dS Spectro Suite Pro



Spectroscopy Fully integrated proton spectroscopy



Spectroscopy XD More precise and more robust spectroscopy

dS Cardiac Suite



П



Cardiac Expert Expand your cardiac MR functionality

Pro



Cardiac Expert extension Fast CMR methods for anatomy, function and more



CardiacQuant Non-invasive assessment of myocardial tissue



k-t BLAST Accelerate your cardiac studies



Coronary Acquisition Perform non-invasive imaging of coronary arteries





3D BreastVIEW View your 3D TSE imaging data in any plane



SmartExam Breast Consistent fat suppression for every patient



4D-THRIVE/BLISS Accelerate dynamic breast imaging

dS Liver Suite Premium



mDIXON Quant Non-invasive liver fat fraction quantification



MR Elastography Non-invasive assessment of liver tissue stiffness



MR Elastography extension More robust liver tissue stiffness map

dS Pediatric Suite



dS Ped NeuroSpine coil High SNR for your pediatric brain and spine studies



dS Ped Torso coil High SNR for your pediatric torso and cardiac studies



Pediatric positioning pack Taking care of the smallest patients





3D PelvisVIEW View your 3D TSE imaging data in any plane



4D FreeBreathing Multi-phase liver studies without breathholds



MEGA Detection of additional metabolites



3D NerveVIEW Review nerve plexus, non-invasively



4D-TRANCE Contrast-free imaging of brain vascular anatomy



MR SpectroView Review metabolite content based on 1H spectroscopy



3D APT Enhanced diagnostic confidence in neuro oncology



CardiacQuant extension Flexibility in creation of T1 maps



MR Permeability Lesion characterization by reviewing vascular leakage



SyntAc Exploring neuro-radiology with Synthetic MR imaging



SENC Identifying asymptomatic patients at risk of heart failure



AV Cartilage Assessment Visualize cartilage structures

Index clinical applications

2K IMAGING	Page
3D APT	Page
3D ASL	Page
3D BrainVIEW	Page
3D BreastVIEW	Page
3D MSK VIEW	Page
3D NerveVIEW	Page
3D Non selective	Page
3D PelvisVIEW	Page
3D SpineVIEW	Page
3D VANE XD	Page
4D Free Breathing	Page
4D-THRIVE/BLISS	Page
4D-TRAK XD	Page
4D-TRANCE	Page
AutoVoice	Page
AV Cartilage Assessment	Page
bFFE XD	Page
Black Blood imaging	Page
Bold	Page
Cardiac Expert	Page
Cardiac Expert extension	Page
CardiacQuant	Page
CardiacQuant extension	Page
ComforTone	Page
Coronary Acquisition	Page
Diffusion Excellence Pack	Page
dS Ped NeuroSpine coil	Page
dS Ped Torso coil	Page
DTI FiberTrak	Page
DWI XD TSE	Page
FiberTrak Extension	Page
IRIS Zoom	Page
k-t BLAST	Page
mDIXON XD FFE	Page
mDIXON XD MultiStation	Page
mDIXON XD TSE	Page
mDIXON Quant	Page
MEGA	Page
MR Advanced Diffusion Analysis	Page
MR Diffusion	Page
MR Echo Accumulation	Page_
MR Elastography	Page_
MR Elastography Extension	Page

MR FiberTrak	Page	/2
MR iViewBold	Page	73
MR MobiView	Page	74
MR Neuro Perfusion	Page	75
MR Permeability	Page	76
MR SpectroView	Page	77
MR Subtraction	Page	78
MR T1 Perfusion	Page	79
MR Workspace	Page	80
Multiband SENSE	Page	81
MultiVane XD	Page	82
NeuroScience	Page	83
NeuroScience extension	Page	84
O-MAR	Page	85
O-MAR XD	Page	86
Pediatric positioning pack	Page	87
ScanTools Pro	Page	88
ScanWise Implant	Page	89
SENC	Page	90
SmartExam Brain	Page	91
SmartExam Breast	Page	92
SmartExam Knee	Page	93
SmartExam Shoulder	Page	94
SmartExam Spine	Page	95
SmartSpeed 3D FreeBreathing	Page	96
SmartSpeed Body	Page	97
SmartSpeed Cardiac	Page	98
SmartSpeed DWI	Page	99
SmartSpeed DWI Body	Page	100
SmartSpeed Implant	Page	101
SmartSpeed MotionFree	Page	102
SmartSpeed MotionFree Body	Page	103
SmartSpeed MSK	Page	104
SmartSpeed Neuro/Spine	Page	105
Spectroscopy	Page	106
Spectroscopy XD	Page	107
SWIp	Page	108
SyntAc	Page	109
VitalScreen	Page	110
VitalEye	Page	111
Whole Body	Page	112
ZOOM Diffusion	Page	113



2k Imaging

High resolution imaging



2k Imaging offers a scan matrix of 2048 x 2048, providing high resolution even with large FOVs, or lower resolution scans with a 2048 matrix reconstruction. Compatible with all imaging methods.



Ultra high resolution imaging



3D APT Enhanced diagnostic confidence in neuro oncology



3D APT (Amide Proton Transfer) is a new unique, contrast-free, brain MR imaging method addressing the need for more confident diagnosis in neuro oncology. 3D APT uses the presence of endogenous cellular proteins, to produce an MR signal that directly correlates with cell proliferation, a marker of tumoral activity. 3D APT can support trained medical professionals in differentiating low grade from high grade gliomas and, in differentiating tumor progression from treatment effect^{1,2}



3D APT image

Additional information:

- 3D APTw images are calculated automatically and displayed as color maps
- Whole glioma coverage can be obtained with a resolution of 2.0 x 2.0 x 5.0 mm

1 Togao et al. (2014) Neuro-Oncology 2 Park KJ et al. (2016) Eur Radiol



3D ASL

Reproducible contrast-free brain perfusion



3D ASL enables you to consistently quantify brain perfusion with an accuracy of 15%¹ in a non-contrast manner with full brain coverage, and better background suppression, compared to 2D pCASL method. 3D ASL includes fully automated calculation of color coded ASL maps.



Quantification of brain perfusion in a non-contrast manner



3D BrainVIEW

View your 3D TSE imaging data in any plane



3D BrainVIEW is an advanced 3D TSE technique that lets you acquire high resolution data in multiple directions, including oblique, in one scan helping you enhance your confidence when diagnosing lesions.



Additional information:

- Isotropic voxel size enabling reformats in any plane without loss of resolution.
- Allows for up to 20% shorter scan times¹.
- Available for a range of contrasts.



3D BreastVIEW



View your 3D TSE imaging data in any plane

3D BreastVIEW is an advanced 3D TSE technique that lets you acquire high resolution data in multiple directions, including oblique, in one scan helping you enhance your confidence when diagnosing lesions.





Data in multiple directions, in one scan

Additional information:

- Isotropic voxel size enabling reformats in any plane without loss of resolution.
- Allows for up to 20% shorter scan times¹.
- Available for a range of contrasts.



3D MSK VIEW

View your 3D TSE imaging data in any plane



3D MSK VIEW is an advanced 3D TSE technique that lets you acquire high resolution data in multiple directions, including oblique, in one scan helping you enhance your confidence when diagnosing lesions.



Additional information:

- Isotropic voxel size enabling reformats in any plane without loss of resolution.
- Allows for up to 20% shorter scan times¹.
- Available for a range of contrasts.



3D NerveVIEW

Review nerve plexus, non-invasively



3D NerveVIEW improves visualization of the brachial and lumbar plexus by providing you with a high resolution T2w TSE acquisition with reduced remaining intra-lumen signal of the veins¹. In addition, the 3D isotropic imaging method allows for reformats in any plane (including oblique) without loss of resolution helping you to save scan time and improve spinal nerve plexus assessment.



Improved visualization of the spinal nerve plexus



3D Non-selective Fast and robust large volume 3D FFE imaging

3D Non-selective enables faster and more robust¹ large volume 3D FFE imaging in brain and cardiac applications. Thanks to shorter TR and TE, 3D Non-selective delivers a 9% faster protocol and improved grey-white matter contrast in Brain 3D TFE.¹ For cardiac applications, 3D Non-selective delivers 3D bFFE with reduced banding artifacts.¹









3D PelvisVIEW

View your 3D TSE imaging data in any plane



3D PelvisVIEW is an advanced 3D TSE technique that lets you acquire high resolution data in multiple directions, including oblique, in one scan helping you enhance your confidence when diagnosing lesions.







Data in multiple directions, in one scan

Additional information:

- Isotropic voxel size enabling reformats in any plane without loss of resolution.
- Allows for up to 20% shorter scan times¹.
- Available for a range of contrasts.


3D SpineVIEW

View your 3D TSE imaging data in any plane



3D SpineVIEW is an advanced 3D TSE technique that lets you acquire high resolution data in multiple directions, including oblique, in one scan helping you enhance your confidence when diagnosing lesions.



Viewing imaging data in oblique directions

Additional information:

- Isotropic voxel size enabling reformats in any plane without loss of resolution.
- Allows for up to 20% shorter scan times¹.
- Available for a range of contrasts.

1 Due to time-efficient, low SAR flip angle sweep technology. Compared to standard 3D TSE.



3D VANE XD

Free breathing abdominal imaging

3D VANE XD supports imaging of the abdomen without the need for the patient to hold their breath, helping you reduce motion artifacts during free breathing¹ and improve patient comfort. With 3D VANE XD, you can now accommodate patients who are unable to hold their breath, including pediatric patients.





Breathhold mDIXON XD (left) versus a free breathing 3D VANE XD (right)



- 3D T1w FFE imaging method.
- Can be combined with fat suppression methods (eTHRIVE, mDIXON XD).



4D FreeBreathing

Multi-phase liver studies without breathholds



Multi-phase contrast-enhanced MRI studies are used to diagnose liver lesions, but many patients are simply not capable of performing multiple breathholds. With 4D FreeBreathing, you can now offer free-breathing MRI liver to a broader population, while improving imaging confidence and the patient experience. This allows you to address patients who have difficulty holding their breath or find it difficult to follow breathing instructions, like the rising elderly segment that has hearing loss, cognitive impairment or respiratory difficulties, as well as children and sedated patients.

4D FreeBreathing allows you to obtain excellent image quality from multi-phase liver studies performed without breathholds. This application is motion robust through its built-in respiratory soft gating and compatibility with high precision external sensors, like VitalEye. As a result, 4D FreeBreathing delivers reliable results that can improve imaging confidence¹. You can easily define variable timings for multiple phases to seamlessly fit 4D FreeBreathing into your current workflow. Real-time reconstruction allows you to view the progress of images as they are acquired to monitor the quality of the results. To provide easy workflow, you can specify the number of temporal phases to be reconstructed. 4D FreeBreathing provides dynamic information with a temporal resolution² down to 3 seconds per phase.



Multi-phase contrast-enhanced MRI liver study



4D-THRIVE/BLISS

Accelerate dynamic breast imaging



4D-THRIVE/BLISS is a time-resolved 3D technique to drastically accelerate dynamic breast imaging through the combination of a keyhole method with CENTRA and SENSE. Combines high spatial resolution with high temporal resolution to facilitate acquisition of multiple dynamic volumetric data sets per breath-hold.



Accelerated, high resolution, sagittal 3D breast imaging



4D-TRAK XD Flexibility in your MR Angiography studies

4D-TRAK XD provides a fast, dynamic contrastenhanced MR Angiography method with flexible sampling of both the arterial- and venous phase, by applying view sharing technique, enabling high spatial and temporal resolution simultaneously.





Fast, dynamic contrast-enhanced MR Angiography



4D-TRANCE

Contrast-free imaging of brain vascular anatomy



4D-TRANCE is a time-resolved technique for noncontrast angiography, promoting patient comfort and enabling you to evaluate the patency of the vascular anatomy in the brain using endogenous contrast with MIP visualization of multiple phases. 4D-TRANCE enables high temporal resolution down to 160 msec.



Non-contrast time-resolved angiography of the brain



AutoVoice

Guiding your patients through the examination

AutoVoice is a fully integrated and automated solution that guides your patients through the MR examination by indicating scan duration, announcing table movements and providing breathhold guidance, helping you enhance patient comfort.

Automated breathhold commands are included, based on your preferred setting:

- Pre-recorded breathhold instructions with timing synchronized to fit patient's respiratory cycle.
- Pre-recorded breathhold instructions with manual timing.
- Manual breathhold instructions with manual timing.



- Announcement:
 - Indication of scan duration.
 - · Announcement of table movements.
- \cdot Customization
 - Record your own text, with local pronunciation.
 - Record your own voice.
- Export
 - Allows export of customized settings to other Philips MR scanners.



AV Cartilage Assessment

Visualize cartilage structures

entres de la conserva de la conserva

AV Cartilage Assessment enables the advanced visualization of cartilage structures and provides tools to support determination of the degradation of the cartilage. Quantitative analysis of T2 relaxation time displayed via integrated color-coded T2 maps. Positioning of cartilage-shaped, layered region of interest is used to assess variation of T2 values across the cartilage depth to determine the degradation of the cartilage.



MR Echo accumulation



bFFE XD Expanding clinical applications of balanced FFE

bFFE XD expands the clinical application of bFFE towards better visualization of fine structures. It delivers robust 3D, high resolution imaging with a spatial resolution up to $0.5 \times 0.5 \times 0.5 \text{ mm}^1$ in less than 6 minutes for inner ear applications, with reduced banding artifacts compared to conventional Philips balanced FFE.







Black Blood imaging Enhance your diagnostic confidence for Brain imaging



Black Blood imaging helps you better differentiate the vessel lumen from the intra lumen blood signal. This enhances your diagnostic confidence by performing your 3D brain imaging with higher and isotropic imaging resolution¹ with a reduction of the intra-lumen brain blood signal² over the complete imaging volume.



Reduction of the intra-lumen brain blood signal

- Fast scan times³ of five minutes.
- 3D isotropic acquisition enables reformats in any plane (including oblique) without loss of resolution.



BOLD

Fast, easy and reliable fMRI

Accurately acquiring fMRI BOLD data during neuro imaging helps visualize task-related areas of activation in the brain. The fMRI paradigms that deliver and control stimuli are fully automated via dedicated ExamCards to make fMRI fast, easy, and reliable. The iView BOLD analysis package provides real-time processing of fMRI BOLD data into functional activation maps.





Visualize task-related areas of activation in the brain



Cardiac Expert

Expand your cardiac MR functionality

Cardiac Expert supports the acquisition of multislice, dynamic tissue studies with T1 weighting and uniform tissue suppression¹ by including Look Locker methods for determining an optimal inversion delay time. Cardiac Expert also provides myocardial tagging² to allow assessment of regional wall motion and allows for real-time interactive planning of challenging cardiac views.





Dynamics



Cardiac Expert extension Fast CMR methods for anatomy, function and more



Cardiac Expert extension is an add-on to the comprehensive Cardiac Expert option. It provides additional techniques for fast black blood imaging, functional imaging and dynamic cardiac MR studies. Cardiac Zoom is a small FOV imaging technique that accelerates black blood TSE of the heart and great vessels. It decreases the required breath hold duration by up to 30% without changing spatial resolution by enabling single beat (shot)imaging, which is challenging for conventional (multi-beat) imaging approaches¹. 3D Non-selective delivers 3D bFFE with reduced banding artifacts compared to Philips 3D Selective 3D FFE imaging. Retrospective EPI combines retrospective triggering with EPI sampling. kt-SENSE is a spatio-temporal acceleration technique that offers all the benefits of k-t BLAST in addition to enhanced image uniformity².



3D Non-selective bFFE in cardiac applications

1 Compared to conventional Philips black blood imaging 2 Compared to regular k-t BLAST



CardiacQuant Non-invasive assessment of myocardial tissue



With CardiacQuant you get access to exciting new applications for cardiology, which can help in the non-invasive assessment of myocardial tissue characteristics by providing you with comprehensive graphs and pixel-based, quantitative information in different regions of the myocardium helping you to make early decisions for therapy.



Quantitative T2*, T2 and T1 maps in a single breathhold scan



★ Ξ

CardiacQuant Extension Flexibility in creation of T1 maps

CardiacQuant Extension is an optional plugin which adds flexibility for the creation of T1 maps. It allows the option of user defined T1 mapping schemes as alternatives for the predefined "native" or "enhanced" schemes as provided by CardiacQuant.





T1 map



ComforTone Reduce acoustic noise and enhance patient experience

Leveraging our years of experience in acoustic noise reduction technologies, our unique ComforTone solution achieves up to 80% acoustic noise reduction¹ with similar image quality and contrast within the same time slot.



T2w TSE

T2w TSE - ComforTone 80% reduction in acoustic noise



- Can be used in routine exams, like brain, spine and MSK.
- · Can be applied with high gradient settings.
- Easy to implement, thanks to ready-to-use ExamCards.



Coronary Acquisition Perform non-invasive imaging of coronary arteries



Coronary Acquisition allows for non-invasive imaging of coronary arteries by displaying good contrast between myocardium and vessels by deploying 3D sequences combined with MotionTrak respiratory navigators for real-time motion correction and T2-preparation.





Non-invasive imaging of coronary arteries



Diffusion Excellence Pack

Step up your diffusion performance



The Diffusion Excellence Pack offers a unique range of innovations to address the common challenges you face to enhance the efficiency and guality of diffusion imaging for areas ranging from oncology to neuro.

SmartShim, delivering reliable fat suppression with simplified workflow due to automating the planning of the shim region and image-based shimming.

Computed DWI, decreasing overall exam time and enhancing clinical workflow¹ by generating synthetic high b-value images.

LOVA ADC, delivering consistent ADC values with up to 95% accuracy² in large field of views.

EPIC Brain, improving the geometrical match between EPI images and MR anatomical brain images, compared to conventional EPI scans.

1 Compared to Philips DWI, by generating synthetic high b-value images. 2 Based on in house testing on phantom at a FOV of 45 cm.





Reliable fat suppression

Synthetic b-value images





Consistent ADC values

Improved geometrical match



dS Ped NeuroSpine coil High SNR for your pediatric brain and spine studies



The dS Ped NeuroSpine 8ch coil is an open-design 8-element coil for high resolution pediatric brain and spine imaging. Specifically designed for neonates, but will accommodate pediatric patients up to 10kg. Open, cradle shaped design enables the operator to position and prepare the patient outside the examination room. Examinations of brain and spine can be performed without having to move the patient.





T2w and T1w Spine imaging

3D T1w Brain imaging



dS Ped Torso coil High SNR for your pediatric torso and cardiac studies

The dS Ped Torso coil is a dedicated 8-element coil designed to provide advanced pediatric torso and cardiac imaging. The coil is optimized for neonates, but will accommodate pediatric patients up to 10kg. Split design allowing the top of the coil to be taken off, enabling easy access to the patient. An insert cradle can be used for additional patient support. A surrounding mattress is available to accommodate larger patients.







T2w abdominal imaging

Cardiac cine imaging



DTI FiberTrak

Fast, easy clinical fiber tracking

Visualize specific white matter fiber tracts in the brain with Diffusion Tensor Imaging (DTI) data and fiber tracking. This package allows you to trace, analyze and process fibers in real-time with minimal mouse clicks. It supports pre-operative surgical planning, post-surgery evaluation, and general evaluation of fiber tracts around tumors and lesions in connection with functional areas. DTI FiberTrak supports up to 32 directions and 16 b-values and includes automatic calculation of Fractional Anisotropy (FA) maps.





Visualization of white matter fiber tracts in the brain



DWI XD TSE Speed up and improve the quality of your diffusion TSE



DWI XD TSE delivers up to 25% faster diffusion TSE imaging with improved resolution due to its multishot approach¹. DWI XD TSE is compatible with MultiVane, contributing to robust suppression of motion artifacts². It also delivers images with less distortion because it is less sensitive to susceptibility differences compared to Philips conventional DWI EPI sequences.



DWI EPI (left) versus robust inner ear DWI XD TSE (right)

1 Compared to Philips DWI TSE (Single-shot) 2 Compared to Philips multi shot DWI TSE



FiberTrak Extension High definition fibertracking in the brain or spine

The FiberTrak Specialist Extension package Allows for diffusion imaging with up to 128 b-vectors and 16 b-values, delivering input for very high definition fiber tracking in the brain or spine.





Visualization of white matter fiber tracts in the brain



IRIS Zoom

Improved small FOV spine diffusion imaging

IRIS Zoom delivers small FOV diffusion imaging with higher resolution, lower distortion and improved fat suppression than Philips Zoom DWI. Higher resolution diffusion imaging in the spine is achieved by employing 2D navigatorbased motion correction integrated into the dS-SENSE framework. IRIS Zoom also delivers higher SNR in spine imaging compared to Philips MultiVane DWI XD TSE.



b800 Zoom DWI (left) versus IRIS Zoom DWI (right)



k-t BLAST

Accelerate your cardiac studies

K-t BLAST provides up to five fold acceleration using an alternative parallel imaging technique employing undersampling in time and space. Suited for dynamic and real-time cardiac studies as well as single breath hold, multi-slice cine studies. Can be combined with most other imaging methods.







mDIXON XD FFE

Improve your fat-free imaging performance



mDIXON XD FFE provides more efficient fat-free imaging in routine scan times. Improve your fat-free imaging over large field-of-views and for high resolution imaging. With up to four image types in one single scan, including with or without fat suppression contrasts, mDIXON XD FFE will enable you to enhance your imaging strategies by simplifying your routine FFE procedures.





Additional information:

- Improved fat-free imaging over large 400-500 mm FOV and for sub-millimetric resolution¹.
- More efficient, faster scanning².
- Increased signal-to-noise ratio².
- Acquire up to four image types in one single scan (water only, in phase, out phase, fat only)

 Compared to the standard mDIXON algorithm, due to unique 7-peak fat model and improved B0 correction.
 Due to the unrestricted echo-time (TE) approach in mDIXON allowing more freedom in protocol optimization.



mDIXON XD MultiStation Non-subtraction peripheral MR Angiography

mDIXON XD MultiStation allows you to perform peripheral MR Angiography with improved vessel-to-background contrast in only one single pass1. You will be able to perform your peripheral MR Angiography acquisitions without the use of a subtraction mask, eliminating artifacts that could arise from misalignment, due to patient motion, between the pre and post contrast scan. Enjoy fast, robust peripheral MR Angiography.



MR Angiography with subtraction (left) and in one single pass (right) with improved vessel-to-background contrast



- Subtraction-less peripheral MR Angiography
- Improved vessel-tobackground contrast by 30-36%¹



mDIXON XD TSE Replace all your FatSat by one single fat-free imaging solution



mDIXON XD TSE brings a new dimension to fat suppression by providing uniform, complete and consistent fat-free imaging, even over large field-of-views and in challenging anatomies. Providing up to four image types in one single scan, including with/without fat suppression contrasts, in routine scan times and resolution simultaneously, you can easily replace your favorite routine TSF scans with it. mDIXON XD TSE will enable you to enhance your imaging strategies by simplifying your routine TSE procedures.



With/without fat suppression contrasts, simultaneously

Additional information:

- 30% faster scanning and up to 30% reduced blurring¹.
- · Increased signal-to-noise ratio².
- Acquire up to four image types in one single scan (water only, in phase, out phase, fat only).

 Due to its unique 2-echo technology, compared to the conventional 3-echo DIXON TSE techniques.
 Compared to a standard non-fat-shift corrected fat-free TSE approach.



mDIXON Quant

Non-invasive liver fat fraction quantification



mDIXON Quant brings a fast and simple 3D procedure for noninvasive liver fat quantification by providing high quality 3D fat fraction maps of the whole liver, even for short T2*, with high accuracy (± 3.5%) and reproducibility (± 1.4%)¹ allowing you to expand your MRI capabilities. T2*/R2* relaxation maps are provided to further help your diagnostic assessment.



Fat fraction maps (left) and T2*/R2* relaxation maps (right)

Additional information:

- Single breathhold acquisition.
- Based on state of the art 6-echo acquisition,
 7-peak fat modeling reconstruction, correction for T2* confounding effect and low flip angle to minimize T1 bias.
- Fat fraction maps are displayed in colors with a quantification bar.

Accuracy and reproducibility were assessed using a reference liver protocol, on fat phantoms [range: 0-100%].
 Reproducibility assessed over systems.



MEGA Detection of additional metabolites



MEGA improves spectroscopy by revealing spectral peaks of interest which would otherwise remain hidden. It also allows detection and relative quantification of J-coupled metabolites such as gamma-aminobutyric acid (GABA) by automatically removing the spectral overlap from other metabolites. (In conventional spectroscopy, removing spectral overlap is only possible with spectral editing.) Frequency-selective RF pulses are included to manipulate the evolution of J-coupled MR signals. In addition, subtraction of on- and off-resonance spectra is used for relative quantification of J-coupled metabolites.



Detection of GABA with single voxel MEGA MR Spectroscopy



MR Advanced Diffusion Analysis Calculate diffusion maps at a b-value of choice



Advanced Diffusion Analyses (ADA) is intended to view, process and analyze MRI diffusion weighted images. It calculates and displays cDWI at a b-value of choice (from 0 to 5,000 s/mm2). The package also supports advanced MR parametric analysis and visualization tools, generating maps of diffusion kurtosis and intravoxel incoherent motion (IVIM) diffusion.



Calculate and generate diffusion maps

- Presents a default diffusion analysis model based on the available original DWI images as well as a selection of alternative models including monoexponential, bi-exponential, simplified IVIM, and kurtosis.
- A 'goodness of fit' value and fitted curve show the fitting quality of the selected model.
- Provides parametric maps of perfusion fraction (f), pseudo diffusivity (D*), Diffusivity (D) and Kurtosis (K).



MR Diffusion Analyze diffusion and anisotropic properties of tissue



MR Diffusion is designed to analyze diffusion and anisotropic properties of tissue. The application evaluates Diffusion Weighted Imaging (DWI) series to generate parametric maps such as ADC and eADC. For Diffusion Tensor Imaging data, additional parametric maps are generated, including fractional anisotropy, axial diffusivity or radial diffusivity. MR Diffusion can help to characterize diffusion properties of tissue in such cases as tumor assessment or stroke.



Additional information:

 The user can make a sub-selection of the acquired b-values for analysis and select preferred color-coding for the parametric maps. **★** Ξ

MR Echo Accumulation Calculation of new images based on series with multiple echoes



MR Echo Accumulation is used to perform pixelwise echo accumulations for imaging series with multiple echoes. MR Echo Accumulation enables the calculation of new images based on the selected sum of echo times of series with multiple echoes. The processing provides interactive update of the results



- MR Echo Accumulation enables the calculation of new images based on the selected sum of echo times.
- Provides the ability to define the accumulation range using an interactive slider.
- Enables preview, save and analysis on the calculated new series.



MR Elastography

Non-invasive assessment of liver tissue stiffness



MR Elastography allows for a non-invasive assessment of differences in tissue stiffness of the liver in a fast breathhold scan providing trained physicians with additional input to help make informed decisions about treatment.





Elastograms reflecting tissue stiffness in kPa

- Image processing is fully integrated at the scanner.
- Automated calculation of Elastograms, reflecting tissue stiffness in kPa.
- Statistical confidence map is provided for reliability assessment.



MR Elastography Extension

More robust liver tissue stiffness maps



The MR Elastography Extension has been developed to enhance diagnostic confidence by delivering faster, more robust stiffness maps across larger regions in the liver. With this new extension. MRE (SE-EPI) can be performed up to 8 times faster¹ and in a single breathhold, while delivering equal or better image confidence. The MR Elastography Extension supports you in getting consistent results for diverse patients. At high field strength the SE based MR Elastography Extension provides more robust stiffness maps, since the technique is more resistant to low T2* signal fall out, compared to FFE based MR Elastography. A noninvasive assessment of tissue stiffness of the liver can be SE based MR Elastography Extension allows for, with a slice coverage of 4-5 slices, in a single breathhold scan. While FFE based MR Elastography requires one breathhold per slice. Image processing is fully integrated at the scanner with automated calculation of elastograms, providing trained physicians with additional input to help make informed decisions about treatment.





MRE, 4 slices in a single breath hold of 16 seconds

MRE, 8 slices in a single breath hold of 16 seconds



MR FiberTrak

DTI tractography reconstruction

Provides visualization and quantification of white matter structure in the brain and spinal tracts using task guidance for generating common or user-defined tracts. This supports advanced anatomical brain assessments via DTI tractography reconstruction.



Visualization of white matter tracts



- The guidance panel suggests which regions of interest and plane are common for identification of certain tracts such as the corticospinal tract.
- The results can be overlaid with other data like fMRI or anatomical series.
- Allows evaluation of fiber tracts around tumors and lesions in combination with functional areas.
- Supports DICOM-based output with merged anatomical tractography information through the Multi Modality Viewer.


MR iView Bold

Functional MRI analysis

MR iView Bold includes the functional MRI (fMRI) analysis package providing real-time processing of fMRI BOLD data into functional activation maps. This helps identify and visualize functional regions of the brain, relying on local metabolic and hemodynamic changes that occur in activated brain areas.



Real-time functional activation maps



- The tool applies a generalized linear regression model to analyze block paradigms, event-related paradigms, and resting state data.
- The paradigms can be user-defined or imported.
- Supports export of functional results including DICOMbased images with co-registered anatomical and fMRI maps.



MR MobiView

Easy viewing of multi-station studies

MR MobiView combines multiple images into a single full-field view to review multi-station acquisitions. This is easily accomplished by fusing images in the head-feet direction. MR MobiView automatically combines "stitches" images from multiple acquisitions of the same examination to create one overall volume.





- Key clinical cases are MRA run-offs, wholebody metastases screening from eye-tothigh, and total spine views to show the complete CNS.
- The resulting image series can be viewed, filmed, and exported using a DICOM compliant tool.



MR Neuro Perfusion

MR parametric brain analysis



MR Neuro Perfusion supports MR parametric brain analyses by evaluating DSC T2* perfusion studies and supporting the visualization and quantification of diffusion perfusion mismatch. MR Neuro Perfusion provides physicians with supporting information for the evaluation of stroke, or assessment and follow-up of brain tumors. The application supports the analysis of T2* Perfusion studies to generate parametric data including TTP, MTT or Tmax.



- Offers several analysis techniques such as leakage correction, which allows to assess the time intensity curves where there is no proper recovery of the baseline after contrast passage, and manual arterial input function (AIF) which enables perfusion-diffusion mismatch if a Diffusion input dataset is available in addition to the Perfusion series.
- The package includes user-selected color coding of the functional data, and maps can be viewed and stored as overlays on anatomical reference images.



MR Permeability Lesion characterization by reviewing vascular leakage



MR Permeability allows to measure the leakage of gadolinium chelates into the extra-vascular extracellular space (EES). Designed to visualize T1 weighted DCE 3D datasets and assist in analyzing the tissue response. The package supports the computation of parametric maps. Its use plays an important role in oncology.



Parametric map helping to analyze tissue response

- Calculates parametric maps such as Ktrans, Kep, ve and vf.
- The application has been validated for prostate and brain cancer.



MR SpectroView Review metabolite maps

MR SpectroView is a task-guided application providing hydrogen single voxel spectra as well as metabolic and ratio maps. It automatically identifies the anatomy to preselect appropriate metabolites or supports user-defined combination of metabolites.



Detection of GABA with single voxel MEGA MR Spectroscopy



- Displays numerical information about metabolites including Peak position and label, SNR, Peak Height, Peak Area, Full Width Half Maximum and Area Ratio of the displayed spectrum.
- Provides metabolite and ratio maps as color overlay on anatomical images or mini spectra on a voxel by voxel basis. Multiple voxels can be selected for spectral comparison.
- Supports automatic and manual phase adjustment as well as a color-coded quality indicator based on field homogeneity and SNR.



MR Subtraction Perform basic calculations between two volumes



MR Subtraction enables basic calculations between two volumes, including addition, subtraction and ratio from within a single dynamic series.



- The application allows to subtract pre-contrast from post-contrast series.
- Weighing factors can be applied to impact the calculation.



MR T1 Perfusion

Review blood supply characteristics



Designed to evaluate time intensity curves of a T1 signal enhancement series, supporting in assessing lesions by reviewing blood supply characteristics. The application produces measurements including relative enhancement, maximum enhancement, time to peak (TTP), and wash-in/ wash-out rates.



- The package includes user-selected color-coding of the parametric maps which can be viewed and stored as overlays on anatomical reference images.
- The opacity of the overlay is user-defined.
- The series can be referenced to any other series such as Diffusion data within the same study.

MR Workspace

MR Workspace supports efficiency and staff satisfaction in the control room through intelligence, guidance and ease of use. Technologists can prepare exams before patients arrive and aim to achieve consistent quality regardless of experience, by using Protocol Assistant, an Al¹-driven solution that learns your protocol preferences and suggests the most appropriate ones based on clinical indication. Advanced visualization includes step-by-step guidance so technologists can begin post-processing, delivering more² diagnostic information to the radiologist, thus saving reporting time. The intuitive interface, large display of clinical images and essential parameter reveal contribute to outstanding ease-of-use.

MR Workspace also supports a large collection of clinical application suites. It is the Philips platform of choice for existing and future clinical innovations.

1 According to the definition of AI from the EU High-Level Expert Group.

2 The addition of step-by-step guidance and automation of routine and complex post-processing applications can now be performed

by the technologist on the console, saved via bookmark functionality, and handed off to the radiologist, which reduces time to results.



MultiBand SENSE High speed, high resolution fMRI and DTI studies

By simultaneously exciting multiple slices, MultiBand SENSE allows you to use state-of-the-art acceleration techniques in functional brain scans to either shorten MRI scan time or enhance diagnostic information. It doubles your capabilities for fMRI. Due to a shorter minimum TR for fMRI, MultiBand SENSE can be used to deliver twice the temporal resolution with virtually no compromise in SNR, or to obtain twice as much anatomical coverage at similar scan times with virtually equal image quality.¹ In your DWI/DTI brain sequences, MultiBand SENSE enables up to 45% acceleration or twice as many diffusion vectors to be acquired with virtually equal image quality.¹ With MultiBand SENSE you can perform fMRI and DTI sequences at high speed and high resolution.²





DTI 16 directions 6:51 min



DTI with MultiBand SENSE 16 directions 3:37 min

1 Compared to Philips DTI/fMRI scans without MultiBand SENSE. 2 Up to a MultiBand SENSE factor of 2.



MultiVane XD

Motion-free imaging in short scan time

MultiVane XD delivers high resolution diagnostic images even in the case of severe patient motion by providing motion correction to a full range of anatomies, in short scan times¹. MultiVane XD works in multiple orientations and for various contrasts (T1w, T2w, FLAIR) helping you to increase your diagnostic confidence.





Diagnostic images, even in the case of severe patient motion



NeuroScience

Explore brain connectivity

NeuroScience comprehensive package helps you to explore brain connectivity by supporting advanced acquisition schemes allowing for high-definition brain fiber tracking, including crossing fibers and advanced fMRI capabilities.



Diffusion acquisition with a b-value of 15.000



- Allows diffusion-weighted multi-shell acquisitions with up to 32 b-values and up to 128 unique diffusion directions
- Easy workflow for user defined gradient direction input
- Perform your fMRI studies with enhanced nyquist ghost stability and extended data storage (up to 64k images)
- Enables monitoring of consistency in longitudinal fMRI studies with a quality assurance tool, in line with fBRIN standards
- Includes BO mapping for offline data correction and image processing
- Easy-to-use export tools in various formats, including NIfTI



NeuroScience extension

Extend your diffusion MRI studies

NeuroScience extension is an add-on to the comprehensive NeuroScience option. The extension brings your multi-shell DTI studies to a higher level. Advanced diffusion gradient control gives the scientific user control of the diffusion encoding gradient duration through selection of multiple diffusion encoding gradient waveforms. Furthermore, 2k DTI provides advanced control over diffusion gradients with up to 2048 independent diffusion encodings (vectors), each with up to 1024 different weightings and 1024 different directions





Multi-shell DTI b4000 (128 directions)



O-MAR

Efficient near-metal soft tissue and bone imaging



O-MAR (Metal Artifact Reduction for Orthopedic implants) allows you to improve visualization of more soft tissue and bone in the near vicinity of MR Conditional Orthopedic implants¹. This allows you to offer postoperative MR imaging to patients with implants who could develop implantrelated conditions.





- Reduction of in-plane susceptibility artifacts² caused by metal implants¹.
- Supports most relevant image contrasts (T1w, T2w, PDw, and STIR).
- Extending MARS (Metal Artifact Reduction Sequence) with the View Angle Tilting (VAT).



O-MAR XD

Efficient near-metal soft tissue and bone imaging



O-MAR XD (Metal Artifact Reduction for Orthopedic implants) allows you to improve visualization of more soft tissue and bone in the near vicinity of MR Conditional Orthopedic implants¹. This allows you to offer postoperative MR imaging to patients with implants who could develop implant-related conditions.





Additional information:

- Reduction of in- and throughplane susceptibility artifacts² caused by metal implants¹.
- Supports most relevant image contrasts (T1w, T2w, PDw, and STIR).
- Extending MARS (Metal Artifact Reduction Sequence) with the View Angle Tilting (VAT) and Slice Encoding for Metal Artifact Correction (SEMAC) techniques.

36 1 Only for use with MR Safe or MR Conditional implants by strictly following the Instructions for Use.

2 Compared to standard high bandwidth spin-echo based techniques.



Pediatric positioning pack

Taking care of the smallest patients



The Pediatric positioning pack for dStream systems consists of a subset of accessories designed to meet the needs of the smallest patients. Included is an anterior coil frame avoiding positioning of the dS Torso coil directly on the patient. Additionally a baby support pad, comfort pad, pediatric knee support and child elevation mattress are included in this package.





ScanTools Pro Clinical imaging performance,enhanced workflow, and ease of use



ScanTools Pro provides comprehensive planning, scanning, processing, viewing, and patient administration tools, as well as data storage and connectivity features. All ScanTools elements support outstanding clinical imaging performance, enhanced workflow, and ease of use. ScanTools Pro contains fast, high resolution imaging methods for the assessment of morphology of all anatomical areas including brain and spine, MSK, body and breast, cardiac, and various blood vessels with or without contrast agents.



★ Ξ

ScanWise Implant A key to confidence withMR Conditional implants

ScanWise Implant software helps you to confidently scan patients within the MR Conditional limits by providing step-by-step guidance to enter the condition values of the implant manufacturer. Your MR system then automatically applies these values for the entire examination helping you to simplify your scanning process for patients with MR Conditional implants. With ScanWise Implant you can increase referrals of patients with MR Conditional implants to your institution by confidently offering MR imaging to this growing patient population.



- Can be used in routine exams, like brain, spine and MSK.
- Step-by-step guidance of scanning parameters as to adhere to the implant manufacturer conditional values. Enter parameters for all scans just once.



SENC Identifying asymptomatic patients at risk of heart failure¹



SENC images are processed through Myocardial Solutions software to measure early and subtle changes in the heart function. The strain-encoded (Fast-SENC) time resolved images provided by Philips are used to extract the quantitative strain information per voxel using third party Myocardial Solutions (MyoStrain) software which generates a clinical report.

With the combination of Philips Fast-SENC and MyoStrain early dysfunction of heart failure can be detected across 48 segments of the heart1 in 10 minutes.



Measure early and subtle changes in the heart function



SmartExam Brain

Standardized exams for consistent MRI results



SmartExam¹ Brain assists in delivering reproducible planning results by using intelligent software which automatically plans the scanning geometries, based on your validated scanning preferences. This enables you to standardize your MRI exam process helping you to enhance consistency in follow-up exams of the same patient and from patient to patient.





Enhanced consistency in follow-up exams

- Dedicated 3D survey scan is included to determine patient positioning.
- Automated planning of the imaging stack is based on anatomic landmarks relating those to a previously defined planning.
- SmartExam planning can be adapted and expanded to fit changing requirements.
- Automated geometry planning can be shared and applied across Philips MRI consoles.



SmartExam Breast

Consistent fat suppression for every patient





SmartExam Breast¹ provides consistent fat suppression for every patient and assists in delivering reproducible planning results by using intelligent software which automatically plans the scanning geometries, based on your validated scanning preferences. This enables you to standardize your MRI exam process helping you to enhance consistency in follow-up exams of the same patient and from patient to patient.









Additional information:

- Dedicated 3D survey scan is included to determine patient positioning.
- Automated planning of the imaging stack is based on anatomic landmarks relating those to a previously defined planning.
- SmartExam planning can be adapted and expanded to fit changing requirements.
- Automated geometry planning can be shared and applied across Philips MRI consoles.

MR Conditional implants.

Consistent fat suppression for every patient



SmartExam Knee

Standardized exams for consistent MRI results



SmartExam Knee¹ assists in delivering reproducible planning results by using intelligent software which automatically plans the scanning geometries, based on your validated scanning preferences. This enables you to standardize your MRI exam process helping you to enhance consistency in follow-up exams of the same patient and from patient to patient.



Consistent reading for any patient



- Dedicated 3D survey scan is included to determine patient positioning.
- Automated planning of the imaging stack is based on anatomic landmarks relating those to a previously defined planning.
- SmartExam planning can be adapted and expanded to fit changing requirements.
- Automated geometry planning can be shared and applied across Philips MRI consoles.



SmartExam Shoulder

Standardized exams for consistent MRI results



SmartExam Shoulder¹ assists in delivering reproducible planning results by using intelligent software which automatically plans the scanning geometries, based on your validated scanning preferences. This enables you to standardize your MRI exam process helping you to enhance consistency in follow-up exams of the same patient and from patient to patient.



Consistent reading for any patient



- Dedicated 3D survey scan is included to determine patient positioning.
- Automated planning of the imaging stack is based on anatomic landmarks relating those to a previously defined planning.
- SmartExam planning can be adapted and expanded to fit changing requirements.
- Automated geometry planning can be shared and applied across Philips MRI consoles.



SmartExam Spine

Standardized exams for consistent MRI results



SmartExam Spine¹ assists in delivering reproducible planning results by using intelligent software which automatically plans the scanning geometries, based on your validated scanning preferences. This enables you to standardize your MRI exam process helping you to enhance consistency in follow-up exams of the same patient and from patient to patient.



Consistent reading for any patient

- Dedicated 3D survey scan is included to determine patient positioning.
- Automated planning of the imaging stack is based on anatomic landmarks relating those to a previously defined planning.
- SmartExam planning can be adapted and expanded to fit changing requirements.
- Includes numbering of the vertebrae and automatically matches the planning of the axial stacks to the disc's orientation.
- Automated geometry planning can be shared and applied across Philips MRI consoles.



SmartSpeed 3D FreeBreathing Image quality and speed without compromise



Philips SmartSpeed delivers fast high-quality imaging for wide range of patients including patients who struggling to hold their breath during the exam. Philips SmartSpeed 3DFreeBreathing helps to acquire ultra-fast, high-quality imaging with reduced artifacts¹. It allows the acquisition of 3D T1-weighted gradient echo scans without the need for breath holding. It is intrinsically robust for motion artifacts that can originate from breathing, peristaltic motion.



Fewer artifacts in breast imaging

- Up to 3 times faster with no loss in image quality²
- Up to 65% higher resolution and improved SNR²
- Improve your productivity with reduced exam time without losing image quality
- Improve signal to noise ratio for challenging anatomies



SmartSpeed Body Image quality and speed without compromise



Philips SmartSpeed delivers image quality and speed at your fingertips. It utilizes the Compressed SENSE speed engine to reduce scan time and an award-winning AI algorithm applied directly at the beginning of the MR reconstruction chain to maximize information and deliver outstanding image quality to body imaging. Philips SmartSpeed can be used in 2D and 3D and for all anatomical contracts. It supports 97% of current clinical MR protocols1 to address the imaging needs of the vast majority of the patients.



High quality 3D MRCP imaging

Additional information:

- Up to 3 times faster with no loss in image quality²
- Up to 65% higher resolution and improved SNR²
- Improve your productivity with reduced exam time without losing image quality
- Improve signal to noise ratio for challenging anatomies

1 On average, measured across a sample of sites from Philips MR installed 2 Compared to Philips SENSE



SmartSpeed Cardiac Image quality and speed without compromise





Fast cardiac cine imaging



Additional information:

- Up to 3 times faster with no loss in image quality²
- Up to 65% higher resolution and improved SNR²
- Improve your productivity with reduced exam time without losing image quality
- Improve signal to noise ratio for challenging anatomies

1 On average, measured across a sample of sites from Philips MR installed 2 Compared to Philips SENSE



SmartSpeed DWI Image quality and speed without compromise

Philips SmartSpeed delivers image quality and speed at your fingertips. By nature, DWI scans have a low signal-to-noise ratio (SNR) as the signal has decayed by the applied diffusion encoding gradients. To have an adequate image quality multiple diffusion directions are involved which can take a considerable amount of time. With Philips SmartSpeed Diffusion the technology of DWI is integrated in the Philips SmartSpeed Engine to reduce the scan time and improve the SNR of individual diffusion-weighted imaging measurements¹.



High quality diffusion TSE imaging



- Up to 3 times faster with no loss in image quality²
- Up to 65% higher resolution and improved SNR²
- Improve your productivity with reduced exam time without losing image quality
- Improve signal to noise ratio for challenging anatomies



SmartSpeed DWI Body Image quality and speed without compromise



Philips SmartSpeed delivers image quality and speed at your fingertips. By nature, DWI scans have a low signal-to-noise ratio (SNR) as the signal has decayed by the applied diffusion encoding gradients. To have an adequate image quality multiple b-values. directions and averages are involved which can take a considerable amount of time. With Philips SmartSpeed Diffusion the technology of DWI is integrated in the Philips SmartSpeed Engine to reduce the scan time and improve the SNR of individual diffusionweighted imaging measurements¹.



High quality DWI breast imaging

- Up to 3 times faster with no loss in image quality²
- Up to 65% higher resolution and improved SNR²
- Improve your productivity with reduced exam time without losing image quality
- Improve signal to noise ratio for challenging anatomies



SmartSpeed Implant Image quality and speed without compromise



Philips SmartSpeed delivers image guality and speed at your fingertips including patients with implants. Performing musculoskeletal (MSK) imaging in patients with implants is challenging. With the introduction of technologies such as Orthopedic Metal Artifact Reduction (O-MAR XD), image guality has substantially improved. However, scan times are longer. With Philips SmartSpeed Implant, the technology of O-MAR XD is integrated in the Philips SmartSpeed Engine to reduce the scan time of the non-Cartesian sequences significantly¹.



Fast imaging in the presence of Conditional implants

- Up to 3 times faster with no loss in image quality²
- Up to 65% higher resolution and improved SNR²
- Improve your productivity with reduced exam time without losing image quality
- Improve signal to noise ratio for challenging anatomies



SmartSpeed MotionFree Image quality and speed without compromise

Philips SmartSpeed delivers fast high-quality imaging for wider range of patients including patients who are in pain, struggling to hold still. Philips SmartSpeed MotionFree utilizes non-Cartesian, motion robust 2D imaging techniques to acquire fast, motion-free images. It reduces gross motion, breathing motion and pulsatility artifacts in over 90% of the cases compared to Cartesian imaging.



Motion free neck imaging



- Up to 3 times faster with no loss in image quality¹
- Up to 65% higher resolution and improved SNR¹
- Improve your productivity with reduced exam time without losing image quality
- Improve signal to noise ratio for challenging anatomies



SmartSpeed MotionFree Body Image quality and speed without compromise



Philips SmartSpeed delivers fast high-quality imaging for wider range of patients including patients who are in pain, struggling to hold still. Philips SmartSpeed MotionFree utilizes non-Cartesian, motion robust 2D imaging techniques to acquire fast, motion-free body images. It reduces gross motion, breathing motion and pulsatility artifacts in over 90% of the cases compared to Cartesian imaging.



High quality pelvis imaging

- Up to 3 times faster with no loss in image quality²
- Up to 65% higher resolution and improved SNR²
- Improve your productivity with reduced exam time without losing image quality
- Improve signal to noise ratio for challenging anatomies



SmartSpeed MSK Image quality and speed without compromise

Philips SmartSpeed delivers image quality and speed at your fingertips. It utilizes the Compressed SENSE speed engine to reduce scan time and an award-winning AI algorithm applied directly at the beginning of the MR reconstruction chain to maximize information and deliver outstanding image quality to MSK imaging. Philips SmartSpeed can be used in 2D and 3D and for all anatomical contracts. It supports 97% of current clinical MR protocols1 to address the imaging needs of the vast majority of the patients.



Fast imaging of the knee



Additional information:

- Up to 3 times faster with no loss in image quality²
- Up to 65% higher resolution and improved SNR²
- Improve your productivity with reduced exam time without losing image quality
- Improve signal to noise ratio for challenging anatomies

1 On average, measured across a sample of sites from Philips MR installed 2 Compared to Philips SENSE



SmartSpeed Neuro/Spine Image quality and speed without compromise

Philips SmartSpeed delivers image quality and speed at your fingertips. It utilizes the Compressed SENSE speed engine to reduce scan time and an award-winning Al algorithm applied directly to maximize information and deliver outstanding image quality to brain and spine imaging. Philips SmartSpeed can be used in 2D and 3D and for all anatomical contracts. It supports 97% of current clinical MR protocols¹ to address the imaging needs of the vast majority of the patients.



Fast imaging of the spine



Additional information:

- Up to 3 times faster with no loss in image quality²
- Up to 65% higher resolution and improved SNR²
- Improve your productivity with reduced exam time without losing image quality
- Improve signal to noise ratio for challenging anatomies

1 On average, measured across a sample of sites from Philips MR installed 2 Compared to Philips SENSE



Spectroscopy

Fully integrated proton spectroscopy



Spectroscopy Specialist provides extra information about the spatial distribution of metabolites in the brain. This package provides a set of single voxel, multi-voxel and multi-slice proton spectroscopy, fully integrated into the MRI console. To reduce scan time, a combination of Turbo Spectroscopic Imaging and dS SENSE can be used. Anisotropic matrix can be used to further reduce acquisition time. Includes SpectroView Analysis package for visualization and processing of all spectroscopic data.



Spectroscopic imaging



Spectroscopy XD More precise and more robust MR brain spectroscopy







No VAPOR, amplitude ~3.5 (top) versus VAPOR, amplitude ~0.35 (bottom)



SWIp

Exquisite susceptibility contrast



SWIp has a high sensitivity to enhance contrast for deoxygenated (venous) blood or calcium deposits and may help you, when used in combination with other clinical information, in the diagnosis of various neurological pathologies. SWIp offers high resolution 3D susceptibility weighted brain imaging allowing you to easily integrate it into your mainstream practice.



3D susceptibility weighted brain imaging, including phase maps

- High signal-to-noise ratio¹.
- Includes detailed phase maps to support advanced diagnosis.
SyntAc Exploring neuroradiology with synthetic MR imaging

SyntAc allows to perform MR imaging with a single quantification scan. The acquisition time can be decreased thanks to compatibility with Philips SmartSpeed and compatibility with Compressed SENSE acceleration technology. The resulting data can be used as input for advanced third party processing software¹ to synthesize MR images with different contrasts, brain parenchyma fraction maps and/or brain segmentation maps.





MR images with different contrasts, retrieved from one single quantification scan



VitalScreen

Guidance at your fingertips

VitalScreen (unilateral) offers guidance and insights on the details of the current patient study. This unilateral 12inch interactive touchscreen provides information on exam duration, which coil to use, patient positioning, physiology signal captors (VCG & respiratory) and – if applicable – contrast usage and breath-hold guidance. Moreover, with VitalScreen (unilateral), the patient position can be easily adjusted and the exam can be started at the patient's side with a single touch. When the door of the exam room is closed, scanning will start instantaneously using SmartStart.



Additional information:

- All information and guidance at a glance without the need to navigate between tabs.
- With only one click, the position of the patient can be switched from head-first to feet-first or vice versa.
- Detailed information about contrast injection protocols is readily available through display on VitalScreen (unilateral).
- As soon as the exam-room door is closed, SmartStart automatically starts the scanning procedure.



VitalEye

Touchless respiratory-triggering

VitalEye touchless patient sensing and AI provide a fast and automated detection of patient's breathing patterns without any operator interaction. Technologist receives a continuous and robust respiratory signal without any interaction, helping to keep a caring eye on the patient. The quality of the physiology signal detected by VitalEye is better than a belt-based approach¹ providing superior image quality, for a broad range of patient sizes.



Additional information:

- Technologist no longer needs to set up an old-fashioned respiratory belt.
- Uses adaptive intelligence to detect even the tiniest microscopic signs of breathing.
- Advanced algorithms analyze up to 50 body locations parallel and in real time.
- Produces a respiratory trace that is more robust compared to traditional respiratory triggering methods¹.



Whole Body Fast head to toe scanning

Whole Body package allows rapid, automated whole body imaging with an effective field of view of 2.1 m (7 ft). Combined with Exam Cards, it delivers complete multi-station head to toe coverage in a single pass by combining all imaging sequences per station. Whole Body Specialist supports whole body oncology and vascular imaging studies. Full coverage with a reduced number of stations due to system's large field of view.





Additional information:

- Fast and comfortable whole body imaging
- Perform complete examinations without needing additional surface coils
- Expands applications and increases referral base
- Feet first whole body scanning supporting patient comfort
- No radiation exposure for patients or clinicians
- Single pass head to toe coverage

Zoom Diffusion Small FOV diffusion imaging for improved image quality



Zoom Diffusion allows you to acquire small FOV imaging, down to 200 x 50 mm, with reduced geometrical distortion, due to reduced EPI echo train length in DWI-EPI compared to conventional full FOV DWI-EPI, and higher spatial resolution, due to smaller acquisition voxel size compared to full FOV DWI-EPI, with same level of geometrical distortion.





Small FOV diffusion imaging with high spatial resolution



© 2023 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

4522 991 79861 * OCT 2023

www.philips.com/mrclinicalapplications