DUNLEE

INSTRUCTIONS FOR USE

X-RAY TUBE ASSEMBLIES
OEM

B Σ YOND SOLUTIONS

Table of contents

English.....	5
Deutsch, German.....	25
Français, French.....	45
Español, Spanish.....	67
简体中文, Chinese simplified.....	89
العربية, Arabic.....	109
Български, Bulgarian.....	129
Hrvatski, Croatian.....	151
Česky, Czech.....	171
Dansk, Danish.....	191
Nederlands, Dutch.....	211
Eesti, Estonian.....	233
Suomi, Finnish.....	253
Ελληνικά, Greek.....	273
Magyar, Hungarian.....	295
Íslenska, Icelandic.....	315
Italiano, Italian.....	335
日本語, Japanese.....	355
한국어, Korean.....	375
Latviešu valodā, Latvian.....	395
Lietuvių, Lithuanian.....	415

Malti, Maltese.....	435
Norsk, Norwegian.....	457
Polski, Polish.....	477
Português, Portuguese.....	497
Română, Romanian.....	517
Русский, Russian.....	537
Srpski, Serbian.....	561
Slovensky, Slovak.....	581
Slovenščina, Slovenian.....	601
Svenska, Swedish.....	621
Türkçe, Turkish.....	641
Українська, Ukrainian.....	661

Instructions for Use

X-Ray Tube Assemblies

English

Table of contents

1 Document Information.....	6
2 Intended Purpose.....	7
3 Safety Information.....	8
Safety Message Used in this Document.....	8
General Safety Information.....	8
Limits for Operation, Storage, and Transportation.....	9
Electrical Safety.....	12
Radiation Protection.....	12
Thermal Safety.....	12
Cooling / Insulating Media.....	13
Electromagnetic Compatibility (EMC).....	13
Decommissioning.....	14
China RoHS Label and Material Declaration Table (China only).....	15
Ukraine RoHS Label and Representative (Ukraine only).....	16
Local Representatives.....	16
.....	16
4 Conformity.....	17
5 Compatibility.....	18
6 Maintenance.....	19
Planned Maintenance by the Organization that is Responsible for the Operation.....	19
Cleaning and Disinfection.....	20
Planned Maintenance Schedule.....	21
7 Product Identification Label.....	22

1 Document Information

Document data

Doc-ID: 300004427482A

Release date: 11/2020

Copyright

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Disclaimer

All rights reserved. Koninklijke Philips N.V. reserves the right to make changes in specifications in accordance with the terms in commercial or other applicable agreements between supplier and buyer, and will not be liable for any consequences resulting from the use of this publication. Reproduction in whole or in part is prohibited without the prior written consent of the copyright owner. The English version of this document is the original version.

Product tradename

The duplicate product identification label which accompanies this Instructions For Use provides the tradename of the X-ray tube assembly.

Contact information of the manufacturer

Philips Medical Systems DMC GmbH

Röntgenstraße 24

22335 Hamburg / Germany

For further help see different contact options on: www.dunlee.com

2 Intended Purpose

Intended Use

The X-ray tube housing assembly is an X-ray generating tube encased in a radiation-shielded housing. The X-ray tube housing assembly is intended to emit X-ray photons to enable the creation of diagnostic images of human patients. The X-ray tube housing assembly is only used as integrated part of an imaging system and has no medical purpose by its own. The intended treatment, duration and treatment parameters are not defined on X-ray tube housing assembly level. This is defined by the manufacturer of the imaging system in accordance with the intended use, the intended purpose and the medical purpose as described in the technical documentation of the imaging system.

Intended User Population

The X-ray tube housing assembly is remotely controlled by the imaging system, which itself is operated by qualified clinical users only. Further details regarding the intended user and patient population are defined by the manufacturer of the imaging system in accordance with the intended use, the intended purpose and the medical purpose as described in the technical documentation of the imaging system. The X-ray tube housing assembly must be installed, repaired and maintained only by technically qualified field service engineers who have received the appropriate instructions and training and who are authorized in particular for the installation and commissioning of radiography equipment in medical applications.

Intended Context

The X-ray tube housing assembly is intended to be used in a non-condensing, climate controlled, indoor, clinical environment and is not intended to be used for homecare or by lay persons. The X-ray tube housing assembly is intended to be used in permanently installed stationary imaging systems. The use frequency is not limited within the expected service lifetime. This X-ray tube housing assembly is intended to be used with multiple patients without any special treatment of the X-ray tube housing assembly between examinations. Additional treatment for use with multiple patients may be required and verified by the manufacturer of the imaging system for specific applications.

Operating Principles

Contraindications are not defined on X-ray tube housing assembly level, because this is defined by the manufacturer of the imaging system in accordance with the intended use, the intended purpose and the medical purpose as described in the technical documentation of the imaging system. Specific warnings are provided within the accompanying documents of the X-ray tube housing assembly. The expected service lifetime is referenced in the risk management file and stated in the accompanying documents of the X-ray tube housing assembly.

Essential performance

The X-ray tube assembly does neither have essential performance per se. Nor is any function of the X-ray tube assembly intended to contribute to the essential performance of the X-ray system in which it is integrated.

3 Safety Information

Warning

These instructions for use are made to make it possible to work safely with the X-ray tube assembly. Operate the X-ray tube assembly only in compliance with the instructions in this document and the operating instructions of the applying X-ray system. Do not use it for purposes other than for that it is intended.

Read me,
understand me,
and obey me!

Safety Message Used in this Document

This safety-relevant message occurs in this instruction for use:

WARNING

shows a hazardous situation for people. If you do not prevent this situation, there is a risk of death or serious injury.

General Safety Information

The for the operation responsible organization of the applying X-ray system is always responsible for the compliance with the regulations that apply to the installation and the operation of this X-ray tube assembly.

- Only technically qualified field service engineers who have received the appropriate instructions and training for the equipment in medical applications and in particular for the installation and commissioning of this X-ray tube assembly are allowed to install and repair this equipment. The installation, calibration and testing of the X-ray tube assembly must be performed in accordance with the respective instructions of the applying X-ray system.
- Only trained and authorized personal are permitted to operate the X-ray tube assembly.
- If the X-ray tube assembly has electrical, mechanical, or functional defects, do not use it. This fact particularly applies to faulty indicators, displays, warnings, and alarms.

- If functional defects or other deviations from the normal operational behavior occur, set the X-ray system to OFF. Tell the company that puts the X-ray system that includes this X-ray tube assembly on the market about the defects or deviations immediately. Continue the operation of the X-ray system only after it is repaired. An operation with defective components increases safety risks.
- Do not install the X-ray tube assembly to other medical devices than this X-ray tube assembly is made for.
- The manufacturer of the X-ray tube assembly is responsible for the safety features of this X-ray tube assembly. The manufacturer rejects the responsibility for the safety features of the X-ray tube assembly if modifications are made.
- For the X-ray tube assembly correct operation and regular, competent maintenance are necessary.
- Operate the X-ray tube assembly correctly. Maintain or let maintain the X-ray tube assembly correctly. In other cases, the manufacturer of the X-ray tube assembly is not responsible for incorrect operation, damages, or injuries that occur.
- Do not make or let make service or maintenance at the X-ray tube assembly while the X-ray tube assembly is in use.
- Any serious incident that occurs in relation to the X-ray tube assembly, must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Limits for Operation, Storage, and Transportation

Safe operation of the X-ray tube assembly is only made sure when you use the X-ray tube assembly within the limits of its specifications.

If you ignore the operation limits, there is the danger of coolant that flows out. The X-ray tube assembly becomes too hot. Component parts that are too hot can blow up as a result of an implosion or explosion. Refer to the respective operation limits of the applying X-ray system.

- Do not operate the X-ray tube assembly in presence of a flammable anesthetic mixture with air or oxygen or a nitrous oxide.
- Do not use detergents and disinfectants, also the detergents and disinfectants used on the patient, that can create explosive gas mixtures.
- The oxygen content of the ambient air during operation must be less than 25 %.


- Obey the specified environmental limits for the X-ray tube assembly:

Environmental limits	Standard values	X-ray tube assembly				
		Type A	Type B	Type C	Type D	
In operation						
Ambient temperature [°C]						
Minimum	+5	+5	+10	+15	+15	
Maximum	+35	+30	---	+35	+35	
Maximum with reduced performance	---	---	+40	---	---	
Ambient humidity [%]						
non-condensing						
Minimum	10	10	20	20	10	
Maximum	90	90	80	80	90	
Ambient atmospheric pressure [kPa]						
Minimum about 3,000 m above the sea level	70	70	70	70	70	
Maximum about -750 m below the sea level	110	110	110	110	110	
For storage and transportation						
Temperature limits [°C]						
Minimum	-30	-30	-25	-25	-25	
Maximum	+70	+70	+70	+70	+70	
Humidity limits [%]						
Minimum	5	5	5	5	5	
Maximum	95	95	95	95	95	

Specific values	X-ray tube assembly	Number	Number	Number	Number
Type A	CTR1740	9890 000 8687x			
Type B	DA 1094 DU 694				
	DU 1750 DA 90				
	DU 2550 DA 90 HS				
	DU 33100 DA 90 HS				
	DU 33100-E				
	DU 2550-E				
	DU 1750-E				
Type C	CTR1725	9890 000 8703x	9890 000 8704x	9890 000 8705x	9890 000 8706x
	CTR1725CSBS /				
	CTR1725GSBS				
	CTR1725S				
	CTR1740	9890 000 8697x	9890 000 8698x	9890 000 8699x	9890 000 8700x
		9890 000 8701x	9890 000 8702x		
	CTR2150	4598 012 2188x	4598 012 2189x	4598 013 0155x	4598 013 0159x
		4598 013 0161x	4598 013 0168x	9890 000 8691x	9890 000 8692x
		9890 000 8693x	9890 000 8694x	9890 000 8695x	9890 000 8696x
	CTR2250	4598 013 2542x	4598 013 2544x	9890 000 8711x	9890 000 8712x
		9890 000 8713x	9890 000 8714x		
	CTR2280KANG				
	CTR2280	4598 013 1433x	4598 013 1434x	4598 013 1437x	4598 013 1438x
		9890 000 8709x	9890 000 8710x		
Type D	XS-100 X-Ray Tube Assembly				
	CT8000				

Electrical Safety

Use this X-ray tube assembly only in medical rooms that obey the requirements of relevant local or transregional standards and laws.

Warning:
Electricity

Radiation Protection

Before you make an X-ray exposure, put all necessary radiation precautions into effect. You can find information about radiation precautions in the instructions for use of the applying X-ray system.

- Before you start with radiation, make sure that the X-ray system supplies sufficient proof of a correct image quality.
- To stop radiation, use the emergency stop switch of the X-ray system!

Warning:
Ionizing radiation

Thermal Safety

To prevent damage due to thermal overloading, make sure that you do not operate the X-ray tube assembly outside its specified load parameters. This prevention decreases the risks to the patient, operating staff, third parties, and the environment.

Warning:
Explosive material

Cooling / Insulating Media

- Do not swallow media that spills.
- Do not let media or its vapors to go into your sewerage system.
- Do not breathe the media vapors.
- Make sure to supply closed rooms with fresh air.
- Do not cause vomiting.
- Remove spilled media. Use a liquid-absorbent material.
- Discard the cleaning materials. Obey the local environmental laws and regulations.

Warning:
General

Electromagnetic Compatibility (EMC)

In accordance with its intended use, this electronic product obeys the law that controls EMC. This law gives the permitted emission levels from electronic equipment and its necessary stability against electromagnetic fields.

Electronic products that obey the EMC requirements have a design so that under normal conditions there is no risk of incorrect operation because of electromagnetic interference. If radio signals from high-frequency transmitters with a high transmitting power occur, the risk of electromagnetic incompatibility when operated closely to medical devices can occur.

In unusual circumstances, unintended functions of the electronic product can occur. Possibly undesirable risks for the patient, operating staff, or third parties occur. For this reason, prevent all kinds of transmission with mobile devices. The risks also apply when the medical device is in STANDBY mode.

The X-ray tube assembly is suitable for professional healthcare facility environment.

Too high EMC disturbances can lead to unwanted X-ray, interrupted X-ray, wrong X-ray, unexpected status changes, error messages, component defects, or loss of installation parameters.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the X-ray tube assembly. Otherwise, degradation of the performance of this equipment could result.

The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The for the operation responsible organization might need to take mitigation measures, such as relocating or re-orienting the equipment.

Warning:
Non-ionizing radiation

Decommissioning

The take-back obligation, correct disposal, and recovery of the X-ray tube assembly refer to the European Waste Electrical and Electronic Equipment (WEEE) directive. It also refers to the requirements of local and transregional legal requirements.

The manufacturer of the X-ray tube assembly assembles state-of-the-art X-ray tube assemblies in terms of safety and environmental protection. If no parts of the X-ray tube assembly are opened and if the X-ray tube assembly is used correctly there are no risks to persons or the environment.

To obey regulations, sometimes it is necessary to use materials that are harmful to the environment. Discard these materials in a correct manner.

This X-ray tube assembly contains materials that are toxic. Do not discard the X-ray tube assembly together with industrial or domestic waste.

Obey the disposal regulations!

The manufacturer

- supports you in the disposal of the X-ray tube assembly in accordance with valid legal requirements.
- takes back the X-ray tube assembly.
- returns re-usable parts to the production cycle via certified disposal companies. Extensive test and quality assurance procedures as well as detailed checks of the components guarantee the same high level concerning quality and functionality that is expected from new materials.
- makes a contribution to the protection of the environment.

In case you have questions concerning safe disposal, please consult the manufacturer in full confidence.

China RoHS Label and Material Declaration Table (China only)

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Ukraine RoHS Label and Representative (Ukraine only)

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Philips Medical Systems DMC GmbH
Röntgenstraße 24, 22335 Hamburg, Germany

Ukraine representative for X-ray tubes, which are marked with the conformity mark of Ukraine RoHS Technical Regulation:

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Local Representatives

Contact information for Australia:

Philips Electronics Australia Ltd
65 Epping Road
North Ryde
NSW 2113
Australia

4 Conformity

The X-ray tube assembly that is shown in these instructions for use puts the provisions of the

- European CE Marking of Conformity, and

Identification num-
ber of the notified
body

- FDA 21 CFR 1020.30

into effect.

If you have further questions that refer to local or transregional legal requirements,
make a contact to:

Philips Medical Systems DMC GmbH
Röntgenstraße 24
22335 Hamburg / Germany

For further help see different contact options on: <http://www.dunlee.com>

5 Compatibility

Only technically qualified field service engineers who have received the appropriate instructions and training for the equipment in medical applications and in particular for the installation and commissioning of this X-ray tube assembly are allowed to install and repair this equipment. The installation, calibration and testing of the X-ray tube assembly must be performed in accordance with the respective instructions of the applying X-ray system.

To make sure that the X-ray tube assembly and the X-ray system are compatible, refer to the technical instructions of the X-ray system.

If necessary, speak to the manufacturer of the applying X-ray system.

The company that puts the X-ray system that includes this X-ray tube assembly on the market is responsible for the installation of the X-ray tube assembly to the X-ray system. The installation must agree with all local and transregional legal requirements.

6 Maintenance

Only technically qualified field service engineers who have received the appropriate instructions and training for the equipment in medical applications and in particular for the maintenance of this X-ray tube assembly are allowed to do planned and corrective maintenance on this equipment. The maintenance inclusive calibration and testing of the X-ray tube assembly must be performed in accordance with the respective instructions of the applying X-ray system.

The manufacturer of this X-ray tube assembly recommends measures for planned maintenance. The recommended measures must agree with the local and transregional legal requirements.

The expected service life time is 10 years. This requires that the X-ray tube assembly is subject to service maintenance according to the requirements of this service instructions and the service instructions of the applying X-ray system.

Planned Maintenance by the Organization that is Responsible for the Operation

The for the operation responsible organization of medical devices is committed to do planned maintenance. This planned maintenance must refer to the regulations for the prevention of industrial accidents, local or transregional medical product laws, and further regulations.

This X-ray tube assembly requires regular planned maintenance:

Make sure that service personnel do planned maintenance on this X-ray tube assembly regularly.

Precautionary planned maintenance measures prevent personal injury and make sure that the for the operation responsible organization puts all commitments into effect.

If functional defects or other deviations from the normal operational behavior occur, set the X-ray system to OFF. Tell the company that puts the X-ray system that includes this X-ray tube assembly on the market about the defects or deviations immediately. Continue the operation of the X-ray system only after it is repaired. An operation with defective components increases safety risks.

Regular maintenance cycles decrease the probability of a spontaneous failure. You cannot prevent spontaneous failures completely. If it is necessary to have the guaranteed availability of the X-ray system, make sure that a redundant standby X-ray system is available.

The structural integrity of the X-ray tube assembly does not wear.

Cleaning and Disinfection

Hazard of overdosed radiation or electrical shock.

- Set the X-ray system to OFF before you start the cleaning procedure of the X-ray tube assembly.
- Do not open the covers of the X-ray system.
- Do not change or remove the covers or other parts of the X-ray tube assembly.
The X-ray tube assembly must be directly accessible.
- Do not change the configuration of the X-ray tube assembly.
- Make sure that no detergents get inside the X-ray equipment.
- If you use cleaning agents such as detergents, make sure that they do not contain explosive substances as they can create explosive gas mixtures.
- Obey the regional legal requirements that refer to the cleaning procedure.
- Obey the manufacturers information and safety instructions of the detergents.

If you do not obey this instruction, there is a risk of death or serious injury.

Cleaning

Only if

- the X-ray tube assembly is not behind covers of the X-ray system,
- the X-ray tube assembly has no covers itself,
- the X-ray tube assembly has direct access,

clean the X-ray tube assembly.

Cleaning procedure

- Do not use detergents or polishes that can cause corrosion, delamination, or abrasion.
 - Do not use caustic, solvent, or abrasive detergents or polishes.
 - Do not use detergents or polishes based on alcohol.
 - Clean enameled parts, chrome parts, and aluminum surfaces only with a moist cloth and a mild detergent or polishes.
 - Wipe with a dry woolen cloth.
-

Disinfection

A disinfection of the X-ray tube assembly is not necessary, unless the system manufacturer requires it.

Planned Maintenance Schedule


- Do checks of the X-ray tube assembly for apparent defects. Do the checks unless they are contrary to the configuration of the X-ray system.

Interval	Scope of work
Always you operate the X-ray system	<p>Do a check if error messages of the X-ray system occur. If the X-ray tube assembly is thermally overloaded and does not operate, wait until a cooling period lets you operate the X-ray tube assembly again.</p> <p>Do a check for visible media leakage or other contamination.</p>
Weekly	Do a check for unusual noises.
Annually	Only if the X-ray tube assembly has free access, clean it.
As per relevant local or transregional standards and laws	Do a constancy test.
As per relevant local or transregional standards and laws	Do a check for image quality.

- If defects occur that refer to the X-ray tube assembly, make a contact to the company that puts the X-ray system that includes this X-ray tube assembly on the market.

7 Product Identification Label

The performance characteristics of the device are defined on the product identification label as part of the accompanying documents.

No.	Label	Designation
01		Registered trade name / trade mark of the manufacturer
02	Philips Medical Systems DMC GmbH Röntgenstraße 24 22335 Hamburg / GERMANY	Address of the legal manufacturer
03		Address of the legal manufacturer
	yyyy-mm	yyyy: year of manufacture mm: month of manufacture
04	X-RAY TUBE HOUSING ASSEMBLY	Type of the medical device: X-ray tube housing assembly
05	Model	Name of the product
06	REF	Type number of the product
07	SN	Serial Number of the product
08	MD	Medical device product
09		Value of the permanent filtration
10		Dimensions of the small focus and the respective standard
11		Dimensions of the large focus and the respective standard
12	Nominal X-ray Tube Voltage	Value of the nominal X-ray tube voltage and the respective standard
13	X-RAY TUBE	X-ray tube
14	US only	Compliance statement, for US only
THIS PRODUCT COMPLIES WITH THE DHHS REQUIREMENT OF 21 CFR SUB-CHAPTER J		

No.	Label	Designation
15	MANUFACTURED: month year	Month and year of manufacture
16	 UDI	The Unique Device Identification [UDI] Data Matrix Code consists of the Global Trade Item Number [GTIN], Serial Number [SN], and Type Number [TN] of the product.
17	(01)xxxxxxxxxxxx	Global Trade Item Number [GTIN]
18	(21)xxxxxxxxxxxx	Serial Number [SN] of the product xxxxxwyyyyy xxxxx : X-ray tube housing w : X-ray tube housing type yyyyy : X-ray tube
19		CSA symbol
20		CE symbol and the number of the notified body
21		Disposal symbol
22		General warning symbol
23		Electronic version of the Instructions For Use. See: www.dunlee.com/IFU
24		China RoHS symbol
25	Duplicate label	Duplicate label for the label overview of the X-ray system
26	(CPO) Refurbished	Certified pre-owned refurbished product
27		Additional label: Instruction manual/booklet must be read
28	 UA.TR.116	Ukraine RoHS Technical Regulation symbol

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All rights reserved. Koninklijke Philips N.V.は、納入業者と購入者の間の商業契約またはその他の適用される契約の条件に従って仕様を変更する権利を留保します。また、この印刷物を使用したことの結果について一切の責任を負いません。著作権保有者から事前に書面による許可を得ることなく本書を複製・転写することは、その全部或いは一部に限らず禁じます。本書の英語版が原本です。