Instructions for Use

X-Ray Tube Housing Assemblies

OEM

Table of contents

English.....	4
Deutsch	25
Français	47
Español.....	69
简体中文.....	91
العربية.....	112
Български.....	133
Hrvatski.....	155
Česky.....	176
Dansk	197
Nederlands.....	219
Eesti	241
Suomi	262
Ελληνικά	284
Magyar	306
Íslenska	329
Italiano.....	350
日本語.....	371
한국어.....	392
Latviešu valodā.....	413
Lietuvių	434

Malti	456
Norsk	478
Polski	499
Português Brasileiro	521
Română	543
Русский	565
Srpski	587
Slovensky	608
Slovenščina	630
Svenska	651
Türkçe	672
Українська	693

Instructions for Use

English

Table of contents

1	Document information	5
2	Intended purpose	6
3	Safety information	7
	Safety message used in this document	7
	General safety information	7
	Limits for operation, storage, and transportation.....	8
	Electrical safety	12
	Radiation protection	12
	Thermal safety.....	12
	Cooling / Insulating media.....	12
	Decommissioning	13
	China RoHS label and material declaration table (China only).....	14
	Ukraine RoHS label and representative (Ukraine only).....	15
	Local representatives	15
4	Conformity	16
5	Compatibility	17
6	Technical data of X-ray tube housing assembly	18
7	Maintenance	19
	Planned maintenance by the organization that is responsible for the operation.....	19
	Cleaning and disinfection	20
	Planned maintenance schedule.....	21
8	Product identification label	22

1 Document information

Document data

Release date: 10/2023

Copyright

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Disclaimer

All rights reserved. Koninklijke Philips N.V. reserves the right to make changes in specifications in accordance with the terms in commercial or other applicable agreements between supplier and buyer, and will not be liable for any consequences resulting from the use of this publication. Reproduction in whole or in part is prohibited without the prior written consent of the copyright owner. The English version of this document is the original version.

Product tradename

The duplicate product identification label which accompanies this Instructions For Use provides the tradename of the X-ray tube housing assembly.

Contact information of the manufacturer

Philips Medical Systems DMC GmbH

Röntgenstraße 24

22335 Hamburg / Germany

For further help see different contact options on: www.dunlee.com

2 Intended purpose

Intended Use

The X-ray tube housing assembly is an X-ray generating tube encased in a radiation-shielded housing. The X-ray tube housing assembly is intended to emit X-ray photons to enable the creation of diagnostic images of human patients. The X-ray tube housing assembly is only used as integrated part of an imaging system and has no medical purpose by its own. The intended treatment, duration and treatment parameters are not defined on X-ray tube housing assembly level. This is defined by the manufacturer of the imaging system in accordance with the intended use, the intended purpose and the medical purpose as described in the technical documentation of the imaging system.

Intended User Population

The X-ray tube housing assembly is remotely controlled by the imaging system, which itself is operated by qualified clinical users only. Further details regarding the intended user and patient population are defined by the manufacturer of the imaging system in accordance with the intended use, the intended purpose and the medical purpose as described in the technical documentation of the imaging system. The X-ray tube housing assembly must be installed, repaired and maintained only by technically qualified field service engineers who have received the appropriate instructions and training and who are authorized in particular for the installation and commissioning of radiography equipment in medical applications.

Intended Context

The X-ray tube housing assembly is intended to be used in a non-condensing, climate controlled, indoor, clinical environment and is not intended to be used for homecare or by lay persons. The X-ray tube housing assembly is intended to be used in permanently installed stationary imaging systems. The use frequency is not limited within the expected service lifetime. This X-ray tube housing assembly is intended to be used with multiple patients without any special treatment of the X-ray tube housing assembly between examinations. Additional treatment for use with multiple patients may be required and verified by the manufacturer of the imaging system for specific applications.

Operating Principles

Contraindications are not defined on X-ray tube housing assembly level, because this is defined by the manufacturer of the imaging system in accordance with the intended use, the intended purpose and the medical purpose as described in the technical documentation of the imaging system. Specific warnings are provided within the accompanying documents of the X-ray tube housing assembly. The expected service lifetime is referenced in the risk management file and stated in the accompanying documents of the X-ray tube housing assembly.

Essential performance

The X-ray tube housing assembly does neither have essential performance per se. Nor is any function of the X-ray tube housing assembly intended to contribute to the essential performance of the X-ray system in which it is integrated.

3 Safety information

Read me, understand me, and obey me!

These instructions for use are made to make it possible to work safely with the X-ray tube housing assembly. Operate the X-ray tube housing assembly only in compliance with the instructions in this document and the operating instructions of the applying X-ray system. Do not use it for purposes other than for that it is intended.

Safety message used in this document

General warnings and safety notices are placed in this instruction for use. Special warnings and safety notices, that applies only to a special X-ray tube housing assembly, are placed in the related Technical Manual of the X-ray tube housing assembly.

This safety-relevant message occurs in this instruction for use:

WARNING

Shows a hazardous situation for people. If you do not prevent this situation, there is a risk of death or serious injury.

General safety information

The for the operation responsible organization of the applying X-ray system is always responsible for the compliance with the regulations that apply to the installation and the operation of this X-ray tube housing assembly.

- Only technically qualified field service engineers who have received the appropriate instructions and training for the equipment in medical applications and in particular for the installation and commissioning of this X-ray tube housing assembly are allowed to install and repair this equipment. The installation, calibration and testing of the X-ray tube housing assembly must be performed in accordance with the respective instructions of the applying X-ray system.
- Only trained and authorized personal are permitted to operate the X-ray tube housing assembly.
- If the X-ray tube housing assembly has electrical, mechanical, or functional defects, do not use it. This fact particularly applies to faulty indicators, displays, warnings, and alarms.
- If functional defects or other deviations from the normal operational behavior occur, set the X-ray system to OFF. Tell the company that puts the X-ray system that includes this X-ray tube housing assembly on the market about the defects or deviations immediately. Continue the operation of the X-ray system only after it is repaired. An operation with defective components increases safety risks.
- Do not install the X-ray tube housing assembly to other medical devices than this X-ray tube housing assembly is made for.

- The manufacturer of the X-ray tube housing assembly is responsible for the safety features of this X-ray tube housing assembly. The manufacturer rejects the responsibility for the safety features of the X-ray tube housing assembly if modifications are made.
- For the X-ray tube housing assembly correct operation and regular, competent maintenance are necessary.
- Operate the X-ray tube housing assembly correctly. Maintain or let maintain the X-ray tube housing assembly correctly. In other cases, the manufacturer of the X-ray tube housing assembly is not responsible for incorrect operation, damages, or injuries that occur.
- Do not make or let make service or maintenance at the X-ray tube housing assembly while the X-ray tube housing assembly is in use.
- Any serious incident that occurs in relation to the X-ray tube housing assembly, must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Limits for operation, storage, and transportation

Safe operation of the X-ray tube housing assembly is only made sure when you use the X-ray tube housing assembly within the limits of its specifications.

If you ignore the operation limits, there is the danger of coolant that flows out. The X-ray tube housing assembly becomes too hot. Component parts that are too hot can blow up as a result of an implosion or explosion. Refer to the respective operation limits of the applying X-ray system.

- Do not operate the X-ray tube housing assembly in presence of a flammable anesthetic mixture with air or oxygen or a nitrous oxide.
- Do not use detergents and disinfectants, also the detergents and disinfectants used on the patient, that can create explosive gas mixtures.
- The oxygen content of the ambient air during operation must be less than 25 %.
- Obey the specified environmental limits for the X-ray tube housing assembly:

Environmental limits	X-ray tube housing assembly				
	Standard values ⁽¹⁾	Specific values			
		Type A	Type B	Type C	Type D
In operation					
Ambient temperature [°C]					
Minimum	+5	+10	+15	+15	+18
Maximum	+35	N/A	+35	+35	+34
Maximum with reduced performance	N/A	+40	N/A	N/A	N/A
Ambient humidity [%]					
non-condensing					
Minimum	10	20	20	10	10
Maximum	90	80	80	90	90
Ambient atmospheric pressure [kPa]					
Minimum (corresponds to a height of about 3,000 m above sea level)	70	70	70	70	70
Maximum	110	110	110	110	110
<p>(1) The values of most X-ray tube housing assemblies (default values) are set in this column. X-ray tube housing assemblies with non-standard limit values are given in the other columns.</p>					

Tab. 1: Environmental limits in operation

300004427486 / D000411792 Rev R / * 10.2023

Dunlee

Environmental limits	X-ray tube housing assembly				
	Standard values ⁽¹⁾	Specific values			
		Type A	Type B	Type C	Type D
For storage and transportation					
Temperature limits [°C]					
Minimum	-30	-25	-25	-25	-25
Maximum	+70	+70	+70	+70	+70
Humidity limits [%]					
Minimum	5	5	5	5	5
Maximum	95	95	95	95	95
Atmospheric pressure [kPa]					
Minimum	70	70	70 / 50 ⁽²⁾	30	70
Maximum	110	110	110	110	110
<p>(1) The values of most X-ray tube housing assemblies (default values) are set in this column. X-ray tube housing assemblies with non-standard limit values are given in the other columns.</p> <p>(2) Only valid for CTR1735, CTR2150 and CTR2150CEPN.</p>					

Tab. 2: Environmental limits for storage and transportation

Specific values	X-ray tube housing assembly			
Type A	DA 1094 DU 694			
	DU 1750 DA 90	DU 1750-E		
	DU 2550 DA 90 HS	DU 2550-E		
	DU 33100 DA 90 HS	DU 33100-E		
Type B	CTR1725	CTR1725CSBS	CTR1725GSBS	CTR1725S
	CTR1740	CTR1740CQNM		
	CTR2140CENM			
	CTR2150	CTR2150CEAN	CTR2150CESV	CTR2150CEUH
	CTR2250			
	CTR2251			
	CTR2280	CTR2280CESV	CTR2280KANG	
	Type C	XS-100 X-Ray Tube Assembly		
CT8000				
Type D	CT4000			
	CT3000			

300004427486 / D000411792 Rev R / * 10.2023

Dunlee

Electrical safety

Warning: Electricity

Use this X-ray tube housing assembly only in medical rooms that obey the requirements of relevant local or transregional standards and laws.

Radiation protection

Warning: Ionizing radiation

Before you make an X-ray exposure, put all necessary radiation precautions into effect. You can find information about radiation precautions in the instructions for use of the applying X-ray system.

- Before you start with radiation, make sure that the X-ray system supplies sufficient proof of a correct image quality.
- To stop radiation, use the emergency stop switch of the X-ray system!

Thermal safety

Warning: Explosive material

To prevent damage due to thermal overloading, make sure that you do not operate the X-ray tube housing assembly outside its specified load parameters. This prevention decreases the risks to the patient, operating staff, third parties, and the environment.

Cooling / Insulating media

Warning: General

- Do not swallow media that spills!
- Do not breathe the media vapors!
- In case of swallowed media:
 - Do not cause vomiting!
 - Immediately call a POISON CENTER or doctor/physician!
- Do not let media or the can go into your sewerage system! Obey the local environmental laws and regulations!

Decommissioning

Obey the disposal regulations!

The take-back obligation, correct disposal, and recovery of the X-ray tube housing assembly refer to the European Waste Electrical and Electronic Equipment (WEEE) directive. It also refers to the requirements of local and transregional legal requirements.

The manufacturer of the X-ray tube housing assembly assembles state-of-the-art X-ray tube housing assemblies in terms of safety and environmental protection. If no parts of the X-ray tube housing assembly are opened and if the X-ray tube housing assembly is used correctly there are no risks to persons or the environment.

To obey regulations, sometimes it is necessary to use materials that are harmful to the environment. Discard these materials in a correct manner.

This X-ray tube housing assembly contains materials that are toxic. Do not discard the X-ray tube housing assembly together with industrial or domestic waste.

The manufacturer

- supports you in the disposal of the X-ray tube housing assembly in accordance with valid legal requirements.
- takes back the X-ray tube housing assembly.
- returns re-usable parts to the production cycle via certified disposal companies. Extensive test and quality assurance procedures as well as detailed checks of the components guarantee the same high level concerning quality and functionality that is expected from new materials.
- makes a contribution to the protection of the environment.

In case you have questions concerning safe disposal, please consult the manufacturer in full confidence.

China RoHS label and material declaration table (China only)

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300004427486 / D000411792 Rev/R / * 10.2023

Dunlee

Ukraine RoHS label and representative (Ukraine only)

Ukraine RoHS Technical Regulation symbol

Address of the legal manufacturer:

Philips Medical Systems DMC GmbH
Röntgenstraße 24, 22335 Hamburg, Germany

Ukraine representative for X-ray tubes, which are marked with the conformity mark of Ukraine RoHS Technical Regulation:

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Local representatives

Contact information for Australia:

Philips Electronics Australia Ltd
65 Epping Road
North Ryde
NSW 2113
Australia

4 Conformity

The X-ray tube housing assembly that is shown in these instructions for use puts the provisions of the

Identification number of the notified body

- European CE Marking of Conformity, and
- FDA 21 CFR 1020.30

into effect.

If you have further questions that refer to local or transregional legal requirements, make a contact to:

Philips Medical Systems DMC GmbH
Röntgenstraße 24
22335 Hamburg / Germany

For further help see different contact options on: <http://www.dunlee.com>

5 Compatibility

Only technically qualified field service engineers who have received the appropriate instructions and training for the equipment in medical applications and in particular for the installation and commissioning of this X-ray tube housing assembly are allowed to install and repair this equipment. The installation, calibration and testing of the X-ray tube housing assembly must be performed in accordance with the respective instructions of the applying X-ray system.

To make sure that the X-ray tube housing assembly and the X-ray system are compatible, refer to the technical instructions of the X-ray system.

If necessary, speak to the manufacturer of the applying X-ray system.

The company that puts the X-ray system that includes this X-ray tube housing assembly on the market is responsible for the installation of the X-ray tube housing assembly to the X-ray system. The installation must agree with all local and transregional legal requirements.

6 Technical data of X-ray tube housing assembly

Technical data of X-ray tube housing assembly for example:

- Single Load Rating
- Serial Load Rating
- Nominal Radiographic Anode Input Power
- Nominal CT Anode Input Power
- Nominal CT Scan Power Index

are specified in the related Technical Data / Technical Manual of the X-ray tube housing assembly, if applicable.

7 Maintenance

Only technically qualified field service engineers who have received the appropriate instructions and training for the equipment in medical applications and in particular for the maintenance of this X-ray tube housing assembly are allowed to do planned and corrective maintenance on this equipment. The maintenance inclusive calibration and testing of the X-ray tube housing assembly must be performed in accordance with the respective instructions of the applying X-ray system.

The manufacturer of this X-ray tube housing assembly recommends measures for planned maintenance. The recommended measures must agree with the local and transregional legal requirements.

The expected service life time is specified in the related Technical Manuals. It is necessary that during this period the X-ray tube housing assembly is subject to the service maintenance. Refer to the requirements of the service instructions and the service instructions of the applying X-ray system.

Planned maintenance by the organization that is responsible for the operation

The organization responsible for the operation of the applying X-ray system is responsible to do planned maintenance. This planned maintenance must refer to the regulations for the prevention of industrial accidents, local or transregional medical product laws, and further regulations.

This X-ray tube housing assembly requires regular planned maintenance:

Make sure that service personnel do planned maintenance on this X-ray tube housing assembly regularly.

Precautionary planned maintenance measures prevent personal injury and make sure that the for the operation responsible organization puts all commitments into effect.

If functional defects or other deviations from the normal operational behavior occur, set the X-ray system to OFF. Tell the company that puts the X-ray system that includes this X-ray tube housing assembly on the market about the defects or deviations immediately. Continue the operation of the X-ray system only after it is repaired. An operation with defective components increases safety risks.

Regular maintenance cycles decrease the probability of a spontaneous failure. You cannot prevent spontaneous failures completely. If it is necessary to have the guaranteed availability of the X-ray system, make sure that a redundant standby X-ray system is available.

The structural integrity of the X-ray tube housing assembly does not wear.

Cleaning and disinfection

WARNING

Hazard of overdosed radiation or electrical shock

- Set the X-ray system to OFF before you start the cleaning procedure of the X-ray tube housing assembly.
- Do not open the covers of the X-ray system.
- Do not change or remove the covers or other parts of the X-ray tube housing assembly.
- Do not change the configuration of the X-ray tube housing assembly.
- Make sure that no detergents get inside the X-ray equipment.
- If you use cleaning agents such as detergents, make sure that they do not contain explosive substances as they can create explosive gas mixtures.
- Obey the regional legal requirements that refer to the cleaning procedure.
- Obey the manufacturers information and safety instructions of the detergents.

If you do not obey this instruction, there is a risk of death or serious injury.

Cleaning

Only if

- the X-ray tube housing assembly is not behind covers of the X-ray system,
- the X-ray tube housing assembly has no covers itself,

clean the X-ray tube housing assembly.

Cleaning procedure

- Do not use detergents or polishes that can cause corrosion, delamination, or abrasion.
- Do not use caustic, solvent, or abrasive detergents or polishes.
- Do not use detergents or polishes based on alcohol.
- Clean enameled parts, chrome parts, and aluminum surfaces only with a moist cloth and a mild detergent or polish.
- Wipe with a dry woolen cloth.

Disinfection

A disinfection of the X-ray tube housing assembly is not necessary.

WARNING

Hazard to device safety

- Do not use direct UV-C light with X-ray tube housing assemblies.
- Do not use disinfection methods that cause the creation of ozone.

Planned maintenance schedule


- Do checks of the X-ray tube housing assembly for apparent defects. Do the checks unless they differ from the corresponding instructions of the X-ray system.

Interval	Scope of work
When you use the X-ray system	Do a check if error messages of the X-ray system occur. If the X-ray tube housing assembly is thermally overloaded and does not operate, wait until a cooling period lets you operate the X-ray tube housing assembly again. Do a check for visible media leakage or other contamination.
Weekly	Do a check for unusual noises.
Annually	Only if the X-ray tube housing assembly has free access, clean it.
As per relevant local or transregional standards and laws	Do a constancy test. Do a check for image quality.

- If defects occur that refer to the X-ray tube housing assembly, make a contact to the company that puts the X-ray system that includes this X-ray tube housing assembly on the market.

8 Product identification label



The performance characteristics of the device are defined on the product identification label as part of the accompanying documents.

No.	Label	Designation
1		Registered trade name / trade mark of the manufacturer
2	Philips Medical Systems DMC GmbH Röntgenstraße 24 22335 Hamburg / GERMANY	Address of the legal manufacturer
3	 yyyy-mm	Address of the legal manufacturer yyyy: year of manufacture (if applicable) mm: month of manufacture (if applicable)
4	X-RAY TUBE HOUSING ASSEMBLY	Type of the medical device: X-ray tube housing assembly
5	Model or #	Name of the product
6	REF	Type number of the product / Reference Number
7	PN	Part Number
8	SN	Serial Number of the product
9	MD	Medical device product
10		Value of the permanent filtration
11		Dimensions of the small focus and the respective standard
12		Together with the symbol  : Dimensions of the smaller focus and the respective standard. Together with the symbol  : Dimensions of the larger focus and the respective standard.
13		Dimensions of the large focus and the respective standard
14	Nominal X-ray Tube Voltage	Value of the nominal X-ray tube voltage and the respective standard
15	X-RAY TUBE	X-ray tube

No.	Label	Designation
16		Compliance statement, for US only
	THIS PRODUCT COMPLIES WITH THE DHHS REQUIREMENT OF 21 CFR SUB-CHAPTER J	
17	MANUFACTURED: month year	Month and year of manufacture
18	 	The Unique Device Identification [UDI] Data Matrix Code consists of the Global Trade Item Number [GTIN], Serial Number [SN], and Type Number [TN] of the product.
19	(01) xxxxxxxxxxxxxxxx	Global Trade Item Number [GTIN]
20	(21) xxxxxxxxxxxx	Serial Number [SN] of the product xxxxxwyyyyyy xxxxxw : X-ray tube housing (w optional) yyyyyy : X-ray tube
21	 or 	CSA symbol CSA symbol with CSA master contract number of the manufacturer
22		CE symbol and the number of the notified body
23		Disposal symbol
24		General warning symbol
25	 	Electronic version of the Instructions For Use. See: www.dunlee.com/IFU
26		China RoHS symbol
27	Duplicate label	Duplicate label for the label overview of the X-ray system
28	(CPO) Refurbished	Certified pre-owned refurbished product

300004427486 / D000411792 Rev R / * 10.2023

Dunlee

No.	Product specific optional additional labels:	Designation
1		Instruction manual/booklet must be read
2		Marking for the functional earth terminal
3		Marking for the protective earth terminal
4		Warning of electricity
5		Warning of toxic substances
6		Warning of hot surface
7		High voltage symbol
8		Safety sign, that indicates a risk of aspiration
9		Polarity sign
10		Marking Anode side
11		Marking Cathode side
12		Focal spot location on the X-ray tube housing assembly
13		UKCA marking
14		Address of manufacturing site/plant
15		Ukraine RoHS Technical Regulation symbol

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All rights reserved. Koninklijke Philips N.V.は、納入業者と購入者の間の商業契約またはその他の適用される契約の条件に従って仕様を変更する権利を留保します。また、この印刷物を使用したことの結果について一切の責任を負いません。著作権保有者から事前に書面による許可を得ることなく本書を複製・転写することは、その全部或いは一部に限らず禁じます。本書の英語版が原本です。
All rights reserved. Koninklijke Philips N.V.는 공급업체와 구매자 간의 상업적 또는 기타 적용 가능한 계약 조건에 따라 사양을 변경할 권한을 보유하며 이 간행물의 사용으로 인한 결과에 어떤 책임도 지지 않습니다. 저작권자의 사전 서면 동의 없이 전체 또는 일부를 복사하는 것은 금지됩니다. 이 문서의 영문 버전이 원본입니다.
All rights reserved. Reproduction or transmission in whole or in part of (the information in) this document, in any form or by any means, electronic, mechanical or otherwise, is prohibited without the prior written consent of the copyright owner.