

PHILIPS

肩頸按摩器

頸部按摩器

PPM3501



Quick Start Guide 快速入門指南

1. Press and hold the power button for 2 seconds to turn on, enter the low-speed massage, soothing mode and the heating function. the indicator light is orange, and the sound of the waves is played.



Indicator light in orange

2. Short press the power button to turn off/on the massage roller.



3. Short press the heating button to turn off the heating function, and the indicator light turns blue; short press again to turn on the heating function, and the indicator light turns orange; long press the heating button to turn off/on the voice prompt.



4. Short press the mode button to switch between soothing mode, vitality mode, strong mode, neck protection mode, and wake-up mode.



5. Short press the sound button to switch the music in sequence: fireplace, wind chime, stream, spring rain, early morning, ocean waves; long press this button to turn off/on the music.



6. The voice setting has a memory function. If it is set to enable voice/music before powering off, the voice/music will be turned on by default after the next power on.



7. Press and hold the power button for 2 seconds to shut down (if the voice is on, it will sound, if the voice is off, it will not prompt).



8. In the operating state, if the battery power is $\leq 10\%$, the blue indicator light of the massager will flash for 5 seconds and shut down; in the shutdown state, if the battery power is $\leq 10\%$, the blue indicator light of the massager will flash for 5 seconds and cannot be started.



Indicator light flashing in blue

- ❗ 1. Please remove personal accessories when using the massager, so as not to affect the massage part.
- 2. If you feel discomfort during massage, please stop operation immediately.
- 3. When the massager is not in use, please turn it off. In case of a malfunction, please turn off the power to avoid accidents.
- 4. When the massager being charged, the blue light ring will be breathing and flashing, when fully charged, the blue light ring will stay on;
- 5. The massager will be automatically turned off after 15 minutes of operation.



請在以下網站連接中了解詳情

www.philips.com

本手冊內容與實物可能存在差異，請以實物為準，如產品有任何更新恕不另行通知。

©2022 AquaShield

All rights reserved

The Philips trademark and the Philips Shield Emblem are registered trademarks of Koninklijke Philips N.V. used under license.

This product has been manufactured by and is sold under the responsibility of HongKong AquaShield Health Technology Company Limited and HongKong AquaShield Health Technology Company Limited is the warrantor in relation to this product.

水護盾是飛利浦淨飲水解決方案、熱水解決方案、衛浴解決方案、水杯及個人按摩器唯一品牌官方授權運營主體。

Rev A JUL 22



1. 長按電源鍵2秒開機，進入低轉速按摩、舒緩模式，加熱功能開啟，指示燈為橙色，播放海浪聲。



指示燈亮橙色

2. 短按電源鍵，關閉/開啟按摩滾輪。



3. 短按加熱鍵關閉加熱功能，此時指示燈亮藍色；再次短按可開啟加熱功能，指示燈亮橙色；長按加熱鍵可關閉/開啟語音提醒。



4. 短按模式鍵，可切換舒緩模式、活力模式、強力模式、護頸模式、喚醒模式。



5. 短按聲音鍵，可依次切換音樂：壁爐，風鈴，溪流，春雨，清晨，海浪；長按此鍵關閉/開啟音樂。



6. 語音設置帶有記憶功能，若開機前設置為開啟語音/音樂，則下一次開啟後默認開啟語音/音樂；若開機前設置為關閉語音/音樂，則下一次開啟後默認關閉語音/音樂。



7. 長按電源鍵2秒開機（若語音開啟狀態下，則響起，若語音處於關閉狀態，則不提示）。



8. 操作狀態下，若電量 $\leq 10\%$ ，按摩器藍色指示燈頻閃5秒並開機；開機狀態下，若電量 $\leq 10\%$ ，按摩器藍色指示燈頻閃5秒且無法啟動。



指示燈頻閃藍色



1. 使用時請除去隨身飾物，以免影響按摩部位。
2. 若在按摩過程中感到不適請立即停止操作。
3. 若不使用按摩器時，請將電源關閉。若發生機器故障，請關掉開關，以免發生意外。
4. 充電時，藍色燈圈呼吸閃爍，充滿後藍色燈圈長亮。
5. 工作15分鐘自動關機。