

Computed Tomography

Whitepaper



Detector-based spectral CT delivers on the promise of increased accuracy for proton therapy Inaccuracy in proton stopping-power ratios (SPR) calculated from conventional HU-based CT images, and the associated proton beam range uncertainty, is a well-known limitation in realizing the full potential of proton therapy. This error has been documented in literature to be approximately 3%. Multi-energy CT has shown promise in reducing SPR inaccuracy and proton range uncertainty. This white paper analyzes peer-reviewed publications that document the accuracy of electron density (ED) and effective atomic number (Z-effective) images from various multi-energy technologies. It then compares the resulting SPR accuracy from these different inputs. The analysis shows a range of quantitative accuracy in the ED and Z-effective images across multi-energy CT technologies, and the resulting SPR values, with dual-layer detector-based spectral CT providing the most accurate results. This analysis demonstrates a spectral SPR uncertainty of less than 1% for detector-based spectral CT.

Over the past decade, the use of proton therapy has been increasing rapidly. The sharp dose cutoff of the proton Bragg peak allows for sharper dose delivery potential than would be possible in external beam photon therapy.^{1,2} As a result, surrounding organs at risk can be more effectively spared, decreasing potential side effects. Alternatively, it enables dose escalation, if necessary.

Locating the proton Bragg peak in the body requires an estimation of the proton stopping-power of tissues. Historically, conventional CT Hounsfield units (HU) are translated to proton stopping-power ratios (SPR) through stoichiometric calibration curves. However, this HU-based SPR assessment is associated with an uncertainty of 3%, which translates to an uncertainty of 3 mm for a target region at 10 cm depth.³ With this, the precision of proton energy deposition is not fully exploited, detracting from the overall potential of proton therapy.

In contrast to conventional HU-based indirect SPR calculations, novel CT techniques can allow for a direct SPR calculation.³⁻¹⁰ These techniques can be categorized as source-based dual-energy CT (DECT) or detector-based spectral CT. The second category consists of dual-layer CT (DLCT) and photon-counting CT (PCCT), employing energy-integrating or photon-counting detectors, respectively. To calculate SPR with the aforementioned multi-energy CT techniques, two datasets are needed: Electron density (ED) and effective atomic number (Z_{eff}). Subsequently, SPR is directly calculated from a formula like the Bethe-Bloch equation. To minimize SPR errors, ED and Z_{eff} values must be highly accurate and quantitative.¹¹

Recent studies have shown that the accuracy and precision of ED and Z_{eff} depend on the technique used.^{12–14} For the source-based DECT techniques (i.e. spin-spin, kVp-switching, split-filter, or dual-source), acquisition of spectral data has to be predefined. The reason for this is the inherent compromise, such as a reduced gantry rotation time, limited reconstruction field-of-view, increased radiation dose, or reduced temporal resolution. On the other hand, for detector-based spectral CT techniques (i.e. DLCT and PCCT), presetting is not required as every scan provides spectral data (in addition to true conventional images).

For DLCT, direct SPR calculation is possible with the use of the Philips MultiModality Simulation Workspace (MM SIM*) software package. For this, ED and Z_{eff} maps can be sent automatically from the CT console to MM SIM*. Here, SPR can be directly calculated using the Bethe Bloch equation based on three different methods (Figure 1).¹⁵⁻¹⁷ Also, the SPR can be calculated for carbon (270.55 MeV), low energy protons (117 MeV) and high energy protons (200 MeV) (Figure 2).

While recent studies have described ED and Z_{eff} accuracies, an analysis on the link of these accuracies with SPR accuracy is lacking. We have, therefore, performed a literature search to identify the accuracies of ED and Z_{eff} on commercially available DECT and spectral CT systems for non-gated scans, to assess how these accuracies translate into SPR uncertainties.

Institution Simulation			•			?
Last Name	First Name	Middle Name	MRN	Date of Birth	Last Modified	
			777 SpectraRS91		2014-04-2511-14:22	
	Brain					
SpectralRSP1 Cardiac	Bruin Precise Cardiac Ition MMSHr "WorkFlow 02, C5, 0927 Siam "WorkFlow 02, C5, 09276, 2	10, 213, 103, 103, 103, 212, 103, 212, 103, 212, 103, 212, 103, 212, 103, 212, 103, 212, 103, 212, 103, 103, 103, 103, 103, 103, 103, 103	Spectral CT 7500 Details		2024-04-25 140222	
Spectral RSP1 Cerclac Cerclac Novalable Images NOV SERVICE SPR via SERVICE SPR via MA	Brain Precise Cardiac Store MM/Stim Tillpod/Sow (D. CT. 0907 MM/Stim Tillpod/Sow (D. CT. 0907) Stim Tillpod/Sow (D. CT. 0907) Stim Tillpod/Sow (D. CT. 0907) Stim Tillpod/Sow (D. CT. 0907) Stim Tillpod/Sow (D. CT. 0907)	V. 228. 103 V. 228. 103 V. 103. 502 V. Modality Shay bo Shay bo Shay bo Shay bo	Spectral CT 700 Denaits of Images e and Time		37	
Spectral ISSP1 Cardiac Source Images Available Images 1025 SEPRICE SPR via SERVICE SPR via too	Bado Precise Cordisc ation MMStm1Nex63xxx 00, C1, 002106, p	V. 213. 102 V. 213. 102 V. 101. VIII. VIII. 101. VIII. 101. VIII. 101. VIII. VIII. 101. VIII. 101. VIIII. 101. VIIII. 101. VIII. 101. VIIII. 101. VIIII. 101. VIIIII. 101. VIIIIIII. 101. VIIIIII. 101. VIIIIIIIIIIIIIII	Spectral CT 700 Details of Images e and Time scripten	990395 CT 219 203 2024-03-15 18 29 Z Effective , Spec	20144-314022 20144-314022 20149-314022	
Sonctra RSP1 Dardas	Bado Precise Circlisc Són MMKSim "Blooking (2), CT, 0917 Scim Blooking (2), CT, 0917 Scim Blooking (2), CT, 0917(5), C	Constant and a second and s	Spectral CT 7800 Databit of Dinapes e and Time scriptice Rimage Crease Options	990095 CT 219 203 2024-03-15.14.39 Z.Effective, Soec Q. Nes ⊕ Mo	201-64-23 140222 201-64-23 140222 201-64-23 140222 201-64-23 140222 201-64-23 140222 201-64-23 140222	
Spectru (KSP) Cardiac Available Images C25507165 550 min UN C25507165 550 min UN	Tanan Preside Cardide Martin Wassfow 82, CT, W20 Martin Wassfow 82, CT, W20 Schwitter Stew 93, CT, W20 President Stew 94, CT, W20	We 200 bits Image Set We 200 bits Anise Michael Set Modelity Start Dis Markie Markie Markie Start Dis Start Dis Sten Dis Sten Dis	Saketark (C 7 700) Decale of Unsupes a and Turne a stratistic 2 Unsupe Casalo Options a Unsupe Crasalo Options be Crasalos Marthood de Crasalos Marthood	990175 CT 229 903 2024-01-15 14 29 Z Ethective , Spec CN End: 4 Sec Set End: 4 Sec Set End: 4 Sec Set End: 4 Sec	22 29 20 20 20 20 20 20 20 20 20 20 20 20 20	
Sovetnakovi Cardiac Provinské Innegen Rozlakáč Koli do Uni Přistovi Koli do Uni	Ternin Preside Cardida: Malacen Teoretica (2), CE, (MP) Malacen Teoretica (2), CE, (MP)	Ye 20 Image Set Image Set Model Model Model Model Sorte 10	Spectral (C1 760) Details a sub Time sortpion 21 maps Creats Options to Crastion Parameter to Transfor University of Creates SPE II	990395 CT 219 903 2024-03-15 1439 2 Effective, Spec SPA Decte Sam SPA Decte Sam SPA Decte Sam SPA Decte Sam SPA Decte Sam SPA Decte Sam SPA Decte Sam	ант на ланова 2014 - 4 - 23 14 62 22 23 23 23 26 26 26 26 26 26 26 26 26 26 26 26 26	

Figure 1 MM SIM*, showing SPR methods

Instrume Instrume Model have Model Data of them Amennica Balencia 2014-64 2014-64 Standing Balencia 2014-64 2014-64 Standing Balencia Balencia 2014-64 Standing Balencia factoria Balencia factoria 2014-64	Last Monthed + 65 100444 66 923833 66 923833 66 923823 66 926323 66 926323 89 56283 39 562802 29 11:14:33 29 1460:22
bi 0 7441 2024 2024 2024 2024 2024 2024 2024 2	65 10 0444 65 93 38 30 65 99 38 30 65 99 39 29 65 98 58 58 38 56 58 8 38 56 58 8 38 56 58 8 29 11 14 33 25 14 60 22
معادمة، فكالمركب (2014) التركية المعادية، معادية، معاد معادية، معادية، معادي معادية، معادية، مع دى معادية، معادية، معادية، معادية، معادية، معادة، معادة، معادية، معادية، معادة، معادية، معادة، معادية، معادى، معادة، معادة	65 993833 65 993259 65 990223 66 985658 38 1658832 29 11:14:33 25 14:02:22
2015月、Adonuma biovom 2015年、20155757575757575757575757575757575757575	65 49:12:59 65 49:02:23 65 48:56:58 38 16:28:02 29 11:14:33 28 14:02:22
2015年、 051 年1 2016年、 2016年 2015年の中国の 2015年の日 2015年1日 2015年 2015年日日 2015年日日 2015年日 2015年年日 2015年年日 2015年年日 2015年年日 2015年年日 2015年年日 2015年年日 2015年年日 2015年年11月11月11月11月11月11月11月11月11月11月11月11月11	65 99:02:23 65 98:50:58 30 16:28:02 29 11:14:33 25 14:02:22
10733 Anorymous Usahom 000035 200-94	05 085058 30 162802 29 11:1433 25 146222
I 00001 2014-11-41 2014-45 Netra 65071 Enric (264-116-12) 2014-44 Perda: Precisi Cardio: 594-274(27.708) 2014-44	38 162812 29 11:14:33 25 14:02:22
antara KRAPI Banka zzz. Sakotra KRAPI 2014-04- antariz Persian Curdas: Sakotra C17300 2014-04- 2014-04- T	29 11:14:33 25 14:02:22
Anticine Cardinace Spectral CT 7:00 2024-64- Spectral CT 7:00 2024-64-	25 140222
Smulation •	
Namber of Images 503	
Shew SFR Image Create Options 📀 Yes 👁 Mo	
SPR Errape Creation Method SPR Bethe Bourque	
SPR Enrope Creation Method SPR Betre Bourgue SPR Enrope Creation Parameter Carbon	
598 Emaige Crastion Marchino 5000, Detro Bioarque 598 Errange Crastions Parameter Secondainy Image Unaid Dia Cavada S97 Errange Secondainy Image Unaid Dia Cavada S97 Errange	- F
39 R Image Cruzisis Michold SP2, Borte Sharevue SP2 Image Cruzisis Planmeter Secondar Image Ubed to Cruzis JP3 Tenes Processing Course	

Figure 2 MM SIM*, showing SPR parameters

Inaccuracy in proton stopping-power ratios (SPR) calculated from conventional HU-based CT images, and the associated proton beam range uncertainty, is a well-known limitation in realizing the full potential of proton therapy. This error has been documented in literature to be approximately 3%. Multi-energy CT has shown promise in reducing SPR inaccuracy and proton range uncertainty.

Key points

- Analyzes peer-reviewed publications that document the accuracy of electron density (ED) and effective atomic number (Z-effective) images from various multi-energy technologies.
- Compares the resulting SPR accuracy from these different inputs.
- The analysis shows a range of quantitative accuracy in the ED and Z-effective images across multi-energy CT technologies, and the resulting SPR values, with dual-layer detector-based spectral CT providing the most accurate results.
- Analysis demonstrates a spectral SPR uncertainty of less than 1% for detector-based spectral CT, which has been further confirmed in published literature.

Utilization of spectral CT is expected to increase for dose calculation in proton therapy because it can reduce this uncertainty to less than 1%.¹⁸



Methods

We have performed a literature search on Pubmed, using the following search query on June 11th 2024: ((photon counting) OR (dual energy)) + ((electron density) or (effective atomic number)) + (phantom).

We have included only peer-reviewed journal articles written in English that describe accuracy measurements of Z-effective and electron density (excluding non-commercial implementations) using the Gammex electron density CT Phantoms 467 (Gammex-RMI, Middleton, WI, USA).

We extracted the reported ED and Z_{eff} accuracies from the tables or figures in the articles and reported them in this white paper. Furthermore, we have used these ED and Z_{eff} accuracies to calculate the uncertainties that these inaccuracies imply on SPR values calculated with the Bethe-Bloch equation, using the polynomial fit proposed by Bourque (see appendix).¹⁹

Results

The Embase search resulted in 154 peer-reviewed journal articles. Out of these, four papers met our inclusion requirements, describing measurements on various DE and spectral CT solutions. One paper reported SPR uncertainty on photon counting CT, but because these results were not based on $Z_{\rm eff}$ and ED maps, this article was excluded.²⁰ Another paper used an incorrect $Z_{\rm eff}$ exponent for the calculation of the reference materials (3.3 and not the system assumed value of 2.94).²¹ Therefore, this article was also excluded. An overview of the included papers can be found in Table 1.

The reported ED and Z_{eff} accuracies based on non-gated CT scans, and the calculated implied SPR accuracies are shown in Table 2, Table 3, and Table 4, respectively. One article did not report ED results and therefore we were not able to calculate the implied SPR uncertainties.

The average absolute reported differences for Z_{eff} ranged from 1.6% (dual-layer) to 10.6% (split-filter) and for ED from 0.4% (dual-layer) to 3.2% (split-filter). The average absolute SPR uncertainties implied by these Z_{eff} and ED uncertainties ranged from 0.3% (dual-layer) to 5.0% (split-filter) (Table 5). Importantly, the maximum SPR uncertainty exceeded 2% for all techniques, apart from the dual-layer solution for which the maximum SPR uncertainty was found to be 0.9%.

Article	Year	Types of multi-energy CT tested	Reported
Goodsitt, et al. ²	2011	KV-switching DECT (KV-DECT)	Z-effective
Almeida, et al. ³	2017	Dual-Source DECT (DS-DECT) and Split-filter DECT (SF-DECT)	Z-effective, electron density
Hua, et al. ⁴	2018	Spectral-detector CT (spectral CT)	Z-effective, electron density

Table 1 Overview of articles included in the literature review analysis

Discussion

The results reported in Table 5 demonstrate that there is a large variation in the accuracy of ED and Z_{eff} between tubebased and detector-based spectral CT solutions. As a result, the implied SPR accuracies varied considerably between the solutions.

One article has reported on the accuracy of dual-layer spectral CT.¹⁴ In this study, Hua et al. compensated for the limited thickness (5 cm) of the Gammex 467 phantom. A reduced pitch and collimation were used to compensate for the assumptions on the body composition of the system's scatter correction algorithms.¹⁴ These adjusted scan settings are not needed in clinical scanning, because normal patient dimensions will fulfill the scatter assumptions.

For PCCT, direct SPR calculation is not possible in clinical practice as both Z_{eff} and ED maps cannot be generated on the currently available system. However, Hu et al.²⁰ reported on the SPR uncertainty while using virtual monoenergetic images (VMI) for an indirect SPR calculation. On the same phantom (Gammex 467), the VMI-based SPR uncertainty was found to be 1.27% for this PCCT.

From this literature study, we can conclude that Hua, et al.¹⁴ were able to report the highest Z_{eff} and ED accuracies. For this, the authors used a dual-layer spectral CT, which showed an implied SPR uncertainty of 0.3%, which is the most optimal implied SPR error in this review.

An SPR uncertainty of 0.3% is in line with both Faller, et al.²² and Longarino, et al.²³ who (using the Gammex 467 Phantom) also reported SPR uncertainties below 1% on the Philips IQon Spectral CT and Spectral CT 7500 (the latest generation detector-based spectral CT from Philips), respectively. The difference is that the calculated value in the current evaluation is a theoretical value, while the SPR uncertainties of Faller et al. (0.6%) and Longarino et al. (0.7%) are end-to-end measurements.

In conclusion, the detector-based spectral CT showed the lowest average absolute differences for both Z_{eff} and ED, which subsequently results in the lowest SPR uncertainty.

The detector-based spectral CT reproducibly delivers on the promise to provide SPR with an uncertainty below 1%, which results in the highest clinically available SPR accuracy.¹⁸

Clinical relevance

To precisely estimate the Bragg peak location, exact estimations of the stopping power of tissues are crucial. Historically, the proton stopping-power has been calculated through a stoichiometric lookup table that converts conventional HU directly into SPR values. This technique is associated with range uncertainties of around 3%, meaning that if the target is 10 cm in the body, the position of the Bragg peak has a location uncertainty of 3 mm in the direction of the beam.

The promise of spectral CT for proton therapy planning is an increase in SPR accuracy. Philips MMSIM offers functionality that can take ED and Z_{eff} results and generate SPR maps according to the Bethe-Bloch equation. However, any error in these ED and Z_{eff} measurements will directly translate into SPR uncertainties. It is therefore crucial that the quantitative accuracy of the spectral measurements is assured.

In this paper we provide an overview of the reported ED and Z_{eff} measured on the Gammex 467. It demonstrates that the accuracy of ED and Z_{eff} results in some DECT implementations are not sufficient to ensure errors under 1%. In some systems the average error even exceeds 3%, meaning that a conventional stoichiometric lookup table would yield better accuracy.

The review shows that detector-based spectral CT reproducibly assures both average and maximum SPR uncertainties under 1%. This conclusion was in line with a recent study by Longarino, et al.²³, that reported an average SPR uncertainty of 0.7% on the Gammex Phantom model 467 on Spectral CT 7500, which uses the latest generation spectral-detector system from Philips.

Gammex 467 insert	Dual-source (Siemens Flash)	Dual-source (Siemens Force)	Split-filter (Siemens Edge)	Dual-layer (Philips IQon)	kVp-switching (GE Discovery CT750 HD)
Lung (LN-300)	-4.2%	-5.3%	-20.4%	7.0%	0.0%
Liver	0.0%	-1.0%	-1.7%	0.7%	-1.1%
Bone (CB2 - 50% mineral)	0.7%	-0.3%	1.4%	-0.2%	-3.9%
Bone (B-200)	0.8%	1.7%	4.2%	0.2%	-3.5%
Cortical bone	-0.2%	0.3%	0.3%	-0.3%	-3.6%
Lung (LN-450)	-7.0%	-2.3%	-14.7%	5.8%	0.0%
Brain	-6.8%	-5.3%	-28.3%	1.2%	15.1%
Adipose	-4.9%	-3.6%	-22.3%	1.6%	-0.6%
Inner bone	-0.3%	1.3%	0.0%	0.0%	-3.4%
Bone (CB2 - 30% mineral)	-0.5%	-1.9%	-5.6%	0.0%	-3.4%
True water	3.1%	0.7%	-11.6%	-2.3%	5.5%
Breast	1.9%	-2.5%	-16.3%	0.2%	3.6%

Table 2 Percentage errors in the effective atomic number (ΔZ_{eff}) measured on dual-source and split-filter DECT as described by Almeida, et al.²⁴, on Dual-layer spectral CT as described by Hua et al.¹⁴, and on kVp-switching DECT as described by Goodsitt, et al.¹²

Gammex 467 insert	Dual-source (Siemens Flash)	Dual-source (Siemens Force)	Split-filter (Siemens Edge)	Dual-layer (Philips IQon)	kVp-switching (GE Discovery CT750 HD)
Lung (LN-300)	0.4%	1.2%	15.3%	0.8%	n/a
Liver	-0.4%	-0.1%	0.7%	0.6%	n/a
Bone (CB2 - 50% mineral)	-0.5%	0.0%	-0.3%	0.5%	n/a
Bone (B-200)	-0.3%	-0.3%	-0.8%	0.1%	n/a
Cortical bone	0.2%	0.0%	-0.5%	-0.1%	n/a
Lung (LN-450)	1.0%	0.2%	7.2%	0.5%	n/a
Brain	1.3%	1.0%	3.1%	1.1%	n/a
Adipose	0.7%	0.5%	2.4%	0.0%	n/a
Inner bone	0.4%	-0.4%	0.6%	0.1%	n/a
Bone (CB2 - 30% mineral)	0.3%	0.5%	2.1%	0.1%	n/a
True water	0.6%	0.6%	3.2%	0.1%	n/a
Breast	0.3%	0.0%	2.1%	0.0%	n/a

Table 3 Percentage errors in electron density (Δ ED) measured on dual-source and split-filter DECT as described by Almeida, et al.²⁴, on detector-based spectral CT as described by Hua et al.¹⁴, and on kVp-switching DECT as described by Goodsitt, et al.¹²

Gammex 467 insert	Dual-source (Siemens Flash)	Dual-source (Siemens Force)	Split-filter (Siemens Edge)	Dual-layer (Philips IQon)	kVp-switching (GE Discovery CT750 HD)
Lung (LN-300)	0.7%	1.6%	18.5%	-2.1%	n/a
Liver	-0.4%	0.0%	0.8%	1.3%	n/a
Bone (CB2 - 50% mineral)	-0.6%	0.1%	-0.6%	0.3%	n/a
Bone (B-200)	-0.4%	-0.5%	-1.1%	0.1%	n/a
Cortical bone	0.2%	-0.1%	-0.6%	0.0%	n/a
Lung (LN-450)	1.6%	0.4%	9.0%	-2.4%	n/a
Brain	2.5%	2.0%	9.4%	-2.1%	n/a
Adipose	1.6%	1.2%	7.1%	0.2%	n/a
Inner bone	0.4%	-0.5%	0.6%	0.1%	n/a
Bone (CB2 - 30% mineral)	0.3%	0.7%	2.6%	0.0%	n/a
True water	0.4%	0.6%	4.5%	-0.4%	n/a
Breast	0.5%	0.3%	4.9%	-0.1%	n/a

Table 4 Implied percentage errors in stopping-power ratio (Δ SPR) based on effective atomic number and electron density data for dual-source and split-filter DECT as described by Almeida, et al.²⁴, on detector-based spectral CT as described by Hua et al.¹⁴, and on kVp-switching DECT as described by Goodsitt, et al.¹² For kVp-switching, Δ SPR could not be calculated, due to missing ED information in the described paper.



		Stopping-power ratio uncertainty			
CT technique category	e CT system	Minimum	Maximum	Absolute mean	
Source-based	Dual-source (Siemens Flash)	-0.6%	2.5%	0.8%	
	Dual-source (Siemens Force)	-0.5%	2.0%	0.7%	
	Split-filter (Siemens Edge)	-1.1%	18.5%	5.0%	
	kVp-switching (GE Discovery CT750 HD)	n/a	n/a	n/a	
Detector-based	Dual-layer (Philips Spectral CT)	-0.3%	0.9%	0.3%	

Table 5 Implied percentage errors in sstopping-power ratio (Δ SPR) over all Gammex 467 inserts (minimum, maximum, and absolute mean values) based on effective atomic number and electron density data for dual-source and split-filter DECT as described by Almeida, et al.24, on detector-based spectral CT as described by Hua et al.¹⁴, and on kVp-switching DECT as described by Goodsitt, et al.¹² For kVp-switching, Δ SPR could not be calculated, due to missing ED information in the described paper.

References

- 1. Bragg WH, Kleeman R. LXXIV. On the ionization curves of radium. The London, Edinburgh, and Dublin Philosophical Magazine and Journal of Science. 1904;8(48):726-738. doi:10.1080/14786440409463246
- 2. Durante M, Orecchia R, Loeffler JS. Charged-particle therapy in cancer: clinical uses and future perspectives. Nat Rev Clin Oncol. 2017;14(8):483-495. doi:10.1038/nrclinonc.2017.30
- 3. Bär E, Lalonde A, Royle G, Lu HM, Bouchard H. The potential of dual-energy CT to reduce proton beam range uncertainties. Med Phys. 2017;44(6):2332-2344. doi:10.1002/mp.12215
- Yang M, Virshup G, Clayton J, Zhu XR, Mohan R, Dong L. Theoretical variance analysis of single- and dual-energy computed tomography methods for calculating proton stopping power ratios of biological tissues. Phys Med Biol. 2010;55(5):1343-1362. doi:10.1088/0031-9155/55/5/006
- 5. Hudobivnik N, Schwarz F, Johnson T, et al. Comparison of proton therapy treatment planning for head tumors with a pencil beam algorithm on dual and single energy CT images. Med Phys. 2016;43(1):495-504. doi:10.1118/1.4939106
- 6. Zhu J, Penfold SN. Dosimetric comparison of stopping power calibration with dual-energy CT and single-energy CT in proton therapy treatment planning. Med Phys. 2016;43(6Part1):2845-2854. doi:10.1118/1.4948683
- 7. Bär E, Lalonde A, Zhang R, et al. Experimental validation of two dual-energy CT methods for proton therapy using heterogeneous tissue samples. Med Phys. 2018;45(1):48-59. doi:10.1002/mp.12666
- 8. Möhler C, Russ T, Wohlfahrt P, et al. Experimental verification of stopping-power prediction from single- and dual-energy computed tomography in biological tissues. Phys Med Biol. 2018;63(2):025001. doi:10.1088/1361-6560/aaa1c9
- Wohlfahrt P, Möhler C, Richter C, Greilich S. Evaluation of Stopping-Power Prediction by Dual- and Single-Energy Computed Tomography in an Anthropomorphic Ground-Truth Phantom. International Journal of Radiation Oncology*Biology*Physics. 2018;100(1):244-253. doi:10.1016/j.ijrobp.2017.09.025
- Xie Y, Ainsley C, Yin L, et al. Ex vivo validation of a stoichiometric dual energy CT proton stopping power ratio calibration. Phys Med Biol. 2018;63(5):055016. doi:10.1088/1361-6560/aaae91
- 11. Paganetti H. Range uncertainties in proton therapy and the role of Monte Carlo simulations. Phys Med Biol. 2012;57(11):R99-R117. doi:10.1088/0031-9155/57/11/R99
- 12. Goodsitt MM, Christodoulou EG, Larson SC. Accuracies of the synthesized monochromatic CT numbers and effective atomic numbers obtained with a rapid kVp switching dual energy CT scanner. Med Phys. 2011;38(4):2222-2232. doi:10.1118/1.3567509
- 13. Almeida IP, Schyns LEJR, Ollers MC. Dual-energy CT quantitative imaging : a comparison study between twin-beam and dual-source CT scanners. Published online 2017:171-179.
- Hua C, Shapira N, Merchant TE, Klahr P, Yagil Y. Accuracy of electron density, effective atomic number, and iodine concentration determination with a dual-layer dual-energy computed tomography system. Med Phys. 2018;45(6):2486-2497. doi:10.1002/mp.12903
- 15. Saito M, Sagara S. Simplified derivation of stopping power ratio in the human body from dual-energy <scp>CT</scp> data. Med Phys. 2017;44(8):4179-4187. doi:10.1002/mp.12386
- Yang M, Virshup G, Clayton J, Zhu XR, Mohan R, Dong L. Theoretical variance analysis of single- and dual-energy computed tomography methods for calculating proton stopping power ratios of biological tissues. Phys Med Biol. 2010;55(5):1343-1362. doi:10.1088/0031-9155/55/5/006
- 17. Bourque AE, Carrier JF, Bouchard H. A stoichiometric calibration method for dual energy computed tomography. Phys Med Biol. 2014;59(8):2059-2088. doi:10.1088/0031-9155/59/8/2059
- Kruis MF. Improving radiation physics, tumor visualisation, and treatment quantification in radiotherapy with spectral or dual-energy CT. J Appl Clin Med Phys. 2022;23(1). doi:10.1002/acm2.13468
- 19. Bourque AE, Carrier JF, Bouchard H. A stoichiometric calibration method for dual energy computed tomography. Phys Med Biol. 2014;59(8):2059-2088. doi:10.1088/0031-9155/59/8/2059
- 20. Hu G, Niepel K, Risch F, et al. Assessment of quantitative information for radiation therapy at a first-generation clinical photon-counting computed tomography scanner. Front Oncol. 2022;12. doi:10.3389/fonc.2022.970299
- Ohira S, Washio H, Yagi M, et al. Estimation of electron density, effective atomic number and stopping power ratio using dual-layer computed tomography for radiotherapy treatment planning. Physica Medica. 2018;56:34-40. doi:10.1016/j. ejmp.2018.11.008
- 22. Faller FK, Mein S, Ackermann B, Debus J, Stiller W, Mairani A. Pre-clinical evaluation of dual-layer spectral computed tomography-based stopping power prediction for particle therapy planning at the Heidelberg Ion Beam Therapy Center. Phys Med Biol. 2020;65(9):095007. doi:10.1088/1361-6560/ab735e
- 23. Longarino FK, Kowalewski A, Tessonnier T, et al. Potential of a Second-Generation Dual-Layer Spectral CT for Dose Calculation in Particle Therapy Treatment Planning. Front Oncol. 2022;12. doi:10.3389/fonc.2022.853495
- 24. Almeida IP, Schyns LEJR, Öllers MC, et al. Dual-energy CT quantitative imaging: a comparison study between twin-beam and dual-source CT scanners. Med Phys. 2017;44(1):171-179. doi:10.1002/mp.12000

Appendix

$$SPR = \frac{\rho_e}{\rho_{e,w}} \bullet \frac{\ln\left(\frac{2m_e c^2 \beta^2}{I_m (1 - \beta^2)}\right) - \beta^2}{\ln\left(\frac{2m_e c^2 \beta^2}{I_w (1 - \beta^2)}\right) - \beta^2}$$

As stated in the introduction, SPR for a medium relative to water is approximated following the Bethe equation: with $\rho_{e,w}$ the electron densities of the medium (input from the spectral reconstruction) and water respectively. Other parameters are:

• β is the velocity of the particle beam relative to the speed of light c (dependent on the energy of the proton beam):

$$3^{2} = 1 - \frac{1}{\left(1 + \frac{E}{m_{o}c^{2}}\right)^{2}}$$

- m_ec² the electron mass (511 keV)
- $m_{\rho}c^{2}$ the proton mass (938 MeV)
- *E* is the proton beam kinetic energy. In this evaluation a relatively low energy of 117.0 MeV is used.
- *I_m* is the mean excitation energy of the medium (computed with the effective atomic number from the spectral reconstruction)
- I_w the mean excitation energy of water (78.73 eV).

$$I_{m} = \begin{cases} a_{1,1}Z + a_{1,0} & Z < 6.26 \\ a_{2,5}Z^{5} + a_{2,4}Z^{4} + a_{2,3}Z^{3} + a_{2,2}Z^{2} + a_{2,1}Z + a_{2,0} & 6.26 \le Z \le 13.52 \\ a_{3,1}Z + a_{3,0} & 13.52 < Z \end{cases}$$

To compute the excitation energy of the medium based on the effective atomic, the polynomial fit proposed by Bourque is used⁷: In these equations, Z is the effective atomic number. The values of the constants are: $a_{1,1}=14.007762$, $a_{1,0}=-24.414214$, $a_{2,5}=-0.005342$, $a_{2,4}=0.207079$, $a_{2,3}=-2.589844$, $a_{2,2}=8.339473$, $a_{2,1}=51.895887$, $a_{2,0}=-219.722173$, $a_{3,1}=11.794847$, $a_{3,0}=-47.707141$.

The error in the stopping power is computed based on the errors in the electron density $\rho_{-}e$ and effective atomic number Z provided:

$$SPR_{error} = \frac{SPR(\rho_e(1 + \rho_{e,error}), Z(1 + Z_{error})) - SPR(\rho_e Z)}{SPR(\rho_e Z)}$$



© 2024 Koninklijke Philips N.V. All rights reserved. Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

How to reach us Please visit www.philips.com

4522 991 87291 * DEC 2024