
A patient is lying in a hospital bed, wearing a blue gown. They are holding a red and yellow stuffed duck. A pulse oximeter is attached to the duck's beak. The patient's hand is visible, holding the duck. The background is blurred, showing the hospital bed and some medical equipment.

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Philips

A young child with dark skin is lying in a hospital bed, wearing a light blue hospital gown. A pulse oximeter is attached to their right index finger, connected by a blue tube. The child is looking towards the camera with a neutral expression. The background is a soft-focus white hospital pillow and blanket.

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Society for Technology in Anesthesia (STA)

- Correlation between Subjective and Objective Pigmentation Measurement Methods for Clinical Research (Corrélation entre méthodes subjectives et objectives de mesure de la pigmentation en recherche clinique)

American Society of Anesthesiologists (ASA)

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World Congress of Anesthesiologists (WCA)

- Are intra-finger pad pigmentation profiles as impactful as pigmentation variation between subjects? (Les profils de pigmentation de la pulpe digitale ont-ils un impact comparable à la variation de pigmentation entre sujets ?)

American Association of Critical Care Nurses (AACN NTI)

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Correlation between Subjective and Objective Pigmentation Measurement Methods for Clinical Research


Co-auteurs : Wim Verkruysse, PhD¹, Michael B. Jaffe, PhD². ¹Philips Hospital Patient Monitoring (HPM), Eindhoven, Pays-Bas ; ²Philips HPM, Cambridge, MA, États-Unis.

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
Are intra-finger pad pigmentation profiles as impactful as pigmentation variation between subjects?

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
Letter Regarding the Potential Improper Use of the Individual Typology Angle in the Context of Disparate Bias in Pulse Oximetry

To the Editor

There has been significant interest recently to develop a better understanding of potential causes underlying an apparent racial bias in pulse oximeters as reported by Spalding et al¹ and others. One area of focus is to better understand how skin tone affects disparate bias in pulse oximeters. The US Food and Drug Administration (FDA) released a discussion paper in late 2023 that proposes expanding requirements for pulse oximeter verification studies, including an increase in the number of included participants and, importantly, the breadth and distribution of their skin pigmentation.² The discussion paper proposes that skin pigmentation should be assessed using both a subjective scale, as well as quantified with an objective scale using the individual typology angle (ITA). The purpose of this letter is to point out some confusion on how the ITA formula is reported in the literature. Novice users computing ITA for the first time may not easily notice a subtle error in the formula propagated in many papers. The resulting error in the actual ITA value is small for medium pigmented skin but is very large for dark and light skin.

The ITA is a metric for characterizing the spectrum of melanin pigmentation in human skin and was originally introduced by Chardon.³ The ITA is typically determined with colorimetry,⁴ measuring color values

L^* , a^* and b^* on skin, where L^* is the perceived lightness, from black (value 0) to white (value 100) and the parameters a^* and b^* represent the chromaticity; a^* represents the balance between magenta (positive values) and green (negative values). Similarly, b^* represents the balance between yellow (positive values) and blue (negative values). The L^* and b^* values are then used to compute the ITA with the following equation:

$$ITA = (180/\pi) \arctan ((L^* - 50)/b^*) \quad (1)$$

Since many people involved in pulse oximetry research, as well as device validation studies (industry, hospitals, governmental institutions, etc) will be new to computing ITA values, we find it useful to note herein that confusion may arise from certain publications, including some from leading authors,^{5,6} and very recent ones.⁷ Namely, in a nonsystematic nor exhaustive search for incorrect formulas in the literature we found at least 33 publications from more than 13 different research groups in which the definition of ITA is printed as follows:

$$ITA = (180/\pi) (\arctan (L^* - 50)/b^*), \quad (2)$$

or variants thereof. The difference between equations 1 and 2 is in the different bracketing which leads to significantly different, incorrect, ITA values for the same set of L^* , b^* values. This is shown in Figure where we use 6 sets of values for L^* , b^* from⁸ to graphically show ITA in Figure A. We also computed ITA using equations 1 and 2. Equation 2 leads to a small range of angles (approximately -10° to $+10^\circ$) for the human skin spectrum, and a large error for light and dark skin, shown in Figure B.


Figure A. Six sets of colorimetric L^* and b^* values (black circles) were estimated from the skin color volume described in showing ITA graphically. The L^* and b^* sets were used to compute angles using equations 1 and 2, giving an error shown in (B). ITA indicates individual typology angle.

Letter Regarding the Potential Improper Use of the Individual Typology Angle in the Context of Disparate Bias in Pulse Oximetry

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https://journals.lww.com/anesthesia-analgesia/citation/2024/05000/letter_regarding_the_potential_improper_use_of_the35.aspx

3. Brass P, Holmich M, Kolditz L, Schick G, Smith AF. Ultrasound guidance versus anatomical landmarks for internal jugular vein catheterization. *Cochrane Database Syst Rev*. 2015;1:CD00962.
4. Wu S, Ling Q, Cao L, Wang L, Xu M, Zeng W. Real-time two-dimensional ultrasound guidance for central venous catheterization: a meta-analysis. *Anesthesiology*. 2013;118:361-375.
5. Evans CS, Mangione MT, Chevit TM, Dinsien JW, Ryker RJ. Eliminating arterial injury during central venous catheterization using manometry. *Anesth Analg*. 2009;109:130-134.
6. Calvoche JA, Rodriguez MS, Troncher A, Klinek M, Stokker RJ, Lissafre E. Incidence of mechanical complications of central venous catheterization using landmark technique do not try more than 3 times. *J Intensive Care Med*. 2016;31:397-402.

DOI: 10.1016/j.ana.2024.06.000

Variability of Printed Monk Skin Tone Scales May Cause Misclassification of Clinical Study Participants: Caveats On Printing

To the Editor

To help better understand health inequities such as racial bias in pulse oximetry,¹ the 10-shade Monk Skin Tone (MST) scale has been recently proposed as an inclusive enrollment tool² to diversify study participant enrollment in pulse oximeter deaeration studies required for regulatory submissions in the United States.³

This letter addresses the challenge of consistency of the printed MST color scale in use. The MST is defined for displays (RGB coordinates) only while a printed version is needed in the actual practice of laboratory studies and clinical studies. Without clear printing guidelines, different color prints may thus lead to different MST classification of people with the same skin. We seek to raise this issue in this matter.


and visualized by mapping them on an L*, b* graph, as is common in describing skin pigmentation.⁴

A strict separation of OP and PP does not exist but price (cut-off about 25\$), print resolution (cut-off about 300 DPI), color quality, calibration frequency (up to yearly for OPs and days for PPs), and operation by skilled personnel are important factors.

OP printing (Figure A) shows large overlaps between different MST grades at the light and dark shades of the scale. Also, the L* and b* values for the colors rendered by the OPs hardly align with the skin tone volume as defined by Del Bino et al.⁵ mostly because the b* values are too low (prints are too grayish) and the L* values too high (prints are not dark enough). The exact cause for this large variation between printers, even of the same model often, may be poor maintenance, but also the fact that they were not calibrated recently. Also, we have observed that unprinted glossy paper is typically somewhat yellower than matte paper, and consequently also tends to have MST colors with slightly larger b* values.

PPs (Figure B) shows a much better separation of MST grades than OPs. More importantly, the rendered colors are more saturated and darker; they now align much better with the skin tone volume. The better alignment is emphasized in Figure C, where the centroids of PP-printed MST grades are shown along with real measured L* and b* values for 21 subjects⁶ as well as selected data measured from a factory-printed skin tone scale (Pantone Skin Tone guide).

This study was not aimed at strictly defining instructions for printing parameters such as calibration frequency, printer type and model, paper type, and weight. In fact, such detailed instructions may be impractical and make it virtually impossible to find such conditions in all areas of the world where MST

Variability of Printed Monk Skin Tone Scales May Cause Misclassification of Clinical Study Participants: Caveats On Printing

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https://journals.lww.com/anesthesia-analgia/citation/2024/06000/variability_of_printed_monk_skin_tone_scales_may.36.aspx

1. FDA, Center for Devices and Radiological Health. Pulse oximeters. Approach for improving the performance evaluation of pulse oximeter devices taking into consideration skin pigmentation, race and ethnicity. Discussion paper and request for feedback. Consulted le 8 décembre 2023; <https://www.fda.gov/media/172005/download>
2. Del Bino S, Sok J, Bessac E, Bernard F. Relationship between skin response to ultraviolet exposure and skin color type. *Pigment Cell Res*. 2005;19:608-614.


Challenges of Subjective Skin Color Scales: The Case for the Use of Objective Pigmentation Measurement Methods in Regulatory Pulse Oximetry Studies

To the Editor

Addressing and understanding skin color-associated bias in pulse oximetry is a patient safety and health equity priority.¹ Updated regulatory requirements and standards are needed to improve population diversity in device validation studies and provide greater assurance concerning performance. Skin pigmentation measurement methods (PMMs) are either subjective or objective, typically using color-matching scales and light-based measurement devices, respectively. The use of new color-matching scales² are being considered for inclusion in a revision to the pulse oximeter International Organisation for Standardization (ISO) standard³ and FDA guidance document⁴ to better ensure the validation. However, we argue that ensuring the diversity of skin pigmentation in device research, development, and regulatory studies must include objective PMM given their high precision and accuracy. We support our opinion by presenting our preliminary results of 2 studies. The first study aimed to quantify operator imprecision of subjective and objective PMM. The second study aimed to quantify potential inaccuracies stemming from the subjective and objective skin tone scales.

CM700d measured CIE color values L^* and b^* (where L^* is the perceived lightness and b^* represents the balance between yellow and blue), from which individual Typology Angle (ITA) was calculated. Color matching was repeated 3 times for each subject, during 3 visits with at least a week between visits. Individual operator scores for each subject, anatomical location, and visit (a total of $12 \times 4 \times 3 = 144$) are plotted against the mean score of 3 operators in the Figure A-D. Different operators typically judged the matching color differently which can lead to an inherent imprecision that spans approximately 20% of the study skin tone range for both PST and VL (± 2 SD, orange bars in Figure A, B). Operator bias can be observed: operators 1 and 3 typically assessed the lighter skin tones as darker and lighter than the mean, respectively (Figure A, B). In contrast, with objective PMM (Figure C, D), the precision is much higher (SD < 5%) and operator bias is negligible.

A second study quantified the potential inaccuracy stemming from different printed specimens of a relatively new scale: the Monk Skin Tone scale (MST).^{1,2} Sixty-eight printed MST prints were acquired from professionally operated high-end printers. "High-end" is loosely defined as costing >\$1 USD and 30k USD for inkjet and laser printers, respectively. Printer operators were asked to print the portable document (pdf file) with embedded hexadecimal color values for each MST shade. We measured CIE color values for each shade with the CM700d and plotted these in the L^* - b^* projection (Figure E), consistent with the definition of the pigmentation metric ITA.⁵ As an illustration, color points for MST shade C are shown for two different printing methods. The error of colorimetric error (MSE) is shown.

Challenges of Subjective Skin Color Scales: The Case for the Use of Objective Pigmentation Measurement Methods in Regulatory Pulse Oximetry Studies

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Citation :

Wim Verkruijsse, Michael B Jaffe, Michael Lipnick, Charifa Zemouri. Challenges of Subjective Skin Color Scales: The Case for the Use of Objective Pigmentation Measurement Methods in Regulatory Pulse Oximetry Studies. *Anesth Analg*. 2024;139(2):e15-E17.

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https://journals.lww.com/anesthesia-analgia/citation/9900/challenges_of_subjective_skin_color_scales_the.823.aspx

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Objective and Perceived Skin Color; Consequences for the Use of Skin-Tone Scales.


To the Editor

Skin color associated bias in pulse oximetry has been a patient safety and health equity priority 1 for some time. Regulatory requirements and standards are being updated to include the use of different types of skin pigmentation measurement methods (PMM) to improve population diversity in device validation studies. In a recent letter² we argued that objective PMMs should be included due to being more precise, more accurate and less subject to operator bias than subjective PMM, e.g. Monk Skin Tone (MST) scale 3. Recent insights, based on data presented herein, reinforce that opinion but also support the use of a subjective PMM as a valuable tool, complementary to an objective PMM. Our data, supported by other data 4, shows that skin color as perceived by humans (Pac) and assessed with a certain class of skin-tone

complementary types of skin tone assessment methodologies.

Osc is assessed with established quantitative metrics used by scientists and cosmetologists in various industries (e.g. textile). A typical use of Osc is to minimize differences between skin and a cosmetics product or a textile. Data from published skin tone references are seen to overlap very nicely with skin data from literature and our own studies (figure A) with the range of L* of Osc between 25 to 75. Notably, Osc data from printed MST swatches as well as the defined MST values have only 3-4 shades that fall within this range.

To assess the differences between Osc and Pac, three skintone scales, PST, Fitzpatrick (FP) and Von Luschan (VL) were compared to Osc as measured with a Konica Minolta CM700d. The Osc for each of the swatches of these scales allows us to quantify the Pac (the swatches selected by assessors). The lightest and darkest swatches of these scales (Figure C) illustrate the differences in L* ranges. We

Objective and perceived skin color; consequences for the use of skin-tone scales

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Citation :

Wim Verkruysse, Monk, E Jr, Michael B Jaffe. Objective and Perceived Skin Color; Consequences for the Use of Skin-Tone Scales.

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Copyright © 2024