**PHILIPS**


**AVENT**

Single/Double  
Electric Breast  
Pump

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


## Introduction

Congratulations on your purchase and welcome to Philips Avent! Unleash your freedom with the Philips Avent Hands-free Electric Breast Pump. The pump that brings the best of both worlds; unique soft comfort combined with an effective\* pump. The unique SkinSense silicone breast shields softly shape around your breast, using your body's natural warmth for a comfortable fit. The ultra-light collection cups are fully transparent for easy nipple placement and milk flow visibility.

Did you know that breast pumps can play a crucial support role in your breastfeeding experience? Breast pumps can enhance milk volume and with a stash of expressed milk in the fridge, partners and caregivers can actively participate in caring for the little one, affording you precious moments of self-care or a smooth return to work. If you are breastfeeding, it is advisable (unless recommended otherwise by your healthcare professional) to wait until your milk supply and breastfeeding schedule are established (normally at least 2 to 4 weeks after giving birth) before you start expressing.

\* 'Effective' is related to the technical performance of the product.

## General description (Fig. 1)

**i** Note: The numbers mentioned below refer to the numbers in the figure on the front foldout page of this user manual.

- |                                |                                |
|--------------------------------|--------------------------------|
| 1 Front cap                    | 9 Collection cup (assembled)   |
| 2 Membrane                     | 10 Motor unit                  |
| 3 Cup                          | 11 Clip (detachable)           |
| 4 Pouring spout                | 12 Silicone tube               |
| 5 Valve                        | 13 Carrying cord               |
| 6 Breast shield (21 mm, 24 mm) | 14 USB cable                   |
| 7 Insert (19 mm)               | 15 Double breast pump (SCF534) |
| 8 Breast shield (dust) cover   | 16 Tube coupler                |

## Intended use and indications for use

The Philips Avent Hands-free Single/Double Electric Breast Pump is intended to express and collect milk from the breast of a lactating woman, and indicated to alleviate the symptoms of engorgement. The device is intended for a single user.

## Important safety information

Only use the product for its intended purpose. Read this important information carefully before you use the product and its batteries and accessories, and save it for future reference. Misuse can lead to hazards or serious injuries.

This user manual can also be found online at: [www.philips.com/support](http://www.philips.com/support).

This product is operated by the lactating woman herself.

You can use this breast pump in an everyday environment.

### **Contraindication**

**Never use the breast pump while you are pregnant, as pumping can induce labor.**

### **Warning**

#### **Warnings to avoid choking, strangulation and injury:**

- Never use the breast pump while driving to avoid accidents due to reduced attention.
- Children and pets should be supervised to ensure that they do not play with the motor unit, packaging materials or accessories.
- Make sure that you use the silicone tube, carrying cord, and USB cable carefully to avoid strangulation.
- This breast pump is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge. Such persons can only use this breast pump if they are supervised by or have received proper instruction concerning use of the appliance by a person responsible for their safety.

#### **Warning to prevent transmitting infectious diseases:**

- If you are a mother who is infected with Hepatitis B, Hepatitis C or Human Immunodeficiency Virus (HIV), expressing breast milk will not reduce or remove the risk of transmitting the virus to your baby through your breast milk.

#### **Warnings to avoid electric shock:**

- If a power adapter (not included) is used for charging, only use a certified double-insulated USB power adapter. Check if the symbol for double insulation is on the adapter: . Use a power source (an adapter or a power bank) with a USB port, of 5 V direct current/1.8 A or 9 W for your breast pump.
- When using a USB power adapter (not included) to charge or operate the breast pump: inspect the adapter for signs of damage before charging. Do not use if the adapter is damaged, or outer casing becomes loose or is separating. Separation of the casing can result in electric shock.

- Inspect the breast pump for signs of damage before each use or if it has been dropped. Do not use the breast pump if it is damaged or does not work properly, or if it has been submerged in water.
- Always ensure that it is easy to disconnect the device from the mains in case of malfunction. The external power supply is considered the disconnection device.

**Warnings to avoid contamination and to ensure hygiene:**

- For hygienic reasons, the breast pump is only intended for repeated use by a single user.
- Clean and disinfect all parts that come into contact with the breast and breast milk before first use and after each use. If you use inserts and/or breast shields of different sizes during the pumping session, make sure to clean and disinfect all such parts.
- Store and transport the collection cups using the breast shield (dust) covers.

**Warnings to avoid breast and nipple problems including pain:**

- Do not use the breast pump if the membrane appears to be damaged or broken. See 'Ordering accessories and spare parts' for information on how to obtain replacement parts.
- Only use accessories and parts recommended by Philips Avent.
- No modification of the breast pump is allowed. If you do so, your warranty becomes invalid.
- There are no user serviceable parts inside the motor unit of the breast pump.
- Never use the breast pump while you are sleepy or drowsy to avoid lack of attention during use.
- If the vacuum is uncomfortable or causes pain, switch off the motor unit and remove the collection cups from your breast.
- Contact your healthcare professional or breastfeeding specialist if you express no or insufficient milk or if your breast is painful, hot, swollen, or red.

**Warnings to prevent damage and malfunction of the breast pump:**

- Portable radio frequency (RF) communications equipment (including mobile phones and peripherals such as antenna cables and external antennas) should not be used closer than 30 cm (12 in) to any part of the breast pump, including external power supply. This could have a negative effect on the performance of the breast pump.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation.

### **Warning to avoid increased electromagnetic emissions or decreased electromagnetic immunity of this device, resulting in improper operation:**

- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

### **Battery safety instructions**

Only use the product for its intended purpose. Read this important information carefully before you use the product and its batteries and accessories, and save it for future reference. Misuse can lead to hazards or serious injuries.

Accessories supplied may vary for different products.

- Charge, use and store the product at the temperature ranges indicated in the technical specifications chapter in this user manual.
- Keep product and batteries away from fire and do not expose them to direct sunlight or high temperatures.
- If the product becomes abnormally hot or smelly, changes color or if charging takes longer than usual, stop using and charging the product and contact Philips.
- Do not place products and their batteries in microwave ovens or on induction cookers.
- Do not open, modify, pierce, damage or dismantle the product or battery to prevent batteries from heating up or releasing toxic or hazardous substances. Do not short-circuit, overcharge or reverse charge batteries.
- Be aware of the risk of terminals of the battery-operated appliance or battery being short-circuited by metal objects.

### **⚠ Caution**

#### **Cautions to prevent damage to and malfunctioning of the breast pump:**

- Prevent the USB power adapter (not included), USB cable, and motor unit from coming into contact with water.
- Keep all parts away from heated surfaces to avoid overheating and deformation of these parts.


### **Before first use**

Disassemble, clean and disinfect all parts that come into contact with the breast and breast milk before you use the breast pump for the first time and after each use.


Fully charge the device before you use it for the first time or after a long period of non-use.

## Cleaning and disinfecting

- Disassemble, clean and disinfect the following parts before first use and after every use: membrane, cup, valve, breast shield, insert (if used), breast shield (dust) cover.

- Clean the front cap, silicone tube, tube coupler (included with the double breast pump), carrying cord, USB cable, and motor unit (including the clip) when needed. Wipe with a clean, damp cloth with water and mild detergent.

### Step 1: Disassembling

Completely disassemble all parts that come into contact with the breast and breast milk. Make sure that you:

- Remove the valve from the cup.
- Separate the membrane and front cap.
- Remove the insert (if used) from the breast shield.

### Step 2: Cleaning

**⚠ Warning: Do not use antibacterial or abrasive cleaning agents to clean the breast pump parts as this may cause damage.**

**⚠ Caution: Never put the motor unit in water or in a dishwasher, as this causes permanent damage to it.**

Clean the following parts before first use and after every use: membrane, cup, valve, breast shield, insert (if used), breast shield (dust) cover. You can clean these parts manually or in the dishwasher.

## Step 2A: Manual cleaning

### Supplies needed:

- Mild dishwashing liquid
- Drinking-quality water
- Soft, clean brush
- Clean tea towel or drying rack
- Clean sink or bowl

**⚠ Caution:** Be careful when you remove and clean the valve. If it gets damaged, your breast pump does not function properly. To clean the valve, rub it gently between your fingers in warm water with some dishwashing liquid. Do not insert objects into the valve, as this may cause damage.

- 1 Rinse the parts under a running tap with lukewarm water.
- 2 Soak the parts for 5 minutes in warm water with some mild dishwashing liquid.
- 3 Clean the parts with a cleaning brush in warm water with some mild dishwashing liquid.
- 4 Thoroughly rinse the parts under a running tap with drinking-quality water.
- 5 Leave the parts to air-dry on a clean tea towel or drying rack.

## Step 2B: Cleaning in the dishwasher

### Supplies needed:

- Mild dishwashing liquid or a dishwashing tablet
- Drinking-quality water

**ⓘ Note:** Food colorings may discolor parts.

- 1 Place all parts on the top rack of the dishwasher. Use a dishwasher basket for the smallest parts.
- 2 Put detergent or a dishwashing tablet in the machine and run a standard program.
- 3 Leave the parts to air-dry on a clean tea towel or drying rack.

## Step 3: Disinfecting

### Supplies needed:

- Household pot
- Drinking-quality water

**⚠ Caution:** During disinfection with boiling water, prevent the parts from touching the sides and bottom of the pot. This can cause irreversible product deformation or damage that Philips cannot be held liable for.

- 1 Fill a household pot with enough water to cover these parts and place all the parts in the pot. Let the water boil for 5 minutes. Make sure the parts do not touch the sides and bottom of the pot.
- 2 Allow the water to cool down.

- 3 Gently remove the parts from the water. Leave the parts to air-dry on a clean tea towel or drying rack.
- 4 Store the dry items safely in a clean, protected area until needed. Ensure that the clean parts have air-dried completely before storing them to help prevent germs and mold from growing.

## Using the breast pump

**i** Note: The figure numbers mentioned below refer to the numbers of the figures on the back foldout pages of this user manual.

## Charging the breast pump

Fully charge the battery before you use the breast pump for the first time and when the battery is low. Fully charging the battery takes up to 1.5 hours. Use a power source (an adapter or a power bank) with a USB port, of 5 V direct current/1.8 A or 9 W for your breast pump. Using a power source with lower specifications might result in slow charging or being unable to pump while connected to a wall socket.

- 1 Switch off the device (Fig. 2).
- 2 Insert the small plug of the USB cable into the motor unit and insert the other end of the USB cable into a USB power source, e.g. adapter or laptop.
- 3 The battery lights light up one by one and then all go out: the battery is charging (Fig. 3).
- 4 All 3 lights are steady on: the battery is fully charged (Fig. 4).
- 5 Remove the USB cable from the USB power source and pull the small plug out of the motor unit.

## Battery indicator light description

The battery indicator lights show the battery capacity.

- Three lights are steady on: the battery is fully charged for approximately 5 pumping sessions.
- Two lights are steady on: the battery is charged sufficiently for 2–3 sessions.
- One light flashes slowly: the battery is charged sufficiently for maximum 1 session. Charging is needed soon.
- First light flashes quickly: the battery is very low, the device will stop operating soon. Charging is needed.

One session is 15 minutes of pumping (double breast pump) or 30 minutes of pumping (single breast pump).

## Sizes of the breast shields and inserts

**⚠ Warning: Using an incorrectly sized breast shield and insert may cause pain, discomfort, and improper nipple positioning during pumping. Prolonged use may lead to nipple or breast tissue injury, inflammation, or blocked milk ducts. Ensure proper fit before use.**

The Philips Avent Hands-free Breast Pump has a soft breast shield that can be used individually or with a flexible silicone insert for a better comfort. There are multiple sizes available.

Measure your nipple to find the most suitable breast shield and, if necessary, the insert for you.

- 1 Use the nipple sizing card that is included in the package or refer to **[www.philips.com/support](http://www.philips.com/support)** to measure your nipples.
- 2 Stimulate your nipple prior to measuring your size.
- 3 Place the nipple sizing card around your nipple. Make sure there is a space of 0.5–2 mm between your nipple and the hole in the nipple sizing card.
- 4 Try various holes to see which opening fits your nipple best. If needed, view fitting in a mirror.
- 5 Select the breast shield (and the insert) that matches the hole that fits your nipple best.
- 6 Breast shield sizes 21 mm, 24 mm, and insert size 19 mm are included in the package.
- 7 Breast shield size 28 mm and insert sizes 15 mm, 17 mm, 26 mm can be purchased separately (see 'Ordering accessories and spare parts').
- 8 Start pumping by using the selected size(s) and closely monitor your comfort. After pumping, the nipple size might slightly change. If pumping becomes uncomfortable, measure your nipple size again and switch to another size, if needed.

**ⓘ Note:** You might need different sizes for your left nipple and your right nipple. Nipples may also change during your pumping journey. For example, the nipple size may change just after giving birth, after several months of breastfeeding or pumping, or may be different again after a second pregnancy.

## Assembling the breast pump

**⚠ Warning: Be careful, when you have disinfected the breast pump parts by boiling them, they can be very hot. To prevent burns, only start assembling the breast pump when disinfected parts have cooled down.**

**⚠ Warning: Wash your hands thoroughly with soap and water before you touch breast pump parts and breasts to prevent contamination. Avoid touching the inside of the cup.**

Make sure you have cleaned and disinfected the parts of the breast pump that come into contact with the breast and breast milk. Gather all parts and place them on a clean surface.

Check the breast pump parts for damage before use and replace them if necessary.

**ⓘ Note:** Correct placement of the breast shield, membrane, insert (if used), valve and silicone tube is essential for the breast pump to form a proper vacuum and to reduce the risk of milk leakage.

- 1 Wash your hands thoroughly with soap and water.
- 2 Push the valve into the cup as far as possible (Fig. 5).
- 3 Attach the membrane to the cup (Fig. 6).
- 4 Attach the front cap (Fig. 7) to the membrane and the cup. You can choose one of 4 possible positions of the front cap, to orient the silicone tube as preferred. Choose the most comfortable position for you. Push down the membrane cap until it is securely in place.
- 5 If you use an insert, first place the insert (Fig. 8) into the breast shield fully (Fig. 9).

**ⓘ Note:** If you use an insert, assemble it to the breast shield prior to attaching the breast shield to the cup, to avoid any gaps.

Make sure that you attach the insert to a breast shield of the corresponding size. The insert sizes 15 mm, 17 mm, 19 mm fit on the breast shield size 21 mm. The insert size 26 mm fits on the breast shield size 28 mm. Breast shield size 24 mm can only be used without an insert.

- 6 Attach the breast shield (including the insert, if any) to the cup (Fig. 10). Push the inner part of the breast shield into the cup as far as possible. Make sure that the rim of the breast shield covers the sides of the cup fully (Fig. 11).
- 7 Connect the silicone tube to the front cap (Fig. 12) and the motor unit (Fig. 13).

- 8 When the breast pump is charged, you can use it cordless. You can also use the breast pump connected to a power source. The battery will not charge while pumping.


**i** Note: Place the breast shield (dust) cover over the breast shield to keep the breast pump clean when you take it with you and while you are preparing for expressing (Fig. 14).

## Wearing the motor unit

The motor unit can be worn in multiple ways:


- 1 Use the provided clip to attach it to your clothes. The clip is pre-attached to the motor unit. Slide the clip over your clothes to fasten it tightly (Fig. 15).
- 2 Use the provided carrying cord to wear it across your body (Fig. 16). Attach the carrying cord to the clip. Make sure that the carrying cord and the silicone tubes do not get entangled. You can adjust the length of the carrying cord.

## Motor unit part description

- 1 USB power inlet
- 2 On/off button with pause function
- 3 Mode selection button
- 4 'Level down' button
- 5 'Level up' button
- 6 LED display
- 7 Clip (detachable)
- 8 Loop
- 9 Tubing port for single breast pump
- 10 Tubing ports for double breast pump

## LED display description



- 1 Battery indicator lights
- 2 Timer (minutes) and vacuum level indication
- 3 Expression icon
- 4 Stimulation icon

## Explanation of modes

The breast pump has two modes. Below you can find an explanation of these modes.

Mode icons	Mode	Explanation	No. of suction levels
	Stimulation mode 	Mode to stimulate the breast to start up the milk flow.	8 suction levels
	Expression mode 	Mode for efficient milk removal after milk has started flowing. After 90 sec. of stimulation, the device automatically switches to the expression mode.	16 suction levels

	The number on the display shows the running time in minutes.
	When you change the suction level, the display shows the selected suction level for a few seconds and then shows the running time again.

If you do not use the buttons for 8 seconds, the LED display dims and then switches off. Pressing any button will re-activate the LED display by showing the running time.

## Pumping instructions

**⚠ Warning: Always switch off the breast pump before you remove the collection cup from your breast, in order to release the vacuum.**

**⚠ Warning: Do not continue pumping for more than five minutes at a time if you do not succeed in expressing any milk. This can lead to pain and damage to your milk ducts. Try to express at another time during the day.**

**⚠ Warning: If the process becomes very uncomfortable or painful, stop using the breast pump and consult your healthcare professional.**

**⚠ Caution: The vacuum level generated by the device is sensitive to the ambient atmospheric pressure. At higher altitudes, the generated vacuum is less intense.**

- 1 Wash your hands thoroughly with soap and water. Clean your breast and nipple prior to placing the collection cup.
- 2 Place the assembled collection cup into your bra, on your breast, with the pouring spout (opening) positioned at the top (Fig. 17).
- 3 Make sure that the nipple is in the center of the breast shield.
- 4 Press and hold the on/off button until the device switches on (Fig. 2). The indicator lights go on to indicate that the device is in stimulation mode.
- 5 You can look through the top of the transparent cup to check whether the collection cup is positioned properly, and whether the nipple is moving. Your nipple should be in the center of the breast shield and should move in and out of the funnel (Fig. 17).
- 6 The breast pump starts in the stimulation mode and slowly increases its vacuum to the last used stimulation suction level (or level 5 when used for the first time) (Fig. 18).
- 7 The display shows the running time (Fig. 19).
- 8 Use the 'level up' and 'level down' buttons to change the suction level to your preference (Fig. 20).
- 9 After 90 seconds, the breast pump smoothly switches to expression mode and slowly increases its vacuum to the last used expression suction level (or level 11 when used for the first time) (Fig. 21).
- 10 If your milk starts flowing earlier, or if you do not have a milk flow after 90 seconds (and want to switch back to the stimulation mode), you can switch modes with the mode selection button (Fig. 22).
- 11 Briefly press the on/off button if you want to pause. To continue, briefly press the on/off button again.

**12** Press and hold the on/off button to switch off (Fig. 2).

**i** Note: You can use the double breast pump as a single breast pump by disconnecting the silicone tube from one collection cup. For your convenience, you can attach the silicone tube to the tube coupler (Fig. 23).

**i** Note: The collection cups can be worn in your regular nursing bra. If you experience discomfort or pain, try switching to a more stretchy bra.

## Using the breast pump when you have symptoms of engorgement

Breast engorgement is the swelling of the breasts, often caused by increased milk production. It occurs when the breasts become overly full and swollen, leading to discomfort. Symptoms may include hard, tight, or warm breasts, difficulty with milk ejection reflex, and redness, tenderness, or pain.

To reduce the risk of overstimulation and engorgement, follow these instructions:

- Gently massage the breasts before and during pumping to encourage milk flow.
- Use the lowest effective suction level and avoid prolonged pumping sessions. High suction levels or extended sessions may cause discomfort or potential breast tissue damage.
- Pump only the amount that your baby needs.

If you experience symptoms of engorgement, consider the following measures:

- Pumping gently for a short duration can help relieve discomfort. Pump just enough to relieve the discomfort as excessive pumping may worsen engorgement.
- If discomfort persists or if you experience redness, warmth, fever, or flu-like symptoms, consult a healthcare professional.

**⚠ Warning: Overstimulation of the breasts through high suction levels or prolonged pumping may lead to discomfort or potential breast tissue damage.**

**⚠ Warning: Excessive or prolonged pumping may increase breast swelling, tenderness, or worsen engorgement. If you experience these symptoms, stop pumping and consult a healthcare professional.**

## After use

**⚠ Warning: Refrigerate or freeze expressed breast milk immediately or keep it at a temperature of 16-29 °C for a maximum of 4 hours before you feed your baby.**

- 1 After completing the pumping session, lift the collection cup out of your bra. Place the collection cup on a clean and flat surface (Fig. 24). Be careful not to spill any milk from the pouring spout.
- 2 To store the breast milk: pour the expressed breast milk from the collection cup into a separate milk storage container (e.g., a bottle or a breast milk storage bag) (Fig. 25). Do not store expressed breast milk in the collection cup.
- 3 Disassemble, clean and disinfect all parts that have been in contact with the breast and breast milk (see 'Cleaning and disinfecting').

## Storing breast milk

Below you find guidelines for storing breast milk:


Storage location	Temperature	Max. storage duration
Counter-top	25 °C or colder (room temperature)	max. 2 hours
Refrigerator	4 °C	max. 24 hours
Freezer	< -18 °C or colder	max. 3 months

**⚠ Warning: Never refreeze thawed breast milk to avoid a decrease in milk quality.**

**⚠ Warning: Never add fresh breast milk to frozen breast milk to avoid a decrease in milk quality and unintended thawing of the frozen milk.**

## Compatibility

The collection cup of the Philips Avent Hands-free Electric Breast Pump is compatible with the following motor unit of Philips Avent breast pump Advanced (SCF391, SCF393), Plus (SCF391, SCF393), Premium (SCF392), Premium Plus (SCF392):

Do not combine it with motor units of other Philips Avent breast pumps or different brands.

## Ordering accessories and spare parts

To buy accessories for convenience, extra collection cups (for example, for using at work and at home), breast shields and inserts of smaller or larger size, or spare parts, go to your Philips dealer. You can also contact the Philips Consumer Care Center in your country.

## Recycling

- This symbol means that electrical products and batteries shall not be disposed of with normal household waste (Fig. 32).
- Follow your country's rules for the separate collection of electrical products and batteries.
- Correct disposal helps prevent negative consequences for the environment and human health.

## Removal of built-in rechargeable battery

The built-in rechargeable battery must only be removed by a qualified professional when the appliance is discarded. Before removing the battery, make sure that the appliance is disconnected from the wall socket and that the battery is completely empty.

---

**Take any necessary safety precautions when you handle tools to open the appliance and when you dispose of the rechargeable battery.**

---

**When you handle batteries, make sure that your hands, the product and the batteries are dry.**

---

**To avoid accidental short-circuiting of batteries after removal, do not let battery terminals come into contact with metal objects (e.g. coins, hairpins, rings). Do not wrap batteries in aluminum foil. Tape battery terminals or put batteries in a plastic bag before you discard them.**

---

**If batteries are damaged or leaking, avoid contact with the skin or eyes. If this occurs, immediately rinse well with water and seek medical care.**

---

- 1 Remove the rubber feet from the bottom of the motor unit (Fig. 26).
- 2 Remove the screw with a torx screwdriver (Fig. 27).
- 3 Open the motor unit by pressing at the sides of the device (Fig. 28).

- 4 Disconnect the connectors from the printed circuit board (Fig. 29).
- 5 Remove the frame from the motor unit (Fig. 30).
- 6 Cut or fold the edges of the plastic battery frame. Remove the battery from the plastic battery frame (Fig. 31).

**(i) Note: Dispose of the rechargeable battery at a separate local collection point. Other parts can be disposed at a specific place for disposal of electric equipment.**

## Support

For all product support, power consumption details and the EU Declaration of Conformity, please visit [www.philips.com/support](http://www.philips.com/support).

## Troubleshooting

This chapter summarizes the most common problems you could encounter with the breast pump. If you are unable to solve the problem with the information below, contact the Consumer Care Center in your country.

Problem	Solution
I experience pain when I use the breast pump.	<ul style="list-style-type: none"> <li>- Check if you use the most suitable breast shield size and, if necessary, the insert size (see 'Sizes of the breast shields and inserts').</li> <li>- Choose a suction level that is comfortable for you (see 'Pumping instructions').</li> <li>- Check if the membrane is not damaged (does not have small cracks, holes etc.).</li> <li>- Make sure that your nipple is placed in the center of the opening of the breast shield.</li> <li>- If the pain persists, stop using the breast pump and consult your healthcare professional.</li> </ul>

Problem	Solution
I do not feel any suction/the suction level is too low.	<ul style="list-style-type: none"> <li>- Check if the valve is not damaged, the opening in its pointed part closes fully, and that the valve is assembled correctly.</li> <li>- Check if the membrane is not damaged (does not have small cracks, holes etc.).</li> <li>- Make sure that all the parts are clean and assembled correctly (see 'Assembling the breast pump').</li> <li>- Check if you use the most suitable breast shield size and, if necessary, the insert size (see 'Sizes of the breast shields and inserts').</li> <li>- Make sure that the collection cup is positioned correctly on the breast to allow a seal between the breast shield and your breast.</li> <li>- If possible, try a higher suction level.</li> <li>- Before pumping, try warming your breast with a warm towel and massaging it to stimulate the areola so the breast is completely unblocked.</li> </ul>
The breast pump is not charging.	<ul style="list-style-type: none"> <li>- Make sure that you use a power source (an adapter or a power bank) with a USB port, of 5 V direct current/1.8 A or 9 W for your breast pump.</li> <li>- The breast pump is too cold or too hot for charging. Let the breast pump acclimatize to room temperature for 30 minutes.</li> </ul>
The breast pump does not switch on (indicator lights do not go on).	<p>Make sure that the battery is charged (see 'Charging the breast pump') or connect the breast pump to a power source. Use a power source (an adapter or a power bank) with a USB port, of 5 V direct current/1.8 A or 9 W for your breast pump.</p>
The suction of the breast pump is too strong.	<ul style="list-style-type: none"> <li>- Before pumping, try warming your breast with a warm towel and massaging it to stimulate the areola so the breast is completely unblocked.</li> <li>- If possible, try a lower suction level. When you express for the first time, start with the default suction level (level 5 in stimulation mode, level 11 in expression mode) and increase/decrease the level when needed. During different expression sessions, different levels may feel most comfortable for you.</li> <li>- Make sure that you only use the breast pump parts provided with this breast pump.</li> <li>- Check if the membrane is not damaged (does not have small cracks, holes etc.).</li> </ul>

## Undesirable side effects

When using the breast pump, the below undesirable side effects may occur. If you experience any of these symptoms, contact a healthcare professional or breastfeeding specialist.

<b>Undesirable side effect</b>	<b>Description</b>
Pain	Pain felt on the breast or nipple, or felt as a response to the suction being above a comfortable level.
Sore nipples	Persistent pain in the nipples at the beginning of the pumping session, or that lasts throughout the entire pumping session, or pain between sessions.
Engorgement	Swelling of the breast. The breast may feel hard, lumpy and tender. May include erythema (redness) of the breast area and fever. Usually occurs during the first days of lactation.
Erythema (redness)	Redness of the skin of the breast.
Bruise, thrombus	A reddish-purple discoloration that does not fade when pressed. When a bruise fades, it may become green and brown.
Nipple trauma	<ul style="list-style-type: none"> <li>- Fissures or cracked nipples.</li> <li>- Skin tissue peeling off the nipple. Normally occurs in combination with cracked nipples and/or blisters.</li> <li>- Nipple tear.</li> <li>- Blister. Looks like small bubbles on the surface of the skin.</li> <li>- Bleeding. Cracked or torn nipples can lead to bleeding of the affected area.</li> </ul>
Clogged mammary ducts	A red, tender lump on the breast. May include erythema (redness) of the breast area. Can lead to mastitis (breast inflammation) and fever if left untreated.
Rash or allergic reaction	Redness or irritation of the skin of the breast, it can be itchy, scaly, painful or blistered.

## Storage

- Keep the breast pump out of direct sunlight as prolonged exposure may cause discoloration. Store the breast pump and its accessories in a safe, clean and dry place.
- If the device has been stored in a hot or cold environment, place it in an environment with a temperature of 20 °C for 30 minutes to let it reach a temperature within the usage conditions (5 °C to 40 °C) before you use it.
- To preserve the battery life, fully charge the battery before storing your breast pump for an extended period of time.

## Electromagnetic fields (EMF)

This Philips device complies with all applicable standards and regulations regarding exposure to electromagnetic fields.

## Technical specifications

Vacuum level	Stimulation: -70 to -225 mbar (-53 to -169 mmHg) Expression: -70 to -360 mbar (-53 to -270 mmHg) Intermittent vacuum
Cycle speed	Stimulation: 105-120 cycles/min, depending on vacuum level Expression: 53-85 cycles/min, depending on vacuum level
Motor unit rated input	5 V direct current / 1.8 A
Battery	Lithium-ion rechargeable battery: 3.6 V, 2600 mAh, 9.4 Wh (internal)
Automatic power-off	After 30 minutes (after the last user interaction), after 15 minutes (in the pause mode)
Power port	USB-C
Operating conditions	Temperature: 5 °C to 40 °C Relative humidity: 15 to 90 % (non-condensing)
Storage and transportation conditions	Temperature: -20 °C to 45 °C Relative humidity: up to 90 % (non-condensing)
Operating and storage pressure	700-1060 hPa (<3000 m altitude)

Net weight	275 g (single breast pump), 355 g (double breast pump) 125 g (assembled collection cup with a breast shield (dust) cover)
External dimensions	Motor unit: 145 mm x 95 mm x 45 mm (length x width x height) Collection cup: 67 mm x 117 mm x 117 mm (length x width x height)
Device classification	Applicable classifications Degree of protection against electric shock: - Collection cup: type BF applied part The motor unit and external power supply (not provided) are a medical electrical system.
Materials	Motor unit: MABS, ABS, silicone, TPE. Collection cup: polypropylene, silicone. Silicone tube: silicone. Clip, tube coupler: ABS. Carrying cord: polyester, polypropylene, ABS.
Service life	2 years
Operational life	500 hours
Mode of operation	Continuous operation
Ingress protection	IP22

## Electromagnetic compatibility – compliance information

This device has been tested and found to comply with IEC 60601-1-2:2020 for electromagnetic compatibility (EMC) as Class B according to CISPR 11:2015 +A1:2016 +A2:2019. Tests levels are listed in the accompanying tables.

<b>Emissions</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11	Group 1	The device does not use RF energy for its intended function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

<b>Emissions</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11	Class B	The device is suitable for use in locations in residential environments and in establishments directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Complies Class A Category	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

### **Guidance and manufacturer's declaration – electromagnetic immunity**

The Philips Avent Hands-free Single/Double Electric Breast Pump is intended for use in the electromagnetic environment specified below. The customer or the user of the breast pump should assure that it is used in such an environment. The essential performance functions are the accuracy of the vacuum pressure and the backflow protection to withstand milk leakage into the motor unit. The Philips Avent Hands-free Single/Double Electric Breast Pump was tested for immunity to electromagnetic disturbances and passed using the criteria of no change in modes, vacuum accuracy ( $\pm 20\%$ ), backflow protection and operation of the breast pump.

<b>Immunity test</b>	<b>Immunity test level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 8$ kV contact $\pm 15$ kV air	Floors can be a source of static electricity. When the air is very dry ( $< 20\%$ RH) the user is advised to prefer to use the device in a room with concrete or wooden floor (instead of using it in a room with a synthetic floor).
Radiated RF EM fields IEC 61000-4-3	10 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	

Immunity test	Immunity test level	Electromagnetic environment - guidance
Proximity fields from RF wireless communications equipment IEC 61000-4-3	385 (27), 450 (28), 710 (9), 745 (9), 780 (9), 810 (28), 870 (28), 930 (28), 1720 (28), 1845 (28), 1970 (28), 2450 (28), 5240 (9), 5500 (9), 5785 (9) MHz (V/m)	The device is immune against any RF wireless communication device such as smart phone, Wi-Fi or Bluetooth device at a distance of 30 cm (12 in) or more. Avoid a closer distance as this may result in improper operation.
RATED power frequency magnetic fields IEC 61000-4-8	30 A/m 50 Hz or 60 Hz	
Electrical fast transient/ burst IEC 61000-4-4	±2 kV 100 kHz repetition frequency	Mains power quality expected of a typical commercial or residential environment.
Surges IEC 61000-4-5	±1 kV Line-to-line ±2 kV Line-to-ground	Mains power quality expected of a typical commercial or residential environment.
Conducted disturbances induced by RF fields IEC 61000-4-6	3 Vrms, before modulation is applied 0.15 MHz – 80 MHz 6 Vrms, before modulation is applied, in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	

Immunity test	Immunity test level	Electromagnetic environment - guidance
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles single phase at 0° 0 % UT; 250/300 cycle	Mains power quality expected of a typical commercial or residential environment. At power interruptions the breast pump will automatically switch-over to battery mode without interruption of the function.
Radiated field in close proximity IEC 61000-4-39	30.00 kHz, CW mod., 8.0 A/m 134.20 kHz, PM mod., 65.0 A/m 13.56 MHz, PM mod., 7.5 A/m	

Note: UT is the AC mains voltage prior to application of the test level. If in the unlikely event the Philips Avent Hands-free Single/Double Electric Breast Pump is locked inoperative in a function by disturbances from other devices beyond the levels stated in the table above, the device will shut-off completely and restart. In that case make sure the distance to the other device causing this disturbance is enlarged.

## Explanation of symbols

The warning signs and symbols are essential to ensure that you use this product safely and correctly and to protect you and others from injury.

Sym- bol	Explanation
	Indicates to follow the instructions for use.
	Indicates important information such as warnings and cautions.

Sym- bol	Explanation
	Indicates usage tips, additional information or a note.
	Indicates the manufacturer.
	Indicates the manufacturing date.
	Indicates the compliance with the general safety and performance requirements of Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.
	Indicates manufacturer's catalog number of the breast pump.
	Indicates the batch number of the medical device.
	Indicates the Global Trade Identification Number and the Lot Code of the breast pump.
	This symbol indicates that the device is a medical device.
	Indicates that the part of the appliance that comes into physical contact with the user (also known as the applied part) is of type BF (Body Floating) according to IEC 60601-1. The applied part is the collection cup.
	Indicates the protection against ingress of solid foreign objects greater than 12.5 mm and against harmful effects due to the ingress of dripping water when tilted at 15 degrees.
	Indicates 'direct current'.

Symbol	Explanation
	Indicates USB.
	Indicates the storage and transportation temperature limits to which the medical device can be safely exposed: -4 °F to 113 °F/-20 °C to 45 °C.
	Indicates the relative humidity upper limits to which the device can be safely exposed: up to 90 %.
	Indicates the operating pressure range: 700 to 1060 hPa.
	Indicates that the breast pump must be kept dry.
	Indicates that the breast pump must be kept out of direct sunlight.
	Indicates 2-year Philips worldwide guarantee.
	Indicates Forest Stewardship Council™. The FSC trademarks enable consumers to choose products that support forest conservation, offer social benefits, and enable the market to provide an incentive for better forest management.
	Indicates that the material used in the product is safe for food contact.
	Identifies the UDI carrier, including the AIDC and human readable information.

Symbol	Explanation
 A recycling symbol consisting of three chasing arrows forming a triangle. The number '20' is in the center, and 'PAP' is written below it.	Indicates that the material used for the packaging is corrugated cardboard/fiberboard.
 A recycling symbol consisting of three chasing arrows forming a triangle. The number '21' is in the center, and 'PAP' is written below it.	Indicates that the material used for the packaging is non-corrugated cardboard/fiberboard.
 A recycling symbol consisting of three chasing arrows forming a triangle. The number '4' is in the center, and 'LDPE' is written below it.	Indicates that the material used for the plastic bag is low density polyethylene.

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Trademarks are the property of Koninklijke Philips N.V.
All rights reserved.
Philips Consumer Lifestyle B.V.