

Clinical Evidence Series

Avatar-based versus conventional patient monitoring with distant vision: a computer-based simulation study



What's different

Avatar-based monitoring converts multiple vital sign parameters into a unified visual representation, allowing clinicians to perceive patient status using visual cues, even at a distance.



Interpretation

When viewed from 8-16 meters, avatar-based monitoring was associated with substantially improved recognition of vital sign abnormalities, suggesting greater efficiency of visual information transfer compared with conventional displays.



Key takeaways

Vital sign recognition was **74% higher at 8 m** and **51% higher at 16 m** with avatar-based monitoring



Background

- Conventional patient monitoring requires clinicians to sequentially process and integrate multiple parameters which may delay recognition of critical changes and contribute to reduced situation awareness
- Avatar-based monitoring is designed to support faster perception and interpretation of patient status, including in situations where clinicians are not positioned directly by the monitor



Study design and methods

- Prospective, single-center, computer-based simulation study conducted at a Swiss academic hospital
- Within-participant comparison of avatar-based vs. conventional monitoring
- Participants viewed monitors at 8-meter and 16-meter distances
- Two predefined scenarios including 11 vital signs per scenario
- Outcomes included the number of correctly recalled vital signs and correct classification (too low/safe/too high)



Objectives

Compare avatar and conventional monitoring for recognition of vital sign changes when viewed from a distance.



Results

Vital sign recognition is significantly higher at a distance with avatar-based monitoring:

- o 74% increase at 8 m ($p < 0.001$)
- o 51% increase at 16 m ($p < 0.001$)