

**Philips Medical Systems  
DICOM Conformance Statement**

**THORAVISION 3.3**

Document Number 4512 130 09631  
16 November 1995

© Copyright Philips Medical Systems Nederland B.V. 1995

All rights reserved

Philips  
Medical  
Systems



**PHILIPS**

**Issued by:**

Philips Medical Systems Nederland B.V.  
Integrated Clinical Solutions  
Building QP-2129  
P.O. Box 10.000  
5680 DA Best  
The Netherlands  
Tel.: +31 40 2763499  
Fax.: +31 40 2762673  
email: dicom@best.ms.philips.com

## Table of Contents

|       |  |   |
|-------|--|---|
| 1     | <b>Introduction</b> .....                              | 1 |
| 1.1   | Scope and field of application .....                   | 1 |
| 1.2   | Intended audience .....                                | 1 |
| 1.3   | Contents and structure .....                           | 1 |
| 1.4   | Used definitions, terms and abbreviations .....        | 1 |
| 1.5   | References .....                                       | 2 |
| 1.6   | Important note to the reader .....                     | 2 |
| 2     | <b>Implementation model</b> .....                      | 3 |
| 2.1   | Application Data Flow Diagram .....                    | 3 |
| 2.2   | Functional definition of Application Entities .....    | 3 |
| 2.3   | Sequencing of Real World Activities .....              | 3 |
| 3     | <b>AE Specifications</b> .....                         | 4 |
| 3.1   | AE THORAVISION DICOM Export Specification .....        | 4 |
| 3.1.1 | Association Establishment Policies .....               | 4 |
| 3.1.2 | Association Initiation Policy .....                    | 4 |
| 3.1.3 | Association Acceptance Policy .....                    | 6 |
| 4     | <b>Communication Profiles</b> .....                    | 7 |
| 4.1   | TCP/IP Stack .....                                     | 7 |
| 4.1.1 | Physical Media Support .....                           | 7 |
| 5     | <b>Extensions/Specializations/Privatizations</b> ..... | 7 |
| 6     | <b>Configuration</b> .....                             | 7 |
| 6.1   | AE Title/Presentation Address mapping .....            | 7 |
| 6.1.1 | Local AE Titles and Presentation Addresses .....       | 7 |
| 6.1.2 | Remote AE Titles and Presentation Addresses .....      | 7 |
| 6.2   | Configurable parameters .....                          | 8 |
| 7     | <b>Support of Extended Character Sets</b> .....        | 8 |
| 8     | <b>Extended Computed Radiography IOD</b> .....         | 9 |



## 1 Introduction

This chapter provides general information about the purpose, scope and contents of this Conformance Statement.

### 1.1 Scope and field of application

The scope of this DICOM Conformance Statement is to facilitate data exchange with equipment of Philips Medical Systems. This document specifies the compliance to the DICOM standard (formally called the NEMA PS 3.X-1993 standards). It contains a short description of the applications involved and provides technical information about the data exchange capabilities of the equipment. The main elements describing these capabilities are: the supported DICOM Service Object Pair (SOP) Classes, Roles, Information Object Definitions (IOD) and Transfer Syntaxes.

The field of application is the integration of the Philips Medical Systems equipment into an environment of medical devices.

This Conformance Statement should be read in conjunction with the DICOM standard and its addenda. The conformance to the DICOM standard is a key element of the Inturis Program (see [INTURIS]).

### 1.2 Intended audience

This Conformance Statement is intended for:

- potential clients,
- system integrators of medical equipment,
- marketing staff interested in system functionality,
- software designers implementing DICOM interfaces.

It is assumed that the reader is familiar with the DICOM standard.

### 1.3 Contents and structure

The DICOM Conformance Statement is contained in chapter 2 through 7 and follows the contents and structuring requirements of DICOM PS 3.2-1993 and Supplement 2 (in case of Media specifications).

Additionally, the chapters following 7 specify the details of the applied IODs.

### 1.4 Used definitions, terms and abbreviations

DICOM definitions, terms and abbreviations are used throughout this Conformance Statement. For a description of these, see NEMA PS 3.3-1993 and PS 3.4-1994.

The word Philips in this document refers to Philips Medical Systems.

## 1.5 References

[DICOM] The Digital Imaging and Communications in Medicine (DICOM) standard:  
NEMA PS 3.X (X refers to the part 1 - 13) and Supplements  
National Electrical Manufacturers Association (NEMA) Publication Sales  
1300 N. 17th Street, Suite 1847  
Rosslyn, Va. 22209, United States of America

[INTURIS] Philips Inturis Program  
Integrated Clinical Solutions  
Philips Medical Systems Nederland B.V. (see address at page ii)

## 1.6 Important note to the reader

This Conformance Statement by itself does not guarantee successful interoperability of Philips equipment with non-Philips equipment. The user (or user's agent) should be aware of the following issues:

- **Interoperability**

Interoperability refers to the ability of application functions, distributed over two or more systems, to work successfully together. The integration of medical devices into a networked environment may require application functions that are not specified within the scope of DICOM. Consequently, using only the information provided by this Conformance Statement does not guarantee interoperability of Philips equipment with non-Philips equipment. It is the user's responsibility to analyse thoroughly the application requirements and to specify a solution that integrates Philips equipment with non-Philips equipment.

- **Validation**

Philips equipment has been carefully tested to assure that the actual implementation of the DICOM interface corresponds with this Conformance Statement.

Where Philips equipment is linked to non-Philips equipment, the first step is to compare the relevant Conformance Statements. If the Conformance Statements indicate that successful information exchange should be possible, additional validation tests will be necessary to ensure the functionality, performance, accuracy and stability of image and image related data. It is the responsibility of the user (or user's agent) to specify the appropriate test suite and to carry out the additional validation tests.

- **New versions of the DICOM Standard**

The DICOM Standard will evolve in future to meet the user's growing requirements and to incorporate new features and technologies. Philips is actively involved in this evolution and plans to adapt its equipment to future versions of the DICOM Standard. In order to do so, Philips reserves the right to make changes to its products or to discontinue its delivery.

The user should ensure that any non-Philips provider linking to Philips equipment, also adapts to future versions of the DICOM Standard. If not, the incorporation of DICOM enhancements into Philips equipment may lead to loss of connectivity (in case of networking) and incompatibility (in case of media).

## 2 Implementation model

THORAVISION 3.3 of Philips Medical Systems is an image generating system. It contains an Export function based on the DICOM Image Storage to transfer image data from the THORAVISION system to a remote system. This DICOM Export function is described in this document.

### 2.1 Application Data Flow Diagram

The THORAVISION DICOM Export transfers a THORAVISION image to a remote DICOM node. This is activated by an operator request or automatically if the system is configured to do so. A remote destination is selected from the user interface, followed by the selection of the image to be transferred.

Post-processed image data can be transferred (pixel value between 0 and 30,000) as an instance of the DICOM Computed Radiography IOD. The images transferred are intended for viewing purpose and VOI LUT transformation (grayscale transformation) only.

The THORAVISION DICOM Export behaves as a single Application Entity. The related Implementation Model is shown in Figure 2-1 on page 3.




Figure 2-1: THORAVISION DICOM Export Implementation Model

### 2.2 Functional definition of Application Entities

The THORAVISION DICOM Export application entity acts as a Service Class User (SCU) of the Storage Service Class. After invoking it will open an association to the remote system. For each image to be transported a retrieve action from the internal THORAVISION storage will take place followed by the conversion to a DICOM message to be transferred to the remote system.

### 2.3 Sequencing of Real World Activities

Not applicable.

### 3 AE Specifications

THORAVISION DICOM Export acts as a single Application Entity.

#### 3.1 AE THORAVISION DICOM Export Specification

The THORAVISION Export Application Entity provides Standard Extended Conformance to the following DICOM 3.0 SOP class as an SCU:

**Table 3-1: Supported SOP class by the THORAVISION Export AE as SCU**

| <i>SOP class Name</i>              | <i>UID</i>                |
|------------------------------------|---------------------------|
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 |

#### 3.1.1 Association Establishment Policies

##### 3.1.1.1 General

THORAVISION Export will offer a configurable maximum PDU size in steps of 1K, 2K, 4K, 8K and 16K (default is 16K = 16384 bytes) on associations initiated by the application itself.

##### 3.1.1.2 Number of Associations

THORAVISION Export will attempt to establish one association at a time.

##### 3.1.1.3 Asynchronous Nature

THORAVISION Export does not support asynchronous operations and will not perform asynchronous window negotiation.

##### 3.1.1.4 Implementation Identifying Information

The Implementation Class UID is: "1.3.46.670589.8.3.3".

The implementation version name is: "THORAVISION 3.3".

#### 3.1.2 Association Initiation Policy

THORAVISION Export initiates associations as a result of the following events:

- the THORAVISION operator requests the export of one THORAVISION image to a remote system,
- the generation of a new image will result in an automatic export of that image when the system is configured in automatic store mode.

##### 3.1.2.1 Export from THORAVISION system

###### 3.1.2.1.1 Associated Real-World Activity

The THORAVISION Export function will be accessible through the THORAVISION User Interface.

After the transfer of an image the association is released. The transferred image will not be deleted from the system. The transfer will be repeated automatically when the transfer was unsuccessful with special response status conditions (e.g. Store SCP down).

### 3.1.2.1.2 Proposed Presentation Contexts

THORAVISION Export will propose the following presentation contexts:

**Table 3-2: Proposed Presentation Contexts for THORAVISION Export**

| <i>Presentation Context table</i>  |                           |                           |                     |             |                             |
|------------------------------------|---------------------------|---------------------------|---------------------|-------------|-----------------------------|
| <i>Abstract Syntax</i>             |                           | <i>Transfer Syntax</i>    |                     | <i>Role</i> | <i>Extended Negotiation</i> |
| <i>Name</i>                        | <i>UID</i>                | <i>Name List</i>          | <i>UID List</i>     |             |                             |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU         | None                        |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU         | None                        |
| Computed Radiography Image Storage | 1.2.840.10008.5.1.4.1.1.1 | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 | SCU         | None                        |

#### 3.1.2.1.2.1 SOP Specific Conformance to Storage SOP Classes

The status of the C-STORE Response (Success, Refused, Error, Warning) will be displayed via the user interface.

Extended negotiation is not supported.

Table 3-3 lists the applied optional and extended modules and attributes of the Extended CR IOD. Conditional attributes Patient Orientation (type 2C), Image Date (type 2C), Image Time (type 2C) are always present.

**Table 3-3: Applied optional/extended Modules and Attributes of the Extended CR IOD**

| <i>IE</i> | <i>Module</i>              | <i>Optional/Extended Attributes</i>   |
|-----------|----------------------------|---|
| Patient   | Patient                    | Other Patient's ID  |
| Study     | General Study              | Referring Physician's Name, Study Description   |
| Series    | General Series             | -   |
|           | CR Series                  | Filter Type, Collimator/grid Name, Focal Spot(s), Plate Type  |
| Equipment | General Equipment          | Institution Name, Station Name, Institutional Department Name, Manufacturer's Model name, Device Serial Number, Software Version(s), Date of Last Calibration, Time of Last Calibration |
| Image     | General Image              | Image Type, Image Comments  |
|           | Image Pixel                | -   |
|           | CR Image (extended module) | KVP, Distance Source to Detector, Exposure Time, Exposure, Generator Power, Sensitivity<br>For Extensions see "Extensions/Specializations/Privatizations" on page 7                     |

**Table 3-3: Applied optional/extended Modules and Attributes of the Extended CR IOD (Continued)**

| <i>IE</i> | <i>Module</i>                        | <i>Optional/Extended Attributes</i>   |
|-----------|--------------------------------------|---|
|           | X-Ray Collimator                     | Collimator Shape, Collimator Left Vertical Edge, Collimator Right Vertical Edge, Collimator Upper Horizontal Edge, Collimator Lower Horizontal Edge |
|           | VOI LUT<br>(applied optional module) | Window Center   |
|           | Stamp<br>(applied private module)    | See "Extensions/Specializations/Privatizations" on page 7   |
|           | SOP Common                           | Specific Character Set  |

**3.1.3 Association Acceptance Policy**

THORAVISION 3.3 does not accept associations.

## 4 Communication Profiles

### 4.1 TCP/IP Stack

THORAVISION Release 3.3 provides DICOM 3.0 TCP/IP Network Communication Support as defined in Part 8 of the DICOM 3.0 Standard.

#### 4.1.1 Physical Media Support

The THORAVISION system supports ISO 8802-3 10BASE5 Ethernet.

## 5 Extensions/Specializations/Privatizations

The applied Standard Extended CR Image Storage SOP Class is extended with additional Standard Type 3 attributes and with one private attribute, as specified in the table below.

**Table 5-1: Applied extended Modules and applied extended Attributes**

| <i>IE</i> | <i>Module</i>              | <i>Extended Attributes</i>  |
|-----------|----------------------------|---|
| Image     | CR Image (extended module) | Preprocessing Function (from standard NM Image Module),<br>Post Processing Function (from standard NM Image Module),<br>Pixel Spacing (from standard Image Plane Module),<br>Image Area Dose Product (from X-Ray Acquisition Module)  |
|           | Stamp (private module)     | Private Creator (private attribute),<br>Stamp Image Sequence containing (Samples per Pixel, Photometric Interpretation, Rows, Columns, Bits Allocated, Bits Stored, High Bit, Pixel Representation, Pixel Data)<br>All these attributes are from standard Image Pixel Module. |

See the specification in Extended Computed Radiography IOD.

## 6 Configuration

The configuration of a THORAVISION system is done by means of updating configuration files. This should be done by Philips service engineers only.

### 6.1 AE Title/Presentation Address mapping

#### 6.1.1 Local AE Titles and Presentation Addresses

The local Application Entity Title and Presentation Address are defined by updating Configuration files.

#### 6.1.2 Remote AE Titles and Presentation Addresses

All remote applications to be selected as export destination (SCP) are configurable for the following items:

- The Application Entity Title of the remote application.
- The Presentation Address at which the remote application should accept association

## Support of Extended Character Sets

requests.

- The Remote Host Name of the system on which the remote application resides.

## 6.2 Configurable parameters

The maximum PDU size is configurable.

The supported Extended Character Set is configurable. Default Character Set is ISO-IR 100 which is the Latin alphabet No. 1, supplementary set.

## 7 Support of Extended Character Sets

THORAVISION Export supports Character Sets ISO-IR 100 (which is default), ISO-IR 101, ISO-IR 109, ISO-IR 110, ISO-IR 144, ISO-IR 127, ISO-IR 126, ISO-IR 138, ISO-IR 148.

## 8 Extended Computed Radiography IOD

The modules selected from the IOD module table of DICOM 3.0 and the extended modules are given in the table below.

**Table 8-1: Applied Modules in the Extended CR IOD**

| <i>IE</i> | <i>Module</i>     |
|-----------|-------------------|
| Patient   | Patient           |
| Study     | General Study     |
| Series    | General Series    |
|           | CR Series         |
| Equipment | General Equipment |
| Image     | General Image     |
|           | Image Pixel       |
|           | Extended CR Image |
|           | X-Ray Collimator  |
|           | VOI LUT           |
|           | Stamp             |
|           | SOP Common        |

The details of these applied modules are given in the tables below. The list of possible values are given (if applicable). The situation that an attribute is present conditionally/optionally or that an attribute may contain a zero length value, is indicated too. Conditions and Defined/Enumerated Values of DICOM 3.0 are applicable but are not shown in the tables.

**Table 8-2: Patient Module**

| <i>Attribute Description</i> | <i>Tag</i> | <i>Note</i>                                      |
|------------------------------|------------|--|
| Patient's Name               | 0010,0010  |  |
| Patient ID                   | 0010,0020  | zero length value if not entered by the operator |
| Patient's Birth Date         | 0010,0030  | zero length value if not entered by the operator |
| Patient's Sex                | 0010,0040  | may have a zero length value                     |
| Other Patient IDs            | 0010,0100  | zero length value if not entered by the operator |

**Table 8-3: General Study Module**

| <i>Attribute Description</i> | <i>Tag</i> | <i>Note</i>                                      |
|------------------------------|------------|--|
| Study Instance UID           | 0020,000D  |  |
| Study Date                   | 0008,0020  |  |
| Study Time                   | 0008,0030  |  |
| Referring Physician's Name   | 0008,0090  | zero length value if not entered by the operator |
| Study ID                     | 0020,0010  |  |
| Accession Number             | 0008,0050  | may have a zero length value                     |
| Study Description            | 0008,1030  | zero length value if not entered by the operator |

**Table 8-4: General Series Module**

| <i>Attribute Description</i> | <i>Tag</i> | <i>Note</i> |
|------------------------------|------------|-------------|
| Modality                     | 0008,0060  | CR          |
| Series Instance UID          | 0020,000E  |             |
| Series Number                | 0020,0011  | 1           |
| Series Date                  | 0008,0021  |             |
| Series Time                  | 0008,0031  |             |

**Table 8-5: CR Series Module**

| <i>Attribute Description</i> | <i>Tag</i> | <i>Note</i>                  |
|------------------------------|------------|------------------------------|
| Body Part Examined           | 0018,0015  | CHEST                        |
| View Position                | 0018,5101  |                              |
| Filter Type                  | 0018,1160  |                              |
| Collimator/grid Name         | 0018,1180  |                              |
| Focal Spot                   | 0018,1190  |                              |
| Plate Type                   | 0018,1260  | may have a zero length value |

**Table 8-6: General Equipment Module**

| <i>Attribute Name</i> | <i>Tag</i> | <i>Note</i>             |
|-----------------------|------------|-------------------------|
| Manufacturer          | 0008,0070  | Philips Medical Systems |

**Table 8-6: General Equipment Module (Continued)**

| <i>Attribute Name</i>         | <i>Tag</i> | <i>Note</i>                                      |
|-------------------------------|------------|--|
| Institution Name              | 0008,0080  | zero length value if not entered by the operator |
| Station Name                  | 0008,1010  |  |
| Institutional Department Name | 0008,1040  |  |
| Manufacturer's Model Name     | 0008,1090  |  |
| Device Serial Number          | 0018,1000  |  |
| Software Versions             | 0018,1020  |  |
| Date of Last Calibration      | 0018,1200  | may have a zero length value                     |
| Time of Last Calibration      | 0018,1201  | may have a zero length value                     |

**Table 8-7: General Image Module**

| <i>Attribute Name</i> | <i>Tag</i> | <i>Note</i>                                      |
|-----------------------|------------|--|
| Image Number          | 0020,0013  |  |
| Patient Orientation   | 0020,0020  |  |
| Image Date            | 0008,0023  |  |
| Image Time            | 0008,0033  |  |
| Image Type            | 0008,0008  | DERIVED\PRIMARY                                  |
| Image Comments        | 0020,4000  | zero length value if not entered by the operator |

**Table 8-8: Image Pixel Module and Extended CR Image Module<sup>a</sup>**

| <i>Attribute Name</i>      | <i>Tag</i> | <i>Note</i>                |
|----------------------------|------------|----------------------------|
| Samples per Pixel          | 0028,0002  | 1                          |
| Photometric Interpretation | 0028,0004  | MONOCHROME1 or MONOCHROME2 |
| Rows                       | 0028,0010  |                            |
| Columns                    | 0028,0011  |                            |
| Bits Allocated             | 0028,0100  | 16                         |
| Bits Stored                | 0028,0101  | 15                         |
| High Bit                   | 0028,0102  | 14                         |
| Pixel Representation       | 0028,0103  | 0000H                      |

**Table 8-8: Image Pixel Module and Extended CR Image Module<sup>a</sup> (Continued)**

| <i>Attribute Name</i>                | <i>Tag</i> | <i>Note</i>  |
|--------------------------------------|------------|--|
| KVP                                  | 0018,0060  |  |
| Distance Source to Detector          | 0018,1110  |  |
| Exposure Time                        | 0018,1150  |  |
| Exposure                             | 0018,1152  |  |
| Generator Power                      | 0018,1170  |  |
| Sensitivity                          | 0018,6000  |  |
| Preprocessing Function <sup>b</sup>  | 0018,5020  | will contain 8 values separated by commas: <ul style="list-style-type: none"> <li>• Scale per Decade: no meaning for DERIVED\PRIMARY image types</li> <li>• Chamber Average: no meaning for DERIVED\PRIMARY image types</li> <li>• Minimum pixel value</li> <li>• Median value, 1% less than minimum pixel value</li> <li>• Median value, 2% less than minimum pixel value</li> <li>• Median value, 50% less than minimum pixel value</li> <li>• Median value, 98% less than minimum pixel value</li> <li>• Maximum pixel value</li> </ul> |
| Postprocessing Function <sup>b</sup> | 0018,5021  |  |
| Pixel Spacing <sup>b</sup>           | 0028,0030  |  |
| Image Area Dose Product <sup>b</sup> | 0018,115E  |  |
| Pixel Data                           | 7FE0,0010  | between 0 and 30,000   |

a. These two modules are combined because some attributes are present in both modules.

b. Additional attribute in the module; classifies the applied SOP Class as Extended SOP Class.

**Table 8-9: X-Ray Collimator Module**

| <i>Attribute Name</i>            | <i>Tag</i> | <i>Note</i> |
|----------------------------------|------------|-------------|
| Collimator Shape                 | 0018,1700  | RECTANGULAR |
| Collimator Left Vertical Edge    | 0018,1702  |             |
| Collimator Right Vertical Edge   | 0018,1704  |             |
| Collimator Upper Horizontal Edge | 0018,1706  |             |
| Collimator Lower Horizontal Edge | 0018,1708  |             |

**Table 8-10: VOI LUT Module**

| <i>Attribute Name</i> | <i>Tag</i> | <i>Note</i> |
|-----------------------|------------|-------------|
| Window Center         | 0028,1050  | 15,000      |
| Window Width          | 0028,1051  | 30,000      |

**Table 8-11: (Private) Stamp Module**

| <i>Attribute Name</i>        | <i>Tag</i> | <i>Note</i>                                   |
|------------------------------|------------|---|
| Private Creator              | 0089,00xx  | VR is LO<br>VM is 1<br>Value is PMS-THORA-3.1 |
| Stamp Image Sequence         | 0089,xx20  | VR is SQ<br>VM is 1                           |
| > Samples per Pixel          | 0028,0002  | 1   |
| > Photometric Interpretation | 0028,0004  | MONOCHROME1 or MONOCHROME2                    |
| > Rows                       | 0028,0010  |   |
| > Columns                    | 0028,0011  |   |
| > Bits Allocated             | 0028,0100  |   |
| > Bits Stored                | 0028,0101  |   |
| > High Bit                   | 0028,0102  |   |
| > Pixel Representation       | 0028,0103  |   |
| > Pixel Data                 | 7FE0,0010  |   |

**Table 8-12: SOP Common Module**

| <i>Attribute Name</i>  | <i>Tag</i> | <i>Note</i>  |
|------------------------|------------|--|
| SOP Class UID          | 0008,0016  | 1.2.840.10008.5.1.4.1.1.1  |
| SOP Instance UID       | 0008,0018  |  |
| Specific Character Set | 0008,0005  | ISO_IR 100 unless configured on other character set (see "Support of Extended Character Sets" on page 8) |