

# DICOM Conformance Statement

CT Big Bore v4.8



*- This page is left intentionally empty -*


## 1. DICOM Conformance Statement Overview

This version of the DICOM Conformance Statement applies to the “Big Bore” scanners running on the iPatient (4.x) platform.

The systems provides the following DICOM data exchange features:

- Store DICOM Images sent from a Workstation or PACS
- Transfer DICOM Images to a Workstation or PACS
- Query/Retrieve a Workstation or PACS for a list of entries representing Series of DICOM Images
- Query/Retrieve support to let a remote system query for a list of entries representing Series of DICOM Images
- Query a HIS/RIS for a MWL
- Update a remote system with information about Performed Procedure Steps (MPPS)
- Store DICOM Images on portable media (CD, CD-RW, DVD+/-R and DVD+/-RW disks)
- Read DICOM Images from a portable media
- Print Images (Grayscale and Color) on a DICOM Printer

**The system is verified as DIN 6862-2 compliant.**



**Figure 1: Big Bore in a workflow**

The following Table presents an overview of all network services and the applicable SOP Classes as provided by the Big Bore scanner, where the first column specifies the used SOP Classes as named in PS 3.6 (Ref PS 3.6 Annex A) of the current DICOM Standard.

Table 1: Network Services

| SOP Class  |                               | User of Service (SCU) | Provider of Service (SCP) | Display |
|--|-------------------------------|-----------------------|---------------------------|---------|
| Name   | UID                           |                       |                           |         |
| Other  |                               |                       |                           |         |
| Verification SOP Class   | 1.2.840.10008.1.1             | Yes                   | Yes                       | N/A     |
| Print Management   |                               |                       |                           |         |
| Basic Grayscale Print Management Meta SOP Class                  | 1.2.840.10008.5.1.1.9         | Yes                   | No                        | N/A     |
| >Basic Film Box SOP Class  | 1.2.840.10008.5.1.1.2         | Yes                   | No                        | N/A     |
| >Basic Film Session SOP Class                                    | 1.2.840.10008.5.1.1.1         | Yes                   | No                        | N/A     |
| >Basic Grayscale Image Box SOP Class                             | 1.2.840.10008.5.1.1.4         | Yes                   | No                        | N/A     |
| Basic Color Print Management Meta SOP Class                      | 1.2.840.10008.5.1.1.18        | Yes                   | No                        | N/A     |
| >Basic Film Session SOP Class                                    | 1.2.840.10008.5.1.1.1         | Yes                   | No                        | N/A     |
| >Basic Film Box SOP Class  | 1.2.840.10008.5.1.1.2         | Yes                   | No                        | N/A     |
| >Basic Color Image Box SOP Class                                 | 1.2.840.10008.5.1.1.4.1       | Yes                   | No                        | N/A     |
| Query/Retrieve   |                               |                       |                           |         |
| Study Root QR Information Model - FIND SOP Class                 | 1.2.840.10008.5.1.4.1.2.2.1   | Yes                   | Yes                       | N/A     |
| Study Root QR Information Model - MOVE SOP Class                 | 1.2.840.10008.5.1.4.1.2.2.2   | Yes                   | Yes                       | N/A     |
| Transfer   |                               |                       |                           |         |
| Computed Radiography Image Storage SOP Class                     | 1.2.840.10008.5.1.4.1.1.1     | Yes                   | Yes                       | Yes     |
| Digital X-Ray Image Storage - For Pres. SOP Class                | 1.2.840.10008.5.1.4.1.1.1.1   | Yes                   | Yes                       | Yes     |
| Digital X-Ray Image Storage - For Proc. SOP Class                | 1.2.840.10008.5.1.4.1.1.1.1.1 | Yes                   | Yes                       | Yes     |
| X-Ray Angiographic Image Storage SOP Class                       | 1.2.840.10008.5.1.4.1.1.12.1  | Yes                   | Yes                       | Yes     |
| X-Ray Radiofluoroscopic Image Storage SOP Class                  | 1.2.840.10008.5.1.4.1.1.12.2  | Yes                   | Yes                       | Yes     |
| CT Image Storage SOP Class                                       | 1.2.840.10008.5.1.4.1.1.2     | Yes                   | Yes                       | Yes     |
| Secondary Capture Image Storage SOP Class                        | 1.2.840.10008.5.1.4.1.1.7     | Yes                   | Yes                       | Yes     |
| X-Ray Radiation Dose SR SOP Class                                | 1.2.840.10008.5.1.4.1.1.88.67 | Yes                   | Yes                       | Yes     |
| General ECG Waveform Storage SOP Class                           | 1.2.840.10008.5.1.4.1.1.9.1.2 | Yes                   | Yes                       | Yes     |
| MR Image Storage SOP Class                                       | 1.2.840.10008.5.1.4.1.1.4     | Yes                   | Yes                       | Yes     |
| Nuclear Medicine Image Storage SOP Class                         | 1.2.840.10008.5.1.4.1.1.20    | Yes                   | Yes                       | Yes     |
| Digital Intraoral X-Ray Image Storage For Presentation SOP Class | 1.2.840.10008.5.1.4.1.1.1.3   | Yes                   | Yes                       | Yes     |
| Digital Intraoral X-Ray Image Storage For Processing SOP Class   | 1.2.840.10008.5.1.4.1.1.1.3.1 | Yes                   | Yes                       | Yes     |
| Grayscale Presentation SOP Class                                 | 1.2.840.10008.5.1.4.1.1.11.1  | Yes                   | Yes                       | Yes     |
| Color Softcopy Presentation State Storage SOP Class              | 1.2.840.10008.5.1.4.1.1.11.2  | Yes                   | Yes                       | Yes     |
| Pseudo-Color Softcopy Presentation State Storage SOP Class       | 1.2.840.10008.5.1.4.1.1.11.3  | Yes                   | Yes                       | Yes     |

| SOP Class   |                               | User of Service (SCU) | Provider of Service (SCP) | Display |
|---|-------------------------------|-----------------------|---------------------------|---------|
| Name  | UID                           |                       |                           |         |
| Blending Softcopy Presentation State Storage SOP Class              | 1.2.840.10008.5.1.4.1.1.11.4  | Yes                   | Yes                       | Yes     |
| Real World Value Mapping Storage SOP Class                          | 1.2.840.10008.5.1.4.1.1.67    | Yes                   | Yes                       | Yes     |
| Raw Data Storage SOP Class  | 1.2.840.10008.5.1.4.1.1.66    | Yes                   | Yes                       | Yes     |
| Key Object Selection Document Storage SOP Class                     | 1.2.840.10008.5.1.4.1.1.88.59 | Yes                   | Yes                       | Yes     |
| Positron Emission Tomography Image Storage SOP Class                | 1.2.840.10008.5.1.4.1.1.128   | Yes                   | Yes                       | Yes     |
| Digital Mammography X-Ray Image Storage For Presentation SOP Class  | 1.2.840.10008.5.1.4.1.1.1.2   | Yes                   | Yes                       | Yes     |
| Digital Mammography X-Ray Image Storage For Processing SOP Class    | 1.2.840.10008.5.1.4.1.1.1.2.1 | Yes                   | Yes                       | Yes     |
| RT Image Storage SOP Class  | 1.2.840.10008.5.1.4.1.1.481.1 | Yes                   | Yes                       | Yes     |
| RT Dose Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.481.2 | Yes                   | Yes                       | Yes     |
| RT Structure Set Storage SOP Class                                  | 1.2.840.10008.5.1.4.1.1.481.3 | Yes                   | Yes                       | Yes     |
| RT Plan Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.481.5 | Yes                   | Yes                       | Yes     |
| Encapsulated PDF Storage SOP Class                                  | 1.2.840.10008.5.1.4.1.1.104.1 | Yes                   | Yes                       | Yes     |
| Spatial Registration Storage SOP Class                              | 1.2.840.10008.5.1.4.1.1.66.1  | Yes                   | Yes                       | Yes     |
| Multiframe Single Bit Secondary Capture Image Storage SOP Class     | 1.2.840.10008.5.1.4.1.1.7.1   | Yes                   | Yes                       | Yes     |
| Multiframe Grayscale Byte Secondary Capture Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7.2   | Yes                   | Yes                       | Yes     |
| Multiframe Grayscale Word Secondary Capture Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7.3   | Yes                   | Yes                       | Yes     |
| Multiframe True Color Secondary Capture Image Storage SOP Class     | 1.2.840.10008.5.1.4.1.1.7.4   | Yes                   | Yes                       | Yes     |
| <b>Workflow Management</b>  |                               |                       |                           |         |
| Modality Worklist Information Model - FIND SOP Class                | 1.2.840.10008.5.1.4.31        | Yes                   | No                        | N/A     |
| Modality Performed Procedure Step SOP Class                         | 1.2.840.10008.3.1.2.3.3       | Yes                   | No                        | N/A     |
| Storage Commitment Push Model SOP Class                             | 1.2.840.10008.1.20.1          | Yes                   | No                        | N/A     |

**Notes:** Normally the system (SCU) requests only supported DICOM objects. All SOP Classes support the default ILE Transfer Syntax. All transfer Syntaxes are configurable in LAN Config.

Not supported is JPEG transfer syntax for all SOP classes if the IODs have no pixel data and for all non storage SOP classes.

For media the BigBore supports:

- FSC service for CD-R, CD-RW, DVD + R, DVD - R, DVD + RW, DVD - RW media
- FSR service for CD-R, CD-RW, DVD + R, DVD - R, DVD + RW, DVD - RW media

After data is written to DVD, the CD/DVD is finalized; the finalized CD/DVD can now be read on mostly every CD/DVD reader.

All the Media Services supported by BigBore are shown in the next table.

Table 2: Media Services

| Media Storage Application Profile         | File-set Creator (FSC) | File-set Updater (FSU) | File-set Reader (FSR) |
|---|------------------------|------------------------|-----------------------|
| <b>Compact Disk-Recordable</b>            |                        |                        |                       |
| CT/MR Studies on CD-R                     | Yes                    | No                     | Yes                   |
| General Purpose CD-R Interchange          | Yes                    | No                     | Yes                   |
| <b>DVD</b>                                |                        |                        |                       |
| CT/MR Studies on DVD Media                | Yes                    | No                     | Yes                   |
| General Purpose DVD Interchange with JPEG | Yes                    | No                     | Yes                   |

## 2. Contents

|  |    |
|--|----|
| 1. DICOM Conformance Statement Overview .....                          | 3  |
| 2. Contents .....  | 7  |
| 3. Introduction.....   | 9  |
| 3.1. Revision History .....  | 9  |
| 3.2. Audience.....   | 9  |
| 3.3. Remarks.....  | 9  |
| 3.4. Definitions, Terms and Abbreviations .....                        | 10 |
| 3.5. References.....   | 11 |
| 4. Networking.....   | 12 |
| 4.1. Implementation model.....   | 12 |
| 4.1.1. Application Data Flow .....                                     | 12 |
| 4.1.2. Functional Definition of AE's .....                             | 13 |
| 4.1.2.1. Functional Definition of DICOM Manager.....                   | 13 |
| 4.1.2.2. Functional Definition of Print Manager.....                   | 15 |
| 4.1.3. Sequencing of Real World Activities .....                       | 15 |
| 4.2. AE Specifications .....   | 16 |
| 4.2.1. DICOM Manager.....  | 16 |
| 4.2.1.1. SOP Classes.....  | 16 |
| 4.2.1.2. Association Policies .....                                    | 17 |
| 4.2.1.3. Association Initiation Policy .....                           | 18 |
| 4.2.1.4. Association Acceptance Policy .....                           | 42 |
| 4.2.2. Print Manager.....  | 52 |
| 4.2.2.1. SOP Classes.....  | 52 |
| 4.2.2.2. Association Policies .....                                    | 52 |
| 4.2.2.3. Association Initiation Policy .....                           | 53 |
| 4.2.2.4. Association Acceptance Policy .....                           | 63 |
| 4.3. Network Interfaces.....   | 64 |
| 4.3.1. Physical Network Interfaces.....                                | 64 |
| 4.3.2. Additional Protocols .....                                      | 64 |
| 4.3.3. IPv4 and IPv6 Support .....                                     | 64 |
| 4.4. Configuration.....  | 64 |
| 4.4.1. AE Title/Presentation Address Mapping .....                     | 64 |
| 4.4.1.1. Local AE Titles.....  | 64 |
| 4.4.1.2. Remote AE Title/Presentation Address Mapping .....            | 64 |
| 4.4.2. Parameters .....  | 65 |
| 5. Media Interchange .....   | 67 |
| 5.1. Implementation model.....   | 67 |
| 5.1.1. Application Data Flow Diagram .....                             | 67 |
| 5.1.2. Functional Definitions of AE's.....                             | 67 |
| 5.1.3. Sequencing of Real World Activities .....                       | 68 |
| 5.1.4. File Meta Information for Implementation Class and Version..... | 68 |
| 5.2. AE Specifications .....   | 68 |
| 5.2.1. Media - Specification.....                                      | 68 |
| 5.2.1.1. DICOM File Meta Information.....                              | 69 |
| 5.2.1.2. Real-World Activities .....                                   | 72 |
| 5.3. Augmented and Private Application Profiles .....                  | 73 |
| 5.4. Media Configuration .....   | 73 |
| 6. Support of Character Sets .....                                     | 74 |

|   |     |
|---|-----|
| 7. Security .....   | 75  |
| 7.1. Security Profiles .....                                    | 75  |
| 7.1.1. Security use Profiles .....                              | 75  |
| 7.1.2. Security Transport Connection Profiles .....             | 75  |
| 7.1.3. Digital Signature Profiles .....                         | 77  |
| 7.1.4. Media Storage Security Profiles .....                    | 77  |
| 7.1.5. Attribute Confidentiality Profiles .....                 | 77  |
| 7.1.6. Network Address Management Profiles .....                | 79  |
| 7.1.7. Time Synchronization Profiles .....                      | 79  |
| 7.1.8. Application Configuration Management Profiles .....      | 79  |
| 7.1.9. Audit Trail Profiles .....                               | 79  |
| 7.2. Association Level Security .....                           | 80  |
| 7.3. Application Level Security .....                           | 80  |
| 8. Annexes of application "Big Bore v4.8" .....                 | 81  |
| 8.1. IOD Contents .....   | 81  |
| 8.1.1. Created SOP Instance .....                               | 81  |
| 8.1.1.1. List of created SOP Classes .....                      | 83  |
| 8.1.2. Usage of Attributes from Received IOD .....              | 102 |
| 8.1.3. Attribute Mapping .....                                  | 102 |
| 8.1.4. Coerced/Modified fields .....                            | 102 |
| 8.2. Data Dictionary of Private Attributes .....                | 103 |
| 8.3. Coded Terminology and Templates .....                      | 103 |
| 8.3.1. Context Groups .....                                     | 103 |
| 8.3.2. Template Specifications .....                            | 103 |
| 8.3.2.1. CT RADIATION DOSE SR IOD TEMPLATES .....               | 103 |
| 8.3.3. Private code definitions .....                           | 108 |
| 8.4. Grayscale Image consistency .....                          | 108 |
| 8.5. Standard Extended/Specialized/Private SOPs .....           | 109 |
| 8.5.1. Standard Extended/Specialized/Private SOP Instance ..... | 109 |
| 8.5.1.1. CT Image Storage SOP Class .....                       | 109 |
| 8.5.1.2. Secondary Capture Image Storage SOP Class .....        | 109 |
| 8.5.1.3. General ECG Waveform Storage SOP Class .....           | 111 |
| 8.5.1.4. X-Ray Radiation Dose SR SOP Class .....                | 111 |
| 8.6. Private Transfer Syntaxes .....                            | 111 |

### 3. Introduction

#### 3.1. Revision History

The revision history provides dates and differences of the different releases.

**Table 3: Revision History**

| Document Version | Date of Issue | Description of change  |
|------------------|---------------|--|
| 01               | 06-Feb-2020   | Final version<br>Version with DIN 6862-2 updates, supported character set  |
| 02               | 26-Apr-2023   | Updated "General Series Module" Table by removing below attributes <ul style="list-style-type: none"><li>○ Performed Procedure Step Start date (0040,0244)</li><li>○ Performed Procedure Step Start time (0040,0245)</li><li>○ Performed Procedure Step ID (0040,0253)</li><li>○ Performed Procedure Step Description (0040,0254)</li></ul> Updated "Table 44: Requested Query Keys for Study Root Information Model" by removing below attribute in study level. <ul style="list-style-type: none"><li>○ Performed Procedure Step Description (0040,0254)</li></ul> |
| 03               | 18-Dec-2024   | Updated document with Security Profiles and Audit Trail Profiles in Section 7. Security  |

#### 3.2. Audience

This Conformance Statement is intended for:

- (Potential) customers
- System integrators of medical equipment
- Marketing staff interested in system functionality
- Software designers implementing DICOM interfaces

It is assumed that the reader is familiar with the DICOM standard.

#### 3.3. Remarks

The DICOM Conformance Statement is contained in chapter 4 through 8 and follows the contents and structuring requirements of DICOM PS 3.2.

This DICOM Conformance Statement by itself does not guarantee successful interoperability of Philips equipment with non-Philips equipment. The user (or user's agent) should be aware of the following issues:

- **Interoperability**

Interoperability refers to the ability of application functions, distributed over two or more systems, to work successfully together. The integration of medical devices into an IT environment may require application functions that are not specified within the scope of DICOM. Consequently, using only the information provided by this Conformance Statement does not guarantee interoperability of Philips equipment with non-Philips equipment. It is the user's responsibility to analyze thoroughly the application requirements and to specify a solution that integrates Philips equipment with non-Philips equipment.

- **Validation**

Philips equipment has been carefully tested to ensure that the actual implementation of the DICOM interface corresponds with this Conformance Statement. Where Philips equipment is linked to non-Philips equipment, the first step is to compare the relevant Conformance Statements. If the Conformance Statements indicate that successful information exchange should be possible, additional

validation tests will be necessary to ensure the functionality, performance, accuracy and stability of image and image related data. It is the responsibility of the user (or user's agent) to specify the appropriate test suite and to carry out the additional validation tests.

- **New versions of the DICOM Standard**

The DICOM Standard will evolve in future to meet the user's growing requirements and to incorporate new features and technologies. Philips is actively involved in this evolution and plans to adapt its equipment to future versions of the DICOM Standard. In order to do so, Philips reserves the right to make changes to its products or to discontinue its delivery. The user should ensure that any non-Philips provider linking to Philips equipment also adapts to future versions of the DICOM Standard. If not, the incorporation of DICOM enhancements into Philips equipment may lead to loss of connectivity (in case of networking) and incompatibility (in case of media).

The user should ensure that any non-Philips provider linking to Philips equipment also adapts to future versions of the DICOM Standard. If not, the incorporation of DICOM enhancements into Philips equipment may lead to loss of connectivity (in case of networking) and incompatibility (in case of media)

### 3.4. Definitions, Terms and Abbreviations

**Table 4: Definitions, Terms and Abbreviations**

| Abbreviation/Term | Explanation                                    |
|-------------------|--|
| AE                | Application Entity                             |
| ANSI              | American National Standard Institute           |
| AP                | Application Profile                            |
| CD                | Compact Disc                                   |
| CD-R              | CD-Recordable                                  |
| CD-M              | CD-Medical                                     |
| CT                | Computed Tomography                            |
| DICOM             | Digital Imaging and Communications in Medicine |
| DIMSE             | DICOM Message Service Element                  |
| DIMSE-C           | DIMSE-Composite                                |
| DIMSE-N           | DIMSE-Normalized                               |
| DIN               | Deutsches Institut für Normung E.V.            |
| DVD               | Digital Versatile Disc                         |
| EBE               | DICOM Explicit VR Big Endian                   |
| ELE               | DICOM Explicit VR Little Endian                |
| FSC               | File-set Creator                               |
| FSR               | File-set Reader                                |
| FSU               | File-set Updater                               |
| GUI               | Graphic User Interface                         |
| HIS               | Hospital Information System                    |
| HL7               | Health Level Seven                             |
| ILE               | DICOM Implicit VR Little Endian                |
| IOD               | Information Object Definition                  |
| MOD               | Magneto-Optical Disk                           |
| MPPS              | Modality Performed Procedure Step              |
| MR                | Magnetic Resonance                             |
| NEMA              | National Electrical Manufacturers Association  |

| Abbreviation/Term | Explanation   |
|-------------------|---|
| NM                | Nuclear Medicine                                    |
| PDU               | Protocol Data Unit                                  |
| P-ELE             | Private CT Transfer Syntax - Explicit Little Endian |
| RF                | X-Ray Radiofluoroscopic                             |
| RIS               | Radiology Information System                        |
| RT                | Radiotherapy  |
| RWA               | Real-World Activity                                 |
| SC                | Secondary Capture                                   |
| SCP               | Service Class Provider                              |
| SCU               | Service Class User                                  |
| SOP               | Service Object Pair                                 |
| TCP/IP            | Transmission Control Protocol/Internet Protocol     |
| UID               | Unique Identifier                                   |

### 3.5. References

[DICOM] Digital Imaging and Communications in Medicine, Parts 1 - 22 (NEMA PS 3.1- PS 3.22),  
National Electrical Manufacturers Association  
1300 North 17th Street  
Suite 900  
Arlington, Virginia 22209  
Internet: <https://www.dicomstandard.org/current>

## 4. Networking

This section contains the networking related services (vs. the media related ones).

### 4.1. Implementation model

The implementation model consists of three sections:

- The application data flow diagram, specifying the relationship between the Application Entities and the "external world" or Real-World Activities,
- A functional description of each Application Entity, and
- The sequencing constraints among them.

#### 4.1.1. Application Data Flow

The system implements and provides DICOM services using the following Application Entities:

- DICOM-Manager
- Print-Manager

The system consists of three Application Entities. The following figure shows the Networking application data flow as a functional overview of the System. As depicted in the Figure, the System incorporates the following functionality.

- After RWA Create and Set Modality Performed Procedure Step, the System as SCU uses the MPPS Service Class to report the modality performed procedure step.
- After operator RWA Request Modality Worklist, the System as SCU uses the Basic Worklist Management Service Class to request the worklist from a DICOM Radiology information system (RIS).
- After RWA Request Verification, the System as SCP provides standard Verification Service Class functionality to the requesting SCU.
- After RWA Import Images, the System as SCP provides standard Storage Service Class functionality to the requesting SCU.
- After RWA Query Local Images/Retrieve Local Images, the System as SCP provides standard Query/Retrieve Service Class functionality to the requesting SCU.
- After RWA Export Images (triggered by either the operator or RWA Retrieve Local Images), the System as SCU uses the Remote SCP Storage Service Class functionality to store Local Images on a Remote Database.
- After operator RWA Find Remote Images, the System as SCU uses the remote SCP Query/Retrieve Service Class functionality to query remote images.
- After operator RWA Move Remote Images, the System as SCU uses the remote SCP Query/Retrieve Service Class functionality to retrieve remote images.
- After operator RWA Request Storage Commitment, the System as SCU uses the remote SCP Storage Commitment Service Class functionality to commit remote images.
- After operator RWA Print Images, the System as SCU uses the remote Print Management Service Class to print local images.




Figure 2: Network Data Flow Diagram

4.1.2. Functional Definition of AE's

This section contains a functional definition for each individual local Application Entity.

4.1.2.1. Functional Definition of DICOM Manager

The DICOM Manager includes the following service classes.

### **Verification Service Class**

The DICOM Manager provides the Verification service as SCU and SCP.

A remote SCU shall request an association with the DICOM Manager for Verification SOP class. After accepting the association the DICOM Manager shall receive and respond to the Verification request and release the association when requested.

The DICOM Manager can request an association to a remote node for Verification SOP class. After receiving the response for the Verification request from the remote SCP system, it releases the association.

### **Basic Worklist Management Service Class**

The DICOM Manager uses the Basic Worklist Management service as SCU.

After initiating a worklist query the DICOM Manager requests an association with the configured remote Basic Worklist Management SCP. After accepting the association the DICOM Manager shall send the find request, wait for response, and then release the association.

The system shall be updated with the query results.

### **Modality Performed Procedure Step Service Class**

The DICOM Manager as SCU uses the Modality Performed Procedure Step service class to report the status of a procedure step to the configured MPPS manager.

As soon as a study is selected on the scanner and the first acquisition is made, a MPPS N-CREATE message is sent with the status IN PROGRESS to the MPPS manager.

After a worklist is finished on the DICOM Manager scanner (indicated by finishing the study), a new association is opened with the MPPS manager and an N-SET message is sent with the status COMPLETED.

### **Storage Service**

When performing a Storage Service Class (SCP), the DICOM Manager will receive images and store them into the system's local database. The same AE may be used (with a configurable different AE title) to access the local CD/DVD or different local hard disk folders.

### **Storage Commitment Service**

The DICOM Manager is responsible to issue and support the storage commitment service as SCU.

The DICOM Manager establishes association with the specified AE title and sends storage commitment (N-ACTION) request using the push model. After that, it may accept storage commitment (N-EVENT-REPORT) requests on the same association or by establishing another association.

### **Query Retrieve Service**

The DICOM Manager waits for another application to connect at the presentation address configured for its AE title. The DICOM Manager will accept associations with Presentation Contexts for Service Object Pair (SOP) classes for:

- Storage Service Classes (C-STORE)
- Query-Retrieve Service Class (C-MOVE and C-FIND only)
- Verification Service Classes.

When performing Query-Retrieve Service Class (C-FIND SCP), the DICOM Manager will query its local database according to the request's parameters, and will send the results to the issuer.

When performing Query-Retrieve Service Class (C-MOVE SCP), the DICOM Manager will issue a C-STORE (SCU) to the target AE for every image found according to the request.

Export Service

When an object is exported from the local database to an external device, the attributes will be preserved unless an Export Converter is applied.

4.1.2.2. Functional Definition of Print Manager

The Print-Manager is a Graphical User Interface (GUI) based application. It enables the user to print predefined images using the DICOM protocol. The user can specify as a printing destination one of several predefined printers. The user can also modify some of the printing parameters such as the film size and format.

4.1.3. Sequencing of Real World Activities

This section contains description of specific sequencing as well as potential constraints of Real-World Activities, including any applicable user interactions, as performed by the DICOM Manager.




Figure 3: Sequencing for Export Images




Figure 4: Sequencing for Retrieve Local Images

## 4.2. AE Specifications

This section in the DICOM Conformance Statement is a set of Application Entity specifications. There are as many of these subsections as there are different AE's in the implementation.

### 4.2.1. DICOM Manager

Details of this specific Application Entity is specified in this section.

#### 4.2.1.1. SOP Classes

This Application Entity provides Standard Conformance to the following SOP Classes.

**Table 5: SOP Classes for DICOM Manager**

| SOP Class Name   | SOP Class UID                 | SCU | SCP |
|--|-------------------------------|-----|-----|
| Verification SOP Class   | 1.2.840.10008.1.1             | Yes | Yes |
| Storage Commitment Push Model SOP Class                            | 1.2.840.10008.1.20.1          | Yes | No  |
| Modality Performed Procedure Step SOP Class                        | 1.2.840.10008.3.1.2.3.3       | Yes | No  |
| Study Root QR Information Model - FIND SOP Class                   | 1.2.840.10008.5.1.4.1.2.2.1   | Yes | Yes |
| Study Root QR Information Model - MOVE SOP Class                   | 1.2.840.10008.5.1.4.1.2.2.2   | Yes | Yes |
| Modality Worklist Information Model - FIND SOP Class               | 1.2.840.10008.5.1.4.31        | Yes | No  |
| Computed Radiography Image Storage SOP Class                       | 1.2.840.10008.5.1.4.1.1.1     | Yes | Yes |
| Digital X-Ray Image Storage - For Pres. SOP Class                  | 1.2.840.10008.5.1.4.1.1.1.1   | Yes | Yes |
| Digital X-Ray Image Storage - For Proc. SOP Class                  | 1.2.840.10008.5.1.4.1.1.1.1.1 | Yes | Yes |
| X-Ray Angiographic Image Storage SOP Class                         | 1.2.840.10008.5.1.4.1.1.12.1  | Yes | Yes |
| X-Ray Radiofluoroscopic Image Storage SOP Class                    | 1.2.840.10008.5.1.4.1.1.12.2  | Yes | Yes |
| CT Image Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.2     | Yes | Yes |
| Secondary Capture Image Storage SOP Class                          | 1.2.840.10008.5.1.4.1.1.7     | Yes | Yes |
| X-Ray Radiation Dose SR SOP Class                                  | 1.2.840.10008.5.1.4.1.1.88.67 | Yes | Yes |
| General ECG Waveform Storage SOP Class                             | 1.2.840.10008.5.1.4.1.1.9.1.2 | Yes | Yes |
| MR Image Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.4     | Yes | Yes |
| Nuclear Medicine Image Storage SOP Class                           | 1.2.840.10008.5.1.4.1.1.20    | Yes | Yes |
| Digital Intraoral X-Ray Image Storage For Presentation SOP Class   | 1.2.840.10008.5.1.4.1.1.1.3   | Yes | Yes |
| Digital Intraoral X-Ray Image Storage For Processing SOP Class     | 1.2.840.10008.5.1.4.1.1.1.3.1 | Yes | Yes |
| Grayscale Presentation SOP Class                                   | 1.2.840.10008.5.1.4.1.1.11.1  | Yes | Yes |
| Color Softcopy Presentation State Storage SOP Class                | 1.2.840.10008.5.1.4.1.1.11.2  | Yes | Yes |
| Pseudo-Color Softcopy Presentation State Storage SOP Class         | 1.2.840.10008.5.1.4.1.1.11.3  | Yes | Yes |
| Blending Softcopy Presentation State Storage SOP Class             | 1.2.840.10008.5.1.4.1.1.11.4  | Yes | Yes |
| Real World Value Mapping Storage SOP Class                         | 1.2.840.10008.5.1.4.1.1.67    | Yes | Yes |
| Raw Data Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.66    | Yes | Yes |
| Key Object Selection Document Storage SOP Class                    | 1.2.840.10008.5.1.4.1.1.88.59 | Yes | Yes |
| Positron Emission Tomography Image Storage SOP Class               | 1.2.840.10008.5.1.4.1.1.128   | Yes | Yes |
| Digital Mammography X-Ray Image Storage For Presentation SOP Class | 1.2.840.10008.5.1.4.1.1.1.2   | Yes | Yes |

| SOP Class Name  | SOP Class UID                 | SCU | SCP |
|---|-------------------------------|-----|-----|
| Digital Mammography X-Ray Image Storage For Processing              | 1.2.840.10008.5.1.4.1.1.1.2.1 | Yes | Yes |
| RT Image Storage SOP Class  | 1.2.840.10008.5.1.4.1.1.481.1 | Yes | Yes |
| RT Dose Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.481.2 | Yes | Yes |
| RT Structure Set Storage SOP Class                                  | 1.2.840.10008.5.1.4.1.1.481.3 | Yes | Yes |
| RT Plan Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.481.5 | Yes | Yes |
| Encapsulated PDF Storage SOP Class                                  | 1.2.840.10008.5.1.4.1.1.104.1 | Yes | Yes |
| Spatial Registration Storage SOP Class                              | 1.2.840.10008.5.1.4.1.1.66.1  | Yes | Yes |
| Multiframe Single Bit Secondary Capture Image Storage               | 1.2.840.10008.5.1.4.1.1.7.1   | Yes | Yes |
| Multiframe Grayscale Byte Secondary Capture Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7.2   | Yes | Yes |
| Multiframe Grayscale Word Secondary Capture Image Storage           | 1.2.840.10008.5.1.4.1.1.7.3   | Yes | Yes |
| Multiframe True Color Secondary Capture Image Storage               | 1.2.840.10008.5.1.4.1.1.7.4   | Yes | Yes |

**Note:** Any SOP specific behavior is documented later in the conformance statement in the applicable SOP specific conformance section.

#### 4.2.1.2. Association Policies

Each AE specification contains a description of the general association establishment and acceptance policies of the AE.

With incoming association requests the system allows acceptance of a range of defined IP addresses which is configurable in the LanConfig application.

##### 4.2.1.2.1 General

The DICOM standard application context is specified below.

**Table 6: DICOM Application Context**

| Description              | Value                 |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |

##### 4.2.1.2.2 Number of Associations

**Table 7: Number of associations as an Association Initiator for this AE**

| Description                                 | Value   |
|---|---|
| Maximum number of simultaneous associations | "Not configurable" (limited to resource availability) |

**Table 8: Number of associations as an Association Acceptor for this AE**

| Description                                 | Value   |
|---|---|
| Maximum number of simultaneous associations | "Not configurable" (limited to resource availability) |

##### 4.2.1.2.3 Asynchronous Nature

The implementation supports negotiation of multiple outstanding transactions, along with the maximum number of outstanding transactions supported.

Table 9: Asynchronous nature as an Association Initiator for this AE

| Description   | Value |
|---|-------|
| Maximum number of outstanding asynchronous transactions | 1     |

#### 4.2.1.2.4 Implementation Identifying Information

The value supplied for Implementation Class UID and version name are documented here.

Table 10: DICOM Implementation Class and Version for DICOM Manager

|                             |                         |
|-----------------------------|-------------------------|
| Implementation Class UID    | 1.3.46.670589.33.103.11 |
| Implementation Version Name | EBW4.8 for CT           |

#### 4.2.1.2.5 Communication Failure Handling

The behavior of the AE during communication failure is summarized in the next table.

Table 11: Communication Failure Behavior

| Exception           | Behavior  | Comment                         |
|---------------------|---|---------------------------------|
| ARTIM Timeout       | The system stops the ARTIM timer and closes the transport connection. | Configurable, minimum value = 1 |
| Association Timeout | A release request is sent in order to close the association.          | Configurable, minimum value = 1 |

#### 4.2.1.3. Association Initiation Policy

The Application Entity will respond to a received Association rejection as shown in the next table.

Table 12: Association Rejection response

| Result                 | Source  | Reason/Diagnosis                           | Behavior  |
|------------------------|---|--|---|
| 1 - rejected-permanent | 1 - DICOM UL service-user                                     | 1 - no-reason-given                        | Association reject is logged and connection is closed |
|                        |   | 2 - application-context-name-not supported | Association reject is logged and connection is closed |
|                        |   | 3 - calling-AE-title-not-recognized        | Association reject is logged and connection is closed |
|                        |   | 7 - called-AE-title-not-recognized         | Association reject is logged and connection is closed |
|                        | 2 - DICOM UL service-provider (ACSE related function)         | 1 - no-reason-given                        | Association reject is logged and connection is closed |
|                        |   | 2 - protocol-version-not-supported         | Association reject is logged and connection is closed |
|                        | 3 - DICOM UL service-provider (Presentation related function) | 1 - temporary-congestion                   | Association reject is logged and connection is closed |
|                        |   | 2 - local-limit-exceeded                   | Association reject is logged and connection is closed |
| 2 - rejected-transient | 1 - DICOM UL service-user                                     | 1 - no-reason-given                        | Association reject is logged and connection is closed |
|                        |   | 2 - application-context-name-not supported | Association reject is logged and connection is closed |

| Result | Source  | Reason/Diagnosis                    | Behavior  |
|--------|---|-------------------------------------|---|
|        |   | 3 - calling-AE-title-not-recognized | Association reject is logged and connection is closed |
|        |   | 7 - called-AE-title-not-recognized  | Association reject is logged and connection is closed |
|        | 2 - DICOM UL service-provider (ACSE related function)         | 1 - no-reason-given                 | Association reject is logged and connection is closed |
|        |   | 2 - protocol-version-not-supported  | Association reject is logged and connection is closed |
|        | 3 - DICOM UL service-provider (Presentation related function) | 1 - temporary-congestion            | Association reject is logged and connection is closed |
|        |   | 2 - local-limit-exceeded            | Association reject is logged and connection is closed |

The behavior of the AE on receiving an Association abort is summarized in the next table.

**Table 13: Association Abort Handling**

| Source  | Reason/Diagnosis                | Behavior                                  |
|---|---------------------------------|---|
| 0 - DICOM UL service-user (initiated abort)     | 0 - reason-not-specified        | Abort is logged and connection is closed. |
| 2 - DICOM UL service-provider (initiated abort) | 0 - reason-not-specified        | Abort is logged and connection is closed. |
|   | 1 - unrecognized-PDU            | Abort is logged and connection is closed. |
|   | 2 - unexpected-PDU              | Abort is logged and connection is closed. |
|   | 4 - unrecognized-PDU-parameter  | Abort is logged and connection is closed. |
|   | 5 - unexpected-PDU-parameter    | Abort is logged and connection is closed. |
|   | 6 - invalid-PDU-parameter-value | Abort is logged and connection is closed. |

The behavior of the AE for sending an association abort is summarized in next table.

**Table 14: DICOM Association Abort Policies**

| Source                        | Reason/Diagnosis         | Behavior  |
|-------------------------------|--------------------------|---|
| 0 - DICOM UL service-user     | 0 - reason-not-specified | When received the system terminates the connection.                                       |
| 2 - DICOM UL service-provider | 1- unrecognized-PDU      | Whenever the system receives unexpected or unrecognized PDUs it terminates the connection |

#### 4.2.1.3.1 (Real-World) Activity – Verification as SCU

##### 4.2.1.3.1.1 Description and Sequencing of Activities




Figure 5: (Real World) Activity - Verification as SCU

DICOM Manager initiates a C-ECHO when the user points to one of the icons in the devices tool-bar in the UI, clicks the right mouse button and selects "Verify Connection" operation.

##### 4.2.1.3.1.2 Proposed Presentation Contexts

Each time an association is initiated, the association initiator proposes a number of presentation contexts supported for the association. The association will be closed immediately upon receiving the response.

The presentation contexts proposed by DICOM Manager for (Real-World) Activity - Verification as SCU are defined in the following table.

Table 15: Proposed Presentation Contexts for (Real-World) Activity – Verification as SCU

| Presentation Context Table |                   |                           |                     |      |                      |
|----------------------------|-------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax            |                   | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name                       | UID               | Name List                 | UID List            |      |                      |
| Verification SOP Class     | 1.2.840.10008.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None                 |
|                            |                   | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |

**Note:** The default supported Transfer Syntax is ILE. ELE has preference over ILE.

##### 4.2.1.3.1.3 SOP Specific Conformance for Verification SOP Class

This section and sub-section includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

DICOM Manager provides standard conformance to the DICOM 3.0.

##### 4.2.1.3.1.3.1 Dataset Specific Conformance for Verification C-ECHO SCU

Detail regarding the Dataset Specific response behavior will be reported in this section.


This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

**Table 16: Status Response**

| Service Status     | Error Code | Further Meaning                  | Behavior   |
|--------------------|------------|----------------------------------|--|
| Success            | 0000       | Success                          | The SCU has successfully send C-ECHO.                |
| Other than Success | <xxxx>     | Problems with sending the C-ECHO | The SCU failed to send the C-ECHO; user is notified. |

#### 4.2.1.3.2 (Real-World) Activity – Modality Worklist as SCU

##### 4.2.1.3.2.1 Description and Sequencing of Activities



**Figure 6: (Real World) Activity - Modality worklist As SCU**

##### 4.2.1.3.2.2 Proposed Presentation Contexts

The presentation contexts are defined in next table.

**Table 17: Proposed Presentation Contexts for (Real-World) Activity – Modality worklist As SCU**

| Presentation Context Table                           |                        |                           |                     |      |                      |
|--|------------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                                      |                        | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name   | UID                    | Name List                 | UID List            |      |                      |
| Modality Worklist Information Model - FIND SOP Class | 1.2.840.10008.5.1.4.31 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None                 |
|  |                        | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |

##### 4.2.1.3.2.3 SOP Specific Conformance for Modality Worklist Information Model - FIND SOP Class

This section and sub-section include the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

##### 4.2.1.3.2.3.1 Dataset Specific Conformance for Modality Worklist Information Model - FIND SOP Class C-FIND-SCU

Detail regarding the Dataset Specific response behavior will be reported in this section.

The table below should be read as follows:

|                   |   |
|-------------------|---|
| Attribute Name:   | Attributes supported to build a Modality Worklist Request Identifier.   |
| Tag:              | DICOM tag for this attribute.   |
| VR:               | DICOM VR for this attribute.  |
| M:                | Matching Keys. An "X" indicates that this attribute is used for (automatic) Worklist Update.  |
| R:                | Return Keys. An "X" will indicate that this attribute will be supplied as a Return Key with zero length for Universal Matching.   |
| Q:                | Interactive Query Key. An "X" will indicate that Worklist attribute can be used as matching key.  |
| D:                | Displayed Keys. An "X" indicates that this Worklist attribute is displayed to the user during a patient registration dialog.  |
| IOD:              | An "X" indicates that this Worklist attribute is included into all object Instances created during performance of the related Procedure Step.                               |
| Type of matching: | The following types of matching exists:<br>Single Value Matching<br>List of UID Matching<br>Wild Card Matching<br>Range Matching<br>Sequence Matching<br>Universal Matching |

Table 18: Worklist Request Identifier (Patient Query)

| Attribute Name   | Tag       | VR | M | R | Q | D | IOD | Type of Matching                  | Comment |
|--|-----------|----|---|---|---|---|-----|-----------------------------------|---------|
| Patient Identification Module                          |           |    |   |   |   |   |     |                                   |         |
| Other Patient IDs                                      | 0010,1000 | LO |   | X |   |   |     |                                   |         |
| Patient ID   | 0010,0020 | LO | X | X | X | X | X   | Single Value, Wildcard, Universal |         |
| Patient's Name   | 0010,0010 | PN | X | X | X | X | X   | Single Value, Wildcard, Universal |         |
| Patient Demographic Module                             |           |    |   |   |   |   |     |                                   |         |
| Confidentiality Constraint on Patient Data Description | 0040,3001 | LO |   | X |   |   |     |                                   |         |
| Ethnic Group   | 0010,2160 | SH |   | X |   |   |     |                                   |         |
| Patient Comments                                       | 0010,4000 | LT |   | X |   |   | X   |                                   |         |
| Patient's Birth Date                                   | 0010,0030 | DA |   | X | X | X | X   | Universal                         |         |
| Patient's Sex  | 0010,0040 | CS |   | X | X | X | X   | Universal                         |         |
| Patient's Size   | 0010,1020 | DS |   | X |   |   |     |                                   |         |
| Patient's Weight                                       | 0010,1030 | DS |   | X |   | X | X   | Universal                         |         |
| Patient's Primary Language Code Sequence               | 0010,0101 | SQ |   | X |   |   |     |                                   |         |
| >Code Value  | 0008,0100 | SH |   | X |   |   |     |                                   |         |
| >Coding Scheme Designator                              | 0008,0102 | SH |   | X |   |   |     |                                   |         |
| >Code Meaning  | 0008,0104 | LO |   | X |   |   |     |                                   |         |

| Attribute Name                          | Tag        | VR | M | R | Q | D | IOD | Type of Matching | Comment |
|---|------------|----|---|---|---|---|-----|------------------|---------|
| Patient Medical Module                  |            |    |   |   |   |   |     |                  |         |
| Additional Patient History              | 0010,21B0  | LT |   | X |   |   | X   |                  |         |
| Allergies                               | 0010,2110  | LO |   | X |   |   |     |                  |         |
| Medical Alerts                          | 0010,2000  | LO |   | X |   |   |     |                  |         |
| Patient State                           | 0038,0500  | LO |   | X |   |   |     |                  |         |
| Pregnancy Status                        | 0010,21C0  | US |   | X |   |   |     |                  |         |
| Smoking Status                          | 0010,21A0  | CS |   | X |   |   |     |                  |         |
| Special Needs                           | 0038,0050  | LO |   | X |   | X |     |                  |         |
| Visit Relationship Module               |            |    |   |   |   |   |     |                  |         |
| Referenced Patient Sequence             | 0008, 1120 | SQ |   | X |   |   |     |                  |         |
| >Referenced SOP Class UID               | 0008, 1150 | UI |   | X |   |   |     |                  |         |
| >Referenced SOP Instance UID            | 0008, 1155 | UI |   | X |   |   |     |                  |         |
| Visit Identification Module             |            |    |   |   |   |   |     |                  |         |
| Admission ID                            | 0038,0010  | LO |   | X |   |   |     |                  |         |
| Visit Status Module                     |            |    |   |   |   |   |     |                  |         |
| Current Patient Location                | 0038, 0300 | LO |   | X |   |   |     | Universal        |         |
| Visit Comments                          | 0038, 4000 | LT |   | X |   |   |     | Universal        |         |
| Visit Admission Module                  |            |    |   |   |   |   |     |                  |         |
| Admitting Diagnoses Description         | 0008,1080  | LO |   | X |   |   |     |                  |         |
| Referring Physician's Address           | 0008,0092  | ST |   | X |   |   |     |                  |         |
| Referring Physician's Telephone Numbers | 0008,0094  | SH |   | X |   |   |     |                  |         |
| Route of Admissions                     | 0038,0016  | LO |   | X |   |   |     |                  |         |
| Admitting Diagnoses Code Sequence       | 0008,1084  | SQ |   | X |   |   |     |                  |         |
| >Code Value                             | 0008,0100  | SH |   | X |   |   |     |                  |         |
| >Coding Scheme Designator               | 0008,0102  | SH |   | X |   |   |     |                  |         |
| >Code Meaning                           | 0008,0104  | LO |   | X |   |   |     |                  |         |
| SOP Common Module                       |            |    |   |   |   |   |     |                  |         |
| Specific Character Set                  | 0008,0005  | CS |   | X |   |   |     |                  |         |
| Scheduled Procedure Step Module         |            |    |   |   |   |   |     |                  |         |
| Scheduled Procedure Step Sequence       | 0040,0100  | SQ |   | X |   |   |     |                  |         |
| >Modality                               | 0008,0060  | CS | X |   | X | X | X   | Single Value     |         |
| >Pre-Medication                         | 0040,0012  | LO |   | X |   |   |     |                  |         |
| >Requested Contrast Agent               | 0032,1070  | LO |   | X |   |   |     |                  |         |

| Attribute Name                          | Tag       | VR | M | R | Q | D | IOD | Type of Matching    | Comment |
|---|-----------|----|---|---|---|---|-----|---------------------|---------|
| >Scheduled Performing Physician's Name  | 0040,0006 | PN |   | X |   |   |     |                     |         |
| >Scheduled Procedure Step Description   | 0040,0007 | LO |   | X |   |   | X   |                     |         |
| >Scheduled Procedure Step ID            | 0040,0009 | SH |   | X |   |   | X   |                     |         |
| >Scheduled Procedure Step Start Date    | 0040,0002 | DA | X | X | X | X |     | Range, Single Value |         |
| >Scheduled Procedure Step Start Time    | 0040,0003 | TM |   | X |   |   |     |                     |         |
| >Scheduled Procedure Step Status        | 0040,0020 | CS |   | X |   |   |     |                     |         |
| >Scheduled Station AE Title             | 0040,0001 | AE | X |   |   |   |     | Single Value        |         |
| >Scheduled Station Name                 | 0040,0010 | SH |   | X |   |   |     |                     |         |
| >Scheduled Protocol Code Sequence       | 0040,0008 | SQ |   | X |   |   | X   |                     |         |
| >>Code Meaning                          | 0008,0104 | LO |   | X |   |   | X   |                     |         |
| >>Code Value                            | 0008,0100 | SH |   | X |   |   | X   |                     |         |
| >>Coding Scheme Designator              | 0008,0102 | SH |   | X |   |   | X   |                     |         |
| Requested Procedure Module              |           |    |   |   |   |   |     |                     |         |
| Names of Intended Recipients of Results | 0040,1010 | PN |   | X |   |   |     |                     |         |
| Patient Transport Arrangements          | 0040,1004 | LO |   | X |   |   |     |                     |         |
| Requested Procedure Comments            | 0040,1400 | LT |   | X |   |   |     |                     |         |
| Requested Procedure Description         | 0032,1060 | LO |   | X |   |   |     |                     |         |
| Requested Procedure ID                  | 0040,1001 | SH | X | X | X |   | X   | Universal           |         |
| Requested Procedure Priority            | 0040,1003 | SH |   | X |   |   |     |                     |         |
| Study Instance UID                      | 0020,000D | UI |   | X |   |   | X   |                     |         |
| Referenced Study Sequence               | 0008,1110 | SQ |   | X |   |   | X   |                     |         |
| >Referenced SOP Class UID               | 0008,1150 | UI |   | X |   |   | X   |                     |         |
| >Referenced SOP Instance UID            | 0008,1155 | UI |   | X |   |   | X   |                     |         |
| Requested Procedure Code Sequence       | 0032,1064 | SQ |   | X |   |   |     |                     |         |

| Attribute Name                   | Tag       | VR | M | R | Q | D | IOD | Type of Matching | Comment |
|----------------------------------|-----------|----|---|---|---|---|-----|------------------|---------|
| >Code Meaning                    | 0008,0104 | LO |   | X |   |   |     |                  |         |
| >Code Value                      | 0008,0100 | SH |   | X |   |   |     |                  |         |
| >Coding Scheme Designator        | 0008,0102 | SH |   | X |   |   |     |                  |         |
| Imaging Service Request Module   |           |    |   |   |   |   |     |                  |         |
| Accession Number                 | 0008,0050 | SH | X | X | X | X | X   | Universal        |         |
| Imaging Service Request Comments | 0040,2400 | LT |   | X |   |   |     |                  |         |
| Referring Physician's Name       | 0008,0090 | PN |   | X | X | X | X   | Universal        |         |
| Requesting Physician             | 0032,1032 | PN |   | X | X | X |     |                  |         |
| Requesting Service               | 0032,1033 | LO |   | X |   |   |     |                  |         |

#### 4.2.1.3.3 (Real-World) Activity – Modality Performed Procedure Step as SCU

##### 4.2.1.3.3.1 Description and Sequencing of Activities




Figure 7: (Real World) Activity - Modality Performed Procedure Step as SCU

##### 4.2.1.3.3.2 Proposed Presentation Contexts

The presentation contexts are defined in next table.

Table 19: Proposed Presentation Contexts for (Real-World) Activity – Modality Performed Procedure Step As SCU

| Presentation Context Table                  |               |                           |                     |      |                      |
|---|---------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                             |               | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name  | UID           | Name List                 | UID List            |      |                      |
| Modality Performed Procedure Step SOP Class | 1.2.840.10008 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None                 |
|   | .3.1.2.3.3    | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |

#### 4.2.1.3.3.3 SOP Specific Conformance for Modality Performed Procedure Step SOP Class

This section and sub-section include the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

##### 4.2.1.3.3.3.1 Dataset Specific Conformance for Modality Performed Procedure Step SOP Class N-CREATE-SCU

Detail regarding the Dataset Specific response behavior will be reported in this section.

Table 20: MPPS Request Identifiers for N-CREATE-RQ

| Attribute Name                               | Tag       | VR | Value      | Comment |
|--|-----------|----|------------|---------|
| SOP Common Module                            |           |    |            |         |
| Specific Character Set                       | 0008,0005 | CS | ISO_IR 100 |         |
| Performed Procedure Step Relationship Module |           |    |            |         |
| Referenced Patient Sequence                  | 0008,1120 | SQ |            |         |
| >Referenced SOP Class UID                    | 0008,1150 | UI |            |         |
| >Referenced SOP Instance UID                 | 0008,1155 | UI |            |         |
| Patient's Name                               | 0010,0010 | PN |            |         |
| Patient ID                                   | 0010,0020 | LO |            |         |
| Patient's Birth Date                         | 0010,0030 | DA |            |         |
| Patient's Sex                                | 0010,0040 | CS |            |         |
| Scheduled Step Attributes Sequence           | 0040,0270 | SQ |            |         |
| >Accession Number                            | 0008,0050 | SH |            |         |
| >Referenced Study Sequence                   | 0008,1110 | SQ |            |         |
| >>Referenced SOP Class UID                   | 0008,1150 | UI |            |         |
| >>Referenced SOP Instance UID                | 0008,1155 | UI |            |         |
| >Study Instance UID                          | 0020,000D | UI |            |         |
| >Requested Procedure description             | 0032,1060 | LO |            |         |
| >Scheduled Procedure Step description        | 0040,0007 | LO |            |         |
| >Scheduled Protocol Code Sequence            | 0040,0008 | SQ |            |         |
| >>Code Value                                 | 0008,0100 | SH |            |         |
| >>Coding Scheme Designator                   | 0008,0102 | SH |            |         |
| >>Coding Scheme Version                      | 0008,0103 | SH |            |         |
| >Scheduled Procedure Step ID                 | 0040,0009 | SH |            |         |
| >Requested Procedure ID                      | 0040,1001 | SH |            |         |
| Performed Procedure Step Information Module  |           |    |            |         |
| Procedure Code Sequence                      | 0008,1032 | SQ |            |         |

| Attribute Name                               | Tag       | VR | Value       | Comment |
|--|-----------|----|-------------|---------|
| >Code Value                                  | 0008,0100 | SH |             |         |
| >Coding Scheme Designator                    | 0008,0102 | SH |             |         |
| >Coding Meaning                              | 0008,0104 | LO |             |         |
| Performed Station AE Title                   | 0040,0241 | AE |             |         |
| Performed Station Name                       | 0040,0242 | SH |             |         |
| Performed Location                           | 0040,0243 | SH |             |         |
| Performed Procedure Step Start Date          | 0040,0244 | DA |             |         |
| Performed Procedure Step Start Time          | 0040,0245 | TM |             |         |
| Performed Procedure Step End Date            | 0040,0250 | DA |             |         |
| Performed Procedure Step End Time            | 0040,0251 | TM |             |         |
| Performed Procedure Step Status              | 0040,0252 | CS | IN PROGRESS |         |
| Performed Procedure Step ID                  | 0040,0253 | SH |             |         |
| Performed Procedure Step Description         | 0040,0254 | LO |             |         |
| Performed Procedure Type Description         | 0040,0255 | LO |             |         |
| <b>Image Acquisition Results Module</b>      |           |    |             |         |
| Modality                                     | 0008,0060 | CS | CT          |         |
| Study ID                                     | 0020,0010 | SH |             |         |
| Performed Protocol Code Sequence             | 0040,0260 | SQ |             |         |
| >Code Meaning                                | 0008,0104 | LO |             |         |
| >Code Value                                  | 0008,0100 | SH |             |         |
| >Coding Scheme Designator                    | 0008,0102 | SH |             |         |
| >Coding Scheme Meaning                       | 0008,0104 | LO |             |         |
| Performed Series Sequence                    | 0040,0340 | SQ |             |         |
| <b>Radiation Dose Module</b>                 |           |    |             |         |
| Total Number of Exposures                    | 0040,0301 | US |             |         |
| Exposure Dose Sequence                       | 0040,030E | SQ |             |         |
| Comments on Radiation Dose                   | 0040,0310 | ST |             |         |
| <b>Extended Dicom and Private attributes</b> |           |    |             |         |
| Specific Character Set                       | 0008,0005 | CS |             |         |
| Implementor ID                               | 00E1,0010 | LO | ELSCINT1    |         |
| Elscint1_DLP total                           | 00E1,1021 | DS |             |         |

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

**Table 21: Status Response**

| Service Status     | Error Code | Further Meaning              | Behavior  |
|--------------------|------------|------------------------------|---|
| Success            | 0000       | Successful operation         | The SCU successfully sends the modality performed procedure step create request. Log entry. |
| Other than Success | 0106       | Invalid attribute value      | Association with problem will be logged. Association will be released.                      |
|                    | 0116       | Attribute Value Out Of Range |   |
|                    | 0105       | No such Attribute            |   |
|                    | 0107       | Attribute List error         |   |

| Service Status | Error Code | Further Meaning               | Behavior |
|----------------|------------|-------------------------------|----------|
|                | 0110       | MPPS may no longer be updated |          |

#### 4.2.1.3.3.3.2 Dataset Specific Conformance for Modality Performed Procedure Step SOP Class N-SET-SCU

Detail regarding the Dataset Specific response behavior will be reported in this section.

**Table 22: MPPS Request Identifiers for N-SET-RQ**

| Attribute Name  | Tag       | VR | Value | Comment   |
|---|-----------|----|-------|---|
| <b>Performed Procedure Step Information Module</b>    |           |    |       |   |
| Procedure Code Sequence                               | 0008,1032 | SQ |       |   |
| >Code Value   | 0008,0100 | SH |       |   |
| >Coding Scheme Designator                             | 0008,0102 | SH |       |   |
| >Code Meaning   | 0008,0104 | LO |       |   |
| Performed Procedure Step End Date                     | 0040,0250 | DA |       |   |
| Performed Procedure Step End Time                     | 0040,0251 | TM |       |   |
| Performed Procedure Step Status                       | 0040,0252 | CS |       | COMPLETED or DISCONTINUED   |
| Performed Procedure Step Description                  | 0040,0254 | LO |       |   |
| Performed Procedure Type Description                  | 0040,0255 | LO |       |   |
| <b>Image Acquisition Results Module</b>               |           |    |       |   |
| Procedure Code Sequence                               | 0040,0260 | SQ |       |   |
| >Code Value   | 0008,0100 | SH |       |   |
| >Coding Scheme Designator                             | 0008,0102 | SH |       |   |
| >Code Meaning   | 0008,0104 | LO |       |   |
| Performed Series Sequence                             | 0040,0340 | SQ |       |   |
| >Retrieve AE Title                                    | 0008,0054 | AE |       |   |
| >Series Description                                   | 0008,103E | LO |       |   |
| >Performing Physician's Name                          | 0008,1050 | PN |       |   |
| >Operators' Name                                      | 0008,1070 | PN |       | Value entered in the new study window for operator name or in case no value entered, the operator name that is logged in the application is used. |
| >Referenced Image Sequence                            | 0008,1140 | SQ |       |   |
| >>Referenced SOP Class UID                            | 0008,1150 | UI |       |   |
| >>Referenced SOP Instance UID                         | 0008,1155 | UI |       |   |
| >Protocol Name  | 0018,1030 | LO |       |   |
| >Series Instance UID                                  | 0020,000E | UI |       |   |
| >Referenced Non-Image Composite SOP Instance Sequence | 0040,0220 | SQ |       |   |
| <b>Radiation Dose Module</b>                          |           |    |       |   |
| Total Number of Exposures                             | 0040,0301 | US |       |   |
| Exposure Dose Sequence                                | 0040,030E | SQ |       |   |
| >KVP  | 0018,0060 | DS |       |   |

| Attribute Name                               | Tag       | VR | Value    | Comment            |
|--|-----------|----|----------|--------------------|
| >Exposure Time                               | 0018,1150 | IS |          |                    |
| >X-ray Tube Current                          | 0018,1151 | IS |          |                    |
| >Radiation Mode                              | 0018,115A | CS |          |                    |
| >CTDIvol                                     | 0018,9345 | FD |          |                    |
| >Comments on Radiation Dose                  | 0040,0310 | ST |          |                    |
| <b>Extended DICOM and Private attributes</b> |           |    |          |                    |
| Implementor ID                               | 00E1,0010 | LO | ELSCINT1 | Implementor ID     |
| Elscint1_DLP total                           | 00E1,1021 | DS |          | Elscint1_DLP total |

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.


**Table 23: Status Response**

| Service Status     | Error Code | Further Meaning                                  | Behavior  |
|--------------------|------------|--|---|
| Success            | 0000       | Message received successfully by SCP             | Association will be released                          |
| Other than Success | XXXX       | Problems with receiving the N-SET Request by SCP | Problem will be logged. Association will be released. |

#### 4.2.1.3.4 (Real-World) Activity – FIND as SCU

##### 4.2.1.3.4.1 Description and Sequencing of Activities

DICOM-Manager initiates an association when the user clicks on one of the icons in the devices tool-bar. The DICOM-Manager searches (C-FIND) by Study Level following by Series level and, optionally (configurable), by Image Level.



**Figure 8: (Real World) Activity - FIND as SCU**

#### 4.2.1.3.4.2 Proposed Presentation Contexts

Each time an association is initiated, the association initiator proposes a number of presentation contexts to be used on that association. In this subsection, the presentation contexts proposed by Archive-Manager for (Real-World) Activity - Find as SCU are defined in the Table below.

**Table 24: Proposed Presentation Contexts for (Real-World) Activity – FIND as SCU**

| Presentation Context Table                       |                             |                           |                     |      |                      |
|--|-----------------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                                  |                             | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name   | UID                         | Name List                 | UID List            |      |                      |
| Study Root QR Information Model - FIND SOP Class | 1.2.840.10008.5.1.4.1.2.2.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None                 |
|  |                             | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |

#### 4.2.1.3.4.3 SOP Specific Conformance for Study Root QR Information Model – FIND SOP Class

This section and sub-section includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

##### 4.2.1.3.4.3.1 Dataset Specific Conformance for Study Root QR Information Model – FIND SOP Class C-FIND-SCU

Detail regarding the Dataset Specific response behavior will be reported in this section.

**Table 25: Supported Query Keys for Study Root Information Model**

| Study Root Information Model      |           |    |   |                       |
|-----------------------------------|-----------|----|---|-----------------------|
| Attribute Name                    | Tag       | VR | Type of Matching                        | Comment               |
| Query/Retrieve Level              | 0008,0052 | CS | Single Value                            | STUDY, SERIES, IMAGES |
| Specific Character Set            | 0008,0005 | CS |   |                       |
| Q/R Study level                   |           |    |   |                       |
| Study Date                        | 0008,0020 | DA | Range, Universal                        |                       |
| Study Time                        | 0008,0030 | TM | Universal                               |                       |
| Accession Number                  | 0008,0050 | SH | Universal                               |                       |
| Modalities in Study               | 0008,0061 | CS | Single Value, Universal, List Of Values |                       |
| Referring Physician's Name        | 0008,0090 | PN | Universal, Wildcard                     |                       |
| Study Description                 | 0008,1030 | LO | Wildcard, Universal                     |                       |
| Patient's Name                    | 0010,0010 | PN | Universal, Wildcard                     |                       |
| Patient ID                        | 0010,0020 | LO | Single Value, Universal, Wildcard       |                       |
| Patient's Birth Date              | 0010,0030 | DA | Universal                               |                       |
| Patient's Birth Time              | 0010,0032 | TM | Universal                               |                       |
| Patient's Sex                     | 0010,0040 | CS | Single Value, Universal                 |                       |
| Study Instance UID                | 0020,000D | UI | Universal                               |                       |
| Study ID                          | 0020,0010 | SH | Single Value, Universal, Wildcard       |                       |
| Number of Study Related Series    | 0020,1206 | IS | Universal                               |                       |
| Number of Study Related Instances | 0020,1208 | IS | Universal                               |                       |

| Study Root Information Model        |           |    |                  |         |
|-------------------------------------|-----------|----|------------------|---------|
| Attribute Name                      | Tag       | VR | Type of Matching | Comment |
| Q/R Series Level                    |           |    |                  |         |
| Series Date                         | 0008,0021 | DA | Universal        |         |
| Series Time                         | 0008,0031 | TM | Universal        |         |
| Modality                            | 0008,0060 | CS | Universal        |         |
| Manufacturer                        | 0008,0070 | LO | Universal        |         |
| Series Description                  | 0008,103E | LO | Universal        |         |
| Body Part Examined                  | 0018,0015 | CS | Universal        |         |
| Protocol Name                       | 0018,1030 | LO | Universal        |         |
| Study Instance UID                  | 0020,000D | UI | Single Value     |         |
| Series Instance UID                 | 0020,000E | UI | Universal        |         |
| Series Number                       | 0020,0011 | IS | Universal        |         |
| Number of Series Related Instances  | 0020,1209 | IS | Universal        |         |
| Performed Procedure Step Start Date | 0040,0244 | DA | Universal        |         |
| Performed Procedure Step Start Time | 0040,0245 | TM | Universal        |         |
| Request Attributes Sequence         | 0040,0275 | SQ | Universal        |         |
| >Scheduled Procedure Step ID        | 0040,0009 | SH | Universal        |         |
| >Requested Procedure ID             | 0040,1001 | SH | Universal        |         |
| Q/R Image Level                     |           |    |                  |         |
| Image Type                          | 0008,0008 | CS | Universal        |         |
| Instance Creation Date              | 0008,0012 | DA | Universal        |         |
| Instance Creation Time              | 0008,0013 | TM | Universal        |         |
| SOP Class UID                       | 0008,0016 | UI | Universal        |         |
| SOP Instance UID                    | 0008,0018 | UI | Universal        |         |
| Contrast/Bolus Agent                | 0018,0010 | LO | Universal        |         |
| Slice Thickness                     | 0018,0050 | DS | Universal        |         |
| KVP                                 | 0018,0060 | DS | Universal        |         |
| Study Instance UID                  | 0020,000D | UI | Single Value     |         |
| Series Instance UID                 | 0020,000E | UI | Single Value     |         |
| Instance Number                     | 0020,0013 | IS | Universal        |         |
| Patient Orientation                 | 0020,0020 | CS | Universal        |         |
| Image Orientation (Patient)         | 0020,0037 | DS | Universal        |         |
| Frame of Reference UID              | 0020,0052 | UI | Universal        |         |
| Slice Location                      | 0020,1041 | DS | Universal        |         |
| Samples per Pixel                   | 0028,0002 | US | Universal        |         |
| Photometric Interpretation          | 0028,0004 | CS | Universal        |         |
| Rows                                | 0028,0010 | US | Universal        |         |
| Columns                             | 0028,0011 | US | Universal        |         |
| Pixel Spacing                       | 0028,0030 | DS | Universal        |         |

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

Table 26: Status Response

| Service Status | Error Code | Further Meaning                               | Behavior   |
|----------------|------------|---|--|
| Success        | 0000       | Success                                       | Matching is successful. Query results are displayed in UI  |
| Cancel         | FE00       | Matching terminated due to Cancel request     | Error code is logged, no query results are displayed in UI |
| Failure        | A700       | Refused - Out of resources                    | Error code is logged, no query results are displayed in UI |
|                | A900       | Failed - Identifier does not match SOP Class  | Error code is logged, no query results are displayed in UI |
|                | C000       | Failed - Unable to process                    | Error code is logged, no query results are displayed in UI |
|                | C100       | Failed - More than one match found            | Error code is logged, no query results are displayed in UI |
|                | C200       | Failed – Unable to support requested template | Error code is logged, no query results are displayed in UI |

#### 4.2.1.3.5 (Real-World) Activity – MOVE as SCU

##### 4.2.1.3.5.1 Description and Sequencing of Activities

The RWA Move Remote Images involves the retrieve of images on a remote system by moving (copying) the matching images from the remote database to another database.

The operator is able to copy the selected images in a patient folder from a remote database to another, local or remote, database by means of the copy tool in the DICOM Manager data handling facility. The DICOM Manager initiates for each copy request an association to the selected peer entity (Remote AE) and uses it to send the Retrieve (C-MOVE) request (and receive the associated responses). The association is released after the final Retrieve (C-MOVE) response for the related request has been received with the status success / failure.




Figure 9: (Real World) Activity - MOVE as SCU

#### 4.2.1.3.5.2 Proposed Presentation Contexts

Each time an association is initiated, the association initiator proposes a number of presentation contexts to be used on that association.

In this subsection, the presentation contexts proposed by DICOM-Manager for (Real-World) Activity - Move as SCU are defined in the following table.

**Table 27: Proposed Presentation Contexts for (Real-World) Activity – MOVE as SCU**

| Presentation Context Table                       |                    |                           |                     |      |                      |
|--|--------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                                  |                    | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name   | UID                | Name List                 | UID List            |      |                      |
| Study Root QR Information Model - MOVE SOP Class | 1.2.840.10008.5.1. | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None                 |
|  | 4.1.2.2.2          | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |

#### 4.2.1.3.5.3 SOP Specific Conformance for Study Root QR Information Model – MOVE SOP Class

This section and sub-section includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

DICOM-Manager provides standard conformance to the DICOM V3.0 Query/Retrieve Service Class as an SCU for the SOP Class Study Root Query/Retrieve Information Model – Move.

##### 4.2.1.3.5.3.1 Dataset Specific Conformance for Study Root QR Information Model – MOVE SOP Class C-MOVE-SCU

Detail regarding the Dataset Specific response behavior will be reported in this section.

**Table 28: Identifiers for MOVE Study Root Information Model as SCU**

| Study Root Information Model |            |    |                       |
|------------------------------|------------|----|-----------------------|
| Attribute Name               | Tag        | VR | Comment               |
| Query/Retrieve Level         | 0008, 0052 | CS | STUDY, SERIES, IMAGES |
| Q/R Study level              |            |    |                       |
| Study Instance UID           | 0020, 000D | UI |                       |
| Q/R Series level             |            |    |                       |
| Series Instance UID          | 0020, 000E | UI |                       |
| Study Instance UID           | 0020, 000D | UI |                       |
| Q/R Image level              |            |    |                       |
| Series Instance UID          | 0020, 000E | UI |                       |
| SOP Instance UID             | 0008, 0018 | UI |                       |
| Study Instance UID           | 0020, 000D | UI |                       |

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

**Table 29: Status Response**

| Service Status | Error Code | Further Meaning | Behavior               |
|----------------|------------|-----------------|------------------------|
| Success        | 0000       | Success         | Storage is successful. |

| Service Status | Error Code | Further Meaning  | Behavior  |
|----------------|------------|--|---|
| Cancel         | FE00       | Sub-operations terminated due to Cancel indication                 | On receiving status 0xFE00 (Sub-operations terminated due to Cancel indication), the SUT will not display any retrieve results. The reason is logged. Retrieve job is failed in the job viewer.                 |
| Failure        | A701       | Refused - Out of Resources – Unable to calculate number of matches | On receiving status 0xA701 (Refused – Out of Resources – Unable to calculate number of matches), the SUT will not display any retrieve results. The reason is logged. Retrieve job is failed in the job viewer. |
|                | A702       | Refused - Out of Resources   | On receiving status 0xA702 (Refused: Out of Resources – Unable to perform sub-operations), the SUT will not display any retrieve results. The reason is logged. Retrieve job is failed in the job viewer.       |
|                | A801       | Refused - Move Destination Unknown                                 | On receiving status 0xA801 (Refused: Move destination unknown) the SUT will not display any retrieve results. The reason is logged. Retrieve job is failed in the job viewer.                                   |
|                | A900       | Error - Identifier Does Not Match SOP Class                        | On receiving status 0xA900 (Identifier does not match SOP Class), the SUT will not display any retrieve results. The reason is logged. Retrieve job is failed in the job viewer.                                |
|                | Cxxx       | Error - Unable to Process  | On receiving status 0xCxxx (Unable to process), the SUT will not display any retrieve results. The reason is logged. Retrieve job is failed in the job viewer.  |
| Warning        | B000       | Sub-operations complete - One or more failures                     | On receiving status 0xB000 (Sub-operations complete - One or more failures), the SUT will not display any retrieve results. The reason is logged. Retrieve job is failed in the job viewer.                     |

#### 4.2.1.3.6 (Real-World) Activity – Image Export

##### 4.2.1.3.6.1 Description and Sequencing of Activities

The associated Real-World Activity is a request for retrieval of images from the disk or save operation from BigBore applications and storage of the images to a remote system using a C-STORE command.




Figure 10: (Real World) Activity – Image Export

#### 4.2.1.3.6.2 Proposed Presentation Contexts

Each time an association is initiated, the association initiator proposes a number of presentation contexts to be used on that association.

The system configuration (by FSE in LanConfig) allows disabling individually each of the supported syntaxes per remote device when establishing and accepting associations. At least one of the syntaxes will be enabled.

When establishing an association to a remote device, all the enabled (for that device) syntaxes will be proposed based on the SOP Class. If more than one Transfer Syntax is accepted by the SCP, the order of selecting the syntax to use is: P-ELE, JPEG, ELE and ILE. However if the proposed SOP class does not contain any pixel data then only ELE and ILE will be proposed.

When a device is newly added in LanConfig, by default the selected transfer syntaxes for the device are ELE & ILE. If the device is a Philips device which is more recent than BigBore then P-ELE will also selected by default. The user will have the option of changing the selected transfer syntax per device using LanConfig.

BigBore supports Level2 DICOM transparency and hence will preserve all the source image data.

The presentation contexts proposed by the DICOM Manager for (Real-World) Activity (C-STORE SCU) are defined in the following table.

Table 30: Proposed Presentation Contexts for (Real-World) Activity – Image Export

| Presentation Context Table                   |                           |   |                        |      |                      |
|--|---------------------------|---|------------------------|------|----------------------|
| Abstract Syntax                              |                           | Transfer Syntax                                   |                        | Role | Extended Negotiation |
| Name   | UID                       | Name List   | UID List               |      |                      |
| Computed Radiography Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.1 | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 | SCU  | None                 |

| Presentation Context Table                        |  |   |                        |      |                      |
|---|--|---|------------------------|------|----------------------|
| Abstract Syntax                                   |  | Transfer Syntax                                   |                        | Role | Extended Negotiation |
| Name  | UID  | Name List   | UID List               |      |                      |
|   |  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    |      |                      |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Digital X-Ray Image Storage - For Pres. SOP Class | 1.2.840.10008.5.1.4.1.1.11.2.840.10008.5.1.4.1.1.1.1 | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Digital X-Ray Image Storage - For Proc. SOP Class | 1.2.840.10008.5.1.4.1.1.1.1.1                        | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| X-Ray Angiographic Image Storage SOP Class        | 1.2.840.10008.5.1.4.1.1.12.1                         | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| X-Ray Radiofluoroscopic Image Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.12.2                         | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| CT Image Storage SOP Class                        | 1.2.840.10008.5.1.4.1.1.2                            | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Secondary Capture Image Storage SOP Class         | 1.2.840.10008.5.1.4.1.1.7                            | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| X-Ray Radiation Dose SR SOP Class                 | 1.2.840.10008.5.1.4.1.1.88.67                        | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| General ECG Waveform Storage SOP Class            | 1.2.840.10008.5.1.4.1.1.9.1.2                        | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| MR Image Storage SOP Class                        | 1.2.840.10008.5.1.4.1.1.4                            | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Nuclear Medicine Image Storage SOP Class          | 1.2.840.10008.5.1.4.1.1.20                           | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |

| Presentation Context Table   |                               |   |                        |      |                      |
|--|-------------------------------|---|------------------------|------|----------------------|
| Abstract Syntax  |                               | Transfer Syntax                                   |                        | Role | Extended Negotiation |
| Name   | UID                           | Name List   | UID List               |      |                      |
| Digital Intraoral X-Ray Image Storage For Presentation SOP Class   | 1.2.840.10008.5.1.4.1.1.1.3   | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|  |                               | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Digital Intraoral X-Ray Image Storage For Processing SOP Class     | 1.2.840.10008.5.1.4.1.1.1.3.1 | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|  |                               | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Grayscale Presentation SOP Class                                   | 1.2.840.10008.5.1.4.1.1.11.1  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Color Softcopy Presentation State Storage SOP Class                | 1.2.840.10008.5.1.4.1.1.11.2  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Pseudo-Color Softcopy Presentation State Storage SOP Class         | 1.2.840.10008.5.1.4.1.1.11.3  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Blending Softcopy Presentation State Storage SOP Class             | 1.2.840.10008.5.1.4.1.1.11.4  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Real World Value Mapping Storage SOP Class                         | 1.2.840.10008.5.1.4.1.1.67    | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Raw Data Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.66    | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Key Object Selection Document Storage SOP Class                    | 1.2.840.10008.5.1.4.1.1.88.59 | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Positron Emission Tomography Image Storage SOP Class               | 1.2.840.10008.5.1.4.1.1.128   | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|  |                               | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Digital Mammography X-Ray Image Storage For Presentation SOP Class | 1.2.840.10008.5.1.4.1.1.1.2   | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|  |                               | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Digital Mammography X-Ray Image Storage For Processing             | 1.2.840.10008.5.1.4.1.1.1.2.1 | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|  |                               | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |

| Presentation Context Table  |  |   |                        |      |                      |
|---|--|---|------------------------|------|----------------------|
| Abstract Syntax   |  | Transfer Syntax                                   |                        | Role | Extended Negotiation |
| Name  | UID  | Name List   | UID List               |      |                      |
| RT Image Storage SOP Class  | 1.2.840.10008.5.1.4.1.1.481.1                                | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| RT Dose Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.481.2                                | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| RT Structure Set Storage SOP Class                                  | 1.2.840.10008.5.1.4.1.1.481.3                                | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| RT Plan Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.481.5                                | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Encapsulated PDF Storage SOP Class                                  | 1.2.840.10008.5.1.4.1.1.104.1                                | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Spatial Registration Storage SOP Class                              | 1.2.840.10008.5.1.4.1.1.66.1                                 | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Multiframe Single Bit Secondary Capture Image Storage               | 1.2.840.10008.5.1.4.1.1.7.1                                  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Multiframe Grayscale Byte Secondary Capture Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7.2                                  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Multiframe Grayscale Word Secondary Capture Image Storage           | 1.2.840.10008.5.1.4.1.1.7.3                                  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Multiframe True Color Secondary Capture Image Storage               | 1.2.840.10008.5.1.4.1.1.7.4<br>1.2.840.10008.5.1.4.1.1.1.3.1 | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCU  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |

"JPEG" here refers to JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14 [Selection Value 1]): Default Transfer Syntax for Lossless JPEG Image Compression.

It will apply the following priorities to the choice of Transfer Syntax which can be configured in LAN Config Tool.

Table 31: Transfer Syntax Priorities

| Transfer Syntax  | UID                    | Comment   |
|--|------------------------|---|
| 1. CT-private-ELE  | 1.3.46.670589.33.1.4.1 | LanConfig. Default for BigBore Refers to RLE compression transfer syntax. |
| 2. DICOM JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 | LanConfig, Transfer Syntax for Lossless JPEG Image Compression (JPEG).    |
| 3. DICOM Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | LanConfig.  |
| 4. DICOM Implicit VR Little Endian                         | 1.2.840.10008.1.2      | LanConfig, default.   |

**Note:** No support of JPEG & P-ELE transfer syntaxes for all SOP classes without pixel data.

#### 4.2.1.3.6.3 SOP Specific Conformance for Storage SOP Class

DICOM Manager AE provides standard conformance to the DICOM V3.0 Storage Service Class as an SCU for SOP Classes mentioned in the previous section.

Any unsuccessful status (error or warning), returned in the C-STORE Response, results in termination of sending further C-STORE requests (if any in the queue) and reporting of the error to the system log file and UI (Queue Manager)

There are two timeouts for the association. One timeout, "Association Timeout" is used to close an idle association. For C-STORE the default is 120 sec and can be configured per remote DICOM node. The other timeout is "Service Timeout" which detects that no data is transmitted over the association and closes it. The default "Service Timeout" for C-STORE is 5 minutes.

#### 4.2.1.3.6.3.1 Dataset Specific Conformance for Study Root QR Information Model – FIND SOP Class C-FIND-SCU

Detail regarding the Dataset Specific response behavior will be reported in this section.

This includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc

**Table 32: Status Response**

| Service Status | Error Code | Further Meaning                     | Behavior  |
|----------------|------------|-------------------------------------|---|
| Success        | 0000       | Success                             | Export is successful.   |
| Failure        | 0212       | Error - Mistyped Argument           | Export will fail in queue manager and the failure reason is logged. |
|                | A701       | Error - Out of Resources            | Export will fail in queue manager and the failure reason is logged. |
|                | 0210       | Error – Duplicate Invocation        | Export will fail in queue manager and the failure reason is logged. |
|                | 0117       | Error - Invalid Object Instance     | Export will fail in queue manager and the failure reason is logged. |
|                | A900       | Error - Data Set does not match SOP | Export will fail in queue manager and the failure reason is logged. |
|                | C000       | Error - Cannot understand           | Export will fail in queue manager and the failure reason is logged. |
|                | A901       | Error - Data Set Does Not Match     | Export will fail in queue manager and the failure reason is logged. |
| Warning        | B000       | Coercion of Data Elements           | Export is success with a warning reason logged.                     |
|                | B006       | Elements Discarded                  | Export is success with a warning reason logged.                     |
|                | B007       | Data Set does not match             | Export is success with a warning reason logged.                     |
|                | 0107       | Attribute List Error                | Export is success with a warning reason logged.                     |

#### 4.2.1.3.7 (Real-World) Activity – Storage Commitment Push Model AS SCU

##### 4.2.1.3.7.1 Description and Sequencing of Activities

DICOM-Manager will attempt to initiate a new association when requested to commit the images that were stored on a remote device, which supports the storage Commitment Service.

The associated real world activity for the N-ACTION is a storage commitment request to the remote storage device.

The associated real world activity for the N-EVENT-REPORT operation is the completion of the storage commitment by the remote device.

This can be as Synchronous storage commitment as the N-EVENT-REPORT-RQ is received inside the configure timeout or as Asynchronous storage commitment after the Release-RQ by the timeout is already sent to the remote system.

DICOM-Manager will issue a failure status if it is unable to properly handle the storage commitment report event.




Figure 11: (Real World) Activity- DICOM Manager (Storage Commitment)

#### 4.2.1.3.7.2 Proposed Presentation Contexts

Each time an association is initiated, the association initiator proposes a number of presentation contexts to be used on that association.

The presentation contexts proposed by DICOM Manager for (Real-World) Activity - Storage Commitment as SCU are defined in the following table.

Table 33: Proposed Presentation Contexts for (Real-World) Activity – Storage Commitment Push Model AS SCU

| Presentation Context Table              |                      |                           |                     |      |                      |
|---|----------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                         |                      | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name                                    | UID                  | Name List                 | UID List            |      |                      |
| Storage Commitment Push Model SOP Class | 1.2.840.10008.1.20.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None                 |
|   |                      | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |

#### 4.2.1.3.7.3 SOP Specific Conformance for Storage Commitment Push Model SOP Class

DICOM-Manager provides standard conformance to the DICOM V3.0 Storage Commitment Service Class using Push Model as an SCU.

Multiple N-ACTION-RQ can be performed over a single association. Multiple N-EVENT-REPORT-RQ can be accepted over a single association. After all N-ACTION-RQ that are waiting in the stack are issued, association will be closed with the timeout which is configurable using LanConfig.

A remote system reports about storage commitment completion using an N-EVENT-REPORT-RQ command. The system accepts the N-EVENT-REPORT-RQ commands over a separate association initiated by the remote system, using reverse role negotiation for the asynchronous behaviour after the earlier connection was timedout.

Storage Commitment for individual images are grouped into large "chunks" and issued as a single Storage Commitment request.

**Table 34: DICOM Command Communication Failure Behavior Storage Commitment**

| Exception                | Behavior   |
|--------------------------|--|
| ARTIM Time-out           | The reason is logged   |
| Reply Time-out           | The association is released. Continues with waiting for storage commitment |
| Association Time-out SCU | The association is released. Continues with waiting for storage commitment |
| Association aborted      | Continues with waiting for storage commitment                              |

#### 4.2.1.3.7.3.1 Dataset Specific Conformance for Storage Commitment Push Model SOP Class N-ACTION-SCU

All details regarding the specific conformance, including response behavior to all status codes, both from an application level and communication errors are provided in following tables for N-ACTION.

**Table 35: Status Response**

| Service Status     | Error Code | Further Meaning                            | Behavior  |
|--------------------|------------|--|---|
| Success            | 0000       | Success                                    | The request for storage commitment is considered successfully stored. |
| Other than Success | <xxxx>     | Problems with sending the N-ACTION Request | The request for storage commitment is marked as failed.               |

#### 4.2.1.3.7.3.2 Dataset Specific Conformance for Storage Commitment Push Model SOP Class N-EVENT-REPORT-SCP

All details regarding the specific conformance, including response behavior to all status codes, both from an application level and communication errors are provided in following tables for N-EVENT-REPORT.

**Table 36: Storage Commitment N-EVENT-REPORT Behavior**

| Event Type Name                                      | Event Type | Behavior  |
|--|------------|---|
| Storage Commitment Request Successful                | 1          | Awaiting Storage Commitment request will be moved to complete in Queue Manager. |
| Storage Commitment Request Complete - Failures Exist | 2          | Awaiting Storage Commitment request will be moved to fail in Queue Manager.     |

**Table 37: Storage Commitment N-EVENT-REPORT Failure Handling Behavior**

| Service Status | Error Code | Further Meaning | Description   |
|----------------|------------|-----------------|---|
| Success        | 0000       | Success         | The storage commitment result has been successfully received. |

#### 4.2.1.4. Association Acceptance Policy

The Application Entity may reject Association attempts as shown in the table below.

**Table 38: DICOM Association Reject Reasons**

| Result                 | Source  | Reason/Diagnosis                           | Behavior  |
|------------------------|---|--|---|
| 1 - rejected-permanent | 1 - DICOM UL service-user                             | 2 - application-context-name-not supported | Association is rejected and the reason is logged. |
|                        |   | 3 - calling-AE-title-not-recognized        | Association is rejected and the reason is logged. |
|                        |   | 7 - called-AE-title-not-recognized         | Association is rejected and the reason is logged. |
|                        | 2 - DICOM UL service-provider (ACSE related function) | 1 - no-reason-given                        | Association is rejected and the reason is logged. |
|                        |   | 2 - protocol-version-not-supported         | Association is rejected and the reason is logged. |
|                        |   |  |   |

The behavior of the AE on DICOM receiving Association Abort Handling is summarized in table below.

**Table 39: DICOM receiving Association Abort Handling**

| Source  | Reason/Diagnosis                | Behavior                                  |
|---|---------------------------------|---|
| 0 - DICOM UL service-user (initiated abort)     | 0 - reason-not-specified        | Abort is logged and connection is closed. |
| 2 - DICOM UL service-provider (initiated abort) | 0 - reason-not-specified        | Abort is logged and connection is closed. |
|   | 1 - unrecognized-PDU            | Abort is logged and connection is closed. |
|   | 2 - unexpected-PDU              | Abort is logged and connection is closed. |
|   | 4 - unrecognized-PDU parameter  | Abort is logged and connection is closed. |
|   | 5 - unexpected-PDU parameter    | Abort is logged and connection is closed. |
|   | 6 - invalid-PDU-parameter value | Abort is logged and connection is closed. |

The behavior of the AE for sending an association abort is summarized in next table.

**Table 40: Association Abort Policies**

| Source  | Reason/Diagnosis         | Behavior   |
|---|--------------------------|--|
| 0 - DICOM UL service-user (initiated abort)     | 0 - reason-not-specified | When an association timeout (configurable per remote device) expired (timeout which determines how long to keep an idle association).When receiving a PDU whose size is bigger than the agreed max PDU size. |
| 2 - DICOM UL service-provider (initiated abort) | 1 - unrecognized-PDU     | Whenever the system receives unexpected or unrecognized PDU (according to the DICOM UPPER LAYER PROTOCOL STATE TRANSITION TABLE in chapter 8 of the DICOM standard).   |

#### 4.2.1.4.1 (Real-World) Activity – Verification as SCP

##### 4.2.1.4.1.1 Description and Sequencing of Activities

A remote system requests verification from DICOM Manager using the C-ECHO command.




Figure 12: (Real World) Activity - Verification as SCP

#### 4.2.1.4.1.2 Accepted Presentation Contexts

The presentation contexts are defined in next table.

Table 41: Accepted Presentation Contexts for (Real-World) Activity – Verification as SCP

| Presentation Context Table |                   |                           |                     |      |                      |
|----------------------------|-------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax            |                   | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name                       | UID               | Name List                 | UID List            |      |                      |
| Verification SOP Class     | 1.2.840.10008.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCP  | None                 |
|                            |                   | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |

#### 4.2.1.4.1.3 SOP Specific Conformance for Verification SOP Class

This section and sub-section includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

DICOM Manager (C-ECHO SCP) provides standard conformance to the DICOM V3.0 verification SOP Class.

##### 4.2.1.4.1.3.1 Dataset Specific Conformance for Verification C-ECHO SCP

Detail regarding the Dataset Specific response behavior will be reported in this section.

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

Table 42: Status Response

| Service Status | Error Code | Further Meaning | Behavior                                |
|----------------|------------|-----------------|---|
| Success        | 0000       | Success         | C-ECHO command was successful received. |

#### 4.2.1.4.2 (Real-World) Activity – FIND as SCP

##### 4.2.1.4.2.1 Description and Sequencing of Activities

The Real World activity associated with the C-FIND-SCP is querying of the local data base based on C-FIND-RQ from the remote DICOM node. DICOM Manager will issue a failure status if it is unable to process the query request.




Figure 13: (Real World) Activity - FIND as SCP

#### 4.2.1.4.2.2 Accepted Presentation Contexts

The presentation contexts are defined in next table.

Table 43: Acceptable Presentation Contexts for (Real-World) Activity – FIND as SCP

| Presentation Context Table                       |                             |                           |                     |      |                      |
|--|-----------------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                                  |                             | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name   | UID                         | Name List                 | UID List            |      |                      |
| Study Root QR Information Model - FIND SOP Class | 1.2.840.10008.5.1.4.1.2.2.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCP  | None                 |
|  |                             | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |

#### 4.2.1.4.2.3 SOP Specific Conformance for Study Root QR Information Model – FIND SOP Class

This section and sub-section includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

DICOM Manager provides standard conformance to the DICOM V3.0 Query/Retrieve Service Class as an SCP for the following SOP Class: Study Root Query/Retrieve Information Model - FIND, UID=1.2.840.10008.5.1.4.1.2.2.1.

##### 4.2.1.4.2.3.1 Dataset Specific Conformance for Study Root QR Information Model – FIND SOP Class C-FIND-SCP

Detail regarding the Dataset Specific response behavior will be reported in this section.

BigBore does not support Relational Search, a query that may contain any combination of keys at any level in the hierarchy. Starting at the top level in the Query/Retrieve Information Model, continuing until the Query/Retrieve level specified in the C-FIND request is reached.

All Required (R) and Unique (U) Study, Series and Image level keys for the Study Root Query/Retrieve Information Model are supported.

Unsupported fields will not be returned in the C-FIND response

**Table 44: Requested Query Keys for Study Root Information Model**

| Study Root Information Model        |           |    |                                   |         |
|-------------------------------------|-----------|----|-----------------------------------|---------|
| Attribute Name                      | Tag       | VR | Type of Matching                  | Comment |
| Query/Retrieve Level                | 0008,0052 | CS | Universal                         |         |
| Q/R Study level                     |           |    |                                   |         |
| Study Date                          | 0008,0020 | DA | Range, Single Value, Universal    |         |
| Study Time                          | 0008,0030 | TM | Universal                         |         |
| Accession Number                    | 0008,0050 | SH | Single Value                      |         |
| Modalities in Study                 | 0008,0061 | CS | Single Value                      |         |
| Referring Physician's Name          | 0008,0090 | PN | Range, Single Value, Universal    |         |
| Study Description                   | 0008,1030 | LO | Universal                         |         |
| Patient's Name                      | 0010,0010 | PN | Single Value                      |         |
| Patient ID                          | 0010,0020 | LO | Single Value                      |         |
| Patient's Birth Date                | 0010,0030 | DA | Single Value                      |         |
| Patient's Birth Time                | 0010,0032 | TM | Single Value, Universal, Wildcard |         |
| Patient's Sex                       | 0010,0040 | CS | Single Value, Universal, Wildcard |         |
| Study Instance UID                  | 0020,000D | UI | Single Value, Universal, Wildcard |         |
| Study ID                            | 0020,0010 | SH | Single Value, Universal, Wildcard |         |
| Number of Study Related Series      | 0020,1206 | IS | Universal                         |         |
| Number of Study Related Instances   | 0020,1208 | IS | Universal                         |         |
| Q/R Series Level                    |           |    |                                   |         |
| Series Date                         | 0008,0021 | DA | Single Value, Universal           |         |
| Series Time                         | 0008,0031 | TM | Universal                         |         |
| Modality                            | 0008,0060 | CS | Universal                         |         |
| Manufacturer                        | 0008,0070 | LO | Universal                         |         |
| Body Part Examined                  | 0018,0015 | CS | Universal                         |         |
| Protocol Name                       | 0018,1030 | LO | Single Value                      |         |
| Series Description                  | 0008,103E | LO | Universal                         |         |
| Series Number                       | 0020,0011 | IS | Universal                         |         |
| Study Instance UID                  | 0020,000D | UI | Single Value, Universal           |         |
| Series Instance UID                 | 0020,000E | UI | Universal                         |         |
| Number of Series Related Instances  | 0020,1209 | IS | Universal                         |         |
| Performed Procedure Step Start Date | 0040,0244 | DA | Universal                         |         |
| Performed Procedure Step Start Time | 0040,0245 | TM | Single Value                      |         |
| Request Attributes Sequence         | 0040,0275 | SQ | Universal                         |         |

| Study Root Information Model |           |    |                  |         |
|------------------------------|-----------|----|------------------|---------|
| Attribute Name               | Tag       | VR | Type of Matching | Comment |
| Q/R Image Level              |           |    |                  |         |
| Image Type                   | 0008,0008 | CS | Single Value     |         |
| Instance Creation Date       | 0008,0012 | DA | Single Value     |         |
| Instance Creation Time       | 0008,0013 | TM | Single Value     |         |
| SOP Class UID                | 0008,0016 | UI | Single Value     |         |
| SOP Instance UID             | 0008,0018 | UI | Universal        |         |
| Contrast/Bolus Agent         | 0018,0010 | LO | Universal        |         |
| Slice Thickness              | 0018,0050 | DS | Universal        |         |
| KVP                          | 0018,0060 | DS | Universal        |         |
| Study Instance UID           | 0020,000D | UI | Single Value     |         |
| Series Instance UID          | 0020,000E | UI | Universal        |         |
| Instance Number              | 0020,0013 | IS | Single Value     |         |
| Patient Orientation          | 0020,0020 | CS | Universal        |         |
| Image Orientation (Patient)  | 0020,0037 | DS | Universal        |         |
| Frame of Reference UID       | 0020,0052 | UI | Single Value     |         |
| Slice Location               | 0020,1041 | DS | Universal        |         |
| Samples per Pixel            | 0028,0002 | US | Single Value     |         |
| Photometric Interpretation   | 0028,0004 | CS | Single Value     |         |
| Rows                         | 0028,0010 | US | Universal        |         |
| Columns                      | 0028,0011 | US | Universal        |         |
| Pixel Spacing                | 0028,0030 | DS | Universal        |         |

Some C-FIND responses will be forwarded before the C-FIND-CANCEL takes effect.

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

**Table 45: Status Response**

| Service Status | Error Code | Further Meaning        | Behavior                                |
|----------------|------------|------------------------|---|
| Success        | 0000       | Matching complete      | Matching successful.                    |
| Failure        | C000       | General failure status | Whenever the find operation failed.     |
| Cancel         | FE00       | Cancel                 | When receiving a cancel C-FIND request. |
| Pending        | FF00       | Pending                | For every C-FIND response received.     |

#### 4.2.1.4.3 (Real-World) Activity – MOVE as SCP

##### 4.2.1.4.3.1 Description and Sequencing of Activities

The Real World activity associated with the C-MOVE command is retrieval of images from the disk and storage of the images to a remote system using a C-STORE command. DICOM Manager will issue a failure status if it is unable to process the transfer request.




Figure 14: (Real World) Activity - MOVE as SCP

#### 4.2.1.4.3.2 Accepted Presentation Contexts

The presentation contexts are defined in next table.

Table 46: Acceptable Presentation Contexts for (Real-World) Activity – MOVE as SCP

| Presentation Context Table                       |                             |                           |                     |      |                      |
|--|-----------------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                                  |                             | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name   | UID                         | Name List                 | UID List            |      |                      |
| Study Root QR Information Model - MOVE SOP Class | 1.2.840.10008.5.1.4.1.2.2.2 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCP  | None                 |
|  |                             | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |

#### 4.2.1.4.3.3 SOP Specific Conformance for Study Root QR Information Model – MOVE SOP Class

This section and sub-section includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

DICOM Manager provides standard conformance to the DICOM V3.0 Query/Retrieve Service Class as an SCP for the following SOP Class: Study Root Query/Retrieve Information Model - MOVE, UID=1.2.840.10008.5.1.4.1.2.2.2. Prioritization of C-MOVE requests is not supported.

##### 4.2.1.4.3.3.1 Dataset Specific Conformance for Study Root QR Information Model – MOVE SOP Class C-MOVE-SCP

Detail regarding the Dataset Specific response behavior will be reported in this section.

DICOM Manager does not support relational C-MOVE requests. All images requested in the C-MOVE will be sent over a single association.

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

All details regarding the specific conformance, including response behavior of all status codes, both from an application level and communication errors are provided in the following table.

Table 47: Status Response

| Service Status | Error Code | Further Meaning                                | Behavior   |
|----------------|------------|--|--|
| Success        | 0000       | Matching complete                              | Move operation is successful.                                |
| Failure        | A801       | Refused - Move Destination Unknown             | Move operation is failed with a failure message in response. |
|                | C000       | Error - Unable to Process                      | Move operation is failed with a failure message in response. |
| Warning        | B000       | Sub-operations Complete - One or more Failures | One of the move operation is successful and other fails.     |
| Cancel         | FE00       | Cancel   | Move operation is terminated with a status code in response. |

#### 4.2.1.4.4 (Real-World) Activity – Image Import

##### 4.2.1.4.4.1 Description and Sequencing of Activities

The real world activity associated with the C-STORE operation is the storage of the image in the memory of the system upon which DICOM Manager is running in order to make it available for immediate processing by applications. DICOM Manager will issue a failure status if it is unable to store the image in the memory.




Figure 15: (Real World) Activity – Image Import

##### 4.2.1.4.4.2 Accepted Presentation Contexts

The presentation contexts are defined in next table.

Table 48: Accepted Presentation Contexts for (Real-World) Activity – Image Import

| Presentation Context Table                   |                           |   |                        |      |                      |
|--|---------------------------|---|------------------------|------|----------------------|
| Abstract Syntax                              |                           | Transfer Syntax                                   |                        | Role | Extended Negotiation |
| Name   | UID                       | Name List   | UID List               |      |                      |
| Computed Radiography Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.1 | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 | SCP  | None                 |
|  |                           | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    |      |                      |
|  |                           | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|  |                           | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |

| Presentation Context Table                        |  |   |                        |      |                      |
|---|--|---|------------------------|------|----------------------|
| Abstract Syntax                                   |  | Transfer Syntax                                   |                        | Role | Extended Negotiation |
| Name  | UID  | Name List   | UID List               |      |                      |
| Digital X-Ray Image Storage - For Pres. SOP Class | 1.2.840.10008.5.1.4.1.1.11.2.840.10008.5.1.4.1.1.1.1 | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Digital X-Ray Image Storage - For Proc. SOP Class | 1.2.840.10008.5.1.4.1.1.1.1                          | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| X-Ray Angiographic Image Storage SOP Class        | 1.2.840.10008.5.1.4.1.1.12.1                         | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| X-Ray Radiofluoroscopic Image Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.12.2                         | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| CT Image Storage SOP Class                        | 1.2.840.10008.5.1.4.1.1.2                            | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Secondary Capture Image Storage SOP Class         | 1.2.840.10008.5.1.4.1.1.7                            | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| X-Ray Radiation Dose SR SOP Class                 | 1.2.840.10008.5.1.4.1.1.88.67                        | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
| General ECG Waveform Storage SOP Class            | 1.2.840.10008.5.1.4.1.1.9.1.2                        | Implicit VR Little Endian                         | 1.2.840.10008.1.2      | SCP  | None                 |
|   |  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    |      |                      |
| MR Image Storage SOP Class                        | 1.2.840.10008.5.1.4.1.1.4                            | Implicit VR Little Endian                         | 1.2.840.10008.1.2      | SCP  | None                 |
|   |  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Nuclear Medicine Image Storage SOP Class          | 1.2.840.10008.5.1.4.1.1.20                           | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Digital Intraoral X-Ray Image Storage For         | 1.2.840.10008.5.1.4.1.1.1.3                          | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |

| Presentation Context Table   |                               |   |                        |      |                      |
|--|-------------------------------|---|------------------------|------|----------------------|
| Abstract Syntax  |                               | Transfer Syntax                                   |                        | Role | Extended Negotiation |
| Name   | UID                           | Name List   | UID List               |      |                      |
| Presentation SOP Class   |                               | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Digital Intraoral X-Ray Image Storage For Processing SOP Class     | 1.2.840.10008.5.1.4.1.1.1.3.1 | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|  |                               | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Grayscale Presentation SOP Class                                   | 1.2.840.10008.5.1.4.1.1.11.1  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Color Softcopy Presentation State Storage SOP Class                | 1.2.840.10008.5.1.4.1.1.11.2  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Pseudo-Color Softcopy Presentation State Storage SOP Class         | 1.2.840.10008.5.1.4.1.1.11.3  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Blending Softcopy Presentation State Storage SOP Class             | 1.2.840.10008.5.1.4.1.1.11.4  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Real World Value Mapping Storage SOP Class                         | 1.2.840.10008.5.1.4.1.1.67    | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Raw Data Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.66    | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Key Object Selection Document Storage SOP Class                    | 1.2.840.10008.5.1.4.1.1.88.59 | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Positron Emission Tomography Image Storage SOP Class               | 1.2.840.10008.5.1.4.1.1.128   | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|  |                               | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Digital Mammography X-Ray Image Storage For Presentation SOP Class | 1.2.840.10008.5.1.4.1.1.1.2   | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|  |                               | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Digital Mammography X-Ray Image Storage For Processing             | 1.2.840.10008.5.1.4.1.1.1.2.1 | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|  |                               | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| RT Image Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.481.1 | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|  |                               | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |

| Presentation Context Table  |  |   |                        |      |                      |
|---|--|---|------------------------|------|----------------------|
| Abstract Syntax   |  | Transfer Syntax                                   |                        | Role | Extended Negotiation |
| Name  | UID  | Name List   | UID List               |      |                      |
| RT Dose Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.481.2                                | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| RT Structure Set Storage SOP Class                                  | 1.2.840.10008.5.1.4.1.1.481.3                                | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| RT Plan Storage SOP Class   | 1.2.840.10008.5.1.4.1.1.481.5                                | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Encapsulated PDF Storage SOP Class                                  | 1.2.840.10008.5.1.4.1.1.104.1                                | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Spatial Registration Storage SOP Class                              | 1.2.840.10008.5.1.4.1.1.66.1                                 | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Multiframe Single Bit Secondary Capture Image Storage               | 1.2.840.10008.5.1.4.1.1.7.1                                  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Multiframe Grayscale Byte Secondary Capture Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7.2                                  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Multiframe Grayscale Word Secondary Capture Image Storage           | 1.2.840.10008.5.1.4.1.1.7.3                                  | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
| Multiframe True Color Secondary Capture Image Storage               | 1.2.840.10008.5.1.4.1.1.7.4<br>1.2.840.10008.5.1.4.1.1.1.3.1 | Explicit VR Little Endian                         | 1.2.840.10008.1.2.1    | SCP  | None                 |
|   |  | Implicit VR Little Endian                         | 1.2.840.10008.1.2      |      |                      |
|   |  | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |

**Note:** The default supported Transfer Syntax is ILE. All Transfer Syntaxes are configurable in LAN Config, in the order Private-ELE, JPEG, ELE, and ILE. JPEG has preference over ELE and ILE.

For all SOP classes without pixel data the JPEG transfer syntax will not supported

#### 4.2.1.4.4.3 SOP Specific Conformance for Storage SOP Class

This section and sub-section includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

DICOM Manager provides standard conformance to the DICOM V3.0 Storage Service Class as a SCP. DICOM Manager conforms to the SOPs of the Storage Service Class at Level 2 (Full - does not discard any data elements). In case of a successful C-STORE, the stored image may be accessed by the processing applications.

##### 4.2.1.4.4.3.1 Dataset Specific Conformance for C-STORE-RSP

Detail regarding the Dataset Specific response behavior will be reported in this section.

This includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc

**Table 49: Status Response**

| Service Status | Error Code | Further Meaning     | Behavior   |
|----------------|------------|---------------------|--|
| Success        | 0000       | Successfully stored | Import is successful.                              |
| Failure        | <xxxx>     | Failed              | Import fails and the reason for failure is logged. |

#### 4.2.2. Print Manager

Details of this specific Application Entity is specified in this section.

##### 4.2.2.1. SOP Classes

This Application Entity provides Standard Conformance to the following SOP Classes.

**Table 50: SOP Classes for Print Manager**

| SOP Class Name                                  | SOP Class UID           | SCU | SCP |
|---|-------------------------|-----|-----|
| Basic Color Print Management Meta SOP Class     | 1.2.840.10008.5.1.1.18  | Yes | No  |
| >Basic Film Session SOP Class                   | 1.2.840.10008.5.1.1.1   | Yes | No  |
| >Basic Film Box SOP Class                       | 1.2.840.10008.5.1.1.2   | Yes | No  |
| >Basic Color Image Box SOP Class                | 1.2.840.10008.5.1.1.4.1 | Yes | No  |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9   | Yes | No  |
| >Basic Film Session SOP Class                   | 1.2.840.10008.5.1.1.1   | Yes | No  |
| >Basic Film Box SOP Class                       | 1.2.840.10008.5.1.1.2   | Yes | No  |
| >Basic Grayscale Image Box SOP Class            | 1.2.840.10008.5.1.1.4   | Yes | No  |

**Note:** Any SOP specific behavior is documented later in the conformance statement in the applicable SOP specific conformance section.

**Note:** The Printer SOP Class is not supported, even though it is included in the table above.

##### 4.2.2.2. Association Policies

Each AE specification contains a description of the general association establishment and acceptance policies of the AE.

###### 4.2.2.2.1 General

The maximum PDU Size that the Print-Manager will use is configurable, with a minimum of 2 Kbytes.

**Table 51: DICOM Application Context**

| Description              | Value                 |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |

###### 4.2.2.2.2 Number of Associations

The number of simultaneous associations that an Application Entity may support as an Initiator or Acceptor is specified.

Print-Manager can have only one open connection at a given time.

**Table 52: Number of associations as an Association Initiator for this AE**

| Description                                 | Value |
|---|-------|
| Maximum number of simultaneous associations | 1     |

#### 4.2.2.2.3 Asynchronous Nature

The implementation supports negotiation of multiple outstanding transactions, along with the maximum number of outstanding transactions supported.

Print-Manager will only allow a single outstanding operation on an association.

**Table 53: Asynchronous nature as an Association Initiator for this AE**

| Description   | Value |
|---|-------|
| Maximum number of outstanding asynchronous transactions | 1     |

#### 4.2.2.2.4 Implementation Identifying Information

The value supplied for Implementation Class UID and version name are documented here.

**Table 54: DICOM Implementation Class and Version for Print Manager**

|                             |                         |
|-----------------------------|-------------------------|
| Implementation Class UID    | 1.3.46.670589.33.103.11 |
| Implementation Version Name | EBW4.8 for CT           |

#### 4.2.2.2.5 Communication Failure Handling

Not applicable. Big Bore System has not implemented communication failure handling for Print Service.

#### 4.2.2.3. Association Initiation Policy

The Application Entity will respond to a received Association rejection as shown in the next table.

**Table 55: Association Rejection response**

| Result                 | Source  | Reason/Diagnosis                           | Behavior                  |
|------------------------|---|--|---------------------------|
| 1 - rejected-permanent | 1 - DICOM UL service-user                                     | 1 - no-reason-given                        | The connection is closed. |
|                        |   | 2 - application-context-name-not supported | The connection is closed. |
|                        |   | 3 - calling-AE-title-not-recognized        | The connection is closed. |
|                        |   | 7 - called-AE-title-not-recognized         | The connection is closed. |
|                        | 2 - DICOM UL service-provider (ACSE related function)         | 1 - no-reason-given                        | The connection is closed. |
|                        |   | 2 - protocol-version-not-supported         | The connection is closed. |
|                        | 3 - DICOM UL service-provider (Presentation related function) | 1 - temporary-congestion                   | The connection is closed. |
|                        |   | 2 - local-limit-exceeded                   | The connection is closed. |
| 2 - rejected-transient | 1 - DICOM UL service-user                                     | 1 - no-reason-given                        | The connection is closed. |
|                        |   | 2 - application-context-name-not supported | The connection is closed. |
|                        |   | 3 - calling-AE-title-not-recognized        | The connection is closed. |
|                        |   | 7 - called-AE-title-not-recognized         | The connection is closed. |
|                        |   | 1 - no-reason-given                        | The connection is closed. |

| Result | Source  | Reason/Diagnosis                   | Behavior                  |
|--------|---|------------------------------------|---------------------------|
|        | 2 - DICOM UL service-providerz (ACSE related function)        | 2 - protocol-version-not-supported | The connection is closed. |
|        | 3 - DICOM UL service-provider (Presentation related function) | 1 - temporary-congestion           | The connection is closed. |
|        |   | 2 - local-limit-exceeded           | The connection is closed. |

The behavior of the AE on receiving an association abort is summarized in next table.


**Table 56: Association Abort Handling**

| Source  | Reason/Diagnosis                | Behavior                  |
|---|---------------------------------|---------------------------|
| 0 - DICOM UL service-user (initiated abort)     | 0 - reason-not-specified        | The connection is closed. |
| 2 - DICOM UL service-provider (initiated abort) | 0 - reason-not-specified        | The connection is closed. |
|   | 1 - unrecognized-PDU            | The connection is closed. |
|   | 2 - unexpected-PDU              | The connection is closed. |
|   | 4 - unrecognized-PDU-parameter  | The connection is closed. |
|   | 5 - unexpected-PDU-parameter    | The connection is closed. |
|   | 6 - invalid-PDU-parameter-value | The connection is closed. |

#### 4.2.2.3.1 (Real-World) Activity – Print Management as SCU

##### 4.2.2.3.1.1 Description and Sequencing of Activities

After selecting the print destination (out of choice list of configured printers) and some print parameters (depending on the configuration and the selected printer; these values can be configured too), Print-Manager initiates an association when a print job is submitted to a DICOM printer (when the user clicks on the print button in the film view). Jobs to different printers are performed simultaneously. Normally, when the job is completed and there are no other jobs to the same printer, the Print manager closes the association with an A-RELEASE request.



**Figure 16: (Real World) Activity – Print Manager Initiates**




Figure 17: (Real World) Activity – Print Management as SCU

#### 4.2.2.3.1.2 Proposed Presentation Contexts

The presentation contexts are defined in next table.

Table 57: Proposed Presentation Contexts for (Real-World) Activity – Print Management As SCU

| Presentation Context Table                      |                         |                           |                     |      |                      |
|---|-------------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                                 |                         | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name  | UID                     | Name List                 | UID List            |      |                      |
| Basic Color Print Management Meta SOP Class     | 1.2.840.10008.5.1.1.18  |                           |                     | SCU  | None                 |
| >Basic Film Session SOP Class                   | 1.2.840.10008.5.1.1.1   | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None                 |
|   |                         | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |
| >Basic Film Box SOP Class                       | 1.2.840.10008.5.1.1.2   | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None                 |
|   |                         | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |
| >Basic Color Image Box SOP Class                | 1.2.840.10008.5.1.1.4.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None                 |
|   |                         | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9   |                           |                     | SCU  | None                 |
|   |                         | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None                 |

| Presentation Context Table           |                       |                           |                     |      |                      |
|--------------------------------------|-----------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                      |                       | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name                                 | UID                   | Name List                 | UID List            |      |                      |
| >Basic Film Session SOP Class        | 1.2.840.10008.5.1.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |
| >Basic Film Box SOP Class            | 1.2.840.10008.5.1.1.2 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None                 |
|                                      |                       | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |
| >Basic Grayscale Image Box SOP Class | 1.2.840.10008.5.1.1.4 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None                 |
|                                      |                       | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |

Note: The only supported Transfer Syntaxes for printing are ELE and ILE.

This section specifies each IOD created (including private IOD's).

Abbreviations used in the Module table for the column "Presence of Value" are:

|        |  |
|--------|--|
| ALWAYS | The attribute is always present with a value   |
| EMPTY  | The attribute is always present without any value (attribute sent zero length)   |
| VNAP   | The attribute is always present and its Value is Not Always Present (attribute sent zero length if no value is present)  |
| ANAP   | The attribute is present under specified condition – if present then it will always have a value   |
| VNAP   | The attribute is present under specified condition – if present then its Value is Not Always Present (attribute sent zero length if condition applies and no value is present) |
| ANAP   | The attribute is present under specified condition – if present then it will not have any value  |

The abbreviations used in the Module table for the column "Source" are:

|          |   |
|----------|---|
| AUTO     | The attribute value is generated automatically                                    |
| CONFIG   | The attribute value source is a configurable parameter                            |
| COPY     | The attribute value source is another SOP instance                                |
| FIXED    | The attribute value is hard-coded in the application                              |
| IMPLICIT | The attribute value source is a user-implicit setting                             |
| MPPS     | The attribute value is the same as that use for Modality Performed Procedure Step |
| MWL      | The attribute value source is a Modality Worklist                                 |
| USER     | The attribute value source is explicit user input                                 |

#### 4.2.2.3.1.3 SOP Specific Conformance for Basic Film Session SOP Class of the Basic Color Print Management Meta SOP Class

This section and sub-section includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

The Printer process conforms to the Basic Film Session SOP Class. The following DIMSE service elements are supported:

- N-CREATE
- N-DELETE

##### 4.2.2.3.1.3.1 Dataset Specific Conformance for Basic Film Session SOP Class N-CREATE-SCU

Details regarding the Dataset Specific response behavior will be reported in this section.

Table 58: Basic Film Session Presentation Module

| Attribute Name   | Tag       | VR | Value                        | Presence of Value | Source       | Comments            |
|------------------|-----------|----|------------------------------|-------------------|--------------|---------------------|
| Number of Copies | 2000,0010 | IS |                              | ALWAYS            | USER         | Value between 1-100 |
| Print Priority   | 2000,0020 | CS | AUTO, HIGH, LOW, MED         | ANAP              | USER         |                     |
| Medium Type      | 2000,0030 | CS | BLUE FILM, CLEAR FILM, PAPER | ALWAYS            | CONFIG, USER |                     |
| Film Destination | 2000,0040 | CS | MAGAZINE, PROCESSOR          | ALWAYS            | CONFIG, USER |                     |

Table 59: Status Response

| Service Status | Error Code | Further Meaning     | Behavior   |
|----------------|------------|---------------------|--|
| Success        | 0000       | Success             | Successful print job status is displayed on UI.      |
| Warning        | <xxxx>     | All warning numbers | Print job continues with a warning code logged.      |
| Failure        | <xxxx>     | All error numbers   | Print job is terminated and error message is logged. |

#### 4.2.2.3.1.3.2 Dataset Specific Conformance for Basic Film Session SOP Class N-DELETE-SCU

Detail regarding the Dataset Specific response behavior will be reported in this section.

Table 60: Status Response

| Service Status | Error Code | Further Meaning                   | Behavior   |
|----------------|------------|-----------------------------------|--|
| Success        | 0000       | Film Session successfully created | Successful print job status is displayed on UI.      |
| Failure        | <xxxx>     | Error code                        | Print job is terminated and error message is logged. |

#### 4.2.2.3.1.4 SOP Specific Conformance for Basic Film Box SOP Class of the Basic Color Print Management Meta SOP Class

This section and sub-sections includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

The Printer process conforms to the Basic Film Box Sop Class. The following DIMSE service elements are supported:

- N-CREATE
- N-ACTION
- N-DELETE

#### 4.2.2.3.1.4.1 Dataset Specific Conformance for Basic Film Box SOP Class N-CREATE-SCU

Table 61: Basic Film Box Presentation Module

| Attribute Name       | Tag       | VR | Value               | Presence of Value | Source       | Comments                          |
|----------------------|-----------|----|---------------------|-------------------|--------------|-----------------------------------|
| Image Display Format | 2010,0010 | ST | STANDARD\1,1        | ALWAYS            | CONFIG       |                                   |
| Film Orientation     | 2010,0040 | CS | LANDSCAPE, PORTRAIT | ALWAYS            | CONFIG, USER |                                   |
| Film Size ID         | 2010,0050 | CS |                     | ALWAYS            | CONFIG, USER | As in printer configuration file. |

| Attribute Name                   | Tag       | VR | Value   | Presence of Value | Source       | Comments                          |
|----------------------------------|-----------|----|---------|-------------------|--------------|-----------------------------------|
| Trim                             | 2010,0140 | CS | NO, YES | ALWAYS            | CONFIG, USER |                                   |
| Magnification Type               | 2010,0060 | CS |         | ALWAYS            | CONFIG       | As in printer configuration file. |
| Min Density                      | 2010,0120 | US |         | ALWAYS            | CONFIG       | As in printer configuration file. |
| Trim                             | 2010,0140 | CS | NO, YES | ALWAYS            | CONFIG, USER |                                   |
| Configuration Information        | 2010,0150 | ST |         | ALWAYS            | CONFIG       | As in printer configuration file. |
| Referenced Film Session Sequence | 2010,0500 | SQ |         | ALWAYS            | AUTO         |                                   |
| >Referenced SOP Class UID        | 0008,1150 | UI |         | ALWAYS            | AUTO         | UID of Parent Film Session.       |
| >Referenced SOP Instance UID     | 0008,1155 | UI |         | ALWAYS            | AUTO         |                                   |

Table 62: Basic Film Box Relationship Module

| Attribute Name                   | Tag       | VR | Value | Presence of Value | Source | Comments                    |
|----------------------------------|-----------|----|-------|-------------------|--------|-----------------------------|
| Referenced Film Session Sequence | 2010,0500 | SQ |       | ALWAYS            | AUTO   |                             |
| >Referenced SOP Class UID        | 0008,1150 | UI |       | ALWAYS            | AUTO   | UID of Parent Film Session. |
| >Referenced SOP Instance UID     | 0008,1155 | UI |       | ALWAYS            | AUTO   |                             |

Table 63: Status Response

| Service Status | Error Code | Further Meaning               | Behavior   |
|----------------|------------|-------------------------------|--|
| Success        | 0000       | Film Box successfully created | Operation is completed successfully                  |
| Warning        | 0116       | Attributes out of range       | Print job continues with a warning code logged.      |
| Failure        | 0106       | Invalid Attribute Value       | Print job is terminated and error message is logged. |
|                | 0105       | No Such Attribute             | Print job is terminated and error message is logged. |

#### 4.2.2.3.1.4.2 Dataset Specific Conformance for Basic Film Box SOP Class N-ACTION-SCU

Table 64: DICOM Command Response Status Handling Behavior for Basic Film Box N-ACTION

| Service Status | Error Code | Further Meaning  | Behavior   |
|----------------|------------|--|--|
| Success        | 0000       | Film accepted for printing   | Operation is completed successfully.                 |
| Failure        | 0119       | Specified SOP Instance is not a member of the specified SOP class            | Print job is terminated and error message is logged. |
| Warning        | B603       | Film Box SOP Instance hierarchy does not contain Image Box SOP Instances     | Print job continues with a warning code logged.      |
|                | B604       | Image Size is larger than Image Box Size - The Image has been de-magnified   | Print job continues with a warning code logged.      |
|                | B609       | Image Size is larger than Image Box Size - The Image has been cropped to fit | Print job continues with a warning code logged.      |

| Service Status | Error Code | Further Meaning   | Behavior  |
|----------------|------------|---|---|
|                | B60A       | Image Size or combined Print Image Size is larger than Image Box Size - The Image or combined Print Image has been decimated to fit | Print job continues with a warning code logged. |

#### 4.2.2.3.1.4.3 Dataset Specific Conformance for Basic Film Box SOP Class N-DELETE-SCU

Table 65: Status Response

| Service Status     | Error Code | Further Meaning                   | Behavior  |
|--------------------|------------|-----------------------------------|---|
| Success            | 0000       | Film Session successfully created | Operation is completed successfully.                      |
| Other than Success | <xxxx>     | Any other status then success     | Print job remains in the queue manager with status failed |

#### 4.2.2.3.1.5 SOP Specific Conformance for Basic Color Image Box SOP Class of the Basic Color Print Management Meta SOP Class

This section and sub-section includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

The Printer process conforms to the Color Image Box Sop Class. The following DIMSE service element is supported: N-SET.

##### 4.2.2.3.1.5.1 Dataset Specific Conformance for Basic Color Image Box SOP Class N-SET-SCU

Details regarding the Dataset Specific response behavior will be reported in this section.

Table 66: Image Box Pixel Presentation Module

| Attribute Name              | Tag       | VR    | Value | Presence of Value | Source | Comments                          |
|-----------------------------|-----------|-------|-------|-------------------|--------|-----------------------------------|
| Image Box Position          | 2020,0010 | US    | 1     | ALWAYS            | AUTO   |                                   |
| Basic Color Image Sequence  | 2020,0111 | SQ    |       | ALWAYS            | AUTO   |                                   |
| >Samples per Pixel          | 0028,0002 | US    | 3     | ALWAYS            | AUTO   |                                   |
| >Photometric Interpretation | 0028,0004 | CS    | RGB   | ALWAYS            | AUTO   |                                   |
| >Planar Configuration       | 0028,0006 | US    | 0, 1  | ALWAYS            | AUTO   |                                   |
| >Rows                       | 0028,0010 | US    |       | ALWAYS            | AUTO   | As in printer configuration file. |
| >Columns                    | 0028,0011 | US    |       | ALWAYS            | AUTO   | As in printer configuration file. |
| >Pixel Aspect Ratio         | 0028,0034 | IS    |       | ALWAYS            | AUTO   | Must be present if not 1/1.       |
| >Bits Allocated             | 0028,0100 | US    | 8     | ALWAYS            | AUTO   |                                   |
| >Bits Stored                | 0028,0101 | US    | 8     | ALWAYS            | AUTO   |                                   |
| >High Bit                   | 0028,0102 | US    | 7     | ALWAYS            | AUTO   |                                   |
| >Pixel Representation       | 0028,0103 | US    | 0     | ALWAYS            | AUTO   |                                   |
| >Pixel Data                 | 7FE0,0010 | OW/OB |       | ALWAYS            | AUTO   |                                   |

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

**Table 67: DICOM Command Response Status Handling Behavior for Basic Color Image Box N-SET**

| Service Status | Error Code | Further Meaning                        | Behavior   |
|----------------|------------|--|--|
| Success        | 0000       | Image successfully stored in Image Box | Operation is completed successfully.                 |
| Warning        | 0107       | Attribute list error                   | Print job continues with a warning code logged.      |
| Failure        | 0106       | Invalid Attribute Value                | Print job is terminated and error message is logged. |
|                | 0105       | No Such Attribute                      | Print job is terminated and error message is logged. |

#### 4.2.2.3.1.6 SOP Specific Conformance for Basic Film Session SOP Class of the Basic Grayscale Print Management Meta SOP Class

This section and sub-section includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

The Printer process conforms to the Basic Film Session SOP Class. The following DIMSE service element is supported:

- N-CREATE
- N-DELETE

##### 4.2.2.3.1.6.1 Dataset Specific Conformance for Basic Film Session SOP Class N-CREATE-SCU

**Table 68: Basic Film Session Presentation Module**

| Attribute Name   | Tag       | VR | Value                        | Presence of Value | Source       | Comments            |
|------------------|-----------|----|------------------------------|-------------------|--------------|---------------------|
| Number of Copies | 2000,0010 | IS |                              | ALWAYS            | CONFIG, USER | Value between 1-100 |
| Print Priority   | 2000,0020 | CS | AUTO, HIGH, LOW, MED         | ALWAYS            | CONFIG, USER |                     |
| Medium Type      | 2000,0030 | CS | BLUE FILM, CLEAR FILM, PAPER | ALWAYS            | USER         |                     |
| Film Destination | 2000,0040 | CS | MAGAZINE, PROCESSOR          | ANAP              | USER         |                     |

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

**Table 69: Status Response**

| Service Status | Error Code | Further Meaning     | Behavior   |
|----------------|------------|---------------------|--|
| Success        | 0000       | Success             | Successful print job status is displayed on UI.      |
| Warning        | <xxxx>     | All warning numbers | Print job continues with a warning code logged.      |
| Failure        | <xxxx>     | All error numbers   | Print job is terminated and error message is logged. |

##### 4.2.2.3.1.6.2 Dataset Specific Conformance for Basic Film Session SOP Class N-DELETE-SCU

Detail regarding the Dataset Specific response behavior will be reported in this section.

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

**Table 70: Status Response**

| Service Status | Error Code | Further Meaning                   | Behavior  |
|----------------|------------|-----------------------------------|---|
| Success        | 0000       | Film Session successfully created | Successful print job status is displayed on UI. |

| Service Status | Error Code | Further Meaning | Behavior   |
|----------------|------------|-----------------|--|
| Failure        | <xxxx>     | Error code      | Print job is terminated and error message is logged. |

#### 4.2.2.3.1.7 SOP Specific Conformance for Basic Film Box SOP Class of the Basic Grayscale Print Management Meta SOP Class

This section and sub-sections includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

The Printer process conforms to the Basic Film Box Sop Class. The following DIMSE service elements are supported: N-CREATE, N-ACTION, and N-DELETE.

##### 4.2.2.3.1.7.1 Dataset Specific Conformance for Basic Film Box SOP Class N-CREATE-SCU

Detail regarding the Dataset Specific response behavior will be reported in this section.

Table 71: Basic Film Box Presentation Module

| Attribute Name            | Tag       | VR | Value               | Presence of Value | Source       | Comment                           |
|---------------------------|-----------|----|---------------------|-------------------|--------------|-----------------------------------|
| Image Display Format      | 2010,0010 | ST | STANDARD\1,1        | ALWAYS            | CONFIG       |                                   |
| Film Orientation          | 2010,0040 | CS | LANDSCAPE, PORTRAIT | ALWAYS            | CONFIG, USER |                                   |
| Film Size ID              | 2010,0050 | CS |                     | ALWAYS            | CONFIG, USER | As in printer configuration file. |
| Magnification Type        | 2010,0060 | CS |                     | ALWAYS            | CONFIG       | As in printer configuration file. |
| Min Density               | 2010,0120 | US |                     | ALWAYS            | CONFIG       | As in printer configuration file. |
| Max Density               | 2010,0130 | US |                     | ALWAYS            | CONFIG       | As in printer configuration file. |
| Trim                      | 2010,0140 | CS | NO, YES             | ALWAYS            | CONFIG, USER |                                   |
| Configuration Information | 2010,0150 | ST |                     | ALWAYS            | CONFIG       | As in printer configuration file. |

Table 72: Basic Film Box Relationship Module

| Attribute Name                   | Tag       | VR | Value | Presence of Value | Source | Comments                    |
|----------------------------------|-----------|----|-------|-------------------|--------|-----------------------------|
| Referenced Film Session Sequence | 2010,0500 | SQ |       | ALWAYS            | AUTO   |                             |
| >Referenced SOP Class UID        | 0008,1150 | UI |       | ALWAYS            | AUTO   | UID of Parent Film Session. |
| >Referenced SOP Instance UID     | 0008,1155 | UI |       | ALWAYS            | AUTO   |                             |

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

Table 73: Status Response

| Service Status | Error Code | Further Meaning               | Behavior  |
|----------------|------------|-------------------------------|---|
| Success        | 0000       | Film Box successfully created | Operation is completed successfully             |
| Warning        | 0116       | Attributes out of range       | Print job continues with a warning code logged. |

| Service Status | Error Code | Further Meaning         | Behavior   |
|----------------|------------|-------------------------|--|
| Failure        | 0106       | Invalid Attribute Value | Print job is terminated and error message is logged. |
|                | 0105       | No Such Attribute       | Print job is terminated and error message is logged. |

#### 4.2.2.3.1.7.2 Dataset Specific Conformance for Basic Film Box SOP Class N-ACTION-SCU

Detail regarding the Dataset Specific response behavior will be reported in this section.

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

Table 74: Status Response

| Service Status | Error Code | Further Meaning   | Behavior   |
|----------------|------------|---|--|
| Success        | 0000       | Film accepted for printing  | Operation is completed successfully.                 |
| Failure        | 0119       | Specified SOP Instance is not a member of the specified SOP class   | Print job is terminated and error message is logged. |
| Warning        | B603       | Film Box SOP Instance hierarchy does not contain Image Box SOP Instances  | Print job continues with a warning code logged.      |
|                | B604       | Image Size is larger than Image Box Size - The Image has been de-magnified  | Print job continues with a warning code logged.      |
|                | B609       | Image Size is larger than Image Box Size - The Image has been cropped to fit  | Print job continues with a warning code logged.      |
|                | B60A       | Image Size or combined Print Image Size is larger than Image Box Size - The Image or combined Print Image has been decimated to fit | Print job continues with a warning code logged.      |

#### 4.2.2.3.1.7.3 Dataset Specific Conformance for Basic Film Box SOP Class N-DELETE-SCU

Detail regarding the Dataset Specific response behavior will be reported in this section.

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

Table 75: Status Response

| Service Status     | Error Code | Further Meaning                   | Behavior  |
|--------------------|------------|-----------------------------------|---|
| Success            | 0000       | Film Session successfully created | Operation is completed successfully.                      |
| Other than Success | <xxxx>     | Any other status then success     | Print job remains in the queue manager with status failed |

#### 4.2.2.3.1.8 SOP Specific Conformance for Basic Grayscale Image Box SOP Class of the Basic Grayscale Print Management Meta SOP Class

This section and sub-section includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

The Printer process conforms to the Color Image Box Sop Class. The following DIMSE service element is supported: N-SET.

##### 4.2.2.3.1.8.1 Dataset Specific Conformance for Basic Grayscale Image Box SOP Class N-SET-SCU

Detail regarding the Dataset Specific response behavior will be reported in this section.

Table 76: Image Box Pixel Presentation Module

| Attribute Name                 | Tag       | VR    | Value       | Presence of Value | Source | Comments                          |
|--------------------------------|-----------|-------|-------------|-------------------|--------|-----------------------------------|
| Image Box Position             | 2020,0010 | US    | 1           | ALWAYS            | AUTO   |                                   |
| Basic Grayscale Image Sequence | 2020,0110 | SQ    |             | ALWAYS            | AUTO   |                                   |
| >Samples per Pixel             | 0028,0002 | US    | 1, 3        | ALWAYS            | AUTO   |                                   |
| >Photometric Interpretation    | 0028,0004 | CS    | MONOCHROME2 | ALWAYS            | AUTO   |                                   |
| >Rows                          | 0028,0010 | US    |             | ALWAYS            | AUTO   | As in printer configuration file. |
| >Columns                       | 0028,0011 | US    |             | ALWAYS            | AUTO   | As in printer configuration file. |
| >Pixel Aspect Ratio            | 0028,0034 | IS    |             | ALWAYS            | AUTO   | Must be present if not 1/1.       |
| >Bits Allocated                | 0028,0100 | US    | 8           | ALWAYS            | AUTO   |                                   |
| >Bits Stored                   | 0028,0101 | US    | 8           | ALWAYS            | AUTO   |                                   |
| >High Bit                      | 0028,0102 | US    | 7           | ALWAYS            | AUTO   |                                   |
| >Pixel Representation          | 0028,0103 | US    | 0           | ALWAYS            | AUTO   |                                   |
| >Pixel Data                    | 7FE0,0010 | OW/OB |             | ALWAYS            | AUTO   |                                   |

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

Table 77: DICOM Command Response Status Handling Behavior for Basic Grayscale Image Box N-SET

| Service Status | Error Code | Further Meaning                        | Behavior   |
|----------------|------------|--|--|
| Success        | 0000       | Image successfully stored in Image Box | Operation is completed successfully.                 |
| Warning        | 0107       | Attribute list error                   | Print job continues with a warning code logged.      |
| Failure        | 0106       | Invalid Attribute Value                | Print job is terminated and error message is logged. |
|                | 0105       | No Such Attribute                      | Print job is terminated and error message is logged. |

#### 4.2.2.3.1.9 SOP Specific Conformance for Printer SOP Class of the Basic Color Print Management Meta SOP Class

Not applicable, Printer SOP Class is not supported.

##### 4.2.2.3.1.9.1 Dataset Specific Conformance for Printer SOP Class N-EVENT-REPORT-SCU

Not applicable, Printer SOP Class is not supported.

#### 4.2.2.3.1.10 SOP Specific Conformance for Printer SOP Class of the Basic Grayscale Print Management Meta SOP Class

Not applicable, Printer SOP Class is not supported.

##### 4.2.2.3.1.10.1 Dataset Specific Conformance for Printer SOP Class N-EVENT-REPORT-SCU

Not applicable, Printer SOP Class is not supported.

#### 4.2.2.4. Association Acceptance Policy

Not applicable, Print Manager AE never accepts an association.

### 4.3. Network Interfaces

#### 4.3.1. Physical Network Interfaces

The System provides only DICOM V3.0 TCP/IP Network Communication Support as defined in PS 3.8 of the DICOM standard.

Big Bore inherits its TCP/IP stack from Windows 10 (i.e. the operating system platform).

Big Bore supports a single network interface Ethernet ISO 8802-3.

With standard supported physical medium include:

IEEE 802.3-1995, 10BASE-T

IEEE 802.3-1995, 100BASE-TX (Fast Ethernet)

IEEE 802.3, 1000BASE-X (Fiber Optic Gigabit Ethernet).

#### 4.3.2. Additional Protocols

Additional protocols such as used for network management are not applicable.

#### 4.3.3. IPv4 and IPv6 Support

Big Bore supports both IPv4 and IPv6 networks communication.

### 4.4. Configuration

Any implementation's DICOM conformance may be dependent upon configuration, which takes place at the time of installation. Issues concerning configuration are addressed in this section.

#### 4.4.1. AE Title/Presentation Address Mapping

The Big Bore system is configured by means of a configuration program. This program is accessible at start-up of the BigBore system. It is password protected and intended to be used by the administrator onsite or Philips Customer Support Engineers only.

Per configuration in LAN Config Tool, the system allows to accept associations from a range of IP addresses. With incoming association requests the system allows acceptance of a range of defined IP addresses which is configurable in the LAN Config application.

An important installation issue is the translation from AE title to presentation address. With incoming association requests the system allows acceptance of a range of defined IP addresses which is configurable in the LAN Config application. How this is performed is described in this section. The system is not IP or AE title sensitive.

##### 4.4.1.1. Local AE Titles

The local AE title mapping and configuration are specified as:

Table 78: AE Title Configuration

| Application Entity | Default AE Title | Default TCP/IP Port |
|--------------------|------------------|---------------------|
| DICOM Manager      | DATABASE         | 104 (configurable)  |
| Print Manager      | DATABASE         | 104 (configurable)  |

##### 4.4.1.2. Remote AE Title/Presentation Address Mapping

Remote AE Title, IP-Address, Port-number, supported DICOM Services and Transfer Syntaxes are freely configurable.

#### 4.4.2. Parameters

The specification of important operational parameters, their default value and range (if configurable) is specified here.

**Table 79: Configuration Local Parameters Table**

| Local Parameter                               | Configurable | Default Value                      |
|---|--------------|------------------------------------|
| <b>Basic Parameters</b>                       |              |                                    |
| Network - Computer Name                       | Yes          | <hostname> Set during installation |
| Network - IP                                  | Yes          | (0.0.0.0)                          |
| Network - NetMask                             | Yes          |                                    |
| Network - Gateway                             | Yes          | (0.0.0.0)                          |
| System Port                                   | Yes          | 104                                |
| System - AE Title                             | Yes          | DATABASE                           |
| MPPS Repair mode                              | Yes          | Checked                            |
| <b>Advanced Parameters</b>                    |              |                                    |
| IP Filter Mode                                | Yes          | Configured Device IPS              |
| Auto Import - Enable                          | Yes          | Unchecked                          |
| Auto Import - Input Folder Name               | Yes          | No Value                           |
| Auto Import - Failed Folder Name              | Yes          | No Value                           |
| Auto Import - Import Type                     | Yes          | DICOM                              |
| Auto Import – Polling Interval In Seconds     | Yes          | 60                                 |
| Verify Service Timeout in Seconds             | Yes          | 60                                 |
| Query Service Timeout in Seconds              | Yes          | 60                                 |
| Storage Service Timeout in Seconds            | Yes          | 300                                |
| Storage Commitment Service Timeout in Seconds | Yes          | 300                                |
| Retrieve Service Timeout in Seconds           | Yes          | 300                                |
| Print Service Timeout in Seconds              | Yes          | 300                                |

**Table 80: Configuration Remote Parameters Table**

| Remote Parameter                      | Configurable | Default Value       |
|---------------------------------------|--------------|---------------------|
| <b>Basic Parameters – Remote Node</b> |              |                     |
| Device Name                           | Yes          | -                   |
| Port                                  | Yes          | -                   |
| AE-Title                              | Yes          | -                   |
| Visible                               | Yes          | Visible             |
| Enable                                | Yes          | Checked             |
| Archived                              | Yes          | Unchecked           |
| CharSet                               | Yes          | Unchecked           |
| Large Archive                         | Yes          | Unchecked           |
| Max PDU                               | Yes          | 65536               |
| ARTimer                               | Yes          | 5                   |
| Transfer Syntax                       | Yes          | ELE and ILE enabled |
| <b>Advanced Parameters - Query</b>    |              |                     |
| Association Timeout In Seconds        | Yes          | 300                 |

| Remote Parameter                               | Configurable | Default Value |
|--|--------------|---------------|
| Lowest Support Level                           | Yes          | Image         |
| Query Response Size                            | Yes          | 100           |
| <b>Advanced Parameters - Store</b>             |              |               |
| Association Timeout In Seconds                 | Yes          | 120           |
| <b>Advanced Parameters - Retrieve</b>          |              |               |
| Association Timeout In Seconds                 | Yes          | 0             |
| <b>Advanced Parameter - Storage Commitment</b> |              |               |
| Association Timeout In Seconds                 | Yes          | 120           |
| <b>Advanced Parameter – Worklist Query</b>     |              |               |
| Association Timeout In Seconds                 | Yes          | 120           |
| Query Response Size                            | Yes          | 100           |
| <b>Advanced Parameter – MPPS</b>               |              |               |
| Association Timeout In Seconds                 | Yes          | 120           |

**Table 81: Configuration General Print Parameters Table**

| General Print Parameter        | Configurable | Default Value |
|--------------------------------|--------------|---------------|
| <b>Advanced Parameter</b>      |              |               |
| Association Timeout In Seconds | Yes          | 150           |

## 5. Media Interchange

### 5.1. Implementation model

The implementation model identifies the DICOM Application Entities for Media in specific implementation and relates the Application Entities to Real-World Activities.

#### 5.1.1. Application Data Flow Diagram

As part of the implementation model, an application data flow diagram is included. The next Figure shows the media interchange application data flow as a functional overview of the Media AE for DICOM CD and DVD.




Figure 18: Media Interchange Application Data Flow Diagram

The Media AE acts as a FSR when reading the directory of the medium. The Media AE acts as a FSC when writing the selected images in a patient folder onto the CD-R / DVD medium.

The Media AE acts as a FSR, for CD-R and DVD, when reading the directory of the medium. The Media AE acts as a FSC for CD-R and DVD, when writing the selected images in a patient folder onto the medium.

#### 5.1.2. Functional Definitions of AE's

This section contains a functional definition for each local Application Entity. It's described in general terms the functions to be performed by the AE, and the DICOM services used to accomplish these functions.

The BigBore system can create and Read CD/DVD.

The Media AE in an BigBore system supports the following functions for CD and DVD as FSR:

- Read the DICOMDIR File from the medium (representing the directory of the DICOM File(s) as recorded on the medium). This information may be displayed as an ordered list of icon images and, if present, with pertinent identifying information (patient name, etc.).
- Read the selected image from the medium and display it on the monitor of the View Station. This information is displayed as an ordered list of frames of the selected image or as a dynamic review of the selected image.

Remark:

1. CD's and DVD's that are created according to the application profiles STD-GEN-CD, STD-CTMR-CD and STD-CTMR-DVD with Explicit VR Little Endian Uncompressed transfer syntax are supported as FSR by the BigBore system.
2. CD's and DVD's that are created according to the application profiles STD-CTMR-CD and STD-CTMR-DVD with JPEG Lossless Process 14 (selection value 1) transfer syntax are not supported as FSR by the BigBore system.

And for CD and DVD as FSC:

- Initialize the medium.
- Write a DICOM File-set onto the medium.
- Create a DICOMDIR File.
- Extend the DICOM File-set and update the DICOMDIR File accordingly. (DICOM Media Storage Service Class).

### 5.1.3. Sequencing of Real World Activities

A Real World Activity of the Media AE is: The user selects a set of objects to write to the CD/DVD. Then the CD /DVD will be created with the selected objects. Once the CD/DVD has been created, the user can read this CD/DVD on the BigBore system or for transport to another device for reading.




Figure 19: (Real World) Activity - Media

Another Real World Activity of the Media AE is: A CD/DVD from another system or previously created CD/DVD can be read by the BigBore system. The BigBore cannot append (FSU) to this created CD/DVD.

**Remark:** DVD's that were created following the definition of the STD-GEN-DVD-JPEG application profile are not supported for reading since BigBore does not support JPEGBaseline, JPEGFOP14 and JPEG2000 compression.

After data is written to DVD, the DVD is finalized; the finalized DVD can now be read on mostly every DVD reader.

### 5.1.4. File Meta Information for Implementation Class and Version

This section contains the values of the file Meta Information that pertain to the Application Entity (see PS 3.10). These are:

- Source Application Entity Title,
- Private Information Creator UID,
- Private Information.

The Application Entity title is registered into the DICOM File Meta Information header and is supported by the CD/DVD-Writer (CD/DVD write option) acting as a FSC.

Table 82: File Meta Information for the Media AE

|                             |                         |
|-----------------------------|-------------------------|
| Implementation Class UID    | 1.3.46.670589.33.103.11 |
| Implementation Version Name | EBW4.8 for CT           |

## 5.2. AE Specifications

This section in the DICOM Conformance Statement specifies a set of Media Application Entities.

### 5.2.1. Media - Specification

The supported Application Profiles, their Roles and the Service Class options, all defined in DICOM terminology, are listed in the following table.

For reading and writing the media AE provides standard conformance to:

- DICOM media Storage Service and File Format ([DICOM] PS 3.10);
- Media Storage Application Profiles ([DICOM] PS 3.11); and
- Media Formats and Physical Media for Media Interchange ([DICOM] PS 3.12) for Reading (FSR) and Writing (FSC).

BigBore system does not support multi-session writes to CD/DVD's.

Supported media:

- CD: CD-R and CD-RW with the profile STD-GEN-CD.
- DVD: DVD-R, DVD+R, DVD-RW and DVD+RW with the profile: STD-GEN-DVD-JPEG as FSR.

The supported Application Profiles, their roles and the Service Class (SC) options, all defined in DICOM terminology, are listed in the next Table.

**Table 83: AE Media AE related Application Profiles, RWA activities and roles**

| Supported Application Profile             | Identifier       | Real-World Activities | Roles |
|---|------------------|-----------------------|-------|
| CT/MR Studies on CD-R                     | STD-CTMR-CD      | Create File-set       | FSC   |
|   |                  | Read File-set         | FSR   |
|   |                  | Display Directory     | DD    |
| General Purpose CD-R Interchange          | STD-GEN-CD       | Create File-set       | FSC   |
|   |                  | Read File-set         | FSR   |
|   |                  | Display Directory     | DD    |
| CT/MR Studies on DVD Media                | STD-CTMR-DVD     | Create File-set       | FSC   |
|   |                  | Read File-set         | FSR   |
|   |                  | Display Directory     | DD    |
| General Purpose DVD Interchange with JPEG | STD-GEN-DVD-JPEG | Create File-set       | FSC   |
|   |                  | Read File-set         | FSR   |
|   |                  | Display Directory     | DD    |

### 5.2.1.1. DICOM File Meta Information

Table below denotes the DICOM file meta attributes included in the DICOMDIR that are created by the Big Bore system.

**Table 84: DICOM Part 10 File Meta Information**

| Attribute Name                     | Tag         | Notes   |
|------------------------------------|-------------|---|
| File Preamble                      | N/A         | All bytes are set to 00H                                |
| DICOM Prefix                       | N/A         | Set to DICOM Prefix "DICM"                              |
| File Meta Information Group Length | (0002,0000) |   |
| File Meta Information Version      | (0002,0001) | Set to 0001H  |
| Media Storage SOP Class UID        | (0002,0002) | 1.2.840.10008.1.3.10                                    |
| Media Storage SOP Instance UID     | (0002,0003) | Big Bore 4.8 system generated UID                       |
| Transfer Syntax UID                | (0002,0010) | Set to Explicit VR Little Endian<br>1.2.840.10008.1.2.1 |
| Implementation Class UID           | (0002,0012) | Set to "1.3.46.670589.33.103.11"                        |
| Implementation Version Name        | (0002,0013) | Set to "EBW4.8 for CT"                                  |

Table 85: DICOMDIR Attributes

| Attribute Name  | Tag       | Notes |
|---|-----------|-------|
| File Set ID   | 0004,1130 |       |
| Offset of the First Directory Record of the Root Directory Entity | 0004,1200 |       |
| Offset of the Last Directory Record of the Root Directory Entity  | 0004,1202 |       |
| File Set Consistency Flag   | 0004,1212 |       |
| Directory Record Sequence   | 0004,1220 |       |
| >Offset of the Next Directory Record                              | 0004,1400 |       |
| >Record In-use Flag   | 0004,1410 |       |
| >Offset of Referenced Lower-Level Directory Entity                | 0004,1420 |       |
| >Directory Record Type  | 0004,1430 |       |
| >Referenced File ID   | 0004,1500 |       |
| >Referenced SOP Class UID in File                                 | 0004,1510 |       |
| >Referenced Transfer Syntax UID in File                           | 0004,1512 |       |
| <b>Patient Level</b>  |           |       |
| Offset Of The Next Dir Record                                     | 0004,1400 |       |
| Record In Use Flag  | 0004,1410 |       |
| Offset of Referenced Lower-Level Directory Entity                 | 0004,1420 |       |
| Directory Record Type   | 0004,1430 |       |
| Specific Character Set  | 0008,0005 |       |
| Patient's Name  | 0010,0010 |       |
| Patient ID  | 0010,0020 |       |
| <b>Study Level</b>  |           |       |
| Offset Of The Next Dir Record                                     | 0004,1400 |       |
| Record In Use Flag  | 0004,1410 |       |
| Offset Of Ref Lower Level Dir Entity                              | 0004,1420 |       |
| Directory Record Type   | 0004,1430 |       |
| Specific Character Set  | 0008,0005 |       |
| Study Date  | 0008,0020 |       |
| Study Time  | 0008,0030 |       |
| Accession Number  | 0008,0050 |       |
| Study Description   | 0008,1030 |       |
| Study Instance UID  | 0020,000D |       |
| Study ID  | 0020,0010 |       |
| Modalities in Study   | 0008,0061 |       |
| Institution Name  | 0008,0080 |       |
| Referring Physician's Name  | 0008,0090 |       |
| Operators Name  | 0008,1070 |       |
| Patient's Birth Date  | 0010,0030 |       |
| Patient's Sex   | 0010,0040 |       |
| Number Of Study Related Series                                    | 0020,1206 |       |
| Number Of Study Related Instances                                 | 0020,1208 |       |
| Requested Procedure ID  | 0040,1001 |       |
| <b>Series Level</b>   |           |       |

| Attribute Name                                    | Tag       | Notes |
|---|-----------|-------|
| Offset Of The Next Directory Record               | 0004,1400 |       |
| Record In-Use Flag                                | 0004,1410 |       |
| Offset of Referenced Lower-Level Directory Entity | 0004,1420 |       |
| Directory Record Type                             | 0004,1430 |       |
| Specific Character Set                            | 0008,0005 |       |
| Series Date                                       | 0008,0021 |       |
| Content Date                                      | 0008,0023 |       |
| Series Time                                       | 0008,0031 |       |
| Content Time                                      | 0008,0033 |       |
| Modality  | 0008,0060 |       |
| Manufacturer                                      | 0008,0070 |       |
| Series Instance UID                               | 0020,000E |       |
| Series Number                                     | 0020,0011 |       |
| Series Description                                | 0008,103E |       |
| Body Part Examined                                | 0018,0015 |       |
| Protocol Name                                     | 0018,1030 |       |
| Number of Series Related Instances                | 0020,1209 |       |
| <b>Image Level</b>                                |           |       |
| Offset Of The Next Directory Record               | 0004,1400 |       |
| Record In-Use Flag                                | 0004,1410 |       |
| Offset Of Referenced Lower-Level Directory Entity | 0004,1420 |       |
| Directory Record Type                             | 0004,1430 |       |
| Referenced File ID                                | 0004,1500 |       |
| Referenced SOP Class UID In File                  | 0004,1510 |       |
| Referenced SOP Instance UID In File               | 0004,1511 |       |
| Referenced Transfer Syntax UID In File            | 0004,1512 |       |
| Specific Character Set                            | 0008,0005 |       |
| Instance Number                                   | 0020,0013 |       |
| Image Type  | 0008,0008 |       |
| Instance Creation Date                            | 0008,0012 |       |
| Instance Creation Time                            | 0008,0013 |       |
| SOP Class UID                                     | 0008,0016 |       |
| SOP Instance UID                                  | 0008,0018 |       |
| Acquisition Date                                  | 0008,0022 |       |
| Acquisition Time                                  | 0008,0032 |       |
| Content Date                                      | 0008,0023 |       |
| Content Time                                      | 0008,0033 |       |
| Acquisition Number                                | 0020,0012 |       |
| Slice Thickness                                   | 0018,0050 |       |
| KVP   | 0018,0060 |       |
| Reconstruction Diameter                           | 0018,1100 |       |
| Gantry/Detector Tilt                              | 0018,1120 |       |
| Angular Position                                  | 0018,1141 |       |
| Image Position (Patient)                          | 0020,0032 |       |
| Image Orientation (Patient)                       | 0020,0037 |       |

| Attribute Name             | Tag       | Notes |
|----------------------------|-----------|-------|
| Frame of Reference UID     | 0020,0052 |       |
| Slice Location             | 0020,1041 |       |
| Image Comments             | 0020,4000 |       |
| Samples per Pixel          | 0028,0002 |       |
| Photometric Interpretation | 0028,0004 |       |
| Rows                       | 0028,0010 |       |
| Columns                    | 0028,0011 |       |
| Pixel Spacing              | 0028,0030 |       |
| Bits Allocated             | 0028,0100 |       |

### 5.2.1.2. Real-World Activities

The AE specification contains a description of the Real-World Activities, which invoke the particular AE.

#### 5.2.1.2.1 RWA - Read File-set

When an image transfer from CD or DVD is initiated then the Media AE acts as an FSR using the interchange option to import SOP Instances from the CD or DVD medium.

##### 5.2.1.2.1.1 Media Storage Application Profile

The media AE supports the RWA Read File-set for the STD-CT/MR studies on CD, the STD-GEN-CD, the General Purpose DVD Interchange with JPEG and CT/MR Studies on DVD Media application profiles.

##### 5.2.1.2.1.1.1 Options

Not applicable. No options are implemented by the system.

#### 5.2.1.2.2 RWA - Create File-set

This Media Application Entity has a File-set Creator functionality which is describe here.

When an image transfer to CD/DVD is initiated then the Media AE acts as an FSC using the interchange option to export SOP Instances from the local database to a CD/DVD medium.

##### 5.2.1.2.2.1 Media Storage Application Profile

The application Profile that is used by this Media Application Entity is specified in this section.

As depicted in the table in section 5.2.1, the Media AE supports the RWA Write Images for the STD-CTMR-CD, STD-GEN-CD, STD-GEN-DVD-JPEG and the STD-CTMR-DVD Application Profile.

The DICOMDIR file will be extended when new images are written. In case some attributes are not present in an image but are specified as mandatory in the DICOMDIR definition in DICOM Media, a generated value will be filled in.

#### Implementation remarks and restrictions

When writing the DICOMDIR records, key values are generated when no value of the corresponding attribute is supplied, according to the following tables.

**Table 86: Generated keys**

| Key               | Tag         | Generated Value                       |
|-------------------|-------------|---------------------------------------|
| <b>Study Keys</b> |             |                                       |
| Study Date        | (0008,0020) | Date on which this Study was created. |

| Key                | Tag         | Generated Value                       |
|--------------------|-------------|---------------------------------------|
| Study Time         | (0008,0030) | Time on which this Study was created. |
| <b>Series Keys</b> |             |                                       |
| Series Number      | (0020,0011) | 1                                     |
| <b>Image Keys</b>  |             |                                       |
| Instance Number    | (0020,0013) | 1 (if empty)                          |

The data selected to write to the media must fit on the currently inserted media. If it does not fit, an error is generated and it is up to the operator to re-select a smaller amount of data to be written to the media. The system will not request additional media or write across multiple media.

#### 5.2.1.2.2.1.1 Options

Not applicable. No options are implemented by the system

#### 5.2.1.2.3 RWA - Display Directory

When a database open action is initiated on the CD/DVD then the Media AE acts as an FSR using the interchange option to read the DICOMDIR of the CD/DVD media.

This will results in an overview of the patients, studies, series and images on the BigBore screen.

#### 5.2.1.2.3.1 Media Storage Application Profile

The media AE supports the RWA Display Directory for the STD-CT/MR studies on CD, the STD-GEN-CD, the General Purpose DVD Interchange with JPEG and CT/MR Studies on DVD Media application profiles.

#### 5.2.1.2.3.1.1 Options

The mandatory DICOMDIR keys are required for the correct display of directory information. The display is structured according the DICOM Composite Information Model: Patient, Study, Series, and Image.

The DICOM standard specifies certain attributes of the DICOMDIR as mandatory. However, these attributes may not be mandatory for the related SOP class IOD. For those attributes the default values apply.

### 5.3. Augmented and Private Application Profiles

Not supported. Big Bore system has not implemented any Augmented and Private Application Profiles.

### 5.4. Media Configuration

Any configuration issues may be found in the Networking Section 4.4 Configuration.

## 6. Support of Character Sets

Any support for character sets in Network and Media services is described here.

**Table 87: Supported DICOM Character Sets**

| Character Set Description | Defined Term | ESC Sequence | ISO Registration Number | Code Element | Character Set                 |
|---------------------------|--------------|--------------|-------------------------|--------------|-------------------------------|
| Latin alphabet No. 1      | ISO_IR 100   | -            | ISO-IR 100              | G1           | Supplementary set of ISO 8859 |
|                           |              | -            | ISO-IR 6                | G0           | ISO 646                       |
| Chinese                   | GB18030      | -            | GB18030                 | -            | -                             |
| Default repertoire        | -            | -            | ISO-IR 6                | G0           | ISO 646                       |

## 7. Security

### 7.1. Security Profiles

Big Bore allows the use of either a conventional (non-secure) DICOM communication or a secure DICOM communication based on the Transport Layer Security (TLS) protocol. If configured Big Bore supports the following security measures:

- Secure authentication of a node
- Integrity and confidentiality of transmitted data
- Generation of audit trail records access control and user authentication

#### 7.1.1. Security use Profiles

Not Applicable. Big Bore system has not implemented Security Use Profiles.

#### 7.1.2. Security Transport Connection Profiles

The TLS Component is a “mode of operation” of Big Bore and will be used for nodes that can authenticate each other before they communicate over sockets. TLS 1.2 can only be used using TCP. Node authentication and encryption are only possible when the node has:

- a “private and public key”;
- a self-signed certificate or certificate signed by a Certificate Authority; and
- a list of certificates with which the system wants to communicate.

Furthermore, the TLS component may communicate using the following Cipher Suites:

TLS\_DHE\_RSA\_WITH\_AES\_128\_GCM\_SHA256  
TLS\_DHE\_RSA\_WITH\_AES\_256\_GCM\_SHA384  
TLS\_ECDHE\_RSA\_WITH\_AES\_128\_GCM\_SHA256  
TLS\_ECDHE\_RSA\_WITH\_AES\_256\_GCM\_SHA384  
TLS\_ECDHE\_RSA\_WITH\_AES\_256\_CBC\_SHA384  
TLS\_ECDHE\_RSA\_WITH\_AES\_256\_CBC\_SHA  
TLS\_RSA\_WITH\_AES\_256\_CBC\_SHA256  
TLS\_RSA\_WITH\_AES\_256\_GCM\_SHA384  
TLS\_ECDHE\_RSA\_WITH\_AES\_128\_CBC\_SHA256  
TLS\_ECDHE\_RSA\_WITH\_AES\_128\_CBC\_SHA  
TLS\_RSA\_WITH\_AES\_128\_CBC\_SHA256  
TLS\_RSA\_WITH\_AES\_128\_GCM\_SHA256  
TLS\_RSA\_WITH\_3DES\_EDE\_CBC\_SHA  
TLS\_RSA\_WITH\_AES\_256\_CBC\_SHA  
TLS\_RSA\_WITH\_AES\_128\_CBC\_SHA  
TLS\_RSA\_WITH\_AES\_128\_GCM\_SHA256  
TLS\_DHE\_RSA\_WITH\_AES\_256\_CBC\_SHA  
TLS\_DHE\_RSA\_WITH\_AES\_128\_CBC\_SHA  
TLS\_RSA\_WITH\_RC4\_128\_SHA  
TLS\_RSA\_WITH\_RC4\_128\_MD5

In case no encryption is used the data is signed and hashed: integrity is present and confidentiality is not present.

#### Certificates

If two systems communicate with each other, one system will be listening on a port (server node) while the other system sets up a connection (client node). The certificate this server node will send to the other client node is the server certificate. The client node initiates the communication and the certificate that the client

node is sending to the server is the client certificate. (Server Client Authentication) The following TLS Certification checks will be done (TLS Handshake). The machine (either server or client) that will send its certificate will choose the certificate according to Common Name (CN) value in the Subject-field. This name is case-sensitive. All present certificates should have unique CN names. **The server verifies:**

- that the client certificate is a valid X.509 certificate;
- that the client certificate is either signed by a CA or is self-signed;
- that the client certificate is in the list of trusted certificates;
- that the client certificate is valid (present time is between "Valid from" and "Valid to" fields of the X.509 certificate);
- that the client certificate has the correct purpose (at least the Client Authenticate purpose).

**The client verifies:**

- that the server certificate is a valid X.509 certificate;
- that the server certificate either is signed by a CA or is self-signed;
- that the server certificate is in the list of trusted certificates;
- that the server certificate is valid (present time is between "Valid from" and "Valid to" fields of the X.509 certificate);
- that the server certificate has the correct purpose (at least Server Authenticate purpose).

In the TLS component no verification is done on:

- revocation of certificates;
- limiting the connection to a limited set of IP-addresses.

Additional information: The value in the Subject-field is determined in the certificate request. The CA will sign the request in case it accepts the values that are present in the request. The CN value can be: IP-number, hostname or hostname. Domain. The value in the CN-field must be equal to the value that is used in making a connection to the server. In case the name is specified as hostname. Domain that same value should be specified during connect. In the ideal situation the name-IP-number translation will be dealt with by the DNS in the hospital. This check is case-insensitive.

Figure below presents the message flow of TLS handshake supported by Big Bore.




Figure 20: Message flow of TLS handshake

### 7.1.3. Digital Signature Profiles

Not Applicable. Big Bore system has not implemented Digital Signature Profiles.

### 7.1.4. Media Storage Security Profiles

Not applicable. Big Bore system has not implemented Media Storage Security Profiles.

### 7.1.5. Attribute Confidentiality Profiles

No instances of the Encrypted Attributes Data Set are created. No Transfer Syntaxes are supported for encoding/decoding of Encrypted Attributes Data Sets.

The table below lists the attributes that are replaced during the identification process.

The terms used to describe the replacement value can be read as below:

- Empty: The attribute will have a value of zero length, is cleared by BigBore system.
- Copied: Attribute has same value as original.

Table 88: Basic Application Level Confidentiality Profile Attributes

| Attribute Name         | Tag       | VR | Replacement Value     |
|------------------------|-----------|----|-----------------------|
| Specific Character Set | 0008,0005 | CS | Copied from original  |
| Instance Creator UID   | 0008,0014 | UI | [empty]               |
| SOP Instance UID       | 0008,0018 | UI | New created UID       |
| Study Date             | 0008,0020 | DT | Replaced by new value |
| Series Date            | 0008,0021 | DT | Replaced by new value |
| Acquisition Date       | 0008,0022 | DT | Replaced by new value |
| Content Date           | 0008,0023 | DT | Replaced by new value |

| Attribute Name                          | Tag       | VR | Replacement Value   |
|---|-----------|----|---|
| Study Time                              | 0008,0030 | TM | Copied from original  |
| Series Time                             | 0008,0031 | TM | Copied from original  |
| Acquisition Time                        | 0008,0032 | TM | Copied from original  |
| Content Time                            | 0008,0033 | TM | Copied from original  |
| Accession Number                        | 0008,0050 | SH | [empty]   |
| Institution Name                        | 0008,0080 | LO | [empty]   |
| Institution Address                     | 0008,0081 | ST | [empty]   |
| Referring Physician's Name              | 0008,0090 | PN | [empty]   |
| Referring Physician's Address           | 0008,0092 | ST | [empty]   |
| Referring Physician's Telephone Numbers | 0008,0094 | SH | [empty]   |
| Station Name                            | 0008,1010 | SH | [empty]   |
| Study Description                       | 0008,1030 | LO | [empty]   |
| Series Description                      | 0008,103E | LO | [empty]   |
| Institutional Department Name           | 0008,1040 | LO | [empty]   |
| Physician(s) of Record                  | 0008,1048 | PN | [empty]   |
| Performing Physicians' Name             | 0008,1050 | PN | [empty]   |
| Name of Physician(s) Reading Study      | 0008,1060 | PN | [empty]   |
| Operators' Name                         | 0008,1070 | PN | [empty]   |
| Admitting Diagnoses Description         | 0008,1080 | LO | [empty]   |
| Manufacturers Model Name                | 0008,1090 | LO | [empty]   |
| Referenced Study Sequence               | 0008,1110 | SQ | [empty]   |
| Referenced Patient Sequence             | 0008,1120 | SQ | [empty]   |
| Referenced Image Sequence               | 0008,1140 | SQ | Replaced by new value   |
| Derivation Description                  | 0008,2111 | ST | [empty]   |
| Patient's Name                          | 0010,0010 | PN | Initials - Last and First Name component can be changed by user |
| Patient ID                              | 0010,0020 | LO | [empty] - Can be changed by user                                |
| Patient's Birth Date                    | 0010,0030 | DA | Replaced by new value   |
| Patient's Birth Time                    | 0010,0032 | TM | Replaced by new value   |
| Patient's Sex                           | 0010,0040 | CS | Copied from original  |
| Other Patient ID's                      | 0010,1000 | LO | [empty]   |
| Other Patient Names                     | 0010,1001 | PN | [empty]   |
| Patient's Age                           | 0010,1010 | AS | [empty]   |
| Patient Size                            | 0010,1020 | DS | Copied from original  |
| Patient Weight                          | 0010,1030 | DS | Copied from original  |
| Medical Record Locator                  | 0010,1090 | LO | [empty]   |
| Medical Alerts                          | 0010,2000 | LO | [empty]   |
| Ethnic Group                            | 0010,2160 | SH | [empty]   |
| Occupation                              | 0010,2180 | SH | [empty]   |
| Additional Patient's History            | 0010,21B0 | LT | [empty]   |
| Pregnancy Status                        | 0010,21C0 | US | [empty]   |
| Patient Comment                         | 0010,4000 | LT | [empty]   |
| Device Serial Number                    | 0018,1000 | LO | [empty]   |
| Protocol Name                           | 0018,1030 | LO | [empty]   |
| Study Instance UID                      | 0020,000D | UI | New created UID   |

| Attribute Name                          | Tag       | VR | Replacement Value    |
|---|-----------|----|----------------------|
| Series Instance UID                     | 0020,000E | UI | New created UID      |
| Study ID                                | 0020,0010 | SH | New created UID      |
| Frame of Reference UID                  | 0020,0052 | UI | New created UID      |
| Synchronization Frame of Reference UID  | 0020,0200 | UI | New created UID      |
| Image Comments                          | 0020,4000 | LT | [empty]              |
| Requesting Physician                    | 0032,1032 | PN | [empty]              |
| Requesting Service                      | 0032,1033 | LO | [empty]              |
| Requested Procedure Description         | 0032,1060 | LO | [empty]              |
| Requested Procedure Code Sequence       | 0032,1064 | SQ | [empty]              |
| Admission ID                            | 0038,0010 | LO | [empty]              |
| Special Needs                           | 0038,0050 | LO | [empty]              |
| Current Patient Location                | 0038,0300 | LO | [empty]              |
| Patient State                           | 0038,0500 | LO | [empty]              |
| Scheduled Procedure Step Sequence       | 0040,0100 | SQ | [empty]              |
| Performed Procedure Step Start Date     | 0040,0244 | DA | [empty]              |
| Performed Procedure Step Start Time     | 0040,0245 | TM | [empty]              |
| Series Description                      | 0008,103E | LO | [empty]              |
| Institutional Department Name           | 0008,1040 | LO | [empty]              |
| Performed Procedure Step ID             | 0040,0253 | SH | Copied from original |
| Performed Procedure Step Description    | 0040,0254 | LO | Copied from original |
| Request Attributes Sequence             | 0040,0275 | SQ | [empty]              |
| Requested Procedure ID                  | 0040,1001 | SH | [empty]              |
| Names of Intended recipients of Results | 0040,1010 | PN | [empty]              |
| Requested Procedure Comments            | 0040,1400 | LT | [empty]              |
| Imaging Service Request Comments        | 0040,2400 | LT | [empty]              |
| UID                                     | 0040,A124 | UI | New created UID      |
| Content Sequence                        | 0040,A730 | SQ | [empty]              |
| Storage Media File Set UID              | 0088,0140 | UI | [empty]              |
| Referenced Frame of Reference UID       | 3006,0024 | UI | [empty]              |
| Related Frame of Reference UID          | 3006,00C2 | UI | [empty]              |

#### 7.1.6. Network Address Management Profiles

Not applicable. Big Bore system has not implemented Network Address Management Profiles.

#### 7.1.7. Time Synchronization Profiles

Big Bore System conforms to the Basic Time Synchronization Profile as NTP Client. Big Bore System does support secure transactions.

#### 7.1.8. Application Configuration Management Profiles

Not applicable. Big Bore system has not implemented Application Configuration Management Profiles.

#### 7.1.9. Audit Trail Profiles

Big Bore creates audit messages according to the IHE ATNA Integration Profile. These messages may contain information that identifies the patient. The following messages will be created and sent to a central Audit Record Repository:

Table 89: List of supported events

| Audit Event Trigger      | Description  | Message<br>DICOM PS 3.15 A.5.3     |
|--------------------------|--|------------------------------------|
| Actor-start-stop         | When Big Bore system starts or shuts down  | Application Activity               |
| Begin-storing-instances  | When an examination is being transferred from the Big Bore system to a remote network node   | Begin Transferring DICOM Instances |
| Instances-Stored         | When an examination is transferred from the Big Bore system to a remote network node   | DICOM Instances Transferred        |
| Study-used               | Study is created, modified, accessed, or deleted   | DICOM Instances Accessed           |
| User Authentication      | When the user logs in or logs out  | User Authentication.               |
| Security Alert           | When an authentication of a secure node during TLS negotiation fails, e.g. due to an invalid certificate and configuration and other changes | Security Alert                     |
| PHI-export               | When printing job is started or export on media e.g CD, DVD  | Export                             |
| PHI-import               | Any import of PHI from removable media.  | Import                             |
| Image Availability Query | A query has been initiated from patient directory to a remote node   | Query                              |
| Audit Log Used           | The audit trail repository has been accessed or modified by something other than the arrival of audit trail messages.                        | Audit Log Used                     |
| Instance deleted         | If SOP Instances are deleted from a specific study. One event covers all instances deleted for the particular study                          | DICOM Study Deleted                |

## 7.2. Association Level Security

Not applicable. Big Bore system has not implemented Association Level Security Profiles.

## 7.3. Application Level Security

Big Bore does not support any specific application level security measures.

- The Application which gives access to Patient records and DICOM communication requires Login with Username and Password.
- The system is used within a secured environment. It is assumed that a secured environment includes at a minimum.
  - The OS is solidified by white-listing applications and files. Not white-listed executable files, libraries, drivers, Java apps, ActiveX controls, scripts, and other code are blocked.
  - Firewall or router protections to ensure that only approved external hosts have network access to Big Bore System.
  - Firewall or router protections to ensure that Big Bore System only has network access to approved external hosts and services.
  - Any communication with external hosts outside the locally secured environment can be configured to use secure network channels.

Other network security procedures such as automated intrusion detection may be appropriate in some environments.

Additional security features may be established by the local security policy and are beyond the scope of this conformance statement.

## 8. Annexes of application "Big Bore v4.8"

### 8.1. IOD Contents

#### 8.1.1. Created SOP Instance

##### General Rules

This section specifies each IOD created by this application.

The Big Bore reflects the fact that the IOD created by the workstation are always based on some source images after the viewing/processing applied and the modified images are saved. Most of the attributes or even the whole modules (Patient, General Study, etc.) are just copied from the source images.

DICOM Overlays are only created for saved Secondary Capture SOP Class images. When images with ROI, Annotations, etc. (group 50xx) are saved as DICOM Secondary Captures, the 50xx groups attributes are converted into DICOM Overlays attributes (group 60xx). In case SCP does not support group 60xx attributes - the Big Bore has a configurable option (in LAN Config) to burn the overlays into the pixel data thus allowing any PACS to display them.

The attribute "Burned In Annotation", for saving displays with multiple images has the value "YES". The attribute "Burned In Annotation" has the value "NO" for derived objects, if saved with "hide titles", (only for secondary capture SOP Class objects).

For Value Representation (VR) equal to Patient's Name (PN), the leading spaces into the Patient's Name will be treated as insignificant for matching purposes.

From the Patient's Name only the first 32 characters are displayed into the Quick View Viewer.

A Time attribute contains a string of characters of the format "hhmmss.frac". The Fractional part is always 3 decimal places.

##### DICOM Origin

The DICOM standard does not define the origin of the patient-based coordinate system used to define voxel coordinates. Following is the definition used by Philips CT scanners:

- X: Center of the CT bore from left to right
- Y: (Center of the CT bore from top to bottom) + Table Height (tag 0018,1130) - 255 mm
- Z: User configurable via the "Zero couch" gantry button.

Note that the image position may be used to determine the coordinates of the center of the CT bore only if the reconstruction offset (tag (01F1,100C), private creator code "ELSCINT1") is "0/0".

Coordinate system polarities are defined as per the DICOM standard.

##### Derived CT Image Attributes

Image Plane module attributes:

- All derived CT images, except curved (panoramic) slab, contain the Image Position (0028,0032) and Image Orientation (0028,0037) attributes.
- All derived CT (including curve slab) images contain the Pixel Spacing (0028,0030) and Slice Thickness (0018,0050) attributes.
- Non-Square pixels are not supported by "Big Bore".

##### Export Converters

A number of configurable export convertors allow the system to modify certain IOD when sent to specific SCP.

### Color-to-Monochrome Converter

The Color (24-bits) SC IOD is converted to a monochrome 8-bit or 12 bit SC IOD, configurable by FSE. A new UID is generated for the converted image.

### Burn Overlays Converter

The overlays on the images are burned into the Pixel data for the image and the DICOM overlay groups (50xx and 60xx) are removed from the IOD. This convertor works only on Secondary Capture IOD's.

### Convert Philips Tags to CCA Cardiac Tags

Private Philips DICOM attributes are converted to allow the IOD to be processed by non-Philips systems.

### Copied modules to the derived IOD's

The following table lists the modules that are always copied from the source images when the created SOP Class IOD is the same as the source SOP Class IOD.

**Table 90: Modules copied to the derived IOD's table**

| Information Entity | Module Name                   |
|--------------------|-------------------------------|
| Patient            | Patient Module                |
|                    | Clinical Trial Subject Module |
| Study              | General Study Module          |
|                    | Patient Study Module          |
|                    | Clinical Trial Study Modules  |
| Series             | General Series Modules        |
|                    | Clinical Trial Series Module  |
| Frame of Reference | Frame of Reference Module     |
| Equipment          | General Equipment Module      |

This section specifies each IOD created (including private IOD's). It should specify the attribute name, tag, VR, and value. The value should specify the range and source (e.g. user input, Modality Worklist, automatically generated, etc.). For content items in templates, the range and source of the concept name and concept values should be specified. Whether the value is always present or not shall be specified.

Abbreviations used in the IOD tables for the column "Presence of Module" are:

ALWAYS            The module is always present  
 CONDITIONAL    The module is used under specified condition

Abbreviations used in the Module table for the column "Presence of Value" are:

ALWAYS            The attribute is always present with a value  
 EMPTY            The attribute is always present without any value (attribute sent zero length)  
 VNAP              The attribute is always present and its Value is Not Always Present  
                       (attribute sent zero length if no value is present)  
 ANAP              The attribute is present under specified condition – if present then it will always have a Value

The abbreviations used in the Module table for the column "Source" are:

AUTO              The attribute value is generated automatically  
 CONFIG           The attribute value source is a configurable parameter  
 COPY              The attribute value source is another SOP instance  
 FIXED             The attribute value is hard-coded in the application

|          |   |
|----------|---|
| IMPLICIT | The attribute value source is a user-implicit setting                             |
| MPPS     | The attribute value is the same as that use for Modality Performed Procedure Step |
| MWL      | The attribute value source is a Modality Worklist                                 |
| USER     | The attribute value source is explicit user input                                 |

Items in the Value and Comment columns in the following tables are filled in where appropriate to further clarify the use or meaning of each attribute beyond the definition provided by the DICOM Standard. All others are left blank for ease of use. See PS3.3 of the DICOM Standard for the complete attribute definitions.

#### 8.1.1.1. List of created SOP Classes

Table 91: List of created SOP Classes

| SOP Class Name                            | SOP Class UID                 |
|---|-------------------------------|
| CT Image Storage SOP Class                | 1.2.840.10008.5.1.4.1.1.2     |
| Secondary Capture Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7     |
| General ECG Waveform Storage SOP Class    | 1.2.840.10008.5.1.4.1.1.9.1.2 |
| X-Ray Radiation Dose SR                   | 1.2.840.10008.5.1.4.1.1.88.67 |

#### 8.1.1.1.1 CT Image Storage SOP Class

Table 92: SOP Class Modules

| Information Entity | Module                     | Presence    |
|--------------------|----------------------------|-------------|
| Patient            | Patient Module             | ALWAYS      |
| Study              | General Study Module       | ALWAYS      |
|                    | Patient Study Module       | CONDITIONAL |
| Series             | General Series Module      | ALWAYS      |
| Frame of Reference | Frame of Reference Module  | ALWAYS      |
| Equipment          | General Equipment Module   | ALWAYS      |
| Acquisition        | General Acquisition Module | ALWAYS      |
| Image              | General Image Module       | ALWAYS      |
|                    | General Reference Module   | ALWAYS      |
|                    | Image Plane Module         | ALWAYS      |
|                    | Image Pixel Module         | ALWAYS      |
|                    | CT Image Module            | ALWAYS      |
|                    | VOI LUT Module             | CONDITIONAL |
|                    | SOP Common Module          | ALWAYS      |

Table 93: Patient Module

| Attribute Name              | Tag       | VR | Value | Presence of Value | Source    | Comment |
|-----------------------------|-----------|----|-------|-------------------|-----------|---------|
| Patient's Name              | 0010,0010 | PN |       | VNAP              | MWL, USER |         |
| Patient ID                  | 0010,0020 | LO |       | VNAP              | MWL, USER |         |
| Patient's Birth Date        | 0010,0030 | DA |       | VNAP              | MWL, USER |         |
| Patient's Sex               | 0010,0040 | CS |       | VNAP              | MWL, USER |         |
| Patient Comments            | 0010,4000 | LT |       | ANAP              | MWL, USER |         |
| Referenced Patient Sequence | 0008,1120 | SQ |       | ANAP              | MWL, USER |         |

|                              |           |    |  |        |      |  |
|------------------------------|-----------|----|--|--------|------|--|
| >Referenced SOP Class UID    | 0008,1150 | UI |  | ALWAYS | AUTO |  |
| >Referenced SOP Instance UID | 0008,1155 | UI |  | ALWAYS | AUTO |  |

Table 94: General Study Module

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source    | Comment                                  |
|------------------------------|-----------|----|-------|-------------------|-----------|--|
| Study Date                   | 0008,0020 | DA |       | VNAP              | AUTO      |  |
| Study Time                   | 0008,0030 | TM |       | VNAP              | AUTO      |  |
| Accession Number             | 0008,0050 | SH |       | VNAP              | MWL, USER |  |
| Referring Physician's Name   | 0008,0090 | PN |       | VNAP              | MWL, USER |  |
| Study Description            | 0008,1030 | LO |       | ANAP              | AUTO      |  |
| Procedure Code Sequence      | 0008,1032 | SQ |       | ANAP              | AUTO      |  |
| >Code Value                  | 0008,0100 | SH |       | ALWAYS            | AUTO      |  |
| >Coding Scheme Designator    | 0008,0102 | SH |       | ALWAYS            | AUTO      |  |
| >Code Meaning                | 0008,0104 | LO |       | ALWAYS            | AUTO      |  |
| Referenced Study Sequence    | 0008,1110 | SQ |       | ANAP              | AUTO      | Shall be present with Zero or more items |
| >Referenced SOP Class UID    | 0008,1150 | UI |       | ALWAYS            | AUTO      |  |
| >Referenced SOP Instance UID | 0008,1155 | UI |       | ALWAYS            | AUTO      |  |
| Study Instance UID           | 0020,000D | UI |       | ALWAYS            | AUTO      |  |
| Study ID                     | 0020,0010 | SH |       | VNAP              | AUTO      |  |

Table 95: Patient Study Module

| Attribute Name                  | Tag       | VR | Value | Presence of Value | Source    | Comment |
|---------------------------------|-----------|----|-------|-------------------|-----------|---------|
| Admitting Diagnoses Description | 0008,1080 | LO |       | ANAP              | MWL, USER |         |
| Patient's Age                   | 0010,1010 | AS |       | ANAP              | MWL, USER |         |
| Patient's Size                  | 0010,1020 | DS |       | ANAP              | MWL, USER |         |
| Patient's Weight                | 0010,1030 | DS |       | ANAP              | MWL, USER |         |

Table 96: General Series Module

| Attribute Name                               | Tag       | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|--------|---------|
| Series Date                                  | 0008,0021 | DA |       | ANAP              | AUTO   |         |
| Series Time                                  | 0008,0031 | TM |       | ANAP              | AUTO   |         |
| Modality                                     | 0008,0060 | CS | CT    | ALWAYS            | AUTO   |         |
| Series Description                           | 0008,103E | LO |       | ANAP              | AUTO   |         |
| Operators' Name                              | 0008,1070 | PN |       | ANAP              | AUTO   |         |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ |       | ANAP              | AUTO   |         |
| >Referenced SOP Class UID                    | 0008,1150 | UI |       | ALWAYS            | AUTO   |         |
| >Referenced SOP Instance                     | 0008,1155 | UI |       | ALWAYS            | AUTO   |         |

| Attribute Name                                | Tag       | VR | Value | Presence of Value | Source | Comment |
|---|-----------|----|-------|-------------------|--------|---------|
| UID   |           |    |       |                   |        |         |
| Body Part Examined                            | 0018,0015 | CS |       | ANAP              | AUTO   |         |
| Protocol Name                                 | 0018,1030 | LO |       | ANAP              | AUTO   |         |
| Patient Position                              | 0018,5100 | CS |       | ANAP              | AUTO   |         |
| Series Instance UID                           | 0020,000E | UI |       | ALWAYS            | AUTO   |         |
| Series Number                                 | 0020,0011 | IS |       | VNAP              | AUTO   |         |
| Performed Procedure Step Start Date           | 0040,0244 | DA |       | ANAP              | AUTO   |         |
| Performed Procedure Step Start Time           | 0040,0245 | TM |       | ANAP              | AUTO   |         |
| Performed Procedure Step ID                   | 0040,0253 | SH |       | ANAP              | AUTO   |         |
| Performed Procedure Step Description          | 0040,0254 | LO |       | ANAP              | AUTO   |         |
| Performed Protocol Code Sequence              | 0040,0260 | SQ |       | ANAP              | AUTO   |         |
| >Code Value                                   | 0008,0100 | SH |       | ALWAYS            | AUTO   |         |
| >Coding Scheme Designator                     | 0008,0102 | SH |       | ALWAYS            | AUTO   |         |
| >Code Meaning                                 | 0008,0104 | LO |       | ALWAYS            | AUTO   |         |
| Request Attributes Sequence                   | 0040,0275 | SQ |       | ANAP              | AUTO   |         |
| >Requested Procedure Description              | 0032,1060 | LO |       | ANAP              | AUTO   |         |
| >Requested Procedure Code Sequence            | 0032,1064 | SQ |       | ANAP              | AUTO   |         |
| >>Code Value                                  | 0008,0100 | SH |       | ALWAYS            | AUTO   |         |
| >>Coding Scheme Designator                    | 0008,0102 | SH |       | ALWAYS            | AUTO   |         |
| >>Code Meaning                                | 0008,0104 | LO |       | ALWAYS            | AUTO   |         |
| >Scheduled Procedure Step Description         | 0040,0007 | LO |       | ANAP              | AUTO   |         |
| >Scheduled Protocol Code Sequence             | 0040,0008 | SQ |       | ANAP              | AUTO   |         |
| >>Code Value                                  | 0008,0100 | SH |       | ALWAYS            | AUTO   |         |
| >>Coding Scheme Designator                    | 0008,0102 | SH |       | ALWAYS            | AUTO   |         |
| >>Code Meaning                                | 0008,0104 | LO |       | ALWAYS            | AUTO   |         |
| >Scheduled Procedure Step ID                  | 0040,0009 | SH |       | ANAP              | AUTO   |         |
| >Requested Procedure ID                       | 0040,1001 | SH |       | ANAP              | AUTO   |         |
| >Reason for the Requested Procedure           | 0040,1002 | LO |       | ALWAYS            | AUTO   |         |
| >Reason for Requested Procedure Code Sequence | 0040,100A | SQ |       | ALWAYS            | AUTO   |         |
| >>Code Value                                  | 0008,0100 | SH |       | ALWAYS            | AUTO   |         |
| >>Coding Scheme                               | 0008,0102 | SH |       | ALWAYS            | AUTO   |         |

| Attribute Name | Tag       | VR | Value | Presence of Value | Source | Comment |
|----------------|-----------|----|-------|-------------------|--------|---------|
| Designator     |           |    |       |                   |        |         |
| >>Code Meaning | 0008,0104 | LO |       | ALWAYS            | AUTO   |         |

Table 97: Frame of Reference Module

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|--------|---------|
| Frame of Reference UID       | 0020,0052 | UI |       | ALWAYS            | AUTO   |         |
| Position Reference Indicator | 0020,1040 | LO |       | VNAP              | AUTO   |         |

Table 98: General Equipment Module

| Attribute Name                | Tag       | VR | Value    | Presence of Value | Source       | Comment |
|-------------------------------|-----------|----|----------|-------------------|--------------|---------|
| Manufacturer                  | 0008,0070 | LO | Philips  | VNAP              | CONFIG       |         |
| Institution Name              | 0008,0080 | LO |          | ANAP              | CONFIG, USER |         |
| Institution Address           | 0008,0081 | ST |          | ANAP              | CONFIG, USER |         |
| Station Name                  | 0008,1010 | SH |          | ANAP              | CONFIG       |         |
| Institutional Department Name | 0008,1040 | LO |          | ANAP              | CONFIG       |         |
| Manufacturer's Model Name     | 0008,1090 | LO | Big Bore | ANAP              | FIXED        |         |
| Device Serial Number          | 0018,1000 | LO |          | ANAP              | AUTO, CONFIG |         |
| Software Version(s)           | 0018,1020 | LO | 4.8      | ANAP              | FIXED        |         |

Table 99: General Acquisition Module

| Attribute Name        | Tag       | VR | Value | Presence of Value | Source | Comment |
|-----------------------|-----------|----|-------|-------------------|--------|---------|
| Acquisition Date      | 0008,0022 | DA |       | ANAP              | AUTO   |         |
| Acquisition Time      | 0008,0032 | TM |       | ANAP              | AUTO   |         |
| Acquisition DateTime  | 0008,002A | DT |       | ANAP              | AUTO   |         |
| Irradiation Event UID | 0008,3010 | UI |       | ANAP              | AUTO   |         |
| Acquisition Number    | 0020,0012 | IS |       | ANAP              | AUTO   |         |

Table 100: General Image Module

| Attribute Name          | Tag       | VR | Value | Presence of Value | Source | Comment |
|-------------------------|-----------|----|-------|-------------------|--------|---------|
| Image Type              | 0008,0008 | CS |       | ANAP              | AUTO   |         |
| Content Date            | 0008,0023 | DA |       | ANAP              | AUTO   |         |
| Content Time            | 0008,0033 | TM |       | ANAP              | AUTO   |         |
| Instance Number         | 0020,0013 | IS |       | VNAP              | AUTO   |         |
| Image Comments          | 0020,4000 | LT |       | ANAP              | AUTO   |         |
| Lossy Image Compression | 0028,2110 | CS |       | ANAP              | AUTO   |         |
| Quality Control Image   | 0028,0300 | CS |       | ANAP              | AUTO   |         |

Table 101: General Reference Module

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source | Comment                             |
|------------------------------|-----------|----|-------|-------------------|--------|-------------------------------------|
| Referenced Image Sequence    | 0008,1140 | SQ |       | ANAP              | AUTO   |                                     |
| >Referenced SOP Class UID    | 0008,1150 | UI |       | ALWAYS            | AUTO   |                                     |
| >Referenced SOP Instance UID | 0008,1155 | UI |       | ALWAYS            | AUTO   |                                     |
| Source Image Sequence        | 0008,2112 | SQ |       | ANAP              | AUTO   | Only present when DIN is configured |
| >Referenced SOP Class UID    | 0008,1150 | UI |       | ALWAYS            | AUTO   |                                     |
| >Referenced SOP Instance UID | 0008,1155 | UI |       | ALWAYS            | AUTO   |                                     |

Table 102: Image Plane Module

| Attribute Name              | Tag       | VR | Value | Presence of Value | Source | Comment |
|-----------------------------|-----------|----|-------|-------------------|--------|---------|
| Slice Thickness             | 0018,0050 | DS |       | VNAP              | AUTO   |         |
| Spacing Between Slices      | 0018,0088 | DS |       | ANAP              | AUTO   |         |
| Image Position (Patient)    | 0020,0032 | DS |       | ALWAYS            | AUTO   |         |
| Image Orientation (Patient) | 0020,0037 | DS |       | ALWAYS            | AUTO   |         |
| Slice Location              | 0020,1041 | DS |       | ANAP              | AUTO   |         |
| Pixel Spacing               | 0028,0030 | DS |       | ALWAYS            | AUTO   |         |

Table 103: Image Pixel Module

| Attribute Name             | Tag       | VR    | Value | Presence of Value | Source | Comment |
|----------------------------|-----------|-------|-------|-------------------|--------|---------|
| Samples per Pixel          | 0028,0002 | US    |       | ALWAYS            | AUTO   |         |
| Photometric Interpretation | 0028,0004 | CS    |       | ALWAYS            | AUTO   |         |
| Rows                       | 0028,0010 | US    |       | ALWAYS            | AUTO   |         |
| Columns                    | 0028,0011 | US    |       | ALWAYS            | AUTO   |         |
| Bits Allocated             | 0028,0100 | US    | 16    | ALWAYS            | AUTO   |         |
| Bits Stored                | 0028,0101 | US    | 12    | ALWAYS            | AUTO   |         |
| High Bit                   | 0028,0102 | US    | 11    | ALWAYS            | AUTO   |         |
| Pixel Representation       | 0028,0103 | US    | 0     | ALWAYS            | FIXED  |         |
| Pixel Data                 | 7FE0,0010 | OW/OB |       | ANAP              | AUTO   |         |

Table 104: CT Image Module

| Attribute Name | Tag       | VR | Value | Presence of Value | Source | Comment |
|----------------|-----------|----|-------|-------------------|--------|---------|
| Image Type     | 0008,0008 | CS |       | ALWAYS            | AUTO   |         |
| Scan Options   | 0018,0022 | CS |       | ANAP              | AUTO   |         |
| KVP            | 0018,0060 | DS |       | VNAP              | AUTO   |         |

|  |           |    |       |        |       |  |
|--|-----------|----|-------|--------|-------|--|
| Filter Type                            | 0018,1160 | SH |       | ALWAYS | AUTO  |  |
| Data Collection Diameter               | 0018,0090 | DS |       | ANAP   | AUTO  |  |
| Reconstruction Diameter                | 0018,1100 | DS |       | ANAP   | AUTO  |  |
| Distance Source to Detector            | 0018,1110 | DS |       | ANAP   | AUTO  |  |
| Distance Source to Patient             | 0018,1111 | DS |       | ANAP   | AUTO  |  |
| Gantry/Detector Tilt                   | 0018,1120 | DS |       | ANAP   | AUTO  |  |
| Table Height                           | 0018,1130 | DS |       | ANAP   | AUTO  |  |
| Exposure Time                          | 0018,1150 | IS |       | ANAP   | AUTO  |  |
| X-ray Tube Current                     | 0018,1151 | IS |       | ANAP   | AUTO  |  |
| Exposure                               | 0018,1152 | IS |       | ANAP   | AUTO  |  |
| Convolution Kernel                     | 0018,1210 | SH |       | ANAP   | AUTO  |  |
| Revolution Time                        | 0018,9305 | FD |       | ANAP   | AUTO  |  |
| Single Collimation Width               | 0018,9306 | FD |       | ANAP   | AUTO  |  |
| Total Collimation Width                | 0018,9307 | FD |       | ANAP   | AUTO  |  |
| Table Speed                            | 0018,9309 | FD |       | ANAP   | AUTO  |  |
| Table Feed per Rotation                | 0018,9310 | FD |       | ANAP   | AUTO  |  |
| Spiral Pitch Factor                    | 0018,9311 | FD |       | ANAP   | AUTO  |  |
| Data Collection Center (Patient)       | 0018,9313 | FD |       | ANAP   | AUTO  |  |
| Reconstruction Target Center (Patient) | 0018,9318 | FD |       | ANAP   | AUTO  |  |
| Exposure Modulation Type               | 0018,9323 | CS |       | ANAP   | AUTO  |  |
| Estimated Dose Saving                  | 0018,9324 | FD |       | ANAP   | AUTO  |  |
| CTDIvol                                | 0018,9345 | FD |       | ANAP   | AUTO  |  |
| CTDI Phantom Type Code Sequence        | 0018,9346 | SQ |       | ANAP   | AUTO  |  |
| >Code Value                            | 0008,0100 | SH |       | ALWAYS | AUTO  |  |
| >Coding Scheme Designator              | 0008,0102 | SH |       | ALWAYS | AUTO  |  |
| >Coding Scheme Version                 | 0008,0103 | SH |       | ALWAYS | AUTO  |  |
| >Code Meaning                          | 0008,0104 | LO |       | ALWAYS | AUTO  |  |
| CT Additional X-Ray Source Sequence    | 0018,9360 | SQ |       | ANAP   | AUTO  |  |
| >X-Ray Tube Current in mA              | 0018,9330 | FD |       | ALWAYS | AUTO  |  |
| Acquisition Number                     | 0020,0012 | IS |       | VNAP   | AUTO  |  |
| Samples per Pixel                      | 0028,0002 | US |       | ALWAYS | AUTO  |  |
| Photometric Interpretation             | 0028,0004 | CS |       | ALWAYS | AUTO  |  |
| Bits Allocated                         | 0028,0100 | US | 16    | ALWAYS | AUTO  |  |
| Bits Stored                            | 0028,0101 | US | 12    | ALWAYS | AUTO  |  |
| High Bit                               | 0028,0102 | US | 11    | ALWAYS | AUTO  |  |
| Rescale Intercept                      | 0028,1052 | DS | -1024 | ALWAYS | FIXED |  |
| Rescale Slope                          | 0028,1053 | DS | 1     | ALWAYS | FIXED |  |

Table 105: VOI LUT Module

| Attribute Name | Tag       | VR | Value | Presence of Value | Source | Comment |
|----------------|-----------|----|-------|-------------------|--------|---------|
| Window Center  | 0028,1050 | DS |       | ANAP              | AUTO   |         |
| Window Width   | 0028,1051 | DS |       | ANAP              | AUTO   |         |

Table 106: SOP Common Module

| Attribute Name                      | Tag       | VR | Value                     | Presence of Value | Source | Comment |
|-------------------------------------|-----------|----|---------------------------|-------------------|--------|---------|
| Specific Character Set              | 0008,0005 | CS |                           | ANAP              | AUTO   |         |
| Instance Creation Date              | 0008,0012 | DA |                           | ANAP              | AUTO   |         |
| Instance Creation Time              | 0008,0013 | TM |                           | ANAP              | AUTO   |         |
| SOP Class UID                       | 0008,0016 | UI | 1.2.840.10008.5.1.4.1.1.2 | ALWAYS            | FIXED  |         |
| SOP Instance UID                    | 0008,0018 | UI |                           | ALWAYS            | AUTO   |         |
| Instance Number                     | 0020,0013 | IS |                           | ANAP              | AUTO   |         |
| Contributing Equipment Sequence     | 0018,A001 | SQ |                           | ANAP              | AUTO   |         |
| >Manufacturer                       | 0008,0070 | LO |                           | ALWAYS            | AUTO   |         |
| >Institution Name                   | 0008,0080 | LO |                           | ANAP              | AUTO   |         |
| >Institution Address                | 0008,0081 | ST |                           | ANAP              | AUTO   |         |
| >Station Name                       | 0008,1010 | SH |                           | ANAP              | AUTO   |         |
| >Institutional Department Name      | 0008,1040 | LO |                           | ANAP              | AUTO   |         |
| >Manufacturer's Model Name          | 0008,1090 | LO |                           | ANAP              | AUTO   |         |
| >Device Serial Number               | 0018,1000 | LO |                           | ANAP              | AUTO   |         |
| >Software Version(s)                | 0018,1020 | LO |                           | ANAP              | AUTO   |         |
| >Purpose of Reference Code Sequence | 0040,A170 | SQ |                           | ALWAYS            | AUTO   |         |
| >>Code Value                        | 0008,0100 | SH |                           | ALWAYS            | AUTO   |         |
| >>Coding Scheme Designator          | 0008,0102 | SH |                           | ALWAYS            | AUTO   |         |
| >>Code Meaning                      | 0008,0104 | LO |                           | ALWAYS            | AUTO   |         |

### 8.1.1.1.2 General ECG Waveform Storage SOP Class

Table 107: SOP Class Modules

| Information Entity | Module                         | Presence    |
|--------------------|--------------------------------|-------------|
| Patient            | Patient Module                 | ALWAYS      |
| Study              | General Study Module           | ALWAYS      |
|                    | Patient Study Module           | CONDITIONAL |
| Series             | General Series Module          | ALWAYS      |
| Equipment          | General Equipment Module       | ALWAYS      |
| Waveform           | Waveform Identification Module | ALWAYS      |
|                    | Waveform Module                | ALWAYS      |
|                    | Acquisition Context Module     | ALWAYS      |
|                    | SOP Common Module              | ALWAYS      |

Table 108: Patient Module

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source    | Comment |
|------------------------------|-----------|----|-------|-------------------|-----------|---------|
| Patient's Name               | 0010,0010 | PN |       | VNAP              | USER      |         |
| Patient ID                   | 0010,0020 | LO |       | VNAP              | USER      |         |
| Patient's Birth Date         | 0010,0030 | DA |       | VNAP              | USER      |         |
| Patient's Sex                | 0010,0040 | CS |       | VNAP              | USER      |         |
| Patient Comments             | 0010,4000 | LT |       | ANAP              | MWL, USER |         |
| Referenced Patient Sequence  | 0008,1120 | SQ |       | ANAP              | MWL, USER |         |
| >Referenced SOP Class UID    | 0008,1150 | UI |       | ALWAYS            | AUTO      |         |
| >Referenced SOP Instance UID | 0008,1155 | UI |       | ALWAYS            | AUTO      |         |

Table 109: General Study Module

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source | Comment                                  |
|------------------------------|-----------|----|-------|-------------------|--------|--|
| Study Date                   | 0008,0020 | DA |       | VNAP              | AUTO   |  |
| Study Time                   | 0008,0030 | TM |       | VNAP              | AUTO   |  |
| Accession Number             | 0008,0050 | SH |       | VNAP              | USER   |  |
| Referring Physician's Name   | 0008,0090 | PN |       | VNAP              | USER   |  |
| Study Description            | 0008,1030 | LO |       | ANAP              | USER   |  |
| Study Instance UID           | 0020,000D | UI |       | ALWAYS            | AUTO   |  |
| Study ID                     | 0020,0010 | SH |       | VNAP              | AUTO   |  |
| Referenced Study Sequence    | 0008,1110 | SQ |       | ANAP              | AUTO   | Shall be present with Zero or more items |
| >Referenced SOP Class UID    | 0008,1150 | UI |       | ALWAYS            | AUTO   |  |
| >Referenced SOP Instance UID | 0008,1155 | UI |       | ALWAYS            | AUTO   |  |

Table 110: Patient Study Module

| Attribute Name                  | Tag       | VR | Value | Presence of Value | Source    | Comment |
|---------------------------------|-----------|----|-------|-------------------|-----------|---------|
| Admitting Diagnoses Description | 0008,1080 | LO |       | ANAP              | MWL, USER |         |
| Patient's Age                   | 0010,1010 | AS |       | ANAP              | MWL, USER |         |
| Patient's Size                  | 0010,1020 | DS |       | ANAP              | MWL, USER |         |
| Patient's Weight                | 0010,1030 | DS |       | ANAP              | MWL, USER |         |

Table 111: General Series Module

| Attribute Name | Tag       | VR | Value | Presence of Value | Source | Comment |
|----------------|-----------|----|-------|-------------------|--------|---------|
| Series Date    | 0008,0021 | DA |       | ANAP              | AUTO   |         |
| Series Time    | 0008,0031 | TM |       | ANAP              | AUTO   |         |

| Attribute Name                               | Tag       | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|--------|---------|
| Modality                                     | 0008,0060 | CS | ECG   | ALWAYS            | FIXED  |         |
| Series Description                           | 0008,103E | LO |       | ANAP              | USER   |         |
| Operators' Name                              | 0008,1070 | PN |       | ANAP              | USER   |         |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ |       | ANAP              | AUTO   |         |
| >Referenced SOP Class UID                    | 0008,1150 | UI |       | ALWAYS            | AUTO   |         |
| >Referenced SOP Instance UID                 | 0008,1155 | UI |       | ALWAYS            | AUTO   |         |
| Body Part Examined                           | 0018,0015 | CS |       | ANAP              | AUTO   |         |
| Protocol Name                                | 0018,1030 | LO |       | ANAP              | AUTO   |         |
| Patient Position                             | 0018,5100 | CS |       | ANAP              | AUTO   |         |
| Series Instance UID                          | 0020,000E | UI |       | ALWAYS            | AUTO   |         |
| Series Number                                | 0020,0011 | IS |       | VNAP              | AUTO   |         |
| Performed Procedure Step Start Date          | 0040,0244 | DA |       | ANAP              | AUTO   |         |
| Performed Procedure Step Start Time          | 0040,0245 | TM |       | ANAP              | AUTO   |         |
| Performed Procedure Step ID                  | 0040,0253 | SH |       | ANAP              | AUTO   |         |
| Performed Procedure Step Description         | 0040,0254 | LO |       | ANAP              | AUTO   |         |
| Performed Protocol Code Sequence             | 0040,0260 | SQ |       | ANAP              | AUTO   |         |
| >Code Value                                  | 0008,0100 | SH |       | ALWAYS            | AUTO   |         |
| >Coding Scheme Designator                    | 0008,0102 | SH |       | ALWAYS            | AUTO   |         |
| >Code Meaning                                | 0008,0104 | LO |       | ALWAYS            | AUTO   |         |
| Request Attributes Sequence                  | 0040,0275 | SQ |       | ANAP              | AUTO   |         |
| >Requested Procedure Description             | 0032,1060 | LO |       | ANAP              | AUTO   |         |
| >Requested Procedure Code Sequence           | 0032,1064 | SQ |       | ANAP              | AUTO   |         |
| >>Code Value                                 | 0008,0100 | SH |       | ALWAYS            | AUTO   |         |
| >>Coding Scheme Designator                   | 0008,0102 | SH |       | ALWAYS            | AUTO   |         |
| >>Code Meaning                               | 0008,0104 | LO |       | ALWAYS            | AUTO   |         |
| >Scheduled Procedure Step Description        | 0040,0007 | LO |       | ANAP              | AUTO   |         |
| >Scheduled Protocol Code Sequence            | 0040,0008 | SQ |       | ANAP              | AUTO   |         |
| >>Code Value                                 | 0008,0100 | SH |       | ALWAYS            | AUTO   |         |
| >>Coding Scheme Designator                   | 0008,0102 | SH |       | ALWAYS            | AUTO   |         |
| >>Code Meaning                               | 0008,0104 | LO |       | ALWAYS            | AUTO   |         |

| Attribute Name                                | Tag       | VR | Value | Presence of Value | Source | Comment |
|---|-----------|----|-------|-------------------|--------|---------|
| >Scheduled Procedure Step ID                  | 0040,0009 | SH |       | ANAP              | AUTO   |         |
| >Requested Procedure ID                       | 0040,1001 | SH |       | ANAP              | AUTO   |         |
| >Reason for the Requested Procedure           | 0040,1002 | LO |       | ALWAYS            | AUTO   |         |
| >Reason for Requested Procedure Code Sequence | 0040,100A | SQ |       | ALWAYS            | AUTO   |         |
| >>Code Value                                  | 0008,0100 | SH |       | ALWAYS            | AUTO   |         |
| >>Coding Scheme Designator                    | 0008,0102 | SH |       | ALWAYS            | AUTO   |         |
| >>Code Meaning                                | 0008,0104 | LO |       | ALWAYS            | AUTO   |         |

Table 112: General Equipment Module

| Attribute Name                | Tag       | VR | Value    | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|----------|-------------------|--------|---------|
| Manufacturer                  | 0008,0070 | LO | Philips  | VNAP              | FIXED  |         |
| Institution Name              | 0008,0080 | LO |          | ANAP              | CONFIG |         |
| Institution Address           | 0008,0081 | ST |          | ANAP              | CONFIG |         |
| Station Name                  | 0008,1010 | SH |          | ANAP              | CONFIG |         |
| Institutional Department Name | 0008,1040 | LO |          | ANAP              | CONFIG |         |
| Manufacturer's Model Name     | 0008,1090 | LO | Big Bore | ANAP              | FIXED  |         |
| Device Serial Number          | 0018,1000 | LO |          | ANAP              | CONFIG |         |
| Software Version(s)           | 0018,1020 | LO | 4.8      | ANAP              | FIXED  |         |

Table 113: Waveform Identification Module

| Attribute Name        | Tag       | VR | Value | Presence of Value | Source | Comment |
|-----------------------|-----------|----|-------|-------------------|--------|---------|
| Content Date          | 0008,0023 | DA |       | ALWAYS            | AUTO   |         |
| Acquisition Date time | 0008,002A | DT |       | ALWAYS            | AUTO   |         |
| Content Time          | 0008,0033 | TM |       | ALWAYS            | AUTO   |         |
| Instance Number       | 0020,0013 | IS |       | ALWAYS            | AUTO   |         |

Table 114: Waveform Module

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|--------|---------|
| Waveform Sequence            | 5400,0100 | SQ |       | ALWAYS            | AUTO   |         |
| >Multiplex Group Time Offset | 0018,1068 | DS |       | ANAP              | AUTO   |         |
| >Waveform Originality        | 003A,0004 | CS |       | ALWAYS            | AUTO   |         |
| >Number of Waveform Channels | 003A,0005 | US |       | ALWAYS            | AUTO   |         |
| >Number of Waveform          | 003A,0010 | UL |       | ALWAYS            | AUTO   |         |

|                                 |           |       |    |        |      |  |
|---------------------------------|-----------|-------|----|--------|------|--|
| Samples                         |           |       |    |        |      |  |
| >Sampling Frequency             | 003A,001A | DS    |    | ALWAYS | AUTO |  |
| >Channel Definition Sequence    | 003A,0200 | SQ    |    | ALWAYS | AUTO |  |
| >>Channel Source Sequence       | 003A,0208 | SQ    |    | ALWAYS | AUTO |  |
| >>>Code Value                   | 0008,0100 | SH    |    | ALWAYS | AUTO |  |
| >>>Coding Scheme Designator     | 0008,0102 | SH    |    | ALWAYS | AUTO |  |
| >>>Code Meaning                 | 0008,0104 | LO    |    | ALWAYS | AUTO |  |
| >>Channel Sample Skew           | 003A,0215 | DS    |    | ANAP   | AUTO |  |
| >>Waveform Bits Stored          | 003A,021A | US    | 8  | ALWAYS | AUTO |  |
| >Waveform Bits Allocated        | 5400,1004 | US    | 16 | ALWAYS | AUTO |  |
| >Waveform Sample Interpretation | 5400,1006 | CS    |    | ALWAYS | AUTO |  |
| >Waveform Data                  | 5400,1010 | OW/OB |    | ALWAYS | AUTO |  |

Table 115: Acquisition Context Module

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|--------|---------|
| Acquisition Context Sequence | 0040,0555 | SQ |       | VNAP              | AUTO   |         |

Table 116: SOP Common Module

| Attribute Name         | Tag       | VR | Value                         | Presence of Value | Source | Comment |
|------------------------|-----------|----|-------------------------------|-------------------|--------|---------|
| Specific Character set | 0008,0005 | CS |                               | ANAP              | COPY   |         |
| SOP Class UID          | 0008,0016 | UI | 1.2.840.10008.5.1.4.1.1.9.1.2 | ALWAYS            | FIXED  |         |
| SOP Instance UID       | 0008,0018 | UI |                               | ALWAYS            | AUTO   |         |
| Instance Number        | 0020,0013 | IS |                               | ALWAYS            | AUTO   |         |

### 8.1.1.1.3 Secondary Capture Image Storage SOP Class

Table 117: SOP Class Modules

| Information Entity | Module                     | Presence    |
|--------------------|----------------------------|-------------|
| Patient            | Patient Module             | ALWAYS      |
| Study              | General Study Module       | ALWAYS      |
|                    | Patient Study Module       | CONDITIONAL |
| Series             | General Series Module      | ALWAYS      |
| Equipment          | SC Equipment Module        | ALWAYS      |
|                    | General Equipment Module   | ALWAYS      |
| Acquisition        | General Acquisition Module | ALWAYS      |
| Image              | General Image Module       | ALWAYS      |
|                    | Image Pixel Module         | ALWAYS      |
|                    | SC Image Module            | ALWAYS      |
|                    | Modality LUT Module        | CONDITIONAL |
|                    | VOI LUT Module             | CONDITIONAL |
|                    | SOP Common Module          | ALWAYS      |

Table 118: Patient Module

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source    | Comment |
|------------------------------|-----------|----|-------|-------------------|-----------|---------|
| Patient's Name               | 0010,0010 | PN |       | VNAP              | COPY      |         |
| Patient ID                   | 0010,0020 | LO |       | VNAP              | COPY      |         |
| Patient's Birth Date         | 0010,0030 | DA |       | VNAP              | COPY      |         |
| Patient's Sex                | 0010,0040 | CS |       | VNAP              | COPY      |         |
| Patient Comments             | 0010,4000 | LT |       | ANAP              | MWL, USER |         |
| Referenced Patient Sequence  | 0008,1120 | SQ |       | ANAP              | MWL, USER |         |
| >Referenced SOP Class UID    | 0008,1150 | UI |       | ALWAYS            | AUTO      |         |
| >Referenced SOP Instance UID | 0008,1155 | UI |       | ALWAYS            | AUTO      |         |

Table 119: General Study Module

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source | Comment                                  |
|------------------------------|-----------|----|-------|-------------------|--------|--|
| Study Date                   | 0008,0020 | DA |       | VNAP              | COPY   |  |
| Study Time                   | 0008,0030 | TM |       | VNAP              | COPY   |  |
| Accession Number             | 0008,0050 | SH |       | VNAP              | COPY   |  |
| Referring Physician's Name   | 0008,0090 | PN |       | VNAP              | COPY   |  |
| Study Description            | 0008,1030 | LO |       | ANAP              | COPY   |  |
| Procedure Code Sequence      | 0008,1032 | SQ |       | ANAP              | AUTO   |  |
| >Code Value                  | 0008,0100 | SH |       | ALWAYS            | AUTO   |  |
| >Coding Scheme Designator    | 0008,0102 | SH |       | ALWAYS            | AUTO   |  |
| >Code Meaning                | 0008,0104 | LO |       | ALWAYS            | AUTO   |  |
| Referenced Study Sequence    | 0008,1110 | SQ |       | ANAP              | AUTO   | Shall be present with Zero or more items |
| >Referenced SOP Class UID    | 0008,1150 | UI |       | ALWAYS            | AUTO   |  |
| >Referenced SOP Instance UID | 0008,1155 | UI |       | ALWAYS            | AUTO   |  |
| Study Instance UID           | 0020,000D | UI |       | ALWAYS            | COPY   |  |
| Study ID                     | 0020,0010 | SH |       | VNAP              | COPY   |  |

Table 120: Patient Study Module

| Attribute Name                  | Tag       | VR | Value | Presence of Value | Source    | Comment |
|---------------------------------|-----------|----|-------|-------------------|-----------|---------|
| Admitting Diagnoses Description | 0008,1080 | LO |       | ANAP              | MWL, USER |         |
| Patient's Age                   | 0010,1010 | AS |       | ANAP              | MWL, USER |         |
| Patient's Size                  | 0010,1020 | DS |       | ANAP              | MWL, USER |         |
| Patient's Weight                | 0010,1030 | DS |       | ANAP              | MWL, USER |         |

Table 121: General Series Module

| Attribute Name                               | Tag       | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|--------|---------|
| Series Date                                  | 0008,0021 | DA |       | ANAP              | AUTO   |         |
| Series Time                                  | 0008,0031 | TM |       | ANAP              | AUTO   |         |
| Series Description                           | 0008,103E | LO |       | ANAP              | AUTO   |         |
| Modality                                     | 0008,0060 | CS | CT    | ALWAYS            | AUTO   |         |
| Operators' Name                              | 0008,1070 | PN |       | ANAP              | AUTO   |         |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ |       | ANAP              | AUTO   |         |
| >Referenced SOP Class UID                    | 0008,1150 | UI |       | ALWAYS            | AUTO   |         |
| >Referenced SOP Instance UID                 | 0008,1155 | UI |       | ALWAYS            | AUTO   |         |
| Body Part Examined                           | 0018,0015 | CS |       | ANAP              | AUTO   |         |
| Protocol Name                                | 0018,1030 | LO |       | ANAP              | AUTO   |         |
| Patient Position                             | 0018,5100 | CS |       | ANAP              | AUTO   |         |
| Series Instance UID                          | 0020,000E | UI |       | ALWAYS            | AUTO   |         |
| Series Number                                | 0020,0011 | IS |       | VNAP              | AUTO   |         |
| Performed Procedure Step Start Date          | 0040,0244 | DA |       | ANAP              | AUTO   |         |
| Performed Procedure Step Start Time          | 0040,0245 | TM |       | ANAP              | AUTO   |         |
| Performed Procedure Step ID                  | 0040,0253 | SH |       | ANAP              | AUTO   |         |
| Performed Procedure Step Description         | 0040,0254 | LO |       | ANAP              | AUTO   |         |
| Performed Protocol Code Sequence             | 0040,0260 | SQ |       | ANAP              | AUTO   |         |
| >Code Value                                  | 0008,0100 | SH |       | ALWAYS            | AUTO   |         |
| >Coding Scheme Designator                    | 0008,0102 | SH |       | ALWAYS            | AUTO   |         |
| >Code Meaning                                | 0008,0104 | LO |       | ALWAYS            | AUTO   |         |
| Request Attributes Sequence                  | 0040,0275 | SQ |       | ANAP              | AUTO   |         |
| >Requested Procedure Description             | 0032,1060 | LO |       | ANAP              | AUTO   |         |
| >Requested Procedure Code Sequence           | 0032,1064 | SQ |       | ANAP              | AUTO   |         |
| >>Code Value                                 | 0008,0100 | SH |       | ALWAYS            | AUTO   |         |
| >>Coding Scheme Designator                   | 0008,0102 | SH |       | ALWAYS            | AUTO   |         |
| >>Code Meaning                               | 0008,0104 | LO |       | ALWAYS            | AUTO   |         |
| >Scheduled Procedure Step Description        | 0040,0007 | LO |       | ANAP              | AUTO   |         |
| >Scheduled Protocol Code Sequence            | 0040,0008 | SQ |       | ANAP              | AUTO   |         |
| >Code Value                                  | 0008,0100 | SH |       | ALWAYS            | AUTO   |         |
| >Coding Scheme Designator                    | 0008,0102 | SH |       | ALWAYS            | AUTO   |         |
| >Code Meaning                                | 0008,0104 | LO |       | ALWAYS            | AUTO   |         |

| Attribute Name                                | Tag       | VR | Value | Presence of Value | Source | Comment |
|---|-----------|----|-------|-------------------|--------|---------|
| >Scheduled Procedure Step ID                  | 0040,0009 | SH |       | ANAP              | AUTO   |         |
| >Requested Procedure ID                       | 0040,1001 | SH |       | ANAP              | AUTO   |         |
| >Reason for the Requested Procedure           | 0040,1002 | LO |       | ALWAYS            | AUTO   |         |
| >Reason for Requested Procedure Code Sequence | 0040,100A | SQ |       | ALWAYS            | AUTO   |         |
| >>Code Value                                  | 0008,0100 | SH |       | ALWAYS            | AUTO   |         |
| >>Coding Scheme Designator                    | 0008,0102 | SH |       | ALWAYS            | AUTO   |         |
| >>Code Meaning                                | 0008,0104 | LO |       | ALWAYS            | AUTO   |         |

Table 122: General Equipment Module

| Attribute Name                | Tag       | VR | Value    | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|----------|-------------------|--------|---------|
| Manufacturer                  | 0008,0070 | LO | Philips  | VNAP              | FIXED  |         |
| Institution Name              | 0008,0080 | LO |          | ANAP              | COPY   |         |
| Institution Address           | 0008,0081 | ST |          | ANAP              | COPY   |         |
| Station Name                  | 0008,1010 | SH |          | ANAP              | COPY   |         |
| Institutional Department Name | 0008,1040 | LO |          | ANAP              | COPY   |         |
| Manufacturer's Model Name     | 0008,1090 | LO | Big Bore | ANAP              | FIXED  |         |
| Device Serial Number          | 0018,1000 | LO |          | ANAP              | COPY   |         |
| Software Version(s)           | 0018,1020 | LO | 4.8      | ANAP              | FIXED  |         |

Table 123: SC Equipment Module

| Attribute Name                                     | Tag       | VR | Value    | Presence of Value | Source | Comment |
|--|-----------|----|----------|-------------------|--------|---------|
| Modality   | 0008,0060 | CS | CT       | ANAP              | FIXED  |         |
| Conversion Type                                    | 0008,0064 | CS | WSD      | ALWAYS            | AUTO   |         |
| Secondary Capture Device Manufacturer              | 0018,1016 | LO | Philips  | ANAP              | FIXED  |         |
| Secondary Capture Device Manufacturer's Model Name | 0018,1018 | LO | Big Bore | ANAP              | FIXED  |         |
| Secondary Capture Device Software Version(s)       | 0018,1019 | LO | 4.8      | ANAP              | FIXED  |         |

Table 124: General Acquisition Module

| Attribute Name       | Tag       | VR | Value | Presence of Value | Source | Comment |
|----------------------|-----------|----|-------|-------------------|--------|---------|
| Acquisition Date     | 0008,0022 | DA |       | ANAP              | AUTO   |         |
| Acquisition Time     | 0008,0032 | TM |       | ANAP              | AUTO   |         |
| Acquisition DateTime | 0008,002A | DT |       | ANAP              | AUTO   |         |

|                       |           |    |  |      |      |  |
|-----------------------|-----------|----|--|------|------|--|
| Irradiation Event UID | 0008,3010 | UI |  | ANAP | AUTO |  |
| Acquisition Number    | 0020,0012 | IS |  | ANAP | AUTO |  |

Table 125: General Image Module

| Attribute Name          | Tag       | VR | Value                                   | Presence of Value | Source | Comment                                     |
|-------------------------|-----------|----|---|-------------------|--------|---|
| Image Type              | 0008,0008 | CS | Value 1: DERIVED,<br>Value 2: SECONDARY | ANAP              | AUTO   |   |
| Content Date            | 0008,0023 | DA |   | ANAP              | AUTO   |   |
| Content Time            | 0008,0033 | TM |   | ANAP              | AUTO   |   |
| Instance Number         | 0020,0013 | IS |   | VNAP              | AUTO   |   |
| Patient Orientation     | 0020,0020 | CS |   | ANAP              | AUTO   |   |
| Image Comments          | 0020,4000 | LT |   | ANAP              | AUTO   |   |
| Lossy Image Compression | 0028,2110 | CS |   | ANAP              | AUTO   |   |
| Quality Control Image   | 0028,0300 | CS |   | ANAP              | AUTO   |   |
| Burned In Annotation    | 0028,0301 | CS |   | ANAP              | AUTO   | For Patient Info<br>"YES" for Other<br>"No" |

Table 126: Image Pixel Module

| Attribute Name             | Tag       | VR    | Value               | Presence of Value | Source | Comment     |
|----------------------------|-----------|-------|---------------------|-------------------|--------|-------------|
| Samples per Pixel          | 0028,0002 | US    | 1, 3                | ALWAYS            | AUTO   |             |
| Photometric Interpretation | 0028,0004 | CS    | MONOCHROME2,<br>RGB | ALWAYS            | AUTO   |             |
| Planar Configuration       | 0028,0006 | US    | 1                   | ANAP              | FIXED  |             |
| Rows                       | 0028,0010 | US    | 512                 | ALWAYS            | AUTO   | 512 or more |
| Columns                    | 0028,0011 | US    | 512                 | ALWAYS            | AUTO   |             |
| Bits Allocated             | 0028,0100 | US    | 16, 8               | ALWAYS            | AUTO   |             |
| Bits Stored                | 0028,0101 | US    | 12, 8               | ALWAYS            | AUTO   |             |
| High Bit                   | 0028,0102 | US    | 11, 7               | ALWAYS            | AUTO   |             |
| Pixel Representation       | 0028,0103 | US    | 0                   | ALWAYS            | FIXED  |             |
| Pixel Data                 | 7FE0,0010 | OW/OB |                     | ANAP              | AUTO   |             |

Table 127: SC Image Module

| Attribute Name            | Tag       | VR | Value | Presence of Value | Source | Comment |
|---------------------------|-----------|----|-------|-------------------|--------|---------|
| Date of Secondary Capture | 0018,1012 | DA |       | ANAP              | AUTO   |         |
| Time of Secondary Capture | 0018,1014 | TM |       | ANAP              | AUTO   |         |
| Pixel Spacing             | 0028,0030 | DS |       | ANAP              | AUTO   |         |

Table 128: Modality LUT Module

| Attribute Name    | Tag       | VR | Value | Presence of Value | Source | Comment |
|-------------------|-----------|----|-------|-------------------|--------|---------|
| Rescale Intercept | 0028,1052 | DS |       | ANAP              | AUTO   |         |

|               |           |    |  |      |      |  |
|---------------|-----------|----|--|------|------|--|
| Rescale Slope | 0028,1053 | DS |  | ANAP | AUTO |  |
| Rescale Type  | 0028,1054 | LO |  | ANAP | AUTO |  |

Table 129: VOI LUT Module

| Attribute Name | Tag       | VR | Value | Presence of Value | Source | Comment |
|----------------|-----------|----|-------|-------------------|--------|---------|
| Window Center  | 0028,1050 | DS |       | ANAP              | AUTO   |         |
| Window Width   | 0028,1051 | DS |       | ANAP              | AUTO   |         |

Table 130: SOP Common Module

| Attribute Name                      | Tag       | VR | Value                     | Presence of Value | Source | Comment |
|-------------------------------------|-----------|----|---------------------------|-------------------|--------|---------|
| Specific Character Set              | 0008,0005 | CS |                           | ANAP              | COPY   |         |
| Instance Creation Date              | 0008,0012 | DA |                           | ANAP              | AUTO   |         |
| Instance Creation Time              | 0008,0013 | TM |                           | ANAP              | AUTO   |         |
| SOP Class UID                       | 0008,0016 | UI | 1.2.840.10008.5.1.4.1.1.7 | ALWAYS            | FIXED  |         |
| SOP Instance UID                    | 0008,0018 | UI |                           | ALWAYS            | AUTO   |         |
| Instance Number                     | 0020,0013 | IS |                           | ANAP              | AUTO   |         |
| Contributing Equipment Sequence     | 0018,A001 | SQ |                           | ANAP              | AUTO   |         |
| >Manufacturer                       | 0008,0070 | LO |                           | ALWAYS            | AUTO   |         |
| >Institution Name                   | 0008,0080 | LO |                           | ANAP              | AUTO   |         |
| >Institution Address                | 0008,0081 | ST |                           | ANAP              | AUTO   |         |
| >Station Name                       | 0008,1010 | SH |                           | ANAP              | AUTO   |         |
| >Institutional Department Name      | 0008,1040 | LO |                           | ANAP              | AUTO   |         |
| >Manufacturer's Model Name          | 0008,1090 | LO |                           | ANAP              | AUTO   |         |
| >Device Serial Number               | 0018,1000 | LO |                           | ANAP              | AUTO   |         |
| >Software Version(s)                | 0018,1020 | LO |                           | ANAP              | AUTO   |         |
| >Purpose of Reference Code Sequence | 0040,A170 | SQ |                           | ALWAYS            | AUTO   |         |
| >>Code Value                        | 0008,0100 | SH |                           | ALWAYS            | AUTO   |         |
| >>Coding Scheme Designator          | 0008,0102 | SH |                           | ALWAYS            | AUTO   |         |
| >>Code Meaning                      | 0008,0104 | LO |                           | ALWAYS            | AUTO   |         |

#### 8.1.1.1.4 X-Ray Radiation Dose SR SOP Class

Table 131: SOP Class Modules

| Information Entity | Module                            | Presence    |
|--------------------|-----------------------------------|-------------|
| Patient            | Patient Module                    | ALWAYS      |
| Study              | General Study Module              | ALWAYS      |
|                    | Patient Study Module              | CONDITIONAL |
| Series             | SR Document Series Module         | ALWAYS      |
| Equipment          | General Equipment Module          | ALWAYS      |
|                    | Enhanced General Equipment Module | ALWAYS      |

|       |                            |        |
|-------|----------------------------|--------|
| Image | SR Document General Module | ALWAYS |
|       | SR Document Content Module | ALWAYS |
|       | SOP Common Module          | ALWAYS |

Table 132: Patient Module

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source    | Comment |
|------------------------------|-----------|----|-------|-------------------|-----------|---------|
| Patient's Name               | 0010,0010 | PN |       | VNAP              | COPY      |         |
| Patient ID                   | 0010,0020 | LO |       | VNAP              | COPY      |         |
| Patient's Birth Date         | 0010,0030 | DA |       | VNAP              | COPY      |         |
| Patient's Sex                | 0010,0040 | CS |       | VNAP              | COPY      |         |
| Patient Comments             | 0010,4000 | LT |       | ANAP              | MWL, USER |         |
| Referenced Patient Sequence  | 0008,1120 | SQ |       | ANAP              | MWL, USER |         |
| >Referenced SOP Class UID    | 0008,1150 | UI |       | ALWAYS            | AUTO      |         |
| >Referenced SOP Instance UID | 0008,1155 | UI |       | ALWAYS            | AUTO      |         |

Table 133: General Study Module

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source | Comment                                  |
|------------------------------|-----------|----|-------|-------------------|--------|--|
| Study Date                   | 0008,0020 | DA |       | VNAP              | COPY   |  |
| Study Time                   | 0008,0030 | TM |       | VNAP              | COPY   |  |
| Accession Number             | 0008,0050 | SH |       | VNAP              | COPY   |  |
| Referring Physician's Name   | 0008,0090 | PN |       | VNAP              | COPY   |  |
| Study Description            | 0008,1030 | LO |       | ANAP              | COPY   |  |
| Study Instance UID           | 0020,000D | UI |       | ALWAYS            | COPY   |  |
| Study ID                     | 0020,0010 | SH |       | VNAP              | COPY   |  |
| Referenced Study Sequence    | 0008,1110 | SQ |       | ANAP              | AUTO   | Shall be present with Zero or more items |
| >Referenced SOP Class UID    | 0008,1150 | UI |       | ALWAYS            | AUTO   |  |
| >Referenced SOP Instance UID | 0008,1155 | UI |       | ALWAYS            | AUTO   |  |

Table 134: Patient Study Module

| Attribute Name                  | Tag       | VR | Value | Presence of Value | Source    | Comment |
|---------------------------------|-----------|----|-------|-------------------|-----------|---------|
| Admitting Diagnoses Description | 0008,1080 | LO |       | ANAP              | MWL, USER |         |
| Patient's Age                   | 0010,1010 | AS |       | ANAP              | MWL, USER |         |
| Patient's Size                  | 0010,1020 | DS |       | ANAP              | MWL, USER |         |
| Patient's Weight                | 0010,1030 | DS |       | ANAP              | MWL, USER |         |

Table 135: SR Document Series Module

| Attribute Name                               | Tag       | VR | Value                      | Presence of Value | Source | Comment |
|--|-----------|----|----------------------------|-------------------|--------|---------|
| Series Date                                  | 0008,0021 | DA |                            | ANAP              | COPY   |         |
| Series Time                                  | 0008,0031 | TM |                            | ANAP              | COPY   |         |
| Modality                                     | 0008,0060 | CS | SR                         | ALWAYS            | FIXED  |         |
| Series Description                           | 0008,103E | LO | Radiation Dose Information | ANAP              | FIXED  |         |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ |                            | VNAP              | AUTO   |         |
| Series Instance UID                          | 0020,000E | UI |                            | ALWAYS            | AUTO   |         |
| Series Number                                | 0020,0011 | IS | 5000                       | ALWAYS            | FIXED  |         |

Table 136: General Equipment Module

| Attribute Name                | Tag       | VR | Value    | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|----------|-------------------|--------|---------|
| Manufacturer                  | 0008,0070 | LO | Philips  | VNAP              | FIXED  |         |
| Institution Name              | 0008,0080 | LO |          | ANAP              | COPY   |         |
| Institution Address           | 0008,0081 | ST |          | ANAP              | COPY   |         |
| Station Name                  | 0008,1010 | SH |          | ANAP              | COPY   |         |
| Institutional Department Name | 0008,1040 | LO |          | ANAP              | COPY   |         |
| Manufacturer's Model Name     | 0008,1090 | LO | Big Bore | ANAP              | FIXED  |         |
| Device Serial Number          | 0018,1000 | LO |          | ANAP              | COPY   |         |
| Software Version(s)           | 0018,1020 | LO | 4.8      | ANAP              | FIXED  |         |

Table 137: Enhanced General Equipment Module

| Attribute Name            | Tag       | VR | Value    | Presence of Value | Source | Comment |
|---------------------------|-----------|----|----------|-------------------|--------|---------|
| Manufacturer              | 0008,0070 | LO | Philips  | ALWAYS            | FIXED  |         |
| Manufacturer's Model Name | 0008,1090 | LO | Big Bore | ALWAYS            | FIXED  |         |
| Device Serial Number      | 0018,1000 | LO |          | ALWAYS            | AUTO   |         |
| Software Version(s)       | 0018,1020 | LO | 4.8      | ALWAYS            | FIXED  |         |

Table 138: SR Document General Module

| Attribute Name                      | Tag       | VR | Value | Presence of Value | Source | Comment |
|-------------------------------------|-----------|----|-------|-------------------|--------|---------|
| Content Date                        | 0008,0023 | DA |       | ALWAYS            | AUTO   |         |
| Content Time                        | 0008,0033 | TM |       | ALWAYS            | AUTO   |         |
| Instance Number                     | 0020,0013 | IS | 1     | ALWAYS            | FIXED  |         |
| Referenced Request Sequence         | 0040,A370 | SQ |       | ALWAYS            | AUTO   |         |
| >Reason for the Requested Procedure | 0040,1002 | LO |       | ALWAYS            | AUTO   |         |
| >Reason for Requested               | 0040,100A | SQ |       | ALWAYS            | AUTO   |         |

| Attribute Name                    | Tag       | VR | Value      | Presence of Value | Source | Comment |
|-----------------------------------|-----------|----|------------|-------------------|--------|---------|
| Procedure Code Sequence           |           |    |            |                   |        |         |
| >>Code Value                      | 0008,0100 | SH |            | ALWAYS            | AUTO   |         |
| >>Coding Scheme Designator        | 0008,0102 | SH |            | ALWAYS            | AUTO   |         |
| >>Code Meaning                    | 0008,0104 | LO |            | ALWAYS            | AUTO   |         |
| Performed Procedure Code Sequence | 0040,A372 | SQ |            | VNAP              | AUTO   |         |
| Completion Flag                   | 0040,A491 | CS | COMPLETE   | ALWAYS            | FIXED  |         |
| Verification Flag                 | 0040,A493 | CS | UNVERIFIED | ALWAYS            | FIXED  |         |

Table 139: SR Document Content Module

| Attribute Name             | Tag       | VR | Value     | Presence of Value | Source | Comment |
|----------------------------|-----------|----|-----------|-------------------|--------|---------|
| Value Type                 | 0040,A040 | CS | CONTAINER | ALWAYS            | FIXED  |         |
| Concept Name Code Sequence | 0040,A043 | SQ |           | ALWAYS            | AUTO   |         |
| >Code Value                | 0008,0100 | SH |           | EMPTY             | AUTO   |         |
| >Coding Scheme Designator  | 0008,0102 | SH |           | EMPTY             | AUTO   |         |
| >Code Meaning              | 0008,0104 | LO |           | EMPTY             | AUTO   |         |
| Continuity Of Content      | 0040,A050 | CS |           | ALWAYS            | AUTO   |         |
| Content Template Sequence  | 0040,A504 | SQ |           | ALWAYS            | FIXED  |         |
| >Mapping Resource          | 0008,0105 | CS | DCMR      | ALWAYS            | FIXED  |         |
| >Template Identifier       | 0040,DB00 | CS | TID 10011 | ALWAYS            | FIXED  |         |
| Content Sequence           | 0040,A730 | SQ |           | ANAP              | AUTO   |         |
| Content Template Sequence  | 0040,A504 | SQ |           | ALWAYS            | FIXED  |         |

Table 140: SOP Common Module

| Attribute Name                  | Tag       | VR | Value                         | Presence of Value | Source | Comment |
|---------------------------------|-----------|----|-------------------------------|-------------------|--------|---------|
| Specific Character Set          | 0008,0005 | CS |                               | ANAP              | COPY   |         |
| Instance Creation Date          | 0008,0012 | DA |                               | ANAP              | COPY   |         |
| Instance Creation Time          | 0008,0013 | TM |                               | ANAP              | COPY   |         |
| SOP Class UID                   | 0008,0016 | UI | 1.2.840.10008.5.1.4.1.1.88.67 | ALWAYS            | FIXED  |         |
| SOP Instance UID                | 0008,0018 | UI |                               | ALWAYS            | COPY   |         |
| Contributing Equipment Sequence | 0018,A001 | SQ |                               | ANAP              |        |         |
| >Manufacturer                   | 0008,0070 | LO |                               | ALWAYS            | AUTO   |         |
| >Institution Name               | 0008,0080 | LO |                               | ANAP              | AUTO   |         |
| >Institution Address            | 0008,0081 | ST |                               | ANAP              | AUTO   |         |
| >Station Name                   | 0008,1010 | SH |                               | ANAP              | AUTO   |         |

| Attribute Name                      | Tag       | VR | Value | Presence of Value | Source | Comment |
|-------------------------------------|-----------|----|-------|-------------------|--------|---------|
| >Institutional Department Name      | 0008,1040 | LO |       | ANAP              | AUTO   |         |
| >Manufacturer's Model Name          | 0008,1090 | LO |       | ANAP              | AUTO   |         |
| >Device Serial Number               | 0018,1000 | LO |       | ANAP              | AUTO   |         |
| >Software Version(s)                | 0018,1020 | LO |       | ANAP              | AUTO   |         |
| >Purpose of Reference Code Sequence | 0040,A170 | SQ |       | ALWAYS            | AUTO   |         |
| >>Code Value                        | 0008,0100 | SH |       | ALWAYS            | AUTO   |         |
| >>Coding Scheme Designator          | 0008,0102 | SH |       | ALWAYS            | AUTO   |         |
| >>Code Meaning                      | 0008,0104 | LO |       | ALWAYS            | AUTO   |         |
| Instance Number                     | 0020,0013 | IS |       | ANAP              | AUTO   |         |

### 8.1.2. Usage of Attributes from Received IOD

The following attributes shall be present in the received IODs in order to be accepted.

#### For all IODs

- SOP Class UID (0008,0016)
- Study Instance UID (0020,000D)
- Series Instance UID (0020,000E)

#### For Image IODs

- Pixel Data (7FE0,0010) - Size may not be 0.
- Bits Allocated (0028,0100)

### 8.1.3. Attribute Mapping

Not Applicable. Attribute Mapping is not claimed by the Big Bore system.

### 8.1.4. Coerced/Modified fields

The Import/Export Transparency of DICOM objects means preserving the attributes' values of the objects imported from an external system (remote or removable), optionally processed and then exported to an external system.

The system complies with Level-2 requirements for Storage SCP as defined in DICOM PS 3.4 Appendix B4.1. In other words, all Type 1, Type 2, and Type 3 Attributes defined in the Information Object Definition (IOD) associated with the SOP Class, as well as any Standard Extended attributes (including Private Attributes) included in the SOP Instance, will be stored and may be accessed.

The system does not coerce any Data Elements, except those defined in the DICOM PS 3.4 Appendix B4.1. In other words, when a DICOM object is imported from another system and later exported, all the attributes values will remain unchanged.

In the received IODs, the following attributes may be modified under certain conditions.

**Table 141: Modified Attributes**

| Attribute      | Tag       | When Modified  |
|----------------|-----------|--|
| Patient's Name | 0010,0010 | If Empty, the Patient's Name will be set to "Unknown". |
| Patient ID     | 0010,0020 | If Empty, the Patient ID will be set to "Unknown".     |

| Attribute        | Tag       | When Modified  |
|------------------|-----------|--|
| Rows             | 0028,0010 | Is Fixed, if rows columns does not match pixel data size.        |
| SOP Instance UID | 0008,0018 | If missing, a new SOP Instance UID will be generate by Big Bore. |

## 8.2. Data Dictionary of Private Attributes

Not Applicable. Big System has not implemented Data Dictionary of Private Attributes.

## 8.3. Coded Terminology and Templates

These are defined in Section 8.5.1 in the specific objects.

### 8.3.1. Context Groups

Not Applicable. Big System has not implemented Context Groups.

### 8.3.2. Template Specifications

#### 8.3.2.1. CT RADIATION DOSE SR IOD TEMPLATES




Figure 21: CT Radiation Dose SR IOD Template Structure

This section describes the content of all the templates used in the CT Radiation Dose Reporting SR.

Table 142: Used Templates for CT Radiation Dose Reporting

| Template Name                          | Template ID |
|--|-------------|
| CT Radiation Dose                      | TID 10011   |
| CT Accumulated Dose Data               | TID 10012   |
| CT Irradiation Event Data              | TID 10013   |
| Scanning Length                        | TID 10014   |
| CT Dose Check Details                  | TID 10015   |
| Observer Context                       | TID 1002    |
| Device Observer Identifying Attributes | TID 1004    |

### 8.3.2.1.1 TID 10011 CT Radiation Dose

Table 143: CT Radiation Dose

| NL | Relation with Parent | Concept Name                                   | VT        | VM  | Presence of Value | Value                                    |
|----|----------------------|--|-----------|-----|-------------------|--|
|    |                      | EV (113701, DCM, X-Ray Radiation Dose Report)  | CONTAINER | 1   | ALWAYS            |  |
| >  | HAS CONCEPT MOD      | EV (121058, DCM, Procedure reported)           | CODE      | 1   | ALWAYS            | P5-08000, SRT, Computed Tomography X-Ray |
| >> | HAS CONCEPT MOD      | EV (G-C0E8, SRT, Has Intent)                   | CODE      | 1   | ALWAYS            | R-408C3, SRT, Diagnostic Intent          |
| >  |                      | DTID (1002) Observer Context                   | INCLUDE   | 1   | ALWAYS            |  |
| >  | HAS OBS CONTEXT      | EV (113809, DCM, Start of X-Ray Irradiation    | DATE/TIME | 1   | ALWAYS            |  |
| >  | HAS OBS CONTEXT      | EV (113810, DCM, End of X-Ray Irradiation)     | DATE/TIME | 1   | ALWAYS            |  |
| >  | HAS OBS CONTEXT      | EV (113705, DCM, Scope of Accumulation)        | CODE      | 1   | ALWAYS            | 113014, DCM, "Study"                     |
| >> | HAS PROPERTIES       | DCID (10001) UID Types                         | UIDREF    | 1   | ALWAYS            | Study Instance UID (0020,000D)           |
| >  | CONTAINS             | DTID (10012) CT Accumulated Dose Data          | INCLUDE   | 1   | ALWAYS            |  |
| >  | CONTAINS             | DTID (10013) CT Irradiation Event Data         | INCLUDE   | 1-n | ALWAYS            |  |
| >  | CONTAINS             | EV (113854, DCM, "Source of Dose Information") | CODE      | 1   | ALWAYS            | 113856, DCM, "Automated Data Collection" |

### 8.3.2.1.2 TID 10012 CT Accumulated Dose

Table 144: CT Accumulated Dose

| NL | Relation with Parent | Concept Name  | VT        | VM | Presence of Value | Value |
|----|----------------------|---|-----------|----|-------------------|-------|
|    |                      | EV (113811, DCM, "CT Accumulated Dose Data")          | CONTAINER | 1  | ALWAYS            |       |
| >  | CONTAINS             | EV (113812, DCM, Total Number of Irradiation Events)] | NUM       | 1  | ALWAYS            |       |
| >  | CONTAINS             | EV (113813, DCM, CT Dose Length Product Total)        | NUM       | 1  | ALWAYS            |       |

### 8.3.2.1.3 TID 10013 CT Irradiation Event Data

Table 145: CT Irradiation Event Data

| NL | Relation with Parent | Concept Name   | VT        | VM | Presence of Value | Value  |
|----|----------------------|--|-----------|----|-------------------|--|
|    |                      | EV (113819, DCM, "CT Acquisition")                   | CONTAINER | 1  | ALWAYS            |  |
| >  | CONTAINS             | EV (125203, DCM, "Acquisition Protocol")             | TEXT      | 1  | ALWAYS            | Protocol Name (0018,1030)  |
| >  | CONTAINS             | EV (123014, DCM, "Target Region")                    | CODE      | 1  | ALWAYS            | Scan Type / Organ Type   |
| >  | CONTAINS             | EV (113820, DCM, "CT Acquisition Type")              | CODE      | 1  | ALWAYS            | Acquisition Type: Constant Angle Acquisition   |
| >  | CONTAINS             | EV (G-C32C, SRT, "Procedure Context")                | CODE      | 1  | ALWAYS            | "Diagnostic radiography with contrast media " or "CT without contrast"   |
| >  | CONTAINS             | EV (113769, DCM, "Irradiation Event UID")            | UIDREF    | 1  | ALWAYS            | Identical to Irradiation Event UID in the images (0008,3010)   |
| >  | CONTAINS             | EV (113822, DCM, "CT Acquisition Parameters")        | CONTAINER | 1  | ALWAYS            |  |
| >> | CONTAINS             | EV (113824, DCM, "Exposure Time")                    | NUM       | 1  | ALWAYS            | Value = Exposure Time (0018, 1150) - per scan. (units = s)   |
| >> | CONTAINS             | DTID 10014 : Scanning Length                         | INCLUDE   | 1  | ALWAYS            | Value = Scan Length (0018,1302) - per scan (units = mm)  |
| >> | CONTAINS             | EV (113826, DCM, "Nominal Single Collimation Width") | NUM       | 1  | ALWAYS            | The width of a single row of acquired data<br>Value = Single Collimation Width(0018,9306) - per scan (units = mm)                                |
| >> | CONTAINS             | EV (113827, DCM, "Nominal Total Collimation Width")  | NUM       | 1  | ALWAYS            | The width of the total collimation over the area of active x-ray detection<br>Value = Total Collimation Width(0018,9307) - per scan (units = mm) |
| >> | CONTAINS             | EV (113828, DCM, "Pitch Factor")                     | NUM       | 1  | CONDITIONAL       | UNITS = EV ({ratio}, UCUM, "ratio")  |
| >> | CONTAINS             | EV (113823, DCM, "Number of X-Ray Sources")          | NUM       | 1  | ALWAYS            | 1  |

|          |                 |  |           |     |             |  |
|----------|-----------------|--|-----------|-----|-------------|--|
| >>       | CONTAINS        | EV (113831, DCM, "CT X-Ray Source Parameters")                                       | CONTAINER | 1   | ALWAYS      |  |
| >>>      | CONTAINS        | EV (113832, DCM, "Identification of the X-Ray Source")                               | TEXT      | 1   | ALWAYS      | A  |
| >>>      | CONTAINS        | EV (113733, DCM, "KVP")  | NUM       | 1   | ALWAYS      | Same as KVP (0018,0060) - per scan<br>Units = kV   |
| >>>      | CONTAINS        | EV (113833, DCM, "Maximum X-Ray Tube Current")                                       | NUM       | 1   | ALWAYS      | Value = Max X-ray Tube Current (01E1, 1052) - per scan.<br>Units = mA  |
| >>>      | CONTAINS        | EV (113734, DCM, "X-Ray Tube Current")   | NUM       | 1   | ALWAYS      | Value = X-ray Tube Current (0018,1151) - per scan<br>Units = mA  |
| >>>      | CONTAINS        | EV (113834, DCM, "Exposure Time per Rotation")                                       | NUM       | 1   | CONDITIONAL | UNITS = EV (s, UCUM, "s")  |
| >        | CONTAINS        | EV (113829, DCM, "CT Dose")  | CONTAINER | 1   | CONDITIONAL | IF EV (113820, DCM, "CT Acquisition Type") does not equal (113805, DCM, "Constant Angle Acquisition")                      |
| >>       | CONTAINS        | EV (113830, DCM, "Mean CTDIvol")   | NUM       | 1   | ALWAYS      | Value = CTDIVol (0018,9345) - per scan<br>Units = mGy  |
| >>       | CONTAINS        | EV (113835, DCM, "CTDIw Phantom Type")   | CODE      | 1   | ALWAYS      | Value = "IEC Head Dosimetry Phantom " for 16cm Head phantom<br>Value = "IEC Body Dosimetry Phantom " for 32cm Body phantom |
| >>       | CONTAINS        | EV (113838, DCM, "DLP")  | NUM       | 1   | ALWAYS      | Value = DLP (00E1, 1021) - per scan (+/- 2% is acceptable).<br>Units = mGy*cm  |
| >>       | CONTAINS        | EV (113930, DCM, "Size Specific Dose Estimation")                                    | NUM       | 1-n | CONDITIONAL | UNITS = EV (mGy, UCUM, "mGy")  |
| >>>      | HAS CONCEPT MOD | EV (G-C036, SRT, "Measurement Method")   | CODE      | 1   | ALWAYS      |  |
| >>><br>> | INFERRED FROM   | EV (113985, DCM, "Series or Instance used for Water Equivalent Diameter estimation") | UIDREF    | 1-n | CONDITIONAL | IF (113984, DCM, "Water Equivalent Diameter From Localizer")   |

|    |          |                                    |         |   |        |  |
|----|----------|------------------------------------|---------|---|--------|--|
| >> | CONTAINS | DTID (10015) CT Dose Check Details | INCLUDE | 1 | ALWAYS |  |
|----|----------|------------------------------------|---------|---|--------|--|

#### 8.3.2.1.4 TID 1002 Observer Context

Table 146: Observer Context

| NL | Relation with Parent | Concept Name                                       | VT      | VM | Presence of Value | Value                 |
|----|----------------------|--|---------|----|-------------------|-----------------------|
| >  | HAS OBS CONTEXT      | EV (121005, DCM, Observer Type)                    | CODE    | 1  | CONDITIONAL       | (121007, DCM, Device) |
| >  | HAS OBS CONTEXT      | DTID (1004) Device observer identifying attributes | INCLUDE | 1  | ALWAYS            |                       |

#### 8.3.2.1.5 TID 1004 Device Observer Identifying Attributes

Table 147: Device Observer Identifying Attributes

| NL | Relation with Parent | Concept Name  | VT     | VM | Presence of Value | Value  |
|----|----------------------|---|--------|----|-------------------|--|
| >  | HAS OBS CONTEXT      | EV (121012, DCM, Device Observer UID)                                 | UIDREF | 1  | ALWAYS            | Concatenation of Philips CT Root (1.3.46.670589.33.1) and Network Board MAC address      |
| >  | HAS OBS CONTEXT      | EV (121013, DCM, Device Observer Name)                                | TEXT   | 1  | CONDITIONAL       | Station Name (0008,1010): Equals the computer name (under Start->System Settings)        |
| >  | HAS OBS CONTEXT      | EV (121014, DCM, Device Observer Manufacturer)                        | TEXT   | 1  | CONDITIONAL       | Manufacturer (0008,0070): Philips  |
| >  | HAS OBS CONTEXT      | EV (121015, DCM, Device Observer Model Name)                          | TEXT   | 1  | CONDITIONAL       | Manufacturer's Model Name (0008,1090): Big Bore  |
| >  | HAS OBS CONTEXT      | EV (121016, DCM, Device Observer Serial Number)                       | TEXT   | 1  | CONDITIONAL       | Device Serial Number (0018,1000) As in Preferences -> Institute -> Product Serial Number |
| >  | HAS OBS CONTEXT      | EV(121017, DCM, Device Observer Physical Location During Observation) | TEXT   | 1  | CONDITIONAL       | Institution Name (0008,0080): As in Preferences -> Institute -> Name                     |

#### 8.3.2.1.6 TID 10014 Scanning Length

Table 148: Scanning Length

| NL | Relation with Parent | Concept Name                        | VT  | VM | Presence of Value | Value          |
|----|----------------------|-------------------------------------|-----|----|-------------------|----------------|
|    | CONTAINS             | EV (113825, DCM, "Scanning Length") | NUM | 1  | ALWAYS            | (mm, UCUM, mm) |

### 8.3.2.1.7 TID 10015 CT Dose Check Details

Table 149: CT Dose Check Details

| NL | Relation with Parent | Concept Name  | VT        | VM | Presence of Value | Value  |
|----|----------------------|---|-----------|----|-------------------|--|
|    |                      | EV (113900, DCM, "Dose Check Alert Details")              | CONTAINER | 1  | ALWAYS            |  |
| >  | CONTAINS             | EV (113901, DCM, DLP Alert Value Configured)              | CODE      | 1  | ALWAYS            | Yes/No<br>Value depends on whether the relevant (head/body) DLP Alert value is set in preferences.     |
| >  | CONTAINS             | EV (113902, DCM, CTDIvol Alert Value Configured)          | CODE      | 1  | ALWAYS            | Yes/No<br>Value depends on whether the relevant (head/body) CTDIvol Alert value is set in preferences. |
| >  | CONTAINS             | EV (113904, DCM, CTDIvol Alert Value)                     | NUM       | 1  | CONDITIONAL       | (mGy, UCUM, mGy)   |
|    |                      | EV (113908, DCM, Dose Check Notification Details)         | CONTAINER | 1  | CONDITIONAL       |  |
| >  | CONTAINS             | EV (113909, DCM, "DLP Notification Value Configured")     | CODE      | 1  | ALWAYS            | Yes/No<br>Value depends on whether a DLP Notification value was set in the protocol step.              |
| >  | CONTAINS             | EV (113910, DCM, "CTDIvol Notification Value Configured") | CODE      | 1  | ALWAYS            | Value depends on whether a CTDIvol Notification value was set in the protocol step.                    |
| >  | CONTAINS             | EV (113912, DCM, "CTDIvol Notification Value")            | NUM       | 1  | CONDITIONAL       |  |

### 8.3.3. Private code definitions

Not Applicable. Big System has not implemented Private code definitions.

### 8.4. Grayscale Image consistency

Not Applicable. Big System has not implemented Grayscale image consistency.

## 8.5. Standard Extended/Specialized/Private SOPs

Table 150: List of created SOP Classes

| SOP Class Name                            | SOP Class UID                 |
|---|-------------------------------|
| CT Image Storage SOP Class                | 1.2.840.10008.5.1.4.1.1.2     |
| Secondary Capture Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7     |
| General ECG Waveform Storage SOP Class    | 1.2.840.10008.5.1.4.1.1.9.1.2 |
| X-Ray Radiation Dose SR                   | 1.2.840.10008.5.1.4.1.1.88.67 |

### 8.5.1. Standard Extended/Specialized/Private SOP Instance

This section describes (particular) Standard Extended SOP Classes, Specialized SOP Classes, or Private SOP Classes that are used.

#### 8.5.1.1. CT Image Storage SOP Class

Table 151: Extended DICOM and private attributes for CT Image Storage SOP Class Instances

| Attribute Name              | Tag       | VR | Value | Presence of Value | Source | Comment |
|-----------------------------|-----------|----|-------|-------------------|--------|---------|
| Angular Position (retired)  | 0018,1141 | DS |       | ALWAYS            | AUTO   |         |
| Acquisition Type            | 0018,9302 | CS |       | ALWAYS            | AUTO   |         |
| Slice Progression Direction | 0054,0500 | CS |       | ALWAYS            | AUTO   |         |
| Units                       | 0054,1001 | CS |       | ALWAYS            | AUTO   |         |

#### 8.5.1.2. Secondary Capture Image Storage SOP Class

Table 152: Extended DICOM and private attributes for Secondary Capture Image Storage SOP Class Instances

| Attribute Name              | Tag       | VR | Value | Presence of Value | Source | Comment |
|-----------------------------|-----------|----|-------|-------------------|--------|---------|
| Medical Alerts              | 0010,2000 | LO |       | VNAP              | COPY   |         |
| Allergies                   | 0010,2110 | LO |       | VNAP              | COPY   |         |
| Pregnancy Status            | 0010,21C0 | US |       | VNAP              | COPY   |         |
| Scan Options                | 0018,0022 | CS |       | ANAP              | COPY   |         |
| Slice Thickness             | 0018,0050 | DS |       | ANAP              | COPY   |         |
| KVP                         | 0018,0060 | DS |       | ANAP              | COPY   |         |
| Data Collection Diameter    | 0018,0090 | DS |       | ANAP              | COPY   |         |
| Reconstruction Diameter     | 0018,1100 | DS |       | ANAP              | COPY   |         |
| Distance Source to Detector | 0018,1110 | DS |       | ANAP              | COPY   |         |
| Distance Source to Patient  | 0018,1111 | DS |       | ANAP              | COPY   |         |
| Gantry/Detector Tilt        | 0018,1120 | DS |       | ANAP              | COPY   |         |
| Table Height                | 0018,1130 | DS |       | ANAP              | COPY   |         |
| Angular Position (retired)  | 0018,1141 | DS |       | ANAP              | COPY   |         |

| Attribute Name              | Tag       | VR | Value | Presence of Value | Source | Comment |
|-----------------------------|-----------|----|-------|-------------------|--------|---------|
| Exposure Time               | 0018,1150 | IS |       | ANAP              | COPY   |         |
| X-ray Tube Current          | 0018,1151 | IS |       | ANAP              | COPY   |         |
| Acquisition Type            | 0018,9302 | CS |       | ANAP              | COPY   |         |
| Single Collimation Width    | 0018,9306 | FD |       | ANAP              | COPY   |         |
| Total Collimation Width     | 0018,9307 | FD |       | ANAP              | COPY   |         |
| Table Speed                 | 0018,9309 | FD |       | ANAP              | COPY   |         |
| Exposure Modulation Type    | 0018,9323 | CS |       | ANAP              | COPY   |         |
| Estimated Dose Saving       | 0018,9324 | FD |       | ANAP              | COPY   |         |
| Requesting Physician        | 0032,1032 | PN |       | VNAP              | COPY   |         |
| Special Needs               | 0038,0050 | LO |       | VNAP              | COPY   |         |
| Patient State               | 0038,0500 | LO |       | VNAP              | COPY   |         |
| Total Number of Exposures   | 0040,0301 | US |       | ANAP              | COPY   |         |
| Exposure Dose Sequence      | 0040,030E | SQ |       | ANAP              | COPY   |         |
| >Acquisition Date time      | 0008,002A | DT |       | ANAP              | COPY   |         |
| >Series Description         | 0008,103E | LO |       | ANAP              | COPY   |         |
| >Contrast/Bolus Agent       | 0018,0010 | LO |       | ANAP              | COPY   |         |
| >Body Part Examined         | 0018,0015 | CS |       | ANAP              | COPY   |         |
| >KVP                        | 0018,0060 | DS |       | ANAP              | COPY   |         |
| >Protocol Name              | 0018,1030 | LO |       | ANAP              | COPY   |         |
| >Exposure Time              | 0018,1150 | IS |       | ANAP              | COPY   |         |
| >X-ray Tube Current         | 0018,1151 | IS |       | ANAP              | COPY   |         |
| >Exposure                   | 0018,1152 | IS |       | ANAP              | COPY   |         |
| >Radiation Mode             | 0018,115A | CS |       | ANAP              | COPY   |         |
| >Scan Length                | 0018,1302 | IS |       | ANAP              | COPY   |         |
| >Acquisition Duration       | 0018,9073 | FD |       | ANAP              | COPY   |         |
| >Acquisition Type           | 0018,9302 | CS |       | ANAP              | COPY   |         |
| >Single Collimation Width   | 0018,9306 | FD |       | ANAP              | COPY   |         |
| >Total Collimation Width    | 0018,9307 | FD |       | ANAP              | COPY   |         |
| >CTDIvol                    | 0018,9345 | FD |       | ANAP              | COPY   |         |
| >Series Number              | 0020,0011 | IS |       | ANAP              | COPY   |         |
| >Comments on Radiation Dose | 0040,0310 | ST |       | ANAP              | COPY   |         |
| Comments on Radiation Dose  | 0040,0310 | ST |       | ANAP              | COPY   |         |

### 8.5.1.3. General ECG Waveform Storage SOP Class

Table 153: Extended DICOM and private attributes for General ECG Waveform Storage SOP Class Instances

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|--------|---------|
| Frame of Reference UID       | 0020,0052 | UI |       | ANAP              | AUTO   |         |
| Position Reference Indicator | 0020,1040 | LO |       | ALWAYS            | AUTO   |         |
| Image Comments               | 0020,4000 | LT |       | ALWAYS            | AUTO   |         |

### 8.5.1.4. X-Ray Radiation Dose SR SOP Class

Extended DICOM and Private Attributes are not applicable for this SOP class instance.

## 8.6. Private Transfer Syntaxes

Supported Private Transfer Syntaxes are shown in the next Table.

Table 154: Supported Private Transfer Syntaxes

| Transfer Syntax Name                                   | Transfer Syntax UID    | Comment              |
|--|------------------------|----------------------|
| Private CT Transfer Syntax - Explicit VR Little Endian | 1.3.46.670589.33.1.4.1 | Private ELE (P-ELE). |

*- This part of the page is left intentionally empty -*

**Issued by:**

Philips Medical Systems Nederland B.V., a Philips Healthcare company,  
P.O. Box 10.000  
5680 DA Best  
The Netherlands

Internet: <https://www.philips.com/dicom>

Doc Id: HSDP-1268145

Date: 18-DEC-2024

