

# DICOM Conformance Statement

## Azurion 2.0



**Issued by:**

Philips Medical Systems Netherlands BV, a Philips Healthcare company,

P.O. Box 10.000  
5680 DA Best  
The Netherlands

Email: [dicom@philips.com](mailto:dicom@philips.com)  
Internet: [www.philips.com/dicom](http://www.philips.com/dicom)

Document Number: ICAP-PF.0035781  
Date: 09-Jun-2018



## 1. DICOM Conformance Statement Overview

Azurion 2.0 system is an interventional X-ray system that is used to perform:

- Image guidance during diagnostic, interventional and minimally invasive surgery procedures for the following clinical application areas: vascular, non-vascular, cardiovascular and neuro procedures.
- Cardiac imaging applications including diagnostics, interventional and minimally invasive surgery procedures.

It provides the following DICOM data exchange features:

- Query the MWL SCP for a Modality Worklist (MWL)
- Update the Modality Performed Procedure Step (MPPS).
- Transfer of DICOM Images and Grayscale Presentation States.
- Query/Retrieve a Workstation or PACS for a list of entries representing Series information of DICOM Images
- Send Storage Commitment to the PACS (for the safe-keeping of the previously transmitted images) and handling the Storage Commitment notifications received from the PACS.
- Print Images on DICOM Printers
- Transfer of X-Ray Radiation Dose Structured Reports to the PACS or Workstation.




Figure 1 : Data Flow of Azurion 2.0 in a DICOM Network

Table 1: Network Services

| SOP Class                                       |                        | User of Service (SCU) | Provider of Service (SCP) |
|---|------------------------|-----------------------|---------------------------|
| Name  | UID                    |                       |                           |
| <b>Other</b>                                    |                        |                       |                           |
| Verification SOP Class                          | 1.2.840.10008.1.1      | Yes                   | Yes                       |
| <b>Print Management</b>                         |                        |                       |                           |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9  | Yes                   | No                        |
| >Basic Film Session SOP Class                   | 1.2.840.10008.5.1.1.1  | Yes                   | No                        |
| >Basic Film Box SOP Class                       | 1.2.840.10008.5.1.1.2  | Yes                   | No                        |
| >Basic Grayscale Image Box SOP Class            | 1.2.840.10008.5.1.1.4  | Yes                   | No                        |
| >Printer SOP Class                              | 1.2.840.10008.5.1.1.16 | Yes                   | No                        |
| Basic Color Print Management Meta SOP Class     | 1.2.840.10008.5.1.1.18 | Yes                   | No                        |
| >Basic Film Session SOP Class                   | 1.2.840.10008.5.1.1.1  | Yes                   | No                        |

Template ID: HSDPCPF-T-30001.09



| SOP Class   |                               | User of Service (SCU) | Provider of Service (SCP) |
|---|-------------------------------|-----------------------|---------------------------|
| Name  | UID                           |                       |                           |
| >Basic Film Box SOP Class                               | 1.2.840.10008.5.1.1.2         | Yes                   | No                        |
| >Basic Color Image Box SOP Class                        | 1.2.840.10008.5.1.1.4.1       | Yes                   | No                        |
| >Printer SOP Class                                      | 1.2.840.10008.5.1.1.16        | Yes                   | No                        |
| Query/Retrieve  |                               |                       |                           |
| Patient Root QR Information Model - FIND SOP Class      | 1.2.840.10008.5.1.4.1.2.1.1   | Yes                   | No                        |
| Study Root QR Information Model - FIND SOP Class        | 1.2.840.10008.5.1.4.1.2.2.1   | Yes                   | No                        |
| Patient Root QR Information Model - MOVE SOP Class      | 1.2.840.10008.5.1.4.1.2.1.2   | Yes                   | No                        |
| Study Root QR Information Model - MOVE SOP Class        | 1.2.840.10008.5.1.4.1.2.2.2   | Yes                   | No                        |
| Transfer  |                               |                       |                           |
| Secondary Capture Image Storage SOP Class               | 1.2.840.10008.5.1.4.1.1.7     | Yes                   | Yes                       |
| Grayscale Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.1  | Yes                   | Yes                       |
| X-Ray Angiographic Image Storage SOP Class              | 1.2.840.10008.5.1.4.1.1.12.1  | Yes                   | Yes                       |
| X-Ray Radiation Dose Structured Report SOP Class        | 1.2.840.10008.5.1.4.1.1.88.67 | Yes                   | No                        |
| Workflow Management                                     |                               |                       |                           |
| Storage Commitment Push Model SOP Class                 | 1.2.840.10008.1.20.1          | Yes                   | No                        |
| Modality Performed Procedure Step SOP Class             | 1.2.840.10008.3.1.2.3.3       | Yes                   | No                        |
| Modality Worklist Information Model - FIND SOP Class    | 1.2.840.10008.5.1.4.31        | Yes                   | No                        |

Table 2: Media Services

| Media Storage Application Profile            | File-set Creator (FSC) | File-set Updater (FSU) | File-set Reader (FSR) |
|--|------------------------|------------------------|-----------------------|
| Compact disk-Recordable                      |                        |                        |                       |
| General Purpose CD-R Interchange             | Yes                    | No                     | Yes                   |
| DVD-RAM                                      |                        |                        |                       |
| General Purpose Interchange on DVD-RAM Media | Yes                    | No                     | Yes                   |
| USB  |                        |                        |                       |
| General Purpose USB Media                    | Yes                    | Yes                    | Yes                   |

## 2. Table of Contents

|  |           |
|--|-----------|
| <b>1. DICOM CONFORMANCE STATEMENT OVERVIEW</b>                             | <b>3</b>  |
| <b>2. TABLE OF CONTENTS</b>  | <b>5</b>  |
| <b>3. INTRODUCTION</b>   | <b>7</b>  |
| 3.1. REVISION HISTORY  | 7         |
| 3.2. AUDIENCE  | 7         |
| 3.3. REMARKS   | 7         |
| 3.4. DEFINITIONS, TERMS AND ABBREVIATIONS                                  | 8         |
| 3.5. REFERENCES  | 9         |
| <b>4. NETWORKING</b>   | <b>10</b> |
| <b>4.1. IMPLEMENTATION MODEL</b>   | <b>10</b> |
| 4.1.1. Application Data Flow   | 10        |
| 4.1.2. Functional Definition of AE's                                       | 11        |
| 4.1.2.1. Functional Definition of Application Entity Azurion 2.0           | 11        |
| 4.1.3. Sequencing of Real World Activities                                 | 12        |
| <b>4.2. AE SPECIFICATIONS</b>  | <b>13</b> |
| 4.2.1. AE Specification of Azurion 2.0                                     | 13        |
| 4.2.1.1. SOP Classes   | 13        |
| 4.2.1.2. Association Policies  | 14        |
| 4.2.1.3. General   | 14        |
| 4.2.1.4. Number of Associations  | 14        |
| 4.2.1.5. Asynchronous Nature   | 14        |
| 4.2.1.6. Implementation Identifying Information                            | 14        |
| 4.2.1.7. Communication Failure Handling                                    | 15        |
| 4.2.1.8. Association Initiation Policy                                     | 15        |
| 4.2.1.9. (Real-World) Activity – Verification as SCU                       | 16        |
| 4.2.1.10. (Real-World) Activity – Modality Worklist as SCU                 | 18        |
| 4.2.1.11. (Real-World) Activity – Modality Performed Procedure Step as SCU | 23        |
| 4.2.1.12. (Real-World) Activity – Image Export                             | 33        |
| 4.2.1.13. (Real-World) Activity – Storage Commitment Push Model as SCU     | 35        |
| 4.2.1.14. (Real-World) Activity – FIND as SCU                              | 38        |
| 4.2.1.15. (Real-World) Activity – MOVE as SCU                              | 44        |
| 4.2.1.16. (Real-World) Activity – Print Management as SCU                  | 48        |
| 4.2.1.17. Association Acceptance Policy                                    | 58        |
| 4.2.1.18. (Real-World) Activity – Verification as SCP                      | 59        |
| 4.2.1.19. (Real-World) Activity – Image Import                             | 60        |
| <b>4.3. NETWORK INTERFACES</b>   | <b>62</b> |
| 4.3.1. Physical Network Interfaces   | 62        |
| 4.3.2. Additional Protocols  | 62        |
| 4.3.2.1. Basic TLS Secure Transport Connection Profile                     | 62        |
| 4.3.2.2. Basic Time Synchronization Profile                                | 62        |
| 4.3.2.3. Basic Application Level Confidentiality Profile                   | 62        |
| 4.3.3. IPv4 and IPv6 Support   | 62        |
| <b>4.4. CONFIGURATION</b>  | <b>62</b> |
| 4.4.1. AE Title/Presentation Address Mapping                               | 62        |
| 4.4.1.1. Local AE Titles   | 62        |
| 4.4.1.2. Remote AE Title/Presentation Address Mapping                      | 63        |
| 4.4.2. Parameters  | 63        |
| <b>5. MEDIA INTERCHANGE</b>  | <b>64</b> |
| <b>5.1. IMPLEMENTATION MODEL</b>   | <b>64</b> |
| 5.1.1. Application Data Flow Diagram                                       | 64        |

|  |            |
|--|------------|
| 5.1.2. Functional Definitions of AE's.....   | 64         |
| 5.1.3. Sequencing of Real World Activities .....                                       | 64         |
| <b>5.2. AE SPECIFICATIONS.....</b>   | <b>64</b>  |
| 5.2.1. Media Storage Media - Specification.....  | 64         |
| 5.2.1.1. File Meta Information for the Media AE .....                                  | 65         |
| 5.2.1.2. Real-World Activities .....   | 65         |
| 5.2.1.3. RWA - Read File-set.....  | 65         |
| 5.2.1.4. RWA - Create File-set.....  | 65         |
| 5.2.1.5. RWA - Update File-set.....  | 65         |
| <b>5.3. AUGMENTED AND PRIVATE APPLICATION PROFILES .....</b>                           | <b>66</b>  |
| <b>5.4. MEDIA CONFIGURATION.....</b>   | <b>66</b>  |
| <b>6. SUPPORT OF CHARACTER SETS .....</b>  | <b>67</b>  |
| <b>7. SECURITY .....</b>   | <b>68</b>  |
| <b>7.1. SECURITY PROFILES .....</b>  | <b>68</b>  |
| 7.1.1. Security use Profiles.....  | 68         |
| 7.1.2. Security Transport Connection Profiles .....                                    | 68         |
| 7.1.3. Digital Signature Profiles .....  | 69         |
| 7.1.4. Media Storage Security Profiles .....   | 69         |
| 7.1.5. Attribute Confidentiality Profiles .....  | 69         |
| 7.1.6. Network Address Management Profiles .....                                       | 70         |
| 7.1.7. Time Synchronization Profiles .....   | 70         |
| 7.1.8. Application Configuration Management Profiles .....                             | 70         |
| 7.1.9. Audit Trail Profiles .....  | 70         |
| <b>7.2. ASSOCIATION LEVEL SECURITY .....</b>   | <b>71</b>  |
| <b>7.3. APPLICATION LEVEL SECURITY.....</b>  | <b>71</b>  |
| <b>8. ANNEXES .....</b>  | <b>72</b>  |
| <b>8.1. IOD CONTENTS .....</b>   | <b>72</b>  |
| 8.1.1. Created SOP Instances.....  | 72         |
| 8.1.1.1. List of Created SOP Classes .....   | 72         |
| 8.1.1.2. Secondary Capture Image Storage SOP Class.....                                | 72         |
| 8.1.1.3. Grayscale Softcopy Presentation State Storage SOP Class .....                 | 77         |
| 8.1.1.4. X-Ray Angiographic Image Storage SOP Class.....                               | 82         |
| 8.1.1.5. X-Ray Radiation Dose SR .....   | 92         |
| 8.1.2. Attribute Mapping .....   | 97         |
| 8.1.3. Coerced/Modified Fields.....  | 97         |
| <b>8.2. DATA DICTIONARY OF PRIVATE ATTRIBUTES .....</b>                                | <b>97</b>  |
| <b>8.3. CODED TERMINOLOGY AND TEMPLATES.....</b>                                       | <b>97</b>  |
| 8.3.1. Context Group .....   | 97         |
| 8.3.2. Template Specifications .....   | 97         |
| 8.3.2.1. TID 10001 Projection X-Ray Radiation Dose.....                                | 98         |
| 8.3.2.2. TID 10002 Accumulated X-Ray Dose .....  | 99         |
| 8.3.2.3. TID 10003 Irradiation Event X-Ray Data.....                                   | 99         |
| 8.3.2.4. TID 10003B Irradiation Event X-Ray Source Data.....                           | 101        |
| 8.3.2.5. TID 10003C Irradiation Event X-Ray Mechanical Data.....                       | 102        |
| 8.3.2.6. TID 10004 Accumulated Fluoroscopy and Acquisition Projection X-Ray Dose ..... | 102        |
| 8.3.2.7. TID 10007 Accumulated Total Projection Radiography Dose.....                  | 102        |
| 8.3.2.8. TID 1002 Observer Context .....   | 103        |
| 8.3.2.9. TID 1004 Device Observer Identifying Attributes .....                         | 103        |
| 8.3.3. Private Code Definitions.....   | 103        |
| <b>8.4. GRAYSCALE IMAGE CONSISTENCY.....</b>   | <b>103</b> |
| <b>8.5. STANDARD/EXTENDED/SPECIALIZED/PRIVATE SOP CLASSES.....</b>                     | <b>103</b> |

## 3. Introduction

### 3.1. Revision History

The revision history provides dates and differences of the different releases.

Table 3: Revision History

| document Version | Date of Issue | Status     | Description                     |
|------------------|---------------|------------|---------------------------------|
| 00               | 09-May-2018   | Authorized | Initial version for Azurion 2.0 |

### 3.2. Audience

This Conformance Statement is intended for:

- (Potential) customers
- System integrators of medical equipment
- Marketing staff interested in system functionality
- Software designers implementing DICOM interfaces

It is assumed that the reader is familiar with the DICOM standard.

### 3.3. Remarks

The DICOM Conformance Statement is contained in chapter 4 through 8 and follows the contents and structuring requirements of DICOM PS 3.2.

This DICOM Conformance Statement by itself does not guarantee successful interoperability of Philips equipment with non-Philips equipment. The user (or user's agent) should be aware of the following issues:

- **Interoperability**  
Interoperability refers to the ability of application functions, distributed over two or more systems, to work successfully together. The integration of medical devices into an IT environment may require application functions that are not specified within the scope of DICOM. Consequently, using only the information provided by this Conformance Statement does not guarantee interoperability of Philips equipment with non-Philips equipment.  
It is the user's responsibility to analyse thoroughly the application requirements and to specify a solution that integrates Philips equipment with non-Philips equipment.
- **Validation**  
Philips equipment has been carefully tested to assure that the actual implementation of the DICOM interface corresponds with this Conformance Statement.  
Where Philips equipment is linked to non-Philips equipment, the first step is to compare the relevant Conformance Statements. If the Conformance Statements indicate that successful information exchange should be possible, additional validation tests will be necessary to ensure the functionality, performance, accuracy and stability of image and image related data. It is the responsibility of the user (or user's agent) to specify the appropriate test suite and to carry out the additional validation tests.
- **New versions of the DICOM Standard**  
The DICOM Standard will evolve in future to meet the user's growing requirements and to incorporate new features and technologies. Philips is actively involved in this evolution and plans to adapt its equipment to future versions of the DICOM Standard. In order to do so, Philips reserves the right to make changes to its products or to discontinue its delivery. The user should ensure that any non-Philips provider linking to Philips equipment also adapts to future versions of the DICOM Standard. If not, the incorporation of DICOM enhancements into Philips equipment may lead to loss of connectivity (in case of networking) and incompatibility (in case of media).

### 3.4. Definitions, Terms and Abbreviations

Table 4: Definitions, Terms and Abbreviations

| Abbreviation/Term | Explanation                                     |
|-------------------|---|
| AE                | Application Entity                              |
| ANSI              | American National Standard Institute            |
| CD                | Compact Disc                                    |
| CD-R              | CD-Recordable                                   |
| CIS               | Clinical Information System                     |
| DICOM             | Digital Imaging and Communications in Medicine  |
| DIMSE             | DICOM Message Service Element                   |
| DX                | Digital X-Ray                                   |
| EBE               | DICOM Explicit VR Big Endian                    |
| ELE               | DICOM Explicit VR Little Endian                 |
| FSC               | File-set Creator                                |
| FSR               | File-set Reader                                 |
| FSU               | File-set Updater                                |
| GUI               | Graphic User Interface                          |
| ILE               | DICOM Implicit VR Little Endian                 |
| IOD               | Information Object Definition                   |
| ISIS              | Information System - Imaging System             |
| MPPS              | Modality Performed Procedure Step               |
| NEMA              | National Electrical Manufacturers Association   |
| PDU               | Protocol Data Unit                              |
| RIS               | Radiology Information Systems                   |
| RWA               | Real-World Activity                             |
| SC                | Secondary Capture                               |
| SCP               | Service Class Provider                          |
| SCU               | Service Class User                              |
| SOP               | Service Object Pair                             |
| SR                | Structured Report                               |
| TCP/IP            | Transmission Control Protocol/Internet Protocol |
| UID               | Unique Identifier                               |
| WLM               | Worklist Management                             |
| XA                | X-Ray Angiographic                              |



### 3.5. References

[DICOM] NEMA PS3 / ISO 12052, Digital Imaging and Communications in Medicine (DICOM) Standard, National Electrical Manufacturers Association, Rosslyn, VA, USA (available free at <https://medical.nema.org/>)

## 4. Networking

### 4.1. Implementation model

The implementation model consists of three sections:

- The application data flow diagram, specifying the relationship between the Application Entities and the "external world" or Real-World Activities,
- A functional description of each Application Entity, and
- The sequencing constraints among them.

#### 4.1.1. Application Data Flow

Azurion 2.0 has a single Application Entity in its implementation, namely Azurion 2.0 Application Entity. Figure 2 shows the relationship between the Local and Remote Real World Activities.

- After RWA Verify Application Level Communication, the Azurion 2.0 as SCU uses the remote Request Verification SCP functionality to verify communication.
- After RWA Modality Worklist, the Azurion 2.0 as SCU uses the remote Modality Worklist Information Model SCP functionality to query for Modality Worklist.
- After RWA Create and Set Modality Performed Procedure Step, the Azurion 2.0 as SCU uses the remote Modality Performed Procedure Step SOP Class functionality to Report Modality Performed Procedure Step.
- After RWA Transfer Images + Presentation States + X-Ray Dose Structured Report, the Azurion 2.0 as SCU uses the remote SCP Storage Service Class functionality to store local images-Ray Dose Structured Report and presentation states in a remote database.
- The Azurion 2.0 as SCP Storage Service class accepts images +presentation states from the remote SCU.
- After RWA Storage Commitment, the Azurion 2.0 as SCU uses the remote SCP Storage Commitment Service Class functionality to commit remote images and presentation states.
- After RWA Print Images, the Azurion 2.0 as SCU uses the remote SCP Print management Service Class functionality to print the images.
- After RWA Query Retrieve, the Azurion 2.0 as SCU uses the remote SCP Query Retrieve images Service Class functionality to Query Retrieve.




Figure 2 : Functional Overview

#### 4.1.2. Functional Definition of AE's

This section contains a functional definition for each individual local Application Entity.

##### 4.1.2.1. Functional Definition of Application Entity Azurion 2.0

###### Verification Service Class

The Azurion 2.0 can perform (only to pre-configured systems) the Verification service as SCU, the Azurion 2.0 shall request an Association. When the association is accepted by the remote system, the Azurion 2.0 shall send the Verification request, receive the Verification response, and request for releasing the association.

The Azurion 2.0 can also perform the verification service as SCP, it shall accept association when the association is requested by the remote system.

###### Storage Service Class

The Azurion 2.0 can perform (only to pre-configured systems) the Storage service as SCU (RWA Transfer Images + Presentation States and Structured Report), triggered by the operator or by an event in the system, e.g. closing of an examination, acquisition of images. The Azurion 2.0 shall request an association with the selected remote SCP for all applicable Storage SOP classes. When the association is accepted, the Azurion 2.0 shall send the Storage requests (including data from local database), receive the Storage responses and act accordingly, and finally request for releasing the association.

The Azurion 2.0 can also perform (only to pre-configured systems) the Storage service as SCP (RWA accept Images +Presentation States). The Azurion 2.0 shall receive an association with the selected remote SCU for all applicable Storage SOP classes, accept the association and receives the Storage requests responses and act accordingly.

Template ID: HSDPCPF-T-30001.09

### ***Print Management Service Class***

The Azurion 2.0 can perform the Print service as SCU (RWA Print Images), triggered by the operator. For each printed sheet, the Azurion 2.0 shall request an association with the selected remote SCP (i.e., a Print Server) for all applicable SOP classes of the applicable Print Management Meta SOP class. When the association is accepted, the Azurion 2.0 shall send the Print requests including data from local database (the N-GET RQ message to get the printer status, the N-CREATE-RQ message to create the Film Session and the Film Box, the N-SET-RQ message to set the Image Box on the printer, finally, the N-ACTION-RQ message to give printer the command to print), receive the Print responses and act accordingly, and finally request for releasing the association. The Azurion 2.0 can perform the Print service as SCU (RWA Get Printer Status), triggered by the operator in the service mode. The Azurion 2.0 shall request an association with the selected remote SCP (Print Server) for the Printer SOP class. When the association is accepted, the Azurion 2.0 shall send the N-GET request, receive the responses from the Print Server and act accordingly, and finally request for releasing the association.

### ***Basic Worklist Management Service Class***

The Azurion 2.0 can perform (only to the pre-configured MWL SCP) the Basic Worklist Management service as SCU (RWA Request Modality Worklist), triggered by the operator. The Azurion 2.0 shall request an association. When the association is accepted, the Azurion 2.0 shall send the Worklist request, receive the Worklist responses, and request for releasing the association.

### ***Study Management Service Class***

The Azurion 2.0 can perform (only to the pre-configured MWL SCP) the Study Management service as SCU (RWA Create and Set Modality Performed Procedure Step), triggered by the start of an examination (triggered by the first x-ray exposure) for acquisition or closing. The Azurion 2.0 shall request an association. When the association is accepted, the Azurion 2.0 shall send Create and Set requests, receive the responses, and request for releasing the association.

### ***Query Retrieve Service Class***

The Azurion 2.0 can perform the Query Retrieve service as SCU. The Azurion 2.0 shall request an association. When the association is accepted, the Azurion 2.0 shall send a Query/Retrieve request, receive the responses, and request for releasing the association.

## **4.1.3. Sequencing of Real World Activities**

The following sequence of Real World activities are supported by Azurion 2.0:

- The clinical user queries the MWL SCP for a (specific) Worklist representing the list of Scheduled Procedure Steps (with demographic information). Based on that query entered at Azurion 2.0, it sends the C-FIND-RQ message with the query criteria. The automatic broad query with configured query criteria shall be triggered whenever user selects scheduled patients page or all patients page.
- The clinical user starts the examination. As a result, Azurion 2.0 notifies the MWL SCP of the start of a new Procedure Step, i.e. it sends the MPPS N-CREATE-RQ message with the "IN PROGRESS" status of the examination.
- The clinical user acquires images with a certain procedure. As a result, if background image transfer is configured, Azurion 2.0 sends automatically the acquired images (and corresponding presentation states) to the PACS and/or the Workstation, i.e., it sends the C-STORE-RQ messages containing the image (and presentation state) information.
- The clinical user completes the examination. As a result, if auto-transfer is configured, Azurion 2.0 sends images (and corresponding presentation states) and dose report to the PACS and/or Workstation (background image transfer), i.e., it sends the C-STORE-RQ messages containing the image (and presentation state) and dose report information.
- When all images (and presentation states), which were to be automatically transferred to the PACS, have been transferred and storage-commit configured, the Azurion 2.0 asks the PACS to take responsibility for the images (and presentation states) that it has stored that originate from the examination, i.e., it sends the N-ACTION-RQ message containing the request for storage commit.
- Azurion 2.0 notifies the MWL SCP of the completion of a Procedure Step, i.e., it sends the N-SET-RQ message with the "COMPLETED" status of the examination. Additionally to the basic flow of activities, the clinical user may also perform the following steps:
  - The clinical user manually transfers images to the PACS and/or Workstation. As a result, Azurion 2.0 sends the C-STORE-RQ messages containing the image information.
  - The clinical user manually prints selected images. As a result, Azurion 2.0 sends the N-GET-RQ message to get the printer status, the N-CREATE-RQ message to create the Film Session and the Film Box, the N-SET-RQ



message to set the Image Box on the printer. Finally, it sends the N-ACTION-RQ message to give printer the command to print.

Figure 3 presents normal scheduled workflow. Other workflow situations (e.g., unscheduled procedure steps) will have other sequencing constraints. For example, printing could equally take place after the acquired images have been stored or after the examination have been closed or could be omitted completely. Query for images could take place before images have been acquired or could be omitted completely.




Figure 3 : Azurion 2.0 Sequence Diagram

## 4.2. AE Specifications

This section in the DICOM Conformance Statement is a set of Application Entity specifications. There are as many of these subsections as there are different AE's in the implementation.

### 4.2.1. AE Specification of Azurion 2.0

Detail of this specific Application Entity is specified in this section.

#### 4.2.1.1. SOP Classes

This Application Entity provides Standard Conformance to the following SOP Classes.

Table 5: SOP Classes for Azurion 2.0

| SOP Class Name                                  | SOP Class UID           | SCU | SCP |
|---|-------------------------|-----|-----|
| Verification SOP Class                          | 1.2.840.10008.1.1       | Yes | Yes |
| Storage Commitment Push Model SOP Class         | 1.2.840.10008.1.20.1    | Yes | No  |
| Modality Performed Procedure Step SOP Class     | 1.2.840.10008.3.1.2.3.3 | Yes | No  |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9   | Yes | No  |
| >Basic Film Session SOP Class                   | 1.2.840.10008.5.1.1.1   | Yes | No  |
| >Basic Film Box SOP Class                       | 1.2.840.10008.5.1.1.2   | Yes | No  |
| >Basic Grayscale Image Box SOP Class            | 1.2.840.10008.5.1.1.4   | Yes | No  |

Template ID: HSDPCPF-T-30001.09

|   |                               |     |     |
|---|-------------------------------|-----|-----|
| >Printer SOP Class                                      | 1.2.840.10008.5.1.1.16        | Yes | No  |
| Basic Color Print Management Meta SOP Class             | 1.2.840.10008.5.1.1.18        | Yes | No  |
| >Basic Film Session SOP Class                           | 1.2.840.10008.5.1.1.1         | Yes | No  |
| >Basic Film Box SOP Class                               | 1.2.840.10008.5.1.1.2         | Yes | No  |
| >Basic Color Image Box SOP Class                        | 1.2.840.10008.5.1.1.4.1       | Yes | No  |
| >Printer SOP Class                                      | 1.2.840.10008.5.1.1.16        | Yes | No  |
| Secondary Capture Image Storage SOP Class               | 1.2.840.10008.5.1.4.1.1.7     | Yes | Yes |
| Grayscale Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.1  | Yes | Yes |
| X-Ray Angiographic Image Storage SOP Class              | 1.2.840.10008.5.1.4.1.1.12.1  | Yes | Yes |
| X-Ray Radiation Dose SR                                 | 1.2.840.10008.5.1.4.1.1.88.67 | Yes | No  |
| Modality Worklist Information Model - FIND SOP Class    | 1.2.840.10008.5.1.4.31        | Yes | No  |
| Patient Root QR Information Model - FIND SOP Class      | 1.2.840.10008.5.1.4.1.2.1.1   | Yes | No  |
| Study Root QR Information Model - FIND SOP Class        | 1.2.840.10008.5.1.4.1.2.2.1   | Yes | No  |
| Patient Root QR Information Model - MOVE SOP Class      | 1.2.840.10008.5.1.4.1.2.1.2   | Yes | No  |
| Study Root QR Information Model - MOVE SOP Class        | 1.2.840.10008.5.1.4.1.2.2.2   | Yes | No  |

#### 4.2.1.2. Association Policies

Each AE specification contains a description of the general association establishment and acceptance policies of the AE.

#### 4.2.1.3. General

DICOM standard application context name for DICOM 3.0 is always proposed:

**Table 6: DICOM Application Context**

| Description              | Value                 |
|--------------------------|-----------------------|
| Application Context Name | 1.2.840.10008.3.1.1.1 |

#### 4.2.1.4. Number of Associations

Azurion 2.0(SCU) can initialize a maximum of one simultaneous associations. The maximum number of simultaneous associations supported by the Azurion 2.0(SCP) is unlimited by default.

**Table 7: Maximum number of associations as an Association Initiator for Azurion 2.0**

| Property                                    | Value |
|---|-------|
| Maximum number of simultaneous associations | 1     |

#### 4.2.1.5. Asynchronous Nature

The Azurion 2.0 does not support asynchronous operations except for storage commitment. After the storage commitment N-ACTION request is transmitted, storage commitment notification may be handled on another association

#### 4.2.1.6. Implementation Identifying Information

The value supplied for Implementation Class UID and version name are documented here.

**Table 8: DICOM Implementation Class and Version for Azurion 2.0**

| Property                    | Value                    |
|-----------------------------|--------------------------|
| Implementation Class UID    | 1.3.46.670589.7.29.2.0.1 |
| Implementation Version Name | Azurion                  |



#### 4.2.1.7. Communication Failure Handling

The behaviour of this application entity during communication failure is summarized in the table below.

**Table 9: Communication Failure Behavior**

| Exception           | Behaviour   |
|---------------------|---|
| Timeout             | The association is aborted using A-ABORT and the command is marked failed. The reason is logged and reported to the user. |
| Association aborted | The ABORT is handled by failing the job and the command is marked failed. The reason is logged and reported to the user.  |

#### 4.2.1.8. Association Initiation Policy

The Application Entity will respond to a received Association rejection as shown in the next table.

**Table 10: Association Rejection response**

| Result               | Source                  | Reason/Diagnosis                         | Behaviour  |
|----------------------|-------------------------|--|------------|
| 1-rejected-permanent | 1-DICOM UL service-user | 1-no-reason-given                        | Log entry. |
|                      |                         | 2-application-context-name-not-supported | Log entry. |
|                      |                         | 3-calling-AE-title-not-recognized        | Log entry. |
|                      |                         | 7-called- AE-title-not-recognized        | Log entry. |

| Result               | Source  | Reason/Diagnosis                         | Behaviour  |
|----------------------|---|--|------------|
|                      | 2-DICOM UL service-provider (ACSE related function)         | 1-no-reason-given                        | Log entry. |
|                      |   | 2-no-reason-given                        | Log entry. |
|                      | 3-DICOM UL service-provider (Presentation related function) | 1-temporary-congestion                   | Log entry. |
|                      |   | 2-local-limit-exceeded                   | Log entry. |
| 2-rejected-transient | 1-DICOM UL service-user                                     | 1-no-reason-given                        | Log entry. |
|                      |   | 2-application-context-name-not-supported | Log entry. |
|                      |   | 3-calling-AE-title-not-recognized        | Log entry. |
|                      |   | 7-called- AE-title-not-recognized        | Log entry. |
|                      | 2-DICOM UL service-provider (ACSE related function)         | 1-no-reason-given                        | Log entry. |
|                      |   | 2-no-reason-given                        | Log entry. |
|                      | 3-DICOM UL service-provider (Presentation related function) | 1-temporary-congestion                   | Log entry. |
|                      |   | 2-local-limit-exceeded                   | Log entry. |

**Table 11: Association Abort Handling**

| Source                        | Reason/Diagnosis                | Behaviour  |
|-------------------------------|---------------------------------|------------|
| 0 - DICOM UL service-user     | 0 - reason-not-specified        | Log entry. |
| 2 - DICOM UL service-provider | 0 - reason-not-specified        | Log entry. |
|                               | 1 - unrecognized-PDU            | Log entry. |
|                               | 2 - unexpected-PDU              | Log entry. |
|                               | 4 - unrecognized-PDU-parameter  | Log entry. |
|                               | 5 - unexpected-PDU-parameter    | Log entry. |
|                               | 6 - invalid-PDU-parameter-value | Log entry. |

#### 4.2.1.9. (Real-World) Activity – Verification as SCU

##### 4.2.1.9.1.1. Description and Sequencing of Activities

For each Verify Application Level Communication request, an association towards the remote system is established and a C-ECHO request is transmitted. Once the response is received, the association is closed.






Figure 4: Sequencing of RWA Verify Application Level Communication

4.2.1.9.1.2. Proposed Presentation Contexts

Each time an association is initiated, the Azurion 2.0 proposes one presentation contexts to be used on that association. The presentation context proposed by the Azurion 2.0 for Verify Application Level Communication is defined in Table below.

Table 12: Proposed Presentation Contexts for (Real-World) Activity – Verification as SCU

| Presentation Context Table |                   |                           |                     |      |           |
|----------------------------|-------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax            |                   | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name                       | UID               | Name                      | UID                 |      |           |
| Verification SOP Class     | 1.2.840.10008.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|                            |                   | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|                            |                   | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |

4.2.1.9.1.3. SOP Specific Conformance for Verification SOP Class

4.2.1.9.1.3.1. Dataset Specific Conformance for Verification SOP Class C-ECHO-SCU

The behavior of the Azurion 2.0 for status codes in a Verification response is summarized in Table 13.

Table 13: Verification C-ECHO Response Status Handling Behavior

| Service Status | Error Code | Further Meaning | Behaviour   |
|----------------|------------|-----------------|---|
| Success        | 0000       | Confirmation    | The SCP has successfully responded to the verification request. |

The behavior of the Azurion 2.0 during communication failure is summarized in Table 14.

Table 14: Verification Communication Failure Behavior


| Service Status | Error Code |
|----------------|------------|
|----------------|------------|

|                      |  |
|----------------------|--|
| Timeout              | The association is aborted using A-ABORT. The reason is logged and reported to the user. |
| Association aborted  | The reason is logged and failure is reported to the user.                                |
| Association rejected | The reason is logged and failure is reported to the user.                                |

**4.2.1.10. (Real-World) Activity – Modality Worklist as SCU**

**4.2.1.10.1.1. Description and Sequencing of Activities**

For each Broad or Specific Worklist request, the Azurion 2.0 opens an association towards the Basic Worklist Management SCP and sends a C-FIND request. After retrieval of all responses containing matching Worklist items, the association is closed (see Figure 5). All returned Worklist items are displayed to the operator who can select an item from the Worklist and perform an examination.



**Figure 5 : Sequencing of RWA Request Modality Worklist**

The clinical user may cancel the query to the MWL SCP. As a result, Azurion 2.0 sends a C-FIND Cancel Request to the MWL SCP.

**4.2.1.10.1.2. Proposed Presentation Contexts**

The presentation contexts are defined in the Table 15.

**Table 15: Proposed Presentation Contexts for (Real-World) Activity – Modality worklist as SCU**

| Presentation Context Table                           |                        |                           |                     |      |           |
|--|------------------------|---------------------------|---------------------|------|-----------|
| Abstract Syntax                                      |                        | Transfer Syntax           |                     | Role | Ext. Neg. |
| Name   | UID                    | Name                      | UID                 |      |           |
| Modality Worklist Information Model - FIND SOP Class | 1.2.840.10008.5.1.4.31 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None      |
|  |                        | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |           |
|  |                        | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |           |

#### 4.2.1.10.1.3. SOP Specific Conformance for Modality Worklist Information Model - FIND SOP Class

Two kinds of queries can be done with the Azurion 2.0: a broad query and a specific query. A broad query for the Worklist is initiated by the operator without filling in any search criteria (the search criteria are based on system configuration). The Matching Keys are presented in Table 16.

A specific Worklist request is initiated by the operator after filling in search criteria in the Graphical User Interface. At least one key should be specified. No verification of query results in relation to the original query criteria is done. The Matching Keys are presented in Table 17.

A received Worklist entry is validated. The entry will be discarded, and an error will be reported when a type-one or type-two attribute is missing, or when the translation of a type-one attribute fails (this includes individual attributes within a sub-sequence).

**Table 16: Broad Query Presentation Contexts for (Real-World) Activity Modality Worklist as SCU**

| Attribute Name                      | Tag         | Matching Key                                 |
|-------------------------------------|-------------|--|
| Scheduled Station AE Title          | (0040,0001) | Single value matching                        |
| Scheduled Procedure Step Start Date | (0040,0002) | Universal matching or range matching         |
| Scheduled Procedure Step Start Time | (0040,0003) | Universal matching or range matching         |
| Modality                            | (0008,0060) | Single value matching and Universal matching |

**Table 17: Matching Table MWL Information Model – Specific Query**

| Attribute Name                      | Tag         | Matching Key  |
|-------------------------------------|-------------|---|
| Scheduled Station AE Title          | (0040,0001) | Universal matching or single value matching                       |
| Scheduled Procedure Step Start Date | (0040,0002) | Universal matching or range matching                              |
| Modality                            | (0008,0060) | Universal matching or single value matching                       |
| Patient's Name                      | (0010,0010) | Universal matching or single value matching or wild card matching |
| Patient ID                          | (0010,0020) | Universal matching or single value matching                       |
| Accession Number                    | (0008,0050) | Universal matching or single value matching                       |
| Requested Procedure ID              | (0040,1001) | Universal matching or single value matching                       |

#### 4.2.1.10.1.3.1. Dataset Specific Conformance for Modality Worklist Information Model - FIND SOP Class C-FIND-SCU

This section specifies the Modality Worklist Request Attributes.

For each attribute in the following information is supplied:

Attribute Name: Attributes supported to build a Modality Worklist Request Identifier.

Tag: DICOM tag for this attribute.

VR: DICOM VR for this attribute.

M: Matching Keys for (automatic) Worklist Update.

R: Return Keys. An "X" will indicate that this attribute as matching key can be used

Q: Interactive Query Key. An "X" will indicate that this attribute as matching key can be used.

D: Displayed Keys. An "X" indicates that this Worklist attribute is displayed on the user during a patient registration dialog.

IOD: An "X" indicates that this Worklist attribute is included into all object Instances created during performance of the related Procedure Step.

Type of Matching: The following types of matching exists:

- Single Value Matching
- List of UID Matching

Template ID: HSDPCPF-T-30001.09



- Wild Card Matching
- Range Matching
- Sequence Matching
- Universal Matching

**Table 18: C-FIND-RQ Dataset Specification.**

| Modality Worklist Information Model - FIND SOP Class |           |    |   |   |   |   |     |  |         |
|--|-----------|----|---|---|---|---|-----|--|---------|
| Attribute Name                                       | Tag       | VR | M | R | Q | D | IOD | Type Of Matching                                   | Comment |
| <b>Scheduled Procedure Step</b>                      |           |    |   |   |   |   |     |  |         |
| Scheduled Procedure Step Sequence                    | 0040,0100 | SQ |   | X |   |   |     | NA   |         |
| >Modality  | 0008,0060 | CS | X | X | X | X | X   | Single Value,<br>Universal,<br>Wild Card           |         |
| >Requested Contrast Agent                            | 0032,1070 | LO | X | X |   |   |     | Universal  |         |
| >Scheduled Station AE Title                          | 0040,0001 | AE | X | X | X | X |     | Single Value,<br>Universal,<br>Wild Card           |         |
| >Scheduled Procedure Step Start Date                 | 0040,0002 | DA | X | X | X | X |     | Single Value,<br>Universal,<br>Wild Card,<br>range |         |
| >Scheduled Procedure Step Start Time                 | 0040,0003 | TM | X | X |   | X |     | Single Value,<br>Range,<br>Universal               |         |
| >Scheduled Performing Physician's Name               | 0040,0006 | PN | X | X | X | X |     | Single Value,<br>Universal,<br>Wild Card           |         |
| >Scheduled Procedure Step Description                | 0040,0007 | LO | X |   |   |   |     | Universal  |         |
| >Scheduled Protocol Code Sequence                    | 0040,0008 | SQ | X |   |   |   |     | Universal  |         |
| >>Code Value   | 0008,0100 | SH | X |   |   |   |     | Universal  |         |
| >>Coding Scheme Designator                           | 0008,0102 | SH | X |   |   |   |     | Universal  |         |
| >>Coding Scheme Version                              | 0008,0103 | SH | X |   |   |   |     | Universal  |         |
| >>Code Meaning                                       | 0008,0104 | LO | X |   |   |   |     | Universal  |         |
| >>Protocol Context Sequence                          | 0040,0440 | SQ | X |   |   |   |     | Universal  |         |
| >>>Measurement Units Code Sequence                   | 0040,08EA | SQ | X |   |   |   |     | Universal  |         |
| >>>>Code Value                                       | 0008,0100 | SH | X |   |   |   |     | Universal  |         |
| >>>>Coding Scheme Designator                         | 0008,0102 | SH | X |   |   |   |     | Universal  |         |
| >>>>Coding Scheme Version                            | 0008,0103 | SH | X |   |   |   |     | Universal  |         |
| >>>>Code Meaning                                     | 0008,0104 | LO | X |   |   |   |     | Universal  |         |
| >>>Value Type  | 0040,A040 | CS | X |   |   |   |     | Universal  |         |
| >>>Concept Name Code Sequence                        | 0040,A043 | SQ | X |   |   |   |     | Universal  |         |
| >>>>Code Value                                       | 0008,0100 | SH | X |   |   |   |     | Universal  |         |
| >>>>Coding Scheme Designator                         | 0008,0102 | SH | X |   |   |   |     | Universal  |         |
| >>>>Coding Scheme Version                            | 0008,0103 | SH | X |   |   |   |     | Universal  |         |
| >>>>Code Meaning                                     | 0008,0104 | LO | X |   |   |   |     | Universal  |         |
| >>>DateTime  | 0040,A120 | DT | X |   |   |   |     | Universal  |         |
| >>>Person Name                                       | 0040,A123 | PN | X |   |   |   |     | Universal  |         |



**Modality Worklist Information Model - FIND SOP Class**

| Attribute Name                     | Tag       | VR | M | R | Q | D | IOD | Type Of Matching                         | Comment |
|------------------------------------|-----------|----|---|---|---|---|-----|--|---------|
| <b>Scheduled Procedure Step</b>    |           |    |   |   |   |   |     |  |         |
| >>>Text Value                      | 0040,A160 | UT | X |   |   |   |     | Universal                                |         |
| >>>Concept Code Sequence           | 0040,A168 | SQ | X |   |   |   |     | Universal                                |         |
| >>>>Code Value                     | 0008,0100 | SH | X |   |   |   |     | Universal                                |         |
| >>>>Coding Scheme Designator       | 0008,0102 | SH | X |   |   |   |     | Universal                                |         |
| >>>>Coding Scheme Version          | 0008,0103 | SH | X |   |   |   |     | Universal                                |         |
| >>>>Code Meaning                   | 0008,0104 | LO | X |   |   |   |     | Universal                                |         |
| >>>Numeric Value                   | 0040,A30A | DS | X |   |   |   |     | Universal                                |         |
| >Scheduled Procedure Step ID       | 0040,0009 | SH | X |   |   |   |     | Universal                                |         |
| >Scheduled Station Name            | 0040,0010 | SH | X |   |   |   |     | Universal                                |         |
| >Scheduled Procedure Step Location | 0040,0011 | SH | X |   |   |   |     | Universal                                |         |
| >Pre-Medication                    | 0040,0012 | LO | X |   |   |   |     | Universal                                |         |
| >Scheduled Procedure Step Status   | 0040,0020 | CS | X |   |   |   |     | Universal                                |         |
| <b>Requested Procedure</b>         |           |    |   |   |   |   |     |  |         |
| Study Instance UID                 | 0020,000D | UI | X | X |   |   |     | Universal                                |         |
| Requested Procedure Code Sequence  | 0032,1064 | SQ | X |   |   |   |     | Universal                                |         |
| >Code Value                        | 0008,0100 | SH | X |   |   |   |     | Universal                                |         |
| >Coding Scheme Designator          | 0008,0102 | SH | X |   |   |   |     | Universal                                |         |
| >Coding Scheme Version             | 0008,0103 | SH | X |   |   |   |     | Universal                                |         |
| >Code Meaning                      | 0008,0104 | LO | X |   |   |   |     | Universal                                |         |
| Requested Procedure ID             | 0040,1001 | SH | X |   | X | X |     | Single Value,<br>Universal,<br>Wild Card |         |
| Patient Transport Arrangements     | 0040,1004 | LO | X |   |   |   |     | Universal                                |         |
| <b>Imaging Service Request</b>     |           |    |   |   |   |   |     |  |         |
| Accession Number                   | 0008,0050 | SH | X |   | X | X | X   | Single Value,<br>Universal,<br>Wild Card |         |
| Referring Physician's Name         | 0008,0090 | PN | X | X |   |   | X   | Universal                                |         |
| Requesting Physician               | 0032,1032 | PN | X |   |   |   |     | Universal                                |         |
| <b>Visit Relationship</b>          |           |    |   |   |   |   |     |  |         |
| Referenced Patient Sequence        | 0008,1120 | SQ |   |   |   |   |     | NA                                       |         |
| >Referenced SOP Class UID          | 0008,1150 | UI | X |   |   |   |     | Universal                                |         |
| >Referenced SOP Instance UID       | 0008,1155 | UI | X |   |   |   |     | Universal                                |         |
| <b>Patient Identification</b>      |           |    |   |   |   |   |     |  |         |
| Patient's Name                     | 0010,0010 | PN | X | X | X | X | X   | Single Value,<br>Universal,<br>Wild Card |         |
| Patient ID                         | 0010,0020 | LO | X | X | X | X | X   | Single Value,<br>Universal,<br>Wild Card |         |
| <b>Patient Demographic</b>         |           |    |   |   |   |   |     |  |         |
| Patients Birth Date                | 0010,0030 | DA | X | X |   |   | X   | Single, Universal,                       |         |

**Modality Worklist Information Model - FIND SOP Class**

| Attribute Name  | Tag       | VR | M | R | Q | D | IOD | Type Of Matching        | Comment    |
|---|-----------|----|---|---|---|---|-----|-------------------------|------------|
| <b>Scheduled Procedure Step</b>                               |           |    |   |   |   |   |     |                         |            |
|   |           |    |   |   |   |   |     | Range, Wild Card        |            |
| Patient's Sex   | 0010,0040 | CS | X | X |   |   | X   | Single Value, Universal |            |
| Patient's Weight  | 0010,1030 | DS | X | X |   |   | X   | Universal               |            |
| <b>Patient Medical</b>  |           |    |   |   |   |   |     |                         |            |
| Medical Alerts  | 0010,2000 | LO | X |   |   |   |     | Universal               |            |
| Allergies   | 0010,2110 | LO | X |   |   |   |     | Universal               |            |
| Pregnancy Status  | 0010,21C0 | US | X | X |   |   |     | Universal               |            |
| Patient State   | 0038,0500 | LO | X | X |   |   |     | Universal               |            |
| <b>Attributes For The Modality Worklist C-Find Identifier</b> |           |    |   |   |   |   |     |                         |            |
| Specific Character Set  | 0008,0005 | CS | X |   |   |   | X   | Universal               | ISO_IR 100 |

The behavior of the Azurion 2.0 for status codes in C-FIND response is summarized in Table 19.

**Table 19: Modality Worklist C-FIND Response Status Handling Behavior**

| Service Status | Error Code            | Further Meaning   | Behaviour  |
|----------------|-----------------------|---|--|
| Success        | 0000                  | Matching is complete - No final Identifier is supplied.   | The result is reported to the user and is logged.  |
| Refused        | A700                  | Out of Resources  | Stops with processing the C-FIND Response(s) from the SCP. No responses displayed to the user.   |
| Failed         | A900                  | Identifier Does Not Match SOP Class   | Stops with processing the C-FIND Response(s) from the SCP. The reason is logged and the failure is reported to the user. No responses displayed to the user. |
|                | C001                  | Unable to process   | Stops with processing the C-FIND Response(s) from the SCP. The reason is logged and the failure is reported to the user. No responses displayed to the user. |
| Cancel         | FE00                  | Matching terminated due to Cancel Match request   | Stops with processing the C-FIND Response(s) from the SCP. No responses displayed to the user.   |
| Pending        | FF00                  | Matches are continuing – Current. Match is supported in the same manner as supplied and any Optional Keys were Required Keys. | Continues with processing of the C-FIND Response(s) from the SCP   |
|                | FF01                  | Matches are continuing – Warning that one or more Optional Keys were not supported for existence for this Identifier.         | Continues with processing of the C-FIND Response(s) from the SCP.  |
| *              | Any other status code | *   | The association is aborted using A-ABORT. The reason is logged and the failure is reported to the user. No responses displayed to the user.                  |

The behavior of the Azurion 2.0 during communication failure is summarized in Table 20.



**Table 20: Modality Worklist Communication Failure Behavior**

| Exception            | Behaviour  |
|----------------------|--|
| Timeout              | The query is marked as failed. The association is aborted using A-ABORT. The reason is logged and reported to the user. The Azurion 2.0 stops processing the C-FIND Response(s) from the SCP.    |
| Association Aborted  | If the association is aborted using A-ABORT, the query is marked as failed. The reason is logged and failure is reported to the user. Stops with processing the C-FIND Response(s) from the SCP. |
| Association Rejected | The query is marked as failed. The reason is logged and failure is reported to the user. No C-FIND request performed.  |

**4.2.1.11. (Real-World) Activity – Modality Performed Procedure Step as SCU**

**4.2.1.11.1. Description and Sequencing of Activities**

For each MPPS Job, a new association towards the Modality Performed Procedure Step is established and closed when the MPPS Job has been transmitted. There are three kinds of MPPS Jobs:

- *MPPS Create Job.* Only an N-CREATE request with status “IN PROGRESS” is transmitted. Once the response is received, the association is closed. The MPPS Create Job is submitted when:
  - A Worklist examination is selected for acquisition;
  - A local examination is selected for acquisition and the system is configured to be connected to an IHE compatible MWL SCP.
  - An already Completed Worklist examination or an already Completed, local examination, is re-selected for acquisition (IHE Append Use Case).
  
- *MPPS Set Job.* Only an N-SET request is transmitted. The status field will respectively be set to “DISCONTINUED” or “COMPLETED”. Once the response is received, the association is closed. The MPPS Set Job is submitted when:
  - An examination is discontinued/deleted/restored to solve patient mixing or closed and the MPPS Create Job is already handled (transmitted).
- *MPPS Create & Set Job.* Over the same association both the N-CREATE request and the N-SET request corresponding to the same examination are transmitted. The MPPS Create & Set Job is submitted when:
  - An examination is discontinued/deleted/restored to solve patient mixing and the MPPS Create Job is not already handled (not transmitted). In such a case the MPPS Create & Set Job replaces the MPPS Create Job.




Figure 6: Sequencing of RWA separate MPPS Create and separate MPPS Set Job






Figure 7: Sequencing of RWA MPPS Create and Set Job

4.2.1.11.1.2. Proposed Presentation Contexts

Each time an association is initiated, the Azurion 2.0 proposes one presentation context to be used on that association. The presentation context proposed by the Azurion 2.0 for Create and Set Modality Performed Procedure Step is defined in table 21.

The implementation chooses Explicit VR Little Endian transfer syntax in case multiple transfer syntaxes are accepted in the association acceptance. All details regarding the specific conformance, including response behavior to all status codes, both from an application level and communication errors are provided in table 23.

Table 21: Proposed Presentation Contexts for (Real-World) Activity – MPPS as SCU

| Presentation Context Table                  |                         |                           |                     |      |                      |
|---|-------------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                             |                         | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name  | UID                     | Name                      | UID                 |      |                      |
| Modality Performed Procedure Step SOP Class | 1.2.840.10008.3.1.2.3.3 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None                 |
|   |                         | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |
|   |                         | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |                      |

4.2.1.11.1.3. SOP Specific Conformance Modality Performed Procedure Step SOP Class

The set of attributes within an N-CREATE and N-SET messages is fixed and it does not depend on configuration settings. In an N-CREATE message, all possible attributes and attribute sequences used in the N-SET are forecasted by defining the attributes

and settings their values to NULL. When an N-SET message is transmitted, it may occur that a forecasted attribute isn't actually used. Table 22 up till Table 24 indicate whether or not an attribute and attribute value is sent during N-CREATE.

Table 24 up till Table 26 indicate whether or not an attribute and attribute value is sent during MPPS N-SET.

**4.2.1.11.1.3.1. Dataset Specific Conformance for Modality Performed Procedure Step SOP Class N-CREATE-SCU**

**Table 22: N-CREATE-RQ Dataset Specification.**

| Modality Performed Procedure Step SOP Class  |           |    |       |                                      |
|--|-----------|----|-------|--------------------------------------|
| Attribute Name                               | Tag       | VR | Value | Comment                              |
| <b>SOP Common Module</b>                     |           |    |       |                                      |
| Specific Character Set                       | 0008,0005 | CS |       | ISO_IR 100                           |
| <b>Performed Procedure Step Relationship</b> |           |    |       |                                      |
| Referenced Patient Sequence                  | 0008,1120 | SQ |       |                                      |
| >Referenced SOP Class UID                    | 0008,1150 | UI |       |                                      |
| >Referenced SOP Instance UID                 | 0008,1155 | UI |       |                                      |
| Patient's Name                               | 0010,0010 | PN |       |                                      |
| Patient ID                                   | 0010,0020 | LO |       |                                      |
| Patient's Birth Date                         | 0010,0030 | DA |       |                                      |
| Patient's Sex                                | 0010,0040 | CS |       |                                      |
| Scheduled Step Attributes Sequence           | 0040,0270 | SQ |       |                                      |
| >Accession Number                            | 0008,0050 | SH |       | Empty in case of an unscheduled exam |
| >Referenced Study Sequence                   | 0008,1110 | SQ |       | Empty in case of an unscheduled exam |
| >>Referenced SOP Class UID                   | 0008,1150 | UI |       |                                      |
| >>Referenced SOP Instance UID                | 0008,1155 | UI |       |                                      |
| >Study Instance UID                          | 0020,000D | UI |       |                                      |
| >Requested Procedure Description             | 0032,1060 | LO |       | Empty in case of an unscheduled exam |
| >>Code Value                                 | 0008,0100 | SH |       |                                      |
| >>Coding Scheme Designator                   | 0008,0102 | SH |       |                                      |
| >>Code Meaning                               | 0008,0104 | LO |       |                                      |
| >>Context Group Extension Flag               | 0008,010B | CS |       |                                      |
| >>Context Identifier                         | 0008,010F | CS |       |                                      |
| >>Context UID                                | 0008,0117 | UI |       |                                      |
| >Scheduled Procedure Step Description        | 0040,0007 | LO |       | Empty in case of an unscheduled exam |
| >>Code Value                                 | 0008,0100 | SH |       |                                      |
| >>Coding Scheme Designator                   | 0008,0102 | SH |       |                                      |
| >>Code Meaning                               | 0008,0104 | LO |       |                                      |
| >>Context Group Extension Flag               | 0008,010B | CS |       |                                      |
| >>Context Identifier                         | 0008,010F | CS |       |                                      |
| >>Context UID                                | 0008,0117 | UI |       |                                      |
| >>Protocol Context Sequence                  | 0040,0440 | SQ |       |                                      |
| >>>Content Item Modifier Sequence            | 0040,0441 | SQ |       |                                      |
| >>>>Value Type                               | 0040,A040 | CS |       |                                      |
| >>>>Concept Name Code Sequence               | 0040,A043 | SQ |       |                                      |
| >>>>>Code Value                              | 0008,0100 | SH |       |                                      |
| >>>>>Coding Scheme Designator                | 0008,0102 | SH |       |                                      |
| >>>>>Code Meaning                            | 0008,0104 | LO |       |                                      |
| >>>>>Context Group Extension Flag            | 0008,010B | CS |       |                                      |

**Modality Performed Procedure Step SOP Class**

| Attribute Name                   | Tag       | VR | Value | Comment                              |
|----------------------------------|-----------|----|-------|--------------------------------------|
| >>>>Context Identifier           | 0008,010F | CS |       |                                      |
| >>>>Context UID                  | 0008,0117 | UI |       |                                      |
| >>>Value Type                    | 0040,A040 | CS |       |                                      |
| >>>Concept Name Code Sequence    | 0040,A043 | SQ |       |                                      |
| >>>>Code Value                   | 0008,0100 | SH |       |                                      |
| >>>>Coding Scheme Designator     | 0008,0102 | SH |       |                                      |
| >>>>Code Meaning                 | 0008,0104 | LO |       |                                      |
| >>>>Context Group Extension Flag | 0008,010B | CS |       |                                      |
| >>>>Context Identifier           | 0008,010F | CS |       |                                      |
| >>>>Context UID                  | 0008,0117 | UI |       |                                      |
| >Scheduled Procedure Step ID     | 0040,0009 | SH |       | Empty in case of an unscheduled exam |
| >Requested Procedure ID          | 0040,1001 | SH |       | Empty in case of an unscheduled exam |

**Performed Procedure Step Information**

|                                      |           |    |  |   |
|--------------------------------------|-----------|----|--|---|
| Procedure Code Sequence              | 0008,1032 | SQ |  |   |
| >Code Value                          | 0008,0100 | SH |  |   |
| >Coding Scheme Designator            | 0008,0102 | SH |  |   |
| >Code Meaning                        | 0008,0104 | LO |  |   |
| >Context Group Extension Flag        | 0008,010B | CS |  |   |
| >Context Identifier                  | 0008,010F | CS |  |   |
| >Context UID                         | 0008,0117 | UI |  |   |
| Performed Station AE Title           | 0040,0241 | AE |  | AE Title as configured by the MWL SCP/CIS unit. |
| Performed Station Name               | 0040,0242 | SH |  |   |
| Performed Location                   | 0040,0243 | SH |  | Always Empty                                    |
| Performed Procedure Step Start Date  | 0040,0244 | DA |  |   |
| Performed Procedure Step Start Time  | 0040,0245 | TM |  |   |
| Performed Procedure Step End Date    | 0040,0250 | DA |  |   |
| Performed Procedure Step End Time    | 0040,0251 | TM |  |   |
| Performed Procedure Step Status      | 0040,0252 | CS |  |   |
| Performed Procedure Step ID          | 0040,0253 | SH |  |   |
| Performed Procedure Step Description | 0040,0254 | LO |  |   |
| Performed Procedure Type Description | 0040,0255 | LO |  |   |

**Image Acquisition Results**

|                                  |           |    |  |  |
|----------------------------------|-----------|----|--|--|
| Modality                         | 0008,0060 | CS |  | Applied Value(s): XA   |
| Study ID                         | 0020,0010 | SH |  | If no Study ID is known, the Accession Number will be used as value. |
| Performed Protocol Code Sequence | 0040,0260 | SQ |  |  |
| >Code Value                      | 0008,0100 | SH |  |  |
| >Coding Scheme Designator        | 0008,0102 | SH |  |  |
| >Coding Scheme Version           | 0008,0103 | SH |  |  |
| >Code Meaning                    | 0008,0104 | LO |  |  |
| >Context Group Extension Flag    | 0008,010B | CS |  |  |
| >Context Identifier              | 0008,010F | CS |  |  |
| >Context UID                     | 0008,0117 | UI |  |  |

**Modality Performed Procedure Step SOP Class**

| Attribute Name                        | Tag       | VR | Value | Comment   |
|---------------------------------------|-----------|----|-------|---|
| >Protocol Context Sequence            | 0040,0440 | SQ |       |   |
| >>Content Item Modifier Sequence      | 0040,0441 | SQ |       |   |
| >>>Value Type                         | 0040,A040 | CS |       |   |
| >>>Concept Name Code Sequence         | 0040,A043 | SQ |       |   |
| >>>>Code Value                        | 0008,0100 | SH |       |   |
| >>>>Coding Scheme Designator          | 0008,0102 | SH |       |   |
| >>>>Code Meaning                      | 0008,0104 | LO |       |   |
| >>>>Context Group Extension Flag      | 0008,010B | CS |       |   |
| >>>>Context Identifier                | 0008,010F | CS |       |   |
| >>>>Context UID                       | 0008,0117 | UI |       |   |
| >>Value Type                          | 0040,A040 | CS |       |   |
| >>Concept Name Code Sequence          | 0040,A043 | SQ |       |   |
| >>>Code Value                         | 0008,0100 | SH |       |   |
| >>>Coding Scheme Designator           | 0008,0102 | SH |       |   |
| >>>Code Meaning                       | 0008,0104 | LO |       |   |
| >>>Context Group Extension Flag       | 0008,010B | CS |       |   |
| >>>Context Identifier                 | 0008,010F | CS |       |   |
| >>>Context UID                        | 0008,0117 | UI |       |   |
| Performed Series Sequence             | 0040,0340 | SQ |       | Sequence will be empty when there are no images to report |
| >Retrieve AE Title                    | 0008,0054 | AE |       |   |
| >Series Description                   | 0008,103E | LO |       |   |
| >Performing Physician's Name          | 0008,1050 | PN |       |   |
| >Operators' Name                      | 0008,1070 | PN |       |   |
| >Operator Identification Sequence     | 0008,1072 | SQ |       |   |
| >>Institution Address                 | 0008,0081 | ST |       |   |
| >>Person Identification Code Sequence | 0040,1101 | SQ |       |   |
| >>>Code Value                         | 0008,0100 | SH |       |   |
| >>>Coding Scheme Designator           | 0008,0102 | SH |       |   |
| >>>Code Meaning                       | 0008,0104 | LO |       |   |
| >>>Context Group Extension Flag       | 0008,010B | CS |       |   |
| >>>Context Identifier                 | 0008,010F | CS |       |   |
| >>>Context UID                        | 0008,0117 | UI |       |   |
| >>Person's Address                    | 0040,1102 | ST |       |   |
| >>Person's Telephone Numbers          | 0040,1103 | LO |       |   |
| >Referenced Image Sequence            | 0008,1140 | SQ |       |   |
| >>Referenced SOP Class UID            | 0008,1150 | UI |       |   |
| >>Referenced SOP Instance UID         | 0008,1155 | UI |       |   |
| >>Container Identifier                | 0040,0512 | LO |       |   |
| >>Specimen Description Sequence       | 0040,0560 | SQ |       |   |
| >>>Specimen Identifier                | 0040,0551 | LO |       |   |
| >>>Specimen UID                       | 0040,0554 | UI |       |   |

**Modality Performed Procedure Step SOP Class**

| Attribute Name  | Tag       | VR | Value | Comment      |
|---|-----------|----|-------|--------------|
| >Protocol Name  | 0018,1030 | LO |       |              |
| >Series Instance UID                                  | 0020,000E | UI |       |              |
| >Referenced Non-Image Composite SOP Instance Sequence | 0040,0220 | SQ |       |              |
| >>Referenced SOP Class UID                            | 0008,1150 | UI |       |              |
| >>Referenced SOP Instance UID                         | 0008,1155 | UI |       |              |
| Radiation Dose  |           |    |       |              |
| Image and Fluoroscopy Area Dose Product               | 0018,115E | DS |       |              |
| Total Time of Fluoroscopy                             | 0040,0300 | US |       |              |
| Total Number of Exposures                             | 0040,0301 | US |       |              |
| Entrance Dose   | 0040,0302 | US |       |              |
| Entrance Dose in mGy                                  | 0040,8302 | DS |       |              |
| Billing And Material Management Codes                 |           |    |       |              |
| Film Consumption Sequence                             | 0040,0321 | SQ |       | Always Empty |
| >Medium Type  | 2000,0030 | CS |       |              |
| >Film Size ID   | 2010,0050 | CS |       |              |
| >Number of Films                                      | 2100,0170 | IS |       |              |

The behavior of the Azurion 2.0 for status codes in an MPPS N-CREATE response and NSET response is presented in Table 23 and Table 24 respectively. In case of the retransmission attempt each message stored in the persistent queue is sent over a separate association.

**Table 23: MPPS N-CREATE Response Status Handling Behavior**

| Service Status | Error Code            | Further Meaning       | Behavior   |
|----------------|-----------------------|-----------------------|--|
| Success        | 0000                  | Successful operation. | The notify status of the related examination is updated (set to in progress). The examination status is not changed (e.g. still in progress).  |
| Failure        | 0213                  | Resource Limitation   | The message contents is made persistent and the message is added to the persistent queue and waits for the next retransmission attempt. The examination status is not changed (e.g. still in progress).  |
| *              | Any other status code | *                     | If the response status is reported during initial transmission the message contents is made persistent and the message is added to the persistent queue. If this response status is the result of the retransmission attempt related examination is updated to the state as if the transmission succeeded. This means that the notify status of the related examination is updated (set to in progress) and if this response status is the result of the retransmission attempt message is removed from the persistent queue. The examination status is not changed (e.g. still in progress).The response status is logged as a warning. |

**Table 24: MPPS Communication Failure Behavior (N-SET, N-CREATE) Exception Behavior**

| Exception            | Behavior   |
|----------------------|--|
| Timeout              | The Association is aborted using A-ABORT. The reason is logged and reported to the user. The message content is made persistent and the message is added to the persistent queue and waits for the next retransmission attempt. The examination status is set to CLOSED. |
| Association Aborted  | The command is marked as failed. The reason is logged and reported to the user. The message content is made persistent and the message is added to the persistent queue and waits for the next retransmission attempt. The examination status is set to CLOSED.          |
| Association Rejected | The command is marked as failed. The reason is logged and reported to the user. The message content is made persistent and the message is added to the persistent queue and waits for the next retransmission attempt. The examination status is set to CLOSED.          |

**4.2.1.11.1.3.2. Dataset Specific Conformance for Modality Performed Procedure Step SOP Class N-SET-SCU**

**Table 25: N-SET-RQ Dataset Specification.**

| Modality Performed Procedure Step SOP Class |           |    |       |  |
|---|-----------|----|-------|--|
| Attribute Name                              | Tag       | VR | Value | Comment                                    |
| <b>Sop Common Module</b>                    |           |    |       |  |
| Specific Character Set                      | 0008,0005 | CS |       | ISO_IR100                                  |
| <b>Performed Procedure Step Information</b> |           |    |       |  |
| Procedure Code Sequence                     | 0008,1032 | SQ |       |  |
| >Code Value                                 | 0008,0100 | SH |       |  |
| >Coding Scheme Designator                   | 0008,0102 | SH |       |  |
| >Code Meaning                               | 0008,0104 | LO |       |  |
| >Context Group Extension Flag               | 0008,010B | CS |       |  |
| >Context Identifier                         | 0008,010F | CS |       |  |
| >Context UID                                | 0008,0117 | UI |       |  |
| Performed Procedure Step End Date           | 0040,0250 | DA |       |  |
| Performed Procedure Step End Time           | 0040,0251 | TM |       |  |
| Performed Procedure Step Status             | 0040,0252 | CS |       | Applied Values: COMPLETED or DISCONTINUED  |
| Performed Procedure Step Description        | 0040,0254 | LO |       | May be Empty by configuration              |
| Performed Procedure Type Description        | 0040,0255 | LO |       |  |
| <b>Image Acquisition Results</b>            |           |    |       |  |
| Performed Protocol Code Sequence            | 0040,0260 | SQ |       | Sequence remains empty                     |
| >Code Value                                 | 0008,0100 | SH |       |  |
| >Coding Scheme Designator                   | 0008,0102 | SH |       |  |
| >Code Meaning                               | 0008,0104 | LO |       |  |
| >Context Group Extension Flag               | 0008,010B | CS |       |  |
| >Context Identifier                         | 0008,010F | CS |       |  |
| >Context UID                                | 0008,0117 | UI |       |  |
| >Protocol Context Sequence                  | 0040,0440 | SQ |       | May be empty when no images to be reported |
| >>Content Item Modifier Sequence            | 0040,0441 | SQ |       |  |
| >>>Value Type                               | 0040,A040 | CS |       |  |
| >>>Concept Name Code Sequence               | 0040,A043 | SQ |       |  |

**Modality Performed Procedure Step SOP Class**

| Attribute Name                                | Tag       | VR | Value | Comment |
|---|-----------|----|-------|---------|
| >>>>Code Value                                | 0008,0100 | SH |       |         |
| >>>>Coding Scheme Designator                  | 0008,0102 | SH |       |         |
| >>>>Code Meaning                              | 0008,0104 | LO |       |         |
| >>>>Context Group Extension Flag              | 0008,010B | CS |       |         |
| >>>>Context Identifier                        | 0008,010F | CS |       |         |
| >>>>Context UID                               | 0008,0117 | UI |       |         |
| >>Value Type                                  | 0040,A040 | CS |       |         |
| >>Concept Name Code Sequence                  | 0040,A043 | SQ |       |         |
| >>>Code Value                                 | 0008,0100 | SH |       |         |
| >>>Coding Scheme Designator                   | 0008,0102 | SH |       |         |
| >>>Code Meaning                               | 0008,0104 | LO |       |         |
| >>>Context Group Extension Flag               | 0008,010B | CS |       |         |
| >>>Context Identifier                         | 0008,010F | CS |       |         |
| >>>Context UID                                | 0008,0117 | UI |       |         |
| Performed Series Sequence                     | 0040,0340 | SQ |       |         |
| >Retrieve AE Title                            | 0008,0054 | AE |       |         |
| >Series Description                           | 0008,103E | LO |       |         |
| >Series Description Code Sequence             | 0008,103F | SQ |       |         |
| >>Code Value                                  | 0008,0100 | SH |       |         |
| >>Coding Scheme Designator                    | 0008,0102 | SH |       |         |
| >>Code Meaning                                | 0008,0104 | LO |       |         |
| >>Context Group Extension Flag                | 0008,010B | CS |       |         |
| >>Context Identifier                          | 0008,010F | CS |       |         |
| >>Context UID                                 | 0008,0117 | UI |       |         |
| >Performing Physician's Name                  | 0008,1050 | PN |       |         |
| >Performing Physician Identification Sequence | 0008,1052 | SQ |       |         |
| >>Institution Address                         | 0008,0081 | ST |       |         |
| >>Person Identification Code Sequence         | 0040,1101 | SQ |       |         |
| >>>Code Value                                 | 0008,0100 | SH |       |         |
| >>>Coding Scheme Designator                   | 0008,0102 | SH |       |         |
| >>>Code Meaning                               | 0008,0104 | LO |       |         |
| >>>Context Group Extension Flag               | 0008,010B | CS |       |         |
| >>>Context Identifier                         | 0008,010F | CS |       |         |
| >>>Context UID                                | 0008,0117 | UI |       |         |
| >>Person's Address                            | 0040,1102 | ST |       |         |
| >>Person's Telephone Numbers                  | 0040,1103 | LO |       |         |
| >Operators' Name                              | 0008,1070 | PN |       |         |
| >Operator Identification Sequence             | 0008,1072 | SQ |       |         |
| >>Institution Address                         | 0008,0081 | ST |       |         |
| >>Person Identification Code Sequence         | 0040,1101 | SQ |       |         |
| >>>Code Value                                 | 0008,0100 | SH |       |         |
| >>>Coding Scheme Designator                   | 0008,0102 | SH |       |         |
| >>>Code Meaning                               | 0008,0104 | LO |       |         |
| >>>Context Group Extension Flag               | 0008,010B | CS |       |         |

**Modality Performed Procedure Step SOP Class**

| Attribute Name  | Tag       | VR | Value | Comment   |
|---|-----------|----|-------|---|
| >>>Context Identifier                                 | 0008,010F | CS |       |   |
| >>>Context UID  | 0008,0117 | UI |       |   |
| >>Person's Address                                    | 0040,1102 | ST |       |   |
| >>Person's Telephone Numbers                          | 0040,1103 | LO |       |   |
| >Referenced Image Sequence                            | 0008,1140 | SQ |       |   |
| >>Referenced SOP Class UID                            | 0008,1150 | UI |       |   |
| >>Referenced SOP Instance UID                         | 0008,1155 | UI |       |   |
| >>Container Identifier                                | 0040,0512 | LO |       |   |
| >>Specimen Description Sequence                       | 0040,0560 | SQ |       |   |
| >>>Specimen Identifier                                | 0040,0551 | LO |       |   |
| >>>Specimen UID                                       | 0040,0554 | UI |       |   |
| >Protocol Name  | 0018,1030 | LO |       |   |
| >Series Instance UID                                  | 0020,000E | UI |       |   |
| >Referenced Non-Image Composite SOP Instance Sequence | 0040,0220 | SQ |       | Refers to Dicom Object that were transferred to the PACS. |
| >>Referenced SOP Class UID                            | 0008,1150 | UI |       |   |
| >>Referenced SOP Instance UID                         | 0008,1155 | UI |       |   |
| >Archive Requested                                    | 0040,A494 | CS |       |   |

**Radiation Dose**

|   |           |    |        |  |
|---|-----------|----|--------|--|
| Image and Fluoroscopy Area Dose Product | 0018,115E | DS |        |  |
| Total Time of Fluoroscopy               | 0040,0300 | US |        |  |
| Total Number of Exposures               | 0040,0301 | US |        |  |
| Entrance Dose                           | 0040,0302 | US |        |  |
| Entrance Dose in mGy                    | 0040,8302 | DS |        |  |
| Exposure Dose Sequence                  | 0040,030E | SQ |        |  |
| >KVP                                    | 0018,0060 | DS |        |  |
| >Exposure Time                          | 0018,1150 | IS |        |  |
| >Radiation Mode                         | 0018,115A | CS | PULSED |  |
| >Filter Type                            | 0018,1160 | SH |        |  |
| >Filter Material                        | 0018,7050 | CS |        |  |
| >X-Ray Tube Current In uA               | 0018,8151 | DS |        |  |
| >Comments on Radiation Dose             | 0040,0310 | ST |        |  |

**Billing And Material Management Codes**

|                           |           |    |  |  |
|---------------------------|-----------|----|--|--|
| Film Consumption Sequence | 0040,0321 | SQ |  |  |
| >Medium Type              | 2000,0030 | CS |  |  |
| >Film Size ID             | 2010,0050 | CS |  |  |
| >Number of Films          | 2100,0170 | IS |  |  |

The behavior of the Azurion 2.0 for status codes in an MPPS N-SET response is presented in Table 26. In case of the retransmission attempt each message stored in the persistent queue is sent over a separate association.

**Table 26: MPPS N-SET Response Status Handling Behavior**

| Service Status | Error Code | Further Meaning       | Behavior  |
|----------------|------------|-----------------------|---|
| Success        | 0000       | Successful operation. | The notify status of the related examination is |






|         |                       |                     |  |
|---------|-----------------------|---------------------|--|
|         |                       |                     | updated (set to notified). The examination status is set to COMPLETED and it is logged.  |
| Failure | 0213                  | Resource Limitation | The message contents is made persistent and the message is added to the persistent queue and waits for the next retransmission attempt. The examination status is set to CLOSED.   |
| *       | Any other status code | *                   | If this response status is reported during initial transmission the message contents is made persistent and the message is added to the persistent queue. If this response status is the result of the retransmission attempt related examination is updated to the state as if the transmission succeeded and the message is removed from the persistent queue. Then the notify status of the related examination is updated (set to notified). The response status is logged as a warning. The examination status is set to COMPLETED. |

#### 4.2.1.12.(Real-World) Activity – Image Export

##### 4.2.1.12.1.1. Description and Sequencing of Activities

The operator can select images (and presentation states) and request them to be sent to (pre-configure) multiple destinations. Each request is forwarded to the job queue and processed as individual request to Transfer Images (and Presentation States). If background image and presentation state transfer is configured, the Azurion 2.0 sends the acquired images and presentation states automatically. It can be configured which instances will be automatically marked and the destinations where the instances are automatically sent to. The background image and presentation state transfer is triggered by the image acquisition event and/or by the close examination event in Azurion 2.0. For each request to Transfer Images (and Presentation States) (i.e., transfer job), one association towards the remote system is established. Within the association, for each image or presentation state, a C-STORE request is transmitted. Once the responses are received, the association is closed. A possible sequence of interactions between the Azurion 2.0 and a remote AE with only one C-STORE request is presented in Figure 8.



**Figure 8 : Sequencing of RWA Image Export**

**4.2.1.12.1.2. Proposed Presentation Contexts**

Each time an association is initiated, the Azurion 2.0 proposes two presentation contexts to be used on that association. The presentation context proposed by the Azurion 2.0 for Transfer Images is defined in Table 27. The implementation proposes each SOP Class only once in the abstract syntax specifying all possible transfer syntaxes for that SOP Class. Due to the fact that the SCP has to react with a chosen transfer syntax, as per SOP Class used transfer syntax is forced by the SCP.

**Table 27: Proposed Presentation Contexts for (Real-World) Activity Image Export**

| Presentation Context Table                              |                               |  |                         |      |                      |
|---|-------------------------------|--|-------------------------|------|----------------------|
| Abstract Syntax   |                               | Transfer Syntax  |                         | Role | Extended Negotiation |
| Name  | UID                           | Name   | UID                     |      |                      |
| Secondary Capture Image Storage SOP Class               | 1.2.840.10008.5.1.4.1.1.7     | Implicit VR Little Endian  | 1.2.840.10008.1.2       | SCU  | None                 |
|   |                               | Explicit VR Little Endian  | 1.2.840.10008.1.2.1     |      |                      |
|   |                               | Explicit VR Big Endian   | 1.2.840.10008.1.2.2     |      |                      |
|   |                               | JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14) | 1.2.840.10008.1.2.4.7.0 |      |                      |
| Grayscale Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.1  | Implicit VR Little Endian  | 1.2.840.10008.1.2       | SCU  | None                 |
|   |                               | Explicit VR Little Endian  | 1.2.840.10008.1.2.1     |      |                      |
|   |                               | Explicit VR Big Endian   | 1.2.840.10008.1.2.2     |      |                      |
| X-Ray Angiographic Image Storage SOP Class              | 1.2.840.10008.5.1.4.1.1.12.1  | Implicit VR Little Endian  | 1.2.840.10008.1.2       | SCU  | None                 |
|   |                               | Explicit VR Little Endian  | 1.2.840.10008.1.2.1     |      |                      |
|   |                               | Explicit VR Big Endian   | 1.2.840.10008.1.2.2     |      |                      |
|   |                               | JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14) | 1.2.840.10008.1.2.4.7.0 |      |                      |
| X-Ray Radiation Dose SR                                 | 1.2.840.10008.5.1.4.1.1.88.67 | Implicit VR Little Endian  | 1.2.840.10008.1.2       | SCU  | None                 |
|   |                               | Explicit VR Little Endian  | 1.2.840.10008.1.2.1     |      |                      |
|   |                               | Explicit VR Big Endian   | 1.2.840.10008.1.2.2     |      |                      |

**4.2.1.12.1.3. SOP Specific Conformance for SOP Classes**

Azurion 2.0 can exchange image data in the following formats:

- X-Ray Angiographic images can either be sent with raw pixel data or processed pixel data.
- Standard Secondary Capture (1.2.840.10008.5.1.4.1.1.7)
- Standard Grayscale Softcopy Presentation State SOP Class (1.2.840.10008.5.1.4.1.1.11.1)

Furthermore, the Azurion 2.0 can exchange non-image data in the following format:

- X-Ray Radiation Dose Structured Report (1.2.840.10008.5.1.4.1.1.88.67)

The behavior of Azurion 2.0 for status codes in a C-STORE response is summarized in Table 28.

**Table 28: Storage C-STORE Response Status Handling Behavior**

| Service Status | Error Code | Further Meaning      | Behavior  |
|----------------|------------|----------------------|---|
| Success        | 0000       | Successful operation | The SCP has successfully stored the SOP Instances. If all SOP Instances in a send job have status success then the job is marked as completed. Success is logged. |

| Service Status | Error Code            | Further Meaning                   | Behavior  |
|----------------|-----------------------|-----------------------------------|---|
| Refused        | A700-A7FF             | Out of Resources                  | The association is aborted using A-ABORT and the send job is marked as failed. The failure reason is logged.  |
| Error          | A900-A9FF             | Data Set does not match SOP Class | The association is aborted using A-ABORT and the send job is marked as failed. The failure reason is logged.  |
|                | C000-CFFF             | Cannot Understand                 | The association is aborted using A-ABORT and the send job is marked as failed. The failure reason is logged.  |
| Warning        | B000                  | Coercion of Data Elements         | The SCP has successfully stored the SOP Instances. If all SOP Instances in a send job have status success then the job is marked as completed. The Warning is logged. |
|                | B006                  | Elements discarded                | The SCP has successfully stored the SOP Instances. If all SOP Instances in a send job have status success then the job is marked as completed. The Warning is logged. |
|                | B007                  | Data set does not match SOP class | The SCP has successfully stored the SOP Instances. If all SOP Instances in a send job have status success then the job is marked as completed. The Warning is logged. |
| *              | Any other status code | *                                 | The association is aborted using A-ABORT and the send job is marked as failed. The failure reason is logged.  |

The behavior of Azurion 2.0 during communication failure is summarized in Table 29.

**Table 29: Storage Communication Failure Behavior**

| Exception            | Behavior   |
|----------------------|--|
| Timeout              | The association is aborted using A-ABORT and the send job is marked as failed. The failure reason is logged. |
| Association aborted  | The association is aborted using A-ABORT and the send job is marked as failed. The failure reason is logged. |
| Association rejected | The association is aborted using A-ABORT and the send job is marked as failed. The failure reason is logged. |

**4.2.1.13. (Real-World) Activity – Storage Commitment Push Model as SCU**

**4.2.1.13.1.1. Description and Sequencing of Activities**

Azurion 2.0 Storage Commitment as a SCU service. It accepts a storage commitment notification (N-EVENT-REPORT) from systems that send them.




Figure 9: Data Flow Diagram – Commit Image (asynchronous)

#### 4.2.1.13.1.2. Proposed Presentation Contexts

The presentation contexts are defined in the Table 30.

Table 30: Proposed Presentation Contexts for (Real-World) Activity – Storage Commitment Push Model as SCU

| Presentation Context Table              |                      |                           |                     |      |                      |
|---|----------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                         |                      | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name                                    | UID                  | Name                      | UID                 |      |                      |
| Storage Commitment Push Model SOP Class | 1.2.840.10008.1.20.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None                 |
|   |                      | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |
|   |                      | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |                      |

#### 4.2.1.13.1.3. SOP Specific Conformance for Storage Commitment Push Model SOP Class

This section and sub-section include the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

##### 4.2.1.13.1.3.1. Dataset Specific Conformance for Storage Commitment Push Model SOP Class N-ACTION-SCU

Details regarding the Dataset Specific response behaviour will be reported in this section.

**Table 31: N-ACTION-RQ Dataset Specification.**

| Storage Commitment Push Model SOP Class |           |    |         |
|---|-----------|----|---------|
| Attribute Name                          | Tag       | VR | Comment |
| <b>Sop Common Module</b>                |           |    |         |
| SOP Class UID                           | 0008,0016 | UI |         |
| SOP Instance UID                        | 0008,0018 | UI |         |
| <b>Storage Commitment Module</b>        |           |    |         |
| Transaction UID                         | 0008,1195 | UI |         |
| Referenced SOP Sequence                 | 0008,1199 | SQ |         |
| >Referenced SOP Class UID               | 0008,1150 | UI |         |
| >Referenced SOP Instance UID            | 0008,1155 | UI |         |

The details regarding the response behaviour to status codes are provided in Table 32.

**Table 32: Storage Commitment N-ACTION Response Status Handling Behavior**

| Service Status | Error Code | Further Meaning      | Behavior   |
|----------------|------------|----------------------|--|
| Success        | 0000       | Successful operation | The storage commitment request has been successfully sent. The storage commitment request job is marked as completed. Success is logged. |

**Table 33: Storage Commitment N-EVENT-REPORT Behavior**

| Event Type   | Event Type ID | Behavior   |
|--|---------------|--|
| Storage Commitment Request Successful                | 1             | The Examination is marked as completed and it becomes a candidate for an automatic deletion from the local database if local resources become scarce.                            |
| Storage Commitment Request Complete - Failures Exist | 2             | The failure is reported to the operator by not marking the examination as completed. The operator may re-transfer the image data (which was previously transferred to the PACS). |

#### 4.2.1.13.1.3.2. Dataset Specific Conformance for Storage Commitment Push Model SOP Class N-EVENT-REPORT-SCP

Details regarding the Dataset Specific response behaviour will be reported in this section.

**Table 34: N-EVENT-REPORT-RSP Dataset Specification.**

| Storage Commitment Push Model SOP Class |           |    |         |
|---|-----------|----|---------|
| Attribute Name                          | Tag       | VR | Comment |
| <b>Sop Common Module</b>                |           |    |         |
| SOP Class UID                           | 0008,0016 | UI |         |
| SOP Instance UID                        | 0008,0018 | UI |         |

The details regarding the response behaviour to status codes are provided in Table 35.

**Table 35: Storage Commitment N-EVENT-REPORT Response Status Handling Behavior**

| Service Status | Error Code | Further Meaning      | Behavior   |
|----------------|------------|----------------------|--|
| Success        | 0000       | Successful operation | The storage commitment result has been successfully received. The SCP has successfully stored the SOP Instances. The examination is marked as completed. |

**Table 36: Storage Commitment N-EVENT-REPORT Response Status**

| Event Type   | Event Type ID | Behavior   |
|--|---------------|--|
| Storage Commitment Request Successful                | 1             | The Examination is marked as completed and it becomes a candidate for an automatic deletion from the local database if local resources become scarce.                            |
| Storage Commitment Request Complete - Failures Exist | 2             | The failure is reported to the operator by not marking the examination as completed. The operator may re-transfer the image data (which was previously transferred to the PACS). |

#### 4.2.1.14. (Real-World) Activity – FIND as SCU

##### 4.2.1.14.1.1. Description and Sequencing of Activities

The operator is able to query a (pre-configured) remote database. The ACP AE initiates an association to the selected Remote AE and uses it to send C-FIND requests (and receive the associated find replies). For each query a number of C-FIND requests is established in one association to the peer entity, which is released when all query results are received. An example sequencing of Activities is presented in Figure 10 and Figure 11.





Figure 10: Sequencing of RWA (Patient Root Q/R Information Model)



**Figure 11: Sequencing of RWA (Study Root Q/R Information Model)**

The clinical user may cancel the query to the PACS or Workstation. As a result, the Azurion 2.0 sends a C-FIND Cancel Request to the PACS or Workstation.

**4.2.1.14.1.2. Proposed Presentation Contexts**

Each time an association is initiated, the ACP AE proposes two presentation contexts to be used on that association. The presentation context proposed by the ACP AE for Find Remote Images is defined in Table 37.

**Table 37: Proposed Presentation Contexts for (Real-World) Activity – FIND As SCU**

| Presentation Context Table                         |                             |                           |                     |      |                      |
|--|-----------------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                                    |                             | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name   | UID                         | Name                      | UID                 |      |                      |
| Patient Root QR Information Model - FIND SOP Class | 1.2.840.10008.5.1.4.1.2.1.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None                 |
|  |                             | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |
|  |                             | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |                      |
| Study Root QR Information Model - FIND SOP Class   | 1.2.840.10008.5.1.4.1.2.2.1 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None                 |
|  |                             | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |
|  |                             | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |                      |

**4.2.1.14.1.3. SOP Specific Conformance for SOP Classes**

The Azurion 2.0 provides standard conformance to this SOP class. The Azurion 2.0 AE does not generate queries containing optional keys and it does not generate relational queries.

**4.2.1.14.1.3.1. Dataset Specific Conformance for Patient Root QR Information Model - FIND SOP Class C-FIND-SCU**

In Table 37 the supported query keys for each query level are described. Universal matching is supported by default.

**Table 38: Supported Query Keys for Patient Root Information Model**

| Patient Root Information Model |           |    |   |                        |
|--------------------------------|-----------|----|---|------------------------|
| Attribute Name                 | Tag       | VR | Type Of Matching  | Comment                |
| Query/Retrieve Level           | 0008,0052 | CS | Single Value  | PATIENT, STUDY, SERIES |
| Specific Character Set         | 0008,0005 | CS |   |                        |
| Q/R Patient Level              |           |    |   |                        |
| Patient's Name                 | 0010,0010 | PN | Single value matching or wild card matching or universal matching |                        |
| Patient ID                     | 0010,0020 | LO | Single value matching or wild card matching or universal matching |                        |
| Patient's Birth Date           | 0010,0030 | DA | Single value matching or universal matching                       |                        |
| Patient's Sex                  | 0010,0040 | CS | Universal matching only   |                        |



| Patient Root Information Model     |            |    |   |         |
|------------------------------------|------------|----|---|---------|
| Attribute Name                     | Tag        | VR | Type Of Matching  | Comment |
| <b>Q/R Study level</b>             |            |    |   |         |
| Study Date                         | 0008,0020  | DA | Range matching or universal matching                              |         |
| Study Time                         | 0008,0030  | TM | Universal matching only   |         |
| Accession Number                   | 0008,0050  | SH | Value matching or wild card matching or universal matching        |         |
| Query/Retrieve Level               | 0008,0052  | CS | Single  |         |
| Modalities in Study                | 0008,0061  | CS | Single  |         |
| Study Description                  | 0008,1030) | LO | Universal matching  |         |
| Patient's Name                     | 0010,0010  | PN | Single value matching or wild card matching or universal matching |         |
| Patient ID                         | 0010,0020  | LO | Single value matching or universal matching                       |         |
| Patient's Birth Date               | 0010,0030  | DA | Single value matching or universal matching                       |         |
| Patient's Sex                      | 0010,0040  | CS | Universal matching only   |         |
| Patient Size                       | 0010,1020  | DS | Universal matching  |         |
| Patient Weight                     | 0010,1030  | DS | Universal matching  |         |
| Patient Comment                    | 0010,4000  | LT | Universal matching  |         |
| Study ID                           | 0020,0010  | SH | Universal matching only   |         |
| Study Instance UID                 | 0020,000D  | UI | Universal matching only   |         |
| Number Of Study Related Series     | 0020,1206  | IS | Universal matching  |         |
| <b>Q/R Series level</b>            |            |    |   |         |
| Modality                           | 0008,0060  | CS | Universal matching only   |         |
| Series Description                 | 0008,103E  | LO | Universal matching  |         |
| Number Of Series Related Instances | 0008,1050  | PN | Universal matching only   |         |
| Study Instance UID                 | 0020,000D  | UI | Single value matching only  |         |
| Series Instance UID                | 0020,000E  | UI | Universal matching only   |         |
| Series Number                      | 0020,0011  | IS | Universal matching only   |         |

The behavior of the Azurion 2.0 for status codes in C-FIND response is summarized in Table 38.

**Table 39: Status Response**

| Response Status | Response Code | Further Meaning   | Behavior   |
|-----------------|---------------|---|--|
| Successful      | 0000          | Matching is complete - No final Identifier is supplied. | Stops with processing the C-Find Response(s) from the SCP. All results are displayed to the operator.                    |
| Refused         | A700          | Out of Resources  | Stops with processing the C-Find Response(s) from the SCP. The reason is logged and the failure is reported to the user. |

| Exception Status | Error Code            | Further Meaning  | Behavior  |
|------------------|-----------------------|--|---|
| Failed           | A900                  | Identifier Does Not Match SOP Class  | Stops with processing the C-Find Response(s) from the SCP. The reason is logged and the failure is reported to the user.            |
|                  | Cxxx                  | Unable to process  | Stops with processing the C-Find Response(s) from the SCP. The reason is logged and the failure is reported to the user.            |
| Cancel           | FE00                  | Matching terminated due to Cancel Match request  | Stops with processing the C-Find Response(s) from the SCP. Results already received up to that point are displayed to the operator. |
| Pending          | FF00                  | Matches are continuing - Current Match is supported in the same manner as supplied and any Optional Keys were Required Keys. | Continues with processing of the C-Find Response(s) from the SCP  |
|                  | FF01                  | Matches are continuing - Warning that one or more Optional Keys were not supported for existence for this Identifier.        | Continues with processing of the C-Find Response(s) from the SCP.   |
| *                | Any other status code | *  | The association is aborted using A-ABORT. The reason is logged and the failure is reported to the user.                             |

The behavior of the Azurion 2.0 during communication failure is summarized in Table 39.

**Table 40: DICOM Command Communication Failure Behavior**

| Exception            | Behavior   |
|----------------------|--|
| Timeout              | The query is marked as failed. The association is aborted using A-ABORT. The reason is logged and reported to the user. Azurion 2.0 stops processing the C-FIND Response(s) from the SCP.        |
| Association aborted  | If the association is aborted using A-ABORT, the query is marked as failed. The reason is logged and failure is reported to the user. Stops with processing the C-FIND Response(s) from the SCP. |
| Association rejected | The query is marked as failed. The reason is logged and failure is reported to the user. No C-FIND request performed.  |

**4.2.1.14.1.4. SOP Specific Conformance for Study Root QR Information Model - FIND SOP Class**

The Azurion 2.0 provides standard conformance to this SOP class. The Azurion 2.0 AE does not generate queries containing optional keys and it does not generate relational queries.

**4.2.1.14.1.4.1. Dataset Specific Conformance for Study Root QR Information Model - FIND SOP Class C-FIND-SCU**

In table 40 the supported query keys for each query level are described. Universal matching is supported as default.

**Table 41: Supported Query Keys for Study Root Information Model**

| Study Root Information Model |           |    |                                      |               |
|------------------------------|-----------|----|--------------------------------------|---------------|
| Attribute Name               | Tag       | VR | Type Of Matching                     | Comment       |
| Query/Retrieve Level         | 0008,0052 | CS | Single Value                         | STUDY, SERIES |
| Specific Character Set       | 0008,0005 | CS |                                      |               |
| Q/R Study level              |           |    |                                      |               |
| Study Date                   | 0008,0020 | DA | Range matching or universal matching |               |
| Study Time                   | 0008,0030 | TM | Universal matching only              |               |



| Study Root Information Model       |            |    |   |         |
|------------------------------------|------------|----|---|---------|
| Attribute Name                     | Tag        | VR | Type Of Matching  | Comment |
| Accession Number                   | 0008,0050  | SH | Value matching or wild card matching or universal matching        |         |
| Query/Retrieve Level               | 0008,0052  | CS | Single  | STUDY   |
| Modalities in Study                | 0008,0061  | CS | Single  |         |
| Study Description                  | 0008,1030) | LO | Universal matching  |         |
| Patient's Name                     | 0010,0010  | PN | Single value matching or wild card matching or universal matching |         |
| Patient ID                         | 0010,0020  | LO | Single value matching or universal matching                       |         |
| Patient's Birth Date               | 0010,0030  | DA | Single value matching or universal matching                       |         |
| Patient's Sex                      | 0010,0040  | CS | Universal matching only   |         |
| Patient Size                       | 0010,1020  | DS | Universal matching  |         |
| Patient Weight                     | 0010,1030  | DS | Universal matching  |         |
| Patient Comment                    | 0010,4000  | LT | Universal matching  |         |
| Study ID                           | 0020,0010  | SH | Universal matching only   |         |
| Study Instance UID                 | 0020,000D  | UI | Universal matching only   |         |
| Number Of Study Related Series     | 0020,1206  | IS | Universal matching  |         |
| Q/R Series level                   |            |    |   |         |
| Modality                           | 0008,0060  | CS | Universal matching only   |         |
| Series Description                 | 0008,103E  | LO | Universal matching  |         |
| Number Of Series Related Instances | 0008,1050  | PN | Universal matching only   |         |
| Study Instance UID                 | 0020,000D  | UI | Single value matching only  |         |
| Series Instance UID                | 0020,000E  | UI | Universal matching only   |         |
| Series Number                      | 0020,0011  | IS | Universal matching only   |         |

The behavior of the Azurion 2.0 for status codes in C-FIND response is summarized in Table 41.

**Table 42: Status Response**

| Response Status | Code | Further Meaning   | Behavior   |
|-----------------|------|---|--|
| Success         | 0000 | Matching is complete - No final Identifier is supplied. | Stops with processing the C-Find Response(s) from the SCP. All results are displayed to the operator.                    |
| Refused         | A700 | Out of Resources  | Stops with processing the C-Find Response(s) from the SCP. The reason is logged and the failure is reported to the user. |
| Failed          | A900 | Identifier Does Not Match SOP Class                     | Stops with processing the C-Find Response(s) from the SCP. The reason is logged and the failure is reported to the user. |
|                 | Cxxx | Unable to process                                       | Stops with processing the C-Find Response(s) from the SCP. The reason is logged and the failure is reported to the user. |

| State Status | Code                  | Further Meaning  | Behavior  |
|--------------|-----------------------|--|---|
| Cancel       | FE00                  | Matching terminated due to Cancel Match request  | Stops with processing the C-Find Response(s) from the SCP. Results already received up to that point are displayed to the operator. |
| Pending      | FF00                  | Matches are continuing - Current Match is supported in the same manner as supplied and any Optional Keys were Required Keys. | Continues with processing of the C-Find Response(s) from the SCP  |
|              | FF01                  | Matches are continuing - Warning that one or more Optional Keys were not supported for existence for this Identifier.        | Continues with processing of the C-Find Response(s) from the SCP.   |
| *            | Any other status code | *  | The association is aborted using A-ABORT. The reason is logged and the failure is reported to the user.                             |

The behavior of the Azurion 2.0 during communication failure is summarized in Table 43.

**Table 43: DICOM Command Communication Failure Behavior**

| Exception            | Behavior   |
|----------------------|--|
| Timeout              | The query is marked as failed. The association is aborted using A-ABORT. The reason is logged and reported to the user. The Azurion 2.0 stops processing the C-FIND Response(s) from the SCP.    |
| Association aborted  | If the association is aborted using A-ABORT, the query is marked as failed. The reason is logged and failure is reported to the user. Stops with processing the C-FIND Response(s) from the SCP. |
| Association rejected | The query is marked as failed. The reason is logged and failure is reported to the user. No C-FIND request performed.  |

#### 4.2.1.15. (Real-World) Activity – MOVE as SCU

##### 4.2.1.15.1.1. Description and Sequencing of Activities

The request to move remote images is forwarded to the job queue. For each move job, one association towards the remote system is established, and C-MOVE requests are transmitted. Once the responses are received, the association is closed. An example of sequencing of activities is presented in Figure C-MOVE requests are done on the series level.




Figure 12: Sequencing of RWA

The clinical user may cancel the move operation. As a result, the Azurion 2.0 Sends a C-MOVE Cancel Request to the PACS or Workstation.

#### 4.2.1.15.1.2. Proposed Presentation Contexts

The presentation contexts for MOVE as SCU are defined in table 44.

Table 44: Proposed Presentation Contexts for (Real-World) Activity – MOVE As SCU

| Presentation Context Table                         |                             |                           |                     |      |                      |  |  |
|--|-----------------------------|---------------------------|---------------------|------|----------------------|--|--|
| Abstract Syntax                                    |                             | Transfer Syntax           |                     | Role | Extended Negotiation |  |  |
| Name   | UID                         | Name List                 | UID List            |      |                      |  |  |
| Patient Root QR Information Model - MOVE SOP Class | 1.2.840.10008.5.1.4.1.2.1.2 | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 | SCU  | None                 |  |  |
|  |                             | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |  |  |
|  |                             | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |  |  |
| Study Root QR Information Model - MOVE SOP Class   | 1.2.840.10008.5.1.4.1.2.2.2 | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 | SCU  | None                 |  |  |
|  |                             | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |  |  |
|  |                             | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |  |  |

The implementation proposes each SOP Class only ones in the abstract syntax specifying all possible transfer syntaxes for that SOP Class. Due to the fact that the SCP has to react with a chosen transfer syntax, the per SOP Class used transfer syntax is forced by the SCP

#### 4.2.1.15.1.3. SOP Specific Conformance for SOP Classes

Selecting a query result can retrieve only whole examinations. It is not possible to retrieve information if Patient ID contains the sign "greater than" or "less than" (> or <).

A C-MOVE can be done with the keys presented in Table 44 or Table 47.

#### 4.2.1.15.1.4. SOP Specific Conformance for Patient Root QR Information Model - MOVE SOP Class

The Azurion 2.0 provides standard conformance for the Patient Root QR Information Model - MOVE SOP Class.

##### 4.2.1.15.1.4.1. Dataset Specific Conformance for Patient Root QR Information Model - MOVE SOP Class C-MOVE-SCU

The behavior of the Identifiers for MOVE is summarized in this section.

**Table 45: Identifiers for MOVE Patient Root Information Model as SCU**

| Patient Root Information Model |           |    |                       |
|--------------------------------|-----------|----|-----------------------|
| Attribute Name                 | Tag       | VR | Comment               |
| Query/Retrieve Level           | 0008,0052 | CS | Applied value: Series |
| Q/R Series level               |           |    |                       |
| Patient ID                     | 0010,0020 | PN |                       |
| Study Instance UID             | 0020,000D | UI |                       |
| Series Instance UID            | 0020,000E | UI |                       |

The DICOM C-MOVE Patient Root Information Model Command Status Response Handling is shown in the Table 46.

**Table 46: Status Response for C-MOVE Patient Root Information Model**

| Service Status | Error Code            | Further Meaning  | Behavior   |
|----------------|-----------------------|--|--|
| Success        | 0000                  | Sub-operations Complete – No Failures                    | The move job is marked as completed. The association is released. Success is logged.                             |
| Refused        | A701                  | Out of Resources – Unable to calculate number of matches | The move job is marked as failed. The association is released. The reason is logged and reported to the user.    |
|                | A702                  | Out of Resources – Unable to perform sub operations      | The move job is marked as failed. The association is released. The reason is logged and reported to the user.    |
|                | A801                  | Move Destination Unknown                                 | The move job is marked as failed. The association is released. The reason is logged and reported to the user.    |
| Failed         | A900                  | Identifier Does Not Match SOP Class                      | The move job is marked as failed. The association is released. The reason is logged and reported to the user.    |
|                | Cxxx                  | Unable to process  | The move job is marked as failed. The association is released. The reason is logged and reported to the user.    |
| Cancel         | FE00                  | Sub-operations terminated due to Cancel Indication       | The move job is marked as cancelled. The association is released. The reason is logged and reported to the user. |
| Warning        | B000                  | Sub-operations Complete – One or more Failures           | The move job is marked as failed. The association is released. The reason is logged and reported to the user.    |
| *              | Any other status code | *  | The association is aborted using A-ABORT. The reason is logged and the failure is reported to the user.          |

The possible Communication Failures during a C-MOVE as SCU are shown in the Table 46.

**Table 47: DICOM Command Communication Failure Behavior for C-MOVE Patient Root Information Model.**

| Exception            | Behavior   |
|----------------------|--|
| Timeout              | The query is marked as failed. The association is aborted using A-ABORT. The reason is logged and reported to the user. The Azurion 2.0 stops processing the C-FIND Response(s) from the SCP.    |
| Association aborted  | If the association is aborted using A-ABORT, the query is marked as failed. The reason is logged and failure is reported to the user. Stops with processing the C-FIND Response(s) from the SCP. |
| Association rejected | The query is marked as failed. The reason is logged and failure is reported to the user. No C-FIND request performed.  |

**4.2.1.15.1.5. SOP Specific Conformance for Study Root QR Information Model - MOVE SOP Class**

The Azurion 2.0 provides standard conformance to this SOP class.

**4.2.1.15.1.5.1. Dataset Specific Conformance for Study Root QR Information Model - MOVE SOP Class C-MOVE-SCU**

The identifiers for C-MOVE as SCU are listed in the Table 47.

**Table 48: Identifiers for MOVE Study Root Information Model as SCU**

| Study Root Information Model |           |    |                       |
|------------------------------|-----------|----|-----------------------|
| Attribute Name               | Tag       | VR | Comment               |
| Query/Retrieve Level         | 0008,0052 | CS | Applied value: SERIES |
| Q/R Series level             |           |    |                       |
| Series Instance UID          | 0020,000E | UI |                       |
| Study Instance UID           | 0020,000D | UI |                       |

The DICOM C-MOVE Study Root Information Model Command Status Response Handling is shown in the Table 49.

**Table 49: Status response for Study Root Information Model C-MOVE-SCU.**

| Service Status | Error Code | Further Meaning  | Behavior   |
|----------------|------------|--|--|
| Success        | 0000       | Sub-operations Complete – No Failures                    | The move job is marked as completed. The association is released. Success is logged.                             |
| Refused        | A701       | Out of Resources – Unable to calculate number of matches | The move job is marked as failed. The association is released. The reason is logged and reported to the user.    |
|                | A702       | Out of Resources – Unable to perform sub operations      | The move job is marked as failed. The association is released. The reason is logged and reported to the user.    |
|                | A801       | Move Destination Unknown                                 | The move job is marked as failed. The association is released. The reason is logged and reported to the user.    |
| Failed         | A900       | Identifier Does Not Match SOP Class                      | The move job is marked as failed. The association is released. The reason is logged and reported to the user.    |
|                | Cxxx       | Unable to process  | The move job is marked as failed. The association is released. The reason is logged and reported to the user.    |
| Cancel         | FE00       | Sub-operations terminated due to Cancel Indication       | The move job is marked as cancelled. The association is released. The reason is logged and reported to the user. |
| Warning        | B000       | Sub-operations Complete – One or more Failures           | The move job is marked as failed. The association is released. The reason is logged and reported to the user.    |

| Service Status | Error Code            | Further Meaning | Behavior  |
|----------------|-----------------------|-----------------|---|
| *              | Any other Status code | *               | The association is aborted using A-ABORT. The reason is logged and the failure is reported to the user. |

The possible Communication Failures for C-MOVE-SCU are shown in the table 49.

**Table 50: DICOM Command Communication Failure Behavior for Study Root Information Model C-MOVE-SCU**

| Exception            | Behavior   |
|----------------------|--|
| Timeout              | The query is marked as failed. The association is aborted using A-ABORT. The reason is logged and reported to the user. The Azurion 2.0 stops processing the C-FIND Response(s) from the SCP.    |
| Association aborted  | If the association is aborted using A-ABORT, the query is marked as failed. The reason is logged and failure is reported to the user. Stops with processing the C-FIND Response(s) from the SCP. |
| Association rejected | The query is marked as failed. The reason is logged and failure is reported to the user. No C-FIND request performed.  |

#### 4.2.1.16.(Real-World) Activity – Print Management as SCU

##### 4.2.1.16.1.1. Description and Sequencing of Activities

The operator can select images and request them to be printed on a printer (out of choice list of configured printers). Each request is forwarded to the job queue and processed as individual request to Print Images. The print job consists of data describing the images and graphics to be printed as well as the requested layout and other parameters. One print job on Azurion 2.0 may result in a number of film sessions with the printer equal to the number of printed film sheets. Each film sheet within the print job is internally processed, converted to a STANDARD/1, 1 page and then an association towards the remote Print Server is established and the page image is sent to that Print Server. Once the transmission of the film sheet is completed, the association is closed. A sequence of interactions between the Azurion 2.0 and a remote AE to print one film sheet is presented in Figure 13.






Figure 13: Sequencing of RWA Print Images

The following implementation remarks are important to achieve successful printing:

- Each film session will be in a single association with one or more film boxes and one or more film sheets
- The number of images per Film Box is one. The images to be printed on one film are rendered by the Azurion 2.0 into one logical image. This logical image is very large, depending on the pixel matrix size (pixels per line, lines per image). A rough indication is 20 Mbytes. One should take this into account when selecting the DICOM printer and the printer configuration (e.g. the amount of memory).
- The Azurion 2.0 will request for releasing the association when the print command is given (i.e. the N-ACTION Request). The association is not kept open for receiving N-EVENT-REPORTs of the Printer SOP Class.

**Table 51: The Applied Order of Print Service Elements**

| Service Element of SOP Class                     | Description   |
|--|---|
| N-GET of the Printer SOP Class                   | Purpose is to retrieve printer information.   |
| N-CREATE of the Basic Film Session SOP Class     | Specifies the DICOM Printer about some general presentation parameters, applicable for all films in the Film Session. Applied attributes are: Number of Copies, Print Priority, Medium Type, Film Destination                                 |
| N-CREATE of the Basic Film Box SOP Class         | Specifies the DICOM Printer about some general presentation parameters, applicable for all images in the Film Box. Applied attributes are: Film Orientation, Film Size ID, Magnification Type, Max. Density, Configuration Information, Trim. |
| N-SET of the Basic Grayscale Image Box SOP Class | Images to be printed. Applied attributes are: Polarity  |
| N-ACTION of the Basic Film Box SOP Class         | Triggers the DICOM Printer to print. This actual print action is done at film box level. No attributes are present.   |

**Table 52: Media Services**

| Service Element of SOP Class            | Description  |
|---|--|
| N-EVENT-REPORT of the Printer SOP Class | When N-EVENT-REPORT is received, no printer status polling on a separate connection is executed. |

**4.2.1.16.1.2. Proposed Presentation Contexts**

Each time an association is initiated, the Azurion 2.0 proposes presentation contexts to be used on that association. The presentation contexts proposed by the Azurion 2.0 for Print Images is defined in Table 53.

The implementation proposes the SOP Class only ones in the abstract syntax specifying all possible transfer syntaxes for that SOP Class. Due to the fact that the SCP has to react with a chosen transfer syntax, as per SOP Class used transfer syntax is forced by the SCP.

**Table 53: Proposed Presentation Contexts for (Real-World) Activity – Print Management as SCU**

| Presentation Context Table                      |                        |                           |                     |      |                      |
|---|------------------------|---------------------------|---------------------|------|----------------------|
| Abstract Syntax                                 |                        | Transfer Syntax           |                     | Role | Extended Negotiation |
| Name  | UID                    | Name                      | UID                 |      |                      |
| Basic Grayscale Print Management Meta SOP Class | 1.2.840.10008.5.1.1.9  | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None                 |
|   |                        | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |
|   |                        | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |                      |
| >Basic Film Session SOP Class                   | 1.2.840.10008.5.1.1.1  | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None                 |
|   |                        | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |
|   |                        | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |                      |
| >Basic Film Box SOP Class                       | 1.2.840.10008.5.1.1.2  | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None                 |
|   |                        | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |
|   |                        | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |                      |
| >Basic Grayscale Image Box SOP Class            | 1.2.840.10008.5.1.1.4  | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None                 |
|   |                        | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |
|   |                        | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 |      |                      |
| >Printer SOP Class                              | 1.2.840.10008.5.1.1.16 | Implicit VR Little Endian | 1.2.840.10008.1.2   | SCU  | None                 |
|   |                        | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |



**Presentation Context Table**

| Abstract Syntax                             |                         | Transfer Syntax           |                     | Role | Extended Negotiation |
|---|-------------------------|---------------------------|---------------------|------|----------------------|
| Name  | UID                     | Name                      | UID                 |      |                      |
| Basic Color Print Management Meta SOP Class | 1.2.840.10008.5.1.1.18  | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 | SCU  | None                 |
|   |                         | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |
|   |                         | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |
| >Basic Film Session SOP Class               | 1.2.840.10008.5.1.1.1   | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 | SCU  | None                 |
|   |                         | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |
|   |                         | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |
| >Basic Film Box SOP Class                   | 1.2.840.10008.5.1.1.2   | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 | SCU  | None                 |
|   |                         | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |
|   |                         | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |
| >Basic Color Image Box SOP Class            | 1.2.840.10008.5.1.1.4.1 | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 | SCU  | None                 |
|   |                         | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |
|   |                         | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |
| >Printer SOP Class                          | 1.2.840.10008.5.1.1.16  | Explicit VR Big Endian    | 1.2.840.10008.1.2.2 | SCU  | None                 |
|   |                         | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |                      |
|   |                         | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |      |                      |

**4.2.1.16.1.3. SOP Specific Conformance for Basic Film Session SOP Class for Basic Grayscale Print Meta**

The Azurion 2.0 provides standard conformance to the Basic Grayscale Print Management Meta SOP Class. A description and the applied optional (i.e. non-mandatory attributes as Print SCU) attributes in these Service Elements are specified as well. Note that the Service Elements order is not specified by the DICOM standard. Azurion 2.0 sends the N-DELETE request for the film session. Overlay, annotation (showing the values of some major identifying attributes) and shutter information is processed in the images sent to the printer, all the processing including annotations will be part of the image.

This section and sub-section include the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

**4.2.1.16.1.3.1. Dataset Specific Conformance for Basic Film Session SOP Class for Basic Grayscale Print Meta N-CREATE-SCU**

**Table 54: N-CREATE-RQ Dataset Specification.**

| Attribute Name     | Tag       | VR | Value   | Presence of Value | Source | Comment |
|--------------------|-----------|----|---|-------------------|--------|---------|
| Number of Copies   | 2000,0010 | IS | Between 1 and 100                                     | ALWAYS            | USER   |         |
| Print Priority     | 2000,0020 | CS | MED   | ALWAYS            | AUTO   |         |
| Medium Type        | 2000,0030 | CS | PAPER, BLUE FILM, CLEAR FILM                          | ALWAYS            | USER   |         |
| Film Destination   | 2000,0040 | CS | MAGAZINE, PROCESSOR                                   | ALWAYS            | AUTO   |         |
| Film Session Label | 2000,0050 | LO | Human readable label that identifies the film session | ANAP              | AUTO   |         |

The details regarding the response behaviour to status codes are provided in Table 65.

**Table 55: Status Response**

| Service Status | Error Code | Further Meaning      | Behavior  |
|----------------|------------|----------------------|---|
| Success        | 0000       | Successful operation | The print job continues.  |
| Warning        | 0116       | Any warning          | The print job continues and the warning is logged.  |
| Error          | 0105       | Any error            | The association is aborted using A-ABORT.<br>The print job will keep resubmitting the failed sheets until the error is solved or the retry timeout is exceeded. |
|                | 0106       | Any error            | The association is aborted using A-ABORT. The print job will keep resubmitting the failed sheets until the error is solved or the retry timeout is exceeded.    |

**4.2.1.16.1.4. SOP Specific Conformance for Basic Film Box SOP Class for Basic Grayscale Print Meta**

**4.2.1.16.1.4.1. Dataset Specific Conformance for Basic Film Box SOP Class for Basic Grayscale Print Meta N-CREATE-SCU**

The behavior of the Azurion 2.0 for status codes in an N-CREATE response is summarized in Table 61.

**Table 56: N-CREATE-RQ Dataset Specification.**

| Attribute Name            | Tag       | VR | Value   | Presence of Value | Source | Comment |
|---------------------------|-----------|----|---|-------------------|--------|---------|
| Image Display Format      | 2010,0010 | ST | STANDARD\1,1  |                   |        |         |
| Workstation Format        | 2010,0010 | ST | STANDARD\C,R, CUSTOMi   | ANAP              | AUTO   |         |
| Film Orientation          | 2010,0040 | CS | PORTRAIT, LANDSCAPE   | ANAP              | USER   |         |
| Film Size ID              | 2010,0050 | CS | DICOM specifies a number of Defined Terms; more values are possible and is print configuration dependent.   | ANAP              | USER   |         |
| Border Density            | 2010,0100 | CS | BLACK   |                   |        |         |
| Magnification Type        | 2010,0060 | CS | Normally sent out, however sometimes send out empty Because some DICOM printers are not able to handle (Value NONE for) this attribute. Applied value(s): NONE                      | ANAP              | AUTO   |         |
| Max Density               | 2010,0130 | US | Maximum density of the images on the film, expressed in hundredths of OD. If Max Density is higher than maximum printer density than Max Density is set to maximum printer density. | ANAP              | AUTO   |         |
| Trim                      | 2010,0140 | CS | NO  | ANAP              | AUTO   |         |
| Configuration Information | 2010,0150 | ST | Contains a vendor specific Lookup-table (LUT); should be applied by the   | ANAP              | AUTO   |         |



| Attribute Name                       | Tag       | VR | Value                                 | Presence of Value | Source | Comment |
|--------------------------------------|-----------|----|---------------------------------------|-------------------|--------|---------|
|                                      |           |    | DICOM printer if LUT data is present. |                   |        |         |
| Illumination                         | 2010,015E | US |                                       | ANAP              |        |         |
| Reflected Ambient Light              | 2010,0160 | US |                                       | ANAP              |        |         |
| Referenced Film Session Sequence     | 2010,0500 | SQ |                                       | ANAP              | AUTO   |         |
| >Referenced SOP Class UID            | 0008,1150 | UI |                                       | ALWAYS            |        |         |
| >Referenced SOP Instance UID         | 0008,1155 | UI |                                       | ALWAYS            |        |         |
| Referenced Presentation LUT Sequence | 2050,0500 | SQ |                                       | ANAP              | AUTO   |         |

The details regarding the response behaviour to status codes are provided in Table 62.

**Table 57: Status Response**

| Service Status | Error Code | Further Meaning      | Behavior   |
|----------------|------------|----------------------|--|
| Success        | 0000       | Successful operation | The print job continues.   |
| Warning        | 0116       | Any warning          | The print job continues and the warning is logged.   |
| Error          | 0105       | Any error            | The association is aborted using A-ABORT. The print job will keep resubmitting the failed sheets until the error is solved or the retry timeout is exceeded. |
|                | 0106       | Any error            | The association is aborted using A-ABORT. The print job will keep resubmitting the failed sheets until the error is solved or the retry timeout is exceeded. |

**4.2.1.16.1.5. SOP Specific Conformance for Basic Grayscale Image Box SOP Class for Basic Grayscale Print Meta**

**4.2.1.16.1.5.1. Dataset Specific Conformance for Basic Grayscale Image Box SOP Class for Basic Grayscale Print Meta N-SET-SCU**

This application entity supports the attributes described in the table 63.

**Table 58: N-SET-RQ Dataset Specification.**

| Attribute Name                 | Tag       | VR | Value                             | Presence of Value | Source | Comment |
|--------------------------------|-----------|----|-----------------------------------|-------------------|--------|---------|
| SOP Class UID                  | 0008,0016 | UI |                                   | ALWAYS            |        |         |
| SOP Instance UID               | 0008,0018 | UI |                                   | ALWAYS            |        |         |
| Image Box Position             | 2020,0010 | US | 1                                 | ANAP              | AUTO   |         |
| Polarity                       | 2020,0020 | CS | NORMAL                            | ANAP              | AUTO   |         |
| Basic Grayscale Image Sequence | 2020,0110 | SQ |                                   | ANAP              |        |         |
| >Samples Per Pixel             | 0028,0002 | US | 1                                 | ANAP              | AUTO   |         |
| >Photometric Interpretation    | 0028,0004 | CS | MONOCHROME2                       | ANAP              | AUTO   |         |
| >Rows                          | 0028,0010 | US | Depending on the selected printer | ANAP              | AUTO   |         |

|                                |           |       |   |        |      |  |
|--------------------------------|-----------|-------|---|--------|------|--|
|                                |           |       | type and film size.                                   |        |      |  |
| >Columns                       | 0028,0011 | US    | Depending on the selected printer type and film size. | ANAP   | AUTO |  |
| >Pixel Aspect Ratio            | 0028,0034 | IS    |   | ANAP   |      |  |
| >Bits Allocated                | 0028,0100 | US    | 8   | ANAP   | AUTO |  |
| >Bits Stored                   | 0028,0101 | US    | 12, 8   | ANAP   | AUTO |  |
| >High Bit                      | 0028,0102 | US    | 11, 7   | ANAP   | AUTO |  |
| >Pixel Representation          | 0028,0103 | US    | 0X0000  | ANAP   | AUTO |  |
| >Pixel Data                    | 7FE0,0010 | OW/OB |   | ANAP   | AUTO |  |
| Image Box Position             | 2020,0010 | US    |   | ALWAYS | AUTO |  |
| Polarity                       | 2020,0020 | CS    |   | ALWAYS | AUTO |  |
| Basic Grayscale Image Sequence | 2020,0110 | SQ    |   | ALWAYS | AUTO |  |

The details regarding the response behaviour to status codes are provided in Table 64.

**Table 59: Status Response**

| Service Status | Error Code | Further Meaning      | Behavior   |
|----------------|------------|----------------------|--|
| Success        | 0000       | Successful operation | The print job continues.   |
| Warning        | 0107       | Any warning          | The print job continues and the warning is logged.   |
| Error          | 0105       | Any error            | The association is aborted using A-ABORT. The print job will keep resubmitting the failed sheets until the error is solved or the retry timeout is exceeded. |
|                | 0106       | Any error            | The association is aborted using A-ABORT. The print job will keep resubmitting the failed sheets until the error is solved or the retry timeout is exceeded. |

#### 4.2.1.16.1.6. SOP Specific Conformance for Basic Grayscale Image Box SOP Class for Basic Grayscale Print Meta

##### 4.2.1.16.1.6.1. Dataset Specific Conformance for Printer SOP Class for Basic Grayscale Print Meta N-GET-SCU

This application entity supports the attributes described in the Table 65.

**Table 60: N-GET-RQ Dataset Specification.**

| Attribute Name        | Tag       | VR | Value | Presence of Value | Source | Comment |
|-----------------------|-----------|----|-------|-------------------|--------|---------|
| SOP Class UID         | 0008,0016 | UI |       | ALWAYS            | AUTO   |         |
| SOP Instance UID      | 0008,0018 | UI |       | ALWAYS            | AUTO   |         |
| Printer Status        | 2110,0010 | CS |       | ANAP              | AUTO   |         |
| Printer Status Info   | 2110,0020 | CS |       | ANAP              | AUTO   |         |
| Print Priority        | 2000,0020 | CS |       | ALWAYS            | AUTO   |         |
| Execution Status      | 2100,0020 | CS |       | ALWAYS            | AUTO   |         |
| Execution Status Info | 2100,0030 | CS |       | ALWAYS            | AUTO   |         |

##### 4.2.1.16.1.6.2. Dataset Specific Conformance for Printer SOP Class for Basic Grayscale Print Meta N-EVENT-REPORT-SCP

This application entity supports the attributes described in the next table.

**Table 61: N-EVENT-REPORT-RSP Dataset Specification.**

| Attribute Name      | Tag       | VR | Value | Presence of Value | Source | Comment |
|---------------------|-----------|----|-------|-------------------|--------|---------|
| SOP Class UID       | 0008,0016 | UI |       | ALWAYS            | AUTO   |         |
| SOP Instance UID    | 0008,0018 | UI |       | ALWAYS            | AUTO   |         |
| Printer Status      | 2110,0010 | CS |       | ANAP              | AUTO   |         |
| Printer Status Info | 2110,0020 | CS |       | ANAP              | AUTO   |         |

**4.2.1.16.1.7. SOP Specific Conformance for Basic Color Print Management Meta SOP Class**

The Azurion 2.0 provides standard conformance to the Basic Color Print Management Meta SOP Class. A description and the applied optional (i.e. non-mandatory attributes as Print SCU) attributes in these Service Elements are specified as well. Note that the Service Elements order is not specified by the DICOM standard. Azurion 2.0 sends the N-DELETE request for the film session. Overlay, annotation (showing the values of some major identifying attributes) and shutter information is processed in the images sent to the printer, all the processing including annotations will be part of the image.

**4.2.1.16.1.7.1. SOP Specific Conformance for Basic Film Session SOP Class for Basic Color Print Meta**

**4.2.1.16.1.7.1.1. Dataset Specific Conformance for Basic Film Session SOP Class for Basic Color Print Meta N-CREATE-SCU**

This application entity supports the attributes described in the Table 67.

**Table 62: N-CREATE-RQ Dataset Specification.**

| Attribute Name     | Tag       | VR | Value   | Presence of Value | Source | Comment |
|--------------------|-----------|----|---|-------------------|--------|---------|
| Number of Copies   | 2000,0010 | IS | Between 1 and 100                                     | ALWAYS            | USER   |         |
| Print Priority     | 2000,0020 | CS | MED   | ALWAYS            | AUTO   |         |
| Medium Type        | 2000,0030 | CS | PAPER, BLUE FILM, CLEAR FILM                          | ALWAYS            | USER   |         |
| Film Destination   | 2000,0040 | CS | MAGAZINE, PROCESSOR                                   | ALWAYS            | AUTO   |         |
| Film Session Label | 2000,0050 | LO | Human readable label that identifies the film session | ANAP              | AUTO   |         |

The details regarding the response behaviour to status codes are provided in Table 68.

**Table 63: Status Response**

| Service Status | Error Code | Further Meaning      | Behavior   |
|----------------|------------|----------------------|--|
| Success        | 0000       | Successful operation | The print job continues.   |
| Warning        | 0107       | Any warning          | The print job continues and the warning is logged.   |
| Error          | 0105       | Any error            | The association is aborted using A-ABORT. The print job will keep resubmitting the failed sheets until the error is solved or the retry timeout is exceeded. |
|                | 0106       | Any error            | The association is aborted using A-ABORT. The print job will keep resubmitting the failed sheets until the error is solved or the retry timeout is exceeded. |

4.2.1.16.1.7.2. SOP Specific Conformance for Basic Film Box SOP Class for Basic Color Print Meta

4.2.1.16.1.7.2.1. Dataset Specific Conformance for Basic Film Box SOP Class for Basic Color Print Meta N-CREATE-SCU

The behavior of the Azurion 2.0 for status codes in an N-CREATE response is summarized in Table 69.

Table 64: N-CREATE-RQ Dataset Specification.

| Attribute Name                   | Tag       | VR | Value   | Presence of Value | Source | Comment |
|----------------------------------|-----------|----|---|-------------------|--------|---------|
| SOP Class UID                    | 0008,0016 | UI |   | ALWAYS            | AUTO   |         |
| SOP Instance UID                 | 0008,0018 | UI |   | ALWAYS            | AUTO   |         |
| Workstation Format               | 2010,0010 | ST | STANDARD\C, R, CUSTOMi  | ALWAYS            | AUTO   |         |
| Film Orientation                 | 2010,0040 | CS | PORTRAIT, LANDSCAPE   | ALWAYS            | USER   |         |
| Film Size ID                     | 2010,0050 | CS | DICOM specifies a number of Defined Terms; more values are possible and is print configuration dependent.   | ALWAYS            | USER   |         |
| Magnification Type               | 2010,0060 | CS | Normally sent out, however sometimes send out empty Because some DICOM printers are not able to handle (Value NONE for) this attribute. Applied value(s): NONE                      | ALWAYS            | AUTO   |         |
| Max Density                      | 2010,0130 | US | Maximum density of the images on the film, expressed in hundredths of OD. If Max Density is higher than maximum printer density than Max Density is set to maximum printer density. | ALWAYS            | AUTO   |         |
| Trim                             | 2010,0140 | CS | NO  | ALWAYS            | AUTO   |         |
| Configuration Information        | 2010,0150 | ST | Contains a vendor specific Lookup-table (LUT); should be applied by the DICOM printer if LUT data is present.   | ALWAYS            | AUTO   |         |
| Reflected Ambient Light          | 2010,0160 | US |   | VNAP              |        |         |
| Referenced Film Session Sequence | 2010,0500 | SQ | Parent Film Session   | ALWAYS            | AUTO   |         |
| >Referenced SOP Class UID        | 0008,1150 | UI | 1.2.840.10008.5.1.1.1   | ALWAYS            | AUTO   |         |
| >Referenced SOP Instance UID     | 0008,1155 | UI | UID of Parent Film Session  | ALWAYS            | AUTO   |         |
| Referenced Image                 | 2010,0510 | SQ |   | ALWAYS            |        |         |





| Attribute Name                           | Tag       | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|--------|---------|
| Box Sequence                             |           |    |       |                   |        |         |
| Referenced Basic Annotation Box Sequence | 2010,0520 | SQ |       | VNAP              |        |         |
| Referenced Presentation LUT Sequence     | 2050,0500 | SQ |       | VNAP              |        |         |

The details regarding the response behaviour to status codes are provided in Table 70.

**Table 65: Status Response**

| Service Status | Error Code | Further Meaning      | Behavior   |
|----------------|------------|----------------------|--|
| Success        | 0000       | Successful operation | The print job continues.   |
| Warning        | 0107       | Any warning          | The print job continues and the warning is logged.   |
| Error          | 0105       | Any error            | The association is aborted using A-ABORT. The print job will keep resubmitting the failed sheets until the error is solved or the retry timeout is exceeded. |
|                | 0106       | Any error            | The association is aborted using A-ABORT. The print job will keep resubmitting the failed sheets until the error is solved or the retry timeout is exceeded. |

**4.2.1.16.1.7.3. Dataset Specific Conformance for Basic Film Box SOP Class for Basic Color Print Meta N-ACTION-SCU**

The behavior of the Azurion 2.0 for status codes in an N-ACTION response is summarized in Table 71.

**Table 66: N-ACTION-RQ Dataset Specification.**

| Attribute Name   | Tag       | VR | Value | Presence of Value | Source | Comment |
|------------------|-----------|----|-------|-------------------|--------|---------|
| SOP Class UID    | 0008,0016 | UI |       | ALWAYS            | AUTO   |         |
| SOP Instance UID | 0008,0018 | UI |       | ALWAYS            | AUTO   |         |

**4.2.1.16.1.7.4. SOP Specific Conformance for Basic Color Image Box SOP Class for Basic Color Print Meta**

**4.2.1.16.1.7.4.1. Dataset Specific Conformance for Basic Color Image Box SOP Class for Basic Color Print Meta N-SET-SCU**

This application entity supports the attributes described in the Table 72.

**Table 67: N-SET-RQ Dataset Specification.**

| Attribute Name     | Tag       | VR | Value | Presence of Value | Source | Comment |
|--------------------|-----------|----|-------|-------------------|--------|---------|
| SOP Class UID      | 0008,0016 | UI |       | ALWAYS            | AUTO   |         |
| SOP Instance UID   | 0008,0018 | UI |       | ALWAYS            | AUTO   |         |
| Image Box Position | 2020,0010 | US |       | ANAP              | AUTO   |         |
| Polarity           | 2020,0020 | CS |       | ANAP              | AUTO   |         |
| Samples Per Pixel  | 0028,0002 | US |       | ANAP              | AUTO   |         |

Template ID: HSDPCPF-T-30001.09



| Attribute Name             | Tag       | VR    | Value | Presence of Value | Source | Comment |
|----------------------------|-----------|-------|-------|-------------------|--------|---------|
| Photometric Interpretation | 0028,0004 | CS    |       | ANAP              | AUTO   |         |
| Rows                       | 0028,0010 | US    |       | ANAP              | AUTO   |         |
| Columns                    | 0028,0011 | US    |       | ANAP              | AUTO   |         |
| Pixel Aspect Ratio         | 0028,0034 | IS    |       | ANAP              | AUTO   |         |
| Bits Allocated             | 0028,0100 | US    |       | ANAP              | AUTO   |         |
| Bits Stored                | 0028,0101 | US    |       | ANAP              | AUTO   |         |
| High Bit                   | 0028,0102 | US    |       | ANAP              | AUTO   |         |
| Pixel Representation       | 0028,0103 | US    |       | ANAP              | AUTO   |         |
| Pixel Data                 | 7FE0,0010 | OW/OB |       | ANAP              | AUTO   |         |

The details regarding the response behaviour to status codes are provided in Table 73.

**Table 68: Status Response**

| Service Status | Error Code | Further Meaning      | Behavior   |
|----------------|------------|----------------------|--|
| Success        | 0000       | Successful operation | The print job continues.   |
| Warning        | 0107       | Any warning          | The print job continues and the warning is logged.   |
| Error          | 0105       | Any error            | The association is aborted using A-ABORT. The print job will keep resubmitting the failed sheets until the error is solved or the retry timeout is exceeded. |
|                | 0106       | Any error            | The association is aborted using A-ABORT. The print job will keep resubmitting the failed sheets until the error is solved or the retry timeout is exceeded. |

#### 4.2.1.16.1.7.5. SOP Specific Conformance for Basic Color Image Box SOP Class for Basic Color Print Meta

##### 4.2.1.16.1.7.5.1. Dataset Specific Conformance for Printer SOP Class for Basic Color Print Meta N-EVENT-REPORT-SCP

This application entity supports the attributes described in the Table 74.

**Table 69: N-EVENT-REPORT-RSP Dataset Specification.**

| Attribute Name   | Tag       | VR | Value | Presence of Value | Source | Comment |
|------------------|-----------|----|-------|-------------------|--------|---------|
| SOP Class UID    | 0008,0016 | UI |       | ALWAYS            | AUTO   |         |
| SOP Instance UID | 0008,0018 | UI |       | ALWAYS            | AUTO   |         |

#### 4.2.1.17. Association Acceptance Policy

The Application Entity will respond to a received Association rejection as shown in the Table 50.

**Table 70: Association Rejection response**

| Result               | Source                  | Reason/Diagnosis                         | Behaviour |
|----------------------|-------------------------|--|-----------|
| 1-rejected-permanent | 1-DICOM UL service-user | 1-no-reason-given                        | Log entry |
|                      |                         | 2-application-context-name-not-supported | Log entry |



| Result               | Source   | Reason/Diagnosis                         | Behaviour |
|----------------------|--|--|-----------|
|                      |  | 3-calling-AE-title-not-recognized        | Log entry |
|                      |  | 7-called- AE-title-not-recognized        | Log entry |
|                      | 2-DICOM UL service-provider<br>(ACSE related function)         | 1-no-reason-given                        | Log entry |
|                      |  | 2-no-reason-given                        | Log entry |
|                      | 3-DICOM UL service-provider<br>(Presentation related function) | 1-temporary-congestion                   | Log entry |
|                      |  | 2-local-limit-exceeded                   | Log entry |
| 2-rejected-transient | 1-DICOM UL service-user  | 1-no-reason-given                        | Log entry |
|                      |  | 2-application-context-name-not-supported | Log entry |
|                      |  | 3-calling-AE-title-not-recognized        | Log entry |
|                      |  | 7-called- AE-title-not-recognized        | Log entry |
|                      | 2-DICOM UL service-provider<br>(ACSE related function)         | 1-no-reason-given                        | Log entry |
|                      |  | 2-no-reason-given                        | Log entry |
|                      | 3-DICOM UL service-provider<br>(Presentation related function) | 1-temporary-congestion                   | Log entry |
|                      |  | 2-local-limit-exceeded                   | Log entry |

The behavior of the Azurion 2.0 during Abort Handling is summarized in Table 51.


**Table 71: Association Abort Handling**

| Source                        | Reason/Diagnosis                | Behaviour |
|-------------------------------|---------------------------------|-----------|
| 0 - DICOM UL service-user     | 0 - reason-not-specified        | Log entry |
| 2 - DICOM UL service-provider | 0 - reason-not-specified        | Log entry |
|                               | 1 - unrecognized-PDU            | Log entry |
|                               | 2 - unexpected-PDU              | Log entry |
|                               | 4 - unrecognized-PDU-parameter  | Log entry |
|                               | 5 - unexpected-PDU-parameter    | Log entry |
|                               | 6 - invalid-PDU-parameter-value | Log entry |

#### 4.2.1.18.(Real-World) Activity – Verification as SCP

##### 4.2.1.18.1.1. Description and Sequencing of Activities

A remote system requests verification from Azurion 2.0 using the C-ECHO command.



**Figure 14: (Real World) Activity - Verification as SCP Accepted Presentation Contexts**

The presentation contexts are defined in next table.

**Table 72: Acceptable Presentation Contexts for (Real-World) Activity – Verification as SCP**

| Presentation Context Table |                   |                           |                     |      |             |
|----------------------------|-------------------|---------------------------|---------------------|------|-------------|
| Abstract Syntax            |                   | Transfer Syntax           |                     | Role | Negotiation |
| Name                       | UID               | Name List                 | UID List            |      |             |
| Verification SOP Class     | 1.2.840.10008.1.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCP  | None        |
|                            |                   | Implicit VR Little Endian | 1.2.840.10008.1.2   |      |             |

**4.2.1.18.1.2. SOP Specific Conformance for Verification SOP Class**

Azurion 2.0 (C-ECHO SCP) provides standard conformance to the DICOM V3.0 verification SOP Class.

**4.2.1.18.1.2.1. Dataset Specific Conformance for Verification C-ECHO SCP**

This part of the section includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

**Table 73: Status Response**

| Service Status | Error Code | Further Meaning | Behavior                                |
|----------------|------------|-----------------|---|
| Success        | 0000       | Success         | C-ECHO command was successful received. |

**4.2.1.19.(Real-World) Activity – Image Import**

**4.2.1.19.1.1. Description and Sequencing of Activities**

The real world activity associated with the C-STORE operation is the storage of the image in the memory of the system upon which Azurion 2.0 is running in order to make it available for immediate processing by applications. Azurion 2.0 will issue a failure status if it is unable to store the image in the memory.




Figure 15: (Real World) Activity - Image Import.

#### 4.2.1.19.1.2. Accepted Presentation Contexts

The presentation contexts are defined in the Table 54.

Table 74: Acceptable Presentation Contexts for (Real-World) Activity – Image Import

| Presentation Context Table                              |                              |  |                        |      |                      |
|---|------------------------------|--|------------------------|------|----------------------|
| Abstract Syntax   |                              | Transfer Syntax  |                        | Role | Extended Negotiation |
| Name  | UID                          | Name List  | UID List               |      |                      |
| Secondary Capture Image Storage SOP Class               | 1.2.840.10008.5.1.4.1.1.7    | Implicit VR Little Endian  | 1.2.840.10008.1.2      | SCP  | None                 |
|   |                              | Explicit VR Little Endian  | 1.2.840.10008.1.2.1    |      |                      |
|   |                              | Explicit VR Big Endian   | 1.2.840.10008.1.2.2    |      |                      |
|   |                              | JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |
| Grayscale Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.1 | Implicit VR Little Endian  | 1.2.840.10008.1.2      | SCP  | None                 |
|   |                              | Explicit VR Little Endian  | 1.2.840.10008.1.2.1    |      |                      |
|   |                              | Explicit VR Big Endian   | 1.2.840.10008.1.2.2    |      |                      |
| X-Ray Angiographic Image Storage SOP Class              | 1.2.840.10008.5.1.4.1.1.12.1 | Implicit VR Little Endian  | 1.2.840.10008.1.2      | SCP  | None                 |
|   |                              | Explicit VR Little Endian  | 1.2.840.10008.1.2.1    |      |                      |
|   |                              | Explicit VR Big Endian   | 1.2.840.10008.1.2.2    |      |                      |
|   |                              | JPEG Lossless, Non-Hierarchical, First-Order Prediction (Process 14) | 1.2.840.10008.1.2.4.70 |      |                      |

#### 4.2.1.19.1.3. SOP Specific Conformance for Storage SOP Classes

This section and sub-section includes the manufacturer SOP and Dataset specific information as well the status codes and their corresponding behavior.

Azurion 2.0 provides standard conformance to the DICOM V3.0 Storage Service Class as a SCP. Azurion 2.0 conforms to the SOPs of the Storage Service Class at Level 2 (Full). In case of a successful C-STORE, the stored image may be accessed by the processing applications.

Template ID: HSDPCPF-T-30001.09

#### 4.2.1.19.1.3.1. Dataset Specific Conformance for C-STORE-RSP

This includes the dataset specific behavior, i.e. error codes, error and exception handling, time-outs, etc.

**Table 75: Status Response**

| Service Status | Error Code | Further Meaning   | Behavior                                |
|----------------|------------|-------------------|---|
| Success        | 0000       | Successful stored | Whenever the store operation succeeded. |
| Failure        | Cxxx       | Failed            | Whenever the store operation failed.    |

### 4.3. Network Interfaces

#### 4.3.1. Physical Network Interfaces

The System provides only DICOM V3.0 TCP/IP Network Communication Support as defined in PS 3.8 of the standard. TCP/IP is the only protocol stack supported.

Supported physical medium include:

- IEEE 802.3-1995, 10BASE-T
- IEEE 802.3-1995, 100BASE-TX (Fast Ethernet)
- IEEE 802.3, 1000BASE-X (Fiber Optic Gigabit Ethernet).

The TCP/IP Stack as supported by the underlying Operating System. The API is the WinSock 2 interface as supported by the underlying Operating System.

#### 4.3.2. Additional Protocols

##### 4.3.2.1. Basic TLS Secure Transport Connection Profile

Azurion 2.0 conforms to the Basic TLS Secure Transport Connection Profile

##### 4.3.2.2. Basic Time Synchronization Profile

Azurion 2.0 conforms to the Basic Time Synchronization Profile as an NTP Client implementing the Maintain Time transaction.

##### 4.3.2.3. Basic Application Level Confidentiality Profile

See Section 7.2.3.

#### 4.3.3. IPv4 and IPv6 Support

Azurion 2.0 Supports IPv4 and IPv6.

### 4.4. Configuration

Any implementation's DICOM conformance may be dependent upon configuration, which takes place at the time of installation. Issues concerning configuration are addressed in this section.

#### 4.4.1. AE Title/Presentation Address Mapping

An important installation issue is the translation from AE title to presentation address. How this is to be performed is described here.

##### 4.4.1.1. Local AE Titles

The local AE title mapping and configuration are specified as:

**Table 76: AE Title Configuration**

| Application Entity | Default AE Title | Default TCP/IP Port |
|--------------------|------------------|---------------------|
| IENGINE_SCU        | IENGINE_SCU      | 29536               |

#### 4.4.1.2. Remote AE Title/Presentation Address Mapping

All relevant remote applications that should be able to setup a DICOM association towards Azurion 2.0 and that should be able to accept a DICOM association from Azurion 2.0 must be configured during Azurion 2.0 configuration time.

#### 4.4.2. Parameters

The specification of important operational parameters, their default value and range (if configurable) are specified here.

**Table 77: Configuration Parameters**

| Parameter            | Configurable | Default Value |
|----------------------|--------------|---------------|
| AE Title             | Yes          | IENGINE_SCU   |
| Port number          | Yes          | 29536         |
| IP host name/address | Yes          | -             |

## 5. Media Interchange

### 5.1. Implementation model

The implementation model identifies the DICOM Application Entities for Media in specific implementation and relates the Application Entities to Real-World Activities.

#### 5.1.1. Application Data Flow Diagram

As part of the implementation model, an application data flow diagram is included. The next Figure shows the media interchange application data flow as a functional overview of the Media AE for DICOM CD and DVD.




Figure 16: Application Data Flow Diagram

The Media AE acts as a FSC for CD-R and DVD, when writing the selected images in a patient folder onto the medium.

#### 5.1.2. Functional Definitions of AE's

The Azurion 2.0 implements the following functions for DICOM media.

- Write a DICOM file-set onto the medium.
- Create a DICOMDIR file.

#### 5.1.3. Sequencing of Real World Activities

Not applicable.

## 5.2. AE Specifications

This section in the DICOM Conformance Statement specifies a set of Media Application Entities.

### 5.2.1. Media Storage Media - Specification

This section contains general policies that apply to all of the Application Entities described in subsequent section. The Azurion 2.0 provides standard conformance to the DICOM interchange option of the media storage service class, and follows the specifications as defined in the DICOM standard – Media Storage and File Format for Data Interchange (PS 3.10) and Media Storage Application Profiles (PS 3.11).

The Azurion 2.0 supports multi-patient and multi-session for CD-R media (both reading and writing). For one or more Application Profiles, the following table shows the Real-World Activities and the roles of each of these Real-World Activities.

Table 78: AE Media Storage related Application Profiles, Real-World Activities and Roles

| Supported Application Profile                | Identifier      | Real-World Activities | Roles           |
|--|-----------------|-----------------------|-----------------|
| General Purpose CD-R Interchange             | STD-GEN-CD      | Create File-set       | FSC/FSR         |
| General Purpose Interchange on DVD-RAM Media | STD-GEN-DVD-RAM | Create File-set       | FSC/FSR         |
| General Purpose USB Media Interchange        | STD-GEN-USB     | Create File-set       | FSC/FSR/FS<br>U |



### 5.2.1.1. File Meta Information for the Media AE

Table 79: Implementation Identifying Information

| Property                     | Value                    |
|------------------------------|--------------------------|
| Implementation Class UID:    | 1.3.46.670589.7.29.2.0.1 |
| Implementation Version Name: | Azurion                  |

### 5.2.1.2. Real-World Activities

The AE specification contains a description of the Real-World Activities, which invoke the particular AE.

#### 5.2.1.3. RWA - Read File-set

Azurion 2.0 can read SC, XA and GSPS objects.

##### 5.2.1.3.1.1. Media Storage Application Profile

Refer to the table in section 5.2.1, Table 77.

##### 5.2.1.3.1.1.1. Options

Not applicable.

#### 5.2.1.4. RWA - Create File-set

This Media Application Entity has a File-set Creator functionality which is described here.

#### Write Images

The Azurion 2.0 acts as an FSC when writing DICOM objects onto DICOM media. The Azurion 2.0 can also store private attributes. When the Azurion 2.0 has to write objects to DICOM media, it can encounter the following situation.

The objects were previously received via C-STORE operations. Some attributes in the received images have a zero-length value (Type 2 attributes). However, the Application Profile specifies some of these attributes as type 1: they must have a value. In such cases the Azurion 2.0 supplies a value for the following attributes (if necessary):

- Patient ID;
- Study ID;
- Series Number;
- Instance number;
- Study Date;
- Study Time.

The mechanism of generating a value for Patient ID is to create a new value (i.e. Study Instance UID) for each new study written to the medium, even if this study belongs to a patient recorded earlier. Study ID is assigned the value of the first Requested Procedure ID (0040, 1001) encountered in the Request Attributes Sequence (0040, 0275).

##### 5.2.1.4.1.1. Media Storage Application Profile

Refer to the table in section 5.2.1, Table 77.

##### 5.2.1.4.1.1.1. Options

Not applicable.

#### 5.2.1.5. RWA - Update File-set

Not applicable.

##### 5.2.1.5.1.1. Media Storage Application Profile

Refer to the table in section 5.2.1, Table 77.

Template ID: HSDPCPF-T-30001.09



**5.2.1.5.1.1.1. Options**

Not applicable.

**5.3. Augmented and Private Application Profiles**

Not applicable

**5.4. Media Configuration**

Not applicable.

## 6. Support of Character Sets

Any support for character sets in Network and Media services is described here.

**Table 80: Supported DICOM Character Sets**

| Character Set Description | Defined Term    | ESC Sequence            | ISO Registration Number | Code Element | Character Set   |
|---------------------------|-----------------|-------------------------|-------------------------|--------------|---|
| Latin alphabet No. 1      | ISO_IR 100      | -                       | ISO-IR 100              | G1           | Supplementary set of ISO 8859                                   |
|                           |                 | -                       | ISO-IR 6                | G0           | ISO 646   |
| Japanese                  | ISO 2022 IR 159 | ESC 02/04<br>02/0804/04 | ISO-IR 159              | G0           | JIS X 0212: Supplementary Kanji set (only for the patient name) |
| Japanese                  | ISO 2022 IR 87  | ESC 02/04 04/02         | ISO-IR 87               | G0           | JIS X 0208: Kanji and Hiragana (only for the patient name)      |
| Japanese                  | ISO 2022IR 13   | ESC 02/09 04/09         | ISO-IR 13               | G1           | JIS X 0201: Katakana (only for the patient name)                |
|                           |                 | ESC 02/08 04/10         | ISO-IR 14               | G0           | JIS X 0201: Romaji (only for the patient name)                  |

## 7. Security

### 7.1. Security Profiles

#### 7.1.1. Security use Profiles

Not applicable

#### 7.1.2. Security Transport Connection Profiles

Secure communication is a "mode of operation" supported by the implementation of the DICOM Basic TLS Secure Transport Connection Profile [DICOM]. This functionality will be used by the DICOM nodes, which can authenticate each other before they exchange DICOM information. For secure communication the TLS protocol v1.0 is used which provides message authentication, integrity, confidentiality, and replay protection. Confidentiality is optional and can be controlled by the encryption settings. The Azurion 2.0 may communicate using the following Cipher Suites:

- TLS\_RSA\_WITH\_NULL\_SHA (Node authentication without encryption)
- TLS\_RSA\_WITH\_3DES\_SHA (Node authentication with encryption)

The Azurion 2.0 supports X.509 certificates. The following TLS Certification checks will be done (TLS Handshake). The machine (either server or client) that will send its certificate will:

- Choose the certificate according to Common Name (CN) value in the Subject-field.
- This name is case-sensitive. All present certificates should have unique CN names.

The server verifies:

- That the client certificate is a X.509 certificate which is not tampered with
- That the client certificate is in the list of trusted certificates
- That the client certificate is not expired (present time is between "Valid From" and "Valid To" fields of the X.509 certificate)
- That the client certificate has the correct purpose (at least the Client Authentication purpose)

The client verifies:

- That the server certificate is a X.509 certificate which is not tampered with
- That the server certificate is in the list of trusted certificates
- That the server certificate is not expired (present time is between "Valid From" and "Valid To" fields of the X.509 certificate)
- That the server certificate has the correct purpose (at least Server Authentication purpose)

No verification is done on:

- Revocation of certificates
- Limiting the connection to a limited set of IP-addresses

Node authentication with or without encryption is only possible when both nodes have:

- An access to their own private keys
- An access to a copy of the certificate of the other node containing its public key

The Azurion 2.0 can only read certificates from the certificate stores of the HKEY\_LOCAL\_MACHINE registry key. It is the responsibility of the Hospital to setup and maintain the certificate stores. This includes the removal of revoked certificates and certificate updates prior to their expiration. Since neither X.500 directories, Lightweight Directory Access Protocol (LDAP) nor Certificate Revocation Lists (CRLs) are supported, the whole certificate chain needs to be replaced after a security breach.

The following figure presents the message flow of TLS handshake supported.




Figure 17: Message flow of TLS handshake

### 7.1.3. Digital Signature Profiles

Not applicable

### 7.1.4. Media Storage Security Profiles

Not applicable

### 7.1.5. Attribute Confidentiality Profiles

Azurion 2.0 conforms to the Basic Application Level Confidentiality Profile as a de-identifier. This functionality is targeted towards creating a special purpose, de-identified version of an already existing Data set.

Table presents all attributes that can be de-identified by the Azurion 2.0. Each Attribute to be protected has its value replaced by a different “replacement value” which does not allow identification of the patient.

Table 81: De-identified Attributes

| Attribute Name       | Tag       | VR | Replacement Value             |
|----------------------|-----------|----|-------------------------------|
| Patient Name         | 0010,0010 | PN | Assign user-specified value   |
| Patient ID           | 0010,0020 | LO | Generate and provide a new ID |
| Patient's Birth Date | 0010,0030 | DA | Make Empty                    |
| Patient's Sex        | 0010,0040 | CS | Make Empty                    |

| Attribute Name                       | Tag       | VR | Replacement Value                |
|--------------------------------------|-----------|----|----------------------------------|
| Other Patient Ids                    | 0010,1000 | LO | Make Empty                       |
| Patient's Size                       | 0010,1020 | DS | Make Empty                       |
| Patient Weight                       | 0010,1030 | DS | Make Empty                       |
| Ethnic Group                         | 0010,2160 | SH | Make Empty                       |
| Additional Patient's History         | 0010,21B0 | LT | Make Empty                       |
| Patient Comments                     | 0010,4000 | LT | Make Empty                       |
| SOP Instance UID                     | 0008,0018 | UI | Generate and provide new ID      |
| Accession Number                     | 0008,0050 | SH | Make Empty                       |
| Institution Name                     | 0008,0080 | LO | Make Empty                       |
| Referring Physician's Name           | 0008,0090 | PN | Make Empty                       |
| Device Serial Number                 | 0008,1000 | LO | Make Empty                       |
| Station Name                         | 0008,1010 | SH | Make Empty                       |
| Institutional Department Name        | 0008,1040 | LO | Make Empty                       |
| Performing Physician's Name          | 0008,1050 | PN | Make Empty                       |
| Operators' Name                      | 0008,1070 | PN | Make Empty                       |
| Referenced SOP Instance UID          | 0008,1155 | UI | Generate and provide a new ID    |
| Protocol Name                        | 0018,1030 | LO | Make Empty                       |
| Study ID                             | 0020,0010 | SH | Make Empty                       |
| Study Instance UID                   | 0020,000D | UI | Generate and provide a new ID    |
| Series Instance UID                  | 0020,000E | UI | Generate and provide a new ID    |
| Performed Procedure Step Description | 0040,0254 | LO | Make Empty                       |
| Request Attributes Sequence          | 0040,0275 | SQ | Generate and provide dummy value |

### 7.1.6. Network Address Management Profiles

Not applicable.

### 7.1.7. Time Synchronization Profiles

Azurion 2.0 conforms to the IHE CT Profile. It is possible to synchronize time with the NTP Timeserver using serviceability. The NTP Timeserver is an element of Hospital Infrastructure.

### 7.1.8. Application Configuration Management Profiles

Not applicable.

### 7.1.9. Audit Trail Profiles

The Audit Trail Component is a component of Azurion 2.0 and can create messages according to the ATNA, IHE defined standard.

Actors are information systems or components of information systems that produce, manage, or act on categories of information required by operational activities in the enterpMWL SCPe. The Audit Trail Component allows security officers in an institution to audit activities, to detect non-compliant behavior in the enterpMWL SCPe, and to facilitate detection of improper creation, access, modification and deletion of Protected Health Information (PHI), where PHI data is considered as information records (Registration, Order, Study/Procedure, Reports and to a lesser degree Images/Presentation States), and not the flow of information between the systems. This includes information exported to and imported from every secured node in the "secured domain".

The messages will be created and sent to a syslog server according to the syslog protocol. The time that is used will be the local time of the system which is synchronized with the NTP Time Server. The timeserver and syslog server are elements of the Hospital infrastructure. The following messages will be created and sent to a central Audit Record Repository

- UserCreated
- UserDeleted
- UserUpdated
- UserGroupMappingChanged



- BeginTransferring
- DataExport
- DataImport
- InstancesTransferred
- StudyCreated
- EmergencyStudyCreated
- StudyUpdated
- StudyDeleted
- SeriesUpdated
- ProcedureOpenWith
- ProcedureStarted
- ProcedureComplete
- ProcedureSuspend
- SecurityAlert
- UserLogin
- UserLoggedOff
- NodeAdded
- NodeRemoved
- PatientDeleted
- ApplicationActivity

## 7.2. Association Level Security

The Azurion 2.0 accepts associations from unknown AEs but only for Storagecommit N-Event-Report, and C-Store as SCP. If Azurion 2.0 is configured to use secure mode, then the incoming associations (for Azurion 2.0 as SCP) should follow secure mode.

## 7.3. Application Level Security

The Azurion 2.0 allows the use of either a conventional (non-secure) DICOM communication or a secure DICOM communication based on the Transport Layer Security (TLS) protocol [TLS]. If configured, the Azurion 2.0 supports security

Measures for:

- Secure authentication of a node
- Integrity and confidentiality of transmitted data
- Replay protection
- Generation of audit trail records
- Access control and user authentication.

## 8. Annexes

### 8.1. IOD Contents

#### 8.1.1. Created SOP Instances

This section specifies each IOD created by this application and specifies the content for each IOD created (including private IODs). For each attribute in the IOD the following information is supplied:

- Attribute name
- Tag
- VR – Value representation
- Value - specifies possible values
- Presence of value - specifies if attribute is always present or only under specific conditions
- Source of value - specifies the source of the value
- Comment - gives additional information on the attribute

Abbreviations used in the IOD tables for the column "Presence of Module" are:

ALWAYS The module is always present.  
CONDITIONAL The module is used under specified condition.

Abbreviations used in the Module table for the column "Presence of Value" are:

ALWAYS The attribute is always present with a value.  
EMPTY The attribute is always present without any value. (attribute sent zero length)  
VNAP The attribute is always present and its Value is Not Always Present. (attribute sent zero length if no value is present)  
ANAP The attribute is present under specified condition – if present then it will always have a value.

The abbreviations used in the Module table for the column "Source" are:

AUTO The attribute value is generated automatically.  
CONFIG The attribute value source is a configurable parameter.  
COPY The attribute value source is another SOP instance.  
FIXED The attribute value is hard-coded in the application.  
IMPLICIT The attribute value source is a user-implicit setting.  
MPPS The attribute value is the same as that use for Modality Performed Procedure Step.  
MWL The attribute value source is a Modality Worklist.  
USER The attribute value source is explicit user input.

#### 8.1.1.1. List of Created SOP Classes

Table 82: List of Created SOP Classes

| SOP Class Name  | SOP Class UID                 |
|---|-------------------------------|
| Secondary Capture Image Storage SOP Class               | 1.2.840.10008.5.1.4.1.1.7     |
| Grayscale Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.1  |
| X-Ray Angiographic Image Storage SOP Class              | 1.2.840.10008.5.1.4.1.1.12.1  |
| X-Ray Radiation Dose Structured Report SOP Class        | 1.2.840.10008.5.1.4.1.1.88.67 |

#### 8.1.1.2. Secondary Capture Image Storage SOP Class

Table 83: SOP Class Modules

| Information Entity | Module         | Presence |
|--------------------|----------------|----------|
| Patient            | Patient Module | ALWAYS   |





|        |                       |                          |
|--------|-----------------------|--------------------------|
| Study  | General Study Module  | ALWAYS                   |
|        | Patient Study Module  | ALWAYS                   |
| Series | General Series Module | ALWAYS                   |
|        | Equipment             | General Equipment Module |
| Image  | SC Equipment Module   | ALWAYS                   |
|        | General Image Module  | ALWAYS                   |
| Image  | Image Pixel Module    | ALWAYS                   |
|        | SC Image Module       | ALWAYS                   |
|        | VOI LUT Module        | ALWAYS                   |
|        | SOP Common Module     | ALWAYS                   |

**Table 84: Patient Module**

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source    | Comment |
|------------------------------|-----------|----|-------|-------------------|-----------|---------|
| Referenced Patient Sequence  | 0008,1120 | SQ |       | ANAP              | MWL       |         |
| >Referenced SOP Class UID    | 0008,1150 | UI |       | ALWAYS            | MWL       |         |
| >Referenced SOP Instance UID | 0008,1155 | UI |       | ALWAYS            | MWL       |         |
| Patient's Name               | 0010,0010 | PN |       | VNAP              | MWL, USER |         |
| Patient ID                   | 0010,0020 | LO |       | VNAP              | MWL, USER |         |
| Patient's Birth Date         | 0010,0030 | DA |       | VNAP              | MWL, USER |         |
| Patient's Sex                | 0010,0040 | CS |       | VNAP              | MWL, USER |         |
| Other Patient IDs            | 0010,1000 | LO |       | ANAP              | MWL, USER |         |
| Ethnic Group                 | 0010,2160 | SH |       | ANAP              | MWL, USER |         |
| Patient Comments             | 0010,4000 | LT |       | ANAP              | MWL       |         |

**Table 85: General Study Module**

| Attribute Name             | Tag       | VR | Value                   | Presence of Value | Source          | Comment |
|----------------------------|-----------|----|-------------------------|-------------------|-----------------|---------|
| Study Date                 | 0008,0020 | DA |                         | VNAP              | AUTO            |         |
| Study Time                 | 0008,0030 | TM |                         | VNAP              | AUTO            |         |
| Accession Number           | 0008,0050 | SH |                         | VNAP              | AUTO, MWL, USER |         |
| Referring Physician's Name | 0008,0090 | PN |                         | VNAP              | MWL             |         |
| Study Description          | 0008,1030 | LO |                         | ANAP              | AUTO, MWL       |         |
| Procedure Code Sequence    | 0008,1032 | SQ |                         | ANAP              | MWL             |         |
| >Code Value                | 0008,0100 | SH |                         | ALWAYS            |                 |         |
| >Coding Scheme Designator  | 0008,0102 | SH |                         | ALWAYS            |                 |         |
| >Code Meaning              | 0008,0104 | LO |                         | ALWAYS            |                 |         |
| Study Instance UID         | 0020,000D | UI |                         | ALWAYS            | AUTO            |         |
| Study ID                   | 0020,0010 | SH | In case the Study ID is | VNAP              | AUTO,           |         |

Template ID: HSDPCPF-T-30001.09



| Attribute Name | Tag | VR | Value   | Presence of Value | Source    | Comment |
|----------------|-----|----|---|-------------------|-----------|---------|
|                |     |    | empty the accession number will be assigned. In case Study ID and Accession Number are the same, the Study ID will be left empty. |                   | MWL, USER |         |

**Table 86: Patient Study Module**

| Attribute Name             | Tag       | VR | Value   | Presence of Value | Source    | Comment |
|----------------------------|-----------|----|---|-------------------|-----------|---------|
| Patient's Size             | 0010,1020 | DS | In meters. When received from the MWL SCP, the value can still be modified    | ANAP              | MWL, USER |         |
| Patient's Weight           | 0010,1030 | DS | In kilograms. When received from the MWL SCP, the value can still be modified | ANAP              | MWL, USER |         |
| Additional Patient History | 0010,21B0 | LT |   | ANAP              | MWL       |         |

**Table 87: General Series Module**

| Attribute Name                               | Tag       | VR | Value   | Presence of Value | Source     | Comment |
|--|-----------|----|---|-------------------|------------|---------|
| Series Date                                  | 0008,0021 | DA | <yyyymmdd>  | ANAP              | AUTO       |         |
| Series Time                                  | 0008,0031 | TM | <hhmmss>  | ANAP              | AUTO       |         |
| Modality                                     | 0008,0060 | CS | OT  | ALWAYS            | AUTO       |         |
| Series Description                           | 0008,103E | LO |   | ANAP              | AUTO       |         |
| Performing Physicians' Name                  | 0008,1050 | PN |   | ANAP              | MWL, USER  |         |
| Operators' Name                              | 0008,1070 | PN |   | ANAP              | MWL, USER  |         |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ | Identifies the MPPS SOP Instance to which this image is related | ANAP              | AUTO, MPPS |         |
| >Referenced SOP Class UID                    | 0008,1150 | UI | MPPS SOP Class UID  | ALWAYS            | MPSS       |         |
| >Referenced SOP Instance UID                 | 0008,1155 | UI | MPPS SOP Instance UID   | ALWAYS            | MPPS       |         |
| Protocol Name                                | 0018,1030 | LO |   | ANAP              | AUTO       |         |



| Attribute Name                       | Tag       | VR | Value | Presence of Value | Source     | Comment |
|--------------------------------------|-----------|----|-------|-------------------|------------|---------|
| Series Instance UID                  | 0020,000E | UI |       | ALWAYS            | AUTO       |         |
| Series Number                        | 0020,0011 | IS |       | VNAP              | AUTO       |         |
| Laterality                           | 0020,0060 | CS |       | ANAP              | AUTO       |         |
| Performed Procedure Step Start Date  | 0040,0244 | DA |       | ANAP              | AUTO, MPPS |         |
| Performed Procedure Step Start Time  | 0040,0245 | TM |       | ANAP              | AUTO, MPPS |         |
| Performed Procedure Step ID          | 0040,0253 | SH |       | ANAP              | AUTO, MPPS |         |
| Performed Procedure Step Description | 0040,0254 | LO |       | ANAP              | MPPS, MWL  |         |
| Performed Protocol Code Sequence     | 0040,0260 | SQ |       | ANAP              | MWL        |         |
| >Code Value                          | 0008,0100 | SH |       | ALWAYS            |            |         |
| >Coding Scheme Designator            | 0008,0102 | SH |       | ALWAYS            |            |         |
| >Code Meaning                        | 0008,0104 | LO |       | ALWAYS            |            |         |

Table 88: General Equipment Module

| Attribute Name                | Tag       | VR | Value   | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|---------|-------------------|--------|---------|
| Manufacturer                  | 0008,0070 | LO | Philips | VNAP              | AUTO   |         |
| Institution Name              | 0008,0080 | LO |         | ANAP              | CONFIG |         |
| Station Name                  | 0008,1010 | SH |         | ANAP              | CONFIG |         |
| Institutional Department Name | 0008,1040 | LO |         | ANAP              | AUTO   |         |
| Manufacturer's Model Name     | 0008,1090 | LO | Azurion | ANAP              | AUTO   |         |
| Device Serial Number          | 0018,1000 | LO |         | ANAP              | AUTO   |         |
| Software Versions             | 0018,1020 | LO | 2.0.0   | ANAP              | AUTO   |         |

Table 89: SC Equipment Module

| Attribute Name  | Tag       | VR | Value | Presence of Value | Source | Comment |
|-----------------|-----------|----|-------|-------------------|--------|---------|
| Modality        | 0008,0060 | CS | OT    | ANAP              | AUTO   |         |
| Conversion Type | 0008,0064 | CS | WSD   | ALWAYS            | AUTO   |         |

Table 90: General Image Module

| Attribute Name              | Tag       | VR | Value           | Presence of Value | Source | Comment |
|-----------------------------|-----------|----|-----------------|-------------------|--------|---------|
| Image Type                  | 0008,0008 | CS | DERIVED\PRIMARY | ALWAYS            | AUTO   |         |
| Content Date                | 0008,0023 | DA | <yyyymmdd>      | ALWAYS            | AUTO   |         |
| Content Time                | 0008,0033 | TM | <hhmmss>        | ALWAYS            | AUTO   |         |
| Instance Number             | 0020,0013 | IS |                 | ALWAYS            | AUTO   |         |
| Patient Orientation         | 0020,0020 | CS |                 | VNAP              | AUTO   |         |
| Lossy Image Compression     | 0028,2110 | CS | 00              | ALWAYS            | AUTO   |         |
| Icon Image Sequence         | 0088,0200 | SQ |                 | ANAP              | AUTO   |         |
| >Samples per Pixel          | 0028,0002 | US | 1               | ALWAYS            | AUTO   |         |
| >Photometric Interpretation | 0028,0004 | CS | PALETTE COLOR   | ALWAYS            | AUTO   |         |

Template ID: HSDPCPF-T-30001.09



|                       |           |       |     |        |      |  |
|-----------------------|-----------|-------|-----|--------|------|--|
| >Rows                 | 0028,0010 | US    | 240 | ALWAYS | AUTO |  |
| >Columns              | 0028,0011 | US    | 240 | ALWAYS | AUTO |  |
| >Bits Allocated       | 0028,0100 | US    | 8   | ALWAYS | AUTO |  |
| >Bits Stored          | 0028,0101 | US    | 8   | ALWAYS | AUTO |  |
| >High Bit             | 0028,0102 | US    | 7   | ALWAYS | AUTO |  |
| >Pixel Representation | 0028,0103 | US    | 0   | ALWAYS | AUTO |  |
| > Pixel Data          | 7FE0,0010 | OW/OB |     | ALWAYS | AUTO |  |

**Table 91: Image Pixel Module**

| Attribute Name             | Tag       | VR    | Value  | Presence of Value | Source | Comment |
|----------------------------|-----------|-------|--|-------------------|--------|---------|
| Samples per Pixel          | 0028,0002 | US    |  | ALWAYS            | AUTO   |         |
| Photometric Interpretation | 0028,0004 | CS    | RGB  | ALWAYS            | AUTO   |         |
| Planar Configuration       | 0028,0006 | US    |  | ANAP              | AUTO   |         |
| Rows                       | 0028,0010 | US    | 1024 or 512  | ALWAYS            | AUTO   |         |
| Columns                    | 0028,0011 | US    | 1024 or 512  | ALWAYS            | AUTO   |         |
| Bits Allocated             | 0028,0100 | US    | 16 Note: For Snapshot function, this number can be 8 | ALWAYS            | AUTO   |         |
| Bits Stored                | 0028,0101 | US    | 8 and 12   | ALWAYS            | AUTO   |         |
| High Bit                   | 0028,0102 | US    |  | ALWAYS            | AUTO   |         |
| Pixel Representation       | 0028,0103 | US    | 0000H  | ALWAYS            | AUTO   |         |
| Smallest Image Pixel Value | 0028,0106 | US/SS |  | ANAP              | AUTO   |         |
| Largest Image Pixel Value  | 0028,0107 | US/SS |  | ANAP              | AUTO   |         |
| Pixel Data                 | 7FE0,0010 | OW/OB |  | ALWAYS            | AUTO   |         |

**Table 92: SC Image Module**

| Attribute Name                 | Tag       | VR | Value | Presence of Value | Source | Comment |
|--------------------------------|-----------|----|-------|-------------------|--------|---------|
| Date of Secondary Capture      | 0018,1012 | DA |       | ANAP              | AUTO   |         |
| Time of Secondary Capture      | 0018,1014 | TM |       | ANAP              | AUTO   |         |
| Nominal Scanned Pixel Spacing  | 0018,2010 | DS |       | ANAP              | AUTO   |         |
| Pixel Spacing Calibration Type | 0028,0A02 | CS |       | ANAP              | AUTO   |         |

**Table 93: VOI LUT Module**

| Attribute Name | Tag       | VR | Value | Presence of Value | Source | Comment |
|----------------|-----------|----|-------|-------------------|--------|---------|
| Window Center  | 0028,1050 | DS |       | ALWAYS            | AUTO   |         |
| Window Width   | 0028,1051 | DS |       | ALWAYS            | AUTO   |         |

**Table 94: SOP Common Module**

| Attribute Name         | Tag       | VR | Value      | Presence of Value | Source | Comment |
|------------------------|-----------|----|------------|-------------------|--------|---------|
| Specific Character Set | 0008,0005 | CS | ISO_IR 100 | ANAP              | AUTO   |         |



|                        |           |    |                           |        |      |  |
|------------------------|-----------|----|---------------------------|--------|------|--|
| Instance Creation Date | 0008,0012 | DA |                           | ANAP   | AUTO |  |
| Instance Creation Time | 0008,0013 | TM |                           | ANAP   | AUTO |  |
| SOP Class UID          | 0008,0016 | UI | 1.2.840.10008.5.1.4.1.1.7 | ALWAYS | AUTO |  |
| SOP Instance UID       | 0008,0018 | UI |                           | ALWAYS | AUTO |  |
| Instance Number        | 0020,0013 | IS |                           | ANAP   | AUTO |  |

### 8.1.1.3. Grayscale Softcopy Presentation State Storage SOP Class

**Table 95: SOP Class Modules**

| Information Entity      | Module                                   | Presence |
|-------------------------|--|----------|
| Patient                 | Patient Module                           | Always   |
| Study                   | General Study Module                     | Always   |
|                         | Patient Study Module                     | Always   |
| Series                  | General Series Module                    | Always   |
|                         | Presentation Series Module               | Always   |
| Equipment               | General Equipment Module                 | Always   |
| Presentation State      | Presentation State Identification Module | Always   |
|                         | Presentation State Relationship Module   | Always   |
|                         | Presentation State Shutter Module        | Always   |
|                         | Display Shutter Module                   | Always   |
|                         | Bitmap Display Shutter Module            | Always   |
|                         | Displayed Area Module                    | Always   |
|                         | Graphic Layer Module                     | Always   |
|                         | Softcopy Presentation LUT Module         | Always   |
|                         | SOP Common Module                        | Always   |
| Softcopy VOI LUT module | Always                                   |          |

**Table 96: Patient Module**

| Attribute Name               | Tag       | VR | Value      | Presence of Value | Source       | Comment |
|------------------------------|-----------|----|------------|-------------------|--------------|---------|
| Referenced Patient Sequence  | 0008,1120 | SQ |            | ANAP              | MWL          |         |
| >Referenced SOP Class UID    | 0008,1150 | UI |            | ALWAYS            | MWL          |         |
| >Referenced SOP Instance UID | 0008,1155 | UI |            | ALWAYS            | MWL          |         |
| Patient's Name               | 0010,0010 | PN |            | VNAP              | MWL,<br>USER |         |
| Patient ID                   | 0010,0020 | LO |            | VNAP              | MWL,<br>USER |         |
| Patient's Birth Date         | 0010,0030 | DA | <yyyymmdd> | VNAP              | MWL,<br>USER |         |
| Patient's Sex                | 0010,0040 | CS |            | VNAP              | MWL,<br>USER |         |
| Other Patient IDs            | 0010,1000 | LO |            | ANAP              | MWL,<br>USER |         |
| Ethnic Group                 | 0010,2160 | SH |            | ANAP              | MWL,<br>USER |         |
| Patient Comments             | 0010,4000 | LT |            | ANAP              | MWL,<br>USER |         |

**Table 97: General Study Module**

| Attribute Name               | Tag       | VR | Value   | Presence of Value | Source                | Comment |
|------------------------------|-----------|----|---|-------------------|-----------------------|---------|
| Study Date                   | 0008,0020 | DA | <yyyymmdd>  | VNAP              | AUTO                  |         |
| Study Time                   | 0008,0030 | TM |   | VNAP              | AUTO                  |         |
| Accession Number             | 0008,0050 | SH |   | VNAP              | AUTO,<br>MWL,<br>USER |         |
| Referring Physician's Name   | 0008,0090 | PN |   | VNAP              | MWL                   |         |
| Study Description            | 0008,1030 | LO |   | ANAP              | AUTO                  |         |
| Procedure Code Sequence      | 0008,1032 | SQ |   | ANAP              | MWL,AUTO              |         |
| >Code Value                  | 0008,0100 | SH |   | ALWAYS            | MWL,AUTO              |         |
| >Coding Scheme Designator    | 0008,0102 | SH |   | ALWAYS            | MWL,AUTO              |         |
| >Code Meaning                | 0008,0104 | LO |   | ALWAYS            | MWL,AUTO              |         |
| Referenced Study Sequence    | 0008,1110 | SQ |   | ANAP              | MWL                   |         |
| >Referenced SOP Class UID    | 0008,1150 | UI |   | ALWAYS            | MWL                   |         |
| >Referenced SOP Instance UID | 0008,1155 | UI |   | ALWAYS            | MWL                   |         |
| Study Instance UID           | 0020,000D | UI |   | ALWAYS            | AUTO                  |         |
| Study ID                     | 0020,0010 | SH | In case the Study ID is empty the accession number will be assigned. In case Study ID and Accession Number are the same, the Study ID will be left empty. | VNAP              | AUTO,<br>MWL,<br>USER |         |

**Table 98: Patient Study Module**

| Attribute Name             | Tag       | VR | Value   | Presence of Value | Source       | Comment |
|----------------------------|-----------|----|---|-------------------|--------------|---------|
| Patient's Size             | 0010,1020 | DS | In meters. When received from the MWL SCP, the value can still be modified    | ANAP              | MWL,<br>USER |         |
| Patient's Weight           | 0010,1030 | DS | In kilograms. When received from the MWL SCP, the value can still be modified | ANAP              | MWL,<br>USER |         |
| Additional Patient History | 0010,21B0 | LT |   | ANAP              | MWL          |         |



Table 99: General Series Module

| Attribute Name                               | Tag       | R  | Value   | Presence of Value | Source           | Comment |
|--|-----------|----|---|-------------------|------------------|---------|
| Series Date                                  | 0008,0021 | DA | <yyyymmdd>  | ANAP              | AUTO             |         |
| Series Time                                  | 0008,0031 | TM | <hhmmss>  | ANAP              | AUTO             |         |
| Modality                                     | 0008,0060 | CS | PR  | ALWAYS            | AUTO             |         |
| Series Description                           | 0008,103E | LO |   | ALWAYS            | AUTO             |         |
| Performing Physicians' Name                  | 0008,1050 | PN |   | ANAP              | MWL, USER        |         |
| Operators' Name                              | 0008,1070 | PN |   | ANAP              | MWL, USER        |         |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ | Identifies the MPPS SOP Instance to which this image is related | ANAP              | AUTO, MPPS       |         |
| >Referenced SOP Class UID                    | 0008,1150 | UI | MPPS SOP Class UID  | ALWAYS            | MWL              |         |
| >Referenced SOP Instance UID                 | 0008,1155 | UI | MPPS SOP Instance UID   | ALWAYS            | MWL              |         |
| Protocol Name                                | 0018,1030 | LO |   | ANAP              | AUTO             |         |
| Series Instance UID                          | 0020,000E | UI |   | ALWAYS            | AUTO             |         |
| Series Number                                | 0020,0011 | IS |   | VNAP              | AUTO             |         |
| Laterality                                   | 0020,0060 | CS |   | ANAP              | AUTO             |         |
| Performed Procedure Step Start Date          | 0040,0244 | DA | <yyyymmdd>  | ALWAYS            | AUTO, MPPS       |         |
| Performed Procedure Step Start Time          | 0040,0245 | TM | <hhmmss>  | ALWAYS            | AUTO, MPPS       |         |
| Performed Procedure Step ID                  | 0040,0253 | SH | Same as MPPS  | ALWAYS            | AUTO, MPPS, USER |         |
| Performed Procedure Step Description         | 0040,0254 | LO |   | ANAP              | MPPS, MWL        |         |
| Request Attributes Sequence                  | 0040,0275 | SQ |   | ANAP              | MWL              |         |
| >Accession Number                            | 0008,0050 | SH |   | ANAP              | MWL              |         |
| >Issuer of Accession Number Sequence         | 0008,0051 | SQ |   | ANAP              | MWL              |         |
| >Referenced Study Sequence                   | 0008,1110 | SQ |   | ANAP              | MWL              |         |
| >>Referenced SOP Class UID                   | 0008,1150 | UI |   | ALWAYS            | MWL              |         |
| >>Referenced SOP Instance UID                | 0008,1155 | UI |   | ALWAYS            | MWL              |         |
| >Study Instance UID                          | 0020,000D | UI |   | ANAP              | MWL              |         |
| >Requested Procedure Description             | 0032,1060 | LO |   | ANAP              | MWL              |         |
| >Requested Procedure Code Sequence           | 0032,1064 | SQ |   | ANAP              | MWL              |         |
| >>Code Value                                 | 0008,0100 | SH |   | ALWAYS            | MWL              |         |
| >>Coding Scheme Designator                   | 0008,0102 | SH |   | ALWAYS            | MWL              |         |
| >>Code Meaning                               | 0008,0104 | LO |   | ALWAYS            | MWL              |         |
| >Scheduled Protocol Code Sequence            | 0040,0008 | SQ |   | ANAP              | MWL              |         |
| >>Code Value                                 | 0008,0100 | SH |   | ALWAYS            | MWL              |         |
| >>Coding Scheme Designator                   | 0008,0102 | SH |   | ALWAYS            | MWL              |         |
| >>Code Meaning                               | 0008,0104 | LO |   | ALWAYS            | MWL              |         |
| >>Protocol Context Sequence                  | 0040,0440 | SQ |   | ANAP              | MWL              |         |
| >>>Value Type                                | 0040,A040 | CS |   | ALWAYS            | MWL              |         |
| >>>Concept Name Code Sequence                | 0040,A043 | SQ |   | ALWAYS            | MWL              |         |

Template ID: HSDPCPF-T-30001.09



| Attribute Name                                | Tag       | R  | Value | Presence of Value | Source | Comment |
|---|-----------|----|-------|-------------------|--------|---------|
| >>>>Code Value                                | 0008,0100 | SH |       | ALWAYS            | MWL    |         |
| >>>>Coding Scheme Designator                  | 0008,0102 | SH |       | ALWAYS            | MWL    |         |
| >>>>Code Meaning                              | 0008,0104 | LO |       | ALWAYS            | MWL    |         |
| >Reason for the Requested Procedure           | 0040,1002 | LO |       | ANAP              | MWL    |         |
| >Reason for Requested Procedure Code Sequence | 0040,100A | SQ |       | ANAP              | MWL    |         |
| >>Code Value                                  | 0008,0100 | SH |       | ALWAYS            | MWL    |         |
| >>Coding Scheme Designator                    | 0008,0102 | SH |       | ALWAYS            | MWL    |         |
| >>Code Meaning                                | 0008,0104 | LO |       | ALWAYS            | MWL    |         |
| >>Context Identifier                          | 0008,010F | CS |       | ANAP              | MWL    |         |
| >>Context UID                                 | 0008,0117 | UI |       | ANAP              | MWL    |         |

Table 100: Presentation Series Module

| Attribute Name | Tag       | VR | Value | Presence of Value | Source | Comment |
|----------------|-----------|----|-------|-------------------|--------|---------|
| Modality       | 0008,0060 | CS | PR    | ALWAYS            | AUTO   |         |

Table 101: General Equipment Module

| Attribute Name                | Tag       | VR | Value   | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|---------|-------------------|--------|---------|
| Manufacturer                  | 0008,0070 | LO | Philips | VNAP              | AUTO   |         |
| Institution Name              | 0008,0080 | LO |         | ANAP              | AUTO   |         |
| Station Name                  | 0008,1010 | SH |         | ANAP              | AUTO   |         |
| Institutional Department Name | 0008,1040 | LO |         | ANAP              | AUTO   |         |
| Manufacturer's Model Name     | 0008,1090 | LO | Azurion | ANAP              | AUTO   |         |
| Device Serial Number          | 0018,1000 | LO |         | ANAP              | AUTO   |         |
| Software Versions             | 0018,1020 | LO | 2.0.0   | ANAP              | AUTO   |         |

Table 102: Presentation State Identification Module

| Attribute Name             | Tag       | VR | Value        | Presence of Value | Source | Comment |
|----------------------------|-----------|----|--------------|-------------------|--------|---------|
| Instance Number            | 0020,0013 | IS |              | ALWAYS            | AUTO   |         |
| Content Label              | 0070,0080 | CS | AS LAST SEEN | ALWAYS            | AUTO   |         |
| Content Description        | 0070,0081 | LO |              | VNAP              | AUTO   |         |
| Presentation Creation Date | 0070,0082 | DA | <yyyymmdd>   | ALWAYS            | AUTO   |         |
| Presentation Creation Time | 0070,0083 | TM | <hhmmss>     | ALWAYS            | AUTO   |         |
| Content Creator's Name     | 0070,0084 | PN |              | VNAP              | AUTO   |         |

Table 103: Presentation State Relationship Module

| Attribute Name             | Tag       | R  | Value | Presence of Value | Source | Comment |
|----------------------------|-----------|----|-------|-------------------|--------|---------|
| Referenced Series Sequence | 0008,1115 | SQ |       | ALWAYS            | AUTO   |         |
| >Referenced Image Sequence | 0008,1140 | SQ |       | ALWAYS            | AUTO   |         |





|                               |           |    |  |        |      |  |
|-------------------------------|-----------|----|--|--------|------|--|
| >>Referenced SOP Class UID    | 0008,1150 | UI |  | ALWAYS | AUTO |  |
| >>Referenced SOP Instance UID | 0008,1155 | UI |  | ALWAYS | AUTO |  |
| >Series Instance UID          | 0020,000E | UI |  | ALWAYS | AUTO |  |

**Table 104: Presentation State Shutter Module**

| Attribute Name             | Tag       | VR | Value | Presence of Value | Source | Comment |
|----------------------------|-----------|----|-------|-------------------|--------|---------|
| Shutter Presentation Value | 0018,1622 | US |       | ANAP              | AUTO   |         |

**Table 105: Display Shutter Module**

| Attribute Name                | Tag       | VR | Value       | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|-------------|-------------------|--------|---------|
| Shutter Shape                 | 0018,1600 | CS | RECTANGULAR | ALWAYS            | AUTO   |         |
| Shutter Left Vertical Edge    | 0018,1602 | IS |             | ANAP              | AUTO   |         |
| Shutter Right Vertical Edge   | 0018,1604 | IS |             | ANAP              | AUTO   |         |
| Shutter Upper Horizontal Edge | 0018,1606 | IS |             | ANAP              | AUTO   |         |
| Shutter Lower Horizontal Edge | 0018,1608 | IS |             | ANAP              | AUTO   |         |
| Shutter Presentation Value    | 0018,1622 | US |             | ANAP              | AUTO   |         |

**Table 106: Displayed Area Module**

| Attribute Name                           | Tag       | VR | Value        | Presence of Value | Source | Comment  |
|--|-----------|----|--------------|-------------------|--------|--|
| Displayed Area Selection Sequence        | 0070,005A | SQ |              | ALWAYS            | AUTO   |  |
| >Referenced Image Sequence               | 0008,1140 | SQ |              | ANAP              | AUTO   |  |
| >>Referenced SOP Class UID               | 0008,1150 | UI |              | ALWAYS            | AUTO   |  |
| >>Referenced SOP Instance UID            | 0008,1155 | UI |              | ALWAYS            | AUTO   |  |
| >Displayed Area Top Left Hand Corner     | 0070,0052 | SL |              | ALWAYS            | AUTO   | calculated by the system but initiated by user i.e. without user zoom and pan displayed area does not have any meaning |
| >Displayed Area Bottom Right Hand Corner | 0070,0053 | SL |              | ALWAYS            | AUTO   |  |
| >Presentation Size Mode                  | 0070,0100 | CS | SCALE TO FIT | ALWAYS            |        |  |
| >Presentation Pixel Aspect Ratio         | 0070,0102 | IS |              | ANAP              | AUTO   |  |

**Table 107: Graphic Layer Module**

| Attribute Name         | Tag       | VR | Value    | Presence of Value | Source | Comment |
|------------------------|-----------|----|----------|-------------------|--------|---------|
| Graphic Layer Sequence | 0070,0060 | SQ |          | ALWAYS            | AUTO   |         |
| >Graphic Layer         | 0070,0002 | CS | GRAPHICS | ALWAYS            | AUTO   |         |
| >Graphic Layer Order   | 0070,0062 | IS |          | ALWAYS            | AUTO   |         |

**Table 108: Softcopy Presentation LUT Module**

| Attribute Name         | Tag       | VR | Value    | Presence of Value | Source | Comment |
|------------------------|-----------|----|----------|-------------------|--------|---------|
| Presentation LUT Shape | 2050,0020 | CS | IDENTITY | ANAP              | AUTO   |         |

**Table 109: SOP Common Module**

| Attribute Name         | Tag       | VR | Value                        | Presence of Value | Source | Comment |
|------------------------|-----------|----|------------------------------|-------------------|--------|---------|
| Specific Character Set | 0008,0005 | CS | ISO_IR 100                   | ANAP              | AUTO   |         |
| Instance Creation Date | 0008,0012 | DA |                              | ANAP              | AUTO   |         |
| Instance Creation Time | 0008,0013 | TM |                              | ANAP              | AUTO   |         |
| SOP Class UID          | 0008,0016 | UI | 1.2.840.10008.5.1.4.1.1.11.1 | ALWAYS            | AUTO   |         |
| SOP Instance UID       | 0008,0018 | UI |                              | ALWAYS            | AUTO   |         |
| Instance Number        | 0020,0013 | IS |                              | ANAP              | AUTO   |         |

**Table 110: Softcopy VOI LUT module**

| Attribute Name                     | Tag       | VR | Value | Presence of Value | Source | Comment |
|------------------------------------|-----------|----|-------|-------------------|--------|---------|
| Softcopy VOI LUT Sequence          | 0028,3110 | SQ |       | ALWAYS            | AUTO   |         |
| >Referenced Image Sequence         | 0008,1140 | SQ |       | ANAP              | AUTO   |         |
| >>Referenced SOP Class UID         | 0008,1150 | UI |       | ALWAYS            | AUTO   |         |
| >>Referenced SOP Instance UID      | 0008,1155 | UI |       | ALWAYS            | AUTO   |         |
| >Window Center                     | 0028,1050 | DS |       | ALWAYS            | AUTO   |         |
| >Window Width                      | 0028,1051 | DS |       | ALWAYS            | AUTO   |         |
| >Window Center & Width Explanation | 0028,1055 | LO |       | ANAP              | AUTO   |         |
| >VOI LUT Function                  | 0028,1056 | CS |       | ANAP              | AUTO   |         |

**8.1.1.4. X-Ray Angiographic Image Storage SOP Class**

**Table 111: SOP Class Modules**

| Information Entity | Module                   | Presence |
|--------------------|--------------------------|----------|
| Patient            | Patient Module           | Always   |
| Study              | General Study Module     | Always   |
|                    | Patient Study Module     | Always   |
| Series             | General Series Module    | Always   |
| Equipment          | General Equipment Module | Always   |
| Image              | General Image Module     | Always   |
|                    | Image Pixel Module       | Always   |
|                    | Contrast/Bolus Module    | Always   |
|                    | Cine Module              | Always   |
|                    | Multi-Frame Module       | Always   |
|                    | Display Shutter Module   | Always   |
|                    | X-Ray Image Module       | Always   |
|                    | X-Ray Acquisition Module | Always   |
| X-Ray Table Module | Always                   |          |



|                      |              |
|----------------------|--------------|
| XA Positioner Module | Always       |
| DX Detector Module   | Always       |
| Modality LUT Module  | Always       |
| VOI LUT Module       | Always       |
| Curve Module         | User Defined |
| SOP Common Module    | Always       |

**Table 112: Patient Module**

| Attribute Name               | Tag       | VR | Value      | Presence of Value | Source    | Comment |
|------------------------------|-----------|----|------------|-------------------|-----------|---------|
| Referenced Patient Sequence  | 0008,1120 | SQ |            | ANAP              | MWL       |         |
| >Referenced SOP Class UID    | 0008,1150 | UI |            | ALWAYS            | MWL       |         |
| >Referenced SOP Instance UID | 0008,1155 | UI |            | ALWAYS            | MWL       |         |
| Patient's Name               | 0010,0010 | PN |            | VNAP              | MWL, USER |         |
| Patient ID                   | 0010,0020 | LO |            | VNAP              | MWL, USER |         |
| Patient's Birth Date         | 0010,0030 | DA | <yyyymmdd> | VNAP              | MWL, USER |         |
| Patient's Sex                | 0010,0040 | CS |            | VNAP              | MWL, USER |         |
| Other Patient IDs            | 0010,1000 | LO |            | ANAP              | MWL, USER |         |
| Ethnic Group                 | 0010,2160 | SH |            | ANAP              | MWL, USER |         |
| Patient Comments             | 0010,4000 | LT |            | ANAP              | MWL, USER |         |

**Table 113: General Study Module**

| Attribute Name             | Tag       | VR | Value   | Presence of Value | Source          | Comment |
|----------------------------|-----------|----|---|-------------------|-----------------|---------|
| Study Date                 | 0008,0020 | DA | <yyyymmdd>  | VNAP              | AUTO            |         |
| Study Time                 | 0008,0030 | TM | <hhmmss>  | VNAP              | AUTO, MWL, USER |         |
| Accession Number           | 0008,0050 | SH |   | VNAP              | MWL             |         |
| Referring Physician's Name | 0008,0090 | PN | Patient's referring physician.  | VNAP              | MWL             |         |
| Study Description          | 0008,1030 | LO | Based on configuration<br>Study Description is:-not exported -<br>based on schedule procedure step description(WLM) -<br>based on requested procedure step description(WLM) -<br>internal generated performed procedure description | ANAP              | AUTO, MWL       |         |
| Procedure Code Sequence    | 0008,1032 | SQ |   | ANAP              | MWL,AUTO        |         |

|                              |           |    |  |        |          |  |
|------------------------------|-----------|----|--|--------|----------|--|
| >Code Value                  | 0008,0100 | SH |  | ALWAYS | MWL,AUTO |  |
| >Coding Scheme Designator    | 0008,0102 | SH |  | ALWAYS | MWL,AUTO |  |
| >Code Meaning                | 0008,0104 | LO |  | ALWAYS | MWL,AUTO |  |
| Referenced Study Sequence    | 0008,1110 | SQ |  | ANAP   | MWL      |  |
| >Referenced SOP Class UID    | 0008,1150 | UI |  | ALWAYS | MWL      |  |
| >Referenced SOP Instance UID | 0008,1155 | UI |  | ALWAYS | MWL      |  |
| Study Instance UID           | 0020,000D | UI |  | ALWAYS | MWL      |  |
| Study ID                     | 0020,0010 | SH |  | VNAP   | MWL      |  |

**Table 114: Patient Study Module**

| Attribute Name             | Tag       | VR | Value  | Presence of Value | Source           | Comment |
|----------------------------|-----------|----|--|-------------------|------------------|---------|
| Patient's Size             | 0010,1020 | DS | In meters. When received from the MWL SCP, the value can still be modified.    | ANAP              | MWL,<br>US<br>ER |         |
| Patient's Weight           | 0010,1030 | DS | In kilograms. When received from the MWL SCP, the value can still be modified. | ANAP              | MWL,<br>US<br>ER |         |
| Additional Patient History | 0010,21B0 | LT |  | ANAP              | MWL              |         |

**Table 115: General Series Module**

| Attribute Name                               | Tag       | VR | Value   | Presence of Value | Source            | Comment |
|--|-----------|----|---|-------------------|-------------------|---------|
| Series Date                                  | 0008,0021 | DA | <yyyymmdd>  | ANAP              | AUTO              |         |
| Series Time                                  | 0008,0031 | TM | <hhmmss>  | ANAP              | AUTO              |         |
| Modality                                     | 0008,0060 | CS | XA  | ALWAYS            | AUTO              |         |
| Series Description                           | 0008,103E | LO |   | ANAP              | AUTO              |         |
| Performing Physicians' Name                  | 0008,1050 | PN |   | ANAP              | MWL,<br>US<br>ER  |         |
| Operators' Name                              | 0008,1070 | PN |   | ANAP              | MWL,<br>US<br>ER  |         |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ | Identifies the MPPS SOP Instance to which this image is related | ANAP              | AUTO,<br>MP<br>PS |         |
| >Referenced SOP Class UID                    | 0008,1150 | UI | MPPS SOP Class UID  | ALWAYS            | MWL               |         |
| >Referenced SOP Instance UID                 | 0008,1155 | UI | MPPS SOP Instance UID   | ALWAYS            | MWL               |         |
| Protocol Name                                | 0018,1030 | LO |   | ANAP              | AUTO              |         |
| Series Instance UID                          | 0020,000E | UI |   | ALWAYS            | AUTO              |         |
| Series Number                                | 0020,0011 | IS |   | VNAP              |                   |         |
| Laterality                                   | 0020,0060 | CS |   | EMPTY             | AUTO              |         |
| Performed Procedure Step Start Date          | 0040,0244 | DA | <yyyymmdd>  | ANAP              | AUTO,<br>MP<br>PS |         |



| Attribute Name                        | Tag       | VR | Value        | Presence of Value | Source                             | Comment |
|---------------------------------------|-----------|----|--------------|-------------------|------------------------------------|---------|
| Performed Procedure Step Start Time   | 0040,0245 | TM |              | ANAP              | AUTO                               |         |
| Performed Procedure Step ID           | 0040,0253 | SH | Same as MPPS | ANAP              | AUTO,<br>MP<br>PS<br>,<br>US<br>ER |         |
| Performed Procedure Step Description  | 0040,0254 | LO | Same as MPPS | ANAP              | AUTO,<br>MP<br>PS<br>,<br>US<br>ER |         |
| Request Attributes Sequence           | 0040,0275 | SQ |              | ANAP              | MWL                                |         |
| >Accession Number                     | 0008,0050 | SH |              | ANAP              | MWL                                |         |
| >Issuer of Accession Number Sequence  | 0008,0051 | SQ |              | ANAP              | MWL                                |         |
| >Referenced Study Sequence            | 0008,1110 | SQ |              | ANAP              | MWL                                |         |
| >>Referenced SOP Class UID            | 0008,1150 | UI |              | ALWAYS            | MWL                                |         |
| >>Referenced SOP Instance UID         | 0008,1155 | UI |              | ALWAYS            | MWL                                |         |
| >Study Instance UID                   | 0020,000D | UI |              | ANAP              | MWL                                |         |
| >Requested Procedure Description      | 0032,1060 | LO |              | ANAP              | MWL                                |         |
| >Requested Procedure Code Sequence    | 0032,1064 | SQ |              | ANAP              | MWL                                |         |
| >>Code Value                          | 0008,0100 | SH |              | ALWAYS            | MWL                                |         |
| >>Coding Scheme Designator            | 0008,0102 | SH |              | ALWAYS            | MWL                                |         |
| >>Code Meaning                        | 0008,0104 | LO |              | ALWAYS            | MWL                                |         |
| >>Context Identifier                  | 0008,010F | CS |              | ANAP              | MWL                                |         |
| >>Context UID                         | 0008,0117 | UI |              | ANAP              | MWL                                |         |
| >Scheduled Procedure Step Description | 0040,0007 | LO |              | ANAP              | MWL                                |         |
| >Scheduled Protocol Code Sequence     | 0040,0008 | SQ |              | ANAP              | MWL                                |         |
| >>Code Value                          | 0008,0100 | SH |              | ALWAYS            | MWL                                |         |
| >>Coding Scheme Designator            | 0008,0102 | SH |              | ALWAYS            | MWL                                |         |
| >>Code Meaning                        | 0008,0104 | LO |              | ALWAYS            | MWL                                |         |
| >>Context Identifier                  | 0008,010F | CS |              | ANAP              | MWL                                |         |
| >>Context UID                         | 0008,0117 | UI |              | ANAP              | MWL                                |         |
| >>Protocol Context Sequence           | 0040,0440 | SQ |              | ANAP              | MWL                                |         |
| >>>Content Item Modifier Sequence     | 0040,0441 | SQ |              | ANAP              | MWL                                |         |
| >>>>Value Type                        | 0040,A040 | CS |              | ALWAYS            | MWL                                |         |
| >>>>Concept Name Code Sequence        | 0040,A043 | SQ |              | ALWAYS            | MWL                                |         |
| >>>>>Code Value                       | 0008,0100 | SH |              | ALWAYS            | MWL                                |         |
| >>>>>Coding Scheme Designator         | 0008,0102 | SH |              | ALWAYS            | MWL                                |         |
| >>>>>Code Meaning                     | 0008,0104 | LO |              | ALWAYS            | MWL                                |         |
| >>>>Value Type                        | 0040,A040 | CS |              | ALWAYS            | MWL                                |         |
| >>>>Concept Name Code Sequence        | 0040,A043 | SQ |              | ALWAYS            | MWL                                |         |
| >>>>>Code Value                       | 0008,0100 | SH |              | ALWAYS            | MWL                                |         |
| >>>>>Coding Scheme Designator         | 0008,0102 | SH |              | ALWAYS            | MWL                                |         |
| >>>>>Code Meaning                     | 0008,0104 | LO |              | ALWAYS            | MWL                                |         |
| >>>>>Context Identifier               | 0008,010F | CS |              | ANAP              | MWL                                |         |
| >>>>>Context UID                      | 0008,0117 | UI |              | ANAP              | MWL                                |         |

| Attribute Name                                | Tag       | VR | Value | Presence of Value | Source | Comment |
|---|-----------|----|-------|-------------------|--------|---------|
| >Scheduled Procedure Step ID                  | 0040,0009 | SH |       | ANAP              | MWL    |         |
| >Reason for Requested Procedure Code Sequence | 0040,100A | SQ |       | ANAP              | MWL    |         |
| >>Code Value                                  | 0008,0100 | SH |       | ALWAYS            | MWL    |         |
| >>Coding Scheme Designator                    | 0008,0102 | SH |       | ALWAYS            | MWL    |         |
| >>Code Meaning                                | 0008,0104 | LO |       | ALWAYS            | MWL    |         |
| >>Context Identifier                          | 0008,010F | CS |       | ANAP              | MWL    |         |
| >>Context UID                                 | 0008,0117 | UI |       | ANAP              | MWL    |         |

**Table 116: General Equipment Module**

| Attribute Name                | Tag       | VR | Value   | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|---------|-------------------|--------|---------|
| Manufacturer                  | 0008,0070 | LO | Philips | VNAP              | AUTO   |         |
| Institution Name              | 0008,0080 | LO |         | ANAP              | AUTO   |         |
| Station Name                  | 0008,1010 | SH |         | ANAP              | AUTO   |         |
| Institutional Department Name | 0008,1040 | LO |         | ANAP              | AUTO   |         |
| Manufacturer's Model Name     | 0008,1090 | LO | Azurion | ANAP              | AUTO   |         |
| Device Serial Number          | 0018,1000 | LO |         | ANAP              | AUTO   |         |
| Software Versions             | 0018,1020 | LO | 2.0.0   | ANAP              | AUTO   |         |

**Table 117: General Image Module**

| Attribute Name   | Tag       | VR | Value   | Presence of Value | Source | Comment |
|------------------|-----------|----|---|-------------------|--------|---------|
| Image Type       | 0008,0008 | CS | Applied value(s): ORIGINAL or DERIVED (if subtraction has been processed into the image)<br>Value 2: PRIMARY<br>Value 3:<br><ul style="list-style-type: none"> <li>SINGLE PLANE (if the image is a single plane acquisition)</li> <li>BIPLANE A (if the image is the first plane of a Bi-plane acquisition)</li> <li>BIPLANE B (if the image is the second plane of a Bi-plane acquisition)</li> </ul> Value 4:<br><ul style="list-style-type: none"> <li>SINGLE A (if the image is derived from plane A of a biplane image and sent as a SINGLE PLANE image)</li> <li>SINGLE B (if the image is derived from plane B of a biplane image and sent as a SINGLE PLANE image)</li> </ul> | ANAP              | AUTO   |         |
| Acquisition Date | 0008,0022 | DA | <yyyymmdd>  | ANAP              | AUTO   |         |



|                              |           |    |            |        |      |  |
|------------------------------|-----------|----|------------|--------|------|--|
| Content Date                 | 0008,0023 | DA | <yyyymmdd> | ANAP   | AUTO |  |
| Acquisition Time             | 0008,0032 | TM | <hhmmss>   | ANAP   | AUTO |  |
| Content Time                 | 0008,0033 | TM | <hhmmss>   | ANAP   |      |  |
| Referenced Image Sequence    | 0008,1140 | SQ |            | ANAP   |      |  |
| >Referenced SOP Class UID    | 0008,1150 | UI |            | ALWAYS |      |  |
| >Referenced SOP Instance UID | 0008,1155 | UI |            | ALWAYS |      |  |
| Derivation Description       | 0008,2111 | ST |            | ANAP   |      |  |
| Instance Number              | 0020,0013 | IS |            | VNAP   |      |  |
| Patient Orientation          | 0020,0020 | CS |            | ANAP   |      |  |
| Lossy Image Compression      | 0028,2110 | CS | 00         | ANAP   | AUTO |  |
| Icon Image Sequence          | 0088,0200 | SQ |            | ANAP   | AUTO |  |
| > Samples per Pixel          | 0028,0002 | US |            | ALWAYS | AUTO |  |
| > Photometric Interpretation | 0028,0004 | CS |            | ALWAYS | AUTO |  |
| > Rows                       | 0028,0010 | US |            | ALWAYS | AUTO |  |
| > Columns                    | 0028,0011 | US |            | ALWAYS | AUTO |  |
| > Bits Allocated             | 0028,0100 | US |            | ALWAYS | AUTO |  |
| > Bits Stored                | 0028,0101 | US |            | ALWAYS | AUTO |  |
| > High Bit                   | 0028,0102 | US |            | ALWAYS | AUTO |  |
| >Pixel Representation        | 0028,0103 | US |            | ALWAYS | AUTO |  |
| >Pixel Data                  | 7FE0,0010 | OW |            | ALWAYS | AUTO |  |
| Presentation LUT Shape       | 2050,0020 | CS |            | ALWAYS |      |  |

Table 118: Image Pixel Module

| Attribute Name             | Tag       | VR        | Value  | Presence of Value | Source | Comment |
|----------------------------|-----------|-----------|--|-------------------|--------|---------|
| Samples per Pixel          | 0028,0002 | US        | 1  | ALWAYS            | AUTO   |         |
| Photometric Interpretation | 0028,0004 | CS        | MONOCHROME2: Upon import, only images with a photometric interpretation MONOCHROME1 or MONOCHROME2 are accepted. | ALWAYS            | AUTO   |         |
| Rows                       | 0028,0010 | US        | 2048, 1024 or 512  | ALWAYS            | AUTO   |         |
| Columns                    | 0028,0011 | US        | 2048, 1024 or 512  | ALWAYS            | AUTO   |         |
| Bits Allocated             | 0028,0100 | US        | 16 or 8  | ALWAYS            | AUTO   |         |
| Bits Stored                | 0028,0101 | US        | 8 and 12   | ALWAYS            | AUTO   |         |
| High Bit                   | 0028,0102 | US        |  | ALWAYS            | AUTO   |         |
| Pixel Representation       | 0028,0103 | US        | 0000H  | ALWAYS            | AUTO   |         |
| Pixel Data                 | 7FE0,0010 | OW/O<br>B |  | ANAP              | AUTO   |         |

Table 119: Contrast/Bolus Module

| Attribute Name       | Tag       | VR | Value | Presence of Value | Source | Comment |
|----------------------|-----------|----|-------|-------------------|--------|---------|
| Contrast/Bolus Agent | 0018,0010 | LO |       | VNAP              | AUTO   |         |

**Table 120: Cine Module**

| Attribute Name                 | Tag       | VR | Value | Presence of Value | Source | Comment |
|--------------------------------|-----------|----|-------|-------------------|--------|---------|
| Recommended Display Frame Rate | 0008,2144 | IS |       | ANAP              | AUTO   |         |
| Cine Rate                      | 0018,0040 | IS |       | ANAP              | AUTO   |         |
| Frame Time                     | 0018,1063 | DS |       | ALWAYS            | AUTO   |         |
| Frame Delay                    | 0018,1066 | DS |       | ANAP              | AUTO   |         |

**Table 121: Multi-Frame Module**

| Attribute Name          | Tag       | VR | Value                       | Presence of Value | Source | Comment |
|-------------------------|-----------|----|-----------------------------|-------------------|--------|---------|
| Number of Frames        | 0028,0008 | IS |                             | ALWAYS            | AUTO   |         |
| Frame Increment Pointer | 0028,0009 | AT | 0x00181065 or<br>0x00181063 | ALWAYS            | AUTO   |         |

**Table 122: Display Shutter Module**

| Attribute Name                | Tag       | VR | Value       | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|-------------|-------------------|--------|---------|
| Shutter Shape                 | 0018,1600 | CS | RECTANGULAR | ALWAYS            | AUTO   |         |
| Shutter Left Vertical Edge    | 0018,1602 | IS |             | ALWAYS            | AUTO   |         |
| Shutter Right Vertical Edge   | 0018,1604 | IS |             | ALWAYS            | AUTO   |         |
| Shutter Upper Horizontal Edge | 0018,1606 | IS |             | ALWAYS            | AUTO   |         |
| Shutter Lower Horizontal Edge | 0018,1608 | IS |             | ALWAYS            | AUTO   |         |

**Table 123: X-Ray Image Module**

| Attribute Name               | Tag       | VR | Value   | Presence of Value | Source | Comment |
|------------------------------|-----------|----|---|-------------------|--------|---------|
| Image Type                   | 0008,0008 | CS |   | ALWAYS            | AUTO   |         |
| Referenced Image Sequence    | 0008,1140 | SQ |   | ALWAYS            | AUTO   |         |
| >Referenced SOP Class UID    | 0008,1150 | UI |   | ALWAYS            | AUTO   |         |
| >Referenced SOP Instance UID | 0008,1155 | UI |   | ALWAYS            | AUTO   |         |
| Samples per Pixel            | 0028,0002 | US | 1   | ALWAYS            | AUTO   |         |
| Photometric Interpretation   | 0028,0004 | CS | MONOCHROME2 Upon import, only images with a photometric interpretation MONOCHROME1 or MONOCHROME2 are accepted. | ALWAYS            | AUTO   |         |
| Frame Increment Pointer      | 0028,0009 | AT |   | ALWAYS            |        |         |
| Bits Allocated               | 0028,0100 | US | 16 or 8   | ALWAYS            | AUTO   |         |
| Bits Stored                  | 0028,0101 | US | 10 or 8. Upon import, only images with bits stored equal to 8 or 10 bits are accepted.                          | ALWAYS            | AUTO   |         |
| High Bit                     | 0028,0102 | US |   | ALWAYS            | AUTO   |         |
| Pixel Representation         | 0028,0103 | US | 0000H   | ALWAYS            | AUTO   |         |





| Attribute Name               | Tag       | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|--------|---------|
| Pixel Intensity Relationship | 0028,1040 | CS |       | ALWAYS            | AUTO   |         |
| Lossy Image Compression      | 0028,2110 | CS |       | ALWAYS            | AUTO   |         |

**Table 124: X-Ray Acquisition Module**

| Attribute Name                | Tag       | VR | Value   | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|---|-------------------|--------|---------|
| KVP                           | 0018,0060 | DS |   | ALWAYS            | AUTO   |         |
| Exposure Time                 | 0018,1150 | IS | Only sent if Exposure (0018, 1152) is not sent.   | ANAP              | AUTO   |         |
| X-Ray Tube Current            | 0018,1151 | IS | Only sent if Exposure (0018, 1152) is not sent.   | ANAP              | AUTO   |         |
| Average Pulse Width           | 0018,1154 | DS |   | ANAP              | AUTO   |         |
| Exposure                      | 0018,1152 | IS | Only sent if Exposure Time (0018, 1150) and X-Ray Tube Current (0018, 1151) are not sent. | ANAP              | AUTO   |         |
| Radiation Setting             | 0018,1155 | CS |   | ALWAYS            | AUTO   |         |
| Imager Pixel Spacing          | 0018,1164 | DS |   | ANAP              | AUTO   |         |
| X-Ray Tube Current in $\mu$ A | 0018,8151 | DS |   | ANAP              |        |         |

**Table 125: X-Ray Table Module**

| Attribute Name               | Tag       | VR | Value  | Presence of Value | Source | Comment |
|------------------------------|-----------|----|--|-------------------|--------|---------|
| Table Motion                 | 0018,1134 | CS | STATIC, DYNAMIC  | VNAP              | AUTO   |         |
| Table Vertical Increment     | 0018,1135 | DS | Incremental change (per frame) in vertical position relatively to the first frame of Multi-Frame image in mm     | ANAP              | AUTO   |         |
| Table Lateral Increment      | 0018,1136 | DS | Incremental change (per frame) in lateral position relatively to the first frame of Multi-Frame image in mm      | ANAP              | AUTO   |         |
| Table Longitudinal Increment | 0018,1137 | DS | Incremental change (per frame) in longitudinal position relatively to the first frame of Multi-Frame image in mm | ANAP              | AUTO   |         |
| Table Angle                  | 0018,1138 | DS | Angle of table plane in  | ANAP              | AUTO   |         |

|  |  |  |   |  |  |  |
|--|--|--|---|--|--|--|
|  |  |  | degrees relative to horizontal plane (gravity plane). Positive values indicate that the head of the table is upwards. |  |  |  |
|--|--|--|---|--|--|--|

**Table 126: XA Positioner Module**

| Attribute Name                       | Tag       | VR | Value   | Presence of Value | Source | Comment |
|--------------------------------------|-----------|----|---|-------------------|--------|---------|
| Distance Source to Detector          | 0018,1110 | DS |   | ALWAYS            | AUTO   |         |
| Distance Source to Patient           | 0018,1111 | DS |   | ALWAYS            | AUTO   |         |
| Positioner Motion                    | 0018,1500 | CS | STATIC, DYNAMIC   | ALWAYS            | AUTO   |         |
| Positioner Primary Angle             | 0018,1510 | DS |   | ALWAYS            | AUTO   |         |
| Positioner Secondary Angle           | 0018,1511 | DS |   | ALWAYS            | AUTO   |         |
| Positioner Primary Angle Increment   | 0018,1520 | DS | An array that contains the Positioner Primary Angle Increments between the n-th frame and the previous frame for a Multi-frame image.   | ANAP              | AUTO   |         |
| Positioner Secondary Angle Increment | 0018,1521 | DS | An array that contains the Positioner Secondary Angle Increments between the n-th frame and the previous frame for a Multi-frame image. | ANAP              | AUTO   |         |

**Table 127: DX Detector Module**

| Attribute Name       | Tag       | VR | Value | Presence of Value | Source | Comment |
|----------------------|-----------|----|-------|-------------------|--------|---------|
| Imager Pixel Spacing | 0018,1164 | DS |       | ALWAYS            | AUTO   |         |
| Detector Type        | 0018,7004 | CS |       | VNAP              | AUTO   |         |

**Table 128: Modality LUT Module**

| Attribute Name        | Tag       | VR    | Value  | Presence of Value | Source | Comment |
|-----------------------|-----------|-------|--|-------------------|--------|---------|
| Modality LUT Sequence | 0028,3000 | SQ    |  | ANAP              | AUTO   |         |
| >LUT Descriptor       | 0028,3002 | US/SS | If 8 bits stored:<br>[value1]=256<br>[value 2]=0<br>[value 3]=8 If 10 bits stored<br>[value 1]=1024<br>[value 2]=0<br>[value 3]=10 | ANAP              | AUTO   |         |
| >Modality LUT Type    | 0028,3004 | LO    | US   | ANAP              | AUTO   |         |
| >LUT Data             | 0028,3006 | UN    |  | ANAP              | AUTO   |         |

**Note:** Modality LUT and pixel intensity “LOG” only when X-Ray Angiographic images are sent with unprocessed pixel data  
-No Modality LUT and pixel intensity “LIN” when X-Ray Angiographic images are sent with processed pixel data.

**Table 129: VOI LUT Module**

| Attribute Name | Tag       | VR | Value | Presence of Value | Source | Comment |
|----------------|-----------|----|-------|-------------------|--------|---------|
| Window Center  | 0028,1050 | DS |       | ALWAYS            | AUTO   |         |
| Window Width   | 0028,1051 | DS |       | ALWAYS            | AUTO   |         |

**Table 130: Curve Module**

| Attribute Name            | Tag       | VR | Value     | Presence of Value | Source | Comment |
|---------------------------|-----------|----|-----------|-------------------|--------|---------|
| Curve Dimensions          | 5000,0005 | US |           | ALWAYS            | AUTO   |         |
| Number Of Points          | 5000,0010 | US |           | ALWAYS            | AUTO   |         |
| Type Of Data              | 5000,0020 | CS |           | ALWAYS            | AUTO   |         |
| Axis Units                | 5000,0030 | SH | DPPS\NONE | ALWAYS            | AUTO   |         |
| Data Value Representation | 5000,0103 | US |           | ALWAYS            | AUTO   |         |
| Minimum Coordinate Value  | 5000,0104 | US |           | ALWAYS            | AUTO   |         |
| Maximum Coordinate Value  | 5000,0105 | US |           | ALWAYS            | AUTO   |         |
| Curve Data Descriptor     | 5000,0110 | US |           | ALWAYS            | AUTO   |         |
| Coordinate Start Value    | 5000,0112 | US |           | ALWAYS            | AUTO   |         |
| Coordinate Step Value     | 5000,0114 | US |           | VNAP              | AUTO   |         |
| Curve Data                | 5000,3000 | OW |           | ALWAYS            | AUTO   |         |

**Table 131: SOP Common Module**

| Attribute Name         | Tag       | VR | Value                        | Presence of Value | Source | Comment |
|------------------------|-----------|----|------------------------------|-------------------|--------|---------|
| Specific Character Set | 0008,0005 | CS | ISO_IR 100                   | ANAP              | AUTO   |         |
| SOP Class UID          | 0008,0016 | UI | 1.2.840.10008.5.1.4.1.1.12.1 | ALWAYS            | AUTO   |         |
| SOP Instance UID       | 0008,0018 | UI | Generated by device          | ALWAYS            | AUTO   |         |

**8.1.1.5. X-Ray Radiation Dose SR**

**Table 132: SOP Class Modules**

| Information Entity | Module                            | Presence |
|--------------------|-----------------------------------|----------|
| Patient            | Patient Module                    | Always   |
| Study              | General Study Module              | Always   |
|                    | Patient Study Module              | Always   |
| Series             | SR Document Series Module         | Always   |
| Equipment          | General Equipment Module          | Always   |
|                    | Enhanced General Equipment Module | Always   |
| Image              | SR Document General Module        | Always   |
|                    | SR Document Content Module        | Always   |



SOP Common Module

Always

**Table 133: Patient Module**

| Attribute Name               | Tag       | VR | Value      | Presence of Value | Source       | Comment |
|------------------------------|-----------|----|------------|-------------------|--------------|---------|
| Referenced Patient Sequence  | 0008,1120 | SQ |            | ANAP              | MWL          |         |
| >Referenced SOP Class UID    | 0008,1150 | UI |            | ALWAYS            | MWL          |         |
| >Referenced SOP Instance UID | 0008,1155 | UI |            | ALWAYS            | MWL          |         |
| Patient's Name               | 0010,0010 | PN |            | ALWAYS            | MWL,<br>USER |         |
| Patient ID                   | 0010,0020 | LO |            | VNAP              | MWL,<br>USER |         |
| Patient's Birth Date         | 0010,0030 | DA | <yyyymmdd> | VNAP              | MWL,<br>USER |         |
| Patient's Sex                | 0010,0040 | CS |            | VNAP              | MWL,<br>USER |         |
| Other Patient IDs            | 0010,1000 | LO |            | ANAP              | MWL,<br>USER |         |
| Ethnic Group                 | 0010,2160 | SH |            | ANAP              | MWL,<br>USER |         |
| Patient Comments             | 0010,4000 | LT |            | ANAP              | MWL,<br>USER |         |

**Table 134: General Study Module**

| Attribute Name             | Tag       | VR | Value  | Presence of Value | Source                | Comment |
|----------------------------|-----------|----|--|-------------------|-----------------------|---------|
| Study Date                 | 0008,0020 | DA | <yyyymmdd>   | ALWAYS            | MWL,AUTO              |         |
| Study Time                 | 0008,0030 | TM | <hhmmss>   | ALWAYS            | MWL,AUTO              |         |
| Accession Number           | 0008,0050 | SH |  | ALWAYS            | AUTO,<br>MWL,<br>USER |         |
| Referring Physician's Name | 0008,0090 | PN | Patient's referring physician  | VNAP              | MWL                   |         |
| Study Description          | 0008,1030 | LO | Based on configuration<br>Study Description is:<br>- not exported<br>- based on schedule<br>procedure step description | ANAP              | AUTO,<br>MWL          |         |

| Attribute Name               | Tag       | VR | Value   | Presence of Value | Source                | Comment |
|------------------------------|-----------|----|---|-------------------|-----------------------|---------|
|                              |           |    | (WLM) - based on requested procedure step description<br>(WLM) - internal generated performed procedure description.                                      |                   |                       |         |
| Procedure Code Sequence      | 0008,1032 | SQ |   | ANAP              | MWL,AUTO              |         |
| Referenced Study Sequence    | 0008,1110 | SQ |   | ANAP              | MWL                   |         |
| >Referenced SOP Class UID    | 0008,1150 | UI |   | ALWAYS            | MWL                   |         |
| >Referenced SOP Instance UID | 0008,1155 | UI |   | ALWAYS            | MWL                   |         |
| Study Instance UID           | 0020,000D | UI |   | ALWAYS            | AUTO                  |         |
| Study ID                     | 0020,0010 | SH | In case the Study ID is empty the accession number will be assigned. In case Study ID and Accession Number are the same, the Study ID will be left empty. | VNAP              | AUTO,<br>MWL,<br>USER |         |

Table 135: Patient Study Module

| Attribute Name             | Tag       | VR | Value  | Presence of Value | Source       | Comment |
|----------------------------|-----------|----|--|-------------------|--------------|---------|
| Patient's Weight           | 0010,1030 | DS | In kilograms. When received from the MWL SCP, the value can still be modified. | ANAP              | MWL,<br>USER |         |
| Patient's Size             | 0010,1020 | DS |  | ANAP              | MWL,USER     |         |
| Additional Patient History | 0010,21B0 | LT |  | ANAP              | MWL          |         |

Table 136: SR Document Series Module

| Attribute Name     | Tag       | VR | Value                        | Presence of Value | Source | Comment |
|--------------------|-----------|----|------------------------------|-------------------|--------|---------|
| Series Date        | 0008,0021 | DA |                              | ANAP              | AUTO   |         |
| Series Time        | 0008,0031 | TM |                              | ANAP              | AUTO   |         |
| Modality           | 0008,0060 | CS | SR                           | ALWAYS            | AUTO   |         |
| Series Description | 0008,103E | LO | "Radiation Dose Information" | ANAP              | AUTO   |         |



|  |           |    |  |        |      |  |
|--|-----------|----|--|--------|------|--|
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ |  | VNAP   | AUTO |  |
| >Referenced SOP Class UID                    | 0008,1150 | UI |  | ALWAYS | AUTO |  |
| >Referenced SOP Instance UID                 | 0008,1155 | UI |  | ALWAYS | AUTO |  |
| Series Instance UID                          | 0020,000E | UI |  | ALWAYS | AUTO |  |
| Series Number                                | 0020,0011 | IS |  | ALWAYS | AUTO |  |

**Table 137: General Equipment Module**

| Attribute Name                | Tag       | VR | Value   | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|---------|-------------------|--------|---------|
| Manufacturer                  | 0008,0070 | LO | Philips | VNAP              | AUTO   |         |
| Institution Name              | 0008,0080 | LO |         | ANAP              | CONFIG |         |
| Station Name                  | 0008,1010 | SH |         | ANAP              | CONFIG |         |
| Institutional Department Name | 0008,1040 | LO |         | ANAP              | CONFIG |         |
| Manufacturer's Model Name     | 0008,1090 | LO | Azurion | ANAP              | AUTO   |         |
| Device Serial Number          | 0018,1000 | LO |         | ANAP              | AUTO   |         |
| Software Versions             | 0018,1020 | LO | 2.0.0   | ANAP              | AUTO   |         |

**Table 138: Enhanced General Equipment Module**

| Attribute Name            | Tag       | VR | Value   | Presence of Value | Source | Comment |
|---------------------------|-----------|----|---------|-------------------|--------|---------|
| Manufacturer              | 0008,0070 | LO | Philips | ALWAYS            | AUTO   |         |
| Manufacturer's Model Name | 0008,1090 | LO | Azurion | ALWAYS            | AUTO   |         |
| Device Serial Number      | 0018,1000 | LO |         | ALWAYS            | AUTO   |         |
| Software Versions         | 0018,1020 | LO | 2.0.0   | ALWAYS            | AUTO   |         |

**Table 139: SR Document General Module**

| Attribute Name                      | Tag       | VR | Value  | Presence of Value | Source | Comment                                   |
|-------------------------------------|-----------|----|--|-------------------|--------|---|
| Content Date                        | 0008,0023 | DA |  | ALWAYS            | AUTO   |   |
| Content Time                        | 0008,0033 | TM |  | ALWAYS            | AUTO   |   |
| Referenced Instance Sequence        | 0008,114A | SQ |  | ANAP              | AUTO   |   |
| >Referenced SOP Class UID           | 0008,1150 | UI | 1.2.840.10008.5.1.4.1.1.12.1                                       | ALWAYS            | AUTO   |   |
| >Referenced SOP Instance UID        | 0008,1155 | UI | SOP Instance UID of AcquiredRun as it will be sent to the Archive. | ALWAYS            | AUTO   | Note: State of Archive should be checked. |
| >Purpose of Reference Code Sequence | 0040,A170 | SQ |  | ALWAYS            |        |   |
| >>Code Value                        | 0008,0100 | SH |  | ALWAYS            |        |   |
| >>Coding Scheme Designator          | 0008,0102 | SH |  | ALWAYS            |        |   |

| Attribute Name                                | Tag       | VR | Value   | Presence of Value | Source | Comment |
|---|-----------|----|---|-------------------|--------|---------|
| >>Code Meaning                                | 0008,0104 | LO |   | ALWAYS            |        |         |
| Instance Number                               | 0020,0013 | IS |   | ALWAYS            | AUTO   |         |
| Referenced Request Sequence                   | 0040,A370 | SQ |   | ANAP              | AUTO   |         |
| >Study Instance UID                           | 0020,000D | UI |   | ALWAYS            | MWL    |         |
| >Requested Procedure Description              | 0032,1060 | LO |   | VNAP              |        |         |
| >Requested Procedure ID                       | 0040,1001 | SH |   | VNAP              |        |         |
| >Reason for the Requested Procedure           | 0040,1002 | LO |   | ANAP              |        |         |
| >Placer Order Number/Imaging Service Request  | 0040,2016 | LO |   | VNAP              |        |         |
| >Filler Order Number/Imaging Service Request  | 0040,2017 | LO |   | VNAP              |        |         |
| Performed Procedure Code Sequence             | 0040,A372 | SQ |   | VNAP              |        |         |
| Current Requested Procedure Evidence Sequence | 0040,A375 | SQ |   | ANAP              |        |         |
| >Referenced Series Sequence                   | 0008,1115 | SQ |   | ALWAYS            |        |         |
| >>Referenced SOP Sequence                     | 0008,1199 | SQ |   | ALWAYS            |        |         |
| >>>Referenced SOP Class UID                   | 0008,1150 | UI |   | ALWAYS            |        |         |
| >>>Referenced SOP Instance UID                | 0008,1155 | UI |   | ALWAYS            |        |         |
| >>Series Instance UID                         | 0020,000E | UI |   | ALWAYS            |        |         |
| >Study Instance UID                           | 0020,000D | UI |   | ALWAYS            |        |         |
| Completion Flag                               | 0040,A491 | CS | COMPLETE  | ALWAYS            | AUTO   |         |
| Completion Flag Description                   | 0040,A492 | LO | "Complete X-Ray Radiation Dose Structured Report" | ANAP              | AUTO   |         |
| Verification Flag                             | 0040,A493 | CS | UNVERIFIED  | ALWAYS            | AUTO   |         |

Table 140: SR Document Content Module

| Attribute Name                | Tag       | VR | Value     | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|-----------|-------------------|--------|---------|
| Value Type                    | 0040,A040 | CS | CONTAINER | ALWAYS            | AUTO   |         |
| Content Template Sequence     | 0040,A504 | SQ |           | ANAP              |        |         |
| Content Sequence              | 0040,A730 | SQ |           | ANAP              |        |         |
| >Referenced SOP Sequence      | 0008,1199 | SQ |           | ALWAYS            |        |         |
| >>Referenced SOP Class UID    | 0008,1150 | UI |           | ALWAYS            |        |         |
| >>Referenced SOP Instance UID | 0008,1155 | UI |           | ALWAYS            |        |         |
| >Relationship Type            | 0040,A010 | CS |           | ALWAYS            |        |         |
| >Value Type                   | 0040,A040 | CS |           | ALWAYS            |        |         |
| >Concept Name Code Sequence   | 0040,A043 | SQ |           | ANAP              |        |         |
| >>Code Value                  | 0008,0100 | SH |           | ALWAYS            |        |         |
| >>Coding Scheme Designator    | 0008,0102 | SH |           | ALWAYS            |        |         |
| >>Code Meaning                | 0008,0104 | LO |           | ALWAYS            |        |         |
| >Temporal Range Type          | 0040,A130 | CS |           | ALWAYS            |        |         |
| >Concept Code Sequence        | 0040,A168 | SQ |           | ALWAYS            |        |         |
| >>Code Value                  | 0008,0100 | SH |           | ALWAYS            |        |         |





|                                    |           |    |  |        |  |  |
|------------------------------------|-----------|----|--|--------|--|--|
| >>Coding Scheme Designator         | 0008,0102 | SH |  | ALWAYS |  |  |
| >>Code Meaning                     | 0008,0104 | LO |  | ALWAYS |  |  |
| >Graphic Data                      | 0070,0022 | FL |  | ALWAYS |  |  |
| >Graphic Type                      | 0070,0023 | CS |  | ALWAYS |  |  |
| >Referenced Frame of Reference UID | 3006,0024 | UI |  | ALWAYS |  |  |

Table 141: SOP Common Module

| Attribute Name         | Tag       | VR | Value                         | Presence of Value | Source | Comment |
|------------------------|-----------|----|-------------------------------|-------------------|--------|---------|
| Specific Character Set | 0008,0005 | CS | ISO_IR 100                    | ANAP              | AUTO   |         |
| Instance Creation Date | 0008,0012 | DA |                               | ANAP              | AUTO   |         |
| Instance Creation Time | 0008,0013 | TM |                               | ANAP              | AUTO   |         |
| SOP Class UID          | 0008,0016 | UI | 1.2.840.10008.5.1.4.1.1.88.67 | ALWAYS            | AUTO   |         |
| SOP Instance UID       | 0008,0018 | UI |                               | ALWAYS            | AUTO   |         |
| Instance Number        | 0020,0013 | IS |                               | ANAP              | AUTO   |         |

### 8.1.2. Attribute Mapping

Not applicable

### 8.1.3. Coerced/Modified Fields

Not applicable

## 8.2. Data Dictionary of Private Attributes

Not applicable

## 8.3. Coded Terminology and Templates

Not applicable

### 8.3.1. Context Group

Not applicable.

### 8.3.2. Template Specifications

#### X-RAY RADIATION DOSE STRUCTURED REPORT IOD TEMPLATES

The templates that comply SCPe the X-Ray Radiation Dose Structured Report are interconnected as indicated in the figure below:




Figure 18: X-Ray Radiation Dose Structured Report IOD Template Structure

This section describes the content of all the templates used in the X-Ray Radiation Dose Reporting SR.

Table 142: Used Templates for X-Ray Radiation Dose Reporting

| Template Name   | Template ID |
|---|-------------|
| Projection X-Ray Radiation Dose                               | TID 10001   |
| Accumulated X-Ray Dose  | TID 10002   |
| Irradiation Event X-Ray Data                                  | TID 10003   |
| Irradiation Event X-Ray Source Data                           | TID 10003B  |
| Irradiation Event X-Ray Mechanical Data                       | TID 10003C  |
| Accumulated Fluoroscopy and Acquisition Projection X-Ray Dose | TID 10004   |
| Accumulated Total Projection Radiography Dose                 | TID 10007   |
| Observer Context  | TID 1002    |
| Device Observer Identifying Attributes                        | TID 1004    |

8.3.2.1. TID 10001 Projection X-Ray Radiation Dose

Table 143: Projection X-Ray Radiation Dose

| NL | Relation with Parent | Concept Name                            | VT        | VM  | Presence of Value | Value   |
|----|----------------------|---|-----------|-----|-------------------|---|
|    |                      | X-Ray Radiation Dose Report(113701,DCM) | CONTAINER | 1   | ALWAYS            |   |
| >  | HAS CONCEPT MOD      | Procedure reported                      | CODE      | 1   | ALWAYS            | Projection X-Ray                              |
| >> | HAS CONCEPT MOD      | Has Intent                              | CODE      | 1   | ALWAYS            | Combined Diagnostic and Therapeutic Procedure |
| >  |                      | DTID 1002 "Observer Context"            | INCLUDE   | 1-n | ALWAYS            |   |
| >  | HAS OBS CONTEXT      | Scope of Accumulation (DCM, 113705)     | CODE      | 1   | ALWAYS            | 113016, DCM, Performed Procedure Step         |

|    |                 |   |         |     |        |  |
|----|-----------------|---|---------|-----|--------|--|
| >> | HAS PROPER TIES | Performed Procedure Step SOP Instance UID (DCM, 121126) | UIDREF  | 1   | ALWAYS |  |
| >  | CONTAINS        | DTID 10002 "Accumulated X-Ray Dose"                     | INCLUDE | 1   | MC     |  |
| >  | CONTAINS        | DTID 10003 "Irradiation Event X-Ray Data"               | INCLUDE | 1-n | ALWAYS |  |
| >  | CONTAINS        | Comment (DCM, 121106)                                   | TEXT    | 1   | ALWAYS | X-Ray Radiation Dose Structured Report related |
| >  | CONTAINS        | Source of Dose Information (DCM, 113854)                | CODE    | 1-n | ALWAYS | (113856, DCM, Automated Data Collection)       |

**8.3.2.2. TID 10002 Accumulated X-Ray Dose**

**Table 144: Accumulated X-Ray Dose**

| NL | Relation with Parent | Concept Name                                | VT        | VM | Presence of Value | Value                        |
|----|----------------------|---|-----------|----|-------------------|------------------------------|
|    |                      | Accumulated X-Ray Dose Data (DCM, 113702)   | CONTAINER | 1  | ALWAYS            |                              |
| >  | HAS CONCEPT MOD      | Acquisition Plane (DCM, 113764)             | CODE      | 1  | ALWAYS            | 113620, DCM, Plane A)        |
| >  | CONTAINS             | Calibration (DCM, 122505)                   | CONTAINER | 1  | ALWAYS            | DCMR, 122505                 |
| >> | CONTAINS             | Calibration Factor (DCM, 122322)            | NUM       | 1  | ALWAYS            | 1 Units: (1, UCUM, no units) |
| >> | CONTAINS             | Calibration DateTime (DCM, 113723)          | DATETIME  | 1  | ALWAYS            |                              |
| >> | CONTAINS             | Calibration Responsible Party (DCM, 113724) | TEXT      | 1  | ALWAYS            |                              |
| >> | CONTAINS             | Calibration Protocol(DCM,113720)            | TEXT      | 1  | ALWAYS            |                              |
| >> | HAS CONCEPT MOD      | Dose Measurement Device (DCM, 113794)       | CODE      | 1  | ALWAYS            | Dosimeter (SRT, A-2C090)     |
| >> | CONTAINS             | Calibration Uncertainty (DCM, 113763)       | NUM       | 1  | ALWAYS            | 1 Units: (% , UCUM, Percent) |

| NL | Relation with Parent | Concept Name                                   | VT  | VM | Presence of Value | Value |
|----|----------------------|--|-----|----|-------------------|-------|
| >  | CONTAINS             | Distance Source to Reference Point(DCM,113737) | NUM | 1  | OPTIONAL          |       |
| >  | CONTAINS             | Height of System (99PHI-IXR-XPER, 001)         | NUM | 1  | ALWAYS            |       |
| >  | CONTAINS             | Focal Spot to ISO Center (99PHI-IXR-XPER, 002) | NUM | 1  | ALWAYS            |       |

**8.3.2.3. TID 10003 Irradiation Event X-Ray Data**

**Table 145: Irradiation Event X-Ray Data**

| NL | Relation with Parent | Concept Name                               | VT        | VM | Presence of Value | Value |
|----|----------------------|--|-----------|----|-------------------|-------|
|    | CONTAINS             | Irradiation Event X-Ray Data (DCM, 113706) | CONTAINER | 1  | ALWAYS            |       |

| NL   | Relation with Parent | Concept Name  | VT        | VM | Presence of Value | Value                                 |
|------|----------------------|---|-----------|----|-------------------|---------------------------------------|
| >    | HAS CONCEPT MOD      | Acquisition Plane (DCM, 113764)                         | CODE      | 1  | ALWAYS            | (113620, DCM, Plane A)                |
| >    | CONTAINS             | DateTime Started (DCM, 111526)                          | DATETIME  | 1  | ALWAYS            |                                       |
| >    | CONTAINS             | Irradiation Event Type (DCM, 113721)                    | CODE      | 1  | ALWAYS            | (113611, DCM, Stationary Acquisition) |
| >    | CONTAINS             | Acquisition Protocol(DCM, 125203)                       | TEXT      | 1  | ALWAYS            |                                       |
| >    | CONTAINS             | Pulse Rate(DCM,113791)                                  | NUM       | 1  | CONDITIONAL       |                                       |
| >    | CONTAINS             | Irradiation Event UID (DCM, 113769)                     | UIDREF    | 1  | ALWAYS            |                                       |
| >    | CONTAINS             | Dose Area Product (DCM, 122130)                         | NUM       | 1  | ALWAYS            | Units: (Gy.m2, UCUM, Gy.m2)           |
| >>>  | CONTAINS             | Collimated Field Height (DCM, 113788)                   | NUM       | 1  | ALWAYS            |                                       |
| >>>  | CONTAINS             | Collimated Field Width (DCM, 113789)                    | NUM       | 1  | ALWAYS            |                                       |
| >    | CONTAINS             | Patient Table Relationship (DCM, 113745)                | CODE      | 1  | ALWAYS            | (F-10470, SRT, headfirst)             |
| >    | CONTAINS             | Patient Orientation (DCM, 113743)                       | CODE      | 1  | ALWAYS            | (F-10450, SRT, recumbent)             |
| >>   | HAS CONCEPT MOD      | Patient Orientation Modifier (DCM, 113744)              | CODE      | 1  | ALWAYS            | (F-10340, SRT, supine)                |
| >    | CONTAINS             | Target Region (DCM, 123014)                             | CODE      | 1  | ALWAYS            | (T-D3000, SRT, Chest)                 |
| >>>  | CONTAINS             | Patient Equivalent Thickness(DCM,111638)                | NUM       | 1  | OPTIONAL          |                                       |
| >    | CONTAINS             | Comment (DCM, 121106)                                   | TEXT      | 1  | ALWAYS            |                                       |
| >    | CONTAINS             | DTID 10003B "Irradiation Event X-Ray Source Data"       | INCLUDE   | 1  |                   |                                       |
| >    | CONTAINS             | DTID 10003C "Irradiation Event X-Ray MechanicalData"    | INCLUDE   | 1  |                   |                                       |
| >>>  | CONTAINS             | Subimages per Frame (99PHI-IXR-XPEN, 004)               | NUM       | 1  | ALWAYS            |                                       |
| >>>  | CONTAINS             | Wedges and Shutters (99PHI-IXR-XPEN, 005)               | CONTAINER | 1  | ALWAYS            | DCMR, 005                             |
| >>>> | CONTAINS             | Bottom Shutter (99PHI-IXR-XPEN, 006)                    | NUM       | 1  | ALWAYS            |                                       |
| >>>> | CONTAINS             | Left Shutter (99PHI-IXR-XPEN, 007)                      | NUM       | 1  | ALWAYS            |                                       |
| >>>> | CONTAINS             | Right Shutter (99PHI-IXR-XPEN, 008)                     | NUM       | 1  | ALWAYS            |                                       |
| >>>> | CONTAINS             | Top Shutter (99PHI-IXR-XPEN, 009)                       | NUM       | 1  | ALWAYS            |                                       |
| >>>> | CONTAINS             | Distance Wedge 1 (99PHI-IXR-XPEN, 010)                  | NUM       | 1  | ALWAYS            |                                       |
| >>>> | CONTAINS             | Distance Wedge 2 (99PHI-IXR-XPEN, 011)                  | NUM       | 1  | ALWAYS            |                                       |
| >>>> | CONTAINS             | Angle Wedge 1 (99PHI-IXR-XPEN, 012)                     | NUM       | 1  | ALWAYS            |                                       |
| >>>> | CONTAINS             | Angle Wedge 2 (99PHI-IXR-XPEN, 013)                     | NUM       | 1  | ALWAYS            |                                       |
| >>>  | CONTAINS             | Beam Position (99PHI-IXR-XPEN, 014)                     | CONTAINER | 1  | ALWAYS            |                                       |
| >>>> | CONTAINS             | Longitudinal Beam Position (99PHI-IXR-XPEN, 015)        | NUM       | 1  | ALWAYS            |                                       |
| >>>> | CONTAINS             | Lateral Beam Position (99PHI-IXR-XPEN, 016)             | NUM       | 1  | ALWAYS            |                                       |
| >>>> | CONTAINS             | Beam Angle (99PHI-IXR-XPEN, 017)                        | NUM       | 1  | ALWAYS            |                                       |
| >>>  | CONTAINS             | Final Distance Source to Detector (99PHI-IXR-XPEN, 018) | NUM       | 1  | ALWAYS            |                                       |

| NL   | Relation with Parent | Concept Name                                   | VT        | VM | Presence of Value | Value     |
|------|----------------------|--|-----------|----|-------------------|-----------|
| >>>  | CONTAINS             | Final Table Cradle Angle (99PHI-IXR-XPER, 019) | NUM       | 1  | ALWAYS            |           |
| >>>  | CONTAINS             | Table Height Position (99PHI-IXR-XPER, 021)    | NUM       | 1  | ALWAYS            |           |
| >>>  | CONTAINS             | Final Table Tilt Angle (99PHI-IXR-XPER, 022)   | NUM       | 1  | ALWAYS            |           |
| >>>  | CONTAINS             | Detector Field Size (99PHI-IXR-XPER, 023)      | CONTAINER | 1  | ALWAYS            | DCMR, 023 |
| >>>> | CONTAINS             | X Side (99PHI-IXR-XPER, 024)                   | NUM       | 1  | ALWAYS            |           |
| >>>> | CONTAINS             | Y Side (99PHI-IXR-XPER, 025)                   | NUM       | 1  | ALWAYS            |           |
| >>>  | CONTAINS             | Object Thickness (99PHI-IXR-XPER, 026)         | NUM       | 1  | ALWAYS            |           |
| >>>  | CONTAINS             | Application Name (99PHI-IXR-XPER, 029)         | TEXT      | 1  | ALWAYS            |           |
| >>>  | CONTAINS             | Fluoro Flavour (99PHI-IXR-XPER, 030)           | TEXT      | 1  | ALWAYS            |           |

**8.3.2.4. TID 10003B Irradiation Event X-Ray Source Data**

**Table 146: Irradiation Event X-Ray Source Data**

| NL | Relation with Parent | Concept Name                                 | VT        | VM  | Presence of Value     | Value   |
|----|----------------------|--|-----------|-----|-----------------------|---|
| >  | CONTAINS             | Dose (RP) (DCM, 113738)                      | NUM       | 1   | ALWAYS                | Units: (Gy, UCUM, Gy)                         |
| >  | CONTAINS             | Reference Point Definition (DCM, 113780)     | TEXT      | 1   | ALWAYS                | 15cm below BeamIsocenter                      |
|    | CONTAINS             | Number of Pulses (DCM, 113768)               | NUM       | 1   | ALWAYS                | Units: (1, UCUM, no units)                    |
|    | CONTAINS             | Pulse Width (DCM, 113793)                    | NUM       | 1   | ALWAYS                | Units: (ms, UCUM, ms)                         |
|    | CONTAINS             | Irradiation Duration (DCM, 113742)           | NUM       | 1   | ALWAYS                | Units: (ms, UCUM, ms)                         |
|    | CONTAINS             | KVP (DCM, 113733)                            | NUM       | 1   | ALWAYS                | Units: (kV, UCUM, kV)                         |
|    | CONTAINS             | X-RayTube Current(DCM,113734)                | NUM       | 1-n | MANDATORY CONDITIONAL | Units: (mA, UCUM, mA)                         |
|    | CONTAINS             | Average X-Ray Tube Current (DCM, 113767)     | NUM       | 1   | ALWAYS                | Units: (mA, UCUM, mA)                         |
|    | CONTAINS             | Exposure (DCM, 113736)                       | NUM       | 1   | MANDATORY CONDITIONAL | Units: (uA.s, UCUM, uA.s)                     |
|    | CONTAINS             | Focal Spot Size(DCM,113766)                  | NUM       | 1   | OPTIONAL              | Units: (mm, UCUM, mm)                         |
|    |                      | X-Ray Filters (DCM, 113771)                  | CONTAINER | 1   | ALWAYS                | DCMR, 10007                                   |
| >  | CONTAINS             | X-Ray Filter Type (DCM, 113772)              | CODE      | 1   | ALWAYS                | (113650, DCM, Strip filter)                   |
| >  | CONTAINS             | X-Ray Filter Material (DCM, 113757)          | CODE      | 1   | ALWAYS                | (C-120F9, SRT, Aluminum or Aluminum compound) |
| >  | CONTAINS             | X-Ray Filter Thickness Minimum (DCM, 113758) | NUM       | 1   | ALWAYS                | Units: (mm, UCUM, mm)                         |
| >  | CONTAINS             | X-Ray Filter Thickness Maximum (DCM, 113773) | NUM       | 1   | ALWAYS                | Units: (mm, UCUM, mm)                         |

|  |  |                                     |     |   |        |                       |
|--|--|-------------------------------------|-----|---|--------|-----------------------|
|  |  | Collimated Field Area (DCM, 113790) | NUM | 1 | ALWAYS | Units: (m2, UCUM, m2) |
|--|--|-------------------------------------|-----|---|--------|-----------------------|

**8.3.2.5. TID 10003C Irradiation Event X-Ray Mechanical Data**

**Table 147: Irradiation Event X-Ray Mechanical Data**

| NL | Relation with Parent | Concept Name                                   | VT  | VM | Presence of Value | Value                   |
|----|----------------------|--|-----|----|-------------------|-------------------------|
|    |                      | Positioner Primary Angle(DCM,112011)           | NUM | 1  | UC                | Units: (deg, UCUM, deg) |
|    |                      | Positioner Secondary Angle(DCM,112012)         | NUM | 1  | UC                | Units: (deg, UCUM, deg) |
|    |                      | Table Head Tilt Angle (DCM, 113754)            | NUM | 1  | ALWAYS            | Units: (deg, UCUM, deg) |
|    |                      | Table Horizontal Rotation Angle(DCM,113755)    | NUM | 1  | ALWAYS            | Units: (deg, UCUM, deg) |
|    |                      | Table Cradle Tilt Angle (DCM, 113756)          | NUM | 1  | ALWAYS            | Units: (deg, UCUM, deg) |
|    |                      | Distance Source to Isocenter (DCM, 113748)     | NUM | 1  | ALWAYS            |                         |
|    |                      | Distance Source to Detector (DCM, 113750)      | NUM | 1  | ALWAYS            |                         |
|    |                      | Table Longitudinal Position (DCM, 113751)      | NUM | 1  | ALWAYS            |                         |
|    |                      | Table Lateral Position (DCM, 113752)           | NUM | 1  | ALWAYS            |                         |
|    |                      | Table Height Position (DCM, 113753)            | NUM | 1  | ALWAYS            |                         |
|    |                      | Distance Source to Reference Point(DCM,113737) | NUM | 1  | ALWAYS            |                         |

**8.3.2.6. TID 10004 Accumulated Fluoroscopy and Acquisition Projection X-Ray Dose**

**Table 148: Accumulated Fluoroscopy and Acquisition Projection X-Ray Dose**

| NL | Relation with Parent | Concept Name                                      | VT  | VM | Presence of Value     | Value                       |
|----|----------------------|---|-----|----|-----------------------|-----------------------------|
| >  | CONTAINS             | Acquisition Dose Area Product Total (DCM, 113727) | NUM | 1  | ALWAYS                | Units: (Gy.m2, UCUM, Gy.m2) |
| >  | CONTAINS             | Acquisition Dose (RP) Total (DCM, 113729)         | NUM | 1  | ALWAYS                | Units: (Gy, UCUM, Gy)       |
| >  | CONTAINS             | TotalFluoroTime,DCM,113730)                       | NUM | 1  | MANDATORY CONDITIONAL |                             |
| >  | CONTAINS             | Total Acquisition Time (DCM, 113855)              | NUM | 1  | ALWAYS                | Units: (s, UCUM, s)         |

**8.3.2.7. TID 10007 Accumulated Total Projection Radiography Dose**

**Table 149: Accumulated Total Projection Radiography Dose**

| NL | Relation with Parent | Concept Name                                      | VT  | VM | Presence of Value | Value                       |
|----|----------------------|---|-----|----|-------------------|-----------------------------|
| >  | CONTAINS             | Dose Area Product Total (DCM, 113722)             | NUM | 1  | ALWAYS            | Units: (Gy.m2, UCUM, Gy.m2) |
| >  | CONTAINS             | Dose (RP) Total (DCM, 113725)                     | NUM | 1  | ALWAYS            | Units: (Gy, UCUM, Gy)       |
| >  | CONTAINS             | Total Number of Radiographic Frames (DCM, 113731) | NUM | 1  | ALWAYS            | Units: (1, UCUM, no units)  |



|   |          |  |      |   |        |                          |
|---|----------|--|------|---|--------|--------------------------|
| > | CONTAINS | Reference Point Definition (DCM, 113780) | TEXT | 1 | ALWAYS | 15cm below BeamIsocenter |
|---|----------|--|------|---|--------|--------------------------|

### 8.3.2.8. TID 1002 Observer Context

Table 150: Observer Context

| NL | Relation with Parent | Concept Name                                       | VT      | VM | Presence of Value | Value                |
|----|----------------------|--|---------|----|-------------------|----------------------|
| >  | HAS OBS CONTEXT      | Observer Type (DCM, 121005)                        | CODE    | 1  | ALWAYS            | Device (DCM, 121007) |
|    | HAS OBS CONTEXT      | DTID 1004 "Device Observer Identifying Attributes" | INCLUDE | 1  | ALWAYS            |                      |

### 8.3.2.9. TID 1004 Device Observer Identifying Attributes

Table 151: Device Observer Identifying Attributes

| NL | Relation with Parent | Concept Name                                | VT     | VM | Presence of Value | Value                |
|----|----------------------|---|--------|----|-------------------|----------------------|
| >  | HAS OBS CONTEXT      | Device Observer Manufacturer (DCM, 121014)  | TEXT   | 1  | ALWAYS            |                      |
| >  | HAS OBS CONTEXT      | Observer Type (DCM, 121005)                 | CODE   | 1  | ALWAYS            | Device (DCM, 121007) |
| >  | HAS OBS CONTEXT      | Device Observer Manufacturer (DCM, 121014)  | TEXT   | 1  | ALWAYS            |                      |
| >  | HAS OBS CONTEXT      | Device Observer Model Name (DCM, 121015)    | TEXT   | 1  | ALWAYS            |                      |
| >  | HAS OBS CONTEXT      | Device Observer Name (DCM, 121013)          | TEXT   | 1  | ALWAYS            |                      |
| >  | HAS OBS CONTEXT      | Device Observer Serial Number (DCM, 121016) | TEXT   | 1  | ALWAYS            |                      |
| >  | HAS OBS CONTEXT      | Device Observer UID (DCM, 121012)           | UIDREF | 1  | ALWAYS            |                      |

### 8.3.3. Private Code Definitions

Not applicable.

## 8.4. Grayscale Image Consistency

The monitors and printers attached to the product can be calibrated by using the Service Application.

## 8.5. Standard/Extended/Specialized/Private SOP Classes

Not applicable.