


DICOM Conformance Statement

Application Annex:

Multi-Modality Applications on Philips IntelliSpace Portal V4.0



Issued by:

Philips Medical Systems Nederland BV, a Philips Healthcare company,

P.O. Box 10.000
5680 DA Best
The Netherlands

Email: dicom@philips.com

Internet: <http://www.healthcare.philips.com/connectivity>

Document Number: 7868707.01

Date: 23-Dec-2011

Table of Contents

| | |
|--|----------|
| 1. INTRODUCTION | 4 |
| 1.1. REVISION HISTORY | 4 |
| 1.2. TERMINOLOGY | 4 |
| 2. DATA SPECIFICATIONS | 5 |
| 2.1. SUPPORTED IOD'S | 5 |
| 2.1.1. Acceptance Criteria | 5 |
| 2.1.2. Contents of Created IOD's | 6 |
| 2.1.2.1. Encapsulated PDF Storage | 7 |
| 2.1.2.2. Grayscale Softcopy Presentation State Storage SOP Class | 9 |
| 2.1.2.3. Positron Emission Tomography Image Storage SOP Class | 15 |
| 2.1.2.4. CT Image Storage SOP Class | 24 |
| 2.1.2.5. Nuclear Medicine Image Storage SOP Class | 32 |
| 2.1.2.6. MR Image Storage SOP Class | 42 |
| 2.1.2.7. RT Structure Set Storage SOP Class | 53 |
| 2.1.2.8. Secondary Capture Image Storage SOP Class | 57 |

1. Introduction

This DICOM Conformance Statement annex is applicable to the Multi-Modality Applications on Philips IntelliSpace Portal V4.0 hosting platform, later referred to as Multi-Modality Applications. In general the Multi-Modality Applications are used for viewing and analyzing the images from multiple modalities.

The following analysis packages are offered by Multi-Modality Applications (licensing might be applicable):

- Multi-Modality Viewer
- Multi-Modality Tumor Tracking

1.1. Revision History

The revision history below provides dates and differences among individual document versions.

Table 1: Revision History

| Document Version | Date of Issue | Status | Description |
|------------------|---------------|---------------|---|
| 00 | 01-Jul-2011 | Final version | Updated after review. |
| 01 | 23-Dec-2011 | Final version | Merged MMTT Annex into this Multi-Modality Annex document |

1.2. Terminology

| | |
|-------|--|
| DICOM | Digital Imaging and Communications in Medicine |
| IOD | Information Object Definition |
| UID | Unique Identifier |
| VR | Value Representation |

2. Data Specifications

2.1. Supported IOD's

This section specifies each IOD accepted and / or created by Multi-Modality Applications.

ACCEPTED The applicable IOD is accepted for storage in the repository of the hosting platform and supported for import in Multi-Modality Applications for viewing and analysis.

CREATED The Multi-Modality Applications supports generation of derived data by using the applicable IOD and are able to store this data in the repository of the hosting platform.

Table 2: Supported IOD's

| IOD | | Support | |
|--|-------------------------------|----------|---------|
| Name | UID | ACCEPTED | CREATED |
| Encapsulated PDF Storage | 1.2.840.10008.5.1.4.1.1.104.1 | No | Yes |
| Grayscale Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.1 | Yes | Yes |
| Positron Emission Tomography Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.128 | Yes | Yes |
| CT Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.2 | Yes | Yes |
| Nuclear Medicine Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.20 | Yes | Yes |
| MR Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.4 | Yes | Yes |
| RT Structure Set Storage SOP Class | 1.2.840.10008.5.1.4.1.1.481.3 | Yes | Yes |
| Secondary Capture Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7 | Yes | Yes |
| Multi-frame Single Bit Secondary Capture Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7.1 | Yes | No |
| Multi-frame Grayscale Byte SC Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7.2 | Yes | No |
| Multi-frame Grayscale Word SC Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7.3 | Yes | No |
| Multi-frame True Color Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.4 | Yes | No |

2.1.1. Acceptance Criteria

This section specifies the acceptance criteria applied by Multi-Modality Applications to which a dataset should adhere before it can be imported into the application. This can be criteria on the highest level (e.g. data from a certain manufacturer or system model) or certain DICOM attributes mandatory to be present into the dataset holding a specific value. In case one or more Philips private attributes are required, then a list of supported Philips system models will be mentioned.

Table 3: Accepted system models

| Manufacturer | Modality | System Model Name(s) |
|----------------|----------------|----------------------|
| Not applicable | Not applicable | Not applicable |

Table 4: Accepted transfer syntaxes per IOD

| IOD | | Transfer Syntax | |
|---|------------------------------|---------------------------|------------------------|
| Name | UID | Name | UID |
| Grayscale Softcopy Presentation State Storage SOP Class | 1.2.840.10008.5.1.4.1.1.11.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| | | Implicit VR Little Endian | 1.2.840.10008.1.2 |
| Positron Emission Tomography Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.128 | CT-private-ELE | 1.3.46.670589.33.1.4.1 |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |

| | | | |
|--|-------------------------------|---|------------------------|
| | | Implicit VR Little Endian | 1.2.840.10008.1.2 |
| | | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |
| CT Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.2 | CT-private-ELE | 1.3.46.670589.33.1.4.1 |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| | | Implicit VR Little Endian | 1.2.840.10008.1.2 |
| | | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |
| Nuclear Medicine Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.20 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| | | Implicit VR Little Endian | 1.2.840.10008.1.2 |
| MR Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.4 | CT-private-ELE | 1.3.46.670589.33.1.4.1 |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| | | Implicit VR Little Endian | 1.2.840.10008.1.2 |
| | | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |
| RT Structure Set Storage SOP Class | 1.2.840.10008.5.1.4.1.1.481.3 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| | | Implicit VR Little Endian | 1.2.840.10008.1.2 |
| Secondary Capture Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7 | CT-private-ELE | 1.3.46.670589.33.1.4.1 |
| | | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| | | Implicit VR Little Endian | 1.2.840.10008.1.2 |
| | | JPEG Lossless, Non-Hierarchical, FOP (Process 14) | 1.2.840.10008.1.2.4.70 |
| Multi-frame Single Bit Secondary Capture Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7.1 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| | | Implicit VR Little Endian | 1.2.840.10008.1.2 |
| Multi-frame Grayscale Byte SC Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7.2 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| | | Implicit VR Little Endian | 1.2.840.10008.1.2 |
| Multi-frame Grayscale Word SC Image Storage SOP Class | 1.2.840.10008.5.1.4.1.1.7.3 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| | | Implicit VR Little Endian | 1.2.840.10008.1.2 |
| Multi-frame True Color Secondary Capture Image Storage | 1.2.840.10008.5.1.4.1.1.7.4 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 |
| | | Implicit VR Little Endian | 1.2.840.10008.1.2 |

2.1.2. Contents of Created IOD's

This section specifies in detail the attribute contents of created data objects. Attributes are grouped together by its corresponding module as specified by DICOM standard. Philips private attributes are excluded for specification.

Abbreviations used in the Module table for the column "Presence of Value" are:

| | |
|--------|--|
| ALWAYS | The attribute is always present with a value |
| EMPTY | The attribute is always present without any value (attribute sent zero length) |
| VNAP | The attribute is always present and its Value is Not Always Present (attribute sent zero length if no value is present) |
| ANAP | The attribute is present under specified condition – if present then it will always have a value |
| ANAPCV | The attribute is present under specified condition – if present then its Value is Not Always Present (attribute sent zero length if condition applies and no value is present) |
| ANAPEV | The attribute is present under specified condition – if present then it will not have any value |

The abbreviations used in the Module table for the column "Source" are:

| | |
|--------|--|
| AUTO | The attribute value is generated automatically |
| CONFIG | The attribute value source is a configurable parameter |
| COPY | The attribute value source is another SOP instance |

| | |
|----------|---|
| FIXED | The attribute value is hard-coded in the application |
| IMPLICIT | The attribute value source is a user-implicit setting |
| MPPS | The attribute value is the same as that use for Modality Performed Procedure Step |
| COPY | The attribute value source is a Modality Worklist |
| USER | The attribute value source is explicit user input |

2.1.2.1. Encapsulated PDF Storage

Table 5: IOD of Created Encapsulated PDF Storage Instances

| Information Entity | Module | Presence Of Module |
|-----------------------|-------------------------------------|--------------------|
| Patient | Patient Module | ALWAYS |
| Study | General Study Module | ALWAYS |
| Series | Encapsulated Document Series Module | ALWAYS |
| Equipment | General Equipment Module | ALWAYS |
| Equipment | SC Equipment Module | ALWAYS |
| Encapsulated Document | Encapsulated Document Module | ALWAYS |
| Encapsulated Document | SOP Common Module | ALWAYS |

Table 6: Patient Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|---------------|---------|
| Referenced Patient Sequence | 0008,1120 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Patient's Name | 0010,0010 | PN | | VNAP | COPY | - |
| Patient ID | 0010,0020 | LO | | VNAP | COPY | - |
| Patient's Birth Date | 0010,0030 | DA | | VNAP | COPY, USER | - |
| Patient's Birth Time | 0010,0032 | TM | | ANAPCV | COPY, USER | - |
| Patient's Sex | 0010,0040 | CS | | VNAP | COPY, USER | - |
| Other Patient IDs | 0010,1000 | LO | | ANAPCV | COPY, USER | - |
| Other Patient Names | 0010,1001 | PN | | ANAPCV | COPY, USER | - |
| Ethnic Group | 0010,2160 | SH | | ANAPCV | COPY | - |
| Patient Comments | 0010,4000 | LT | | ANAPCV | COPY | - |
| Issuer of Patient ID | 0010,0021 | LO | | ANAPCV | COPY | - |

Table 7: General Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---|-----------|----|-------|-------------------|---------------|---------|
| Study Date | 0008,0020 | DA | | VNAP | COPY | - |
| Study Time | 0008,0030 | TM | | VNAP | COPY | - |
| Accession Number | 0008,0050 | SH | | VNAP | COPY, USER | - |
| Referring Physician's Name | 0008,0090 | PN | | VNAP | COPY, USER | - |
| Referring Physician Identification Sequence | 0008,0096 | SQ | | ANAPCV | COPY | - |
| Study Description | 0008,1030 | LO | | ANAPCV | COPY, USER | - |

| | | | | | | |
|--|-----------|----|--|--------|------------|---|
| Procedure Code Sequence | 0008,1032 | SQ | | ANAPCV | COPY | - |
| Physician(s) of Record | 0008,1048 | PN | | ANAPCV | COPY | - |
| Physician(s) of Record Identification Sequence | 0008,1049 | SQ | | ANAP | COPY | - |
| Name of Physician(s) Reading Study | 0008,1060 | PN | | ANAPCV | COPY | - |
| Physician(s) Reading Study Identification Sequence | 0008,1062 | SQ | | ANAPCV | COPY | - |
| Referenced Study Sequence | 0008,1110 | SQ | | ANAPCV | COPY | - |
| Study Instance UID | 0020,000D | UI | | ALWAYS | AUTO, COPY | - |
| Study ID | 0020,0010 | SH | | VNAP | AUTO, COPY | - |

Table 8: Encapsulated Document Series Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---------------------|-----------|----|---------------------|-------------------|--------|---|
| Modality | 0008,0060 | CS | | ALWAYS | AUTO | - |
| Series Description | 0008,103E | LO | QC date-time-status | ANAPCV | AUTO | Describes date and time of QC procedure and pass/fail status. |
| Series Instance UID | 0020,000E | UI | | ALWAYS | AUTO | - |
| Series Number | 0020,0011 | IS | | ALWAYS | AUTO | - |

Table 9: General Equipment Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------|-----------|-------|-------|-------------------|------------|---------|
| Manufacturer | 0008,0070 | LO | | ALWAYS | COPY | - |
| Institution Name | 0008,0080 | LO | | ANAPCV | COPY, USER | - |
| Institution Address | 0008,0081 | ST | | ANAPCV | COPY, USER | - |
| Station Name | 0008,1010 | SH | | ALWAYS | COPY | - |
| Institutional Department Name | 0008,1040 | LO | | ANAPCV | COPY | - |
| Manufacturer's Model Name | 0008,1090 | LO | | ALWAYS | COPY | - |
| Device Serial Number | 0018,1000 | LO | | ANAPCV | COPY | - |
| Software Version(s) | 0018,1020 | LO | 4.0 | ALWAYS | COPY | 4.0 |
| Spatial Resolution | 0018,1050 | DS | | ANAPCV | COPY | - |
| Date of Last Calibration | 0018,1200 | DA | | ANAPCV | COPY | - |
| Time of Last Calibration | 0018,1201 | TM | | ANAPCV | COPY | - |
| Pixel Padding Value | 0028,0120 | US/SS | | ANAPCV | COPY | - |

Table 10: SC Equipment Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------|-----------|----|-------|-------------------|--------|---------|
| Modality | 0008,0060 | CS | | ALWAYS | AUTO | - |
| Conversion Type | 0008,0064 | CS | WSD | ALWAYS | AUTO | - |

Table 11: Encapsulated Document Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|----------------------|-----------|----|-------|-------------------|--------|---------|
| Content Date | 0008,0023 | DA | | ANAP | AUTO | - |
| Acquisition Datetime | 0008,002A | DT | | ANAP | AUTO | - |
| Content Time | 0008,0033 | TM | | ANAP | AUTO | - |
| Instance Number | 0020,0013 | IS | 1 | ALWAYS | AUTO | - |

| | | | | | | |
|------------------------------------|-----------|----|-------------------|--------|------|---|
| Burned In Annotation | 0028,0301 | CS | YES | ALWAYS | AUTO | - |
| Concept Name Code Sequence | 0040,A043 | SQ | | EMPTY | AUTO | - |
| Document Title | 0042,0010 | ST | | EMPTY | AUTO | - |
| Encapsulated Document | 0042,0011 | OB | | ALWAYS | AUTO | - |
| MIME Type of Encapsulated Document | 0042,0012 | LO | application - PDF | ALWAYS | AUTO | - |

Table 12: SOP Common Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------------|-----------|----|----------------------|-------------------|------------|----------------------|
| Specific Character Set | 0008,0005 | CS | | ANAP | AUTO, COPY | - |
| Instance Creation Date | 0008,0012 | DA | | ALWAYS | AUTO | - |
| Instance Creation Time | 0008,0013 | TM | | ALWAYS | AUTO | - |
| Instance Creator UID | 0008,0014 | UI | | ANAPCV | AUTO | - |
| SOP Class UID | 0008,0016 | UI | | ALWAYS | COPY | - |
| SOP Instance UID | 0008,0018 | UI | | ALWAYS | AUTO | - |
| Timezone Offset From UTC | 0008,0201 | SH | | ANAPCV | COPY | - |
| Contributing Equipment Sequence | 0018,A001 | SQ | | ANAP | AUTO | - |
| >Manufacturer | 0008,0070 | LO | Philips | ANAP | FIXED | PHILIPS. |
| >Institution Address | 0008,0081 | ST | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Station Name | 0008,1010 | SH | IntelliSpace Portal | ANAP | CONFIG | IntelliSpace Portal. |
| >Institutional Department Name | 0008,1040 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Software Version(s) | 0018,1020 | LO | 4.0 | ANAP | FIXED | 4.0 |
| >Device Serial Number | 0018,1000 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Institution Name | 0008,0080 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Manufacturer's Model Name | 0008,1090 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Purpose of Reference Code Sequence | 0040,A170 | SQ | | ANAP | FIXED | - |
| >>Code Value | 0008,0100 | SH | 109102 | ANAP | FIXED | 109102 |
| >>Coding Scheme Designator | 0008,0102 | SH | DCM | ANAP | FIXED | DCM |
| >>Code Meaning | 0008,0104 | LO | Processing Equipment | ANAP | FIXED | Processing Equipment |

2.1.2.2. Grayscale Softcopy Presentation State Storage SOP Class

Table 13: IOD of Created Grayscale Softcopy Presentation State Storage SOP Class Instances

| Information Entity | Module | Presence Of Module |
|--------------------|--|--------------------|
| Patient | Patient Module | ALWAYS |
| Study | General Study Module | ALWAYS |
| Study | Patient Study Module | CONDITIONAL |
| Series | General Series Module | ALWAYS |
| Series | Presentation Series Module | ALWAYS |
| Equipment | General Equipment Module | ALWAYS |
| Presentation State | Presentation State Identification Module | ALWAYS |
| Presentation State | Presentation State Relationship Module | ALWAYS |
| Presentation State | Display Shutter Module | CONDITIONAL |
| Presentation State | Overlay Plane Module | CONDITIONAL |
| Presentation State | Overlay Activation Module | CONDITIONAL |
| Presentation State | Displayed Area Module | CONDITIONAL |
| Presentation State | Graphic Annotation Module | CONDITIONAL |

| | | |
|--------------------|----------------------------------|-------------|
| Presentation State | Spatial Transformation Module | CONDITIONAL |
| Presentation State | Graphic Layer Module | CONDITIONAL |
| Presentation State | Modality LUT Module | CONDITIONAL |
| Presentation State | Softcopy VOI LUT Module | CONDITIONAL |
| Presentation State | Softcopy Presentation LUT Module | ALWAYS |
| Presentation State | SOP Common Module | ALWAYS |

Table 14: Patient Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|--------|---------|
| Referenced Patient Sequence | 0008,1120 | SQ | | ANAP | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Patient's Name | 0010,0010 | PN | | VNAP | COPY | - |
| Patient ID | 0010,0020 | LO | | VNAP | COPY | - |
| Patient's Birth Date | 0010,0030 | DA | | VNAP | COPY | - |
| Patient's Birth Time | 0010,0032 | TM | | ANAPCV | COPY | - |
| Patient's Sex | 0010,0040 | CS | | VNAP | COPY | - |
| Other Patient IDs | 0010,1000 | LO | | ANAPCV | COPY | - |
| Other Patient Names | 0010,1001 | PN | | ANAPCV | COPY | - |
| Ethnic Group | 0010,2160 | SH | | ANAPCV | COPY | - |
| Patient Comments | 0010,4000 | LT | | ANAPCV | COPY | - |
| Issuer of Patient ID | 0010,0021 | LO | | ANAPCV | COPY | - |

Table 15: General Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|--------|---------|
| Study Date | 0008,0020 | DA | | VNAP | COPY | - |
| Study Time | 0008,0030 | TM | | VNAP | COPY | - |
| Accession Number | 0008,0050 | SH | | VNAP | COPY | - |
| Referring Physician's Name | 0008,0090 | PN | | VNAP | COPY | - |
| Referring Physician Identification Sequence | 0008,0096 | SQ | | ANAP | COPY | - |
| Study Description | 0008,1030 | LO | | ANAPCV | COPY | - |
| Procedure Code Sequence | 0008,1032 | SQ | | ANAP | COPY | - |
| Physician(s) of Record | 0008,1048 | PN | | ANAPCV | COPY | - |
| Physician(s) of Record Identification Sequence | 0008,1049 | SQ | | ANAP | COPY | - |
| Name of Physician(s) Reading Study | 0008,1060 | PN | | ANAPCV | COPY | - |
| Physician(s) Reading Study Identification Sequence | 0008,1062 | SQ | | ANAP | COPY | - |
| Referenced Study Sequence | 0008,1110 | SQ | | ANAPCV | COPY | - |
| Study Instance UID | 0020,000D | UI | | ALWAYS | COPY | - |
| Study ID | 0020,0010 | SH | | VNAP | COPY | - |

Table 16: Patient Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------------|-----------|----|-------|-------------------|--------|---------|
| Admitting Diagnoses Description | 0008,1080 | LO | | ANAPCV | COPY | - |
| Admitting Diagnoses Code Sequence | 0008,1084 | SQ | | ANAPCV | COPY | - |

| | | | | | | |
|----------------------------|-----------|----|--|--------|------|---|
| Patient's Age | 0010,1010 | AS | | ANAPCV | COPY | - |
| Patient's Size | 0010,1020 | DS | | ANAPCV | COPY | - |
| Patient's Weight | 0010,1030 | DS | | ANAPCV | COPY | - |
| Occupation | 0010,2180 | SH | | ANAPCV | COPY | - |
| Additional Patient History | 0010,21B0 | LT | | ANAPCV | COPY | - |

Table 17: General Series Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|-------|-------|-------------------|------------|---------|
| Series Date | 0008,0021 | DA | | ANAPCV | AUTO | - |
| Series Time | 0008,0031 | TM | | ANAPCV | AUTO | - |
| Modality | 0008,0060 | CS | | ALWAYS | COPY | - |
| Series Description | 0008,103E | LO | | ANAPCV | COPY, USER | - |
| Performing Physician's Name | 0008,1050 | PN | | ANAPCV | COPY | - |
| Operators' Name | 0008,1070 | PN | | ANAPCV | COPY | - |
| Operator Identification Sequence | 0008,1072 | SQ | | ANAPCV | COPY | - |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Body Part Examined | 0018,0015 | CS | | ANAPCV | COPY | - |
| Protocol Name | 0018,1030 | LO | | ANAPCV | COPY | - |
| Patient Position | 0018,5100 | CS | | VNAP | COPY | - |
| Series Instance UID | 0020,000E | UI | | ALWAYS | AUTO | - |
| Series Number | 0020,0011 | IS | | VNAP | AUTO | - |
| Laterality | 0020,0060 | CS | | VNAP | COPY | - |
| Smallest Pixel Value in Series | 0028,0108 | US/SS | | ANAPCV | AUTO | - |
| Largest Pixel Value in Series | 0028,0109 | US/SS | | ANAPCV | AUTO | - |
| Request Attributes Sequence | 0040,0275 | SQ | | ANAPCV | COPY | - |
| Performed Procedure Step Start Date | 0040,0244 | DA | | ANAPCV | COPY | - |
| Performed Procedure Step Start Time | 0040,0245 | TM | | ANAPCV | COPY | - |
| Performed Procedure Step Description | 0040,0254 | LO | | ANAPCV | COPY | - |
| Performed Protocol Code Sequence | 0040,0260 | SQ | | ANAPCV | COPY | - |
| Comments on the Performed Procedure Step | 0040,0280 | ST | | ANAPCV | COPY | - |

Table 18: Presentation Series Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|----------------|-----------|----|-------|-------------------|--------|---------|
| Modality | 0008,0060 | CS | | ALWAYS | AUTO | - |

Table 19: General Equipment Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|-------|-------------------|--------|---|
| Manufacturer | 0008,0070 | LO | | ALWAYS | AUTO | Applied value: Philips medical systems. |
| Institution Name | 0008,0080 | LO | | ALWAYS | CONFIG | - |
| Station Name | 0008,1010 | SH | | ALWAYS | AUTO | Same as the host Name. |
| Institutional Department Name | 0008,1040 | LO | | ALWAYS | CONFIG | - |

| | | | | | | |
|---------------------------|-----------|----|-----|--------|------|-----|
| Manufacturer's Model Name | 0008,1090 | LO | | ALWAYS | AUTO | - |
| Device Serial Number | 0018,1000 | LO | | ALWAYS | AUTO | - |
| Software Version(s) | 0018,1020 | LO | 4.0 | ALWAYS | AUTO | 4.0 |

Table 20: Presentation State Identification Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|----------------------------|-----------|----|--------------------------------|-------------------|--------|---|
| Presentation Creation Date | 0070,0082 | DA | | ALWAYS | AUTO | - |
| Presentation Creation Time | 0070,0083 | TM | | ALWAYS | AUTO | - |
| Instance Number | 0020,0013 | IS | | ALWAYS | AUTO | - |
| Content Label | 0070,0080 | CS | AS LAST SEEN, NEW AT IMPORT | ALWAYS | AUTO | Applied values: AS LAST SEEN, NEW AT IMPORT |
| Content Description | 0070,0081 | LO | | VNAP | AUTO | - |
| Content Creator's Name | 0070,0084 | PN | | VNAP | AUTO | - |

Table 21: Presentation State Relationship Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|-------|-------------------|--------|---------|
| Referenced Series Sequence | 0008,1115 | SQ | | ALWAYS | AUTO | - |
| >Referenced Image Sequence | 0008,1140 | SQ | | ANAP | AUTO | - |
| >>Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | AUTO | - |
| >>Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | AUTO | - |
| >Series Instance UID | 0020,000E | UI | | ALWAYS | AUTO | - |

Table 22: Display Shutter Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|----------------------------|-----------|----|-------|-------------------|--------|------------------|
| Shutter Presentation Value | 0018,1622 | US | | ANAP | AUTO | Applied value: 0 |

Table 23: Overlay Plane Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------|-----------|-------|-------|-------------------|--------|---------|
| Overlay Rows | 6000,0010 | US | | ALWAYS | AUTO | - |
| Overlay Columns | 6000,0011 | US | | ALWAYS | AUTO | - |
| Overlay Description | 6000,0022 | LO | | ANAPEV | AUTO | - |
| Overlay Type | 6000,0040 | CS | | ALWAYS | AUTO | - |
| Overlay Subtype | 6000,0045 | LO | | ANAPEV | AUTO | - |
| Overlay Origin | 6000,0050 | SS | | ALWAYS | AUTO | - |
| Overlay Bits Allocated | 6000,0100 | US | | ALWAYS | AUTO | - |
| Overlay Bit Position | 6000,0102 | US | | ALWAYS | AUTO | - |
| ROI Area | 6000,1301 | IS | | ANAPEV | AUTO | - |
| ROI Mean | 6000,1302 | DS | | ANAPEV | AUTO | - |
| ROI Standard Deviation | 6000,1303 | DS | | ANAPEV | AUTO | - |
| Overlay Label | 6000,1500 | LO | | ANAPCV | AUTO | - |
| Overlay Data | 6000,3000 | OW/OB | | ALWAYS | AUTO | - |

Table 24: Overlay Activation Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--------------------------|-----------|----|-------|-------------------|--------|------------------|
| Overlay Activation Layer | 6000,1001 | CS | | ANAP | AUTO | Applied value: 1 |

Table 25: Displayed Area Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-----------------------|-------------------|----------|---------------------------------------|
| Displayed Area Selection Sequence | 0070,005A | SQ | | ALWAYS | IMPLICIT | - |
| >Referenced Image Sequence | 0008,1140 | SQ | | ANAP | IMPLICIT | - |
| >>Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | IMPLICIT | - |
| >>Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | IMPLICIT | - |
| >Displayed Area Top Left Hand Corner | 0070,0052 | SL | | ALWAYS | IMPLICIT | - |
| >Displayed Area Bottom Right Hand Corner | 0070,0053 | SL | | ALWAYS | IMPLICIT | - |
| >Presentation Size Mode | 0070,0100 | CS | MAGNIFY, SCALE TO FIT | ALWAYS | IMPLICIT | Applied values: MAGNIFY, SCALE TO FIT |
| >Presentation Pixel Spacing | 0070,0101 | DS | | ANAP | IMPLICIT | Applied values: (0.0, 0.0) |
| >Presentation Pixel Aspect Ratio | 0070,0102 | IS | | ANAP | IMPLICIT | - |
| >Presentation Pixel Magnification Ratio | 0070,0103 | FL | | ANAP | IMPLICIT | Applied value: 1.0 |
| >Zoom Mode | 2001,103F | CS | | VNAP | IMPLICIT | - |

Table 26: Graphic Annotation Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---------------------------------|-----------|----|-------|-------------------|----------|---------|
| Graphic Annotation Sequence | 0070,0001 | SQ | | ALWAYS | IMPLICIT | - |
| >Referenced Image Sequence | 0008,1140 | SQ | | ANAP | IMPLICIT | - |
| >>Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | IMPLICIT | - |
| >>Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | IMPLICIT | - |
| >Graphic Layer | 0070,0002 | CS | | ALWAYS | IMPLICIT | - |
| >Text Object Sequence | 0070,0008 | SQ | | ANAP | IMPLICIT | - |
| >>Anchor Point Annotation Units | 0070,0004 | CS | | ANAP | IMPLICIT | - |
| >>Unformatted Text Value | 0070,0006 | ST | | ALWAYS | IMPLICIT | - |
| >>Anchor Point | 0070,0014 | FL | | ALWAYS | IMPLICIT | - |
| >>Anchor Point Visibility | 0070,0015 | CS | | ALWAYS | IMPLICIT | - |
| >Graphic Object Sequence | 0070,0009 | SQ | | ANAP | AUTO | - |
| >>Graphic Annotation Units | 0070,0005 | CS | | ALWAYS | IMPLICIT | - |
| >>Graphic Dimensions | 0070,0020 | US | | ALWAYS | IMPLICIT | - |
| >>Number of Graphic Points | 0070,0021 | US | | ALWAYS | IMPLICIT | - |
| >>Graphic Data | 0070,0022 | FL | | ALWAYS | IMPLICIT | - |
| >>Graphic Type | 0070,0023 | CS | | ALWAYS | IMPLICIT | - |
| >>Graphic Filled | 0070,0024 | CS | | ANAP | IMPLICIT | - |

Table 27: Spatial Transformation Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------|-----------|----|-------|-------------------|----------|---------|
| Image Horizontal Flip | 0070,0041 | CS | | ALWAYS | IMPLICIT | - |
| Image Rotation | 0070,0042 | US | | ALWAYS | IMPLICIT | - |

Table 28: Graphic Layer Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------|-----------|----|-------|-------------------|----------|---------|
| Graphic Layer Sequence | 0070,0060 | SQ | | ALWAYS | IMPLICIT | - |
| >Graphic Layer | 0070,0002 | CS | | ALWAYS | IMPLICIT | - |
| >Graphic Layer Order | 0070,0062 | IS | | ALWAYS | IMPLICIT | - |

Table 29: Modality LUT Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------|-----------|----|--|-------------------|--------|---|
| Rescale Intercept | 0028,1052 | DS | | ALWAYS | COPY | - |
| Rescale Slope | 0028,1053 | DS | | ALWAYS | COPY | - |
| Rescale Type | 0028,1054 | LO | cm/sec, millirads, milliseconds, mm ² /sec, normalized, seconds, US | ALWAYS | COPY | Applied Value(s): cm/sec, milliradians, milliseconds, mm ² /sec, normalized, seconds, US |

Table 30: Softcopy VOI LUT Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|-------|-------------------|--------|---------|
| Softcopy VOI LUT Sequence | 0028,3110 | SQ | | ALWAYS | AUTO | - |
| >Referenced Image Sequence | 0008,1140 | SQ | | ANAP | AUTO | - |
| >>Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | AUTO | - |
| >>Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | AUTO | - |
| >Window Center | 0028,1050 | DS | | ALWAYS | AUTO | - |
| >Window Width | 0028,1051 | DS | | ALWAYS | AUTO | - |

Table 31: Softcopy Presentation LUT Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---------------------------|-----------|-------|-------|-------------------|--------|--|
| Presentation LUT Sequence | 2050,0010 | SQ | | ANAP | AUTO | Present if Presentation LUT Shape not present. |
| >LUT Descriptor | 0028,3002 | US/SS | | ALWAYS | AUTO | - |
| >LUT Data | 0028,3006 | US/OW | | ALWAYS | AUTO | - |
| Presentation LUT Shape | 2050,0020 | CS | | ALWAYS | AUTO | - |

Table 32: SOP Common Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------|-----------|----|-------|-------------------|--------|--|
| Specific Character Set | 0008,0005 | CS | | ANAP | AUTO | - |
| Instance Creation Date | 0008,0012 | DA | | ALWAYS | AUTO | - |
| Instance Creator UID | 0008,0014 | UI | | ALWAYS | AUTO | - |
| SOP Class UID | 0008,0016 | UI | | ALWAYS | AUTO | Applied Value(s): 1.2.840.10008.5.1.4.1.1.11.1 |
| SOP Instance UID | 0008,0018 | UI | | ALWAYS | AUTO | - |

2.1.2.3. Positron Emission Tomography Image Storage SOP Class

Table 33: IOD of Created Positron Emission Tomography Image Storage SOP Class Instances

| Information Entity | Module | Presence Of Module |
|--------------------|---------------------------------------|--------------------|
| Patient | Patient Module | ALWAYS |
| Study | General Study Module | ALWAYS |
| Study | Patient Study Module | CONDITIONAL |
| Series | General Series Module | ALWAYS |
| Series | PET Series Module | ALWAYS |
| Frame of Reference | Frame of Reference Module | ALWAYS |
| Equipment | General Equipment Module | ALWAYS |
| Image | General Image Module | ALWAYS |
| Image | Image Plane Module | ALWAYS |
| Image | Image Pixel Module | ALWAYS |
| Image | NM/PET Patient Orientation Module | ALWAYS |
| Image | PET Isotope Module | ALWAYS |
| Image | PET Multi-Gated Acquisition Module | CONDITIONAL |
| Image | PET Image Module | ALWAYS |
| Image | Overlay Plane Module | CONDITIONAL |
| Image | VOI LUT Module | CONDITIONAL |
| Image | SOP Common Module | ALWAYS |
| | Extended Dicom and Private attributes | CONDITIONAL |

Table 34: Patient Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|---------------|---------|
| Referenced Patient Sequence | 0008,1120 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Patient's Name | 0010,0010 | PN | | VNAP | COPY | - |
| Patient ID | 0010,0020 | LO | | VNAP | COPY | - |
| Patient's Birth Date | 0010,0030 | DA | | VNAP | COPY, USER | - |
| Patient's Birth Time | 0010,0032 | TM | | ANAPCV | COPY, USER | - |
| Patient's Sex | 0010,0040 | CS | | VNAP | COPY, USER | - |
| Other Patient IDs | 0010,1000 | LO | | ANAPCV | COPY, USER | - |
| Other Patient Names | 0010,1001 | PN | | ANAPCV | COPY, USER | - |
| Ethnic Group | 0010,2160 | SH | | ANAPCV | COPY | - |
| Patient Comments | 0010,4000 | LT | | ANAPCV | COPY | - |
| Issuer of Patient ID | 0010,0021 | LO | | ANAPCV | COPY, USER | - |

Table 35: General Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------|-----------|----|-------|-------------------|---------------|---------|
| Study Date | 0008,0020 | DA | | VNAP | COPY | - |
| Study Time | 0008,0030 | TM | | VNAP | COPY | - |
| Accession Number | 0008,0050 | SH | | VNAP | COPY, USER | - |

| | | | | | | |
|--|-----------|----|--|--------|------------|---|
| Referring Physician's Name | 0008,0090 | PN | | VNAP | COPY, USER | - |
| Referring Physician Identification Sequence | 0008,0096 | SQ | | ANAPCV | COPY | - |
| Study Description | 0008,1030 | LO | | ANAPCV | COPY, USER | - |
| Procedure Code Sequence | 0008,1032 | SQ | | ANAPCV | COPY | - |
| Physician(s) of Record | 0008,1048 | PN | | ANAPCV | COPY | - |
| Physician(s) of Record Identification Sequence | 0008,1049 | SQ | | ANAPCV | COPY | - |
| Name of Physician(s) Reading Study | 0008,1060 | PN | | ANAPCV | COPY | - |
| Physician(s) Reading Study Identification Sequence | 0008,1062 | SQ | | ANAPCV | COPY | - |
| Referenced Study Sequence | 0008,1110 | SQ | | ANAPCV | COPY | - |
| Study Instance UID | 0020,000D | UI | | ALWAYS | AUTO, COPY | - |
| Study ID | 0020,0010 | SH | | VNAP | AUTO, COPY | - |

Table 36: Patient Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------------|-----------|----|-------|-------------------|------------|---------|
| Admitting Diagnoses Description | 0008,1080 | LO | | ANAPCV | COPY | - |
| Admitting Diagnoses Code Sequence | 0008,1084 | SQ | | ANAPCV | COPY | - |
| Patient's Age | 0010,1010 | AS | | ANAPCV | COPY, USER | - |
| Patient's Size | 0010,1020 | DS | | ANAPCV | COPY, USER | - |
| Patient's Weight | 0010,1030 | DS | | ANAPCV | COPY, USER | - |
| Occupation | 0010,2180 | SH | | ANAPCV | COPY | - |
| Additional Patient History | 0010,21B0 | LT | | ANAPCV | COPY | - |

Table 37: General Series Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|------------|---------|
| Modality | 0008,0060 | CS | | ALWAYS | COPY | - |
| Series Description | 0008,103E | LO | | ANAPCV | COPY, USER | - |
| Performing Physician's Name | 0008,1050 | PN | | ANAPCV | COPY | - |
| Operators' Name | 0008,1070 | PN | | ANAPCV | COPY, USER | - |
| Operator Identification Sequence | 0008,1072 | SQ | | ANAPCV | COPY | - |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Body Part Examined | 0018,0015 | CS | | ANAPCV | COPY | - |
| Protocol Name | 0018,1030 | LO | | ANAPCV | COPY | - |
| Patient Position | 0018,5100 | CS | | ANAPCV | COPY | - |
| Series Instance UID | 0020,000E | UI | | ALWAYS | AUTO | - |

| | | | | | | |
|--|-----------|-------|--|--------|------|---|
| Series Number | 0020,0011 | IS | | VNAP | AUTO | - |
| Laterality | 0020,0060 | CS | | ANAPCV | COPY | - |
| Smallest Pixel Value in Series | 0028,0108 | US/SS | | ANAPCV | AUTO | - |
| Largest Pixel Value in Series | 0028,0109 | US/SS | | ANAPCV | AUTO | - |
| Request Attributes Sequence | 0040,0275 | SQ | | ANAPCV | COPY | - |
| Performed Procedure Step Start Date | 0040,0244 | DA | | ANAPCV | COPY | - |
| Performed Procedure Step Start Time | 0040,0245 | TM | | ANAPCV | COPY | - |
| Performed Procedure Step Description | 0040,0254 | LO | | ANAPCV | COPY | - |
| Performed Protocol Code Sequence | 0040,0260 | SQ | | ANAPCV | COPY | - |
| Comments on the Performed Procedure Step | 0040,0280 | ST | | ANAPCV | COPY | - |

Table 38: PET Series Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|--------|---------|
| Series Date | 0008,0021 | DA | | ALWAYS | AUTO | - |
| Series Time | 0008,0031 | TM | | ALWAYS | AUTO | - |
| Series Type | 0054,1000 | CS | | ALWAYS | AUTO | - |
| Acquisition Start Condition | 0018,0073 | CS | | ANAPEV | COPY | - |
| Acquisition Start Condition Data | 0018,0074 | IS | | ANAPEV | COPY | - |
| Acquisition Termination Condition | 0018,0071 | CS | | ANAPEV | COPY | - |
| Acquisition Termination Condition Data | 0018,0075 | IS | | ANAPEV | COPY | - |
| Attenuation Correction Method | 0054,1101 | LO | | ANAPEV | COPY | - |
| Axial Mash | 0054,1201 | IS | | ANAPEV | COPY | - |
| Coincidence Window Width | 0054,1210 | DS | | ANAPEV | COPY | - |
| Collimator Type | 0018,1181 | CS | | VNAP | COPY | - |
| Corrected Image | 0028,0051 | CS | | VNAP | COPY | - |
| Counts Source | 0054,1002 | CS | | ALWAYS | COPY | - |
| Decay Correction | 0054,1102 | CS | | ALWAYS | COPY | - |
| Field of View Dimension(s) | 0018,1149 | IS | | ANAPEV | COPY | - |
| Field of View Shape | 0018,1147 | CS | | ANAPEV | COPY | - |
| Gantry/Detector Tilt | 0018,1120 | DS | | ANAPEV | COPY | - |
| Number of Slices | 0054,0081 | US | | ALWAYS | COPY | - |
| Number of Time Slices | 0054,0101 | US | | ANAP | COPY | - |

| | | | | | | |
|---------------------------------|-----------|----|--|--------|------|---|
| Randoms Correction Method | 0054,1100 | CS | | ANAPEV | COPY | - |
| Reconstruction Method | 0054,1103 | LO | | ANAPEV | COPY | - |
| Scatter Correction Method | 0054,1105 | LO | | ANAPEV | COPY | - |
| Transverse Mash | 0054,1202 | IS | | ANAPEV | COPY | - |
| Type of Detector Motion | 0054,0202 | CS | | ANAPEV | COPY | - |
| Units | 0054,1001 | CS | | ALWAYS | COPY | - |
| Energy Window Range Sequence | 0054,0013 | SQ | | VNAP | COPY | - |
| >Energy Window Lower Limit | 0054,0014 | DS | | VNAP | COPY | - |
| >Energy Window Upper Limit | 0054,0015 | DS | | VNAP | COPY | - |
| Axial Acceptance | 0054,1200 | DS | | ANAPEV | COPY | - |
| Reprojection Method | 0054,1004 | CS | | ANAPCV | COPY | - |
| Number or R R Intervals | 0054,0061 | US | | ANAP | COPY | - |
| Reconstruction Diameter | 0018,1100 | DS | | ANAPEV | COPY | - |
| Convolution Kernel | 0018,1210 | SH | | ANAPEV | COPY | - |
| Detector lines of Response used | 0054,1104 | LO | | ANAPEV | COPY | - |
| Gantry/Detector Slew | 0018,1121 | DS | | ANAPEV | COPY | - |
| Collimator/grid Name | 0018,1180 | SH | | ANAPEV | COPY | - |
| Detector Element Size | 0054,1203 | DS | | ANAPEV | COPY | - |
| Secondary Counts Type | 0054,1220 | CS | | ANAPEV | COPY | - |

Table 39: Frame of Reference Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|--------|---------|
| Frame of Reference UID | 0020,0052 | UI | | ALWAYS | AUTO | - |
| Position Reference Indicator | 0020,1040 | LO | | VNAP | AUTO | - |

Table 40: General Equipment Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---------------------|-----------|----|--------------------|-------------------|------------------|---------------------|
| Manufacturer | 0008,0070 | LO | Philips Healthcare | VNAP | AUTO | Philips Healthcare. |
| Institution Name | 0008,0080 | LO | | ANAPCV | AUTO, COPY, USER | - |
| Institution Address | 0008,0081 | ST | | ANAPCV | AUTO, COPY, USER | - |

| | | | | | | |
|-------------------------------|-----------|-------|-----|--------|------------|-----|
| Institutional Department Name | 0008,1040 | LO | | ANAPCV | AUTO, COPY | - |
| Manufacturer's Model Name | 0008,1090 | LO | | VNAP | AUTO, COPY | - |
| Software Version (s) | 0018,1020 | LO | 4.0 | VNAP | COPY | 4.0 |
| Station Name | 0008,1010 | SH | | ALWAYS | COPY | - |
| Spatial Resolution | 0018,1050 | DS | | ANAPCV | COPY | - |
| Date of Last Calibration | 0018,1200 | DA | | ANAPCV | COPY | - |
| Time of Last Calibration | 0018,1201 | TM | | ANAPCV | COPY | - |
| Pixel Padding Value | 0028,0120 | US/SS | | ANAPCV | COPY | - |

Table 41: General Image Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|-------|-------|-------------------|------------|---------|
| Acquisition Date | 0008,0022 | DA | | ANAPCV | COPY | - |
| Content Date | 0008,0023 | DA | | ANAPCV | COPY | - |
| Acquisition Datetime | 0008,002A | DT | | ANAPCV | COPY | - |
| Acquisition Time | 0008,0032 | TM | | ANAPCV | COPY | - |
| Content Time | 0008,0033 | TM | | ANAPCV | COPY | - |
| Referenced Image Sequence | 0008,1140 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ANAPCV | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ANAPCV | COPY | - |
| Derivation Description | 0008,2111 | ST | | ANAPCV | COPY | - |
| Source Image Sequence | 0008,2112 | SQ | | ANAPCV | COPY | - |
| Derivation Code Sequence | 0008,9215 | SQ | | ANAPCV | COPY | - |
| Instance Number | 0020,0013 | IS | | VNAP | AUTO, COPY | - |
| Patient Orientation | 0020,0020 | CS | | ANAPCV | COPY | - |
| Images in Acquisition | 0020,1002 | IS | | ANAPCV | COPY | - |
| Quality Control Image | 0028,0300 | CS | | ANAPCV | COPY | - |
| Burned In Annotation | 0028,0301 | CS | | ANAPCV | COPY | - |
| Lossy Image Compression | 0028,2110 | CS | | ANAPCV | COPY | - |
| Lossy Image Compression Ratio | 0028,2112 | DS | | ANAPCV | COPY | - |
| Icon Image Sequence | 0088,0200 | SQ | | ANAPCV | COPY | - |
| >Samples per Pixel | 0028,0002 | US | | ALWAYS | COPY | - |
| >Photometric Interpretation | 0028,0004 | CS | | ALWAYS | COPY | - |
| >Planar Configuration | 0028,0006 | US | | ANAP | COPY | - |
| >Rows | 0028,0010 | US | | ALWAYS | COPY | - |
| >Columns | 0028,0011 | US | | ALWAYS | COPY | - |
| >Pixel Aspect Ratio | 0028,0034 | IS | | ANAP | COPY | - |
| >Bits Allocated | 0028,0100 | US | | ALWAYS | COPY | - |
| >Bits Stored | 0028,0101 | US | | ALWAYS | COPY | - |
| >High Bit | 0028,0102 | US | | ALWAYS | COPY | - |
| >Pixel Representation | 0028,0103 | US | | ALWAYS | COPY | - |
| >Smallest Image Pixel Value | 0028,0106 | US/SS | | ANAPCV | COPY | - |
| >Largest Image Pixel Value | 0028,0107 | US/SS | | ANAPCV | COPY | - |
| >Red Palette Color Lookup Table Descriptor | 0028,1101 | US/SS | | ANAP | COPY | - |
| >Green Palette Color Lookup Table Descriptor | 0028,1102 | US/SS | | ANAP | COPY | - |
| >Blue Palette Color Lookup Table Descriptor | 0028,1103 | US/SS | | ANAP | COPY | - |

| | | | | | | |
|--|-----------|-------|--|--------|------|---|
| >Red Palette Color Lookup Table Data | 0028,1201 | OW | | ANAP | COPY | - |
| >Green Palette Color Lookup Table Data | 0028,1202 | OW | | ANAP | COPY | - |
| >Blue Palette Color Lookup Table Data | 0028,1203 | OW | | ANAP | COPY | - |
| >Pixel Data | 7FE0,0010 | OW/OB | | ANAP | COPY | - |
| Presentation LUT Shape | 2050,0020 | CS | | ANAPCV | COPY | - |

Table 42: Image Plane Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------|-----------|----|-------|-------------------|--------|---------|
| Slice Thickness | 0018,0050 | DS | | VNAP | AUTO | - |
| Image Position (Patient) | 0020,0032 | DS | | ALWAYS | AUTO | - |
| Image Orientation (Patient) | 0020,0037 | DS | | ALWAYS | AUTO | - |
| Pixel Spacing | 0028,0030 | DS | | ALWAYS | AUTO | - |
| Slice Location | 0020,1041 | DS | | ANAPEV | COPY | - |

Table 43: Image Pixel Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---|-----------|-------|-------|-------------------|------------|---------|
| Planar Configuration | 0028,0006 | US | | ANAP | COPY | - |
| Rows | 0028,0010 | US | | ALWAYS | AUTO, COPY | - |
| Columns | 0028,0011 | US | | ALWAYS | AUTO, COPY | - |
| Pixel Aspect Ratio | 0028,0034 | IS | | ANAP | COPY | - |
| Pixel Representation | 0028,0103 | US | | ALWAYS | AUTO, COPY | - |
| Smallest Image Pixel Value | 0028,0106 | US/SS | | ANAPCV | COPY | - |
| Largest Image Pixel Value | 0028,0107 | US/SS | | ANAPCV | COPY | - |
| Red Palette Color Lookup Table Descriptor | 0028,1101 | US/SS | | ANAP | COPY | - |
| Green Palette Color Lookup Table Descriptor | 0028,1102 | US/SS | | ANAP | COPY | - |
| Blue Palette Color Lookup Table Descriptor | 0028,1103 | US/SS | | ANAP | COPY | - |
| Red Palette Color Lookup Table Data | 0028,1201 | OW | | ANAP | COPY | - |
| Green Palette Color Lookup Table Data | 0028,1202 | OW | | ANAP | COPY | - |
| Blue Palette Color Lookup Table Data | 0028,1203 | OW | | ANAP | COPY | - |
| Pixel Data | 7FE0,0010 | OW/OB | | ALWAYS | AUTO, COPY | - |

Table 44: NM/PET Patient Orientation Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---|-----------|----|-------|-------------------|--------|---------|
| Patient Gantry Relationship Code Sequence | 0054,0414 | SQ | | VNAP | COPY | - |
| >Code Meaning | 0008,0104 | LO | | ALWAYS | COPY | - |

| | | | | | | |
|---|-----------|----|--|--------|------|---|
| >Code Value | 0008,0100 | SH | | ALWAYS | COPY | - |
| >Coding Scheme Designator | 0008,0102 | SH | | ALWAYS | COPY | - |
| Patient Orientation Code Sequence | 0054,0410 | SQ | | VNAP | COPY | - |
| >Patient Orientation Modifier Code Sequence | 0054,0412 | SQ | | ALWAYS | COPY | - |
| >>Code Meaning | 0008,0104 | LO | | ALWAYS | COPY | - |
| >>Code Value | 0008,0100 | SH | | ALWAYS | COPY | - |
| >>Coding Scheme Designator | 0008,0102 | SH | | ALWAYS | COPY | - |
| >Code Meaning | 0008,0104 | LO | | ALWAYS | COPY | - |
| >Code Value | 0008,0100 | SH | | ALWAYS | COPY | - |
| >Coding Scheme Designator | 0008,0102 | SH | | ALWAYS | COPY | - |

Table 45: PET Isotope Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|--------|---------|
| Radiopharmaceutical Information Sequence | 0054,0016 | SQ | | VNAP | COPY | - |
| >Radionuclide Half Life | 0018,1075 | DS | | ANAPEV | COPY | - |
| >Radionuclide Positron Fraction | 0018,1076 | DS | | ANAPEV | COPY | - |
| >Radionuclide Total Dose | 0018,1074 | DS | | ANAPEV | COPY | - |
| >Radiopharmaceutical | 0018,0031 | LO | | ANAPEV | COPY | - |
| >Radiopharmaceutical Start Time | 0018,1072 | TM | | ANAPEV | COPY | - |
| >Radiopharmaceutical Volume | 0018,1071 | DS | | ANAPEV | COPY | - |
| >Radionuclide Code Sequence | 0054,0300 | SQ | | ANAPEV | COPY | - |
| >>Code Meaning | 0008,0104 | LO | | ALWAYS | COPY | - |
| >>Code Value | 0008,0100 | SH | | ALWAYS | COPY | - |
| >>Coding Scheme Designator | 0008,0102 | SH | | ALWAYS | COPY | - |
| >Radiopharmaceutical Code Sequence | 0054,0304 | SQ | | ANAPEV | COPY | - |
| >>Code Meaning | 0008,0104 | LO | | VNAP | COPY | - |
| >>Code Value | 0008,0100 | SH | | VNAP | COPY | - |

| | | | | | | |
|--|-----------|----|--|--------|------|---|
| >>Coding Scheme Designator | 0008,0102 | SH | | VNAP | COPY | - |
| Interventional Drug Information Sequence | 0018,0026 | SQ | | ANAPEV | COPY | - |

Table 46: PET Multi-Gated Acquisition Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------|-----------|----|-------|-------------------|--------|---------|
| Beat Rejection Flag | 0018,1080 | CS | | VNAP | COPY | - |
| Trigger Source or Type | 0018,1061 | LO | | ANAPCV | COPY | - |
| PVC Rejection | 0018,1085 | LO | | ANAPCV | COPY | - |
| Skip Beats | 0018,1086 | IS | | ANAPCV | COPY | - |
| Heart Rate | 0018,1088 | IS | | ANAPCV | COPY | - |
| Cardiac Framing Type | 0018,1064 | LO | | ANAPCV | COPY | - |

Table 47: PET Image Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--------------------------------------|-----------|----|--------------------------------|-------------------|--------|---------|
| Image Type | 0008,0008 | CS | Value 1: DERIVED, SECONDARY | ALWAYS | AUTO | - |
| Samples per Pixel | 0028,0002 | US | | ALWAYS | AUTO | - |
| Photometric Interpretation | 0028,0004 | CS | | ALWAYS | AUTO | - |
| Bits Allocated | 0028,0100 | US | | ALWAYS | AUTO | - |
| Bits Stored | 0028,0101 | US | | ALWAYS | AUTO | - |
| High Bit | 0028,0102 | US | | ALWAYS | AUTO | - |
| Rescale Intercept | 0028,1052 | DS | | ALWAYS | AUTO | - |
| Rescale Slope | 0028,1053 | DS | | ALWAYS | AUTO | - |
| Acquisition Date | 0008,0022 | DA | | VNAP | COPY | - |
| Acquisition Time | 0008,0032 | TM | | VNAP | COPY | - |
| Actual Frame Duration | 0018,1242 | IS | | VNAP | COPY | - |
| Decay Factor | 0054,1321 | DS | | ANAP | COPY | - |
| Dose Calibration Factor | 0054,1322 | DS | | ANAPEV | COPY | - |
| Frame Reference Time | 0054,1300 | DS | | ALWAYS | COPY | - |
| Trigger Time | 0018,1060 | DS | | ANAP | COPY | - |
| Frame Time | 0018,1063 | DS | | ANAP | COPY | - |
| Low R R Value | 0018,1081 | IS | | ANAP | COPY | - |
| High R R Value | 0018,1082 | IS | | ANAP | COPY | - |
| Lossy Image Compression | 0028,2110 | CS | | ANAP | COPY | - |
| Image Index | 0054,1330 | US | | ALWAYS | COPY | - |
| Nominal Interval | 0018,1062 | IS | | ANAPEV | COPY | - |
| Intervals Acquired | 0018,1083 | IS | | ANAPEV | COPY | - |
| Intervals Rejected | 0018,1084 | IS | | ANAPEV | COPY | - |
| Primary (Prompts) Counts Accumulated | 0054,1310 | IS | | ANAPEV | COPY | - |
| Secondary Counts Accumulated | 0054,1311 | IS | | ANAPEV | COPY | - |
| Slice Sensitivity Factor | 0054,1320 | DS | | ANAPEV | COPY | - |
| Scatter Fraction Factor | 0054,1323 | DS | | ANAPEV | COPY | - |
| Dead Time Factor | 0054,1324 | DS | | ANAPEV | COPY | - |
| View Code Sequence | 0054,0220 | SQ | | ANAPEV | COPY | - |
| Slice Progression Direction | 0054,0500 | CS | | ANAPEV | COPY | - |
| Anatomic Region Sequence | 0008,2218 | SQ | | ANAPEV | COPY | - |

Table 48: Overlay Plane Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------|-----------|-------|-------|-------------------|--------|---------|
| Overlay Rows | 6000,0010 | US | | ALWAYS | AUTO | - |
| Overlay Columns | 6000,0011 | US | | ALWAYS | AUTO | - |
| Overlay Description | 6000,0022 | LO | | ANAPCV | AUTO | - |
| Overlay Type | 6000,0040 | CS | | ALWAYS | AUTO | - |
| Overlay Subtype | 6000,0045 | LO | | ANAPCV | AUTO | - |
| Overlay Origin | 6000,0050 | SS | | ALWAYS | AUTO | - |
| Overlay Bits Allocated | 6000,0100 | US | | ALWAYS | AUTO | - |
| Overlay Bit Position | 6000,0102 | US | | ALWAYS | AUTO | - |
| ROI Area | 6000,1301 | IS | | ANAPCV | AUTO | - |
| ROI Mean | 6000,1302 | DS | | ANAPCV | AUTO | - |
| ROI Standard Deviation | 6000,1303 | DS | | ANAPCV | AUTO | - |
| Overlay Label | 6000,1500 | LO | | ANAPCV | AUTO | - |
| Overlay Data | 6000,3000 | OW/OB | | ANAP | AUTO | - |

Table 49: VOI LUT Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------------|-----------|----|-------|-------------------|--------|---------|
| Window Center | 0028,1050 | DS | | ALWAYS | AUTO | - |
| Window Width | 0028,1051 | DS | | ANAP | AUTO | - |
| Window Center & Width Explanation | 0028,1055 | LO | | ANAPCV | AUTO | - |
| VOI LUT Sequence | 0028,3010 | SQ | | ANAPCV | AUTO | - |

Table 50: SOP Common Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---------------------------------|-----------|----|-----------------------------|-------------------|------------|----------------------|
| Instance Creation Date | 0008,0012 | DA | | VNAP | AUTO, COPY | - |
| Instance Creation Time | 0008,0013 | TM | | VNAP | AUTO, COPY | - |
| Instance Creator UID | 0008,0014 | UI | | VNAP | AUTO, COPY | - |
| Specific Character Set | 0008,0005 | CS | | ANAP | AUTO, COPY | - |
| SOP Class UID | 0008,0016 | UI | 1.2.840.10008.5.1.4.1.1.128 | ALWAYS | AUTO, COPY | - |
| SOP Instance UID | 0008,0018 | UI | | ALWAYS | AUTO | - |
| Timezone Offset From UTC | 0008,0201 | SH | | ANAPCV | COPY | - |
| Contributing Equipment Sequence | 0018,A001 | SQ | | ANAP | AUTO | - |
| >Manufacturer | 0008,0070 | LO | Philips | ANAP | FIXED | PHILIPS. |
| >Institution Address | 0008,0081 | ST | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Station Name | 0008,1010 | SH | IntelliSpace Portal | ANAP | CONFIG | IntelliSpace Portal. |
| >Institutional Department Name | 0008,1040 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Software Version(s) | 0018,1020 | LO | 4.0 | ANAP | FIXED | 4.0 |
| >Device Serial Number | 0018,1000 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Institution Name | 0008,0080 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Manufacturer's Model Name | 0008,1090 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |

| | | | | | | |
|-------------------------------------|-----------|----|----------------------|------|-------|----------------------|
| >Purpose of Reference Code Sequence | 0040,A170 | SQ | | ANAP | FIXED | - |
| >>Code Value | 0008,0100 | SH | 109102 | ANAP | FIXED | 109102 |
| >>Coding Scheme Designator | 0008,0102 | SH | DCM | ANAP | FIXED | DCM |
| >>Code Meaning | 0008,0104 | LO | Processing Equipment | ANAP | FIXED | Processing Equipment |

Table 51: Extended DICOM and private attributes for Positron Emission Tomography Image Storage SOP Class Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------------|-----------|-------|-------|-------------------|--------|---------|
| Frame Number | 2001,107C | IS | | ALWAYS | AUTO | - |
| Slice Index | 2001,10DE | DS | | ALWAYS | AUTO | - |
| Referenced Overlay Sequence | 0008,1130 | SQ | | ANAPEV | COPY | - |
| Reference Curve Sequence | 0008,1145 | SQ | | ANAPEV | COPY | - |
| Primary Anatomic Structure Sequence | 0008,2228 | SQ | | ANAPEV | COPY | - |
| Curve Dimensions | 50xx,0005 | US | | ALWAYS | AUTO | - |
| Number of Points | 50xx,0010 | US | | ALWAYS | AUTO | - |
| Type of Data | 50xx,0020 | CS | | ALWAYS | AUTO | - |
| Data Value Representation | 50xx,0103 | US | | ALWAYS | AUTO | - |
| Curve Data | 50xx,3000 | OB/OW | | ALWAYS | AUTO | - |
| Curve Description | 50xx,0022 | LO | | ANAPCV | AUTO | - |
| Axis Units | 50xx,0030 | SH | | ANAPCV | AUTO | - |
| Axis Labels | 50xx,0040 | SH | | ANAPCV | AUTO | - |
| Minimum Coordinate Value | 50xx,0104 | US | | ANAPCV | AUTO | - |
| Maximum Coordinate Value | 50xx,0105 | US | | ANAPCV | AUTO | - |
| Curve Range | 50xx,0106 | SH | | ANAPCV | AUTO | - |
| Curve Data Descriptor | 50xx,0110 | US | | ANAP | AUTO | - |
| Coordinate Start Value | 50xx,0112 | US | | ANAP | AUTO | - |
| Coordinate Step Value | 50xx,0114 | US | | ANAP | AUTO | - |
| Curve Label | 50xx,2500 | LO | | ANAPCV | AUTO | - |
| Curve Referenced Overlay Sequence | 50xx,2600 | SQ | | ANAPCV | AUTO | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ANAP | AUTO | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ANAP | AUTO | - |
| >Referenced Overlay Group | 50xx,2610 | US | | ANAP | AUTO | - |

2.1.2.4. CT Image Storage SOP Class

Table 52: IOD of Created CT Image Storage SOP Class Instances

| Information Entity | Module | Presence Of Module |
|--------------------|---------------------------|--------------------|
| Patient | Patient Module | ALWAYS |
| Study | General Study Module | ALWAYS |
| Study | Patient Study Module | CONDITIONAL |
| Series | General Series Module | ALWAYS |
| Frame of Reference | Frame of Reference Module | ALWAYS |
| Equipment | General Equipment Module | ALWAYS |
| Image | General Image Module | ALWAYS |
| Image | Image Plane Module | ALWAYS |
| Image | Image Pixel Module | ALWAYS |
| Image | Contrast/Bolus Module | CONDITIONAL |

| | | |
|-------|---------------------------------------|-------------|
| Image | Device Module | CONDITIONAL |
| Image | CT Image Module | ALWAYS |
| Image | Overlay Plane Module | CONDITIONAL |
| Image | VOI LUT Module | CONDITIONAL |
| Image | SOP Common Module | ALWAYS |
| | Extended Dicom and Private attributes | CONDITIONAL |

Table 53: Patient Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|---------------|---------|
| Referenced Patient Sequence | 0008,1120 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Patient's Name | 0010,0010 | PN | | VNAP | COPY | - |
| Patient ID | 0010,0020 | LO | | VNAP | COPY | - |
| Patient's Birth Date | 0010,0030 | DA | | VNAP | COPY, USER | - |
| Patient's Birth Time | 0010,0032 | TM | | ANAPCV | COPY, USER | - |
| Patient's Sex | 0010,0040 | CS | | VNAP | COPY, USER | - |
| Other Patient IDs | 0010,1000 | LO | | ANAPCV | COPY, USER | - |
| Other Patient Names | 0010,1001 | PN | | ANAPCV | COPY, USER | - |
| Ethnic Group | 0010,2160 | SH | | ANAPCV | COPY | - |
| Patient Comments | 0010,4000 | LT | | ANAPCV | COPY | - |
| Issuer of Patient ID | 0010,0021 | LO | | ANAPCV | COPY, USER | - |

Table 54: General Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|---------------|---------|
| Study Date | 0008,0020 | DA | | VNAP | COPY | - |
| Study Time | 0008,0030 | TM | | VNAP | COPY | - |
| Accession Number | 0008,0050 | SH | | VNAP | COPY, USER | - |
| Referring Physician's Name | 0008,0090 | PN | | VNAP | COPY, USER | - |
| Referring Physician Identification Sequence | 0008,0096 | SQ | | ANAPCV | COPY | - |
| Study Description | 0008,1030 | LO | | ANAPCV | COPY, USER | - |
| Procedure Code Sequence | 0008,1032 | SQ | | ANAPCV | COPY | - |
| Physician(s) of Record | 0008,1048 | PN | | ANAPCV | COPY | - |
| Physician(s) of Record Identification Sequence | 0008,1049 | SQ | | ANAPCV | COPY | - |
| Name of Physician(s) Reading Study | 0008,1060 | PN | | ANAPCV | COPY | - |
| Physician(s) Reading Study Identification Sequence | 0008,1062 | SQ | | ANAPCV | COPY | - |
| Referenced Study Sequence | 0008,1110 | SQ | | ANAPCV | COPY | - |
| Study Instance UID | 0020,000D | UI | | ALWAYS | AUTO, COPY | - |

| | | | | | | |
|----------|-----------|----|--|------|---------------|---|
| Study ID | 0020,0010 | SH | | VNAP | AUTO, COPY | - |
|----------|-----------|----|--|------|---------------|---|

Table 55: Patient Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------------|-----------|----|-------|-------------------|---------------|---------|
| Admitting Diagnoses Description | 0008,1080 | LO | | ANAPCV | COPY | - |
| Admitting Diagnoses Code Sequence | 0008,1084 | SQ | | ANAPCV | COPY | - |
| Patient's Age | 0010,1010 | AS | | ANAPCV | COPY, USER | - |
| Patient's Size | 0010,1020 | DS | | ANAPCV | COPY, USER | - |
| Patient's Weight | 0010,1030 | DS | | ANAPCV | COPY, USER | - |
| Occupation | 0010,2180 | SH | | ANAPCV | COPY | - |
| Additional Patient History | 0010,21B0 | LT | | ANAPCV | COPY | - |

Table 56: General Series Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|-------|-------|-------------------|---------------|---------|
| Series Date | 0008,0021 | DA | | ANAPCV | AUTO | - |
| Series Time | 0008,0031 | TM | | ANAPCV | AUTO | - |
| Modality | 0008,0060 | CS | | ALWAYS | COPY | - |
| Series Description | 0008,103E | LO | | ANAPCV | COPY, USER | - |
| Performing Physician's Name | 0008,1050 | PN | | ANAPCV | COPY | - |
| Operators' Name | 0008,1070 | PN | | ANAPCV | COPY, USER | - |
| Operator Identification Sequence | 0008,1072 | SQ | | ANAPCV | COPY | - |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Body Part Examined | 0018,0015 | CS | | ANAPCV | COPY | - |
| Protocol Name | 0018,1030 | LO | | ANAPCV | COPY | - |
| Patient Position | 0018,5100 | CS | | ANAPCV | COPY | - |
| Series Instance UID | 0020,000E | UI | | ALWAYS | AUTO | - |
| Series Number | 0020,0011 | IS | | VNAP | AUTO | - |
| Laterality | 0020,0060 | CS | | ANAPCV | COPY | - |
| Smallest Pixel Value in Series | 0028,0108 | US/SS | | ANAPCV | AUTO, COPY | - |
| Largest Pixel Value in Series | 0028,0109 | US/SS | | ANAPCV | AUTO, COPY | - |
| Request Attributes Sequence | 0040,0275 | SQ | | ANAPCV | COPY | - |
| Performed Procedure Step Start Date | 0040,0244 | DA | | ANAPCV | COPY | - |
| Performed Procedure Step Start Time | 0040,0245 | TM | | ANAPCV | COPY | - |
| Performed Procedure Step Description | 0040,0254 | LO | | ANAPCV | COPY | - |
| Performed Protocol Code Sequence | 0040,0260 | SQ | | ANAPCV | COPY | - |

| | | | | | | |
|--|-----------|----|--|--------|------|---|
| Comments on the Performed Procedure Step | 0040,0280 | ST | | ANAPCV | COPY | - |
|--|-----------|----|--|--------|------|---|

Table 57: Frame of Reference Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|--------|---------|
| Frame of Reference UID | 0020,0052 | UI | | ALWAYS | COPY | - |
| Position Reference Indicator | 0020,1040 | LO | | VNAP | COPY | - |

Table 58: General Equipment Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------|-----------|-------|-------|-------------------|------------|---------|
| Manufacturer | 0008,0070 | LO | | VNAP | COPY | - |
| Institution Name | 0008,0080 | LO | | ANAPCV | COPY, USER | - |
| Institution Address | 0008,0081 | ST | | ANAPCV | COPY, USER | - |
| Station Name | 0008,1010 | SH | | ANAPCV | COPY | - |
| Institutional Department Name | 0008,1040 | LO | | ANAPCV | COPY | - |
| Manufacturer's Model Name | 0008,1090 | LO | | ANAPCV | COPY | - |
| Device Serial Number | 0018,1000 | LO | | ANAPCV | COPY | - |
| Software Version(s) | 0018,1020 | LO | 4.0 | ANAPCV | COPY | 4.0 |
| Spatial Resolution | 0018,1050 | DS | | ANAPCV | COPY | - |
| Date of Last Calibration | 0018,1200 | DA | | ANAPCV | COPY | - |
| Time of Last Calibration | 0018,1201 | TM | | ANAPCV | COPY | - |
| Pixel Padding Value | 0028,0120 | US/SS | | ANAPCV | COPY | - |

Table 59: General Image Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|-------|-------------------|------------|---------|
| Acquisition Date | 0008,0022 | DA | | ANAPCV | COPY | - |
| Content Date | 0008,0023 | DA | | ANAPCV | COPY | - |
| Acquisition Datetime | 0008,002A | DT | | ANAPCV | COPY | - |
| Acquisition Time | 0008,0032 | TM | | ANAPCV | COPY | - |
| Content Time | 0008,0033 | TM | | ANAPCV | COPY | - |
| Referenced Image Sequence | 0008,1140 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ANAPCV | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ANAPCV | COPY | - |
| Derivation Description | 0008,2111 | ST | | ANAPCV | COPY | - |
| Source Image Sequence | 0008,2112 | SQ | | ANAPCV | COPY | - |
| Derivation Code Sequence | 0008,9215 | SQ | | ANAPCV | COPY | - |
| Instance Number | 0020,0013 | IS | | VNAP | AUTO, COPY | - |
| Patient Orientation | 0020,0020 | CS | | ANAPCV | COPY | - |
| Images in Acquisition | 0020,1002 | IS | | ANAPCV | COPY | - |
| Quality Control Image | 0028,0300 | CS | | ANAPCV | COPY | - |
| Burned In Annotation | 0028,0301 | CS | | ANAPCV | COPY | - |
| Lossy Image Compression | 0028,2110 | CS | | ANAPCV | COPY | - |
| Lossy Image Compression Ratio | 0028,2112 | DS | | ANAPCV | COPY | - |
| Icon Image Sequence | 0088,0200 | SQ | | ANAPCV | COPY | - |
| >Samples per Pixel | 0028,0002 | US | | ALWAYS | COPY | - |
| >Photometric Interpretation | 0028,0004 | CS | | ALWAYS | COPY | - |

| | | | | | | |
|--|-----------|-------|--|--------|------|---|
| >Planar Configuration | 0028,0006 | US | | ANAP | COPY | - |
| >Rows | 0028,0010 | US | | ALWAYS | COPY | - |
| >Columns | 0028,0011 | US | | ALWAYS | COPY | - |
| >Pixel Aspect Ratio | 0028,0034 | IS | | ANAP | COPY | - |
| >Bits Allocated | 0028,0100 | US | | ALWAYS | COPY | - |
| >Bits Stored | 0028,0101 | US | | ALWAYS | COPY | - |
| >High Bit | 0028,0102 | US | | ALWAYS | COPY | - |
| >Pixel Representation | 0028,0103 | US | | ALWAYS | COPY | - |
| >Smallest Image Pixel Value | 0028,0106 | US/SS | | ANAPCV | COPY | - |
| >Largest Image Pixel Value | 0028,0107 | US/SS | | ANAPCV | COPY | - |
| >Red Palette Color Lookup Table Descriptor | 0028,1101 | US/SS | | ANAP | COPY | - |
| >Green Palette Color Lookup Table Descriptor | 0028,1102 | US/SS | | ANAP | COPY | - |
| >Blue Palette Color Lookup Table Descriptor | 0028,1103 | US/SS | | ANAP | COPY | - |
| >Red Palette Color Lookup Table Data | 0028,1201 | OW | | ANAP | COPY | - |
| >Green Palette Color Lookup Table Data | 0028,1202 | OW | | ANAP | COPY | - |
| >Blue Palette Color Lookup Table Data | 0028,1203 | OW | | ANAP | COPY | - |
| >Pixel Data | 7FE0,0010 | OW/OB | | ANAP | COPY | - |
| Presentation LUT Shape | 2050,0020 | CS | | ANAPCV | COPY | - |

Table 60: Image Plane Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------|-----------|----|-------|-------------------|--------|---------|
| Slice Thickness | 0018,0050 | DS | | VNAP | COPY | - |
| Image Position (Patient) | 0020,0032 | DS | | ANAP | COPY | - |
| Image Orientation (Patient) | 0020,0037 | DS | | ANAP | COPY | - |
| Slice Location | 0020,1041 | DS | | ANAPCV | COPY | - |
| Pixel Spacing | 0028,0030 | DS | | ANAP | COPY | - |

Table 61: Image Pixel Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---|-----------|-------|-------|-------------------|------------|---------|
| Planar Configuration | 0028,0006 | US | | ANAP | COPY | - |
| Rows | 0028,0010 | US | | ALWAYS | AUTO, COPY | - |
| Columns | 0028,0011 | US | | ALWAYS | AUTO, COPY | - |
| Pixel Aspect Ratio | 0028,0034 | IS | | ANAP | COPY | - |
| Pixel Representation | 0028,0103 | US | | ALWAYS | COPY | - |
| Smallest Image Pixel Value | 0028,0106 | US/SS | | ANAPCV | COPY | - |
| Largest Image Pixel Value | 0028,0107 | US/SS | | ANAPCV | COPY | - |
| Red Palette Color Lookup Table Descriptor | 0028,1101 | US/SS | | ANAP | COPY | - |
| Green Palette Color Lookup Table Descriptor | 0028,1102 | US/SS | | ANAP | COPY | - |
| Blue Palette Color Lookup Table Descriptor | 0028,1103 | US/SS | | ANAP | COPY | - |
| Red Palette Color Lookup Table Data | 0028,1201 | OW | | ANAP | COPY | - |

| | | | | | | |
|---------------------------------------|-----------|-------|--|------|------|---|
| Green Palette Color Lookup Table Data | 0028,1202 | OW | | ANAP | COPY | - |
| Blue Palette Color Lookup Table Data | 0028,1203 | OW | | ANAP | COPY | - |
| Pixel Data | 7FE0,0010 | OW/OB | | ANAP | COPY | - |

Table 62: Contrast/Bolus Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|--------|---------|
| Contrast/Bolus Agent | 0018,0010 | LO | | VNAP | COPY | - |
| Contrast/Bolus Agent Sequence | 0018,0012 | SQ | | ANAPCV | COPY | - |
| Contrast/Bolus Administration Route Sequence | 0018,0014 | SQ | | ANAPCV | COPY | - |
| >Additional Drug Sequence | 0018,002A | SQ | | ANAP | COPY | - |
| Contrast/Bolus Route | 0018,1040 | LO | | ANAPCV | COPY | - |
| Contrast/Bolus Volume | 0018,1041 | DS | | ANAPCV | COPY | - |
| Contrast/Bolus Start Time | 0018,1042 | TM | | ANAPCV | COPY | - |
| Contrast/Bolus Stop Time | 0018,1043 | TM | | ANAPCV | COPY | - |
| Contrast/Bolus Total Dose | 0018,1044 | DS | | ANAPCV | COPY | - |
| Contrast Flow Rate | 0018,1046 | DS | | ANAPCV | COPY | - |
| Contrast Flow Duration | 0018,1047 | DS | | ANAPCV | COPY | - |
| Contrast/Bolus Ingredient | 0018,1048 | CS | | ANAPCV | COPY | - |
| Contrast/Bolus Ingredient Concentration | 0018,1049 | DS | | ANAPCV | COPY | - |

Table 63: Device Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--------------------------------------|-----------|----|-------|-------------------|--------|---------|
| Device Sequence | 0050,0010 | SQ | | ANAPCV | COPY | - |
| >Manufacturer | 0008,0070 | LO | | ANAPCV | COPY | - |
| >Manufacturer's Model Name | 0008,1090 | LO | | ANAPCV | COPY | - |
| >Device Serial Number | 0018,1000 | LO | | ANAPCV | COPY | - |
| >Device ID | 0018,1003 | LO | | ANAPCV | COPY | - |
| >Device Length | 0050,0014 | DS | | ANAPCV | COPY | - |
| >Device Diameter | 0050,0016 | DS | | ANAPCV | COPY | - |
| >Device Diameter Units | 0050,0017 | CS | | ANAPCV | COPY | - |
| >Device Volume | 0050,0018 | DS | | ANAPCV | COPY | - |
| >Intermarker Distance | 0050,0019 | DS | | ANAPCV | COPY | - |
| >Device Description | 0050,0020 | LO | | ANAPCV | COPY | - |
| >Code Value | 0008,0100 | SH | | ALWAYS | COPY | - |
| >Coding Scheme Designator | 0008,0102 | SH | | ALWAYS | COPY | - |
| >Coding Scheme Version | 0008,0103 | SH | | ANAP | COPY | - |
| >Code Meaning | 0008,0104 | LO | | ALWAYS | COPY | - |
| >Mapping Resource | 0008,0105 | CS | | ANAP | COPY | - |
| >Context Group Version | 0008,0106 | DT | | ANAP | COPY | - |
| >Context Group Local Version | 0008,0107 | DT | | ANAP | COPY | - |
| >Context Group Extension Flag | 0008,010B | CS | | ANAPCV | COPY | - |
| >Context Group Extension Creator UID | 0008,010D | UI | | ANAP | COPY | - |
| >Context Identifier | 0008,010F | CS | | ANAPCV | COPY | - |

Table 64: CT Image Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------|-----------|----|-------|-------------------|------------|---------|
| Image Type | 0008,0008 | CS | | ALWAYS | AUTO | - |
| Scan Options | 0018,0022 | CS | | ANAPCV | AUTO, COPY | - |
| KVP | 0018,0060 | DS | | VNAP | COPY | - |
| Data Collection Diameter | 0018,0090 | DS | | ANAPCV | COPY | - |
| Reconstruction Diameter | 0018,1100 | DS | | ANAPCV | COPY | - |
| Distance Source to Detector | 0018,1110 | DS | | ANAPCV | COPY | - |
| Distance Source to Patient | 0018,1111 | DS | | ANAPCV | COPY | - |
| Gantry/Detector Tilt | 0018,1120 | DS | | ANAPCV | COPY | - |
| Table Height | 0018,1130 | DS | | ANAPCV | COPY | - |
| Rotation Direction | 0018,1140 | CS | | ANAPCV | COPY | - |
| Exposure Time | 0018,1150 | IS | | ANAPCV | COPY | - |
| X-ray Tube Current | 0018,1151 | IS | | ANAPCV | COPY | - |
| Exposure | 0018,1152 | IS | | ANAPCV | COPY | - |
| Exposure in μ As | 0018,1153 | IS | | ANAPCV | COPY | - |
| Filter Type | 0018,1160 | SH | | ANAPCV | COPY | - |
| Generator Power | 0018,1170 | IS | | ANAPCV | COPY | - |
| Focal Spot(s) | 0018,1190 | DS | | ANAPCV | COPY | - |
| Convolution Kernel | 0018,1210 | SH | | ANAPCV | COPY | - |
| Revolution Time | 0018,9305 | FD | | ANAPCV | COPY | - |
| Single Collimation Width | 0018,9306 | FD | | ANAPCV | COPY | - |
| Total Collimation Width | 0018,9307 | FD | | ANAPCV | COPY | - |
| Table Speed | 0018,9309 | FD | | ANAPCV | COPY | - |
| Table Feed per Rotation | 0018,9310 | FD | | ANAPCV | COPY | - |
| Spiral Pitch Factor | 0018,9311 | FD | | ANAPCV | COPY | - |
| Acquisition Number | 0020,0012 | IS | | VNAP | COPY | - |
| Samples per Pixel | 0028,0002 | US | | ALWAYS | COPY | - |
| Photometric Interpretation | 0028,0004 | CS | | ALWAYS | COPY | - |
| Bits Allocated | 0028,0100 | US | | ALWAYS | COPY | - |
| Bits Stored | 0028,0101 | US | | ALWAYS | COPY | - |
| High Bit | 0028,0102 | US | | ALWAYS | COPY | - |
| Rescale Intercept | 0028,1052 | DS | | ALWAYS | COPY | - |
| Rescale Slope | 0028,1053 | DS | | ALWAYS | COPY | - |

Table 65: Overlay Plane Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------|-----------|----|-------|-------------------|--------|---------|
| Overlay Rows | 6000,0010 | US | | ALWAYS | AUTO | - |
| Overlay Columns | 6000,0011 | US | | ALWAYS | AUTO | - |
| Overlay Description | 6000,0022 | LO | | ANAPCV | AUTO | - |
| Overlay Type | 6000,0040 | CS | | ALWAYS | AUTO | - |
| Overlay Subtype | 6000,0045 | LO | | ALWAYS | AUTO | - |
| Overlay Origin | 6000,0050 | SS | | ALWAYS | AUTO | - |
| Overlay Bits Allocated | 6000,0100 | US | | ALWAYS | AUTO | - |
| Overlay Bit Position | 6000,0102 | US | | ALWAYS | AUTO | - |
| ROI Area | 6000,1301 | IS | | ANAPCV | AUTO | - |
| ROI Mean | 6000,1302 | DS | | ANAPCV | AUTO | - |
| ROI Standard Deviation | 6000,1303 | DS | | ANAPCV | AUTO | - |

| | | | | | | |
|---------------|-----------|-------|--|--------|------|---|
| Overlay Label | 6000,1500 | LO | | ANAPCV | AUTO | - |
| Overlay Data | 6000,3000 | OW/OB | | ANAP | AUTO | - |

Table 66: VOI LUT Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------------|-----------|----|-------|-------------------|--------|---------|
| Window Center | 0028,1050 | DS | | ANAPCV | AUTO | - |
| Window Width | 0028,1051 | DS | | ANAPCV | AUTO | - |
| Window Center & Width Explanation | 0028,1055 | LO | | ANAPCV | AUTO | - |
| VOI LUT Sequence | 0028,3010 | SQ | | ANAPCV | AUTO | - |

Table 67: SOP Common Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------------|-----------|----|----------------------|-------------------|------------|----------------------|
| Specific Character Set | 0008,0005 | CS | | ANAP | AUTO, COPY | - |
| Instance Creation Date | 0008,0012 | DA | | ANAPCV | AUTO | - |
| Instance Creation Time | 0008,0013 | TM | | ANAPCV | AUTO | - |
| Instance Creator UID | 0008,0014 | UI | | ANAPCV | AUTO | - |
| SOP Class UID | 0008,0016 | UI | | ALWAYS | COPY | - |
| SOP Instance UID | 0008,0018 | UI | | ALWAYS | AUTO | - |
| Timezone Offset From UTC | 0008,0201 | SH | | ANAPCV | COPY | - |
| Contributing Equipment Sequence | 0018,A001 | SQ | | ANAP | AUTO | - |
| >Manufacturer | 0008,0070 | LO | Philips | ANAP | FIXED | PHILIPS. |
| >Institution Address | 0008,0081 | ST | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Station Name | 0008,1010 | SH | IntelliSpace Portal | ANAP | CONFIG | IntelliSpace Portal. |
| >Institutional Department Name | 0008,1040 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Software Version(s) | 0018,1020 | LO | 4.0 | ANAP | FIXED | 4.0 |
| >Device Serial Number | 0018,1000 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Institution Name | 0008,0080 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Manufacturer's Model Name | 0008,1090 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Purpose of Reference Code Sequence | 0040,A170 | SQ | | ANAP | FIXED | - |
| >>Code Value | 0008,0100 | SH | 109102 | ANAP | FIXED | 109102 |
| >>Coding Scheme Designator | 0008,0102 | SH | DCM | ANAP | FIXED | DCM |
| >>Code Meaning | 0008,0104 | LO | Processing Equipment | ANAP | FIXED | Processing Equipment |

Table 68: Extended DICOM and private attributes for CT Image Storage SOP Class Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---------------------------|-----------|----|-------|-------------------|--------|---------|
| Modalities in Study | 0008,0061 | CS | | ANAPCV | COPY | - |
| Medical Alerts | 0010,2000 | LO | | ANAPCV | COPY | - |
| Allergies | 0010,2110 | LO | | ANAPCV | COPY | - |
| Pregnancy Status | 0010,21C0 | US | | ANAPCV | COPY | - |
| Spacing Between Slices | 0018,0088 | DS | | ANAPCV | COPY | - |
| Heart Rate | 0018,1088 | IS | | ANAPCV | COPY | - |
| Angular Position | 0018,1141 | DS | | ANAPCV | COPY | - |
| Scan Arc | 0018,1143 | DS | | ANAPCV | COPY | - |
| CT Exposure Sequence | 0018,9321 | SQ | | ANAPCV | COPY | - |
| >Exposure Modulation Type | 0018,9323 | CS | | ANAPCV | COPY | - |

| | | | | | | |
|--|-----------|-------|--|--------|------|---|
| >Estimated Dose Saving | 0018,9324 | FD | | ANAPCV | COPY | - |
| >Exposure Time in ms | 0018,9328 | FD | | ANAPCV | COPY | - |
| >X-Ray Tube Current in mA | 0018,9330 | FD | | ANAPCV | COPY | - |
| >Exposure in mAs | 0018,9332 | FD | | ANAPCV | COPY | - |
| Number of Study Related Series | 0020,1206 | IS | | ANAPCV | COPY | - |
| Number of Study Related Instances | 0020,1208 | IS | | ANAPCV | COPY | - |
| Requesting Physician | 0032,1032 | PN | | ANAPCV | COPY | - |
| Requesting Service | 0032,1033 | LO | | ANAPCV | COPY | - |
| Requested Procedure Description | 0032,1060 | LO | | ANAPCV | COPY | - |
| Requested Contrast Agent | 0032,1070 | LO | | ANAPCV | COPY | - |
| Admission ID | 0038,0010 | LO | | ANAPCV | COPY | - |
| Special Needs | 0038,0050 | LO | | ANAPCV | COPY | - |
| Patient State | 0038,0500 | LO | | ANAPCV | COPY | - |
| Confidentiality Constraint on Patient Data Description | 0040,3001 | LO | | ANAPCV | COPY | - |
| DateTime | 0040,A120 | DT | | ANAPCV | COPY | - |
| Curve Dimensions | 50xx,0005 | US | | ALWAYS | AUTO | - |
| Number of Points | 50xx,0010 | US | | ALWAYS | AUTO | - |
| Type of Data | 50xx,0020 | CS | | ALWAYS | AUTO | - |
| Data Value Representation | 50xx,0103 | US | | ALWAYS | AUTO | - |
| Curve Data | 50xx,3000 | OB/OW | | ALWAYS | AUTO | - |
| Curve Description | 50xx,0022 | LO | | ANAPCV | AUTO | - |
| Axis Units | 50xx,0030 | SH | | ANAPCV | AUTO | - |
| Axis Labels | 50xx,0040 | SH | | ANAPCV | AUTO | - |
| Minimum Coordinate Value | 50xx,0104 | US | | ANAPCV | AUTO | - |
| Maximum Coordinate Value | 50xx,0105 | US | | ANAPCV | AUTO | - |
| Curve Range | 50xx,0106 | SH | | ANAPCV | AUTO | - |
| Curve Data Descriptor | 50xx,0110 | US | | ANAP | AUTO | - |
| Coordinate Start Value | 50xx,0112 | US | | ANAP | AUTO | - |
| Coordinate Step Value | 50xx,0114 | US | | ANAP | AUTO | - |
| Curve Label | 50xx,2500 | LO | | ANAPCV | AUTO | - |
| Referenced Overlay Sequence | 50xx,2600 | SQ | | ANAPCV | AUTO | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ANAP | AUTO | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ANAP | AUTO | - |
| >Referenced Overlay Group | 50xx,2610 | US | | ANAP | AUTO | - |

2.1.2.5. Nuclear Medicine Image Storage SOP Class

Table 69: IOD of Created Nuclear Medicine Image Storage SOP Class Instances

| Information Entity | Module | Presence Of Module |
|--------------------|---------------------------|--------------------|
| Patient | Patient Module | ALWAYS |
| Study | General Study Module | ALWAYS |
| Study | Patient Study Module | CONDITIONAL |
| Series | General Series Module | ALWAYS |
| Frame of Reference | Frame of Reference Module | CONDITIONAL |
| Equipment | General Equipment Module | ALWAYS |
| Image | General Image Module | ALWAYS |
| Image | Image Pixel Module | ALWAYS |
| Image | Multi-Frame Module | ALWAYS |

| | | |
|-------|---------------------------------------|-------------|
| Image | NM/PET Patient Orientation Module | ALWAYS |
| Image | NM Image Pixel Module | ALWAYS |
| Image | NM Multi-Frame Module | ALWAYS |
| Image | NM Image Module | ALWAYS |
| Image | NM Isotope Module | ALWAYS |
| Image | NM Detector Module | ALWAYS |
| Image | NM Tomo Acquisition Module | CONDITIONAL |
| Image | NM Multi-gated Acquisition Module | CONDITIONAL |
| Image | NM Reconstruction Module | CONDITIONAL |
| Image | Overlay Plane Module | CONDITIONAL |
| Image | VOI LUT Module | CONDITIONAL |
| Image | SOP Common Module | ALWAYS |
| Image | NM Phase Module | CONDITIONAL |
| | Extended Dicom and Private attributes | CONDITIONAL |

Table 70: Patient Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|---------------|---------|
| Referenced Patient Sequence | 0008,1120 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Patient's Name | 0010,0010 | PN | | VNAP | COPY | - |
| Patient ID | 0010,0020 | LO | | VNAP | COPY | - |
| Patient's Birth Date | 0010,0030 | DA | | VNAP | COPY, USER | - |
| Patient's Birth Time | 0010,0032 | TM | | ANAPCV | COPY, USER | - |
| Patient's Sex | 0010,0040 | CS | | VNAP | COPY, USER | - |
| Other Patient IDs | 0010,1000 | LO | | ANAPCV | COPY, USER | - |
| Other Patient Names | 0010,1001 | PN | | ANAPCV | COPY, USER | - |
| Ethnic Group | 0010,2160 | SH | | ANAPCV | COPY | - |
| Patient Comments | 0010,4000 | LT | | ANAPCV | COPY | - |
| Issuer of Patient ID | 0010,0021 | LO | | ANAPCV | COPY, USER | - |

Table 71: General Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---|-----------|----|-------|-------------------|---------------|---------|
| Study Date | 0008,0020 | DA | | VNAP | COPY | - |
| Study Time | 0008,0030 | TM | | VNAP | COPY | - |
| Accession Number | 0008,0050 | SH | | VNAP | COPY, USER | - |
| Referring Physician's Name | 0008,0090 | PN | | VNAP | COPY, USER | - |
| Referring Physician Identification Sequence | 0008,0096 | SQ | | ANAPCV | COPY | - |
| Study Description | 0008,1030 | LO | | ANAPCV | COPY, USER | - |
| Procedure Code Sequence | 0008,1032 | SQ | | ANAPCV | COPY | - |
| Physician(s) of Record | 0008,1048 | PN | | ANAPCV | COPY | - |

| | | | | | | |
|--|-----------|----|--|--------|------------|---|
| Physician(s) of Record Identification Sequence | 0008,1049 | SQ | | ANAPCV | COPY | - |
| Name of Physician(s) Reading Study | 0008,1060 | PN | | ANAPCV | COPY | - |
| Physician(s) Reading Study Identification Sequence | 0008,1062 | SQ | | ANAPCV | COPY | - |
| Referenced Study Sequence | 0008,1110 | SQ | | ANAPCV | COPY | - |
| Study Instance UID | 0020,000D | UI | | ALWAYS | AUTO, COPY | - |
| Study ID | 0020,0010 | SH | | VNAP | AUTO, COPY | - |

Table 72: Patient Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------------|-----------|----|-------|-------------------|------------|---------|
| Admitting Diagnoses Description | 0008,1080 | LO | | ANAPCV | COPY | - |
| Admitting Diagnoses Code Sequence | 0008,1084 | SQ | | ANAPCV | COPY | - |
| Patient's Age | 0010,1010 | AS | | ANAPCV | COPY, USER | - |
| Patient's Size | 0010,1020 | DS | | ANAPCV | COPY, USER | - |
| Patient's Weight | 0010,1030 | DS | | ANAPCV | COPY, USER | - |
| Occupation | 0010,2180 | SH | | ANAPCV | COPY | - |
| Additional Patient History | 0010,21B0 | LT | | ANAPCV | COPY | - |

Table 73: General Series Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|-------|-------|-------------------|------------|---------|
| Series Date | 0008,0021 | DA | | ANAPCV | AUTO | - |
| Series Time | 0008,0031 | TM | | ANAPCV | AUTO | - |
| Modality | 0008,0060 | CS | | ALWAYS | COPY | - |
| Series Description | 0008,103E | LO | | ANAPCV | COPY, USER | - |
| Performing Physician's Name | 0008,1050 | PN | | ANAPCV | COPY | - |
| Operators' Name | 0008,1070 | PN | | ANAPCV | COPY, USER | - |
| Operator Identification Sequence | 0008,1072 | SQ | | ANAPCV | COPY | - |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Body Part Examined | 0018,0015 | CS | | ANAPCV | COPY | - |
| Protocol Name | 0018,1030 | LO | | ANAPCV | COPY | - |
| Patient Position | 0018,5100 | CS | | ANAPCV | COPY | - |
| Series Instance UID | 0020,000E | UI | | ALWAYS | AUTO | - |
| Series Number | 0020,0011 | IS | | VNAP | AUTO | - |
| Laterality | 0020,0060 | CS | | ANAPCV | COPY | - |
| Smallest Pixel Value in Series | 0028,0108 | US/SS | | ANAPCV | AUTO, COPY | - |
| Largest Pixel Value in Series | 0028,0109 | US/SS | | ANAPCV | AUTO, COPY | - |

| | | | | | | |
|--|-----------|----|--|--------|------|---|
| Request Attributes Sequence | 0040,0275 | SQ | | ANAPCV | COPY | - |
| Performed Procedure Step Start Date | 0040,0244 | DA | | ANAPCV | COPY | - |
| Performed Procedure Step Start Time | 0040,0245 | TM | | ANAPCV | COPY | - |
| Performed Procedure Step Description | 0040,0254 | LO | | ANAPCV | COPY | - |
| Performed Protocol Code Sequence | 0040,0260 | SQ | | ANAPCV | COPY | - |
| Comments on the Performed Procedure Step | 0040,0280 | ST | | ANAPCV | COPY | - |

Table 74: Frame of Reference Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|--------|---------|
| Frame of Reference UID | 0020,0052 | UI | | ALWAYS | AUTO | - |
| Position Reference Indicator | 0020,1040 | LO | | VNAP | AUTO | - |

Table 75: General Equipment Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------|-----------|-------|--------------------|-------------------|------------------|---------------------|
| Manufacturer | 0008,0070 | LO | Philips Healthcare | VNAP | AUTO, COPY | Philips Healthcare. |
| Institution Name | 0008,0080 | LO | | ANAPCV | AUTO, COPY, USER | - |
| Institution Address | 0008,0081 | ST | | ANAPCV | AUTO, COPY, USER | - |
| Institutional Department Name | 0008,1040 | LO | | ANAPCV | AUTO, COPY | - |
| Manufacturer's Model Name | 0008,1090 | LO | | VNAP | AUTO, COPY | - |
| Device Serial Number | 0018,1000 | LO | | ANAPCV | COPY | - |
| Software Version (s) | 0018,1020 | LO | 4.0 | ALWAYS | COPY | 4.0 |
| Station Name | 0008,1010 | SH | | ALWAYS | COPY | - |
| Spatial Resolution | 0018,1050 | DS | | ANAPCV | COPY | - |
| Date of Last Calibration | 0018,1200 | DA | | ANAPCV | COPY | - |
| Time of Last Calibration | 0018,1201 | TM | | ANAPCV | COPY | - |
| Pixel Padding Value | 0028,0120 | US/SS | | ANAPCV | COPY | - |

Table 76: General Image Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|--------|---------|
| Acquisition Date | 0008,0022 | DA | | ANAPCV | COPY | - |
| Acquisition Datetime | 0008,002A | DT | | ANAPCV | COPY | - |
| Acquisition Time | 0008,0032 | TM | | ANAPCV | COPY | - |
| Content Date | 0008,0023 | DA | | ANAPCV | COPY | - |
| Content Time | 0008,0033 | TM | | ANAPCV | COPY | - |
| Image Comments | 0020,4000 | LT | | VNAP | COPY | - |
| Referenced Image Sequence | 0008,1140 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ANAPCV | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ANAPCV | COPY | - |
| Derivation Description | 0008,2111 | ST | | ANAPCV | COPY | - |

| | | | | | | |
|--|-----------|-------|--|--------|---------------|---|
| Source Image Sequence | 0008,2112 | SQ | | ANAPCV | COPY | - |
| Derivation Code Sequence | 0008,9215 | SQ | | ANAPCV | COPY | - |
| Instance Number | 0020,0013 | IS | | ALWAYS | AUTO, COPY | - |
| Patient Orientation | 0020,0020 | CS | | ANAPCV | COPY | - |
| Images in Acquisition | 0020,1002 | IS | | ANAPCV | COPY | - |
| Quality Control Image | 0028,0300 | CS | | ANAPCV | COPY | - |
| Burned In Annotation | 0028,0301 | CS | | ANAPCV | COPY | - |
| Lossy Image Compression Ratio | 0028,2112 | DS | | ANAPCV | COPY | - |
| Icon Image Sequence | 0088,0200 | SQ | | ANAPCV | COPY | - |
| >Samples per Pixel | 0028,0002 | US | | ALWAYS | COPY | - |
| >Photometric Interpretation | 0028,0004 | CS | | ALWAYS | COPY | - |
| >Planar Configuration | 0028,0006 | US | | ANAP | COPY | - |
| >Rows | 0028,0010 | US | | ALWAYS | COPY | - |
| >Columns | 0028,0011 | US | | ALWAYS | COPY | - |
| >Pixel Aspect Ratio | 0028,0034 | IS | | ANAP | COPY | - |
| >Bits Allocated | 0028,0100 | US | | ALWAYS | COPY | - |
| >Bits Stored | 0028,0101 | US | | ALWAYS | COPY | - |
| >High Bit | 0028,0102 | US | | ALWAYS | COPY | - |
| >Pixel Representation | 0028,0103 | US | | ALWAYS | COPY | - |
| >Smallest Image Pixel Value | 0028,0106 | US/SS | | ANAPCV | COPY | - |
| >Largest Image Pixel Value | 0028,0107 | US/SS | | ANAPCV | COPY | - |
| >Red Palette Color Lookup Table Descriptor | 0028,1101 | US/SS | | ANAP | COPY | - |
| >Green Palette Color Lookup Table Descriptor | 0028,1102 | US/SS | | ANAP | COPY | - |
| >Blue Palette Color Lookup Table Descriptor | 0028,1103 | US/SS | | ANAP | COPY | - |
| >Red Palette Color Lookup Table Data | 0028,1201 | OW | | ANAP | COPY | - |
| >Green Palette Color Lookup Table Data | 0028,1202 | OW | | ANAP | COPY | - |
| >Blue Palette Color Lookup Table Data | 0028,1203 | OW | | ANAP | COPY | - |
| >Pixel Data | 7FE0,0010 | OW/OB | | ANAP | COPY | - |
| Presentation LUT Shape | 2050,0020 | CS | | ANAPCV | COPY | - |

Table 77: Image Pixel Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---|-----------|-------|-------|-------------------|---------------|---------|
| Planar Configuration | 0028,0006 | US | | ANAP | COPY | - |
| Rows | 0028,0010 | US | | ALWAYS | AUTO, COPY | - |
| Columns | 0028,0011 | US | | ALWAYS | AUTO, COPY | - |
| Pixel Aspect Ratio | 0028,0034 | IS | | ANAP | COPY | - |
| Pixel Representation | 0028,0103 | US | | ALWAYS | AUTO, COPY | - |
| Largest Image Pixel Value | 0028,0107 | US/SS | | ANAPCV | COPY | - |
| Smallest Image Pixel Value | 0028,0106 | US/SS | | ANAPCV | COPY | - |
| Red Palette Color Lookup Table Descriptor | 0028,1101 | US/SS | | ANAP | COPY | - |

| | | | | | | |
|---|-----------|-------|--|--------|------|---|
| Green Palette Color Lookup Table Descriptor | 0028,1102 | US/SS | | ANAP | COPY | - |
| Blue Palette Color Lookup Table Descriptor | 0028,1103 | US/SS | | ANAP | COPY | - |
| Red Palette Color Lookup Table Data | 0028,1201 | OW | | ANAP | COPY | - |
| Green Palette Color Lookup Table Data | 0028,1202 | OW | | ANAP | COPY | - |
| Blue Palette Color Lookup Table Data | 0028,1203 | OW | | ANAP | COPY | - |
| Pixel Data | 7FE0,0010 | OW/OB | | ALWAYS | AUTO | - |

Table 78: Multi-Frame Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------|-----------|----|-------|-------------------|------------|---------|
| Number of Frames | 0028,0008 | IS | | ALWAYS | AUTO, COPY | - |

Table 79: NM/PET Patient Orientation Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---|-----------|----|-------|-------------------|--------|---------|
| Patient Gantry Relationship Code Sequence | 0054,0414 | SQ | | ALWAYS | COPY | - |
| >Code Meaning | 0008,0104 | LO | | ALWAYS | COPY | - |
| >Code Value | 0008,0100 | SH | | ALWAYS | COPY | - |
| >Coding Scheme Designator | 0008,0102 | SH | | ALWAYS | COPY | - |
| Patient Orientation Code Sequence | 0054,0410 | SQ | | ALWAYS | COPY | - |
| >Patient Orientation Modifier Code Sequence | 0054,0412 | SQ | | ALWAYS | COPY | - |
| >>Code Meaning | 0008,0104 | LO | | ALWAYS | COPY | - |
| >>Code Value | 0008,0100 | SH | | ALWAYS | COPY | - |
| >>Code Scheme Designator | 0008,0102 | SH | | ALWAYS | COPY | - |
| >Code Meaning | 0008,0104 | LO | | ALWAYS | COPY | - |
| >Code Value | 0008,0100 | SH | | ALWAYS | COPY | - |
| >Coding Scheme Designator | 0008,0102 | SH | | ALWAYS | COPY | - |

Table 80: NM Image Pixel Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|----------------------------|-----------|----|-------|-------------------|--------|---------|
| Samples per Pixel | 0028,0002 | US | | ALWAYS | AUTO | - |
| Photometric Interpretation | 0028,0004 | CS | | ALWAYS | AUTO | - |
| Pixel Spacing | 0028,0030 | DS | | VNAP | AUTO | - |
| Bits Allocated | 0028,0100 | US | | ALWAYS | AUTO | - |
| Bits Stored | 0028,0101 | US | | ALWAYS | AUTO | - |
| High Bit | 0028,0102 | US | | ALWAYS | AUTO | - |

Table 81: NM Multi-Frame Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------|-----------|----|-------|-------------------|--------|---------|
| Detector Vector | 0054,0020 | US | | ANAP | COPY | - |
| Energy Window Vector | 0054,0010 | US | | ANAP | COPY | - |
| Frame Increment Pointer | 0028,0009 | AT | | ALWAYS | COPY | - |
| Number of Detectors | 0054,0021 | US | | ANAP | COPY | - |

| | | | | | | |
|--------------------------|-----------|----|--|--------|------|---|
| Number of Energy Windows | 0054,0011 | US | | ALWAYS | COPY | - |
| Number of R-R Intervals | 0054,0061 | US | | ANAP | COPY | - |
| Number of Rotations | 0054,0051 | US | | ANAP | COPY | - |
| Number of Slices | 0054,0081 | US | | ANAP | COPY | - |
| Number of Time Slots | 0054,0071 | US | | ANAP | COPY | - |
| R-R Interval Vector | 0054,0060 | US | | ANAP | COPY | - |
| Rotation Vector | 0054,0050 | US | | ANAP | COPY | - |
| Slice Vector | 0054,0080 | US | | ANAP | COPY | - |
| Time Slot Vector | 0054,0070 | US | | ANAP | COPY | - |

Table 82: NM Image Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------------|-----------|----|---|-------------------|--------|---------|
| Image Type | 0008,0008 | CS | Value 1: DERIVED, Value 2: SECONDARY | ALWAYS | AUTO | - |
| Acquisition Termination Condition | 0018,0071 | CS | | ANAPEV | COPY | - |
| Corrected Image | 0028,0051 | CS | | ANAPEV | COPY | - |
| Counts Accumulated | 0018,0070 | IS | | VNAP | COPY | - |
| Image ID | 0054,0400 | SH | | ANAPEV | COPY | - |
| Lossy Image Compression | 0028,2110 | CS | | ANAP | COPY | - |
| Table Height | 0018,1130 | DS | | ANAPEV | COPY | - |
| Table Traverse | 0018,1131 | DS | | ANAPEV | COPY | - |
| Actual Frame Duration | 0018,1242 | IS | | ANAP | COPY | - |
| Count Rate | 0018,1243 | IS | | ANAPEV | COPY | - |
| Processing Function | 0018,5020 | LO | | ANAPEV | COPY | - |
| Whole Body Technique | 0018,1301 | CS | | ANAPEV | COPY | - |
| Scan Velocity | 0018,1300 | DS | | ANAPCV | COPY | - |
| Scan Length | 0018,1302 | IS | | ANAPCV | COPY | - |
| Trigger Source OR Type | 0018,1061 | LO | | ANAPEV | COPY | - |
| Anatomic Region Sequence | 0008,2218 | SQ | | ANAPEV | COPY | - |

Table 83: NM Isotope Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|--------|---------|
| Energy Window Information Sequence | 0054,0012 | SQ | | VNAP | COPY | - |
| >Energy Window Name | 0054,0018 | SH | | VNAP | COPY | - |
| >Energy Window Range Sequence | 0054,0013 | SQ | | VNAP | COPY | - |
| >>Energy Window Lower Limit | 0054,0014 | DS | | VNAP | COPY | - |
| >>Energy Window Upper Limit | 0054,0015 | DS | | VNAP | COPY | - |
| Radiopharmaceutical Information Sequence | 0054,0016 | SQ | | VNAP | COPY | - |
| >Radionuclide Total Dose | 0018,1074 | DS | | ALWAYS | COPY | - |
| >Radiopharmaceutical | 0018,0031 | LO | | VNAP | COPY | - |

| | | | | | | |
|---------------------------------|-----------|----|--|--------|------|---|
| >Radiopharmaceutical Start Time | 0018,1072 | TM | | ALWAYS | COPY | - |
| >Radiopharmaceutical Volume | 0018,1071 | DS | | VNAP | COPY | - |
| >Radionuclide Code Sequence | 0054,0300 | SQ | | VNAP | COPY | - |
| >>Code Meaning | 0008,0104 | LO | | ALWAYS | COPY | - |
| >>Code Value | 0008,0100 | SH | | ALWAYS | COPY | - |
| >>Coding Scheme Designator | 0008,0102 | SH | | ALWAYS | COPY | - |
| Intervention Drug Sequence | 0018,0028 | SQ | | ANAPEV | COPY | - |

Table 84: NM Detector Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|-------|-------------------|--------|---------|
| Detector Information Sequence | 0054,0022 | SQ | | VNAP | COPY | - |
| >Image Orientation (Patient) | 0020,0037 | DS | | ALWAYS | COPY | - |
| >Zoom Center | 0028,0032 | DS | | ALWAYS | COPY | - |
| >Collimator Type | 0018,1181 | CS | | ALWAYS | COPY | - |
| >Collimator/grid Name | 0018,1180 | SH | | VNAP | COPY | - |
| >Focal Distance | 0018,1182 | IS | | VNAP | COPY | - |
| >Image Position (Patient) | 0020,0032 | DS | | ALWAYS | COPY | - |
| >View Code Sequence | 0054,0220 | SQ | | VNAP | COPY | - |
| >>Code Meaning | 0008,0104 | LO | | VNAP | COPY | - |
| >>Code Value | 0008,0100 | SH | | VNAP | COPY | - |
| >>Coding Scheme Designator | 0008,0102 | SH | | VNAP | COPY | - |

Table 85: NM Tomo Acquisition Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|-------|-------------------|--------|---------|
| Rotation Information Sequence | 0054,0052 | SQ | | VNAP | COPY | - |
| >Actual Frame Duration | 0018,1242 | IS | | VNAP | COPY | - |
| >Angular Step | 0018,1144 | DS | | VNAP | COPY | - |
| >Number of Frames in Rotation | 0054,0053 | US | | VNAP | COPY | - |
| >Rotation Direction | 0018,1140 | CS | | VNAP | COPY | - |
| >Scan Arc | 0018,1143 | DS | | VNAP | COPY | - |
| >Start Angle | 0054,0200 | DS | | VNAP | COPY | - |

| | | | | | | |
|-------------------------|-----------|----|--|--------|------|---|
| >Table Height | 0018,1130 | DS | | VNAP | COPY | - |
| >Table Traverse | 0018,1131 | DS | | VNAP | COPY | - |
| Type of Detector Motion | 0054,0202 | CS | | ANAPEV | COPY | - |

Table 86: NM Multi-gated Acquisition Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|----------------------------------|-----------|----|-------|-------------------|--------|---------|
| Gated Information Sequence | 0054,0062 | SQ | | ANAPCV | COPY | - |
| >Data Information Sequence | 0054,0063 | SQ | | VNAP | COPY | - |
| >>Frame Time | 0018,1063 | DS | | VNAP | COPY | - |
| >>Time Slot Information Sequence | 0054,0072 | SQ | | VNAP | COPY | - |
| >>>Time Slot Time | 0054,0073 | DS | | ANAPCV | COPY | - |
| Beat Rejection Flag | 0018,1080 | CS | | ANAPEV | COPY | - |
| PVC Rejection | 0018,1085 | LO | | ANAPEV | COPY | - |
| Skip Beats | 0018,1086 | IS | | ANAPEV | COPY | - |
| Heart Rate | 0018,1088 | IS | | ANAPEV | COPY | - |

Table 87: NM Reconstruction Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------|-----------|----|-------|-------------------|--------|---------|
| Slice Thickness | 0018,0050 | DS | | VNAP | AUTO | - |
| Spacing Between Slices | 0018,0088 | DS | | VNAP | COPY | - |
| Reconstruction Diameter | 0018,1100 | DS | | ANAPEV | COPY | - |
| Convolution Kernel | 0018,1210 | SH | | ANAPEV | COPY | - |
| Slice Location | 0020,1041 | DS | | ANAPEV | COPY | - |
| Slice Progression Direction | 0054,0500 | CS | | ANAPEV | COPY | - |

Table 88: Overlay Plane Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------|-----------|----|-------|-------------------|--------|---------|
| Overlay Rows | 6000,0010 | US | | ALWAYS | AUTO | - |
| Overlay Columns | 6000,0011 | US | | ALWAYS | AUTO | - |
| Overlay Description | 6000,0022 | LO | | ANAPCV | AUTO | - |
| Overlay Type | 6000,0040 | CS | | ALWAYS | AUTO | - |
| Overlay Subtype | 6000,0045 | LO | | ANAPCV | AUTO | - |
| Overlay Origin | 6000,0050 | SS | | ALWAYS | AUTO | - |
| Overlay Bits Allocated | 6000,0100 | US | | ALWAYS | AUTO | - |
| Overlay Bit Position | 6000,0102 | US | | ALWAYS | AUTO | - |
| ROI Area | 6000,1301 | IS | | ANAPCV | AUTO | - |
| ROI Mean | 6000,1302 | DS | | ANAPCV | AUTO | - |
| ROI Standard Deviation | 6000,1303 | DS | | ANAPCV | AUTO | - |
| Overlay Label | 6000,1500 | LO | | ANAPCV | AUTO | - |

| | | | | | | |
|--------------|-----------|-------|--|------|------|---|
| Overlay Data | 6000,3000 | OW/OB | | ANAP | AUTO | - |
|--------------|-----------|-------|--|------|------|---|

Table 89: VOI LUT Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------------|-----------|----|-------|-------------------|--------|---------|
| Window Center | 0028,1050 | DS | | ALWAYS | AUTO | - |
| Window Width | 0028,1051 | DS | | ANAP | AUTO | - |
| Window Center & Width Explanation | 0028,1055 | LO | | ANAPCV | AUTO | - |
| VOI LUT Sequence | 0028,3010 | SQ | | ANAPCV | AUTO | - |

Table 90: SOP Common Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------------|-----------|----|----------------------------|-------------------|------------|----------------------|
| Specific Character Set | 0008,0005 | CS | | ANAP | AUTO, COPY | - |
| Instance Creation Date | 0008,0012 | DA | | VNAP | AUTO, COPY | - |
| Instance Creation Time | 0008,0013 | TM | | VNAP | AUTO, COPY | - |
| SOP Class UID | 0008,0016 | UI | 1.2.840.10008.5.1.4.1.1.20 | ALWAYS | COPY, AUTO | - |
| SOP Instance UID | 0008,0018 | UI | | ALWAYS | AUTO | - |
| Timezone Offset From UTC | 0008,0201 | SH | | ANAPCV | COPY | - |
| Contributing Equipment Sequence | 0018,A001 | SQ | | ANAP | AUTO | - |
| >Manufacturer | 0008,0070 | LO | Philips | ANAP | FIXED | PHILIPS. |
| >Institution Address | 0008,0081 | ST | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Station Name | 0008,1010 | SH | IntelliSpace Portal | ANAP | CONFIG | IntelliSpace Portal. |
| >Institutional Department Name | 0008,1040 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Software Version(s) | 0018,1020 | LO | 4.0 | ANAP | FIXED | 4.0 |
| >Device Serial Number | 0018,1000 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Institution Name | 0008,0080 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Manufacturer's Model Name | 0008,1090 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Purpose of Reference Code Sequence | 0040,A170 | SQ | | ANAP | FIXED | - |
| >>Code Value | 0008,0100 | SH | 109102 | ANAP | FIXED | 109102 |
| >>Coding Scheme Designator | 0008,0102 | SH | DCM | ANAP | FIXED | DCM |
| >>Code Meaning | 0008,0104 | LO | Processing Equipment | ANAP | FIXED | Processing Equipment |

Table 91:NM Phase Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|--------|---------|
| Phase Information Sequence | 0054,0032 | | | ANAPCV | COPY | - |
| >Phase Delay | 0054,0036 | | | ALWAYS | COPY | - |
| >Actual Frame Duration | 0018,1242 | | | ALWAYS | COPY | - |
| >Pause Between Frames | 0054,0038 | | | ALWAYS | COPY | - |
| >Number of Frames in Phase | 0054,0033 | | | ALWAYS | COPY | - |
| >Trigger Vector | 0054,0210 | | | ANAPEV | COPY | - |
| >Number of Triggers in Phase | 0054,0211 | | | ANAP | COPY | - |
| >Phase Description | 0054,0039 | | | ANAPEV | COPY | - |

Table 92: Extended DICOM and private attributes for Nuclear Medicine Image Storage SOP Class Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------|-----------|----|-------|-------------------|--------|---------|
| Rescale Intercept | 0028,1052 | DS | | ALWAYS | AUTO | - |
| Rescale Slope | 0028,1053 | DS | | ALWAYS | AUTO | - |

| | | | | | | |
|-------------------------------------|-----------|-------|------------|--------|------|---|
| Series Type (PIIM) | 2001,106E | SH | | ALWAYS | AUTO | - |
| Frame Number | 2001,107C | IS | Value 1: 1 | ALWAYS | AUTO | - |
| Frame Pixel Data | 2001,107D | OW/OB | | ALWAYS | AUTO | - |
| Frame Index Pointer | 2001,10D0 | DS | | ALWAYS | AUTO | - |
| Index Number Vector | 2001,10D1 | DS | | ALWAYS | AUTO | - |
| Frame Slice Number | 2001,10D7 | DS | 1 | ALWAYS | AUTO | - |
| Curve Dimensions | 50xx,0005 | US | | ALWAYS | AUTO | - |
| Number of Points | 50xx,0010 | US | | ALWAYS | AUTO | - |
| Type of Data | 50xx,0020 | CS | | ALWAYS | AUTO | - |
| Data Value Representation | 50xx,0103 | US | | ALWAYS | AUTO | - |
| Curve Data | 50xx,3000 | OB/OW | | ALWAYS | AUTO | - |
| Curve Description | 50xx,0022 | LO | | ANAPCV | AUTO | - |
| Axis Units | 50xx,0030 | SH | | ANAPCV | AUTO | - |
| Axis Labels | 50xx,0040 | SH | | ANAPCV | AUTO | - |
| Minimum Coordinate Value | 50xx,0104 | US | | ANAPCV | AUTO | - |
| Maximum Coordinate Value | 50xx,0105 | US | | ANAPCV | AUTO | - |
| Curve Range | 50xx,0106 | SH | | ANAPCV | AUTO | - |
| Curve Data Descriptor | 50xx,0110 | US | | ANAP | AUTO | - |
| Coordinate Start Value | 50xx,0112 | US | | ANAP | AUTO | - |
| Coordinate Step Value | 50xx,0114 | US | | ANAP | AUTO | - |
| Curve Label | 50xx,2500 | LO | | ANAPCV | AUTO | - |
| Curve Referenced Overlay Sequence | 50xx,2600 | SQ | | ANAPCV | AUTO | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ANAP | AUTO | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ANAP | AUTO | - |
| >Referenced Overlay Group | 50xx,2610 | US | | ANAP | AUTO | - |
| Referenced Overlay Sequence | 0008,1130 | SQ | | ANAPEV | COPY | - |
| Reference Curve Sequence | 0008,1145 | SQ | | ANAPEV | COPY | - |
| Primary Anatomic Structure Sequence | 0008,2228 | SQ | | ANAPEV | COPY | - |

2.1.2.6. MR Image Storage SOP Class

Table 93: IOD of Created MR Image Storage SOP Class Instances

| Information Entity | Module | Presence Of Module |
|--------------------|---------------------------------------|--------------------|
| Patient | Patient Module | ALWAYS |
| Study | General Study Module | ALWAYS |
| Study | Patient Study Module | CONDITIONAL |
| Series | General Series Module | ALWAYS |
| Frame of Reference | Frame of Reference Module | ALWAYS |
| Equipment | General Equipment Module | ALWAYS |
| Image | General Image Module | ALWAYS |
| Image | Image Plane Module | ALWAYS |
| Image | Image Pixel Module | ALWAYS |
| Image | MR Image Module | ALWAYS |
| Image | Overlay Plane Module | CONDITIONAL |
| Image | VOI LUT Module | CONDITIONAL |
| Image | SOP Common Module | ALWAYS |
| Image | Contrast/Bolus Module | CONDITIONAL |
| | Extended Dicom and Private attributes | CONDITIONAL |

Table 94: Patient Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|---------------|---------|
| Referenced Patient Sequence | 0008,1120 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Patient's Name | 0010,0010 | PN | | VNAP | COPY | - |
| Patient ID | 0010,0020 | LO | | VNAP | COPY | - |
| Patient's Birth Date | 0010,0030 | DA | | VNAP | COPY, USER | - |
| Patient's Birth Time | 0010,0032 | TM | | ANAPCV | COPY, USER | - |
| Patient's Sex | 0010,0040 | CS | | VNAP | COPY, USER | - |
| Other Patient IDs | 0010,1000 | LO | | ANAPCV | COPY, USER | - |
| Other Patient Names | 0010,1001 | PN | | ANAPCV | COPY, USER | - |
| Ethnic Group | 0010,2160 | SH | | ANAPCV | COPY | - |
| Patient Comments | 0010,4000 | LT | | ANAPCV | COPY | - |
| Issuer of Patient ID | 0010,0021 | LO | | ANAPCV | COPY, USER | - |

Table 95: General Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|---------------|---------|
| Study Date | 0008,0020 | DA | | VNAP | COPY | - |
| Study Time | 0008,0030 | TM | | VNAP | COPY | - |
| Accession Number | 0008,0050 | SH | | VNAP | COPY, USER | - |
| Referring Physician's Name | 0008,0090 | PN | | VNAP | COPY, USER | - |
| Referring Physician Identification Sequence | 0008,0096 | SQ | | ANAPCV | COPY | - |
| Study Description | 0008,1030 | LO | | ANAPCV | COPY, USER | - |
| Procedure Code Sequence | 0008,1032 | SQ | | ANAPCV | COPY | - |
| Physician(s) of Record | 0008,1048 | PN | | ANAPCV | COPY | - |
| Physician(s) of Record Identification Sequence | 0008,1049 | SQ | | ANAPCV | COPY | - |
| Name of Physician(s) Reading Study | 0008,1060 | PN | | ANAPCV | COPY | - |
| Physician(s) Reading Study Identification Sequence | 0008,1062 | SQ | | ANAPCV | COPY | - |
| Referenced Study Sequence | 0008,1110 | SQ | | ANAPCV | COPY | - |
| Study Instance UID | 0020,000D | UI | | ALWAYS | AUTO, COPY | - |
| Study ID | 0020,0010 | SH | | VNAP | AUTO, COPY | - |

Table 96: Patient Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------------|-----------|----|-------|-------------------|------------|---------|
| Admitting Diagnoses Description | 0008,1080 | LO | | ANAPCV | COPY | - |
| Admitting Diagnoses Code Sequence | 0008,1084 | SQ | | ANAPCV | COPY | - |
| Patient's Age | 0010,1010 | AS | | ANAPCV | COPY, USER | - |
| Patient's Size | 0010,1020 | DS | | ANAPCV | COPY, USER | - |
| Patient's Weight | 0010,1030 | DS | | ANAPCV | COPY, USER | - |
| Occupation | 0010,2180 | SH | | ANAPCV | COPY | - |
| Additional Patient History | 0010,21B0 | LT | | ANAPCV | COPY | - |

Table 97: General Series Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|-------|-------|-------------------|------------|--------------------------------------|
| Series Date | 0008,0021 | DA | | ANAPCV | AUTO | - |
| Series Time | 0008,0031 | TM | | ANAPCV | AUTO | - |
| Modality | 0008,0060 | CS | | ALWAYS | COPY | - |
| Series Description | 0008,103E | LO | | ANAPCV | COPY, USER | - |
| Performing Physician's Name | 0008,1050 | PN | | ANAPCV | COPY | - |
| Operators' Name | 0008,1070 | PN | | ANAPCV | COPY, USER | - |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Body Part Examined | 0018,0015 | CS | | ANAPCV | COPY | - |
| Protocol Name | 0018,1030 | LO | | ANAPCV | COPY | - |
| Patient Position | 0018,5100 | CS | | ANAPCV | COPY | - |
| Series Instance UID | 0020,000E | UI | | ALWAYS | AUTO | Generated by IntelliSpace Portal 4.0 |
| Series Number | 0020,0011 | IS | | ALWAYS | AUTO | - |
| Laterality | 0020,0060 | CS | | ANAPCV | COPY | - |
| Smallest Pixel Value in Series | 0028,0108 | US/SS | | ANAPCV | COPY | - |
| Largest Pixel Value in Series | 0028,0109 | US/SS | | ANAPCV | COPY | - |
| Request Attributes Sequence | 0040,0275 | SQ | | VNAP | COPY | - |
| >Scheduled Procedure Step Description | 0040,0007 | LO | | VNAP | COPY | - |
| >Scheduled Protocol Code Sequence | 0040,0008 | SQ | | ANAPCV | COPY | - |
| >Scheduled Procedure Step ID | 0040,0009 | SH | | ALWAYS | COPY | - |
| >Requested Procedure ID | 0040,1001 | SH | | ALWAYS | COPY | - |
| Performed Procedure Step Start Date | 0040,0244 | DA | | ALWAYS | COPY | - |
| Performed Procedure Step Start Time | 0040,0245 | TM | | ALWAYS | COPY | - |
| Performed Procedure Step ID | 0040,0253 | SH | | ALWAYS | COPY | - |
| Performed Procedure Step Description | 0040,0254 | LO | | ANAPCV | COPY | - |

| | | | | | | |
|--|-----------|----|--|--------|------|---|
| Performed Protocol Code Sequence | 0040,0260 | SQ | | ANAPCV | COPY | - |
| >Code Value | 0008,0100 | SH | | ALWAYS | COPY | - |
| >Coding Scheme Designator | 0008,0102 | SH | | ALWAYS | COPY | - |
| >Coding Scheme Version | 0008,0103 | SH | | ANAP | COPY | - |
| >Mapping Resource | 0008,0105 | CS | | ANAP | COPY | - |
| >Context Group Version | 0008,0106 | DT | | ANAP | COPY | - |
| >Context Group Local Version | 0008,0107 | DT | | ANAP | COPY | - |
| >Context Group Extension Flag | 0008,010B | CS | | ALWAYS | COPY | - |
| >Context Group Extension Creator UID | 0008,010D | UI | | ANAP | COPY | - |
| >Context Identifier | 0008,010F | CS | | VNAP | COPY | - |
| Comments on the Performed Procedure Step | 0040,0280 | ST | | ANAPCV | COPY | - |

Table 98: Frame of Reference Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|-------------|---------|
| Frame of Reference UID | 0020,0052 | UI | | ALWAYS | COPY | - |
| Position Reference Indicator | 0020,1040 | LO | | EMPTY | COPY, FIXED | - |

Table 99: General Equipment Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------|-----------|-------|-------|-------------------|------------|-----------------------|
| Manufacturer | 0008,0070 | LO | | ALWAYS | COPY | - |
| Institution Name | 0008,0080 | LO | | ANAPCV | COPY, USER | - |
| Institution Address | 0008,0081 | ST | | ANAPCV | COPY, USER | - |
| Station Name | 0008,1010 | SH | | ANAPCV | COPY | - |
| Institutional Department Name | 0008,1040 | LO | | ANAPCV | COPY | - |
| Manufacturer's Model Name | 0008,1090 | LO | | ANAPCV | COPY | - |
| Device Serial Number | 0018,1000 | LO | | ANAPCV | FIXED | System serial number. |
| Software Version(s) | 0018,1020 | LO | 4.0 | ANAPCV | COPY | 4.0 |
| Spatial Resolution | 0018,1050 | DS | | ANAPCV | COPY | - |
| Date of Last Calibration | 0018,1200 | DA | | ANAPCV | COPY | - |
| Time of Last Calibration | 0018,1201 | TM | | ANAPCV | COPY | - |
| Pixel Padding Value | 0028,0120 | US/SS | | ANAPCV | COPY | - |

Table 100: General Image Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|----------|---------|
| Acquisition Date | 0008,0022 | DA | | ANAPCV | COPY | - |
| Content Date | 0008,0023 | DA | | ANAPCV | COPY | - |
| Acquisition Datetime | 0008,002A | DT | | ANAPCV | COPY | - |
| Acquisition Time | 0008,0032 | TM | | ANAPCV | COPY | - |
| Content Time | 0008,0033 | TM | | ANAPCV | COPY | - |
| Referenced Image Sequence | 0008,1140 | SQ | | ANAPCV | AUTO | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ANAPCV | IMPLICIT | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ANAPCV | IMPLICIT | - |
| Acquisition Number | 0020,0012 | IS | | ALWAYS | COPY | - |
| Patient Orientation | 0020,0020 | CS | | ANAPCV | COPY | - |

| | | | | | | |
|--|-----------|-------|--|--------|----------------|-------------------------|
| Derivation Description | 0008,2111 | ST | | ANAPCV | COPY | - |
| Source Image Sequence | 0008,2112 | SQ | | ANAPCV | COPY | - |
| Derivation Code Sequence | 0008,9215 | SQ | | ANAPCV | COPY | - |
| Instance Number | 0020,0013 | IS | | ALWAYS | AUTO, COPY | - |
| Patient Orientation | 0020,0020 | CS | | ANAPCV | COPY | - |
| Images in Acquisition | 0020,1002 | IS | | ANAPCV | COPY | - |
| Quality Control Image | 0028,0300 | CS | | ANAPCV | COPY | - |
| Burned In Annotation | 0028,0301 | CS | | ANAPCV | COPY | - |
| Lossy Image Compression | 0028,2110 | CS | | ANAPCV | COPY, FIXED | Applied value: 00 |
| Lossy Image Compression Ratio | 0028,2112 | DS | | ANAPCV | COPY | - |
| Icon Image Sequence | 0088,0200 | SQ | | ANAPCV | COPY | - |
| >Samples per Pixel | 0028,0002 | US | | ALWAYS | COPY | - |
| >Photometric Interpretation | 0028,0004 | CS | | ALWAYS | COPY | - |
| >Planar Configuration | 0028,0006 | US | | ANAP | COPY | - |
| >Rows | 0028,0010 | US | | ALWAYS | COPY | - |
| >Columns | 0028,0011 | US | | ALWAYS | COPY | - |
| >Pixel Aspect Ratio | 0028,0034 | IS | | ANAP | COPY | - |
| >Bits Allocated | 0028,0100 | US | | ALWAYS | COPY | - |
| >Bits Stored | 0028,0101 | US | | ALWAYS | COPY | - |
| >High Bit | 0028,0102 | US | | ALWAYS | COPY | - |
| >Pixel Representation | 0028,0103 | US | | ALWAYS | COPY | - |
| >Smallest Image Pixel Value | 0028,0106 | US/SS | | ANAPCV | COPY | - |
| >Largest Image Pixel Value | 0028,0107 | US/SS | | ANAPCV | COPY | - |
| >Red Palette Color Lookup Table Descriptor | 0028,1101 | US/SS | | ANAP | COPY | - |
| >Green Palette Color Lookup Table Descriptor | 0028,1102 | US/SS | | ANAP | COPY | - |
| >Blue Palette Color Lookup Table Descriptor | 0028,1103 | US/SS | | ANAP | COPY | - |
| >Red Palette Color Lookup Table Data | 0028,1201 | OW | | ANAP | COPY | - |
| >Green Palette Color Lookup Table Data | 0028,1202 | OW | | ANAP | COPY | - |
| >Blue Palette Color Lookup Table Data | 0028,1203 | OW | | ANAP | COPY | - |
| >Pixel Data | 7FE0,0010 | OW/OB | | ANAP | COPY | - |
| Presentation LUT Shape | 2050,0020 | CS | | ANAPCV | COPY, FIXED | Applied value: IDENTITY |

Table 101: Image Plane Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------|-----------|----|-------|-------------------|---------------|---|
| Slice Thickness | 0018,0050 | DS | | VNAP | AUTO, COPY | - |
| Image Position (Patient) | 0020,0032 | DS | | ALWAYS | AUTO, COPY | - |
| Image Orientation (Patient) | 0020,0037 | DS | | ALWAYS | AUTO, COPY | - |
| Slice Location | 0020,1041 | DS | | ANAPCV | AUTO, COPY | Value is shortest distance between TLHC of the images from TLHC of first image of the series. |

| | | | | | | |
|---------------|-----------|----|--|--------|---------------|---|
| Pixel Spacing | 0028,0030 | DS | | ALWAYS | AUTO, COPY | - |
|---------------|-----------|----|--|--------|---------------|---|

Table 102: Image Pixel Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---|-----------|-------|-------|-------------------|----------------|-------------------------------------|
| Planar Configuration | 0028,0006 | US | | ANAP | AUTO, COPY | - |
| Rows | 0028,0010 | US | | ALWAYS | IMPLICIT | Applied values: min: 64 - max: 2048 |
| Columns | 0028,0011 | US | | ALWAYS | IMPLICIT | Applied values: min: 64 - max: 2048 |
| Pixel Aspect Ratio | 0028,0034 | IS | | ANAP | COPY, FIXED | Applied value: (1\1) |
| Bits Stored | 0028,0101 | US | | ALWAYS | COPY | - |
| High Bit | 0028,0102 | US | | ALWAYS | COPY | - |
| Pixel Representation | 0028,0103 | US | | ALWAYS | COPY | - |
| Smallest Image Pixel Value | 0028,0106 | US/SS | | ANAPCV | COPY | - |
| Largest Image Pixel Value | 0028,0107 | US/SS | | ANAPCV | COPY | - |
| Red Palette Color Lookup Table Descriptor | 0028,1101 | US/SS | | ANAP | COPY | - |
| Green Palette Color Lookup Table Descriptor | 0028,1102 | US/SS | | ANAP | COPY | - |
| Blue Palette Color Lookup Table Descriptor | 0028,1103 | US/SS | | ANAP | COPY | - |
| Red Palette Color Lookup Table Data | 0028,1201 | OW | | ANAP | COPY | - |
| Green Palette Color Lookup Table Data | 0028,1202 | OW | | ANAP | COPY | - |
| Blue Palette Color Lookup Table Data | 0028,1203 | OW | | ANAP | COPY | - |
| Pixel Data | 7FE0,0010 | OW/OB | | ALWAYS | IMPLICIT | - |

Table 103: MR Image Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---------------------|-----------|----|-------|-------------------|---------------|--|
| Image Type | 0008,0008 | CS | | ALWAYS | AUTO | Applied values: {{ORIGINAL, DERIVED}}, PRIMARY, {T2_STAR, R2, R2_STAR, W, F, IP, OP, KTRANS, KEP, VE, VP}, {KTRANS, KEP, VE, VP}, {UNSPECIFIED}} |
| Scanning Sequence | 0018,0020 | CS | | ALWAYS | AUTO, COPY | - |
| Sequence Variant | 0018,0021 | CS | | ALWAYS | AUTO, COPY | - |
| Scan Options | 0018,0022 | CS | | VNAP | IMPLICIT | - |
| MR Acquisition Type | 0018,0023 | CS | | ANAP | AUTO, COPY | - |
| Sequence Name | 0018,0024 | SH | | ANAPCV | AUTO, COPY | - |
| Angio Flag | 0018,0025 | CS | | ANAPCV | AUTO, COPY | - |
| Repetition Time | 0018,0080 | DS | | ANAPCV | COPY | - |
| Echo Time | 0018,0081 | DS | | ANAP | COPY | - |
| Inversion Time | 0018,0082 | DS | | ANAPCV | COPY | - |
| Number of Averages | 0018,0083 | DS | | ANAPCV | COPY | - |
| Imaging Frequency | 0018,0084 | DS | | ANAPCV | COPY | - |

| | | | | | | |
|-----------------------------------|-----------|----|--|--------|------|---|
| Imaged Nucleus | 0018,0085 | SH | | ANAPCV | COPY | - |
| Echo Number(s) | 0018,0086 | IS | | ANAPCV | COPY | - |
| Magnetic Field Strength | 0018,0087 | DS | | ANAPCV | COPY | - |
| Spacing Between Slices | 0018,0088 | DS | | ANAPCV | COPY | - |
| Number of Phase Encoding Steps | 0018,0089 | IS | | ANAPCV | COPY | - |
| Echo Train Length | 0018,0091 | IS | | VNAP | COPY | - |
| Percent Sampling | 0018,0093 | DS | | ANAPCV | COPY | - |
| Percent Phase Field of View | 0018,0094 | DS | | ANAPCV | COPY | - |
| Pixel Bandwidth | 0018,0095 | DS | | ANAPCV | COPY | - |
| Trigger Time | 0018,1060 | DS | | ANAPCV | COPY | Will only have a value if Dynamic Series (2001,1012) Equals 1 |
| Nominal Interval | 0018,1062 | IS | | ANAPCV | COPY | - |
| Beat Rejection Flag | 0018,1080 | CS | | ANAPCV | COPY | - |
| Low R-R Value | 0018,1081 | IS | | ANAPCV | COPY | - |
| High R-R Value | 0018,1082 | IS | | ANAPCV | COPY | - |
| Intervals Acquired | 0018,1083 | IS | | ANAPCV | COPY | - |
| Intervals Rejected | 0018,1084 | IS | | ANAPCV | COPY | - |
| PVC Rejection | 0018,1085 | LO | | ANAPCV | COPY | - |
| Skip Beats | 0018,1086 | IS | | ANAPCV | COPY | - |
| Heart Rate | 0018,1088 | IS | | ANAPCV | COPY | - |
| Cardiac Number of Images | 0018,1090 | IS | | ANAPCV | COPY | - |
| Trigger Window | 0018,1094 | IS | | ANAPCV | COPY | - |
| Reconstruction Diameter | 0018,1100 | DS | | ANAPCV | COPY | Value is a copy of the largest value of the Field of View. |
| Receive Coil Name | 0018,1250 | SH | | ANAPCV | COPY | - |
| Transmit Coil Name | 0018,1251 | SH | | ANAPCV | COPY | - |
| Acquisition Matrix | 0018,1310 | US | | ANAPCV | COPY | - |
| In-plane Phase Encoding Direction | 0018,1312 | CS | | ANAPCV | COPY | - |
| Flip Angle | 0018,1314 | DS | | ANAPCV | COPY | - |
| Variable Flip Angle Flag | 0018,1315 | CS | | ANAPCV | COPY | - |
| SAR | 0018,1316 | DS | | ANAPCV | COPY | - |
| dB/dt | 0018,1318 | DS | | ANAPCV | COPY | - |
| Temporal Position Identifier | 0020,0100 | IS | | ANAPCV | COPY | - |
| Number of Temporal Positions | 0020,0105 | IS | | ANAPCV | COPY | - |
| Temporal Resolution | 0020,0110 | DS | | ALWAYS | COPY | - |
| Samples per Pixel | 0028,0002 | US | | ALWAYS | COPY | - |
| Photometric Interpretation | 0028,0004 | CS | | ALWAYS | COPY | - |
| Bits Allocated | 0028,0100 | US | | ANAPCV | COPY | - |

Table 104: Overlay Plane Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------|-----------|----|-------|-------------------|--------|---------|
| Overlay Rows | 6000,0010 | US | | ALWAYS | AUTO | - |
| Overlay Columns | 6000,0011 | US | | ALWAYS | AUTO | - |
| Overlay Description | 6000,0022 | LO | | ANAPCV | AUTO | - |
| Overlay Type | 6000,0040 | CS | | ALWAYS | AUTO | - |
| Overlay Subtype | 6000,0045 | LO | | ANAPCV | AUTO | - |
| Overlay Origin | 6000,0050 | SS | | ALWAYS | AUTO | - |
| Overlay Bits Allocated | 6000,0100 | US | | ALWAYS | AUTO | - |
| Overlay Bit Position | 6000,0102 | US | | ALWAYS | AUTO | - |

| | | | | | | |
|------------------------|-----------|-------|--|--------|------|---|
| ROI Area | 6000,1301 | IS | | ANAPCV | AUTO | - |
| ROI Mean | 6000,1302 | DS | | ANAPCV | AUTO | - |
| ROI Standard Deviation | 6000,1303 | DS | | ANAPCV | AUTO | - |
| Overlay Label | 6000,1500 | LO | | ANAPCV | AUTO | - |
| Overlay Data | 6000,3000 | OW/OB | | ANAP | AUTO | - |

Table 105: VOI LUT Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------------|-----------|-------|--|-------------------|--------|---------|
| Window Center | 0028,1050 | DS | | ALWAYS | AUTO | - |
| Window Width | 0028,1051 | DS | | ANAP | AUTO | - |
| Window Center & Width Explanation | 0028,1055 | LO | | ANAPCV | AUTO | - |
| VOI LUT Sequence | 0028,3010 | SQ | | ANAPCV | AUTO | - |
| >LUT Descriptor | 0028,3002 | US/SS | Value 1: 0, Value 2: 0, Value 3: 0 | ALWAYS | AUTO | - |
| >LUT Explanation | 0028,3003 | LO | | EMPTY | AUTO | - |
| >LUT Data | 0028,3006 | US/OW | Value 1: 0 | ALWAYS | AUTO | - |

Table 106: SOP Common Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------------|-----------|----|----------------------|-------------------|-------------|----------------------|
| Specific Character Set | 0008,0005 | CS | | ANAP | AUTO, COPY | - |
| Instance Creation Date | 0008,0012 | DA | | ALWAYS | AUTO | - |
| Instance Creation Time | 0008,0013 | TM | | ALWAYS | AUTO | - |
| Instance Creator UID | 0008,0014 | UI | | ANAPCV | AUTO | - |
| SOP Class UID | 0008,0016 | UI | | ALWAYS | COPY, FIXED | - |
| SOP Instance UID | 0008,0018 | UI | | ALWAYS | AUTO | - |
| Timezone Offset From UTC | 0008,0201 | SH | | ANAPCV | COPY | - |
| Contributing Equipment Sequence | 0018,A001 | SQ | | ANAP | AUTO | - |
| >Manufacturer | 0008,0070 | LO | Philips | ANAP | FIXED | PHILIPS. |
| >Institution Address | 0008,0081 | ST | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Station Name | 0008,1010 | SH | IntelliSpace Portal | ANAP | CONFIG | IntelliSpace Portal. |
| >Institutional Department Name | 0008,1040 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Software Version(s) | 0018,1020 | LO | 4.0 | ANAP | FIXED | 4.0 |
| >Device Serial Number | 0018,1000 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Institution Name | 0008,0080 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Manufacturer's Model Name | 0008,1090 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Purpose of Reference Code Sequence | 0040,A170 | SQ | | ANAP | FIXED | - |
| >>Code Value | 0008,0100 | SH | 109102 | ANAP | FIXED | 109102 |
| >>Coding Scheme Designator | 0008,0102 | SH | DCM | ANAP | FIXED | DCM |
| >>Code Meaning | 0008,0104 | LO | Processing Equipment | ANAP | FIXED | Processing Equipment |

Table 107: Contrast/Bolus Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------|-----------|----|-------|-------------------|--------|---------|
| Contrast/Bolus Agent | 0018,0010 | LO | | VNAP | COPY | - |
| Contrast/Bolus Agent Sequence | 0018,0012 | SQ | | ANAPCV | COPY | - |

| | | | | | | |
|--|-----------|----|--|--------|------|---|
| Contrast/Bolus Administration Route Sequence | 0018,0014 | SQ | | ANAPCV | COPY | - |
| >Additional Drug Sequence | 0018,002A | SQ | | ANAP | COPY | - |
| Contrast/Bolus Route | 0018,1040 | LO | | ANAPCV | COPY | - |
| Contrast/Bolus Volume | 0018,1041 | DS | | ANAPCV | COPY | - |
| Contrast/Bolus Start Time | 0018,1042 | TM | | ANAPCV | COPY | - |
| Contrast/Bolus Stop Time | 0018,1043 | TM | | ANAPCV | COPY | - |
| Contrast/Bolus Total Dose | 0018,1044 | DS | | ANAPCV | COPY | - |
| Contrast Flow Rate | 0018,1046 | DS | | ANAPCV | COPY | - |
| Contrast Flow Duration | 0018,1047 | DS | | ANAPCV | COPY | - |
| Contrast/Bolus Ingredient | 0018,1048 | CS | | ANAPCV | COPY | - |
| Contrast/Bolus Ingredient Concentration | 0018,1049 | DS | | ANAPCV | COPY | - |

Table 108: Extended DICOM and private attributes for MR Image Storage SOP Class Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---------------------------------------|-----------|----|-------|-------------------|--------|---------------|
| Code Value | 0008,0100 | SH | | VNAP | AUTO | - |
| Coding Scheme Designator | 0008,0102 | SH | | ALWAYS | AUTO | - |
| Coding Scheme Version | 0008,0103 | SH | | ALWAYS | AUTO | - |
| Code Meaning | 0008,0104 | LO | | VNAP | AUTO | - |
| Derivation Image Sequence | 0008,9124 | SQ | | ANAP | AUTO | - |
| >Source Image Sequence | 0008,2112 | SQ | | ALWAYS | AUTO | - |
| >>Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | AUTO | - |
| >>Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | AUTO | - |
| >Derivation Code Sequence | 0008,9215 | SQ | | ALWAYS | AUTO | - |
| >>Code Value | 0008,0100 | SH | | ALWAYS | AUTO | From CID 7203 |
| >>Coding Scheme Designator | 0008,0102 | SH | | ALWAYS | AUTO | From CID 7203 |
| >>Coding Scheme Version | 0008,0103 | SH | | ANAP | AUTO | From CID 7203 |
| >>Code Meaning | 0008,0104 | LO | | ALWAYS | AUTO | From CID 7203 |
| Medical Alerts | 0010,2000 | LO | | ANAPCV | COPY | - |
| Allergies | 0010,2110 | LO | | ANAPCV | COPY | - |
| Pregnancy Status | 0010,21C0 | US | | VNAP | COPY | - |
| Acquisition Duration | 0018,9073 | FD | | ANAP | COPY | - |
| Diffusion b-value | 0018,9087 | FD | | ANAP | AUTO | - |
| Diffusion Gradient Orientation | 0018,9089 | FD | | ANAP | AUTO | - |
| Rescale Intercept | 0028,1052 | DS | | ALWAYS | COPY | - |
| Rescale Slope | 0028,1053 | DS | | ALWAYS | COPY | - |
| Rescale Type | 0028,1054 | LO | | ALWAYS | COPY | - |
| Requesting Physician | 0032,1032 | PN | | VNAP | COPY | - |
| Requesting Service | 0032,1033 | LO | | VNAP | COPY | - |
| Requested Procedure Description | 0032,1060 | LO | | VNAP | COPY | - |
| Requested Contrast Agent | 0032,1070 | LO | | VNAP | COPY | - |
| Study Comments (retired) | 0032,4000 | LT | | VNAP | COPY | - |
| Special Needs | 0038,0050 | LO | | ANAPCV | COPY | - |
| Patient State | 0038,0500 | LO | | ANAPCV | COPY | - |
| Scheduled Performing Physician's Name | 0040,0006 | PN | | VNAP | COPY | - |
| Performed Station AE Title | 0040,0241 | AE | | ANAPCV | COPY | - |
| Performed Station Name | 0040,0242 | SH | | ANAPCV | COPY | - |

| | | | | | | |
|--|-----------|----|--|--------|----------|---|
| Performed Location | 0040,0243 | SH | | ANAPCV | COPY | - |
| Performed Procedure Step End Date | 0040,0250 | DA | | VNAP | COPY | - |
| Performed Procedure Step End Time | 0040,0251 | TM | | VNAP | COPY | - |
| Performed Procedure Step Status | 0040,0252 | CS | | ANAP | COPY | - |
| Film Consumption Sequence | 0040,0321 | SQ | | EMPTY | AUTO | - |
| Requested Procedure ID | 0040,1001 | SH | | VNAP | COPY | - |
| Reason for the Requested Procedure | 0040,1002 | LO | | VNAP | COPY | - |
| Requested Procedure Priority | 0040,1003 | SH | | VNAP | COPY | - |
| Patient Transport Arrangements | 0040,1004 | LO | | VNAP | COPY | - |
| Requested Procedure Location | 0040,1005 | LO | | VNAP | COPY | - |
| Requested Procedure Comments | 0040,1400 | LT | | VNAP | COPY | - |
| Reason for the Imaging Service Request (retired) | 0040,2001 | LO | | VNAP | COPY | - |
| Issue Date of Imaging Service Request | 0040,2004 | DA | | VNAP | COPY | - |
| Issue Time of Imaging Service Request | 0040,2005 | TM | | VNAP | COPY | - |
| Order Enterer's Location | 0040,2009 | SH | | VNAP | COPY | - |
| Order Callback Phone Number | 0040,2010 | SH | | VNAP | COPY | - |
| Imaging Service Request Comments | 0040,2400 | LT | | VNAP | COPY | - |
| Real World Value Mapping Sequence | 0040,9096 | SQ | | ALWAYS | AUTO | - |
| >Real World Value Intercept | 0040,9224 | FD | | ALWAYS | AUTO | - |
| >Real World Value Slope | 0040,9225 | FD | | ALWAYS | AUTO | - |
| Private Creator Group 2001 | 2001,0010 | LO | | ALWAYS | FIXED | - |
| Chemical Shift | 2001,1001 | FL | | ANAPCV | COPY | - |
| Chemical Shift Number MR | 2001,1002 | IS | | ANAPCV | COPY | - |
| Diffusion B-Factor | 2001,1003 | FL | | ANAPCV | COPY | - |
| Diffusion Direction | 2001,1004 | CS | | ANAPCV | COPY | - |
| Image Enhanced | 2001,1006 | CS | | VNAP | COPY | - |
| Image Type ED ES | 2001,1007 | CS | | VNAP | COPY | - |
| Phase Number | 2001,1008 | IS | | VNAP | COPY | - |
| Image Prepulse Delay | 2001,1009 | FL | | ALWAYS | COPY | - |
| Slice Number MR | 2001,100A | IS | | VNAP | IMPLICIT | - |
| Slice Orientation | 2001,100B | CS | | ALWAYS | IMPLICIT | - |
| Arrhythmia Rejection | 2001,100C | CS | | ALWAYS | COPY | - |
| Cardiac Cycled | 2001,100E | CS | | ALWAYS | COPY | - |
| Cardiac Gate Width | 2001,100F | SS | | ALWAYS | COPY | - |
| Cardiac Sync | 2001,1010 | CS | | ALWAYS | COPY | - |
| Diffusion Echo Time | 2001,1011 | FL | | ANAPCV | IMPLICIT | - |
| Dynamic Series | 2001,1012 | CS | | VNAP | COPY | - |
| EPI Factor | 2001,1013 | SL | | ALWAYS | COPY | - |
| Number of Echoes | 2001,1014 | SL | | VNAP | COPY | - |
| Number of Locations | 2001,1015 | SS | | VNAP | COPY | - |
| Number of PC Directions | 2001,1016 | SS | | VNAP | COPY | - |
| Number of Phases MR | 2001,1017 | SL | | VNAP | COPY | - |

| | | | | | | |
|---------------------------------|-----------|-------|--|--------|-------------------|---|
| Number of Slices MR | 2001,1018 | SL | | VNAP | IMPLICIT, USER | - |
| Partial Matrix Scanned | 2001,1019 | CS | | VNAP | COPY | - |
| PC Velocity | 2001,101A | FL | | ALWAYS | COPY | - |
| Prepulse Delay | 2001,101B | FL | | VNAP | COPY | - |
| Prepulse Type | 2001,101C | CS | | VNAP | COPY | - |
| Reconstruction Number MR | 2001,101D | IS | | VNAP | IMPLICIT | - |
| Respiration Sync | 2001,101F | CS | | VNAP | COPY | - |
| Scanning Technique | 2001,1020 | LO | | ALWAYS | COPY | - |
| SPIR | 2001,1021 | CS | | VNAP | COPY | - |
| Water Fat Shift | 2001,1022 | FL | | VNAP | COPY | - |
| Flip Angle Philips | 2001,1023 | DS | | ALWAYS | COPY | - |
| Series is Interactive | 2001,1024 | CS | | VNAP | COPY | - |
| Echo Time Display | 2001,1025 | SH | | VNAP | COPY | - |
| Contrast Transfer Taste | 2001,1058 | UL | | ALWAYS | COPY | - |
| Number of Stacks | 2001,1060 | SL | | VNAP | IMPLICIT | - |
| Series Transmitted | 2001,1061 | CS | | VNAP | AUTO | - |
| Acquisition Number | 2001,107B | IS | | ALWAYS | COPY | - |
| Number of Dynamic Scans | 2001,1081 | IS | | VNAP | COPY | - |
| IsRawImage | 2001,10A1 | CS | | ANAP | COPY | - |
| Prospective Motion Correction | 2001,10F1 | FL | | ANAP | COPY | - |
| Retrospective Motion Correction | 2001,10F2 | FL | | ANAP | COPY | - |
| Private Creator Group 2005 | 2005,0010 | LO | | ALWAYS | FIXED | - |
| Number of Chemical Shift | 2005,1020 | SL | | ANAPCV | COPY | - |
| Syncra Scan Type | 2005,10A1 | CS | | ANAPCV | COPY | - |
| Diffusion Direction RL | 2005,10B0 | FL | | ANAP | COPY | - |
| Diffusion Direction AP | 2005,10B1 | FL | | ANAP | COPY | - |
| Diffusion Direction FH | 2005,10B2 | FL | | ANAP | COPY | - |
| Curve Dimensions | 50xx,0005 | US | | ALWAYS | AUTO | - |
| Number of Points | 50xx,0010 | US | | ALWAYS | AUTO | - |
| Type of Data | 50xx,0020 | CS | | ALWAYS | AUTO | - |
| Data Value Representation | 50xx,0103 | US | | ALWAYS | AUTO | - |
| Curve Data | 50xx,3000 | OB/OW | | ALWAYS | AUTO | - |
| Curve Description | 50xx,0022 | LO | | ANAPCV | AUTO | - |
| Axis Units | 50xx,0030 | SH | | ANAPCV | AUTO | - |
| Axis Labels | 50xx,0040 | SH | | ANAPCV | AUTO | - |
| Minimum Coordinate Value | 50xx,0104 | US | | ANAPCV | AUTO | - |
| Maximum Coordinate Value | 50xx,0105 | US | | ANAPCV | AUTO | - |
| Curve Range | 50xx,0106 | SH | | ANAPCV | AUTO | - |
| Curve Data Descriptor | 50xx,0110 | US | | ANAP | AUTO | - |
| Coordinate Start Value | 50xx,0112 | US | | ANAP | AUTO | - |
| Coordinate Step Value | 50xx,0114 | US | | ANAP | AUTO | - |
| Curve Label | 50xx,2500 | LO | | ANAPCV | AUTO | - |
| Referenced Overlay Sequence | 50xx,2600 | SQ | | ANAPCV | AUTO | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ANAP | AUTO | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ANAP | AUTO | - |
| >Referenced Overlay Group | 50xx,2610 | US | | ANAP | AUTO | - |

2.1.2.7. RT Structure Set Storage SOP Class

Table 109: IOD of Created RT Structure Set Storage SOP Class Instances

| Information Entity | Module | Presence Of Module |
|--------------------|---------------------------------------|--------------------|
| Patient | Patient Module | ALWAYS |
| Study | General Study Module | ALWAYS |
| Study | Patient Study Module | CONDITIONAL |
| Series | RT Series Module | ALWAYS |
| Equipment | General Equipment Module | ALWAYS |
| Structure Set | Structure Set Module | ALWAYS |
| Structure Set | ROI Contour Module | ALWAYS |
| Structure Set | RT ROI Observations Module | ALWAYS |
| Structure Set | SOP Common Module | ALWAYS |
| | Extended Dicom and Private attributes | CONDITIONAL |

Table 110: Patient Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|---------------|---------|
| Referenced Patient Sequence | 0008,1120 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Patient's Name | 0010,0010 | PN | | ALWAYS | COPY | - |
| Patient ID | 0010,0020 | LO | | ALWAYS | COPY | - |
| Patient's Birth Date | 0010,0030 | DA | | VNAP | COPY, USER | - |
| Patient's Birth Time | 0010,0032 | TM | | ANAPCV | COPY, USER | - |
| Patient's Sex | 0010,0040 | CS | | VNAP | COPY, USER | - |
| Other Patient IDs | 0010,1000 | LO | | ANAPCV | COPY, USER | - |
| Other Patient Names | 0010,1001 | PN | | ANAPCV | COPY, USER | - |
| Ethnic Group | 0010,2160 | SH | | ANAPCV | COPY | - |
| Patient Comments | 0010,4000 | LT | | ANAPCV | COPY | - |
| Issuer of Patient ID | 0010,0021 | LO | | ANAPCV | COPY, USER | - |

Table 111: General Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---|-----------|----|-------|-------------------|---------------|---------|
| Study Date | 0008,0020 | DA | | VNAP | COPY | - |
| Study Time | 0008,0030 | TM | | VNAP | COPY | - |
| Accession Number | 0008,0050 | SH | | VNAP | COPY, USER | - |
| Referring Physician's Name | 0008,0090 | PN | | VNAP | COPY, USER | - |
| Referring Physician Identification Sequence | 0008,0096 | SQ | | ANAPCV | COPY | - |
| Study Description | 0008,1030 | LO | | ANAPCV | COPY, USER | - |
| Procedure Code Sequence | 0008,1032 | SQ | | ANAPCV | COPY | - |
| Physician(s) of Record | 0008,1048 | PN | | ANAPCV | COPY | - |

| | | | | | | |
|--|-----------|----|--|--------|------|---|
| Physician(s) of Record Identification Sequence | 0008,1049 | SQ | | ANAPCV | COPY | - |
| Name of Physician(s) Reading Study | 0008,1060 | PN | | ANAPCV | COPY | - |
| Physician(s) Reading Study Identification Sequence | 0008,1062 | SQ | | ANAPCV | COPY | - |
| Referenced Study Sequence | 0008,1110 | SQ | | ANAPCV | COPY | - |
| Study Instance UID | 0020,000D | UI | | ALWAYS | AUTO | - |
| Study ID | 0020,0010 | SH | | VNAP | AUTO | - |

Table 112: Patient Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------------|-----------|----|-------|-------------------|------------|---------|
| Admitting Diagnoses Description | 0008,1080 | LO | | ANAPCV | COPY | - |
| Admitting Diagnoses Code Sequence | 0008,1084 | SQ | | ANAPCV | COPY | - |
| Patient's Age | 0010,1010 | AS | | ANAPCV | COPY, USER | - |
| Patient's Size | 0010,1020 | DS | | ANAPCV | COPY, USER | - |
| Patient's Weight | 0010,1030 | DS | | ANAPCV | COPY, USER | - |
| Occupation | 0010,2180 | SH | | ANAPCV | COPY | - |
| Additional Patient History | 0010,21B0 | LT | | ANAPCV | COPY | - |

Table 113: RT Series Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|------------|---------|
| Modality | 0008,0060 | CS | | ALWAYS | AUTO | - |
| Series Description | 0008,103E | LO | | ANAPEV | AUTO/ USER | - |
| Series Description Code Sequence | 0008,103F | SQ | | EMPTY | AUTO | - |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ | | ANAPEV | COPY | - |
| Series Instance UID | 0020,000E | UI | | ALWAYS | AUTO | - |
| Series Number | 0020,0011 | IS | | VNAP | AUTO | - |
| Request Attributes Sequence | 0040,0275 | SQ | | ANAPEV | COPY | - |
| Performed Procedure Step Start Date | 0040,0244 | DA | | ANAPEV | COPY | - |
| Performed Procedure Step Start Time | 0040,0245 | TM | | ANAPEV | COPY | - |
| Performed Procedure Step ID | 0040,0253 | SH | | ANAPEV | COPY | - |
| Performed Procedure Step Description | 0040,0254 | LO | | ANAPEV | COPY | - |
| Performed Protocol Code Sequence | 0040,0260 | SQ | | ANAPEV | COPY | - |
| Comments on the Performed Procedure Step | 0040,0280 | ST | | ANAPEV | COPY | - |

Table 114: General Equipment Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|----------------|-----------|----|-------|-------------------|--------|---------|
| Manufacturer | 0008,0070 | LO | | ALWAYS | COPY | - |

| | | | | | | |
|-------------------------------|-----------|-------|-----|--------|---------------|-----|
| Institution Name | 0008,0080 | LO | | ANAPCV | COPY, USER | - |
| Institution Address | 0008,0081 | ST | | ANAPCV | COPY, USER | - |
| Station Name | 0008,1010 | SH | | ALWAYS | COPY | - |
| Institutional Department Name | 0008,1040 | LO | | ANAPCV | COPY | - |
| Manufacturer's Model Name | 0008,1090 | LO | | ALWAYS | COPY | - |
| Device Serial Number | 0018,1000 | LO | | ANAPCV | COPY | - |
| Software Version(s) | 0018,1020 | LO | 4.0 | ALWAYS | COPY | 4.0 |
| Spatial Resolution | 0018,1050 | DS | | ANAPCV | COPY | - |
| Date of Last Calibration | 0018,1200 | DA | | ANAPCV | COPY | - |
| Time of Last Calibration | 0018,1201 | TM | | ANAPCV | COPY | - |
| Pixel Padding Value | 0028,0120 | US/SS | | ANAPCV | COPY | - |

Table 115: Structure Set Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---|-----------|----|-------|-------------------|---------------|---------|
| Instance Number | 0020,0013 | IS | | ANAPEV | AUTO | - |
| Structure Set Label | 3006,0002 | SH | | ALWAYS | AUTO/ USER | - |
| Structure Set Name | 3006,0004 | LO | | ANAPEV | AUTO/ USER | - |
| Structure Set Description | 3006,0006 | ST | | ANAPEV | AUTO/ USER | - |
| Structure Set Date | 3006,0008 | DA | | VNAP | AUTO | - |
| Structure Set Time | 3006,0009 | TM | | VNAP | AUTO | - |
| Referenced Frame of Reference Sequence | 3006,0010 | SQ | | ANAPEV | AUTO | - |
| >Frame of Reference UID | 0020,0052 | UI | | ALWAYS | AUTO | - |
| >RT Referenced Study Sequence | 3006,0012 | SQ | | ANAPEV | AUTO | - |
| >>RT Referenced Series Sequence | 3006,0014 | SQ | | ALWAYS | AUTO | - |
| >>>Series Instance UID | 0020,000E | UI | | ALWAYS | AUTO | - |
| >>>Contour Image Sequence | 3006,0016 | SQ | | ALWAYS | AUTO | - |
| >>>>Referenced Frame Number | 0008,1160 | IS | | ANAP | AUTO | - |
| >>>>Referenced Segment Number | 0062,000B | US | | ANAP | AUTO | - |
| >>>>Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | AUTO | - |
| >>>>Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | AUTO | - |
| >>Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | AUTO | - |
| >>Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | AUTO | - |
| >Frame of Reference Relationship Sequence | 3006,00C0 | SQ | | ANAPEV | AUTO | - |
| Structure Set ROI Sequence | 3006,0020 | SQ | | ANAPEV | AUTO | - |
| >Derivation Code Sequence | 0008,9215 | SQ | | ANAPEV | AUTO | - |
| >ROI Number | 3006,0022 | IS | | ALWAYS | AUTO | - |
| >Referenced Frame of Reference UID | 3006,0024 | UI | | ALWAYS | AUTO | - |
| >ROI Name | 3006,0026 | LO | | VNAP | AUTO | - |

| | | | | | | |
|-----------------------------|-----------|----|--|--------|------|---|
| >ROI Description | 3006,0028 | ST | | ANAPEV | AUTO | - |
| >ROI Volume | 3006,002C | DS | | ANAPEV | AUTO | - |
| >ROI Generation Algorithm | 3006,0036 | CS | | VNAP | AUTO | - |
| >ROI Generation Description | 3006,0038 | LO | | ANAPEV | AUTO | - |

Table 116: ROI Contour Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--------------------------------|-----------|----|-------|-------------------|--------|---------|
| ROI Contour Sequence | 3006,0039 | SQ | | ALWAYS | AUTO | - |
| >ROI Display Color | 3006,002A | IS | | ANAPEV | AUTO | - |
| >Contour Sequence | 3006,0040 | SQ | | ANAPEV | AUTO | - |
| >>Contour Image Sequence | 3006,0016 | SQ | | ANAPEV | AUTO | - |
| >>>Referenced Frame Number | 0008,1160 | IS | | ANAP | AUTO | - |
| >>>Referenced Segment Number | 0062,000B | US | | ANAP | AUTO | - |
| >>>Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | AUTO | - |
| >>>Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | AUTO | - |
| >>Contour Geometric Type | 3006,0042 | CS | | ALWAYS | AUTO | - |
| >>Contour Slab Thickness | 3006,0044 | DS | | ANAPEV | AUTO | - |
| >>Contour Offset Vector | 3006,0045 | DS | | ANAPEV | AUTO | - |
| >>Number of Contour Points | 3006,0046 | IS | | ALWAYS | AUTO | - |
| >>Contour Number | 3006,0048 | IS | | ANAPEV | AUTO | - |
| >>Attached Contours | 3006,0049 | IS | | ANAPEV | AUTO | - |
| >>Contour Data | 3006,0050 | DS | | ALWAYS | AUTO | - |
| >Referenced ROI Number | 3006,0084 | IS | | ALWAYS | AUTO | - |

Table 117: RT ROI Observations Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|--------|---------|
| RT ROI Observations Sequence | 3006,0080 | SQ | | ALWAYS | AUTO | - |
| >RT Related ROI Sequence | 3006,0030 | SQ | | EMPTY | AUTO | - |
| >> Referenced ROI Number | 3006,0084 | | | ALWAYS | AUTO | - |
| >>RT ROI Relationship | 3006,0033 | | | ANAPEV | AUTO | - |
| >Observation Number | 3006,0082 | IS | | ALWAYS | AUTO | - |
| >Referenced ROI Number | 3006,0084 | IS | | ALWAYS | AUTO | - |
| >ROI Observation Label | 3006,0085 | SH | | ANAPEV | AUTO | - |
| >RT ROI Identification Code Sequence | 3006,0086 | SQ | | ANAPEV | AUTO | - |
| >ROI Observation Description | 3006,0088 | ST | | ANAPEV | AUTO | - |
| >Related RT ROI Observations Sequence | 3006,00A0 | SQ | | ANAPEV | AUTO | - |
| >> Observation Number | 3006,0082 | | | ALWAYS | AUTO | - |
| >RT ROI Interpreted Type | 3006,00A4 | CS | | VNAP | AUTO | - |
| >ROI Interpreter | 3006,00A6 | PN | | VNAP | AUTO | - |
| >Material ID | 300A,00E1 | SH | | ANAPEV | AUTO | - |
| >ROI Physical Properties Sequence | 3006,00B0 | SQ | | ANAPEV | AUTO | - |
| >>ROI Physical Property | 3006,00B2 | CS | | ALWAYS | AUTO | - |
| >>ROI Elemental Composition Sequence | 3006,00B6 | SQ | | ANAP | AUTO | - |
| >>>ROI Elemental Composition Atomic Number | 3006,00B7 | US | | ALWAYS | AUTO | - |

| | | | | | | |
|---|-----------|----|--|--------|------|---|
| >>>ROI Elemental Composition Atomic Mass Fraction | 3006,00B8 | FL | | ALWAYS | AUTO | - |
| >>ROI Physical Property Value | 3006,00B4 | DS | | ALWAYS | AUTO | - |

Table 118: SOP Common Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------------|-----------|----|----------------------|-------------------|------------|----------------------|
| Specific Character Set | 0008,0005 | CS | | ANAP | AUTO, COPY | - |
| Instance Creation Date | 0008,0012 | DA | | ALWAYS | AUTO | - |
| Instance Creation Time | 0008,0013 | TM | | ALWAYS | AUTO | - |
| Instance Creator UID | 0008,0014 | UI | | ANAPCV | AUTO | - |
| SOP Class UID | 0008,0016 | UI | | ALWAYS | COPY | - |
| SOP Instance UID | 0008,0018 | UI | | ALWAYS | AUTO | - |
| Timezone Offset From UTC | 0008,0201 | SH | | ANAPCV | COPY | - |
| Contributing Equipment Sequence | 0018,A001 | SQ | | ANAP | AUTO | - |
| >Manufacturer | 0008,0070 | LO | Philips | ANAP | FIXED | PHILIPS. |
| >Institution Address | 0008,0081 | ST | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Station Name | 0008,1010 | SH | IntelliSpace Portal | ANAP | CONFIG | IntelliSpace Portal. |
| >Institutional Department Name | 0008,1040 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Software Version(s) | 0018,1020 | LO | 4.0 | ANAP | FIXED | 4.0 |
| >Device Serial Number | 0018,1000 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Institution Name | 0008,0080 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Manufacturer's Model Name | 0008,1090 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Purpose of Reference Code Sequence | 0040,A170 | SQ | | ANAP | FIXED | - |
| >>Code Value | 0008,0100 | SH | 109102 | ANAP | FIXED | 109102 |
| >>Coding Scheme Designator | 0008,0102 | SH | DCM | ANAP | FIXED | DCM |
| >>Code Meaning | 0008,0104 | LO | Processing Equipment | ANAP | FIXED | Processing Equipment |

Table 119: Extended DICOM and private attributes for RT Structure Set Storage SOP Class Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---------------------|-----------|----|-------|-------------------|--------|---------|
| Modalities in Study | 0008,0061 | CS | | ANAPEV | AUTO | - |
| Patient Orientation | 0020,0020 | CS | | ANAPEV | AUTO | - |

2.1.2.8. Secondary Capture Image Storage SOP Class

Table 120: IOD of Created Secondary Capture Image Storage SOP Class Instances

| Information Entity | Module | Presence Of Module |
|--------------------|--------------------------|--------------------|
| Patient | Patient Module | ALWAYS |
| Study | General Study Module | ALWAYS |
| Study | Patient Study Module | CONDITIONAL |
| Series | General Series Module | ALWAYS |
| Equipment | General Equipment Module | ALWAYS |
| Equipment | SC Equipment Module | ALWAYS |
| Image | General Image Module | ALWAYS |
| Image | Image Pixel Module | ALWAYS |
| Image | Device Module | CONDITIONAL |
| Image | SC Image Module | ALWAYS |

| | | |
|-------|---------------------------------------|-------------|
| Image | Overlay Plane Module | CONDITIONAL |
| Image | Modality LUT Module | CONDITIONAL |
| Image | VOI LUT Module | CONDITIONAL |
| Image | SOP Common Module | ALWAYS |
| | Extended Dicom and Private attributes | CONDITIONAL |

Table 121: Patient Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|---------------|---------|
| Referenced Patient Sequence | 0008,1120 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Patient's Name | 0010,0010 | PN | | VNAP | COPY | - |
| Patient ID | 0010,0020 | LO | | VNAP | COPY | - |
| Patient's Birth Date | 0010,0030 | DA | | VNAP | COPY, USER | - |
| Patient's Birth Time | 0010,0032 | TM | | ANAPCV | COPY, USER | - |
| Patient's Sex | 0010,0040 | CS | | VNAP | COPY, USER | - |
| Other Patient IDs | 0010,1000 | LO | | ANAPCV | COPY, USER | - |
| Other Patient Names | 0010,1001 | PN | | ANAPCV | COPY, USER | - |
| Ethnic Group | 0010,2160 | SH | | ANAPCV | COPY | - |
| Patient Comments | 0010,4000 | LT | | ANAPCV | COPY | - |
| Issuer of Patient ID | 0010,0021 | LO | | ANAPCV | COPY, USER | - |

Table 122: General Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|---------------|---------|
| Study Date | 0008,0020 | DA | | VNAP | COPY | - |
| Study Time | 0008,0030 | TM | | VNAP | COPY | - |
| Accession Number | 0008,0050 | SH | | VNAP | COPY, USER | - |
| Referring Physician's Name | 0008,0090 | PN | | VNAP | COPY, USER | - |
| Referring Physician Identification Sequence | 0008,0096 | SQ | | ANAPCV | COPY | - |
| Study Description | 0008,1030 | LO | | ANAPCV | COPY, USER | - |
| Procedure Code Sequence | 0008,1032 | SQ | | ANAPCV | COPY | - |
| Physician(s) of Record | 0008,1048 | PN | | ANAPCV | COPY | - |
| Physician(s) of Record Identification Sequence | 0008,1049 | SQ | | ANAPCV | COPY | - |
| Name of Physician(s) Reading Study | 0008,1060 | PN | | ANAPCV | COPY | - |
| Physician(s) Reading Study Identification Sequence | 0008,1062 | SQ | | ANAPCV | COPY | - |
| Referenced Study Sequence | 0008,1110 | SQ | | ANAPCV | COPY | - |
| Study Instance UID | 0020,000D | UI | | ALWAYS | AUTO, COPY | - |

| | | | | | | |
|----------|-----------|----|--|------|---------------|---|
| Study ID | 0020,0010 | SH | | VNAP | AUTO, COPY | - |
|----------|-----------|----|--|------|---------------|---|

Table 123: Patient Study Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------------------|-----------|----|-------|-------------------|---------------|---------|
| Admitting Diagnoses Description | 0008,1080 | LO | | ANAPCV | COPY | - |
| Admitting Diagnoses Code Sequence | 0008,1084 | SQ | | ANAPCV | COPY | - |
| Patient's Age | 0010,1010 | AS | | ANAPCV | COPY, USER | - |
| Patient's Size | 0010,1020 | DS | | ANAPCV | COPY, USER | - |
| Patient's Weight | 0010,1030 | DS | | ANAPCV | COPY, USER | - |
| Occupation | 0010,2180 | SH | | ANAPCV | COPY | - |
| Additional Patient History | 0010,21B0 | LT | | ANAPCV | COPY | - |

Table 124: General Series Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|-------|-------|-------------------|---------------|---------|
| Series Date | 0008,0021 | DA | | ANAPCV | AUTO | - |
| Series Time | 0008,0031 | TM | | ANAPCV | AUTO | - |
| Modality | 0008,0060 | CS | | ALWAYS | COPY | - |
| Series Description | 0008,103E | LO | | ANAPCV | COPY, USER | - |
| Performing Physician's Name | 0008,1050 | PN | | ANAPCV | COPY | - |
| Operators' Name | 0008,1070 | PN | | ANAPCV | COPY, USER | - |
| Operator Identification Sequence | 0008,1072 | SQ | | ANAPCV | COPY | - |
| Referenced Performed Procedure Step Sequence | 0008,1111 | SQ | | ANAPCV | COPY | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ALWAYS | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ALWAYS | COPY | - |
| Body Part Examined | 0018,0015 | CS | | ANAPCV | COPY | - |
| Protocol Name | 0018,1030 | LO | | ANAPCV | COPY | - |
| Patient Position | 0018,5100 | CS | | ANAPCV | COPY | - |
| Series Instance UID | 0020,000E | UI | | ALWAYS | AUTO | - |
| Series Number | 0020,0011 | IS | | VNAP | AUTO | - |
| Laterality | 0020,0060 | CS | | ANAPCV | COPY | - |
| Smallest Pixel Value in Series | 0028,0108 | US/SS | | ANAPCV | AUTO, COPY | - |
| Largest Pixel Value in Series | 0028,0109 | US/SS | | ANAPCV | AUTO, COPY | - |
| Request Attributes Sequence | 0040,0275 | SQ | | ANAPCV | COPY | - |
| Performed Procedure Step Start Date | 0040,0244 | DA | | ANAPCV | COPY | - |
| Performed Procedure Step Start Time | 0040,0245 | TM | | ANAPCV | COPY | - |
| Performed Procedure Step Description | 0040,0254 | LO | | ANAPCV | COPY | - |
| Performed Protocol Code Sequence | 0040,0260 | SQ | | ANAPCV | COPY | - |

| | | | | | | |
|--|-----------|----|--|--------|------|---|
| Comments on the Performed Procedure Step | 0040,0280 | ST | | ANAPCV | COPY | - |
|--|-----------|----|--|--------|------|---|

Table 125: General Equipment Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------|-----------|-------|-------------------------|-------------------|------------------|-------------------------|
| Manufacturer | 0008,0070 | LO | | VNAP | COPY | Philips Healthcare. |
| Institution Name | 0008,0080 | LO | | ANAPCV | AUTO, COPY, USER | - |
| Institution Address | 0008,0081 | ST | | ANAPCV | AUTO, COPY, USER | - |
| Station Name | 0008,1010 | SH | | ANAPCV | COPY | - |
| Institutional Department Name | 0008,1040 | LO | | ANAPCV | AUTO | - |
| Manufacturer's Model Name | 0008,1090 | LO | | ANAPCV | COPY | - |
| Device Serial Number | 0018,1000 | LO | | ANAPCV | COPY | - |
| Software Version(s) | 0018,1020 | LO | IntelliSpace Portal 4.0 | ANAPCV | COPY | IntelliSpace Portal 4.0 |
| Spatial Resolution | 0018,1050 | DS | | ANAPCV | COPY | - |
| Date of Last Calibration | 0018,1200 | DA | | ANAPCV | COPY | - |
| Time of Last Calibration | 0018,1201 | TM | | ANAPCV | COPY | - |
| Pixel Padding Value | 0028,0120 | US/SS | | ANAPCV | COPY | - |

Table 126: SC Equipment Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|--|-----------|----|-------|-------------------|-------------|---------|
| Modality | 0008,0060 | CS | CT | ANAPCV | AUTO, FIXED | - |
| Conversion Type | 0008,0064 | CS | WSD | ALWAYS | AUTO | - |
| Secondary Capture Device ID | 0018,1010 | LO | | ANAPCV | AUTO | - |
| Secondary Capture Device Manufacturer | 0018,1016 | LO | | ANAPCV | AUTO | - |
| Secondary Capture Device Manufacturer's Model Name | 0018,1018 | LO | | ANAPCV | AUTO | - |
| Secondary Capture Device Software Version(s) | 0018,1019 | LO | | ANAPCV | AUTO | - |
| Video Image Format Acquired | 0018,1022 | SH | | ANAPCV | AUTO | - |
| Digital Image Format Acquired | 0018,1023 | LO | | ANAPCV | AUTO | - |

Table 127: General Image Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|--------------------------------------|-------------------|------------|---------|
| Image Type | 0008,0008 | CS | Value 1: DERIVED, Value 2: SECONDARY | ANAPCV | AUTO | - |
| Acquisition Date | 0008,0022 | DA | | ANAPCV | AUTO, COPY | - |
| Content Date | 0008,0023 | DA | | ANAPCV | AUTO, COPY | - |
| Acquisition Datetime | 0008,002A | DT | | ANAPCV | COPY | - |
| Acquisition Time | 0008,0032 | TM | | ANAPCV | COPY | - |
| Content Time | 0008,0033 | TM | | ANAPCV | COPY | - |
| Referenced Waveform Sequence | 0008,113A | SQ | | ANAP | AUTO | - |
| Referenced Image Sequence | 0008,1140 | SQ | | ANAPCV | AUTO, COPY | - |

| | | | | | | |
|--|-----------|-------|-----|--------|---------------|---|
| >Referenced SOP Class UID | 0008,1150 | UI | | ANAPCV | COPY | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ANAPCV | COPY | - |
| Derivation Description | 0008,2111 | ST | | ANAPCV | AUTO, COPY | - |
| Source Image Sequence | 0008,2112 | SQ | | ANAPCV | AUTO, COPY | - |
| Derivation Code Sequence | 0008,9215 | SQ | | ANAPCV | AUTO, COPY | - |
| Acquisition Number | 0020,0012 | IS | | ANAPCV | COPY | - |
| Instance Number | 0020,0013 | IS | | ANAPCV | AUTO, COPY | - |
| Patient Orientation | 0020,0020 | CS | | ANAPCV | AUTO, COPY | - |
| Images in Acquisition | 0020,1002 | IS | | ANAPCV | AUTO, COPY | - |
| Image Comments | 0020,4000 | LT | | ANAPCV | AUTO, COPY | - |
| Quality Control Image | 0028,0300 | CS | | ANAPCV | AUTO, COPY | - |
| Burned In Annotation | 0028,0301 | CS | YES | ANAPCV | AUTO, COPY | - |
| Lossy Image Compression | 0028,2110 | CS | | ANAPCV | AUTO, COPY | - |
| Lossy Image Compression Ratio | 0028,2112 | DS | | ANAPCV | AUTO, COPY | - |
| Lossy Image Compression Method | 0028,2114 | CS | | ANAPCV | AUTO | - |
| Icon Image Sequence | 0088,0200 | SQ | | ANAPCV | AUTO, COPY | - |
| >Samples per Pixel | 0028,0002 | US | | ALWAYS | COPY | - |
| >Photometric Interpretation | 0028,0004 | CS | | ALWAYS | COPY | - |
| >Planar Configuration | 0028,0006 | US | | ANAP | COPY | - |
| >Rows | 0028,0010 | US | | ALWAYS | COPY | - |
| >Columns | 0028,0011 | US | | ALWAYS | COPY | - |
| >Pixel Aspect Ratio | 0028,0034 | IS | | ANAP | COPY | - |
| >Bits Allocated | 0028,0100 | US | | ALWAYS | COPY | - |
| >Bits Stored | 0028,0101 | US | | ALWAYS | COPY | - |
| >High Bit | 0028,0102 | US | | ALWAYS | COPY | - |
| >Pixel Representation | 0028,0103 | US | | ALWAYS | COPY | - |
| >Smallest Image Pixel Value | 0028,0106 | US/SS | | ANAPCV | COPY | - |
| >Largest Image Pixel Value | 0028,0107 | US/SS | | ANAPCV | COPY | - |
| >Red Palette Color Lookup Table Descriptor | 0028,1101 | US/SS | | ANAP | COPY | - |
| >Green Palette Color Lookup Table Descriptor | 0028,1102 | US/SS | | ANAP | COPY | - |
| >Blue Palette Color Lookup Table Descriptor | 0028,1103 | US/SS | | ANAP | COPY | - |
| >Red Palette Color Lookup Table Data | 0028,1201 | OW | | ANAP | COPY | - |
| >Green Palette Color Lookup Table Data | 0028,1202 | OW | | ANAP | COPY | - |
| >Blue Palette Color Lookup Table Data | 0028,1203 | OW | | ANAP | COPY | - |
| >Pixel Data | 7FE0,0010 | OW/OB | | ANAP | COPY | - |

| | | | | | | |
|------------------------|-----------|----|--|--------|---------------|---|
| Presentation LUT Shape | 2050,0020 | CS | | ANAPCV | AUTO, COPY | - |
|------------------------|-----------|----|--|--------|---------------|---|

Table 128: Image Pixel Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---|-----------|-------|-------|-------------------|---------------|---------|
| Samples per Pixel | 0028,0002 | US | 3 | ALWAYS | AUTO | - |
| Photometric Interpretation | 0028,0004 | CS | | ALWAYS | AUTO | - |
| Planar Configuration | 0028,0006 | US | | ANAP | AUTO, COPY | - |
| Rows | 0028,0010 | US | | ALWAYS | AUTO, COPY | - |
| Columns | 0028,0011 | US | | ALWAYS | AUTO, COPY | - |
| Pixel Aspect Ratio | 0028,0034 | IS | | ANAP | AUTO, COPY | - |
| Bits Allocated | 0028,0100 | US | 8 | ALWAYS | AUTO | - |
| Bits Stored | 0028,0101 | US | 8 | ALWAYS | AUTO | - |
| High Bit | 0028,0102 | US | 7 | ALWAYS | AUTO | - |
| Pixel Representation | 0028,0103 | US | 0 | ALWAYS | AUTO, COPY | - |
| Smallest Image Pixel Value | 0028,0106 | US/SS | | ANAPCV | AUTO, COPY | - |
| Largest Image Pixel Value | 0028,0107 | US/SS | | ANAPCV | AUTO, COPY | - |
| Red Palette Color Lookup Table Descriptor | 0028,1101 | US/SS | | ANAP | AUTO, COPY | - |
| Green Palette Color Lookup Table Descriptor | 0028,1102 | US/SS | | ANAP | AUTO, COPY | - |
| Blue Palette Color Lookup Table Descriptor | 0028,1103 | US/SS | | ANAP | AUTO, COPY | - |
| Red Palette Color Lookup Table Data | 0028,1201 | OW | | ANAP | AUTO, COPY | - |
| Green Palette Color Lookup Table Data | 0028,1202 | OW | | ANAP | AUTO, COPY | - |
| Blue Palette Color Lookup Table Data | 0028,1203 | OW | | ANAP | AUTO, COPY | - |
| Pixel Data | 7FE0,0010 | OW/OB | | ANAP | AUTO, COPY | - |

Table 129: Device Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|----------------------------|-----------|----|-------|-------------------|--------|---------|
| Device Sequence | 0050,0010 | SQ | | ANAPCV | AUTO | - |
| >Manufacturer | 0008,0070 | LO | | ANAPCV | AUTO | - |
| >Manufacturer's Model Name | 0008,1090 | LO | | ANAPCV | AUTO | - |
| >Device Serial Number | 0018,1000 | LO | | ANAPCV | AUTO | - |
| >Device ID | 0018,1003 | LO | | ANAPCV | AUTO | - |
| >Device Length | 0050,0014 | DS | | ANAPCV | AUTO | - |
| >Device Diameter | 0050,0016 | DS | | ANAPCV | AUTO | - |
| >Device Diameter Units | 0050,0017 | CS | | ANAPCV | AUTO | - |
| >Device Volume | 0050,0018 | DS | | ANAPCV | AUTO | - |
| >Intermarker Distance | 0050,0019 | DS | | ANAPCV | AUTO | - |
| >Device Description | 0050,0020 | LO | | ANAPCV | AUTO | - |

| | | | | | | |
|--------------------------------------|-----------|----|--|--------|------|---|
| >Code Value | 0008,0100 | SH | | ALWAYS | AUTO | - |
| >Coding Scheme Designator | 0008,0102 | SH | | ALWAYS | AUTO | - |
| >Coding Scheme Version | 0008,0103 | SH | | ANAP | AUTO | - |
| >Code Meaning | 0008,0104 | LO | | ALWAYS | AUTO | - |
| >Mapping Resource | 0008,0105 | CS | | ANAP | AUTO | - |
| >Context Group Version | 0008,0106 | DT | | ANAP | AUTO | - |
| >Context Group Local Version | 0008,0107 | DT | | ANAP | AUTO | - |
| >Context Group Extension Flag | 0008,010B | CS | | ANAPCV | AUTO | - |
| >Context Group Extension Creator UID | 0008,010D | UI | | ANAP | AUTO | - |
| >Context Identifier | 0008,010F | CS | | ANAPCV | AUTO | - |

Table 130: SC Image Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---------------------------|-----------|----|-------|-------------------|--------|---------|
| Date of Secondary Capture | 0018,1012 | DA | | ANAPCV | AUTO | - |
| Time of Secondary Capture | 0018,1014 | TM | | ANAPCV | AUTO | - |

Table 131: Overlay Plane Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|------------------------|-----------|-------|-------|-------------------|--------|---------|
| Overlay Rows | 6000,0010 | US | | ALWAYS | AUTO | - |
| Overlay Columns | 6000,0011 | US | | ALWAYS | AUTO | - |
| Overlay Description | 6000,0022 | LO | | ANAPCV | AUTO | - |
| Overlay Type | 6000,0040 | CS | | ALWAYS | AUTO | - |
| Overlay Subtype | 6000,0045 | LO | | ANAPCV | AUTO | - |
| Overlay Origin | 6000,0050 | SS | | ALWAYS | AUTO | - |
| Overlay Bits Allocated | 6000,0100 | US | | ALWAYS | AUTO | - |
| Overlay Bit Position | 6000,0102 | US | | ALWAYS | AUTO | - |
| ROI Area | 6000,1301 | IS | | ANAPCV | AUTO | - |
| ROI Mean | 6000,1302 | DS | | ANAPCV | AUTO | - |
| ROI Standard Deviation | 6000,1303 | DS | | ANAPCV | AUTO | - |
| Overlay Label | 6000,1500 | LO | | ANAPCV | AUTO | - |
| Overlay Data | 6000,3000 | OW/OB | | ANAP | AUTO | - |

Table 132: Modality LUT Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-----------------------|-----------|-------|-------|-------------------|--------|---------|
| Rescale Intercept | 0028,1052 | DS | | ANAP | AUTO | - |
| Rescale Slope | 0028,1053 | DS | | ANAP | AUTO | - |
| Rescale Type | 0028,1054 | LO | | ANAP | AUTO | - |
| Modality LUT Sequence | 0028,3000 | SQ | | ANAP | AUTO | - |
| >LUT Descriptor | 0028,3002 | US/SS | | ANAP | AUTO | - |
| >LUT Explanation | 0028,3003 | LO | | ANAP | AUTO | - |
| >Modality LUT Type | 0028,3004 | LO | | ANAP | AUTO | - |
| >LUT Data | 0028,3006 | US/OW | | ANAP | AUTO | - |

Table 133: VOI LUT Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|----------------|-----------|----|-------|-------------------|------------|---------|
| Window Center | 0028,1050 | DS | | ANAPCV | AUTO, COPY | - |

| | | | | | | |
|-----------------------------------|-----------|----|--|--------|------------|---|
| Window Width | 0028,1051 | DS | | ANAP | AUTO, COPY | - |
| Window Center & Width Explanation | 0028,1055 | LO | | ANAPCV | AUTO, COPY | - |
| VOI LUT Sequence | 0028,3010 | SQ | | ANAPCV | AUTO, COPY | - |

Table 134: SOP Common Module

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|-------------------------------------|-----------|----|---------------------------|-------------------|-------------|----------------------|
| Specific Character Set | 0008,0005 | CS | | ANAP | AUTO, COPY | - |
| Instance Creation Date | 0008,0012 | DA | | ANAPCV | AUTO | - |
| Instance Creation Time | 0008,0013 | TM | | ANAPCV | AUTO | - |
| Instance Creator UID | 0008,0014 | UI | | ANAPCV | AUTO | - |
| SOP Class UID | 0008,0016 | UI | 1.2.840.10008.5.1.4.1.1.7 | ALWAYS | FIXED, COPY | - |
| SOP Instance UID | 0008,0018 | UI | | ALWAYS | AUTO | - |
| Timezone Offset From UTC | 0008,0201 | SH | | ANAPCV | COPY | - |
| Contributing Equipment Sequence | 0018,A001 | SQ | | ANAP | AUTO | - |
| >Manufacturer | 0008,0070 | LO | Philips | ANAP | FIXED | PHILIPS. |
| >Institution Address | 0008,0081 | ST | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Station Name | 0008,1010 | SH | IntelliSpace Portal | ANAP | CONFIG | IntelliSpace Portal. |
| >Institutional Department Name | 0008,1040 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Software Version(s) | 0018,1020 | LO | 4.0 | ANAP | FIXED | 4.0 |
| >Device Serial Number | 0018,1000 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Institution Name | 0008,0080 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Manufacturer's Model Name | 0008,1090 | LO | IntelliSpace Portal | ANAP | FIXED | IntelliSpace Portal. |
| >Purpose of Reference Code Sequence | 0040,A170 | SQ | | ANAP | FIXED | - |
| >>Code Value | 0008,0100 | SH | 109102 | ANAP | FIXED | 109102 |
| >>Coding Scheme Designator | 0008,0102 | SH | DCM | ANAP | FIXED | DCM |
| >>Code Meaning | 0008,0104 | LO | Processing Equipment | ANAP | FIXED | Processing Equipment |

Table 135: Extended DICOM and private attributes for Secondary Capture Image Storage SOP Class Instances

| Attribute Name | Tag | VR | Value | Presence of Value | Source | Comment |
|---------------------------|-----------|-------|----------|-------------------|--------|----------|
| Series Type (PIIM) | 2001,106E | SH | SCSERIES | ALWAYS | AUTO | SCSERIES |
| Curve Dimensions | 50xx,0005 | US | | ALWAYS | AUTO | - |
| Number of Points | 50xx,0010 | US | | ALWAYS | AUTO | - |
| Type of Data | 50xx,0020 | CS | | ALWAYS | AUTO | - |
| Data Value Representation | 50xx,0103 | US | | ALWAYS | AUTO | - |
| Curve Data | 50xx,3000 | OB/OW | | ALWAYS | AUTO | - |
| Curve Description | 50xx,0022 | LO | | ANAPCV | AUTO | - |
| Axis Units | 50xx,0030 | SH | | ANAPCV | AUTO | - |
| Axis Labels | 50xx,0040 | SH | | ANAPCV | AUTO | - |
| Minimum Coordinate Value | 50xx,0104 | US | | ANAPCV | AUTO | - |
| Maximum Coordinate Value | 50xx,0105 | US | | ANAPCV | AUTO | - |
| Curve Range | 50xx,0106 | SH | | ANAPCV | AUTO | - |
| Curve Data Descriptor | 50xx,0110 | US | | ANAP | AUTO | - |
| Coordinate Start Value | 50xx,0112 | US | | ANAP | AUTO | - |

| | | | | | | |
|------------------------------|-----------|----|--|--------|------|---|
| Coordinate Step Value | 50xx,0114 | US | | ANAP | AUTO | - |
| Curve Label | 50xx,2500 | LO | | ANAPCV | AUTO | - |
| Referenced Overlay Sequence | 50xx,2600 | SQ | | ANAPCV | AUTO | - |
| >Referenced SOP Class UID | 0008,1150 | UI | | ANAP | AUTO | - |
| >Referenced SOP Instance UID | 0008,1155 | UI | | ANAP | AUTO | - |
| >Referenced Overlay Group | 50xx,2610 | US | | ANAP | AUTO | - |