Autoclavable Temperature Probes

Instructions for Use

REF 989803203561 (Adult)

Esophageal/Rectal 12 Fr
Temperature Probe

**REF 989803203571 (Pediatric/
Adult)**

Esophageal/Rectal 10 Fr
Temperature Probe

**REF 989803203581 (Adult/
Pediatric)**

Skin Temperature Probe

PHILIPS

Contents

(BG) Указания за употреба.....	1	(LV) Lietošanas instrukcija	87
(CS) Návod k použití	6	(MK) Упатства за користење	92
(DA) Brugervejledning	11	(NL) Gebruiksaanwijzing	98
(DE) Gebrauchsanweisung	16	(NO) Brukerhåndbok	103
(EL) Οδηγίες χρήσης	21	(PL) Instrukcja obsługi.....	108
(EN) Instructions for Use	27	(PT-BR) Instruções de uso	113
(ES) Instrucciones de uso	32	(RO) Instrucțiuni de utilizare	118
(ET) Kasutusjuhend	37	(RU) Инструкция по эксплуатации.....	123
(FI) Käyttöohjeet.....	42	(SK) Návod na použitie.....	129
(FR) Manuel d'utilisation	47	(SL) Navodila za uporabo	134
(HR) Upute za uporabu.....	52	(SR) Uputstvo za upotrebu.....	139
(HU) Használati útmutató.....	57	(SV) Bruksanvisning.....	144
(IT) Istruzioni d'uso	62	(TR) Kullanım Talimatları	149
(JA) ユーザーズガイド	67	(UK) Інструкції з використання	154
(KK) Пайдалану нұсқаулығы	72	(VI) Hướng Dẫn Sử Dụng.....	160
(KO) 사용 설명서.....	77	(ZH-CN) 使用说明	165
(LT) Naudojimo nurodymai.....	82	(ZH-TW) 使用說明	169

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Intended Use

The Philips reusable, autoclavable temperature probes are intended for use when continuous monitoring of human esophageal, rectal, or skin temperature is required. These probes should be applied and used only by licensed and professionally trained clinicians, such as nurses and physicians, in healthcare facilities. These probes, connected only to compatible patient monitors (see Compatibility below), capture the patient's temperature, providing the user with the patient's continuous temperature.

Indications for Use

The Esophageal/Rectal and Skin temperature probes are indicated for continuous body temperature monitoring. Continuous temperature monitoring provides the clinician information to help diagnose conditions that cause changes in body temperature and guide a patient's temperature management.

The Philips 989803203561 Esophageal/Rectal probe is indicated for adult patients, the 989803203571 Esophageal/Rectal probe is indicated primarily for pediatric patients and secondarily for adult patients, and the 989803203581 Skin probe is indicated for adult and pediatric patients.

Contraindication

There are no known contraindications.

Adverse Reactions During Patient Application

- Adverse reactions reported during application of esophageal/rectal probes include: accidental trachea or bronchial intubation accompanied with airway obstruction, esophageal or rectal abrasion and/or perforation, pharyngeal abrasion and tissue burns due to aberrant electrosurgical RF current pathways.
- Adverse reactions reported during application of skin probes include skin irritation and tissue burns due to aberrant electrosurgical RF current pathways.

Measurement Range and Accuracy

The Philips reusable, autoclavable temperature probes are direct mode thermometers. Total system accuracy is the sum of the temperature probe and monitor accuracies. Check the monitor instructions to determine total measurement system accuracy.

Temperature accuracy of these probes at the site of measurement is:

$\pm 0.1^{\circ}\text{C}$ within the temperature range of 25°C to 45°C .

The minimum measurement time for accurate temperature readings is three minutes.

Compatibility

These temperature probes connect to Philips IntelliVue monitor (model MX100). These probes are intended for use with Philips monitors designed to accept 400 series temperature probes. For compatibility with other patient monitors and defibrillators, refer to the instrument's Instructions for Use.

Note: Check compatibility of monitoring equipment and probe prior to use. Incompatible components may result in performance issues.

Warnings

- Be sure to read and understand all warnings listed below.
- These probes are not designed for use involving patient temperature control in incubators.
- These probes are sold non-sterile. Before use, sterilize the Esophageal/Rectal temperature probe. (Refer to Autoclave Sterilization Section below.)
- Before use, inspect and replace damaged or deteriorated temperature probes.
- Connect temperature probes only to Philips monitors that use 400 series probe technology.
- Do NOT apply these temperature probes to patients who are undergoing Magnetic Resonance Imaging (MRI) procedures.
- Use caution when applying, inserting, or removing a temperature probe from a patient.
- Assure probe is not positioned where it could become entangled around patient, or cause choking, strangulation, or inhibit circulation in extremities.
- Carefully position probe to avoid it becoming entangled around Intravenous (IV) tubes and restrict IV flow.
- Mishandling of the probes could result in damage to the internal wires and the loss of electrical isolation or improper temperature readings.
- The use of the Esophageal/Rectal Temperature Probe may be cautioned in neonates and small infants undergoing tracheostomy or internal jugular vein catheterization and laser surgical procedures.
- The use of the Skin Temperature Probe may be cautioned in neonates and small infants whose skin may be more sensitive to adhesives and placement of surface probes.
- All wire-lead patient-connected transducer assemblies are subject to reading error, local heating and possible damage from high-intensity sources of RF energy. Electrosurgical equipment represents one such source since currents may seek alternate paths to ground through probe cables and associated instruments and patient burns may result. If possible, remove the temperature probe from patient contact before activating the electrosurgical equipment or other RF source to minimize the possibility of localized tissue burns. If temperature probes must be used simultaneously with electrosurgical equipment, use currently acceptable procedures to locate both active and ground electrodes of the electrosurgical equipment in close proximity to each other, and locate the temperature probe outside the RF current field as far as possible.

Getting Started

*Note: Read and understand all **Warnings** provided before applying your esophageal/rectal and skin temperature probe.*

Prior to patient application visually inspect the probe. Check that the probe and cable are not damaged, connector is not bent/damaged, and probe tip is not torn. Discard if any of these conditions exist.

Using the Esophageal/Rectal Probes

1. Ensure the probe is compatible with the monitoring equipment. Incompatible components may result in incorrect readings.
2. If patient is being intubated with an E.T. tube, perform intubation prior to inserting the temperature probe into the esophagus.
3. Lubricate probe prior to insertion, and place the probe in position in accordance with currently accepted medical procedures.
4. Verify the probe's location in accordance with acceptable medical procedures.
5. Place the connector end of the probe into the monitor. Check to make sure the mating ends are secure.
6. Three minutes of equilibration time are necessary prior to the initiation of monitoring.
7. When the temperature probe ceases to function, dispose per local medical waste regulations.

Using the Skin Probe

1. Ensure the probe is compatible with the monitoring equipment. Incompatible components may result in incorrect readings.
2. Select site for monitoring skin temperature according to currently accepted medical practices.
3. Place the connector end of the probe into the monitor. Check to make sure the mating ends are secure.
4. For best probe adhesion use tape or adhesive foam, clean and dry site selected for monitoring.
5. Three minutes of equilibration time are necessary prior to the initiation of monitoring.
6. When the temperature probe ceases to function, dispose per local medical waste regulations.

Electromagnetic Interference

If electromagnetic interference (EMI) is encountered, a number of things can be done to mitigate the problem:

- Eliminate the source. Possible sources of EMI can be turned off or moved away to reduce their strength.
- Attenuate the coupling. If coupling is through accessory leads, the interference may be reduced by moving and/or rearranging leads.
- Add external attenuators. If EMI becomes an unusually difficult problem, external devices may be of help.

Cleaning, Disinfection, and Sterilization Instructions for Reusable Temperature Probes

To reduce risk of cross-infection always clean all probes (disinfect/sterilize accordingly) before applying to a different patient. Esophageal/Rectal Probes must be cleaned and sterilized between patient use. Skin probes can be cleaned and disinfected or cleaned and sterilized between patient use. Use only the approved agents listed in the instructions below. Other agents may damage the probe or its connecting wires, shorten the product lifetime, or cause safety hazards.

Note: These temperature probes are capable of withstanding up to 100 cleaning cycles, and up to 100 disinfection cycles.

Caution: Do not expose the metallic connector pins to any cleaning or disinfection solution. Exposing the conductive components could cause the device to malfunction and/or give incorrect results.

Cleaning Instructions

1. Prepare an enzymatic detergent solution, (for example, Prolystica® 2x Concentrate Presoak and Cleaner, or Enzol®) per manufacturer's instructions.
2. Ensure the temperature probe is disconnected from the monitor, and then soak the temperature probe as shown in Figure 1 below (do not submerge the connector) in the prepared enzymatic detergent solution per manufacturer's instructions.
3. Remove the temperature probe from the prepared detergent solution.
4. Dampen a lint-free cloth with the detergent solution and then wipe the entire temperature probe including the exterior surface of the connector. Do not expose the interior electrical components of the connector to moisture. When wiping clean, hold the probe in one hand at the sensing tip and wipe the probe toward the connector end. Be careful not to apply excessive pressure that could stretch the cable jacket and break the internal wires. Wipe all surfaces a minimum of one time.

5. Use a nylon brush dampened in the prepared detergent solution to remove debris from connector (refer to Figure 2 below). Brush the outside of the connector for a minimum of 30 seconds and until visibly clean.

Figure 1


Figure 2

6. Thoroughly rinse the temperature probe under running tap water (20-25°C) for a minimum of 30 seconds and until any detergent or chemical residue is no longer visible. Do not expose the interior electrical components of the connector to moisture.
7. Perform a secondary rinse by submerging the temperature probe into Deionized water (DI) for a minimum of one (1) minute. Again, do not allow the electrical components of the connector to get wet.
8. Remove the device from the DI water and gently wipe the entire temperature probe cable including the exterior surface of the connector with a lint-free cloth. Allow to air dry completely.
9. Visually inspect the device to ensure all surfaces are clean. If visible debris is observed, repeat the cleaning steps above.

Disinfection Instructions for Skin Probes

1. Disinfect probe by immersing it into 70% Isopropyl Alcohol, being careful to avoid immersing the connector as shown in Figure below. Immerse the probe for 10 minutes so that all external surfaces of the probe (excluding the connector) are in contact with the solution.

2. Disinfect connector by wrapping the connector in a dampened lint-free cloth with 70% Isopropyl Alcohol for 10 minutes.

Caution: Prolonged immersion could cause a loss in cable flexibility.

3. Allow surfaces to air dry completely.

Autoclave Sterilization for all Probes

To sterilize a probe, perform the Philips Validated Autoclave Procedure below.

Note: These temperature probes are capable of withstanding up to 100 autoclave cycles.

Caution: To avoid damage to the probe, never exceed 138°C (280°F) maximum saturated steam temperature.

Autoclave Warning

Never exceed recommended maximum steam temperature of 138°C (280°F) during any autoclave cycle.

Philips Validated Autoclave Procedure

1. Coil the cleaned temperature probe assembly into an approximately 4-inch diameter coil.
2. Enclose the coiled probe in a breathable sterilization wrap.
3. Using a pre-vacuum autoclave expose the wrapped probe to a saturated steam temperature of 132°C (270°F) for 4 minutes and a dry time of 30 minutes.
4. Or a Saturated Steam temperature of 134°C/273°F (maximum) for 3 minutes with a dry time for 20 - 30 minutes.

Environmental Requirements

Below are the temperature, humidity, and atmospheric pressure ranges that must be maintained for these probes. Use, store, and transport probes within the following ranges:

Temperature Requirements

- Operating Temperature Range: 0°C to 45°C (32°F to 113°F)
- Storage Temperature Range: 0°C to 40°C (32°F to 104°F)
- Transport Temperature Range: -30°C to 60°C (-22°F to 140°F)

Humidity Requirements

- Operating Humidity Range: 5% RH to 95%RH
- Storage Humidity Range: 5% RH to 75% RH
- Transport Humidity Range: 15% RH to 85%RH

Atmospheric Pressure Requirements

- Operating Atmospheric Pressure Range: 540hPa to 1075hPa
- Storage Atmospheric Pressure Range: 540hPa to 1075hPa
- Transport Atmospheric Pressure Range: 540hPa to 1075hPa

Temperature Probe Disposal

Discard any temperature probe that becomes damaged or deteriorated. Follow approved medical waste disposal methods specified by your facility or local regulations.

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Copyright © 2019
Böblingen GmbH