
English

Introduction

Congratulations on your purchase and welcome to Philips Avent! To fully benefit from the support that Philips Avent offers, register your breast pump at www.philips.com/mybreastpump.

General description (Fig. 1)

1 Pump cover	5 Silicone diaphragm
2 Massage cushion	6 Handle
3 Pump body	7 Valve
4 Stem	8 Bottle (125 ml)

Intended use

The Philips Avent Breast Pump Manual Essential is intended to express and to collect milk from the breast of a lactating woman. The device is intended for a single user.

Important safety information

Read this user manual carefully before you use the breast pump and save it for future reference.

Contraindication

Never use the breast pump while you are pregnant, as pumping can induce labor.

WARNINGS!

Warnings to avoid choking, strangulation and injury:

- **Do not allow children or pets to play with the breast pump parts, packaging materials or accessories.**
- Inspect the breast pump for signs of damage before each use. Do not use the breast pump if it is damaged or if it does not work properly.

Warnings to avoid contamination and to ensure hygiene:

- For hygienic reasons, the breast pump is only intended for repeated use by a single user.
- Clean and disinfect all parts before first use and after each use.
- Do not use antibacterial or abrasive cleaning agents to clean the breast pump parts as this may cause damage.
- Wash your hands thoroughly with soap and water before you touch breast pump parts and breasts to prevent contamination. Avoid touching the inside of containers or lids.

Warnings to avoid breast and nipple problems including pain:

- Always release the handle before you remove the pump body from your breast in order to release the vacuum.
- If the vacuum is uncomfortable or causes pain, stop pumping and remove the breast pump from your breast.
- Do not attempt pumping for more than five minutes at a time if you do not succeed in expressing any milk. Try to express at another time during the day.
- If the process becomes very uncomfortable or painful, stop using the pump and consult your healthcare professional.
- Only use accessories and parts recommended by Philips Avent.
- No modification of the breast pump is allowed.

Warning: If you are infected with Human Immunodeficiency Virus (HIV) or Hepatitis B or C, using the breast pump will not reduce or remove the risk of transmitting the virus to your baby through breast milk.

Caution

Caution to prevent damage to and malfunctioning of the breast pump:

- Keep the parts away from heated surfaces to avoid deformation of these parts.

Preparing for use

Disassemble, clean and disinfect all parts before you use the breast pump for the first time and after every use.

Caution: Be careful when you remove and clean the white valve. Do not insert objects into the white valve, as this may cause damage.

Cleaning and disinfecting

Step 1: Disassembling

Disassemble the breast pump and bottle completely. Also remove the valve from the pump body and the stem from the silicone diaphragm.

Step 2: Cleaning

The parts can be cleaned manually or in the dishwasher.

Step 2A: Manual cleaning

- 1 Rinse all parts under a running tap with lukewarm water.
- 2 Clean all parts with a cleaning brush in warm water with some mild dishwashing liquid.
- 3 Thoroughly rinse all parts under a running tap with cold clear water.
- 4 Leave all parts to air-dry on a clean tea towel or drying rack.

Step 2B: Cleaning in the dishwasher

Note: Food colorings may discolor parts.

- 1 Place all parts on the top rack of the dishwasher.
- 2 Put dishwashing liquid or a tablet in the machine and run a standard program.
- 3 Leave all parts to air-dry on a clean tea towel or drying rack.

Step 3: Disinfecting

Caution: During disinfecting with boiling water, prevent the bottle or other parts from touching the side of the pot. This can cause irreversible product deformation or damage that Philips cannot be held liable for.

- 1 Fill a household pot with enough water to cover all parts and place all parts in the pot. Let the water boil for 5 minutes. Make sure the parts do not touch the side of the pot.
- 2 Allow the water to cool down.
- 3 Gently remove the parts from the water. Place the parts on a clean surface, a tea towel or drying rack to air-dry.

Tip: You can also disinfect the parts that come into contact with milk in a Philips Avent sterilizer.

Assembling the breast pump

Make sure you have cleaned and disinfected the parts of the breast pump that come into contact with milk.

Warning: Be careful, when you have disinfected the breast pump parts by boiling them, they can be very hot. To prevent burns, only start assembling the breast pump when disinfected parts have cooled down.

Check the breast pump kit parts for wear or damage before use and replace them if necessary.

Note: Correct placement of the massage cushion and the silicone diaphragm is essential for the breast pump to form a proper vacuum.

- 1 Insert the valve into the pump body from underneath. Push the valve in as far as possible (Fig. 2).
- 2 Screw the pump body clockwise onto the bottle (Fig. 3) until it is securely fixed.
- 3 Push the stem into the silicone diaphragm fully to the end (Fig. 4).
- 4 Insert the silicone diaphragm into the pump body from above. Make sure it fits securely around the rim by pressing down with your fingers to ensure a perfect seal (Fig. 5).
- 5 Attach the handle onto the diaphragm with stem by hooking the hole in the handle over the end of the stem. Push down the handle onto the pump body until it clicks into place (Fig. 6).
- 6 Place the massage cushion into the funnel-shaped section of the pump body. Push in the inner part of the massage cushion as far as possible and make sure that it is perfectly sealed all around the rim of the pump body. Press in between the petals to remove any trapped air (Fig. 7).

Note: Place the cover over the massage cushion to keep the breast pump clean while you prepare for expressing.

Using the breast pump

Warning: Always use the breast pump with the massage cushion.

Note: Contact your healthcare professional in case of questions or if you experience side effects.

Operating the breast pump

- 1 Wash your hands thoroughly with soap and water and make sure that your breast is centered, so that the massage cushion creates an airtight seal.
- 2 Gently start to press down the handle until you feel the suction of the pump (Fig. 8). Then allow the handle to return to its resting position. Repeat this step rapidly 5 or 6 times to initiate the "let-down" reflex.
- Note:** You do not have to press down the handle fully. Press it down as much as is comfortable for you. Your milk will soon start flowing, even though you are not using all the suction the pump can generate.
- 4 When the milk starts flowing, adopt a slower rhythm by pressing down the handle and keeping it pressed down for a moment before you let it return to its resting position. Continue with this rhythm while your milk is flowing. If your hand becomes tired, try to use the other hand or pump from the other breast.
- Note:** Do not worry if your milk does not flow immediately. Relax and continue pumping. Repositioning the breast pump on your breast from time to time can help stimulate milk flow.
- 5 Continue pumping until you feel your breast is empty.
- 6 When you have finished expressing, carefully remove the breast pump from your breast and unscrew the bottle from the pump body. Clean the other used parts of the breast pump (see "Cleaning and disinfecting").
- 7
- 8

After use

Warning: Only store breast milk collected with a clean and disinfected pump.

If you intend to feed your baby with the expressed breast milk within 4 days, you can store the breast milk in the refrigerator in an assembled Philips Avent bottle, storage cup, or milk bag. Otherwise, you can store breast milk in the freezer for up to 6 months.

Note: For safe use of bottles, storage cups and milk bags, please read the separately provided user manuals.

Warning: Refrigerate or freeze expressed breast milk immediately or keep it at a temperature of 16-29 °C (60-85 °F) for a maximum of 4 hours before you feed your baby. Never refreeze thawed breast milk to avoid a decrease in milk quality. Never add fresh breast milk to frozen breast milk to avoid a decrease in milk quality and unintended thawing of the frozen milk.

Note: For more information about the breast pump and tips for expressing milk, visit www.philips.com/avent.

Storage

Keep the breast pump out of direct sunlight as prolonged exposure may cause discoloration. Store the breast pump and its accessories in a safe, clean and dry place.

Compatibility

The Philips Avent Breast Pump Manual Essential is compatible with all baby bottles in the Philips Avent Natural range, Grow range (*) and Classic/ Anti-colic range, and Philips Avent storage cups with an adapter for cup. We do not recommend using the Philips Avent Natural glass bottles with the breast pump.

*only available in India

Technical information

- Operating conditions:
 - Temperature 5 °C to 40 °C (41 °F to 104 °F)
 - Relative humidity 15 to 90 % (non-condensing)
 - Pressure 700–1060 hPa (<2000 m/9842 ft. altitude)
- Storage and transportation conditions:
 - Temperature -25 °C to 70 °C (-13 °F to 158 °F)
 - Relative humidity up to 90 % (non-condensing)
- Materials: polypropylene, silicone

Troubleshooting

This chapter summarizes the most common problems you could encounter with the manual breast pump. If you are unable to solve the problem with the information below, visit www.philips.com/support for a list of frequently asked questions or contact the Consumer Care Center in your country.

Problem	Solution
I experience pain when I use the breast pump.	- Make sure that your nipple is placed in the center of the massage cushion. - Press the handle down to a vacuum level that is comfortable for you. - If the pain persists, stop using the breast pump and consult your healthcare professional.
I do not feel any suction.	- Make sure you have assembled the breast pump correctly. - See "Assembling the breast pump". - Make sure that the breast pump is positioned correctly on the breast to allow a seal between breast pump and breast.
The breast pump is damaged.	In case of regular use, some light scratching is normal and will not cause problems. However, if a part of the breast pump is badly scratched or cracked, stop using the breast pump and contact the Philips Consumer Care Center or visit www.philips.com/support to obtain a replacement part.

Supplemental information

When using the breast pump, some common breastfeeding-related conditions may occur, including pain sensation, sore nipples, engorgement (swelling of the breast), bruises, thrombus resulting into redness, swelling and/or pain, injured tissue on nipple (nipple trauma), or clogged mammary ducts. If you experience any of these symptoms, contact your healthcare professional.

Explanation of symbols

The warning signs and symbols are essential to ensure that you use this product safely and correctly and to protect you and others from injury.

This symbol indicates:

	to follow the instructions for use.
	important information such as warnings and cautions.
	usage tips, additional information or a note.
	the manufacturer.
	the batch number of the device.
	that the material used for the packaging is corrugated cardboard/fiberoast.

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