Amara View minimal contact full-face mask |

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1109350 R10
RWS 3/2/2018
© Koninklijke Philips N.V., 2018. All rights reserved.

Amara View minimal contact full-face mask

EN	001	INSTRUCTIONS FOR USE	LV	198	LIETOŠANAS INSTRUKCIJA
FR	014	MODE D'EMPLOI	LT	209	NAUDOJIMO INSTRUKCIJOS
IT	025	ISTRUZIONI PER L'USO	ET	220	KASUTUSJUHEND
DE	038	GEBRAUCHSANWEISUNG	RO	231	INSTRUCȚIUNI DE UTILIZARE
ES	050	INSTRUCCIONES DE USO	BG	242	ИНСТРУКЦИИ ЗА УПОТРЕБА
PT	062	INSTRUÇÕES DE UTILIZAÇÃO	UK	254	ІНСТРУКЦІЇ ІЗ ЗАСТОСУВАННЯ
EL	074	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	HR	266	UPUTE ZA UPORABU
DA	086	BRUGERVEJLEDNING	ZH-TW	277	使用說明
NL	097	GEBRUIKSAANWIJZING	ID	287	PETUNJUK PENGGUNAAN
SV	109	BRUKSANVISNING	MS	299	ARAHAN UNTUK PENGGUNAAN
FI	120	KÄYTTÖOHJEET	RU	310	ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ
NO	131	BRUKSANVISNING	KK	322	ПАЙДАЛАНУ НҮСҚАУЛАРЫ
PL	142	INSTRUKCJA OBSŁUGI	SR	335	UPUTSTVO ZA UPOTREBU
CS	154	NÁVOD K POUŽITÍ	HE	346	הוראות שימוש
HU	165	HASZNÁLATI UTASÍTÁS	TR	358	KULLANIM TALİMATLARI
SK	176	NÁVOD NA POUŽITIE	KO	369	사용 지침
SL	187	NAVODILA ZA UPORABO	AR	381	تعليمات للاستخدام

Amara View minimal contact full-face mask

Instructions for Use - English

Intended Use

This mask is intended to provide an interface for application of CPAP or bi-level therapy to patients. The mask is for single patient use in the home or multi-patient use in the hospital/institutional environment. The mask is to be used on patients (>66lbs/30kg) for whom CPAP or bi-level therapy has been prescribed.

Note:

- An exhalation port is built into this mask so a separate exhalation port is not required.
 - This mask is not made with natural rubber latex or DEHP.
-

Symbols

 Warning or Caution

 Note

 Tip

 Consult Instructions for Use

 Not Made with Natural Rubber Latex

X1 System One Resistance Control

Warnings

- This mask is not suitable for providing life support ventilation.
- This mask is designed for use with CPAP or bi-level systems recommended by your health care professional or respiratory therapist. Do not wear this mask unless the CPAP or bi-level system is turned on and operating properly. Do not block or try to seal the exhalation port. **Explanation of the Warning:** CPAP systems are intended to be used with special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. This warning applies to most models of CPAP systems.
- If oxygen is used with the device, the oxygen flow must be turned off when the device is not operating. **Explanation of the Warning:** When the device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the device enclosure. Oxygen accumulated in the device enclosure will create a risk of fire.

- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate. This warning applies to most types of CPAP and bi-level machines.
- At low CPAP or EPAP pressures the flow through the exhalation port may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- Some users may experience skin redness, irritation, or discomfort. If this happens, discontinue use and contact your healthcare professional.
- Consult a physician if any of the following symptoms occur: Unusual chest discomfort, shortness of breath, stomach distension, belching, severe headache, blurred vision, drying of the eyes, eye pain or eye infections.
- This mask should not be used on patients who are uncooperative, obtunded, unresponsive, or unable to remove the mask.
- This mask is not recommended if the patient is taking a prescription drug that may cause vomiting.
- Consult a physician or dentist if you encounter tooth, gum, or jaw soreness. Use of a mask may aggravate an existing dental condition.
- A minimum of 3cm H₂O (hPa) must be maintained when using this mask.
- Attaching an exhalation device requires therapy pressure level adjustment to compensate for increased leak.
- Do not overtighten the headgear straps. Watch for signs of overtightening, such as excessive redness, sores or bulging skin around the edges of the mask. Loosen the headgear straps to alleviate symptoms.
- Do not block or seal the anti-asphyxia valve or exhalation ports.

If your mask uses magnetic headgear clips, the following warnings apply:

- The mask assembly contains magnets. Some medical devices may be affected by magnetic fields. The magnetic clips in this mask should be kept at least 2 in. (50 mm) away from any active medical device with special attention to implanted devices such as pacemakers, defibrillators and cochlear implants.
- Do not use in or near magnetic resonance imaging (MRI) equipment.
- Keep unassembled magnetic headgear clips out of reach of children. In case of accidental swallowing, seek medical assistance immediately.



Figure ①

- ① Mask Cushion with Exhalation Ports
(Do not block)
- ② Mask Frame
- ③ Elbow with Anti-Asphyxia Valve
(Do not block)
- ④ Quick Release Tabs
- ⑤ Quick Release Tube with swivel
- ⑥ Headgear Crown Strap
- ⑦ Headgear Top Side Strap
- ⑧ Headgear Bottom Side Strap
- ⑨ Mask Frame Top Headgear Slot
- ⑩ Magnetic Headgear Clip Receiver
- ⑪ Magnetic Headgear Clip


Optional:

- ⑫ Talon Headgear Clip

Verify the Anti-Asphyxia Valve

The anti-asphyxia valve consists of an air inlet and a flapper ②. With the airflow turned off, verify that the elbow with the flapper is lying flat **a** so that room air can flow in and out through the air inlet. Next, with the airflow on, the flapper should now cover the air inlet and air from the CPAP or bi-level device should flow into the mask **b**. If the flapper does not close or does not function properly, replace the mask.

⚠Warning: Do not block or seal the anti-asphyxia valve.

Contraindications

This mask may not be suitable for use on patients with the following conditions: recent eye surgery or dry eyes, hiatal hernia, excessive reflux, impaired cough reflex, impaired cardiac sphincter function; or on patients unable to remove the mask by themselves.

Before Use Read and Understand the Instructions Completely.

- Hand wash the mask.
- Wash your face.

⚠ Caution: Do not use moisturizer/lotion on your hands or face.

- Inspect the mask and replace if the cushion has hardened or is torn, or if any parts are broken.
- Verify that the therapy device, i.e., ventilator, including the alarms and safety systems, has been validated prior to use.
- Verify therapy device pressure(s).

Cleaning Instructions

Hand wash the mask before first use and daily. The headgear should be hand washed weekly. The headgear does not need to be removed for daily cleaning of the mask.

1. Hand wash mask and headgear in warm water with liquid dishwashing detergent.
2. Rinse thoroughly. Air dry completely before use. Make sure the mask is dry before use. Lay the headgear flat or line dry.

⚠Caution: Do not place the headgear into the dryer.

⚠Caution: Do not use bleach, alcohol, cleaning solutions containing bleach or alcohol, or cleaning solutions containing conditioners or moisturizers.

⚠Caution: Any deviation from these instructions may damage the product.

⚠Caution: Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.

Dishwashing Instructions

In addition to hand washing, the mask may be cleaned in the dishwasher once a week.

⚠Caution: Use a mild liquid dish washing detergent only to wash the mask.

1. Remove the headgear. Do not wash the headgear in the dishwasher.
2. Separate the cushion and the faceplate and wash in the top shelf of the dishwasher.
3. Air dry. Make sure the mask is dry before use.

⚠Caution: Inspect the mask for damage or wear (cracking, crazing, tears, etc). Discard and replace any components as necessary.


Institutional Disinfection

Clean and disinfect the mask daily and between patients when in use in a clinical setting. Refer to the Disinfection Guide for Professional Users to reprocess between patients in a clinical setting. Access the latest version of the Disinfection Guide at www.healthcare.philips.com or by contacting customer service at 1-800-345-6443 (USA or Canada only) or +1-724-387-4000.

Achieving the Right Fit

Sizing the Mask



The mask should rest comfortably against your face. The cushion nostril opening goes under your nostrils (3). The headgear should be away from your eyes. The bottom of the mask cushion should rest just above your chin with your mouth slightly open.

Before putting on the mask:

1. Headgear - Loosen to a large setting.
2. Magnetic Headgear Clips - Grip and twist outward to disconnect from the magnetic headgear clip receiver (4a).

Optional: If your mask uses talon headgear clips, grip and twist away from the frame to disconnect. (4b)

Adjusting the mask:

1. Hold the mask against your face. Pull the headgear over your head (5).
2. Mask cushion - Place your nose over the cushion nostril opening. Push the mask cushion up.

Push the mask cushion up.


Note: Do not place your nose in the cushion nostril opening.

3. Magnetic Headgear Clips - Connect the magnetic headgear clip to the magnetic headgear clip receiver. (4a)

Optional: if your mask uses talon headgear clips, press into place to connect. (4b)

4. Headgear - Pull back the tabs to adjust the top and bottom side straps. The headgear should fit loosely against your face. Do not overtighten (6)(7).

5. Headgear - Pull the crown strap tab over to reduce leaks around your nose (8).

6. Headgear - Make final adjustments by pulling back the tabs to evenly adjust the top and bottom side straps. The headgear should feel snug, but not uncomfortably tight (6)(7).

Note: You should be able to fit one of your fingers underneath the top and bottom side straps.

Using the mask:

1. Connect the flexible tubing (included with your CPAP or bi-level device) to the quick release tube.

2. Lie down. Turn the therapy device on. Breathe normally. Keep your mouth slightly open.

3. Assume different sleeping positions. Move around until comfortable. If there are any excessive air leaks, make final adjustments. If there are any excessive air leaks around the nose, adjust the crown strap first. This will move the mask closer or further away from the bottom of your nose. Some air leaking is normal.

✱ Comfort Tips

- The most common mistake is overtightening the headgear. The headgear should fit loose and comfortable. If your skin bulges around the mask or if you see red marks on your face, loosen the headgear.
- Adjust the crown headgear strap to reduce leaks around your nose and to lift the straps off of your ears (8).
- Adjust the top and bottom side straps to reduce leaks at the sides or bottom of the mask cushion (6)(7).

Removing the Mask

Disconnect the magnetic headgear clips and pull over the top of your head.

Optional: If your mask uses talon headgear clips, disconnect the headgear clips and pull over the top of your head.

Disassembly


1. Magnetic Headgear Clips - Grip and twist away from the magnetic headgear clip receiver to disconnect.

Optional: If your mask uses talon headgear clips, grip and twist away from the frame to disconnect.


2. Headgear - Undo the straps. Pull out of headgear clips and mask frame top headgear slots.

3. Mask Cushion - Grasp the top and the bottom of the mask cushion, and pull it off the mask frame (9).

⚠ Caution: Do not pull the mask cushion by the nostril cushion opening.

4. Quick Release Tube - Squeeze the quick release tabs. Pull the tube away from the mask frame (10).

Assembly

1. Headgear - Thread the headgear straps through the mask frame top headgear slots and magnetic headgear clips.
2. Magnetic Headgear Clips - Connect the headgear clip to the magnetic headgear clip receiver.

Optional: If your mask uses talon headgear clips, push the headgear clips onto the mask frame.

3. Mask Cushion - Grasp both sides and push it onto the mask frame.
4. Quick Release Tube - Push onto the elbow until the quick release tabs click into place.

Philips Respironics System One Resistance Control


Your mask when combined with a Philips Respironics System One device, provides optimal resistance compensation. This mask value is X1 and should be set by your provider. Notes:

- Compare the mask to the device. See your device manual if the values do not match.
- System One is not compatible with masks that use a separate exhalation device.

Specifications

⚠ Warning: The technical specifications of the mask are provided for your healthcare professional to determine if it is compatible with your CPAP or bi-level therapy device. If used outside these specifications, or if used with incompatible devices, the mask may be uncomfortable, the seal of the mask may not be effective, optimum therapy may not be achieved, disconnect alarms may not be detected, and leak, or variation in the rate of leak, may affect device function. The pressure flow curve shown below is an approximation of expected performance. Exact measurements may vary.

Pressure Flow Curve

Resistance with Anti-Asphyxia Valve Closed to Atmosphere Drop in Pressure at

	50 SLPM	100 SLPM
All sizes	0.8 cm H ₂ O	2.0 cm H ₂ O

Closed to Atmosphere: A minimum of 3 cm H₂O must be maintained to ensure that the anti-asphyxia valve closes.

Open to Atmosphere: In the absence of positive airway pressure, the anti-asphyxia valve will open.

Inspiratory and Expiratory Resistance with Anti-Asphyxia Valve Open to Atmosphere at 50 SLPM

Inspiratory Resistance	1.0 cm H ₂ O
Expiratory Resistance	1.5 cm H ₂ O

Deadspace

S	150.4 mL
M	148.8 mL
L	167.8 mL

Sound Levels

A-weighted Sound Power Level	32.5 dBA
A-weighted Sound Pressure Level at 1m	24.2 dBA

Disposal

Dispose of in accordance with local regulations.

Dispose of mask and headgear annually.

Storage Conditions

Temperature: -4° F to 140° F (-20° C to 60° C)

Relative Humidity: 15% to 95%, non-condensing

Limited Warranty

Respironics, Inc., a Philips company, warrants that the system shall be free from defects of workmanship and materials and will perform in accordance with the product specifications for a period of ninety (90) days from the date of sale by Respironics, Inc. to the dealer. If the product fails to perform in accordance with the product specifications, Respironics, Inc. will repair or replace – at its option – the defective material or part. Respironics, Inc. will pay customary freight charges from Respironics, Inc. to the dealer location only. This warranty does not cover damage caused by accident, misuse, abuse, alteration, water ingress, and other defects not related to material or workmanship.. The Respironics, Inc. Service department shall examine any devices returned for service, and Respironics, Inc. reserves the right to charge an evaluation fee for any returned device as to which no problem is found after investigation by Respironics, Inc. Service.

This warranty is non-transferable by unauthorized distributors of Respironics, Inc. products and Respironics, Inc. reserves the right to charge dealers for warranty service of failed product not purchased directly from Respironics, Inc. or authorized distributors.

Respironics, Inc. disclaims all liability for economic loss, loss of profits, overhead, or consequential damages which may be claimed to arise from any sale or use of this product. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

This warranty is given in lieu of all other express warranties. In addition, any implied warranties – including any warranty of merchantability or fitness for the particular purpose – are limited to two years. Some states do not allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you. This warranty gives you specific legal rights, and you may also have other rights which vary from state to state. To exercise your rights under this warranty, contact your local authorized Respironics, Inc. dealer or contact Respironics, Inc. at:

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Limited Warranty For Australia and New Zealand

Note: For Australian and New Zealand customers this warranty replaces the warranty contained above.

1. Respironics, Inc., a Philips company ("Philips Respironics"), warrants that the Products shall be free from defects of workmanship and materials and will perform in accordance with the Product specifications. 2. This warranty is valid for a period of 30 days from the date of purchase from an authorized Philips Respironics dealer. 3 If the Product is found to contain a defect of workmanship or materials or fails to perform in accordance with the Product specifications during the applicable warranty period, Philips Respironics will repair or replace, at its option, the defective material or part. 4. The customer is responsible for returning the product to an authorised Philips Respironics dealer, and collecting the product from the authorised Philips Respironics dealer after repair or replacement, at its own cost. Philips Respironics is responsible only for the freight cost of transporting the product between the authorised Philips Respironics dealer and Philips Respironics. Philips Respironics reserves the right to charge an evaluation and postage fee for any

returned Product as to which no problem is found following investigation. 5. This warranty does not cover damage caused by accident, misuse, abuse, alteration, and other defects not related to materials or workmanship. 6. The warranty provided by Philips Respironics herein is not transferrable by the Buyer in the event of any sale or transfer of Products purchased by the Buyer from an authorized Philips Respironics dealer. 7. To the extent permitted by law, where the Buyer has the benefit of an implied guarantee under the Australian Consumer Law, but the Product is not of a kind ordinarily acquired for personal, domestic or household use or consumption Philips Respironics' liability shall be limited, at the option of Philips Respironics, to the replacement or repair of the Product or the supply of an equivalent Product. 8. To exercise your rights under this warranty, contact your local authorised Philips Respironics dealer. A list of all authorised dealers is available at the following link: <http://www.philips.com.au/healthcare>.

Alternatively, you can make a claim under this warranty by contacting Philips Respironics directly at: Philips Electronics Australia Limited, 65 Epping Road, North Ryde NSW 2113, Australia. Tel: 1300 766 488, Email: prcontact@philips.com. 9. The following statement is provided to a Buyer who is a "consumer" under the Australian Consumer Law: Our goods come with guarantees that cannot be excluded under the Australian Consumer Law. You are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the good repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major failure. 10. The following statement is provided to a Buyer who is a "consumer" under the Consumer Guarantees Act 1993, New Zealand: Our goods come with guarantees that cannot be excluded under the Consumer Guarantees Act 1993. This guarantee applies in addition to the conditions and guarantees implied by that legislation.

AUSTRALIAN SPONSOR DETAILS:

Philips Electronics Australia Ltd.
65 Epping Road, North Ryde, NSW 2113
Australia

Amara View minimal contact full face mask

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Respironics Deutschland GmbH & Co. KG