Autoclavable Temperature Probes

Instructions for Use

Probes with 2-pin connectors

REF 21075A

12 Fr Esophageal/Rectal Probe

REF 21076A

10 Fr Esophageal/Rectal Probe

REF 21078A

Skin Probe

Probes with 1/4 inch phone-plug connectors

REF 989803162621

12 Fr Esophageal/Rectal Probe

REF 989803162631

10 Fr Esophageal/Rectal Probe

REF 989803162641

Skin Probe

PHILIPS

About This Edition

Notice

Copyright © 2009

Koninklijke Philips Electronics N.V.
All rights are reserved.

Printing History

Current Updated Edition/Revision:

March 2009
Part Number 453564146211-L-IFU Rev A

Previous Edition/Revision:

None

Contents

Instructions for Use.....	1	Návod na použitie	71
Manuel d'utilisation	6	Navodila za uporabo	76
Gebrauchsanweisung	11	Instrucțiuni de folosire	81
Gebruiksaanwijzing	16	Указания за употреба.....	86
Instrucciones de uso.....	21	Kasutusjuhend.....	91
Instruções de Uso	26	Naudojimo instrukcija.....	96
Istruzioni d'uso	31	Lietošanas instrukcija.....	101
Brugerhåndbok.....	36	Инструкция по эксплуатации....	106
Bruksanvisning	41	Οδηγίες Χρήσης.....	111
Käyttöopas	46	Kullanım Talimatları.....	116
Brugsanvisning	51	ユ－ザ－ズ・ガイド.....	121
Használati útmutató	56	使用説明.....	126
Instrukcja obsługi.....	61	使用说明书.....	131
Návod k obsluze.....	66		

Product Overview

Intended Use

Philips autoclavable temperature probes are reusable devices used to continuously monitor patient body temperature during clinical care. The probes are available in several configurations that can be applied to the patient orally, rectally, or attached to the skin. (Note: core temperature readings are normally different than skin temperature readings in the same patient.)

These probes should be applied to a patient only under direct supervision of a licensed physician or healthcare provider, and connected only to compatible patient monitors (see *Compatibility* below). The probes are shipped non-sterile. It is recommended that each probe be cleaned, disinfected and inspected for damage before use.

Measurement Range and Accuracy

Temperature accuracy of these probes at the site of measurement is:

- $\pm 0.2^{\circ}\text{C}$ within the temperature range of 0°C to 60°C
- $\pm 0.1^{\circ}\text{C}$ within the temperature range of 25°C to 45°C .

Compatibility

Six temperature probe models are available (see table below). All are compatible with patient monitors that use 400-Series technology.

Autoclavable Temperature Probes

Probe Model No.	Patient Application	Length	Probe Connector
21075A	12 Fr Esophageal/Rectal	10 ft	 <i>2-pin plugs into Philips monitors</i>
21076A	10 Fr Esophageal/Rectal	10 ft	
21078A	Skin	10 ft	
989803162621	12 Fr Esophageal/Rectal	10 ft	 <i>1/4-inch phone plug plugs into non-Philips monitors</i>
989803162631	10 Fr Esophageal/Rectal	10 ft	
989803162641	Skin	10 ft	

Optional Connector Cables

Two optional reusable connector cables, orderable separately, offer adaptability between the temperature probes and 400-Series patient monitors (see table below):

Optional Connector Cables

Optional Cable Model No.	Cable Length	Temperature Probe to Instrument Adaptability
989803163271	10 ft	Adapts a 2-pin probe to a 1/4-inch phone plug patient monitor.
989803163281	10 ft	Adapts a 1/4-inch phone plug probe to a 2-pin patient monitor.

Warning: DO NOT autoclave these optional connector cables

Definition of Product Symbols

Symbols identified below appear on probe product packaging.

R_x only Prescription use only	 Refer to instruction manual/booklet	 Contains No Latex	 Contains No PVC	 Defibrillator Proof
----------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------	----------------------------------------------------------------------------------------------------------

Probe Disposal

Discard any temperature probe or optional connector cable that becomes damaged or deteriorated. Follow approved medical waste disposal methods specified by your facility or local regulations.

Warnings

Be sure to read and understand all warnings listed below.

Patient Care Warnings

- Use temperature probes only under direct supervision of a licensed and trained healthcare professional.
- Before use, clean and disinfect temperature probes and optional connector cables.
- Before use, inspect and replace damaged or deteriorated temperature probes and optional connector cables.
- When using these probes, follow standard application practices recommended by your medical facility.
- Connect temperature probes only to monitors that use 400-Series probe technology.
- To avoid patient injury:
 1. Do NOT apply these temperature probes to patients who are undergoing Magnetic Resonance Imaging (MRI) procedures.
 2. Always use caution when applying, inserting, or removing a temperature probe from a patient.
 3. Assure that probe and optional connector cable are not positioned where they could become entangled around patient, or cause choking, strangulation, or inhibit circulation in extremities.
 4. Assure that probe and optional connector cable are not positioned where they could become entangled around IV tubes and restrict IV flow.

RF Energy and Electrosurgical (ESU) Warnings

All wire-lead patient-connected sensor assemblies (including temperature probes) are subject to reading error, local heating and possible damage from high intensity sources of RF energy. Electrosurgical (ESU) equipment represents one such source. Capacitively coupled currents may seek

alternate paths to ground through probe cables and associated instruments. Patient burns may result.

If at all possible, remove the probe from patient contact before activating the surgical unit or other RF source. If probes must be used simultaneously with electrosurgical apparatus, the instruments attached to the probes should be checked for adequate isolation from RF grounds. Hazards can be reduced by selecting a temperature monitoring point located away from the expected RF current path between the active probe and the return pad.

Cleaning and Disinfection

Cleaning and disinfection of a probe, and optional connector cable if used, is recommended after every patient procedure. Proceed as follows:

1. Clean with a soft cloth moistened with 70% isopropanol, activated dialdehyde (Cidex) or sodium hypochlorite (1 part bleach diluted with 10 parts water minimum).
2. Wipe from one end to the other, being careful not to apply excessive pressure that could stretch the cable jacket and break the internal wires.
3. Disinfect by briefly immersing the probe and optional connector cable into the same cleaning solution, being careful to *avoid immersing the connectors*. Immersion of the cable for 2 minutes or less is not harmful, however prolonged immersion could cause a loss in cable flexibility.
4. Rinse with clean water and dry thoroughly.

Note: Philips Medical Systems does not make any claims as to the efficacy of these chemicals for infection control. Please consult your medical facility's Infection Control Officer for the applicable disinfection policies.

Autoclave Sterilization

Cleaned and disinfected temperature probes are capable of being sterilized prior to use. To sterilize a probe, perform either the *Philips Validated Autoclave Procedure*, or the *Alternative Autoclave Procedure*, as outlined below.

Autoclave Warnings

- Never autoclave optional connector cables.
- Never autoclave any gray color Philips temperature probe. Only tan color Philips probes are designed to withstand steam autoclaving.
- Never exceed maximum steam temperature recommendations of 138°C during any autoclave cycle.

Philips Validated Autoclave Procedure

1. Coil the cleaned and disinfected temperature probe assembly into a four-inch diameter coil.
2. Enclose the coiled probe in a breathable sterilization wrap.
3. Using a pre-vacuum single pulse process, achieve 27-inches Hg vacuum.
4. Expose the wrapped probe to a saturated steam temperature of 132°C (minimum) - 138°C (maximum) for four minutes.

Alternative Autoclave Procedure

Alternatively, you may prepare and expose the temperature probe to the autoclave sterilization cycle recommended by your medical facility and autoclave equipment manufacturer. To avoid damage to the probe, never exceed 138°C maximum saturated steam temperature.

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For more information, please call the local Philips sales office listed in your telephone directory; or contact one of the regional sales offices listed below for the location of your nearest Philips sales office.

United States:

Philips Medical Systems
North America Corporation
3000 Minuteman Road
Andover, MA 01810-1099
+1 800 225-0230 toll free

Canada:

Philips Medical Systems Canada
281 Hillmount Road
Markham, ON
L6C 2S3
(800) 291-6743

Latin America Headquarters:

Philips Argentina S.A.
Depto Sistemas Medicos, Vedia 3892,
Casilla
3479, 1429 Buenos Aires
Tel: +54 11 4546 7699
Fax: +54 11 4546 7697

0123 This product complies with the requirements of the Council Directive 93/42/EEC June 1993 (Medical Device Directive).

Made in USA

US Status: Distributed by
Philips Medical Systems

Manufacturing Address:
Philips Medical Systems
3000 Minuteman Road
Andover, MA 01810-1099 USA
Tel: (800) 225-0230

Europe, Middle East and Africa:

Philips Medizin Systeme
Böblingen GmbH
Hewlett-Packard Str. 2
71034 Böblingen, Germany
Tel: (+49) 7031 463 2254
Fax: (+49) 7031 463 1552

Asia Pacific Headquarters:

Philips Medical Systems
30/F Hopewell Centre
17 Kennedy Road
Wanchai
Hong Kong
Tel: (852) 2821 5888
Fax: (852) 2527 6727

Contains no Latex

<http://philips.com/healthcarestore>

Subject to modification

© 2009
Koninklijke Philips Electronics N.V.
Printed in USA March 2009

Part No. 453564146211