

# DICOM Conformance Statement

Application Annex:  
RTgo



**Issued by:**

Philips Medical Systems Netherlands BV, a Philips Healthcare company,

P.O. Box 10.000  
5680 DA Best  
The Netherlands

Email: [dicom@philips.com](mailto:dicom@philips.com)  
Internet: [www.philips.com/dicom](http://www.philips.com/dicom)

Document Number: ICAP-PF.0028233

Date: 12-Apr-2017

## Table of Contents

|   |          |
|---|----------|
| <b>1. INTRODUCTION</b> .....                          | <b>4</b> |
| <b>1.1. REVISION HISTORY</b> .....                    | <b>4</b> |
| <b>1.2. TERMINOLOGY</b> .....                         | <b>4</b> |
| <b>2. DATA SPECIFICATIONS</b> .....                   | <b>5</b> |
| <b>2.1. SUPPORTED IOD'S</b> .....                     | <b>5</b> |
| 2.1.1. Acceptance Criteria .....                      | 5        |
| 2.1.2. Contents of Created IOD's .....                | 5        |
| 2.1.2.1. Implementation Identifying Information ..... | 6        |
| 2.1.2.2. CT Image Storage SOP Class .....             | 6        |
| 2.1.2.3. RT Structure Set Storage SOP Class .....     | 9        |

## 1. Introduction

This DICOM Conformance Statement annex is applicable to the RTgo plug-in that can be installed on top of MR software.

MR RT Plugin SW is interfaced with MR SW. This Plugin will include the MRCAT feature (previously this was part of MR SW) and new features for MR RT product, AutoContouring and Sim-to-treat SW options. The MRCAT images are transferred as CT objects and autocontouring is transferred as DICOM RT structures. This data is used by treatment planning systems.

The plug-in enables following applications,

- MR-only sim
  - o MR images can be processed into CT-like MRCAT (MR for Calculating Attenuation) images, to provide the information that is needed for assigning Hounsfield units to various tissue types, facilitating MR-only treatment planning for radiation therapy
  - o Images can be exported as CT images
- Auto-Contouring
  - o A set of anatomical structures can be automatically generated based on MR images.
  - o Contours can be exported as RT structure sets

### 1.1. Revision History

The revision history below provides dates and differences among individual document versions.

**Table 1: Revision History**

| Document Version | Date of Issue | Status     | Description    |
|------------------|---------------|------------|----------------|
| 00               | 12-April-2017 | Authorized | Final version. |

### 1.2. Terminology

|       |  |
|-------|--|
| DICOM | Digital Imaging and Communications in Medicine |
| IOD   | Information Object Definition                  |
| UID   | Unique Identifier                              |
| VR    | Value Representation                           |
| RT    | Radio Therapy                                  |



The abbreviations used in the Module table for the column "Source" are:

|          |   |
|----------|---|
| AUTO     | The attribute value is generated automatically                                    |
| CONFIG   | The attribute value source is a configurable parameter                            |
| COPY     | The attribute value source is another SOP instance                                |
| FIXED    | The attribute value is hard-coded in the application                              |
| IMPLICIT | The attribute value source is a user-implicit setting                             |
| MPPS     | The attribute value is the same as that use for Modality Performed Procedure Step |
| MWL      | The attribute value source is a Modality Worklist                                 |
| USER     | The attribute value source is explicit user input                                 |

### 2.1.2.1. Implementation Identifying Information

The value supplied for Implementation Class UID and version name are specified in table 4.

**Table 4: DICOM Implementation Class and Version for Created Instances**

|                             |                         |
|-----------------------------|-------------------------|
| Implementation Class UID    | 1.3.46.670589.54.2.16.4 |
| Implementation Version Name | 16.4.0.0                |

Note: The Implementation Class and Implementation Version for RTgo Plug-in DICOM export refers to different DICOM library and are therefore has different identification compared to the standard system.

**Table 5: List of created SOP Classes**

| SOP Class Name                     | SOP Class UID                 |
|------------------------------------|-------------------------------|
| CT Image Storage                   | 1.2.840.10008.5.1.4.1.1.2     |
| RT Structure Set Storage SOP Class | 1.2.840.10008.5.1.4.1.1.481.3 |

### 2.1.2.2. CT Image Storage SOP Class

**Table 6: IOD of Created CT Image Storage SOP Class Instances**

| Information Entity | Module                                | Presence Of Module |
|--------------------|---------------------------------------|--------------------|
| Patient            | Patient Module                        | ALWAYS             |
| Study              | General Study Module                  | ALWAYS             |
| Study              | Patient Study Module                  | CONDITIONAL        |
| Series             | General Series Module                 | ALWAYS             |
| Frame of Reference | Frame of Reference Module             | ALWAYS             |
| Equipment          | General Equipment Module              | ALWAYS             |
| Image              | General Image Module                  | ALWAYS             |
| Image              | Image Plane Module                    | ALWAYS             |
| Image              | Image Pixel Module                    | ALWAYS             |
| Image              | CT Image Module                       | ALWAYS             |
| Image              | VOI LUT Module                        | ALWAYS             |
| Image              | SOP Common Module                     | ALWAYS             |
|                    | Extended DICOM and Private attributes | ALWAYS             |

**Table 7: Patient Module**

| Attribute Name       | Tag       | VR | Value | Presence of Value | Source    | Comment |
|----------------------|-----------|----|-------|-------------------|-----------|---------|
| Patient's Name       | 0010,0010 | PN |       | ALWAYS            | MWL, USER | -       |
| Patient ID           | 0010,0020 | LO |       | ALWAYS            | MWL, USER | -       |
| Patient's Birth Date | 0010,0030 | DA |       | ALWAYS            | MWL, USER | -       |

|                   |           |    |  |        |           |   |
|-------------------|-----------|----|--|--------|-----------|---|
| Patient's Sex     | 0010,0040 | CS |  | ALWAYS | MWL, USER | - |
| Other Patient IDs | 0010,1000 | LO |  | ANAP   | MWL, USER |   |

Table 8: General Study Module

| Attribute Name             | Tag       | VR | Value | Presence of Value | Source          | Comment |
|----------------------------|-----------|----|-------|-------------------|-----------------|---------|
| Study Date                 | 0008,0020 | DA |       | ALWAYS            | AUTO, MWL       | -       |
| Study Time                 | 0008,0030 | TM |       | ALWAYS            | AUTO, MWL       | -       |
| Accession Number           | 0008,0050 | SH |       | ALWAYS            | AUTO, MWL, USER | -       |
| Referring Physician's Name | 0008,0090 | PN |       | VNAP              | MWL, USER       | -       |
| Study Description          | 0008,1030 | LO |       | VNAP              | MWL, USER       | -       |
| Study Instance UID         | 0020,000D | UI |       | ALWAYS            | AUTO, MWL       | -       |
| Study ID                   | 0020,0010 | SH |       | ALWAYS            | AUTO            | -       |

Table 9: Patient Study Module

| Attribute Name   | Tag       | VR | Value | Presence of Value | Source    | Comment |
|------------------|-----------|----|-------|-------------------|-----------|---------|
| Patient's Weight | 0010,1030 | DS |       | ALWAYS            | MWL, USER | -       |
| Medical Alerts   | 0010,2000 | LO |       | ANAP              | MWL, USER | -       |
| Allergies        | 0010,2110 | LO |       | ANAP              | MWL, USER | -       |
| Pregnancy Status | 0010,21C0 | US |       | ANAP              | MWL, USER | -       |

Table 10: General Series Module

| Attribute Name              | Tag       | VR | Value | Presence of Value | Source     | Comment                             |
|-----------------------------|-----------|----|-------|-------------------|------------|-------------------------------------|
| Series Date                 | 0008,0021 | DA |       | ALWAYS            | AUTO       | Copy from MRCAT mDixon source scan. |
| Series Time                 | 0008,0031 | TM |       | ALWAYS            | AUTO       | Copy from MRCAT mDixon source scan. |
| Modality                    | 0008,0060 | CS | CT    | ALWAYS            | FIXED      |                                     |
| Series Description          | 0008,103E | LO |       | ANAP              | AUTO, USER |                                     |
| Performing Physicians' Name | 0008,1050 | PN |       | ANAP              | AUTO, USER |                                     |
| Body Part Examined          | 0018,0015 | CS |       | ANAP              | AUTO       |                                     |
| Protocol Name               | 0018,1030 | LO |       | ALWAYS            | USER       |                                     |
| Patient Position            | 0018,5100 | CS |       | ALWAYS            | AUTO       |                                     |
| Series Instance UID         | 0020,000E | UI |       | ALWAYS            | AUTO       |                                     |
| Series Number               | 0020,0011 | IS |       | ALWAYS            | AUTO       |                                     |
| Laterality                  | 0020,0060 | CS |       | ANAP              | USER       |                                     |

Table 11: Frame of Reference Module

| Attribute Name               | Tag       | VR | Value | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-------|-------------------|--------|---------|
| Frame of Reference UID       | 0020,0052 | UI |       | ALWAYS            | AUTO   | -       |
| Position Reference Indicator | 0020,1040 | LO |       | EMPTY             | FIXED  | -       |

Table 12: General Equipment Module

| Attribute Name                | Tag       | VR | Value | Presence of Value | Source | Comment                   |
|-------------------------------|-----------|----|-------|-------------------|--------|---------------------------|
| Manufacturer                  | 0008,0070 | LO |       | ALWAYS            | FIXED  |                           |
| Institution Name              | 0008,0080 | LO |       | ALWAYS            | CONFIG | Configured on the system. |
| Institution Address           | 0008,0081 | ST |       | ANAP              |        |                           |
| Station Name                  | 0008,1010 | SH |       | ALWAYS            | CONFIG | Same as the Host Name.    |
| Institutional Department Name | 0008,1040 | LO |       | ALWAYS            | CONFIG | -                         |

|                           |           |    |  |        |              |   |
|---------------------------|-----------|----|--|--------|--------------|---|
| Manufacturer's Model Name | 0008,1090 | LO |  | ALWAYS | FIXED        |   |
| Device Serial Number      | 0018,1000 | LO |  | ALWAYS | AUTO         | System serial number.                   |
| Software Version(s)       | 0018,1020 | LO |  | ALWAYS | AUTO or COPY | The release text of the original Image. |

Table 13: General Image Module

| Attribute Name     | Tag       | VR | Value   | Presence of Value | Source | Comment |
|--------------------|-----------|----|---|-------------------|--------|---------|
| Image Type         | 0008,0008 | CS | DERIVED\PRIMARY\REFORMATTED\REFORMATTED\DERIVED | ANAP              | AUTO   | -       |
| Acquisition Date   | 0008,0022 | DA |   | ANAP              | AUTO   | -       |
| Content Date       | 0008,0023 | DA |   | ANAP              | AUTO   | -       |
| Acquisition Time   | 0008,0032 | TM |   | ANAP              | AUTO   | -       |
| Content Time       | 0008,0033 | TM |   | ANAP              | AUTO   | -       |
| Acquisition Number | 0020,0012 | IS |   | ANAP              | AUTO   |         |
| Instance Number    | 0020,0013 | IS |   | VNAP              | AUTO   |         |

Table 14: Image Plane Module

| Attribute Name              | Tag       | VR | Value | Presence of Value | Source | Comment |
|-----------------------------|-----------|----|-------|-------------------|--------|---------|
| Slice Thickness             | 0018,0050 | DS |       | VNAP              | AUTO   | -       |
| Image Position (Patient)    | 0020,0032 | DS |       | ALWAYS            | AUTO   | -       |
| Image Orientation (Patient) | 0020,0037 | DS |       | ALWAYS            | AUTO   | -       |
| Pixel Spacing               | 0028,0030 | DS |       | ALWAYS            | AUTO   | -       |

Table 15: Image Pixel Module

| Attribute Name             | Tag       | VR       | Value | Presence of Value | Source   | Comment |
|----------------------------|-----------|----------|-------|-------------------|----------|---------|
| Samples per Pixel          | 0028,0002 | US       | 1     | ALWAYS            | FIXED    |         |
| Photometric Interpretation | 0028,0004 | CS       |       | ALWAYS            | FIXED    |         |
| Rows                       | 0028,0010 | US       |       | ALWAYS            | IMPLICIT |         |
| Columns                    | 0028,0011 | US       |       | ALWAYS            | IMPLICIT |         |
| Bits Allocated             | 0028,0100 | US       | 16    | ALWAYS            | FIXED    |         |
| Bits Stored                | 0028,0101 | US       | 12    | ALWAYS            | IMPLICIT |         |
| High Bit                   | 0028,0102 | US       | 11    | ALWAYS            | IMPLICIT |         |
| Pixel Representation       | 0028,0103 | US       | 0     | ALWAYS            | IMPLICIT |         |
| Pixel Data                 | 7FE0,0010 | OB or OW |       | ANAP              | IMPLICIT |         |

Table 16: CT Image Module

| Attribute Name | Tag       | VR | Value   | Presence of Value | Source | Comment                         |
|----------------|-----------|----|---|-------------------|--------|---------------------------------|
| Image Type     | 0008,0008 | CS | DERIVED\PRIMARY\REFORMATTED\REFORMATTED\DERIVED | ALWAYS            | AUTO   |                                 |
| KVP            | 0018,0060 | DS | 0   | VNAP              | FIXED  | Derived from MR image so no KVP |



|                            |           |    |                 |        |       |  |
|----------------------------|-----------|----|-----------------|--------|-------|--|
| Acquisition Number         | 0020,0012 | IS |                 | VNAP   | AUTO  |  |
| Samples per Pixel          | 0028,0002 | US |                 | ALWAYS | AUTO  |  |
| Photometric Interpretation | 0028,0004 | CS | MONOC<br>HROME2 | ALWAYS | AUTO  |  |
| Bits Allocated             | 0028,0100 | US |                 | ALWAYS | AUTO  |  |
| Bits Stored                | 0028,0101 | US |                 | ALWAYS | AUTO  |  |
| High Bit                   | 0028,0102 | US |                 | ALWAYS | AUTO  |  |
| Rescale Intercept          | 0028,1052 | DS |                 | ALWAYS | AUTO  |  |
| Rescale Slope              | 0028,1053 | DS | 1               | ALWAYS | AUTO  |  |
| Rescale Type               | 0028,1054 | LO | HU              | ANAP   | FIXED |  |

Table 17: VOI LUT Module

| Attribute Name | Tag       | VR | Value | Presence of Value | Source | Comment |
|----------------|-----------|----|-------|-------------------|--------|---------|
| Window Center  | 0028,1050 | DS |       | ANAP              | AUTO   |         |
| Window Width   | 0028,1051 | DS |       | ANAP              | AUTO   |         |

Table 18: SOP Common Module

| Attribute Name         | Tag       | VR | Value | Presence of Value | Source | Comment |
|------------------------|-----------|----|-------|-------------------|--------|---------|
| Specific Character Set | 0008,0005 | CS |       | ANAP              | AUTO   |         |
| Instance Creation Date | 0008,0012 | DA |       | ANAP              | AUTO   |         |
| Instance Creation Time | 0008,0013 | TM |       | ANAP              | AUTO   |         |
| SOP Class UID          | 0008,0016 | UI |       | ALWAYS            | AUTO   |         |
| SOP Instance UID       | 0008,0018 | UI |       | ALWAYS            | AUTO   |         |
| Instance Number        | 0020,0013 | IS |       | ANAP              | AUTO   |         |

Table 19: Private Tags

| Attribute Name             | Tag       | VR | Value                     | Presence of Value | Source | Comment  |
|----------------------------|-----------|----|---------------------------|-------------------|--------|--|
| Private Creator Group 2005 | 2005,0010 | LO | Philips MR Imaging DD 006 | ANAP              | AUTO   |  |
| RT blob data               | 2005,1588 | LO |                           | ANAP              | AUTO   |  |
| RT Generated Series        | 2005,1589 | LO |                           | ANAP              | AUTO   | This is a string identifying the version of the plug-in. |

### 2.1.2.3. RT Structure Set Storage SOP Class

Table 20: IOD of Created RT Structure Set Storage SOP Class Instances

| Information Entity | Module                                | Presence Of Module |
|--------------------|---------------------------------------|--------------------|
| Patient            | Patient                               | ALWAYS             |
| Study              | General Study                         | ALWAYS             |
| Study              | Patient Study                         | ALWAYS             |
| Series             | RT Series                             | ALWAYS             |
| Equipment          | General Equipment                     | ALWAYS             |
| Structure Set      | Structure Set                         | ALWAYS             |
| Structure Set      | ROI Contour                           | ALWAYS             |
| Structure Set      | RT ROI Observations                   | ALWAYS             |
| Structure Set      | SOP Common                            | ALWAYS             |
|                    | Extended DICOM and Private attributes | CONDITIONAL        |

Table 21: Patient Module

| Attribute Name       | Tag       | VR | Value | Presence of Value | Source    | Comment |
|----------------------|-----------|----|-------|-------------------|-----------|---------|
| Patient's Name       | 0010,0010 | PN |       | ALWAYS            | MWL, USER |         |
| Patient ID           | 0010,0020 | LO |       | ALWAYS            | MWL, USER |         |
| Patient's Birth Date | 0010,0030 | DA |       | ALWAYS            | MWL, USER |         |
| Patient's Sex        | 0010,0040 | CS |       | ALWAYS            | MWL, USER |         |
| Other Patient IDs    | 0010,1000 | LO |       | VNAP              | MWL, USER |         |

Table 22: General Study Module

| Attribute Name             | Tag       | VR | Value | Presence of Value | Source          | Comment |
|----------------------------|-----------|----|-------|-------------------|-----------------|---------|
| Study Date                 | 0008,0020 | DA |       | ALWAYS            | AUTO, MWL       |         |
| Study Time                 | 0008,0030 | TM |       | ALWAYS            | MWL             |         |
| Accession Number           | 0008,0050 | SH |       | ALWAYS            | AUTO, MWL, USER |         |
| Referring Physician's Name | 0008,0090 | PN |       | VNAP              | MWL, USER       |         |
| Study Description          | 0008,1030 | LO |       | VNAP              | MWL, USER       |         |
| Study Instance UID         | 0020,000D | UI |       | ALWAYS            | AUTO, MWL       |         |
| Study ID                   | 0020,0010 | SH |       | ALWAYS            | AUTO            |         |

Table 23: Patient Study Module

| Attribute Name   | Tag       | VR | Value | Presence of Value | Source    | Comment |
|------------------|-----------|----|-------|-------------------|-----------|---------|
| Patient's Weight | 0010,1030 | DS |       | ALWAYS            | MWL, USER |         |
| Medical Alerts   | 0010,2000 | LO |       | ANAP              | MWL, USER | -       |
| Allergies        | 0010,2110 | LO |       | ANAP              | MWL, USER | -       |
| Pregnancy Status | 0010,21C0 | US |       | ANAP              | MWL, USER | -       |

Table 24: RT Series Module

| Attribute Name      | Tag       | VR | Value    | Presence of Value | Source     | Comment                    |
|---------------------|-----------|----|----------|-------------------|------------|----------------------------|
| Modality            | 0008,0060 | CS | RTSTRUCT | ALWAYS            | FIXED      | Applied value:<br>RTSTRUCT |
| Series Description  | 0008,103E | LO |          | ANAP              | AUTO, USER |                            |
| Series Instance UID | 0020,000E | UI |          | ALWAYS            | AUTO       |                            |
| Series Number       | 0020,0011 | IS |          | ALWAYS            | AUTO       |                            |
| Operator's Name     | 0008,1070 | PN |          | EMPTY             | FIXED      |                            |

Table 25: General Equipment Module

| Attribute Name                | Tag       | VR | Value | Presence of Value | Source | Comment                |
|-------------------------------|-----------|----|-------|-------------------|--------|------------------------|
| Manufacturer                  | 0008,0070 | LO |       | ALWAYS            | FIXED  |                        |
| Institution Name              | 0008,0080 | LO |       | ALWAYS            | CONFIG |                        |
| Institution Address           | 0008,0081 | ST |       | ANAPCV            | CONFIG |                        |
| Station Name                  | 0008,1010 | SH |       | ALWAYS            | CONFIG | Same as the Host Name. |
| Institutional Department Name | 0008,1040 | LO |       | ALWAYS            | CONFIG |                        |
| Manufacturer's Model Name     | 0008,1090 | LO |       | ALWAYS            | FIXED  |                        |
| Device Serial Number          | 0018,1000 | LO |       | ALWAYS            | FIXED  |                        |
| Software Version(s)           | 0018,1020 | LO |       | ALWAYS            | FIXED  |                        |

Table 26: Frame of Reference Module

| Attribute Name         | Tag       | VR | Value   | Presence of Value | Source | Comment |
|------------------------|-----------|----|---|-------------------|--------|---------|
| Frame Of Reference UID | 0020,0052 | UI | 1.3.46.670589.11.0.5.0.4<br>952.2017013010032136<br>000 | ALWAYS            | AUTO   |         |

Table 27: Structure Set Module

| Attribute Name                         | Tag       | VR | Value                     | Presence of Value | Source    | Comment                                     |
|--|-----------|----|---------------------------|-------------------|-----------|---|
| Structure Set Label                    | 3006,0002 | SH | MR-RT                     | ALWAYS            | FIXED     |   |
| Structure Set Name                     | 3006,0004 | LO | MR-RT                     | ALWAYS            | FIXED     |   |
| Structure Set Description              | 3006,0006 | ST | MR-RT AutoContouring      | ALWAYS            | FIXED     |   |
| Structure Set Date                     | 3006,0008 | DA |                           | ALWAYS            | AUTO      |   |
| Structure Set Time                     | 3006,0009 | TM |                           | ALWAYS            | AUTO      |   |
| Referenced Frame of Reference Sequence | 3006,0010 | SQ |                           | ALWAYS            | AUTO      |   |
| >Frame of Reference UID                | 0020,0052 | UI |                           | ALWAYS            | AUTO      |   |
| >RT Referenced Study Sequence          | 3006,0012 | SQ |                           | ALWAYS            | AUTO, MWL |   |
| >>Referenced SOP Class UID             | 0008,1150 | UI | 1.2.840.10008.5.1.4.1.1.2 | ALWAYS            | AUTO      | Applied value:<br>1.2.840.10008.5.1.4.1.1.2 |
| >>Referenced SOP Instance UID          | 0008,1155 | UI |                           | ALWAYS            | AUTO      |   |
| >>RT Referenced Series Sequence        | 3006,0014 | SQ |                           | ALWAYS            | AUTO      |   |
| >>>Series Instance UID                 | 0020,000E | UI |                           | ALWAYS            | AUTO      |   |
| >>>Contour Image Sequence              | 3006,0016 | SQ |                           | ALWAYS            | AUTO      |   |
| >>>>Referenced SOP Class UID           | 0008,1150 | UI |                           | ALWAYS            | AUTO      | Applied value:<br>1.2.840.10008.5.1.4.1.1.2 |
| >>>>Referenced SOP Instance UID        | 0008,1155 | UI |                           | ALWAYS            | AUTO      |   |
| Structure Set ROI Sequence             | 3006,0020 | SQ |                           | ALWAYS            | AUTO      |   |
| >ROI Number                            | 3006,0022 | IS |                           | ALWAYS            | AUTO      |   |
| >Referenced Frame of Reference UID     | 3006,0024 | UI |                           | ALWAYS            | AUTO      |   |
| >ROI Name                              | 3006,0026 | LO |                           | ALWAYS            | CONFIG    |   |
| >ROI Generation Algorithm              | 3006,0036 | CS |                           | ALWAYS            | FIXED     | Applied value:<br>AUTOMATIC                 |

Table 23: ROI Contour Module

| Attribute Name                 | Tag       | VR | Value | Presence of Value | Source | Comment   |
|--------------------------------|-----------|----|-------|-------------------|--------|---|
| ROI Contour Sequence           | 3006,0039 | SQ |       | ALWAYS            | AUTO   |   |
| >ROI Display Color             | 3006,002A | IS |       | ALWAYS            | CONFIG |   |
| >Contour Sequence              | 3006,0040 | SQ |       | VNAP              | AUTO   | Sequence may be empty or zero length if ROI does not contain any contours |
| >>Contour Image Sequence       | 3006,0016 | SQ |       | ANAP              | AUTO   |   |
| >>>Referenced SOP Class UID    | 0008,1150 | UI |       | ANAP              | AUTO   | Applied value:<br>1.2.840.10008.5.1.4.1.1.2                               |
| >>>Referenced SOP Instance UID | 0008,1155 | UI |       | ANAP              | AUTO   |   |
| >>Contour Geometric Type       | 3006,0042 | CS |       | ANAP              | AUTO   |   |
| >>Number of Contour Points     | 3006,0046 | IS |       | ANAP              | AUTO   |   |

| Attribute Name         | Tag       | VR | Value | Presence of Value | Source | Comment |
|------------------------|-----------|----|-------|-------------------|--------|---------|
| >>Contour Data         | 3006,0050 | DS |       | ANAP              | AUTO   |         |
| >Referenced ROI Number | 3006,0084 | IS |       | ALWAYS            | AUTO   |         |

Table 24: RT ROI Observations Module

| Attribute Name               | Tag       | VR | Value           | Presence of Value | Source | Comment |
|------------------------------|-----------|----|-----------------|-------------------|--------|---------|
| RT ROI Observations Sequence | 3006,0080 | SQ |                 | ALWAYS            | AUTO   |         |
| >Observation Number          | 3006,0082 | IS |                 | ALWAYS            | AUTO   |         |
| >Referenced ROI Number       | 3006,0084 | IS |                 | ALWAYS            | AUTO   |         |
| >RT ROI Interpreted Type     | 3006,00A4 | CS | EXTERNAL, ORGAN | ALWAYS            | AUTO   |         |
| >ROI Interpreter             | 3006,00A6 | PN |                 | ALWAYS            | FIXED  |         |

Table 30: SOP Common Module

| Attribute Name         | Tag       | VR | Value | Presence of Value | Source | Comment   |
|------------------------|-----------|----|-------|-------------------|--------|---|
| Specific Character Set | 0008,0005 | CS |       | ALWAYS            | AUTO   | Default: ISO_IR 100.                            |
| Instance Creation Date | 0008,0012 | DA |       | ALWAYS            | AUTO   |   |
| Instance Creation Time | 0008,0013 | TM |       | ALWAYS            | AUTO   |   |
| SOP Class UID          | 0008,0016 | UI |       | ALWAYS            | FIXED  | Applied value:<br>1.2.840.10008.5.1.4.1.1.481.3 |
| SOP Instance UID       | 0008,0018 | UI |       | ALWAYS            | AUTO   |   |

Table 31: Private Tags

| Attribute Name             | Tag       | VR | Value                     | Presence of Value | Source | Comment  |
|----------------------------|-----------|----|---------------------------|-------------------|--------|--|
| Private Creator Group 2005 | 2005,0010 | LO | Philips MR Imaging DD 006 | ANAP              | AUTO   |  |
| RT blob data               | 2005,1588 | LO |                           | ANAP              | AUTO   |  |
| RT Generated Series        | 2005,1589 | LO |                           | ANAP              | AUTO   | This is a string identifying the version of the plug-in. |