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(APACEC)	AB	ATRIAL-PACED COMPLEXES	other complexes also detected
(APACED)	AB	A-PACED COMPLEXES WITH SOME INHIBITION	non-paced complexes also detected
(APACE)	AB	ATRIAL-PACED RHYTHM	
(VPACEC)	AB	VENTRICULAR-PACED COMPLEXES	other complexes also detected
(VPACCF)	AB	AFIB/FLUT AND V-PACED COMPLEXES	other complexes, A-rate>240
(VPACCD)	AB	V-PACED COMPLEXES WITH SOME INHIBITION	non-paced complexes also detected
(VPACFD)	AB	AFIB/FLUT, V-PACED COMPLEXES WITH INHIBITION	non-paced complexes, A-rate>240
(VPACE)	AB	VENTRICULAR-PACED RHYTHM	
(ASVPC)	AB	ATRIAL-SENSED VENTRICULAR-PACED COMPLEXES	other complexes also detected
(ASVP)	AB	ATRIAL-SENSED VENTRICULAR-PACED RHYTHM	ventricular pacing tracks p-waves
(VPACEF)	AB	AFIB/FLUTTER AND VENTRICULAR-PACED RHYTHM	V-paced rhythm, A-rate>240
(AVDPC)	AB	ATRIAL-VENTRICULAR DUAL-PACED COMPLEXES	other complexes also detected

(AVDPCF)	AB	DUAL-PACEMAKER W/ A-NONCAPT DUE TO AFIB/FLUT other complexes and A-rate>240
(AVDP)	AB	ATRIAL-VENTRICULAR DUAL-PACED RHYTHM
(AVDPF)	AB	DUAL-PACEMAKER W/ A-NONCAPT DUE TO AFIB/FLUT dual pacing with A-rate>240
(PCMMC)	AB	A-V DUAL-PACED COMPLEXES W/ SOME INHIBITION other complexes also detected
(PCMM)	AB	A-V DUAL-PACED RHYTHM WITH SOME INHIBITION atrial and/or vent inhibition
(BVPACE)	AB	BIVENTRICULAR PACED RHYTHM non-simultaneous bi-vent pacing
(ABVPC)	AB	ATRIAL- BIVENTRICULAR PACED RHYTHM non-simultaneous bi-vent pacing
(PACENC)	AB	PACEMAKER FAILURE TO CAPTURE APPROPRIATELY
(PACENS)	AB	PACEMAKER FAILURE TO SENSE APPROPRIATELY
(PCNSNC)	AB	PACEMAKER FAILURE TO CAPTURE AND SENSE
(PACEM)	AB	FAILURE TO SENSE AND/OR CAPTURE (?MAGNET) fixed pacing with async rhythm
(AOO)	NS	RHYTHM CONSISTENT WITH AOO PACING
(VOO)	NS	RHYTHM CONSISTENT WITH VOO PACING
(DOO)	NS	RHYTHM CONSISTENT WITH DOO PACING
(AAI)	NS	RHYTHM CONSISTENT WITH AAI PACING
(VVI)	NS	RHYTHM CONSISTENT WITH VVI PACING

(DVI)	NS	RHYTHM CONSISTENT WITH DVI PACING	
(DDI)	NS	RHYTHM CONSISTENT WITH DDI PACING	
(VDD)	NS	RHYTHM CONSISTENT WITH VDD PACING	
(DDD)	NS	RHYTHM CONSISTENT WITH DDD PACING	
(UNKRM)	NS	UNDETERMINED RHYTHM:REVIEW	rhythm measurements incomplete
(PSAR)	AB	PACEMAKER SPIKES OR ARTIFACTS	timing non-diagnostic
(NFRA)	NS	NO FURTHER RHYTHM ANALYSIS ATTEMPTED DUE TO PACED RHYTHM	
(NFAD)	NS	NO FURTHER ANALYSIS ATTEMPTED DUE TO PACED RHYTHM	

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(SR)	NO	SINUS RHYTHM	normal P axis, V-rate ***_***
(SB)	ON	SINUS BRADYCARDIA	V-rate<***
(ST)	ON	SINUS TACHYCARDIA	V-rate>***
(SEAR)	ON	SINUS OR ECTOPIC ATRIAL RHYTHM	P axis (-45,135)
(SEAB)	ON	SINUS OR ECTOPIC ATRIAL BRADYCARDIA	P axis (-45,135), V-rate<***
(SEAT)	ON	SINUS OR ECTOPIC ATRIAL TACHYCARDIA	P axis (-45,135), V-rate>***

(EAR)	BO	ECTOPIC ATRIAL RHYTHM	abnormal P axis, normal rate
(EAB)	BO	ECTOPIC ATRIAL BRADYCARDIA	abnormal P axis, V-rate<***
(EAT)	AB	ECTOPIC ATRIAL TACHYCARDIA	abnormal P axis, V-rate>***
(LLAR)	NS	LOW LEFT ATRIAL RHYTHM	
(HLAR)	NS	HIGH LEFT ATRIAL RHYTHM	
(LRAR)	NS	LOW RIGHT ATRIAL RHYTHM	
(HRAR)	NS	HIGH RIGHT ATRIAL RHYTHM	
(JERA)	AB	ACCELERATED JUNCTIONAL ESCAPE RHYTHM	absent P waves, V-rate 50-70
(JER)	AB	JUNCTIONAL ESCAPE RHYTHM	absent P waves, slow V-rate
(JRA)	AB	ACCELERATED JUNCTIONAL RHYTHM	absent P waves, accele'd V-rate
(JT)	AB	JUNCTIONAL TACHYCARDIA	absent P waves, rapid V-rate
(RVAR)	BO	UNKNOWN RHYTHM, IRREGULAR RATE ***-***	V-rate variation >10%
(BWRV)	BO	BRADYCARDIA WITH IRREGULAR RATE ***-***	mean V-rate<***, variation>8%
(TWRV)	BO	SINUS TACHYCARDIA WITH IRREGULAR RATE ***-***	V-rate>***, variation>10%

(SA)	ON	SINUS ARRHYTHMIA, RATE ***_***	V-rate variation >10%
(SAB)	ON	SLOW SINUS ARRHYTHMIA, RATE ***_***	varied V-rate, mean<***
(SAT)	ON	FAST SINUS ARRHYTHMIA, RATE ***_***	varied V-rate, mean>***
(WPACE)	BO	WANDERING PACEMAKER	varying PR interval & P axis
(AVDIS)	AB	AV DISSOCIATION	PR variation>15%
(ETACH)	AB	EXTREME TACHYCARDIA	V-rate >(220-age)
(SVT)	AB	SUPRAVENTRICULAR TACHYCARDIA	V-rate>(220-age), QRSd<***
(AFIBT)	AB	ATRIAL FIBRILLATION WITH RAPID V-RATE	A-rate>240, V-rate>(180-age)
(TACHW)	AB	WIDE COMPLEX TACHYCARDIA	V-rate>***, QRSd>***
(VTACH)	AB	EXTREME TACHYCARDIA WITH WIDE COMPLEX, NO FURTHER RHYTHM ANALYSIS ATTEMPTED	
(ARYP)	AB	POSSIBLE ATRIAL ARRHYTHMIA, A-RATE ***	multiple Ps
(FLFIB)	AB	ATRIAL FLUTTER/FIBRILLATION, A-RATE ***	multiple Ps
(AFIB0)	AB	ATRIAL FIBRILLATION	? Atrial activity

(AFIB)	AB	ATRIAL FIBRILLATION, V-RATE ***-***	var'd rate, irreg atrial activity
(AFLT)	AB	ATRIAL FLUTTER, A-RATE ***	A-rate 220-340
(AFLT2)	AB	A-FLUTTER W/ PREDOM 2:1 AV BLOCK, A-RATE ***	A-rate 220-340, multiple Ps
(AFL2)	AB	ATRIAL FLUTTER WITH 2:1 AV BLOCK	A-rate 220-340, V-rate>***
(AFLT3)	AB	A-FLUTTER W/ PREDOM 3:1 AV BLOCK, A-RATE ***	A-rate 220-340, multiple Ps
(AFLT4)	AB	A-FLUTTER W/ PREDOM 4:1 AV BLOCK, A-RATE ***	A-rate 220-340, multiple Ps
(AFLTV)	AB	A-FLUTTER W/ VARIED AV BLOCK, A-RATE ***	A-rate 220-340, var'd AV conduc'n
(2AVB)	AB	SECOND DEGREE AV BLOCK	multiple P waves
(2AVB2)	AB	PREDOMINANT 2:1 AV BLOCK	most complexes 2 Ps
(2AVB3)	AB	PREDOMINANT 3:1 AV BLOCK	most complexes 3 Ps
(2AVB4)	AB	PREDOMINANT 4:1 AV BLOCK	most complexes 4 Ps
(2AVBV)	AB	VARYING SECOND DEGREE AV BLOCK	multiple Ps, varied AV conduction
(3AVB)	AB	COMPLETE AV BLOCK, A-RATE ***	V-rate<45, AV dissociation

(3AVBIR) AB COMPLETE AV BLOCK WITH WIDE QRS COMPLEX
V-rate<***, QRSd>***, dissociation

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(UNKBIG) NS BIGEMINY PATTERN, UNCERTAIN MECHANISM

(UNKTRI) NS TRIGEMINY PATTERN, UNCERTAIN MECHANISM

(SVTRI) NS SUPRAVENTRICULAR TRIGEMINY

(JBIG) NS JUNCTIONAL RHYTHM WITH VPC'S IN A BIGEMINY PATTERN

(JTRI) NS JUNCTIONAL RHYTHM WITH VPC'S IN A TRIGEMINY PATTERN

(ABAPC) NS ABERRANTLY CONDUCTED ATRIAL PREMATURE COMPLEXES

(UNKPC) NS PREMATURE COMPLEX, UNCERTAIN MECHANISM

(VSVPC) NS PREMATURE COMPLEX, VENT OR ABERRANT SUPRAVENT

(APC) ON ATRIAL PREMATURE COMPLEX

SV complex w/ short R-R interval

(JPC) ON JUNCTIONAL PREMATURE COMPLEX

SV complex w/ short R-R, absent P

(MAPC) AB MULTIPLE ATRIAL PREMATURE COMPLEXES

SV complexes w/ short R-R intervals

(VPC) ON VENTRICULAR PREMATURE COMPLEX

V complex w/ short R-R interval

(MVPC)	AB	MULTIPLE VENTRICULAR PREMATURE COMPLEXES	V complexes w/ short R-R intervals
(MVSPC)	AB	MULTIPLE PREMATURE COMPLEXES, VENT & SUPRAVEN	V and SV complexes w/ short R-R
(SVBIG)	AB	SUPRAVENTRICULAR BIGEMINY	bigeminy string>4 w/ SV complexes
(VBIG)	AB	VENTRICULAR BIGEMINY	bigeminy string>4 w/ V complexes
(VTRI)	AB	VENTRICULAR TRIGEMINY	trigeminy string>6 w/ V complexes
(MFVPC)	AB	MULTIFORM VENTRICULAR PREMATURE COMPLEXES	short R-R, variable morphology
(PVPC)	AB	PAIRED VENTRICULAR PREMATURE COMPLEXES	sequence of 2 V complexes
(RVPC)	AB	RUN OF VENTRICULAR PREMATURE COMPLEXES	sequence of 3 or more V complexes
(MFPVPC)	AB	PAIRED MULTIFORM VENTRICULAR COMPLEXES	sequence of 2 V complexes
(MFRVPC)	AB	RUN OF MULTIFORM VENTRICULAR COMPLEXES	sequence of 3 or more V complexes
(LRRV)	BO	LONG R-R WITH VENTRICULAR ESCAPE	R-R>175% of normal, wide QRS
(SARV)	AB	SINUS PAUSE/ARREST WITH VENTRICULAR ESCAPE	long R-R interval, wide QRS
(WENCK)	AB	MOBITZ I AV BLOCK (WENCKEBACH)	PR lengthens & dropped complexes

(RECA)	NS	RETROGRADE ATRIAL CAPTURE	
(VIC)	ON	VENTRICULAR INTERPOLATED COMPLEX	interpolated complex, wide QRS
(MVIC)	AB	MULTIPLE VENTRICULAR INTERPOLATED COMPLEXES	interpolated complexes, wide QRS
(IVPC)	ON	INTERPOLATED VENTRICULAR PREMATURE COMPLEX	interpolated complex, wide QRS
(MIVPC)	AB	MULT INTERPOLATED VENT PREMATURE COMPLEXES	interpolated complexes, wide QRSd
(ABC)	ON	ABERRANT COMPLEX	small R-R variation, aberrant QRS
(ABCS)	ON	ABERRANT COMPLEX, POSSIBLY SUPRAVENTRICULAR	aberrant shape, PR 80-220

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(SPRB)	ON	BORDERLINE SHORT PR INTERVAL	PR int <*** mS
(SPR)	BO	SHORT PR INTERVAL, ACCELERATED AV CONDUCTION	PR <*** mS
(BAVCD)	BO	BORDERLINE AV CONDUCTION DELAY	PR >***, V-rate ***-***
(1AVB)	AB	FIRST DEGREE AV BLOCK	PR >***, V-rate ***-***
(2AVBA)	NS	ADVANCED SECOND DEGREE AV BLOCK	

(SARSV)	AB	SINUS PAUSE/ARREST W/ SUPRAVENTRICULAR ESCAPE	long R-R interval, normal QRSd
(SARN)	AB	SINUS PAUSE/ARREST WITH JUNCTIONAL ESCAPE	long R-R, normal QRSd, absent P
(SARA)	AB	SINUS PAUSE/ARREST WITH ATRIAL ESCAPE	long R-R, normal QRSd, normal P
(I2AVB)	AB	INTERMITTENT SECOND DEGREE AV BLOCK	long R-R with multiple Ps
(MOBII)	AB	MOBITZ II AV BLOCK	dropped ventricular complex
(A2AVB)	AB	ALTERNATING SECOND DEGREE AV BLOCK	alternating long R-R, multiple Ps

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(VPELP)	NS	VENTRICULAR PREEXCITATION, A LEFT POSTEROSEPTAL ACCESSORY PATHWAY	
(VPERP)	NS	VENTRICULAR PREEXCITATION, A RIGHT POSTEROSEPTAL ACCESSORY PATHWAY	
(VPERA)	NS	VENTRICULAR PREEXCITATION, A RIGHT ANTEROSEPTAL ACCESSORY PATHWAY	
(VPELA)	NS	VENTRICULAR PREEXCITATION, A LEFT ANTEROSEPTAL ACCESSORY PATHWAY	
(VPELL)	NS	VENTRICULAR PREEXCITATION, A LEFT LATERAL ACCESSORY PATHWAY	
(VPERL)	NS	VENTRICULAR PREEXCITATION, A RIGHT LATERAL ACCESSORY PATHWAY	
(VPE)	AB	VENTRICULAR PREEXCITATION	Delta waves

(VPEL) AB VENT PREEXCITATION, LEFT ACCESSORY PATHWAY
Delta wave & initial axis(30,120)

(VPER) AB VENT PREEXCITATION, RIGHT ACCESSORY PATHWAY
Delta wave & initial axis(-60,29)

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(DEXC) AB CONSIDER DEXTROCARDIA
P, QRS axis rightward

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(RAE) NS RIGHT ATRIAL ENLARGEMENT

(CRAA) ON CONSIDER RIGHT ATRIAL ABNORMALITY
P >0.24 mV limb lead

(PRAA) ON PROBABLE RIGHT ATRIAL ABNORMALITY
biphasic P >0.20 mV in V1

(RAA) AB RIGHT ATRIAL ABNORMALITY
P >0.25 mV 2 lds or <-0.24 mV aVR/aVL

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(LAE) NS LEFT ATRIAL ENLARGEMENT

(CLAA) ON CONSIDER LEFT ATRIAL ABNORMALITY
wide or notched P waves

(PLAA) BO PROBABLE LEFT ATRIAL ABNORMALITY
P >50 mS, <-0.10 mV V1

(PPND) BO PROMINENT P WAVES, NONDIAGNOSTIC
wide/notched/biphasic P waves

(LAA) AB LEFT ATRIAL ABNORMALITY
P, P' > 60 mS, < -0.15 mV V1

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(LAACB) AB LAA, CONSIDER BIATRIAL ABNORMALITIES
P > 80 mS < -0.15 mV V1 & > 0.25 mV limb lds

(RAACB) AB RAA, CONSIDER BIATRIAL ABNORMALITIES
P > 0.30 mV 2 lds & < -0.30 mV aVR/aVL

(BAA) AB BIATRIAL ABNORMALITIES
P > 80 mS, < -0.15 mV V1 & > 0.30 mV 2 lds

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(AXR) ON BORDERLINE RIGHT AXIS DEVIATION
QRS axis (***, ***)

(RAD) ON RIGHT AXIS DEVIATION
QRS axis (***, ***)

(AXL) ON BORDERLINE LEFT AXIS DEVIATION
QRS axis (***, ***)

(LAD) ON LEFT AXIS DEVIATION
QRS axis (***, ***)

(AXSUP) ON SUPERIOR QRS AXIS
QRS axis (-91, 240)

(AXIND) ON INDETERMINATE QRS AXIS
QRS axis indeterminate

(S123) ON S1, S2, S3 PATTERN
S > 30 mS & > 0.2 mV, I II III

(AXPST) BO MARKEDLY POSTERIOR QRS AXIS
late V-lead transition

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(IVCD) NS INTRAVENTRICULAR CONDUCTION DELAY

(BIVCD) ON BORDERLINE INTRAVENTRICULAR CONDUCTION DELAY
QRSd >*** mS

(BIVCDL) BO BORDERLINE IVCD WITH LAD
QRSd >*** mS, axis(-90,-30)

(NIVCD) AB NONSPECIFIC INTRAVENTRICULAR CONDUCTION DELAY
QRSd >*** mS, not LBBB/RBBB

(NIVCDL) AB NONSPECIFIC IVCD WITH LAD
QRSd >*** mS & LAD

(IRBBB) AB INCOMPLETE RIGHT BUNDLE BRANCH BLOCK
QRSd >***, terminal axis (90,270)

(ARBBB) AB IVCD, CONSIDER ATYPICAL RBBB
QRSd>120 mS, terminal axis(90,270)

(CLAFB) AB LAD, CONSIDER LEFT ANTERIOR FASCICULAR BLOCK
axis(240,-40), S>R II III aVF

(LAFB) AB LEFT ANTERIOR FASCICULAR BLOCK
axis(240,-40), init forces inf

(CAFBI) AB LAD, CONSIDER LAFB OR INFERIOR INFARCT
axis(240,-30), Q&R II III aVF

(IRAFB) AB INCOMPLETE RBBB AND LAFB
axis(240,-40), S>R II III aVF

(LPFB)	AB	LEFT POSTERIOR FASCICULAR BLOCK	trm axis(110,210), init force sup
(IRPFB)	AB	IRBBB AND LPFB	RAD, QRSd>120, term axis(90,270)
(RBBB)	AB	RIGHT BUNDLE BRANCH BLOCK	QRSd>120, terminal axis(90,270)
(RLAFB)	AB	RBBB AND LAFB	QRSd >120 mS, axis(-40,240)
(RLPFB)	AB	RBBB AND LPFB	QRSd >120 mS, axis(90,210)
(ILBBB)	AB	INCOMPLETE LEFT BUNDLE BRANCH BLOCK	QRSd>110 mS, terminal axis(-90,-1)
(ALBBB)	AB	IVCD, CONSIDER ATYPICAL LBBB	QRSd>***, notch/slur R I aVL V5-6
(LBBB)	AB	LEFT BUNDLE BRANCH BLOCK	QRSd>***, broad/notched R

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(RSR1)	ON	RSR' IN V1 OR V2, PROBABLY NORMAL VARIANT	small R' only
(LT)	ON	LATE PRECORDIAL R/S TRANSITION	QRS area negative in V5/V6
(ET)	ON	EARLY PRECORDIAL R/S TRANSITION	QRS area positive in V2
(ETRSR1)	ON	RSR' IN V1 OR V2, RIGHT VCD OR RVH	QRS area positive & R' V1/V2

(CRHPI)	BO	CONSIDER RVH OR POSTERIOR INFARCT	large R in V1
(CRHPIR)	BO	CONSIDER RVH OR PMI W/ SEC REPOL ABNORMALITY	large R V1, repol abnormality
(CRVH)	BO	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	large R or R' V1/V2
(CRVHR)	AB	CONSIDER RVH W/ SECONDARY REPOL ABNORMALITY	large R in V1/V2 & repol abnrm
(PRVH)	AB	PROBABLE RIGHT VENTRICULAR HYPERTROPHY	prominent R or R' w/ RAD or RAA
(PRVHR)	AB	PROBABLE RVH W/ SECONDARY REPOL ABNORMALITY	prominent R or R' & repol abnrm
(RVH)	AB	RIGHT VENTRICULAR HYPERTROPHY	prominent R or R' w/ RAD or RAA
(RVHR)	AB	RVH WITH SECONDARY REPOLARIZATION ABNORMALITY	prom R/R', RAD/RAA & repol abnrm

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(LVHST)	NS	LVH WITH SECONDARY REPOLARIZATION CHANGES	
(HVOLT)	NS	HIGH QRS VOLTAGE	
(LVHV)	BO	LVH BY VOLTAGE	R >*** in aVL
(LVHR56)	BO	LVH BY VOLTAGE	R >*** mV in V5 or V6
(LVHRSI)	BO	LVH BY VOLTAGE	(R I+S III) >*** mV

(LVHSR)	AB	CONSIDER LEFT VENTRICULAR HYPERTROPHY (S V1/V2+R V5/V6) >*** mV
(LVHCNV)	AB	CONSIDER LEFT VENTRICULAR HYPERTROPHY (R aVL+S V3) >*** mV
(LVHC)	AB	CONSIDER LEFT VENTRICULAR HYPERTROPHY R5/6/aVL, RISIII, S12R56, S3RaVL
(LVHVP)	AB	PROBABLE LEFT VENTRICULAR HYPERTROPHY R56L/RISIII/S12R56/S3RL & LAA/LAD
(LVHCNP)	AB	PROBABLE LEFT VENTRICULAR HYPERTROPHY (RaVL+SV3) x QRSd >***
(LVHPRE)	AB	PROBABLE LVH WITH SECONDARY REPOL ABNRM R56L/RISIII/S12R56/S3RL & rep abn
(LVH)	AB	LEFT VENTRICULAR HYPERTROPHY (SV1+RV5)>3.5/(RaVL+SV3)>***
(LVH1)	AB	LEFT VENTRICULAR HYPERTROPHY R56L/RISIII/S12R56/S3RL & LAA/LAD
(LVHREP)	AB	LVH WITH SECONDARY REPOLARIZATION ABNORMALITY R56L/RISIII/S12R56/S3RL & rep abn
(LVHCO)	AB	LVH WITH IVCD AND SECONDARY REPOL ABNRM RISIII/S12R56, wQRSd, rep abnrm
(LVHCOL)	AB	LVH WITH IVCD, LAD AND SECONDARY REPOL ABNRM RISIII/S12R56, wQRS, LAD, rep abn
(BVH)	AB	BIVENTRICULAR HYPERTROPHY R/R'1 & R56L/RISIII/S12R56/S3RaVL

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(LVOLFB)	ON	BORDERLINE LOW VOLTAGE IN FRONTAL LEADS	all frontal leads <0.6 mV
(LVOLF)	ON	LOW VOLTAGE IN FRONTAL LEADS	all frontal leads <0.5 mV
(LVOLT)	BO	LOW VOLTAGE THROUGHOUT	frontal<0.5 mV, precordial<1.0 mV
(LVORAD)	BO	LOW VOLTAGE WITH RIGHT AXIS DEVIATION	low voltage, RAD
(CPDP)	BO	CHRONIC PULMONARY DISEASE PATTERN	P rightward, QRS small & vertical
(CPDLV)	BO	LOW VOLTAGE CONSISTENT WITH COPD	low voltage and Dx COPD

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(IMI)	NS	INFERIOR INFARCT	
(IMI3)	BO	BORDERLINE INFERIOR Q WAVES	Qs add to 80 mS in II III aVF
(IMI4)	BO	CONSIDER LAFB OR INFERIOR INFARCT	Qs add to 65 mS II III aVF & LAD
(IMI10)	BO	CONSIDER INFERIOR INFARCT	Q >35 mS in II III aVF
(IMI12)	BO	CONSIDER INFERIOR INFARCT	Q >25 mS, initial axis(240,-30)
(IMI18)	BO	INFERIOR Q WAVES, PROBABLY NORMAL VARIATION	Q >30 mS, age<31 male, <40 female

(IMI20)	AB	PROBABLE INFERIOR INFARCT, AGE INDETERMINATE	Q>35 mS, II III aVF
(IMI22)	AB	PROBABLE INFERIOR INFARCT, AGE INDETERMINATE	Q >35 mS, initial axis(240,-30)
(IMI26)	AB	PROBABLE INFERIOR INFARCT, AGE INDETERMINATE	Q >35 mS, T neg, II III aVF
(IMI24)	AB	PROBABLE INFERIOR INFARCT, OLD	Q>35 mS, abnormal ST-T, II III aVF
(IMI30)	AB	PROBABLE INFEROLATERAL INFARCT, AGE INDETERM	Q >30 mS in V5 V6 & IMI
(IMI49M)	AB	PROBABLE INFERIOR INFARCT, POSSIBLY RECENT	Q>35 mS, ST>0.1mV, T neg, II-aVF
(PINJI)	AB	ST ELEVATION, PROBABLE INFERIOR INJURY	inf ST>0.1 mV, lat ST<-0.05 mV
(IMI50)	AB	PROBABLE INFERIOR INFARCT, ACUTE	Q>25 mS, ST>0.10 mV, II III aVF
(IMI54)	AB	PROBABLE INFERIOR INFARCT, RECENT	Q>25 mS, ST>0.07 mV, T neg, II-aVF
(IMIQ)	AB	INFERIOR INFARCT, AGE INDETERMINATE	Q>35 mS, II III aVF
(IMI62)	AB	INFERIOR INFARCT, AGE INDETERMINATE	Q >35 mS, initial axis(240,-30)
(IMI64)	AB	INFERIOR INFARCT, OLD	Q >35 mS, flat T, II III aVF
(IMI66)	AB	INFERIOR INFARCT, AGE INDETERMINATE	Q >35 mS, T neg, II III aVF

(IMI67)	AB	INFERIOR INFARCT, POSSIBLY ACUTE	Q>35 mS, ST>0.10 mV, II III aVF
(IMIEA)	AB	INFERIOR INJURY, PROBABLE EARLY ACUTE INFARCT	ST>0.15 mV, II III aVF
(IMI80)	AB	INFERIOR Q WAVES, POSSIBLY DUE TO LBBB	Q>35 mS, II III aVF & LBBB
(IMI81)	AB	INFERIOR ST ELEVATION, POSSIBLY DUE TO LBBB	ST>0.15 mV, II III aVF & LBBB
(IMI82)	AB	PROBABLE INFERIOR INFARCT WITH LBBB	Q>35 mS, II III aVF & LBBB
(IMI74)	AB	INFERIOR INFARCT, RECENT	Q>35 mS, ST>0.07 mV, T neg, II-aVF
(IMIA)	AB	INFERIOR INFARCT, ACUTE	Q>35 mS, ST>0.10 mV, II III aVF

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(LMI)	NS	LATERAL INFARCT	
(ILMI)	NS	INFEROLATERAL INFARCT	
(ILMIQ)	NS	INFEROLATERAL INFARCT, AGE INDETERMINATE	
(ILMIA)	NS	INFEROLATERAL INFARCT, ACUTE	
(LMI10)	BO	BORDERLINE LATERAL Q WAVES	Q >35 mS, I aVL V5 V6
(LMI20)	AB	PROBABLE LATERAL INFARCT, AGE INDETERMINATE	Q >35 mS, I aVL V5 V6

(LMI26)	AB	PROBABLE LATERAL INFARCT, AGE INDETERMINATE Q>35 mS, T neg, I aVL V5 V6
(LMI24)	BO	PROBABLE LATERAL INFARCT, OLD Q>35 mS, abnormal ST-T, I aVL V5-6
(LMI28)	BO	LATERAL Q WAVES, PROBABLY DUE TO LVH Q >35 mS, I aVL V5 V6 & LVH
(LMI40)	AB	LATERAL INFARCT, AGE INDETERMINATE Q >35 mS, I aVL V5 V6
(LMI44)	AB	LATERAL INFARCT, OLD Q>35 mS, abnormal ST-T, I aVL V5 V6
(LMI46)	AB	LATERAL INFARCT, AGE INDETERMINATE Q>35 mS, T neg, I aVL V5 V6
(LMI49)	ON	LATERAL Q WAVES, PROBABLY NORMAL VARIATION Q >35 mS, age<31 male, <40 female
(LMI54)	AB	PROBABLE LATERAL INFARCT, RECENT Q>35 mS, ST>.07 mV, T neg, I aVL V5-6
(LMI50)	AB	PROBABLE LATERAL INFARCT, ACUTE Q >25 mS, ST>0.10 mV, I aVL V5 V6
(LMIQ)	AB	LATERAL INFARCT, AGE INDETERMINATE Q >35 mS, I aVL V5 V6
(LMI64)	AB	LATERAL INFARCT, OLD Q>35 mS, flat T, I aVL V5 V6
(LMI66)	AB	LATERAL INFARCT, AGE INDETERMINATE Q>35 mS, T neg, I aVL V5 V6
(LMI67)	AB	LATERAL INFARCT, POSSIBLY ACUTE Q >35 mS, ST >0.07 mV, I aVL V5 V6

(PINJL)	AB	ST ELEVATION, PROBABLE LATERAL INJURY	ST >0.08 mV, I aVL V5 V6
(LMIEA)	AB	LATERAL INJURY, PROBABLE EARLY ACUTE INFARCT	ST >0.10 mV, I aVL V5 V6
(LMI74)	AB	LATERAL INFARCT, RECENT	ST>.07 mV, T neg, Q>35, I aVL V5-6
(LMIA)	AB	LATERAL INFARCT, ACUTE	ST >.20 mV, Q >35 mS, I aVL V5 V6

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(AMI)	NS	ANTERIOR INFARCT	
(ASMI)	NS	ANTEROSEPTAL INFARCT	
(ASMIQ)	NS	ANTEROSEPTAL INFARCT, AGE INDETERMINATE	
(AMI1)	BO	BORDERLINE R WAVE PROGRESSION, ANTERIOR LEADS	R < 0.15 mV
(AMI3)	BO	Q WAVE IN V1	Q >15 mS in V1
(AMI4)	AB	ABNRM R PROG, CONSIDER ASMI OR LEAD PLACEMENT	Q >30 mS, diminished R, V2
(AMI8)	AB	CONSIDER ANTEROSEPTAL INFARCT	Q >30 mS, V1 V2
(AMI10)	AB	CONSIDER ANTEROSEPTAL INFARCT, POSSIBLY ACUTE	Q >30 mS, dimin R, ST>0.15 mV, V1-V3
(AMI12)	AB	CONSIDER ANT-SEPT INFARCT, POSSIBLY RECENT	Q, dim R, ST>0.15 mV, T neg, V1-V3

(AMI14)	AB	PROBABLE ANTEROSEPTAL INFARCT, OLD	Q >30 mS, V1 V2
(AMI16)	AB	ANTERIOR Q WAVES, POSSIBLY DUE TO ILBBB	Q >30 mS, V1 V2 & ILBBB
(AMI17)	AB	ANTERIOR Q WAVES, POSSIBLY DUE TO LVH	Q >30 mS, V1 V2 & LVH
(AMI20)	AB	PROBABLE ANTEROSEPTAL INFARCT, OLD	Q >30 mS & abn ST-T, V1-V3
(AMI21)	AB	PROBABLE ANTEROSEPTAL INFARCT, AGE INDETERM	Q >30 mS, T neg, V1-V3
(AMI21A)	AB	PROBABLE ANTEROSEPTAL INFARCT, ACUTE	Q >30 mS, ST >0.15 mV, V1-V3
(AMI22)	AB	ANT-SEPT INJURY, PROBABLE EARLY ACUTE INFARCT	ST >0.40 mV V1-V3
(ASMA)	AB	ANTEROSEPTAL INFARCT, ACUTE	Q >30 mS, ST >0.25 mV, V1-V3
(AMI26)	AB	ANTEROSEPTAL INFARCT, RECENT	Q >30 mS, ST >0.15 mS, T neg, V1-V3
(AMI30)	AB	PROBABLE ANTERIOR INFARCT, ACUTE	Q >30 mS, ST >0.15 mV, V1-V4
(AMI32)	AB	ANTERIOR INFARCT, ACUTE	Q >30 mS, ST >0.25 mV, V1-V4
(AMI34)	AB	PROBABLE ANTERIOR INFARCT, RECENT	Q >30 mS, ST >0.15 mV, T neg, V2-V4
(AMI36)	AB	ANTERIOR INFARCT, RECENT	Q >30 mS, ST >0.15 mV, T neg, V1-V4

(AMI41)	BO	CONSIDER ANTERIOR INFARCT	diminished R <0.15 mV V3
(AMI44)	BO	CONSIDER ANTERIOR INFARCT	Q >30 mS in V3
(AMI48)	BO	CONSIDER ANTERIOR INFARCT	diminished R <0.15 mV in V4
(AMI49)	BO	CONSIDER ANTERIOR INFARCT	Q >30 mS in V4
(AMI50)	AB	PROBABLE ANTERIOR INFARCT, ACUTE	Q >30 mS, dim R, ST >0.15 mV, T neg
(AMI52)	AB	PROBABLE ANTERIOR INFARCT, RECENT	Q >30 mS, dim R, ST >0.15 mV, T neg
(AMI54)	AB	ANTERIOR INFARCT, AGE INDETERMINATE	Q >30 mS in V2 V3
(AMI57)	AB	ANTERIOR Q WAVES, POSSIBLY DUE TO LVH	Q >30 mS in V1-V3 & LVH
(AMIQ)	AB	ANTERIOR INFARCT, AGE INDETERMINATE	Q >30 mS in V2-V5
(AMI60)	AB	ANTERIOR INFARCT, OLD	Q >30 mS, abnormal ST-T, V2-V5
(AMI61)	AB	ANTERIOR INFARCT, AGE INDETERMINATE	Q >30 mS, T neg, V2-V5
(AMI61A)	AB	ANTERIOR INFARCT, POSSIBLY ACUTE	Q >30 mS, ST >0.15 mV, V1-V5
(PINJA)	AB	ST ELEVATION, PROBABLE ANTERIOR INJURY	ST >0.25 mV in V1-V5

(AMIEA)	AB	ANTERIOR INJURY, EARLY ACUTE INFARCT	ST >0.35 mV in V1-V5
(AMI66)	AB	ANTERIOR INFARCT, RECENT	Q >30 mS, ST >0.15 mV, T neg, V1-V5
(AMIA)	AB	ANTERIOR INFARCT, ACUTE	ST >0.25 mV, T neg, V1-V5

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(NDSTD)	NS	NONDIAGNOSTIC ST DEPRESSION	
(SDJ)	ON	JUNCTIONAL ST DEPRESSION	ST <-0.10 mV any 3 leads
(SDM)	ON	MINIMAL ST DEPRESSION	ST <-0.05 mV in 2 leads
(SDCU)	ON	MINIMAL ST DEPRESSION	ST concave upward
(SDONS)	ON	MINIMAL ST DEPRESSION	ST <-0.03 mV, T neg, any 2 leads
(SD0AN)	ON	MINIMAL ST DEPRESSION, ANTERIOR LEADS	ST <-0.03 mV, V2-V4
(SD0LA)	ON	MINIMAL ST DEPRESSION, LATERAL LEADS	ST <-0.03 mV, I aVL V5 V6
(SD0AL)	ON	MINIMAL ST DEPRESSION, ANTEROLATERAL LEADS	ST <-0.03 mV, I aVL V2-V6
(SD0IN)	ON	MINIMAL ST DEPRESSION, INFERIOR LEADS	ST <-0.03 mV, II III aVF

(SD0DI)	ON	MINIMAL ST DEPRESSION, DIFFUSE LEADS	ST <-0.03 mV, ant/lat/inf
(SD1AN)	BO	BORDERLINE ST DEPRESSION, ANTERIOR LEADS	ST <-0.07 mV, V2-V4
(SD1LA)	BO	BORDERLINE ST DEPRESSION, LATERAL LEADS	ST <-0.07 mV, I aVL V5 V6
(SD1AL)	BO	BORDERLINE ST DEPRESSION, ANTEROLATERAL LEADS	ST <-0.07 mV, I aVL V2-V6
(SD1IN)	BO	BORDERLINE ST DEPRESSION, INFERIOR LEADS	ST <-0.07 mV, II III aVF
(SD1DI)	BO	BORDERLINE ST DEPRESSION, DIFFUSE LEADS	ST <-0.07 mV, ant/lat/inf
(SD15NS)	AB	NONSPECIFIC ST DEPRESSION	ST <-0.10 mV any 2 leads
(SD15AN)	AB	NONSPECIFIC ST DEPRESSION, ANTERIOR LEADS	ST <-0.10 mV, V2-V4
(SD15LA)	AB	NONSPECIFIC ST DEPRESSION, LATERAL LEADS	ST <-0.10 mV, I aVL V5 V6
(SD15AL)	AB	NONSPECIFIC ST DEPRESSION, ANT-LAT LEADS	ST <-0.10 mV, I aVL V2-V6
(SD15IN)	AB	NONSPECIFIC ST DEPRESSION, INFERIOR LEADS	ST <-0.10 mV, II III aVF
(SD15WI)	AB	NONSPECIFIC ST DEPRESSION, DIFFUSE LEADS	ST <-0.10 mV, ant/lat/inf
(SD2NS)	AB	NONSPECIFIC ST DEPRESSION	ST <-0.10 mV, any 2 leads

(SD2AN)	AB	ST DEPRESSION, CONSIDER ISCHEMIA, ANT LEADS ST <-0.10 mV, V2-V4
(SD2LA)	AB	ST DEPRESSION, CONSIDER ISCHEMIA, LAT LEADS ST <-0.10 mV, I aVL V5 V6
(SD2AL)	AB	ST DEPRESSION, CONSIDER ISCHEMIA, ANT-LAT LDS ST <-0.10 mV, I aVL V2-V6
(SD2IN)	AB	ST DEPRESSION, CONSIDER ISCHEMIA, INF LEADS ST <-0.10 mV, II III aVF
(SD2WI)	AB	ST DEPRESSION, CONSIDER ISCHEMIA, DIFFUSE LDS ST <-0.10 mV, ant/lat/inf
(SDPRR)	AB	ST DEPRESSION, PROBABLY RATE RELATED ST <-0.10 mV & extreme tachycardia

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(PUW)	NS	PROMINENT U WAVES
(TALVH)	BO	ABNORMAL T, PROBABLY DUE TO LVH, ANT-LAT LDS LVH & T neg, I aVL V2-V6
(LOWT)	BO	BORDERLINE T WAVE ABNORMALITIES flat T
(TAXAB)	BO	BORDERLINE T WAVE ABNORMALITIES T axis not between (-10,100)
(TAXQT)	BO	BORDERLINE T WAVE ABNORMALITIES QRS-T axis angle (91,180)
(TONS)	BO	BORDERLINE T WAVE ABNORMALITIES T/QRS ratio < 1/20 or flat T

(T0AN)	BO	BORDERLINE T ABNORMALITIES, ANTERIOR LEADS	T flat or neg, V2-V4
(T0LA)	BO	BORDERLINE T ABNORMALITIES, LATERAL LEADS	T flat/neg, I aVL V5 V6
(T0AL)	BO	BORDERLINE T ABNORMALITIES, ANT-LAT LEADS	T flat/neg, I aVL V2-V6
(T0IN)	BO	BORDERLINE T ABNORMALITIES, INFERIOR LEADS	T flat/neg, II III aVF
(T0DI)	BO	BORDERLINE T ABNORMALITIES, DIFFUSE LEADS	T flat/neg
(T1AN)	AB	NONSPECIFIC T ABNORMALITIES, ANTERIOR LEADS	T <-0.10 mV, V2-V4
(T1LA)	AB	NONSPECIFIC T ABNORMALITIES, LATERAL LEADS	T <-0.10 mV, I aVL V5 V6
(T1AL)	AB	NONSPECIFIC T ABNORMALITIES, ANT-LAT LEADS	T <-0.10 mV, I aVL V2-V6
(T1IN)	AB	NONSPECIFIC T ABNORMALITIES, INFERIOR LEADS	T <-0.10 mV, II III aVF
(T1DI)	AB	NONSPECIFIC T ABNORMALITIES, DIFFUSE LEADS	T <-0.10 mV, ant/lat/inf
(T3AN)	AB	ABNORMAL T, CONSIDER ISCHEMIA, ANTERIOR LEADS	T <-0.25 mV, V2-V4
(T1ALVH)	AB	LVH W/ REPOL ABNORMALITIES, POSSIBLE ISCHEMIA	T <-0.25 mV, V1-V3 & LVH
(T3LA)	AB	ABNORMAL T, CONSIDER ISCHEMIA, LATERAL LEADS	T <-0.25 mV, I aVL V5 V6

(T3AL)	AB	ABNORMAL T, CONSIDER ISCHEMIA, ANT-LAT LEADS T <-0.25 mV, I aVL V2-V6
(T3IN)	AB	ABNORMAL T, CONSIDER ISCHEMIA, INFERIOR LEADS T <-0.20 mV, II III aVF
(T3WI)	AB	ABNORMAL T, CONSIDER ISCHEMIA, DIFFUSE LEADS T <-0.20 mV, ant/lat/inf
(T6AN)	AB	ABNORMAL T, PROBABLE ISCHEMIA, ANTERIOR LEADS T <-0.50 mV, V2-V4
(T6LA)	AB	ABNORMAL T, PROBABLE ISCHEMIA, LATERAL LEADS T <-0.50 mV, I aVL V5 V6
(T6AL)	AB	ABNORMAL T, PROBABLE ISCHEMIA, ANT-LAT LEADS T <-0.50 mV, I aVL V2-V6
(T6IN)	AB	ABNORMAL T, PROBABLE ISCHEMIA, INFERIOR LEADS T <-0.40 mV, II III aVF
(T6IL)	AB	ABNORMAL T, PROBABLE ISCHEMIA, INFEROLATERAL T <-0.40 mV, I-III aVL aVF V5-6
(T6WI)	AB	ABNORMAL T, PROBABLE ISCHEMIA, WIDESPREAD T <-0.50 mV, ant/lat/inf

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(REPB)	BO	BORDERLINE REPOLARIZATION ABNORMALITY	ST dep & abnormal T
(REPBAN)	BO	BORDERLINE REPOL ABNORMALITY, ANT LEADS	ST dep, T flat/neg, V2-V4
(REPBLA)	BO	BORDERLINE REPOL ABNORMALITY, LATERAL LEADS	ST dep, T flat/neg, I aVL V5 V6
(REPBAL)	BO	BORDERLINE REPOL ABNORMALITY, ANT-LAT LEADS	ST dep, T flat/neg, I aVL V2-V6
(REPBIN)	BO	BORDERLINE REPOL ABNORMALITY, INFERIOR LEADS	ST dep, T flat/neg, II III aVF
(REPBIL)	BO	BORDERLINE REPOL ABNORMALITY, INF-LAT LEADS	ST dep, T flat/neg, inf/lat
(REPBDI)	BO	BORDERLINE REPOL ABNORMALITY, DIFFUSE LEADS	ST dep, T flat/neg, ant/lat/inf
(REPNS)	AB	NONSPECIFIC REPOLARIZATION ABNORMALITIES	ST dep, T neg, 2-3 leads
(REPAN)	AB	NONSPECIFIC REPOL ABNORMALITY, ANTERIOR LEADS	ST dep, T neg, V2-V4
(REPLA)	AB	NONSPECIFIC REPOL ABNORMALITY, LATERAL LEADS	ST dep, T neg, I aVL V5 V6
(REPAL)	AB	NONSPECIFIC REPOL ABNORMALITY, ANT-LAT LEADS	ST dep, T neg, I aVL V2-V6
(REPLVH)	AB	REPOL ABNORMALITY PROBABLY SECONDARY TO LVH	ST dep, T neg, I aVL V2-V6
(REPIN)	AB	NONSPECIFIC REPOL ABNORMALITY, INFERIOR LEADS	ST dep, T neg, II III aVF

(REPIL)	AB	NONSPECIFIC REPOL ABNORMALITY, INF-LAT LEADS ST dep, T neg, I-III aVL aVF V5-6
(REPIDI)	AB	NONSPECIFIC REPOL ABNORMALITY, DIFFUSE LEADS ST dep, T flat/neg, ant/lat/inf
(REPIA)	AB	REPOL ABNRM SUGGESTS ISCHEMIA, ANTERIOR LEADS ST dep, T neg, V2-V4
(REPILA)	AB	REPOL ABNRM SUGGESTS ISCHEMIA, LATERAL LEADS ST dep, T neg, I aVL V5 V6
(REPIAL)	AB	REPOL ABNRM SUGGESTS ISCHEMIA, ANT-LAT LEADS ST dep, T neg, I aVL V2-V6
(REPII)	AB	REPOL ABNRM SUGGESTS ISCHEMIA, INFERIOR LEADS ST dep, T neg, II III aVF
(REPIIL)	AB	REPOL ABNRM SUGGESTS ISCHEMIA, INFEROLATERAL ST dep, T neg, I-III aVL aVF V5-6
(REPIDI)	AB	REPOL ABNRM SUGGESTS ISCHEMIA, DIFFUSE LEADS ST-T neg, ant/lat/inf
(REPPAN)	AB	REPOL ABNRM, PROBABLE ISCHEMIA, ANTERIOR LDS ST dep, T neg, V2-V4
(REPPLA)	AB	REPOL ABNRM, PROBABLE ISCHEMIA, LATERAL LEADS ST dep, T neg, I aVL V5 V6
(REPPAL)	AB	REPOL ABNRM, PROBABLE ISCHEMIA, ANT-LAT LEADS ST dep, T neg, I aVL V2-V6
(REPPIN)	AB	REPOL ABNRM, PROBABLE ISCHEMIA, INFERIOR LDS ST dep, T neg, II III aVF
(REPPIL)	AB	REPOL ABNRM, PROBABLE ISCHEMIA, INF-LAT LDS ST dep, T neg, I-III aVL aVF V5-6

(REPPWI)	AB	REPOL ABNRM, PROBABLE ISCHEMIA, DIFFUSE LEADS	ST dep, T neg, ant/lat/inf
(REPRR)	AB	REPOLARIZATION ABNORMALITY, PROB RATE RELATED	ST dep, T neg, tachycardia
(LLINV)	AB	LATERAL LEADS ARE ALSO INVOLVED	lat Q or ST-T abnormalities

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(STEND)	NS	NONDIAGNOSTIC ST ELEVATION	
(STE)	NS	ST ELEVATION, SUBEPICARDIAL INJURY	
(MSTEA)	ON	MINIMAL ST ELEVATION, ANTERIOR LEADS	ST >0.08 mV, V1-V4
(MSTEL)	ON	MINIMAL ST ELEVATION, LATERAL LEADS	ST >0.07 mV, I aVL V5 V6
(MSTEAL)	ON	MINIMAL ST ELEVATION, ANTEROLATERAL LEADS	ST >0.06 mV, I aVL V2-V6
(MSTEI)	ON	MINIMAL ST ELEVATION, INFERIOR LEADS	ST >0.06 mV, II III aVF
(MSTED)	ON	MINIMAL ST ELEVATION, DIFFUSE LEADS	ST >0.10 mV, ant/lat/inf
(BSTE)	BO	BORDERLINE ST ELEVATION	ST >0.10 mV in 2 leads
(BSTEAL)	BO	BORDERLINE ST ELEVATION, ANTERIOR LEADS	ST >0.10 mV in V1-V4

(STELVH)	BO	ANTERIOR ST ELEVATION, PROBABLY DUE TO LVH ST >0.20 mV in V1-V4 & LVH
(BSTEL)	BO	BORDERLINE ST ELEVATION, LATERAL LEADS ST >0.06 mV, I aVL V5 V6
(BSTEAL)	BO	BORDERLINE ST ELEVATION, ANTEROLATERAL LEADS ST >0.06 mV, I aVL V2-V6
(BSTEI)	BO	BORDERLINE ST ELEVATION, INFERIOR LEADS ST >0.06 mV, II III aVF
(PERI)	AB	ST ELEVATION SUGGESTS PERICARDITIS ST >0.06 mV, ant/lat/inf
(CINJI)	AB	ST ELEVATION, CONSIDER INFERIOR INJURY ST >0.08 mV, II III aVF
(CINJA)	AB	ST ELEVATION, CONSIDER ANTERIOR INJURY ST >0.15 mV, V1-V5
(CINJL)	AB	ST ELEVATION, CONSIDER LATERAL INJURY ST >0.10 mV, I aVL V5 V6
(CINJAL)	AB	ST ELEVATION, CONSIDER ANTEROLATERAL INJURY ST >0.15 mV, I aVL V2-V6
(EREPOL)	NO	ST ELEV, PROBABLE NORMAL EARLY REPOL PATTERN ST elevation, age<55
(PERI1)	AB	ST ELEVATION SUGGESTS PERICARDITIS ST >0.10 mV, ant/lat/inf

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(TTW)	NS	TALL T WAVES
-------	----	--------------

(TTW10)	BO	TALL T, CONSIDER METABOLIC/ISCHEMIC ABNRM	T >1.2 mV
(TTW20)	BO	TALL T WAVES, CONSIDER HYPERKALEMIA	widespread tall T
(TTW30)	ON	TALL T WAVES, PROBABLY NORMAL VARIANT	T >1.2 mV, age 16-30

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(SQT)	ON	SHORT QT INTERVAL	QTc <340 mS
(HPRCA)	BO	SHORT QT INTERVAL, CONSIDER HYPERCALCEMIA	QTc <310 mS
(LQTB)	BO	BORDERLINE PROLONGED QT INTERVAL	QTc >*** mS
(LQTS)	AB	PROLONGED QT, PROBABLY SECONDARY TO WIDE QRS	QTc >*** mS w/ VCD/RVH/LVH
(LQT)	AB	PROLONGED QT INTERVAL	QTc >*** mS
(HPOCA)	AB	PROLONGED QT INTERVAL, CONSIDER HYPOCALCEMIA	QTc >520 mS
(HPOK)	AB	PROLONGED QT INTERVAL, CONSIDER HYPOKALEMIA	QTc >520 mS & ST-T abnormalities
(DIG1)	AB	REPOL ABNORMALITY, CONSIDER DIGITALIS EFFECT	short QTc & negative ST
(DIG2)	AB	REPOL ABNORMALITIES C/W DIGITALIS EFFECT	ST concave upward & digitalis

(DIG3) AB REPOL ABNORMALITIES C/W DIGITALIS EFFECT
ST-T negative & digitalis

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(DEXC) AB CONSIDER DEXTROCARDIA
P, QRS axis rightward

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(RAE) NS RIGHT ATRIAL ENLARGEMENT

(CRAA) ON CONSIDER RIGHT ATRIAL ABNORMALITY
P >0.24 mV limb lead

(PRAA) ON PROBABLE RIGHT ATRIAL ABNORMALITY
biphasic P >0.20 mV in V1

(RAA) AB RIGHT ATRIAL ABNORMALITY
P >0.25 mV 2 lds or <-0.24 mV aVR/aVL

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(LAE) NS LEFT ATRIAL ENLARGEMENT

(CLAA) ON CONSIDER LEFT ATRIAL ABNORMALITY
wide or notched P waves

(PLAA) BO PROBABLE LEFT ATRIAL ABNORMALITY
P >50 mS, <-0.10 mV V1

(PPND) BO PROMINENT P WAVES, NONDIAGNOSTIC
wide/notched/biphasic P waves

(LAA) AB LEFT ATRIAL ABNORMALITY
P, P' >60 mS, <-0.15 mV V1

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(LAACB)	AB	LAA, CONSIDER BIATRIAL ABNORMALITIES P>80mS <-.15 mV V1 & >.25 mV limb lds
(RAACB)	AB	RAA, CONSIDER BIATRIAL ABNORMALITIES P>0.30 mV 2 lds & <-0.30 mV aVR/aVL
(BAA)	AB	BIATRIAL ABNORMALITIES P>80 mS, <-0.15 mV V1 & >0.30 mV 2 lds

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(AXR)	ON	BORDERLINE RIGHT AXIS DEVIATION QRS axis ***-***
(RAD)	ON	RIGHT AXIS DEVIATION QRS axis ***-***
(AXL)	ON	BORDERLINE LEFT AXIS DEVIATION QRS axis ***-***
(LAD)	ON	LEFT AXIS DEVIATION QRS axis ***-***
(AXSUP)	ON	SUPERIOR QRS AXIS QRS axis (-91,240)
(AXIND)	ON	INDETERMINATE QRS AXIS QRS axis indeterminate
(S123)	ON	S1,S2,S3 PATTERN S >30 mS & >0.2 mV, I II III

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(IVCDP)	AB	NONSPECIFIC INTRAVENTRICULAR CONDUCTION DELAY QRS >*** mS
---------	----	--

(LAFBP)	AB	LEFT ANTERIOR FASCICULAR BLOCK	QRS axis (-60,-90)
(LBBBP)	AB	LEFT BUNDLE BRANCH BLOCK	QRSd >*** mS, late forces leftward
(IRBBTA)	BO	INCOMPLETE RIGHT BUNDLE BRANCH BLOCK	RSR' in V1, late forces anterior
(IRBBBP)	BO	INCOMPLETE RIGHT BUNDLE BRANCH BLOCK	QRSd >***, RSR' or pure R
(RBBBP)	AB	RIGHT BUNDLE BRANCH BLOCK	QRSd >***, RSR' or pure R or QR
(RBBBM)	AB	MARKED RIGHT BUNDLE BRANCH BLOCK	QRSd >160 mS
(RLAFBP)	AB	RBBB AND LAFB	QRSd >90, QRS(-60,-90)
(IVCD)	NS	INTRAVENTRICULAR CONDUCTION DELAY	

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(RSRNV)	NO	RSR' IN V1, NORMAL VARIATION	term-vector post-rightward
(IRBBRV)	BO	IRBBB, THE RSR' PATTERN MAY ALSO REFLECT RVH	IRBBB, R or R' >0.5 mV in V1-V3
(RVHS6)	BO	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	S <*** mV in V6
(RVHS5)	BO	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	S <*** mV in V5

(RVHRS6)	BO	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	R/S <*** in V6
(RVHTA)	AB	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	late forces posterior rightward
(RVHA)	AB	RIGHT AXIS DEVIATION, CONSIDER RVH	frontal & init-horiz'l axis right
(RVHRP1)	AB	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	R' >0.5 mV in V1
(RVHRS)	AB	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	R V1 + S V5 >*** mV
(RVHR1)	AB	PROBABLE RIGHT VENTRICULAR HYPERTROPHY	prominent R>*** V1 or *** V2
(RVHPR1)	AB	PROBABLE RIGHT VENTRICULAR HYPERTROPHY	pure R>*** mV in V1
(RVHT1)	AB	UPRIGHT T IN V1 OR V2, PROBABLE RVH	T >0.10 V1, 3d-9y
(RVHRD)	AB	PROBABLE RIGHT VENTRICULAR HYPERTROPHY	RAD & 1 of R/R'1/2, S5/6, R1S5, T1
(RVHQR)	AB	PROBABLE RIGHT VENTRICULAR HYPERTROPHY	QR pattern V1, 0h-2d
(RVH2V)	AB	RIGHT VENTRICULAR HYPERTROPHY	2 of R/R'V1/2, SV5/6, RV1SV5, TV1
(RVHAT)	AB	RIGHT VENTRICULAR HYPERTROPHY	RAD & upright T
(RVHVT)	AB	RIGHT VENTRICULAR HYPERTROPHY	TV1 & 1 of R/R'V1/2, SV5/6, R1S5

(RVHQRV) AB RIGHT VENTRICULAR HYPERTROPHY
QRV1 & 1 of R/R'V1/2, SV5/6, R1S5

(RVHQR3) AB RIGHT VENTRICULAR HYPERTROPHY
QR pattern V1, 3d-15y

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(LSHC) AB PROMINENT Q, CONSIDER LEFT SEPTAL HYPERTROPHY
deep Q in V5-6

(LSH) AB LEFT SEPTAL HYPERTROPHY
deep Q in V5-6, tall R in V1

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(LVHQ) BO CONSIDER LEFT VENTRICULAR HYPERTROPHY
deep Q in V5-6 or II III aVF

(LVHTA) BO CONSIDER LEFT VENTRICULAR HYPERTROPHY
prominent leftward forces

(LVHR6) BO LVH BY VOLTAGE
R >*** mV in V6

(LVHS12) BO LVH BY VOLTAGE
S <*** in V1 or *** in V2

(LVHRS) BO CONSIDER LEFT VENTRICULAR HYPERTROPHY
RV6+SV1 >*** mV

(LVHQR) AB PROBABLE LEFT VENTRICULAR HYPERTROPHY
Q>0.4 & R >*** in V5 or *** in V6

(LVHQV) AB PROBABLE LEFT VENTRICULAR HYPERTROPHY
Q56/II-aVF & 1 of S1/2, R6, S1R6

(LVHSTE)	AB	REPolarization Abnormality suggests LVH ST>0.1 mV, T>1.0 mV, I aVL V4-6
(LVHSTD)	AB	REPolarization Abnormality suggests LVH ST<-0.01 mV, T<-0.05, I aVL V4-6
(LVHR)	AB	REPolarization Abnormality suggests LVH ST dep, T neg, I aVL V4-V6
(LVHVA)	AB	PROBABLE LEFT VENTRICULAR HYPERTROPHY LAD & 1 of SV1/2, RV6, SV1+RV6
(LVHP)	AB	PROBABLE LVH W/ SECONDARY REPOL ABNORMALITIES LAD, S1/2, R6, S1R6 & repol abnrm
(LVHEV)	AB	LEFT VENTRICULAR HYPERTROPHY extreme leftward forces
(LVHVAQ)	AB	LEFT VENTRICULAR HYPERTROPHY LAD, Q or 1 of SV1/2, RV6, SV1RV6
(LVHRE)	AB	LVH W/ SECONDARY REPolarization Abnormalities LAD, Q/SV12/RV6/S1R6, repol abnrm

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(LCRVH)	AB	LVH BY VOLTAGE, ALSO CONSIDER RVH R >1 V1 & 1 of SV1/2, RV6, SV1RV6
(RCLVH)	AB	RVH, CONSIDER ASSOCIATED LVH RVH & Q<-.07 mV, R >1 mV V6
(BVHVC)	AB	CONSIDER BIVENTRICULAR HYPERTROPHY LVH & 1 of R/R'1/2, S5/6, R1+S5, T1
(BVHC)	AB	CONSIDER BIVENTRICULAR HYPERTROPHY R + S >6.0 mV in 2 of V2-V4

(BVHPED) AB BIVENTRICULAR HYPERTROPHY
R/R'1/2, S5/6, R1S5 & S1/2, R6, S1R6

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(LVOLFB) ON BORDERLINE LOW VOLTAGE IN FRONTAL LEADS
all frontal leads <0.6 mV

(LVOLF) ON LOW VOLTAGE IN FRONTAL LEADS
all frontal leads <0.5 mV

(LVOLT) BO LOW VOLTAGE THROUGHOUT
frontal<0.5 mV, precordial<1.0 mV

(LVORAD) BO LOW VOLTAGE WITH RIGHT AXIS DEVIATION
low voltage, RAD

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(PQIN) BO BORDERLINE Q WAVES IN INFERIOR LEADS
Qs add to 80 mS in II III aVF

(PQLA) BO BORDERLINE Q WAVES IN LATERAL LEADS
Q >35 mS in I aVL V5 V6

(PQAN) BO BORDERLINE Q WAVE IN ANTERIOR LEADS
Q >30 mS in V2-V5

(PQAL) BO BORDERLINE Q WAVE IN ANTEROLATERAL LEADS
Q >35 mS, I aVL V3-V6

(PIMI) AB ABNORMAL Q SUGGESTS INFERIOR INFARCT
Q >35 mS in II III aVF

(PLMI) AB ABNORMAL Q SUGGESTS LATERAL INFARCT
Q >35 mS in I aVL V5 V6

(PASMI)	AB	ABNORMAL Q SUGGESTS ANTEROSEPTAL INFARCT	Q >30 mS in V1 V2
(PAMI)	AB	ABNORMAL Q SUGGESTS ANTERIOR INFARCT	Q >30 mS in V2-V4
(PALMI)	AB	ABNORMAL Q SUGGESTS ANTEROLATERAL INFARCT	Q >30 mS I aVL V4-V6

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(NDSTD)	NS	NONDIAGNOSTIC ST DEPRESSION	
(SDANP)	BO	NONSPECIFIC ST DEPRESSION, ANTERIOR LEADS	ST <-0.10 mV, V2-V5
(SDINP)	BO	NONSPECIFIC ST DEPRESSION, INFERIOR LEADS	ST <-0.10 mV, II III aVF
(SDALP)	BO	NONSPECIFIC ST DEPRESSION, ANTEROLATERAL LDS	ST <-0.10 mV, I aVL V2-V6
(SDPRR)	AB	ST DEPRESSION, PROBABLY RATE RELATED	ST <-0.10 mV & extreme tachycardia

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(PUW)	NS	PROMINENT U WAVES	
(TIN1)	AB	ABNORMAL T WAVES, INFERIOR LEADS	T neg, II III aVF
(TAS1)	AB	ABNORMAL T WAVES, ANTEROSEPTAL LEADS	T neg, V1 V2 V3
(TARVH)	BO	ABNORMAL T, PROB SECONDARY TO RVH, ANT LEADS	RVH & T neg, V1-V3

(TAN1)	AB	ABNORMAL T WAVES, ANTERIOR LEADS	T neg, V1-V5
(TLA1)	AB	ABNORMAL T WAVES, LATERAL LEADS	T neg, I aVL V5-V6
(TAL1)	AB	ABNORMAL T WAVES, ANTEROLATERAL LEADS	T neg, I aVL V2-V6

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(ISCAS)	NS	REPolarization ABNORMALITIES SUGGEST ANTEROSEPTAL ISCHEMIA	
(ISCIL)	NS	REPolarization ABNORMALITIES SUGGEST INFEROLATERAL ISCHEMIA	
(ISCPS)	NS	REPolarization ABNORMALITIES SUGGEST POSTERIOR ISCHEMIA	
(REPB)	BO	BORDERLINE REPolarization ABNORMALITY	ST dep & abnormal T
(REPBAN)	BO	BORDERLINE REPOL ABNORMALITY, ANT LEADS	ST dep, T flat/neg, V2-V4
(REPBLA)	BO	BORDERLINE REPOL ABNORMALITY, LATERAL LEADS	ST dep, T flat/neg, I aVL V5 V6
(REPBAL)	BO	BORDERLINE REPOL ABNORMALITY, ANT-LAT LEADS	ST dep, T flat/neg, I aVL V2-V6
(REPBIN)	BO	BORDERLINE REPOL ABNORMALITY, INFERIOR LEADS	ST dep, T flat/neg, II III aVF
(REPBIL)	BO	BORDERLINE REPOL ABNORMALITY, INF-LAT LEADS	ST dep, T flat/neg, inf/lat

(REPBDI)	BO	BORDERLINE REPOL ABNORMALITY, DIFFUSE LEADS ST dep, T flat/neg, ant/lat/inf
(REPNS)	AB	NONSPECIFIC REPOLARIZATION ABNORMALITIES ST dep, T neg, 2-3 leads
(REPAN)	AB	NONSPECIFIC REPOL ABNORMALITY, ANTERIOR LEADS ST dep, T neg, V2-V4
(REPLA)	AB	NONSPECIFIC REPOL ABNORMALITY, LATERAL LEADS ST dep, T neg, I aVL V5 V6
(REPAL)	AB	NONSPECIFIC REPOL ABNORMALITY, ANT-LAT LEADS ST dep, T neg, I aVL V2-V6
(REPLVH)	AB	REPOL ABNORMALITY PROBABLY SECONDARY TO LVH ST dep, T neg, I aVL V2-V6
(REPIN)	AB	NONSPECIFIC REPOL ABNORMALITY, INFERIOR LEADS ST dep, T neg, II III aVF
(REPLI)	AB	NONSPECIFIC REPOL ABNORMALITY, INF-LAT LEADS ST dep, T neg, I-III aVL aVF V5-6
(REPDI)	AB	NONSPECIFIC REPOL ABNORMALITY, DIFFUSE LEADS ST dep, T flat/neg, ant/lat/inf

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(STEND)	NS	NONDIAGNOSTIC ST ELEVATION
(SEANP)	NO	ST ELEV, PROBABLY NORMAL VARIATION, ANT LEADS ST>0.15 mV, V2-V5
(SEINP)	NO	ST ELEVATION, PROBABLY NORMAL VARIATION, INF ST>0.15 mV, II III aVF

(SEALP)	NO	ST ELEVATION, PROB NORMAL VARIATION, ANT-LAT ST >0.15 mV, I aVL V2-V6
(PERI)	AB	ST ELEVATION SUGGESTS PERICARDITIS ST >0.06 mV, ant/lat/inf
(EREPOL)	NO	ST ELEV, PROBABLE NORMAL EARLY REPOL PATTERN ST elevation, age<55

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(TTW)	NS	TALL T WAVES
(TTW1)	ON	TALL T, PROBABLY NORMAL VARIANT, ANT-LAT LDS T >1.0 mV, I aVL V2-V6

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(SQT)	ON	SHORT QT INTERVAL QTc <340 mS
(HPRCA)	BO	SHORT QT INTERVAL, CONSIDER HYPERCALCEMIA QTc <310 mS
(LQTB)	BO	BORDERLINE PROLONGED QT INTERVAL QTc >*** mS
(LQTS)	AB	PROLONGED QT, PROBABLY SECONDARY TO WIDE QRS QTc >*** mS w/ VCD/RVH/LVH
(LQT)	AB	PROLONGED QT INTERVAL QTc >*** mS
(HPOCA)	AB	PROLONGED QT INTERVAL, CONSIDER HYPOCALCEMIA QTc >520 mS
(HPOK)	AB	PROLONGED QT INTERVAL, CONSIDER HYPOKALEMIA QTc >520 mS & ST-T abnormalities

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(ARVO)	NS	ACUTE RIGHT VENTRICULAR OVERLOAD	
(ACP)	NS	ACUTE COR PULMONALE	
(ASD)	NS	ATRIAL SEPTAL DEFECT	
(AVSD)	NS	ATRIOVENTRICULAR SEPTAL DEFECT	
(CM)	NS	CARDIOMYOPATHY	
(CTA)	AB	CONSIDER TRICUSPID ATRESIA	RAA, LAD & LVH
(CECD)	AB	CONSIDER ENDOCARDIAL CUSHION DEFECT	axis(-30,-170), RVH or RBBB
(CASD)	AB	CONSIDER ATRIAL SEPTAL DEFECT, SEPTUM SECUNDUM	QRS(1,180), RSR' in V1
(CAOCA)	AB	PROB ANT-LAT INFARCT, CONS ANOM ORIG OF C. A.	T <-0.1 mV, Q >30 mS, I aVL V4-6
(CEA)	AB	CONSIDER EBSTEIN ANOMALY	RAA, RBBB, R' <1.0 mV in V1

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(AGEUNK)	NS	AGE IS NOT ENTERED, ASSUMED TO BE AGE *** YEARS OLD FOR PURPOSE OF ECG INTERPRETATION
(PLMP)	NS	PRECARDIAL LEADS MISPLACED
(PLRV)	NS	PRECARDIAL LEAD REVERSAL
(LALLV)	NS	LEFT ARM AND LEFT LEG LEAD REVERSAL

(ECGSIM)	DE	INPUT FROM ECG SIMULATOR OR V1-V4 SHORTED ... NO ANALYSIS PERFORMED
(NAPHF)	DE	NO ANALYSIS DUE TO POSSIBLE HARDWARE FAILURE Channels 1, 2, 3 identical
(TPT)	DE	TECHNICALLY POOR TRACING - PLEASE REPEAT ECG!
(RALARV)	DE	RIGHT AND LEFT ARM LEADS REVERSED, PLEASE REPEAT ECG
(12ML)	DE	ALL 12 LEADS ARE MISSING
(MISLDS)	NS	INCOMPLETE ANALYSIS DUE TO MISSING DATA IN PRECARDIAL LEAD(S)
(MAGNET)	NS	ECG ACQUIRED WITH MAGNET IN PLACE
(NFAMLD)	DE	NO FURTHER ANALYSIS ATTEMPTED FOR THIS ECG - NOT ENOUGH LEADS COULD BE MEASURED
(QMA04)	NS	LEAD(S) *LEAD* WERE NOT USED FOR MORPHOLOGY ANALYSIS
(QMART)	NS	ARTIFACT IN LEAD(S) *LEAD*
(QMBW)	NS	BASELINE WANDER IN LEAD(S) *LEAD*
(QMAB)	NS	ARTIFACT IN LEAD(S) *LEAD* AND BASELINE WANDER IN LEAD(S) *LEAD*
(QMMLD)	NS	MISSING LEAD(S) *LEAD*
(PSREC)	NS	RECONSTRUCTED PACER SPIKES IN LD(S) *LEAD*

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A

(A2AVB)	AB	ALTERNATING SECOND DEGREE AV BLOCK alternating long R-R, multiple Ps
(AAI)	NS	RHYTHM CONSISTENT WITH AAI PACING
(ABAPC)	NS	ABERRANTLY CONDUCTED ATRIAL PREMATURE COMPLXS
(ABC)	ON	ABERRANT COMPLEX small R-R variation, aberrant QRS
(ABCS)	ON	ABERRANT COMPLEX, POSSIBLY SUPRAVENTRICULAR aberrant shape, PR 80-220
(ABVPC)	AB	ATRIAL- BIVENTRICULAR PACED RHYTHM non-simultaneous bi-vent pacing
(ACP)	NS	ACUTE COR PULMONALE
(AFIB)	AB	ATRIAL FIBRILLATION RATE, V RATE ***-*** var'd rate, irreg atrial activity
(AFIB0)	AB	ATRIAL FIBRILLATION ? Atrial activity
(AFIBT)	AB	ATRIAL FIBRILLATION WITH RAPID V-RATE A-rate>240, V-rate>(180-age)

(AFL2)	AB	ATRIAL FLUTTER WITH 2:1 AV BLOCK	A-rate 220-340, V-rate>***
(AFLT)	AB	ATRIAL FLUTTER, A-RATE ***	A-rate 220-340
(AFLT2)	AB	A-FLUTTER W/ PREDOM 2:1 AV BLOCK, A-RATE ***	A-rate 220-340, multiple Ps
(AFLT3)	AB	A-FLUTTER W/ PREDOM 3:1 AV BLOCK, A-RATE ***	A-rate 220-340, multiple Ps
(AFLT4)	AB	A-FLUTTER W/ PREDOM 4:1 AV BLOCK, A-RATE ***	A-rate 220-340, multiple Ps
(AFLTV)	AB	A-FLUTTER W/ VARIED AV BLOCK, A-RATE ***	A-rate 220-340, var'd AV conduc'n
(ALBBB)	AB	IVCD, CONSIDER ATYPICAL LBBB	QRSd>***, notch/slur R I aVL V5-6
(ALI)	NS	ANTEROLATERAL INFARCT	
(ALI10)	AB	CONSIDER ANTEROLATERAL INFARCT	Q >30 mS, I aVL V3-V6
(ALI20)	AB	PROBABLE ANTEROLATERAL INFARCT, AGE INDETERM	Q >30 mS, V3-V6
(ALI24)	AB	PROBABLE ANTEROLATERAL INFARCT, OLD	Q >30 mS, abnormal ST-T, V2-V6
(ALI26)	AB	PROBABLE ANTEROLATERAL INFARCT, AGE INDETERM	Q >30 mS, T neg, V2-V6
(ALI40)	AB	ANTEROLATERAL INFARCT, AGE INDETERMINATE	Q >35 mS, V4-V6

(ALI44)	AB	ANTEROLATERAL INFARCT, OLD	Q >35 mS, abnormal ST-T, V2-V6
(ALI46)	AB	ANTEROLATERAL INFARCT, AGE INDETERMINATE	Q >35 mS, T neg, V2-V6
(ALI48)	BO	ANTEROLATERAL Q WAVES, PROBABLY DUE TO LVH	Q >35 mS in V4-V6 & LVH
(ALI49)	BO	ANTEROLATERAL Q WAVE, PROBABLY NORMAL FOR AGE	Q >35 mS, age <31 male, <40 female
(ALI50)	AB	PROBABLE ANTEROLATERAL INFARCT, ACUTE	ST >0.15 mV, Q >30 mS, V2-V5
(ALI54)	AB	PROBABLE ANTEROLATERAL INFARCT, RECENT	Q >30 mS, ST >0.07 mV, T neg, V2-V6
(ALI64)	AB	ANTEROLATERAL INFARCT, OLD	Q >35 mS & >.10 mV, abnrm ST-T, V3-V6
(ALI66)	AB	ANTEROLATERAL INFARCT, AGE INDETERMINATE	Q >35 mS & >0.10 mV, T neg, V3-V6
(ALI67)	AB	ANTEROLATERAL INFARCT, POSSIBLY ACUTE	Q >35 mS, ST >0.15 mV, V2-V6
(ALI86)	AB	EXTENSIVE ANTERIOR INFARCT, AGE INDETERMINATE	Q >35 mS, flat/neg T, V1-V6
(ALI88)	AB	EXTENSIVE ANTERIOR INFARCT, POSSIBLY ACUTE	Q >35 mS, ST >0.15 mV, V1-V6
(ALI94)	AB	EXTENSIVE ANTERIOR INFARCT, RECENT	Q >35 mS, ST >0.07 mV, T neg, V1-V6
(ALIA)	AB	ANTEROLATERAL INFARCT, ACUTE	Q >35 mS, ST >0.20 mV, V2-V6

(ALIEA)	AB	ANTEROLATERAL INJURY, EARLY ACUTE INFARCT ST >0.15 mV, I aVL V2-V6
(ALIQ)	AB	ANTEROLATERAL INFARCT, AGE INDETERMINATE Q >35 mS & >0.10 mV in V3-V6
(ALIR)	AB	ANTEROLATERAL INFARCT, RECENT Q >35 mS, ST >0.07 mV, T neg, V2-V6
(AMI)	NS	ANTERIOR INFARCT
(AMI1)	BO	BORDERLINE R WAVE PROGRESSION, ANTERIOR LEADS R < 0.15 mV
(AMI3)	BO	Q WAVE IN V1 Q >15 mS in V1
(AMI4)	AB	ABNRM R PROG, CONSIDER ASMI OR LEAD PLACEMENT Q >30 mS, diminished R, V2
(AMI8)	AB	CONSIDER ANTEROSEPTAL INFARCT Q >30 mS, V1 V2
(AMI10)	AB	CONSIDER ANTEROSEPTAL INFARCT, POSSIBLY ACUTE Q >30 mS, dimin R, ST >0.15 mV, V1-V3
(AMI12)	AB	CONSIDER ANT-SEPT INFARCT, POSSIBLY RECENT Q, dim R, ST >0.15 mV, T neg, V1-V3
(AMI14)	AB	PROBABLE ANTEROSEPTAL INFARCT, OLD Q >30 mS, V1 V2
(AMI16)	AB	ANTERIOR Q WAVES, POSSIBLY DUE TO ILBBB Q >30 mS, V1 V2 & ILBBB
(AMI17)	AB	ANTERIOR Q WAVES, POSSIBLY DUE TO LVH Q >30 mS, V1 V2 & LVH

(AMI20)	AB	PROBABLE ANTEROSEPTAL INFARCT, OLD	Q >30 mS & abn ST-T, V1-V3
(AMI21)	AB	PROBABLE ANTEROSEPTAL INFARCT, AGE INDETERM	Q >30 mS, T neg, V1-V3
(AMI21A)	AB	PROBABLE ANTEROSEPTAL INFARCT, ACUTE	Q >30 mS, ST>0.15 mV, V1-V3
(AMI22)	AB	ANT-SEPT INJURY, PROBABLE EARLY ACUTE INFARCT	ST >0.40 mV V1-V3
(AMI26)	AB	ANTEROSEPTAL INFARCT, RECENT	Q >30 mS, ST >0.15 mS, T neg, V1-V3
(AMI30)	AB	PROBABLE ANTERIOR INFARCT, ACUTE	Q >30 mS, ST >0.15 mV, V1-V4
(AMI32)	AB	ANTERIOR INFARCT, ACUTE	Q >30 mS, ST >0.25 mV, V1-V4
(AMI34)	AB	PROBABLE ANTERIOR INFARCT, RECENT	Q >30 mS, ST >0.15 mV, T neg, V2-V4
(AMI36)	AB	ANTERIOR INFARCT, RECENT	Q >30 mS, ST >0.15 mV, T neg, V1-V4
(AMI41)	BO	CONSIDER ANTERIOR INFARCT	diminished R <0.15 mV V3
(AMI44)	BO	CONSIDER ANTERIOR INFARCT	Q >30 mS in V3
(AMI48)	BO	CONSIDER ANTERIOR INFARCT	diminished R <0.15 mV in V4
(AMI49)	BO	CONSIDER ANTERIOR INFARCT	Q >30 mS in V4

(AMI50)	AB	PROBABLE ANTERIOR INFARCT, ACUTE	Q >30 mS, dim R, ST >0.15 mV, T neg
(AMI52)	AB	PROBABLE ANTERIOR INFARCT, RECENT	Q >30 mS, dim R, ST >0.15 mV, T neg
(AMI54)	AB	ANTERIOR INFARCT, AGE INDETERMINATE	Q >30 mS in V2 V3
(AMI57)	AB	ANTERIOR Q WAVES, POSSIBLY DUE TO LVH	Q >30 mS in V1-V3 & LVH
(AMI60)	AB	ANTERIOR INFARCT, OLD	Q >30 mS, abnormal ST-T, V2-V5
(AMI61)	AB	ANTERIOR INFARCT, AGE INDETERMINATE	Q >30 mS, T neg, V2-V5
(AMI61A)	AB	ANTERIOR INFARCT, POSSIBLY ACUTE	Q >30 mS, ST >0.15 mV, V1-V5
(AMI66)	AB	ANTERIOR INFARCT, RECENT	Q >30 mS, ST >0.15 mV, T neg, V1-V5
(AMIA)	AB	ANTERIOR INFARCT, ACUTE	ST >0.25 mV, T neg, V1-V5
(AMIEA)	AB	ANTERIOR INJURY, EARLY ACUTE INFARCT	ST >0.35 mV in V1-V5
(AMIQ)	AB	ANTERIOR INFARCT, AGE INDETERMINATE	Q >30 mS in V2-V5
(AOO)	NS	RHYTHM CONSISTENT WITH AOO PACING	
(APACE)	AB	ATRIAL-PACED RHYTHM	

(APACEC)	AB	ATRIAL-PACED COMPLEXES	other complexes also detected
(APACED)	AB	A-PACED COMPLEXES WITH SOME INHIBITION	non-paced complexes also detected
(APC)	ON	ATRIAL PREMATURE COMPLEX	SV complex w/ short R-R interval
(ARBBB)	AB	IVCD, CONSIDER ATYPICAL RBBB	QRSd>120 mS, terminal axis(90,270)
(ARVO)	NS	ACUTE RIGHT VENTRICULAR OVERLOAD	
(ARYP)	AB	POSSIBLE ATRIAL ARRHYTHMIA, A-RATE ***	multiple Ps
(ASD)	NS	ATRIAL SEPTAL DEFECT	
(ASMI)	NS	ANTEROSEPTAL INFARCT	
(ASMIA)	AB	ANTEROSEPTAL INFARCT, ACUTE	Q >30 mS, ST >0.25 mV, V1-V3
(ASMIQ)	NS	ANTEROSEPTAL INFARCT, AGE INDETERMINATE	
(ASVP)	AB	ATRIAL-SENSED VENTRICULAR-PACED RHYTHM	ventricular pacing tracks p-waves
(ASVPC)	AB	ATRIAL-SENSED VENTRICULAR-PACED COMPLEXES	other complexes also detected
(AVDIS)	AB	AV DISSOCIATION	PR variation>15%
(AVDP)	AB	ATRIAL-VENTRICULAR DUAL-PACED RHYTHM	

(AVDPC)	AB	ATRIAL-VENTRICULAR DUAL-PACED COMPLEXES	other complexes also detected
(AVDPCF)	AB	DUAL-PACEMAKER W/ A-NONCAPT DUE TO AFIB/FLUT	other complexes and A-rate>240
(AVDPF)	AB	DUAL-PACEMAKER W/ A-NONCAPT DUE TO AFIB/FLUT	dual pacing with A-rate>240
(AVSD)	NS	ATRIOVENTRICULAR SEPTAL DEFECT	
(AXIND)	ON	INDETERMINATE QRS AXIS	QRS axis indeterminate
(AXL)	ON	BORDERLINE LEFT AXIS DEVIATION	QRS axis (***_***)
(AXPST)	BO	MARKEDLY POSTERIOR QRS AXIS	late V-lead transition
(AXR)	ON	BORDERLINE RIGHT AXIS DEVIATION	QRS axis (***_***)
(AXSUP)	ON	SUPERIOR QRS AXIS	QRS axis (-91,240)

B

(BAA)	AB	BIATRIAL ABNORMALITIES	P>80 mS, <-0.15 mV V1 & >0.30 mV 2 lds
(BAVCD)	BO	BORDERLINE AV CONDUCTION DELAY	PR >***, V-rate ***_***
(BIVCD)	ON	BORDERLINE INTRAVENTRICULAR CONDUCTION DELAY	QRSd >*** mS

(BIVCDL)	BO	BORDERLINE IVCD WITH LAD	QRSd >*** mS, axis(-90,-30)
(BSTE)	BO	BORDERLINE ST Elevation	ST >0.10 mV in 2 leads
(BSTE A)	BO	BORDERLINE ST ELEVATION, ANTERIOR LEADS	ST >0.10 mV in V1-V4
(BSTEAL)	BO	BORDERLINE ST ELEVATION, ANTEROLATERAL LEADS	ST >0.06 mV, I aVL V2-V6
(BSTE I)	BO	BORDERLINE ST ELEVATION, INFERIOR LEADS	ST >0.06 mV, II III aVF
(BSTE L)	BO	BORDERLINE ST ELEVATION, LATERAL LEADS	ST >0.06 mV, I aVL V5 V6
(BVH)	AB	BIVENTRICULAR HYPERTROPHY	R/R'1 & R5/6L/RISIII/S12R5/6/S3RaVL
(BVHC)	AB	CONSIDER BIVENTRICULAR HYPERTROPHY	R + S >6.0 mV in 2 of V2-V4
(BVHPED)	AB	BIVENTRICULAR HYPERTROPHY	R/R'1/2, S5/6, R1S5 & S1/2, R6, S1R6
(BVHVC)	AB	CONSIDER BIVENTRICULAR HYPERTROPHY	LVH & 1 of R/R'1/2, S5/6, R1+S5, T1
(BVPACE)	AB	BIVENTRICULAR PACED RHYTHM	non-simultaneous bi-vent pacing
(BWRV)	BO	BRADYCARDIA WITH IRREGULAR RATE ***-***	mean V-rate<***, variation>8%

C

(CAFBI)	AB	LAD, CONSIDER LAFB OR INFERIOR INFARCT axis(240,-30), Q&R II III aVF
(CAOCA)	AB	PROB ANT-LAT INFARCT, CONS ANOM ORIG OF C. A. T <-0.1 mV, Q >30 mS, I aVL V4-6
(CASD)	AB	CONSIDER ATRIAL SEPTAL DEFECT, SEPTUM SECUNDUM QRS(1,180), RSR' in V1
(CEA)	AB	CONSIDER EBSTEIN ANOMALY RAA, RBBB, R' <1.0 mV in V1
(CECD)	AB	CONSIDER ENDOCARDIAL CUSHION DEFECT axis(-30,-170), RVH or RBBB
(CINJA)	AB	ST ELEVATION, CONSIDER ANTERIOR INJURY ST >0.15 mV, V1-V5
(CINJAL)	AB	ST ELEVATION, CONSIDER ANTEROLATERAL INJURY ST >0.15 mV, I aVL V2-V6
(CINJI)	AB	ST ELEVATION, CONSIDER INFERIOR INJURY ST >0.08 mV, II III aVF
(CINJL)	AB	ST ELEVATION, CONSIDER LATERAL INJURY ST >0.10 mV, I aVL V5 V6
(CIPMI)	AB	CONSIDER INFEROPOSTERIOR INFARCT inf Q, ant R or ST dep V1-3
(CLAA)	ON	CONSIDER LEFT ATRIAL ABNORMALITY wide or notched P waves
(CLAFB)	AB	LAD, CONSIDER LEFT ANTERIOR FASCICULAR BLOCK axis(240,-40), S >R II III aVF
(CM)	NS	CARDIOMYOPATHY

(CPDLV)	BO	LOW VOLTAGE CONSISTENT WITH COPD	low voltage and Dx COPD
(CPDP)	BO	CHRONIC PULMONARY DISEASE PATTERN	P rightward, QRS small & vertical
(CPMI)	AB	CONSIDER POSTERIOR INFARCT	prom R & T in V1 V2
(CRAA)	ON	CONSIDER RIGHT ATRIAL ABNORMALITY	P >0.24 mV limb lead
(CRHPI)	BO	CONSIDER RVH OR POSTERIOR INFARCT	large R in V1
(CRHPIR)	BO	CONSIDER RVH OR PMI W/ SEC REPOL ABNORMALITY	large R V1, repol abnormality
(CRPMI)	BO	TALL R WAVE IN V2, CONSIDER RVH OR PMI	R/S ratio >3, T >0.30 mV V1 V2
(CPWI)	AB	CONSIDER POSTERIOR WALL INVOLVEMENT	prominent R T in V1 V2
(CRVH)	BO	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	large R or R' V1/V2
(CRVHR)	AB	CONSIDER RVH W/ SECONDARY REPOL ABNORMALITY	large R in V1/V2 & repol abnrm
(CTA)	AB	CONSIDER TRICUSPID ATRESIA	RAA, LAD & LVH

D

(DDD)	NS	RHYTHM CONSISTENT WITH DDD PACING	
(DDI)	NS	RHYTHM CONSISTENT WITH DDI PACING	
(DEXC)	AB	CONSIDER DEXTROCARDIA	P, QRS axis rightward
(DIG1)	AB	REPOL ABNORMALITY, CONSIDER DIGITALIS EFFECT	short QTc & negative ST
(DIG2)	AB	REPOL ABNORMALITIES C/W DIGITALIS EFFECT	ST concave upward & digitalis
(DIG3)	AB	REPOL ABNORMALITIES C/W DIGITALIS EFFECT	ST-T negative & digitalis
(DOO)	NS	RHYTHM CONSISTENT WITH DOO PACING	
(DVI)	NS	RHYTHM CONSISTENT WITH DVI PACING	

E

(EAB)	BO	ECTOPIC ATRIAL BRADYCARDIA	abnormal P axis, V-rate<***
(EAMI)	NS	EXTENSIVE ANTERIOR INFARCT	
(EAMIA)	AB	EXTENSIVE ANTERIOR INFARCT, ACUTE	Q >35 mS, ST >0.15 mV, V1-V6
(EAMIQ)	AB	EXTENSIVE ANTERIOR INFARCT, AGE INDETERMINATE	Q >35 mS, V1-V6

(EAR)	BO	ECTOPIC ATRIAL RHYTHM	abnormal P axis, normal rate
(EAT)	AB	ECTOPIC ATRIAL TACHYCARDIA	abnormal P axis, V-rate>***
(EREPOL)	NO	ST ELEV, PROBABLE NORMAL EARLY REPOL PATTERN	ST elevation, age<55
(ET)	ON	EARLY PRECORDIAL R/S TRANSITION	QRS area positive in V2
(ETACH)	AB	EXTREME TACHYCARDIA	V-rate >(220-age)
(ETRSR1)	ON	RSR' IN V1 OR V2, RIGHT VCD OR RVH	QRS area positive & R' V1/V2

F

(FLFIB)	AB	ATRIAL FLUTTER/FIBRILLATION, A-RATE ***	multiple Ps
---------	----	---	-------------

H

(HLAR)	NS	HIGH LEFT ATRIAL RHYTHM	
(HPOCA)	AB	PROLONGED QT INTERVAL, CONSIDER HYPOCALCEMIA	QTc >520 mS
(HPOK)	AB	PROLONGED QT INTERVAL, CONSIDER HYPOKALEMIA	QTc >520 mS & ST-T abnormalities
(HPRCA)	BO	SHORT QT INTERVAL, CONSIDER HYPERCALCEMIA	QTc <310 mS

(HRAR) NS HIGH RIGHT ATRIAL RHYTHM

(HVOLT) NS HIGH QRS VOLTAGE

I

(I2AVB) AB INTERMITTENT SECOND DEGREE AV BLOCK
long R-R with multiple Ps

(ILBBB) AB INCOMPLETE LEFT BUNDLE BRANCH BLOCK
QRSd>110 mS, terminal axis(-90,-1)

(ILMI) NS INFEROLATERAL INFARCT

(ILMIA) NS INFEROLATERAL INFARCT, ACUTE

(ILMIQ) NS INFEROLATERAL INFARCT, AGE INDETERMINATE

(IMI) NS INFERIOR INFARCT

(IMI3) BO BORDERLINE INFERIOR Q WAVES
Qs add to 80 mS in II III aVF

(IMI4) BO CONSIDER LAFB OR INFERIOR INFARCT
Qs add to 65 mS II III aVF & LAD

(IMI10) BO CONSIDER INFERIOR INFARCT
Q >35 mS in II III aVF

(IMI12) BO CONSIDER INFERIOR INFARCT
Q >25 mS, initial axis(240,-30)

(IMI18) BO INFERIOR Q WAVES, PROBABLY NORMAL VARIATION
Q >30 mS, age<31 male, <40 female

(IMI20) AB PROBABLE INFERIOR INFARCT, AGE INDETERMINATE
Q>35 mS, II III aVF

(IMI22)	AB	PROBABLE INFERIOR INFARCT, AGE INDETERMINATE Q >35 mS, initial axis(240,-30)
(IMI24)	AB	PROBABLE INFERIOR INFARCT, OLD Q>35 mS, abnormal ST-T, II III aVF
(IMI26)	AB	PROBABLE INFERIOR INFARCT, AGE INDETERMINATE Q >35 mS, T neg, II III aVF
(IMI30)	AB	PROBABLE INFEROLATERAL INFARCT, AGE INDETERM Q >30 mS in V5 V6 & IMI
(IMI49M)	AB	PROBABLE INFERIOR INFARCT, POSSIBLY RECENT Q>35 mS, ST>0.1 mV, T neg, II-aVF
(IMI50)	AB	PROBABLE INFERIOR INFARCT, ACUTE Q>25 mS, ST>0.10 mV, II III aVF
(IMI54)	AB	PROBABLE INFERIOR INFARCT, RECENT Q>25 mS, ST>0.07 mV, T neg, II-aVF
(IMI62)	AB	INFERIOR INFARCT, AGE INDETERMINATE Q >35 mS, initial axis(240,-30)
(IMI64)	AB	INFERIOR INFARCT, OLD Q >35 mS, flat T, II III aVF
(IMI66)	AB	INFERIOR INFARCT, AGE INDETERMINATE Q >35 mS, T neg, II III aVF
(IMI67)	AB	INFERIOR INFARCT, POSSIBLY ACUTE Q>35 mS, ST>0.10 mV, II III aVF
(IMI74)	AB	INFERIOR INFARCT, RECENT Q>35 mS, ST>0.07 mV, T neg, II-aVF
(IMI80)	AB	INFERIOR Q WAVES, POSSIBLY DUE TO LBBB Q>35 mS, II III aVF & LBBB

(IMI81)	AB	INFERIOR ST ELEVATION, POSSIBLY DUE TO LBBB ST>0.15 mV, II III aVF & LBBB
(IMI82)	AB	PROBABLE INFERIOR INFARCT WITH LBBB Q>35 mS, II III aVF & LBBB
(IMIA)	AB	INFERIOR INFARCT, ACUTE Q>35 mS, ST>0.10 mV, II III aVF
(IMIEA)	AB	INFERIOR INJURY, PROBABLE EARLY ACUTE INFARCT ST>0.15 mV, II III aVF
(IMIQ)	AB	INFERIOR INFARCT, AGE INDETERMINATE Q>35 mS, II III aVF
(IPMI)	AB	INFEROPOSTERIOR INFARCT inf Q & prom R T, ST dep V1-V3
(IPMIA)	AB	INFEROPOSTERIOR INFARCT, ACUTE ST >.10 II III aVF, <-.05 V1-V4
(IRAFB)	AB	INCOMPLETE RBBB AND LAFB axis(240,-40), S>R II III aVF
(IRBBB)	AB	INCOMPLETE RIGHT BUNDLE BRANCH BLOCK QRSd >***, terminal axis (90,270)
(IRBBBP)	BO	INCOMPLETE RIGHT BUNDLE BRANCH BLOCK QRSd >***, RSR' or pure R
(IRBBRV)	BO	IRBBB, THE RSR' PATTERN MAY ALSO REFLECT RVH IRBBB, R or R' >0.5 mV in V1-V3
(IRBBTA)	BO	INCOMPLETE RIGHT BUNDLE BRANCH BLOCK RSR' in V1, late forces anterior
(IRPFB)	AB	IRBBB AND LPFB RAD, QRSd>120, term axis(90,270)

(ISCAS)	NS	REPOLARIZATION ABNORMALITIES SUGGEST ANTEROSEPTAL ISCHEMIA	
(ISCIL)	NS	REPOLARIZATION ABNORMALITIES SUGGEST INFEROLATERAL ISCHEMIA	
(ISCPs)	NS	REPOLARIZATION ABNORMALITIES SUGGEST POSTERIOR ISCHEMIA	
(IVCD)	NS	INTRAVENTRICULAR CONDUCTION DELAY	
(IVCDP)	AB	NONSPECIFIC INTRAVENTRICULAR CONDUCTION DELAY	QRS >*** mS
(IVPC)	ON	INTERPOLATED VENTRICULAR PREMATURE COMPLEX	interpolated complex, wide QRS

J

(JBIG)	NS	JUNCTIONAL RHYTHM WITH VPC'S IN A BIGEMINY PATTERN	
(JER)	AB	JUNCTIONAL ESCAPE RHYTHM	absent P waves, slow V-rate
(JERA)	AB	ACCELERATED JUNCTIONAL ESCAPE RHYTHM	absent P waves, V-rate 50-70
(JPC)	ON	JUNCTIONAL PREMATURE COMPLEX	SV complex w/ short R-R, absent P
(JRA)	AB	ACCELERATED JUNCTIONAL RHYTHM	absent P waves, accele'd V-rate
(JT)	AB	JUNCTIONAL TACHYCARDIA	absent P waves, rapid V-rate
(JTRI)	NS	JUNCTIONAL RHYTHM WITH VPC'S IN A TRIGEMINY PATTERN	

L

(LAA)	AB	LEFT ATRIAL ABNORMALITY	P,P'>60 mS, <-0.15mV V1
(LAACB)	AB	LAA, CONSIDER BIATRIAL ABNORMALITIES	P>80 mS <-.15 mV V1&>.25 mV limb lds
(LAD)	ON	LEFT AXIS DEVIATION	QRS axis ***,***
(LAE)	NS	LEFT ATRIAL ENLARGEMENT	
(LAFB)	AB	LEFT ANTERIOR FASCICULAR BLOCK	axis(240,-40), init forces inf
(LAFBP)	AB	LEFT ANTERIOR FASCICULAR BLOCK	QRS axis (-60,-90)
(LBBB)	AB	LEFT BUNDLE BRANCH BLOCK	QRSd>***, broad/notched R
(LBBBP)	AB	LEFT BUNDLE BRANCH BLOCK	QRSd>*** mS, late forces leftward
(LCRVH)	AB	LVH BY VOLTAGE, ALSO CONSIDER RVH	R >1 V1 & 1 of SV1/2, RV6, SV1RV6
(LLAR)	NS	LOW LEFT ATRIAL RHYTHM	
(LLINV)	AB	LATERAL LEADS ARE ALSO INVOLVED	lat Q or ST-T abnormalities
(LMI)	NS	LATERAL INFARCT	
(LMI10)	BO	BORDERLINE LATERAL Q WAVES	Q >35 mS, I aVL V5 V6

(LMI20)	AB	PROBABLE LATERAL INFARCT, AGE INDETERMINATE	Q >35 mS, I aVL V5 V6
(LMI24)	BO	PROBABLE LATERAL INFARCT, OLD	Q>35mS, abnormal ST-T, I aVL V5-6
(LMI26)	AB	PROBABLE LATERAL INFARCT, AGE INDETERMINATE	Q>35 mS, T neg, I aVL V5 V6
(LMI28)	BO	LATERAL Q WAVES, PROBABLY DUE TO LVH	Q >35 mS, I aVL V5 V6 & LVH
(LMI40)	AB	LATERAL INFARCT, AGE INDETERMINATE	Q >35 mS, I aVL V5 V6
(LMI44)	AB	LATERAL INFARCT, OLD	Q>35 mS, abnormal ST-T, I aVL V5 V6
(LMI46)	AB	LATERAL INFARCT, AGE INDETERMINATE	Q>35 mS, T neg, I aVL V5 V6
(LMI49)	ON	LATERAL Q WAVES, PROBABLY NORMAL VARIATION	Q >35 mS, age<31 male, <40 female
(LMI50)	AB	PROBABLE LATERAL INFARCT, ACUTE	Q >25 mS, ST>0.10 mV, I aVL V5 V6
(LMI54)	AB	PROBABLE LATERAL INFARCT, RECENT	Q>35 mS, ST>.07 mV,T neg,I aVL V5-6
(LMI64)	AB	LATERAL INFARCT, OLD	Q>35 mS, flat T, I aVL V5 V6
(LMI66)	AB	LATERAL INFARCT, AGE INDETERMINATE	Q>35 mS, T neg, I aVL V5 V6
(LMI67)	AB	LATERAL INFARCT, POSSIBLY ACUTE	Q >35 mS, ST >0.07 mV, I aVL V5 V6

(LMI74)	AB	LATERAL INFARCT, RECENT	ST>.07 mV, T neg, Q>35, I aVL V5-6
(LMIA)	AB	LATERAL INFARCT, ACUTE	ST >.20 mV, Q >35 mS, I aVL V5 V6
(LMIEA)	AB	LATERAL INJURY, PROBABLE EARLY ACUTE INFARCT	ST >0.10 mV, I aVL V5 V6
(LMIQ)	AB	LATERAL INFARCT, AGE INDETERMINATE	Q >35 mS, I aVL V5 V6
(LOWT)	BO	BORDERLINE T WAVE ABNORMALITIES	flat T
(LPFB)	AB	LEFT POSTERIOR FASCICULAR BLOCK	trm axis(110,210), init force sup
(LQT)	AB	PROLONGED QT INTERVAL	QTc >*** mS
(LQTB)	BO	BORDERLINE PROLONGED QT INTERVAL	QTc >*** mS
(LQTS)	AB	PROLONGED QT, PROBABLY SECONDARY TO WIDE QRS	QTc >*** mS w/ VCD/RVH/LVH
(LRAR)	NS	LOW RIGHT ATRIAL RHYTHM	
(LRRV)	BO	LONG R-R WITH VENTRICULAR ESCAPE	R-R>175% of normal, wide QRS
(LSH)	AB	LEFT SEPTAL HYPERTROPHY	deep Q in V5-6, tall R in V1
(LSHC)	AB	PROMINENT Q, CONSIDER LEFT SEPTAL HYPERTROPHY	deep Q in V5-6

(LT)	ON	LATE PRECORDIAL R/S TRANSITION	QRS area negative in V5/V6
(LVH)	AB	LEFT VENTRICULAR HYPERTROPHY	$(SV1+RV5)>3.5/(RaVL+SV3)>***$
(LVH1)	AB	LEFT VENTRICULAR HYPERTROPHY	R56L/RISIII/S12R56/S3RL & LAA/LAD
(LVHC)	AB	CONSIDER LEFT VENTRICULAR HYPERTROPHY	R5/6/aVL, RISIII, S12R56, S3RaVL
(LVHCNP)	AB	PROBABLE LEFT VENTRICULAR HYPERTROPHY	$(RaVL+SV3) \times QRSd >***$
(LVHCNV)	AB	CONSIDER LEFT VENTRICULAR HYPERTROPHY	$(R aVL+S V3) >*** mV$
(LVHCO)	AB	LVH WITH IVCD AND SECONDARY REPOL ABNRM	RISIII/S12R56, wQRSd, repol abnrm
(LVHCOL)	AB	LVH WITH IVCD, LAD AND SECONDARY REPOL ABNRM	RISIII/S12R56, wQRS, LAD, rep abn
(LVHEV)		AB LEFT VENTRICULAR HYPERTROPHY	extreme leftward forces
(LVHP)		AB PROBABLE LVH W/ SECONDARY REPOL ABNORMALITIES	LAD, S1/2, R6, S1R6 & repol abnrm
(LVHPRE)	AB	PROBABLE LVH WITH SECONDARY REPOL ABNRM	R56L/RISIII/S12R56/S3RL & rep abn
(LVHQ)	BO	CONSIDER LEFT VENTRICULAR HYPERTROPHY	deep Q in V5-6 or II, III, aVF
(LVHQR)		AB PROBABLE LEFT VENTRICULAR HYPERTROPHY	$Q>0.4$ & $R >***$ in V5 or $***$ in V6

(LVHQV)		AB PROBABLE LEFT VENTRICULAR HYPERTROPHY	Q56/II-aVF & 1 of S1/2, R6, S1R6
(LVHR)	AB	REPolarIZATION ABNORMALITY SUGGESTS LVH	ST dep, T neg, I aVL V4-V6
(LVHR56)	BO	LVH BY VOLTAGE	R >*** mV in V5 or V6
(LVHR6)	BO	LVH BY VOLTAGE	R >*** in V6
(LVHRE)	AB	LVH W/ SECONDARY REPolarIZATION ABNORMALITIES	LAD, Q/SV12/RV6/S1R6, repol abnrm
(LVHREP)	AB	LVH WITH SECONDARY REPolarIZATION ABNORMALITY	R56L/RISIII/S12R56/S3RL & rep abn
(LVHRS)	BO	CONSIDER LEFT VENTRICULAR HYPERTROPHY	RV6+SV1 >***
(LVHRSD)	BO	LVH BY VOLTAGE	(R I+S III) >*** mV
(LVHS12)	BO	LVH BY VOLTAGE	S <*** in V1 or *** in V2
(LVHSR)	AB	CONSIDER LEFT VENTRICULAR HYPERTROPHY	(S V1/V2+R V5/V6) >*** mV
(LVHST)	NS	LVH WITH SECONDARY REPolarIZATION CHANGES	
(LVHSTD)	AB	REPolarIZATION ABNORMALITY SUGGESTS LVH	ST<-0.01 mV, T<-0.05, I aVL V4-6
(LVHSTE)	AB	REPolarIZATION ABNORMALITY SUGGESTS LVH	ST>0.1 mV, T>1.0 mV, I aVL V4-6

(LVHTA)	BO	CONSIDER LEFT VENTRICULAR HYPERTROPHY	prominent leftward forces
(LVHV)	BO	LVH BY VOLTAGE	R >*** in aVL
(LVHVA)		AB PROBABLE LEFT VENTRICULAR HYPERTROPHY	LAD & 1 of SV1/2, RV6, SV1+RV6
(LVHVAQ)	AB	LEFT VENTRICULAR HYPERTROPHY	LAD, Q or 1 of SV1/2, RV6, SV1RV6
(LVHVP)	AB	PROBABLE LEFT VENTRICULAR HYPERTROPHY	R56L/RISIII/S12R56/S3RL & LAA/LAD
(LVOLF)	ON	LOW VOLTAGE IN FRONTAL LEADS	all frontal leads <0.5 mV
(LVOLFB)	ON	BORDERLINE LOW VOLTAGE IN FRONTAL LEADS	all frontal leads <0.6 mV
(LVOLT)	BO	LOW VOLTAGE THROUGHOUT	frontal<0.5 mV, precordial<1.0 mV
(LVORAD)	BO	LOW VOLTAGE WITH RIGHT AXIS DEVIATION	low voltage, RAD

M

(MAPC)	AB	MULTIPLE ATRIAL PREMATURE COMPLEXES	SV complexes w/ short R-R intvls
(MFPVPC)	AB	PAIRED MULTIFORM VENTRICULAR COMPLEXES	sequence of 2 V complexes
(MFRVPC)	AB	RUN OF MULTIFORM VENTRICULAR COMPLEXES	sequence of 3 or more V complexes

(MFVPC)	AB	MULTIFORM VENTRICULAR PREMATURE COMPLEXES	short R-R, variable morphology
(MIVPC)	AB	MULT INTERPOLATED VENT PREMATURE COMPLEXES	interpolated complexes, wide QRSd
(MOBII)	AB	MOBITZ II AV BLOCK	dropped ventricular complex
(MSTEA)	ON	MINIMAL ST ELEVATION, ANTERIOR LEADS	ST >0.08 mV, V1-V4
(MSTEAL)	ON	MINIMAL ST ELEVATION, ANTEROLATERAL LEADS	ST >0.06 mV, I aVL V2-V6
(MSTED)	ON	MINIMAL ST ELEVATION, DIFFUSE LEADS	ST >0.10 mV, ant/lat/inf
(MSTEI)	ON	MINIMAL ST ELEVATION, INFERIOR LEADS	ST >0.06 mV, II III aVF
(MSTEL)	ON	MINIMAL ST ELEVATION, LATERAL LEADS	ST >0.07 mV, I aVL V5 V6
(MVIC)	AB	MULTIPLE VENTRICULAR INTERPOLATED COMPLEXES	interpolated complexes, wide QRS
(MVPC)	AB	MULTIPLE VENTRICULAR PREMATURE COMPLEXES	V complexes w/ short R-R intervals
(MVSPC)	AB	MULTIPLE PREMATURE COMPLEXES, VENT & SUPRAVEN	V and SV complexes w/ short R-R

N

(NAPHF)	DE	NO ANALYSIS DUE TO POSSIBLE HARDWARE FAILURE Channels 1, 2, 3 identical
(NDSTD)	NS	NONDIAGNOSTIC ST DEPRESSION
(NFAD)	NS	NO FURTHER ANALYSIS ATTEMPTED DUE TO PACED RHYTHM
(NFRA)	NS	NO FURTHER RHYTHM ANALYSIS ATTEMPTED DUE TO PACED RHYTHM
(NIVCD)	AB	NONSPECIFIC INTRAVENTRICULAR CONDUCTION DELAY QRSd >*** mS, not LBBB/RBBB
(NIVCDL)	AB	NONSPECIFIC IVCD WITH LAD QRSd >*** mS & LAD

P

(PACEM)	AB	FAILURE TO SENSE AND/OR CAPTURE (?MAGNET) fixed pacing with async rhythm
(PACENC)	AB	PACEMAKER FAILURE TO CAPTURE APPROPRIATELY
(PACENS)	AB	PACEMAKER FAILURE TO SENSE APPROPRIATELY
(PALMI)	AB	ABNORMAL Q SUGGESTS ANTEROLATERAL INFARCT Q>30 mS I aVL V4-V6
(PAMI)	AB	ABNORMAL Q SUGGESTS ANTERIOR INFARCT Q >30 mS in V2-V4
(PASMI)	AB	ABNORMAL Q SUGGESTS ANTEROSEPTAL INFARCT Q >30 mS in V1 V2

(PCMM)	AB	A-V DUAL-PACED RHYTHM WITH SOME INHIBITION	atrial and/or vent inhibition
(PCMMC)	AB	A-V DUAL-PACED COMPLEXES W/ SOME INHIBITION	other complexes also detected
(PCNSNC)	AB	PACEMAKER FAILURE TO CAPTURE AND SENSE	
(PERI)	AB	ST ELEVATION SUGGESTS PERICARDITIS	ST >0.06 mV, ant/lat/inf
(PERI1)	AB	ST ELEVATION SUGGESTS PERICARDITIS	ST >0.10 mV, ant/lat/inf
(PIMI)	AB	ABNORMAL Q SUGGESTS INFERIOR INFARCT	Q >35 mS in II III aVF
(PINJA)	AB	ST ELEVATION, PROBABLE ANTERIOR INJURY	ST >0.25 mV in V1-V5
(PINJAL)	AB	ST ELEVATION, PROBABLE ANTEROLATERAL INJURY	ST >0.15 mV, I aVL V2-V6
(PINJI)	AB	ST ELEVATION, PROBABLE INFERIOR INJURY	inf ST >0.1 mV, lat ST <-0.05 mV
(PINJL)	AB	ST ELEVATION, PROBABLE LATERAL INJURY	ST >0.08 mV, I aVL V5 V6
(PIPMI)	AB	PROBABLE INFEROPOSTERIOR INFARCT	IMI, R>S V1-2 or ST dep V1-V3
(PLAA)	BO	PROBABLE LEFT ATRIAL ABNORMALITY	P >50 mS, <-0.10 mV V1
(PLMI)	AB	ABNORMAL Q SUGGESTS LATERAL INFARCT	Q >35 mS in I aVL V5 V6

(PLMP)	NS	PRECARDIAL LEADS MISPLACED	
(PLRV)	NS	PRECARDIAL LEAD REVERSAL	
(PMI)	AB	POSTERIOR INFARCT	prominent R T, ST dep V1-V3
(PMIA)	AB	POSTERIOR INFARCT, ACUTE	prominent R T, ST <-.05 V1-V4
(PMIQ)	NS	POSTERIOR INFARCT, AGE INDETERMINATE	
(PPMI)	AB	PROBABLE POSTERIOR INFARCT	prominent R T & ST dep V1-V3
(PPMIA)	AB	PROBABLE POSTERIOR INFARCT, ACUTE	prominent R T, ST <-.05 V1-V3
(PPND)	BO	PROMINENT P WAVES, NONDIAGNOSTIC	wide/notched/biphasic P waves
(PQAL)	BO	BORDERLINE Q WAVE IN ANTEROLATERAL LEADS	Q >35 mS, I aVL V3-V6
(PQAN)	BO	BORDERLINE Q WAVE IN ANTERIOR LEADS	Q >30 mS in V2-V5
(PQIN)	BO	BORDERLINE Q WAVES IN INFERIOR LEADS	Qs add to 80 mS in II III aVF
(PQLA)	BO	BORDERLINE Q WAVES IN LATERAL LEADS	Q >35 mS in I aVL V5 V6
(PRAA)		ON PROBABLE RIGHT ATRIAL ABNORMALITY	biphasic P >0.20 mV in V1
(PRVH)	AB	PROBABLE RIGHT VENTRICULAR HYPERTROPHY	prominent R or R' w/ RAD or RAA

(PRVHR)	AB	PROBABLE RVH W/ SECONDARY REPOL ABNORMALITY	prominent R or R' & repol abnrm
(PSAR)	AB	PACEMAKER SPIKES OR ARTIFACTS	timing non-diagnostic
(PSREC)	NS	RECONSTRUCTED PACER SPIKES IN LD(S) ***	bedside recording
(PUW)	NS	PROMINENT U WAVES	
(PVPC)	AB	PAIRED VENTRICULAR PREMATURE COMPLEXES	sequence of 2 V complexes

R

(RAA)	AB	RIGHT ATRIAL ABNORMALITY	P>0.25 mV 2 lds or <-0.24 mV aVR/aVL
(RAACB)	AB	RAA, CONSIDER BIATRIAL ABNORMALITIES	P>0.30 mV 2 lds & <-0.30 mV aVR/aVL
(RAD)	ON	RIGHT AXIS DEVIATION	QRS axis (***,***)
(RAE)	NS	RIGHT ATRIAL ENLARGEMENT	
(RBBB)	AB	RIGHT BUNDLE BRANCH BLOCK	QRSd>120, terminal axis(90,270)
(RBBBM)	AB	MARKED RIGHT BUNDLE BRANCH BLOCK	QRSd >160 mS
(RBBBP)	AB	RIGHT BUNDLE BRANCH BLOCK	QRSd >***, RSR' or pure R or QR
(RCLVH)	AB	RVH, CONSIDER ASSOCIATED LVH	RVH & Q<-.07 mV, R >1 mV V6

(RECA)	NS	RETROGRADE ATRIAL CAPTURE	
(REPAL)	AB	NONSPECIFIC REPOL ABNORMALITY, ANT-LAT LEADS	ST dep, T neg, I aVL V2-V6
(REPAN)	AB	NONSPECIFIC REPOL ABNORMALITY, ANTERIOR LEADS	ST dep, T neg, V2-V4
(REPB)	BO	BORDERLINE REPOLARIZATION ABNORMALITY	ST dep & abnormal T
(REPBAL)	BO	BORDERLINE REPOL ABNORMALITY, ANT-LAT LEADS	ST dep, T flat/neg, I aVL V2-V6
(REPBAN)	BO	BORDERLINE REPOL ABNORMALITY, ANT LEADS	ST dep, T flat/neg, V2-V4
(REPBDI)	BO	BORDERLINE REPOL ABNORMALITY, DIFFUSE LEADS	ST dep, T flat/neg, ant/lat/inf
(REPBIL)	BO	BORDERLINE REPOL ABNORMALITY, INF-LAT LEADS	ST dep, T flat/neg, inf/lat
(REPBIN)	BO	BORDERLINE REPOL ABNORMALITY, INFERIOR LEADS	ST dep, T flat/neg, II III aVF
(REPBLA)	BO	BORDERLINE REPOL ABNORMALITY, LATERAL LEADS	ST dep, T flat/neg, I aVL V5 V6
(REPDI)	AB	NONSPECIFIC REPOL ABNORMALITY, DIFFUSE LEADS	ST dep, T flat/neg, ant/lat/inf
(REPIA)	AB	REPOL ABNRM SUGGESTS ISCHEMIA, ANTERIOR LEADS	ST dep, T neg, V2-V4
(REPIAL)	AB	REPOL ABNRM SUGGESTS ISCHEMIA, ANT-LAT LEADS	ST dep, T neg, I aVL V2-V6

(REPPLA)	AB	REPOL ABNRM, PROBABLE ISCHEMIA, LATERAL LEADS	ST dep, T neg, I aVL V5 V6
(REPPWI)	AB	REPOL ABNRM, PROBABLE ISCHEMIA, DIFFUSE LEADS	ST dep, T neg, ant/lat/inf
(REPRR)	AB	REPOLARIZATION ABNORMALITY, PROB RATE RELATED	ST dep, T neg, tachycardia
(RLAFB)	AB	RBBB AND LAFB	QRSd >120 mS, axis(-40,240)
(RLAFBP)	AB	RBBB AND LAFB	QRSd>90, QRS(-60,-90)
(RLPFB)	AB	RBBB AND LFPB	QRSd >120 mS, axis(90,210)
(RSRNV)	NO	RSR' IN V1, NORMAL VARIATION	term-vector post-rightward
(RSR1)	ON	RSR' IN V1 OR V2, PROBABLY NORMAL VARIANT	small R' only
(RSRNV)	NO	RSR' IN V1, NORMAL VARIATION	term-vector post-rightward
(RVAR)	BO	UNKNOWN RHYTHM, IRREGULAR RATE ***-***	V-rate variation >10%
(RVH)	AB	RIGHT VENTRICULAR HYPERTROPHY	prominent R or R' w/ RAD or RAA
(RVH2V)	AB	RIGHT VENTRICULAR HYPERTROPHY	2 of R/R'V1/2, SV5/6, RV1SV5, TV1
(RVHA)	AB	RIGHT AXIS DEVIATION, CONSIDER RVH	frontal & init-horiz'l axis right

(RVHAT)	AB	RIGHT VENTRICULAR HYPERTROPHY	RAD & upright T
(RVHPR1)	AB	PROBABLE RIGHT VENTRICULAR HYPERTROPHY	pure R>*** mV in V1
(RVHQR)	AB	PROBABLE RIGHT VENTRICULAR HYPERTROPHY	QR pattern V1, 0h-2d
(RVHQR3)	AB	RIGHT VENTRICULAR HYPERTROPHY	QR pattern V1, 3d-15y
(RVHQRV)	AB	RIGHT VENTRICULAR HYPERTROPHY	QRV1 & 1 of R/R'V1/2, SV5/6, R1S5
(RVHR)	AB	RVH WITH SECONDARY REPOLARIZATION ABNORMALITY	prom R/R', RAD/RAA & repol abnrm
(RVHR1)		AB PROBABLE RIGHT VENTRICULAR HYPERTROPHY	prominent R>*** V1 or *** V2
(RVHRD)		AB PROBABLE RIGHT VENTRICULAR HYPERTROPHY	RAD & 1 of R/R'1/2,S5/6,R1S5,T1
(RVHRP1)	AB	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	R' >0.5mV in V1
(RVHRS)	AB	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	R V1 + S V5 >*** mV
(RVHRS6)	BO	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	R/S <*** in V6
(RVHS5)	BO	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	S <*** mV in V5
(RVHS6)	BO	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	S <*** mV in V6

(RVHT1)	AB	UPRIGHT T IN V1 OR V2, PROBABLE RVH	T >0.10 V1, 3d-9y
(RVHTA)	AB	CONSIDER RIGHT VENTRICULAR HYPERTROPHY	late forces posterior rightward
(RVHVT)	AB	RIGHT VENTRICULAR HYPERTROPHY	TV1 & 1 of R/R'V1/2, SV5/6, R1S5
(RVPC)	AB	RUN OF VENTRICULAR PREMATURE COMPLEXES	sequence of 3 or more V complexes

S

(S123)	ON	S1,S2,S3 PATTERN	S >30 mS & >0.2 mV, I II III
(SA)	ON	SINUS ARRHYTHMIA, RATE ***_***	V-rate variation >10%
(SAB)	ON	SLOW SINUS ARRHYTHMIA, RATE ***_***	varied V-rate, mean<***
(SARA)	AB	SINUS PAUSE/ARREST WITH ATRIAL ESCAPE	long R-R, normal QRSd, normal P
(SARN)	AB	SINUS PAUSE/ARREST WITH JUNCTIONAL ESCAPE	long R-R, normal QRSd, absent P
(SARSV)	AB	SINUS PAUSE/ARREST W/ SUPRAVENTRICULAR ESCAPE	long R-R interval, normal QRSd
(SARV)	AB	SINUS PAUSE/ARREST WITH VENTRICULAR ESCAPE	long R-R interval, wide QRS
(SAT)	ON	FAST SINUS ARRHYTHMIA, RATE ***_***	varied V-rate, mean>***

(SB)	ON	SINUS BRADYCARDIA	V-rate<***
(SDALP)	BO	NONSPECIFIC ST DEPRESSION, ANTEROLATERAL LDS	ST <-0.10 mV, I aVL V2-V6
(SDANP)	BO	NONSPECIFIC ST DEPRESSION, ANTERIOR LEADS	ST <-0.10 mV, V2-V5
(SDINP)	BO	NONSPECIFIC ST DEPRESSION, INFERIOR LEADS	ST <-0.10 mV, II III aVF
(SD0AL)	ON	MINIMAL ST DEPRESSION, ANTEROLATERAL LEADS	ST <-0.03 mV, I aVL V2-V6
(SD0AN)	ON	MINIMAL ST DEPRESSION, ANTERIOR LEADS	ST <-0.03 mV, V2-V4
(SD0DI)	ON	MINIMAL ST DEPRESSION, DIFFUSE LEADS	ST <-0.03 mV, ant/lat/inf
(SD0IN)	ON	MINIMAL ST DEPRESSION, INFERIOR LEADS	ST <-0.03 mV, II III aVF
(SD0LA)	ON	MINIMAL ST DEPRESSION, LATERAL LEADS	ST <-0.03 mV, I aVL V5 V6
(SD0NS)	ON	MINIMAL ST DEPRESSION	ST <-0.03 mV, T neg, any 2 leads
(SDPRR)	AB	ST DEPRESSION, PROBABLY RATE RELATED	ST <-0.10 mV & extreme tachycardia
(SD1AL)	BO	BORDERLINE ST DEPRESSION, ANTEROLATERAL LEADS	ST <-0.07 mV, I aVL V2-V6
(SD1AN)	BO	BORDERLINE ST DEPRESSION, ANTERIOR LEADS	ST <-0.07 mV, V2-V4"

(SD1DI)	BO	BORDERLINE ST DEPRESSION, DIFFUSE LEADS	ST <-0.07 mV, ant/lat/inf
(SD1IN)	BO	BORDERLINE ST DEPRESSION, INFERIOR LEADS	ST <-0.07 mV, II III aVF
(SD1LA)	BO	BORDERLINE ST DEPRESSION, LATERAL LEADS	ST <-0.07 mV, I aVL V5 V6
(SD15AL)	AB	NONSPECIFIC ST DEPRESSION, ANT-LAT LEADS	ST <-0.10 mV, I aVL V2-V6
(SD15AN)	AB	NONSPECIFIC ST DEPRESSION, ANTERIOR LEADS	ST <-0.10 mV, V2-V4
(SD15LA)	AB	NONSPECIFIC ST DEPRESSION, LATERAL LEADS	ST <-0.10 mV, I aVL V5 V6
(SD15IN)	AB	NONSPECIFIC ST DEPRESSION, INFERIOR LEADS	ST <-0.10 mV, II III aVF
(SD15NS)	AB	NONSPECIFIC ST DEPRESSION	ST <-0.10 mV any 2 leads
(SD15WI)	AB	NONSPECIFIC ST DEPRESSION, DIFFUSE LEADS	ST <-0.10 mV, ant/lat/inf
(SD2AL)	AB	ST DEPRESSION, CONSIDER ISCHEMIA, ANT-LAT LDS	ST <-0.10 mV, I aVL V2-V6
(SD2AN)	AB	ST DEPRESSION, CONSIDER ISCHEMIA, ANT LEADS	ST <-0.10 mV, V2-V4
(SD2IN)	AB	ST DEPRESSION, CONSIDER ISCHEMIA, INF LEADS	ST <-0.10 mV, II III aVF
(SD2LA)	AB	ST DEPRESSION, CONSIDER ISCHEMIA, LAT LEADS	ST <-0.10 mV, I aVL V5 V6

(SD2NS)	AB	NONSPECIFIC ST DEPRESSION	ST <-0.10 mV, any 2 leads
(SD2WI)	AB	ST DEPRESSION, CONSIDER ISCHEMIA, DIFFUSE LDS	ST <-0.10 mV, ant/lat/inf
(SDCU)	ON	MINIMAL ST DEPRESSION	ST concave upward
(SDJ)	ON	JUNCTIONAL ST DEPRESSION	ST <-0.10 mV any 3 leads
(SDM)	ON	MINIMAL ST DEPRESSION	ST <-0.05 mV in 2 leads
(SDPRR)	AB	ST DEPRESSION, PROBABLY RATE RELATED	ST <-0.10 mV & extreme tachycardia
(SEAB)	ON	SINUS OR ECTOPIC ATRIAL BRADYCARDIA	P axis (-45,135), V-rate<***
(SEALP)	NO	ST ELEVATION, PROB NORMAL VARIATION, ANT-LAT	ST >0.15 mV, I aVL V2-V6
(SEANP)	NO	ST ELEV, PROBABLY NORMAL VARIATION, ANT LEADS	ST >0.15 mV, V2-V5
(SEAR)	ON	SINUS OR ECTOPIC ATRIAL RHYTHM	P axis (-45,135)
(SEAT)	ON	SINUS OR ECTOPIC ATRIAL TACHYCARDIA	P axis (-45,135), V-rate>***
(SEINP)	NO	ST ELEVATION, PROBABLY NORMAL VARIATION, INF	ST >0.15 mV, II III aVF
(SPR)	BO	SHORT PR INTERVAL, ACCELERATED AV CONDUCTION	PR <*** mS

(SPRB)	ON	BORDERLINE SHORT PR INTERVAL	PR int <*** mS
(SQT)	ON	SHORT QT INTERVAL	QTc <340 mS
(SR)	NO	SINUS RHYTHM	normal P axis, V-rate ***-***
(ST)	ON	SINUS TACHYCARDIA	V-rate>***
(STE)	NS	ST ELEVATION, SUBEPICARDIAL INJURY	
(STELVH)	BO	ANTERIOR ST ELEVATION, PROBABLY DUE TO LVH	ST >0.20 mV in V1-V4 & LVH
(STEND)	NS	NONDIAGNOSTIC ST ELEVATION	
(SVBIG)	AB	SUPRAVENTRICULAR BIGEMINY	bigeminy string>4 w/ SV complexes
(SVT)	AB	SUPRAVENTRICULAR TACHYCARDIA	V-rate>(220-age), QRSd<***
(SVTRI)	NS	SUPRAVENTRICULAR TRIGEMINY	

T

(T0AL)	BO	BORDERLINE T ABNORMALITIES, ANT-LAT LEADS	T flat/neg, I aVL V2-V6
(T0AN)	BO	BORDERLINE T ABNORMALITIES, ANTERIOR LEADS	T flat or neg, V2-V4
(T0DI)	BO	BORDERLINE T ABNORMALITIES, DIFFUSE LEADS	T flat/neg

(T0IN)	BO	BORDERLINE T ABNORMALITIES, INFERIOR LEADS	T flat/neg, II III aVF
(T0LA)	BO	BORDERLINE T ABNORMALITIES, LATERAL LEADS	T flat/neg, I aVL V5 V6
(T0NS)	BO	BORDERLINE T WAVE ABNORMALITIES	T/QRS ratio < 1/20 or flat T
(T1AL)	AB	NONSPECIFIC T ABNORMALITIES, ANT-LAT LEADS	T < -0.10 mV, I aVL V2-V6
(T1AN)	AB	NONSPECIFIC T ABNORMALITIES, ANTERIOR LEADS	T < -0.10 mV, V2-V4
(T1DI)	AB	NONSPECIFIC T ABNORMALITIES, DIFFUSE LEADS	T < -0.10 mV, ant/lat/inf
(T1IN)	AB	NONSPECIFIC T ABNORMALITIES, INFERIOR LEADS	T < -0.10 mV, II III aVF
(T1LA)	AB	NONSPECIFIC T ABNORMALITIES, LATERAL LEADS	T < -0.10 mV, I aVL V5 V6
(T3AL)	AB	ABNORMAL T, CONSIDER ISCHEMIA, ANT-LAT LEADS	T < -0.25 mV, I aVL V2-V6
(T3AN)	AB	ABNORMAL T, CONSIDER ISCHEMIA, ANTERIOR LEADS	T < -0.25 mV, V2-V4
(T3IN)	AB	ABNORMAL T, CONSIDER ISCHEMIA, INFERIOR LEADS	T < -0.20 mV, II III aVF
(T3LA)	AB	ABNORMAL T, CONSIDER ISCHEMIA, LATERAL LEADS	T < -0.25 mV, I aVL V5 V6
(T3WI)	AB	ABNORMAL T, CONSIDER ISCHEMIA, DIFFUSE LEADS	T < -0.20 mV, ant/lat/inf

(T6AL)	AB	ABNORMAL T, PROBABLE ISCHEMIA, ANT-LAT LEADS T <-0.50 mV, I aVL V2-V6
(T6AN)	AB	ABNORMAL T, PROBABLE ISCHEMIA, ANTERIOR LEADS T <-0.50 mV, V2-V4
(T6IL)	AB	ABNORMAL T, PROBABLE ISCHEMIA, INFEROLATERAL T <-0.40 mV, I-III aVL aVF V5-6
(T6IN)	AB	ABNORMAL T, PROBABLE ISCHEMIA, INFERIOR LEADS T <-0.40 mV, II III aVF
(T6LA)	AB	ABNORMAL T, PROBABLE ISCHEMIA, LATERAL LEADS T <-0.50 mV, I aVL V5 V6
(T6WI)	AB	ABNORMAL T, PROBABLE ISCHEMIA, WIDESPREAD T <-0.50 mV, ant/lat/inf
(TACHW)	AB	WIDE COMPLEX TACHYCARDIA V-rate>***, QRSd>***
(TAL1)	AB	ABNORMAL T WAVES, ANTEROLATERAL LEADS T neg, I aVL V2-V6
(TALVH)	BO	ABNORMAL T, PROBABLY DUE TO LVH, ANT-LAT LDS LVH & T neg, I aVL V2-V6
(TAN1)	AB	ABNORMAL T WAVES, ANTERIOR LEADS T neg, V1-V5
(TAS1)	AB	ABNORMAL T WAVES, ANTEROSEPTAL LEADS T neg, V1 V2 V3
(TAXAB)	BO	BORDERLINE T WAVE ABNORMALITIES T axis not between (-10,100)
(TAXQT)	BO	BORDERLINE T WAVE ABNORMALITIES QRS-T axis angle (91,180)

(TIALVH)	AB	LVH W/ REPOL ABNORMALITIES, POSSIBLE ISCHEMIA T <-0.25 mV, V1-V3 & LVH
(TIN1)	AB	ABNORMAL T WAVES, INFERIOR LEADS T neg, II III aVF
(TLA1)	AB	ABNORMAL T WAVES, LATERAL LEADS T neg, I aVL V5-V6
(TPT)	DE	TECHNICALLY POOR TRACING - PLEASE REPEAT ECG!
(TTW)	NS	TALL T WAVES
(TTW1)	ON	TALL T, PROBABLY NORMAL VARIANT, ANT-LAT LDS T >1.0 mV, I aVL V2-V6
(TTW10)	BO	TALL T, CONSIDER METABOLIC/ISCHEMIC ABNRM T >1.2 mV
(TTW20)	BO	TALL T WAVES, CONSIDER HYPERKALEMIA widespread tall T
(TTW30)	ON	TALL T WAVES, PROBABLY NORMAL VARIANT T >1.2 mV, age 16-30
(TWRV)	BO	SINUS TACHYCARDIA WITH IRREGULAR RATE ***_*** V-rate>***,variation>10%

U

(UNKBIG)	NS	BIGEMINY PATTERN, UNCERTAIN MECHANISM
(UNKPC)	NS	PREMATURE COMPLEX, UNCERTAIN MECHANISM
(UNKRM)	NS	UNDETERMINED RHYTHM:REVIEW rhythm measurements incomplete
(UNKTRI)	NS	TRIGEMINY PATTERN, UNCERTAIN MECHANISM

V

(VBIG)	AB	VENTRICULAR BIGEMINY	bigeminy string>4 w/ V complexes
(VDD)	NS	RHYTHM CONSISTENT WITH VDD PACING	
(VIC)	ON	VENTRICULAR INTERPOLATED COMPLEX	interpolated complex, wide QRS
(VOO)	NS	RHYTHM CONSISTENT WITH VOO PACING	
(VPACCD)	AB	V-PACED COMPLEXES WITH SOME INHIBITION	non-paced complexes also detected
(VPACCF)	AB	AFIB/FLUT AND V-PACED COMPLEXES	other complexes, A-rate>240
(VPACE)	AB	VENTRICULAR-PACED RHYTHM	
(VPACEC)	AB	VENTRICULAR-PACED COMPLEXES	other complexes also detected
(VPACEF)	AB	AFIB/FLUTTER AND VENTRICULAR-PACED RHYTHM	V-paced rhythm, A-rate>240
(VPACFD)	AB	AFIB/FLUT, V-PACED COMPLEXES WITH INHIBITION	non-paced complexes, A-rate>240
(VPC)	ON	VENTRICULAR PREMATURE COMPLEX	V complex w/ short R-R interval
(VPE)	AB	VENTRICULAR PREEXCITATION	Delta waves
(VPEL)	AB	VENT PREEXCITATION, LEFT ACCESSORY PATHWAY	Delta wave & initial axis(30,120)

(VPELA)	NS	VENTRICULAR PREEXCITATION, A LEFT ANTEROSEPTAL ACCESSORY PATHWAY	
(VPELL)	NS	VENTRICULAR PREEXCITATION, A LEFT LATERAL ACCESSORY PATHWAY	
(VPELP)	NS	VENTRICULAR PREEXCITATION, A LEFT POSTEROSEPTAL ACCESSORY PATHWAY	
(VPER)	AB	VENT PREEXCITATION, RIGHT ACCESSORY PATHWAY	Delta wave & initial axis(-60,29)
(VPERA)	NS	VENTRICULAR PREEXCITATION, A RIGHT ANTEROSEPTAL ACCESSORY PATHWAY	
(VPERL)	NS	VENTRICULAR PREEXCITATION, A RIGHT LATERAL ACCESSORY PATHWAY	
(VPERP)	NS	VENTRICULAR PREEXCITATION, A RIGHT POSTEROSEPTAL ACCESSORY PATHWAY	
(VSVPC)	NS	PREMATURE COMPLEX, VENT OR ABERRANT SUPRAVENT	
(VTACH)	AB	EXTREME TACHYCARDIA WITH WIDE COMPLEX, NO FURTHER RHYTHM ANALYSIS ATTEMPTED	
(VTRI)	AB	VENTRICULAR TRIGEMINY	trigeminy string>6 w/ V complexes
(VVI)	NS	RHYTHM CONSISTENT WITH VVI PACING	

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