Always here to help you

Register your product and get support at
www.philips.com/welcome

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

Introduction

Congratulations on your purchase and welcome to Philips! To fully benefit from the support that Philips offers, register your product at www.philips.com/welcome.

Light-based hair removal is one of the most effective methods to treat unwanted body hair with longer-lasting results. It is different from current home hair removal methods because it significantly reduces hair regrowth. Philips has been investigating this method together with leading dermatologists for over 15 years. You can now use this innovative technology comfortably in the privacy of your own home, whenever you want.

This appliance has been specifically designed and optimised for use by women, but can be used by men as well.

Philips Lumea Essential enables treatment of unwanted hair on all areas below the cheekbones, including the face (upper lip, chin and cheeks), using two tailored attachments.

Philips Lumea Essential is designed for individual use only. For hygienic reasons and to ensure a longer life time, we advise you not to share the appliance with others.

For more information about Philips Lumea Essential, see www.philips.com/lumea. Check our homepage regularly for updates.

Benefits

This appliance provides the following benefits:

Intense Pulsed Light (IPL) technology for home use

Philips Lumea Essential uses a light-based technology called 'Intense Pulsed Light' (IPL). IPL is also used in the professional beauty market for hair removal and has been in use for over fifteen years. Philips Lumea Essential now brings this innovative technology to the comfort of your home.

Effective hair regrowth prevention for everyday smoothness

This method provides you with long-lasting hair regrowth prevention and silky-smooth skin. Our studies have shown a significant hair reduction after only two treatments. Optimal results were generally achieved after four to five treatments. A reduction in hair density is visible from treatment to treatment. How fast the effect becomes visible and how long the results last varies from person to person. For more information, see section 'How to achieve optimal results' in chapter 'Using the appliance'.

Gentle treatment, even on sensitive body areas


Philips Lumea Essential has been developed in close cooperation with licensed dermatologists for an effective and gentle treatment, even on the most sensitive body areas. The precision attachment even allows you to safely use the appliance on your face.

No replacement lamps needed, no hidden costs

Philips Lumea Essential is equipped with a high-performance lamp which does not require replacement. This lamp can generate over 100,000 flashes and should enable you to continuously maintain results for more than 5 years* while saving money on expensive replacement lamps.

*Based on average usage scenario and recommended treatment regime for one user.

How Philips' Intense Pulsed Light technology works

Hair growth

Hair growth differs from person to person, depending on age, metabolism and other factors. However, everyone's hair normally grows in 3 phases.

1 Growing phase (anagen phase)



The hair grows actively from the root. In this phase the concentration of melanin is highest. Melanin is responsible for the colour of the hair: Melanin concentration is also the determining factor for the efficacy of the Intense Pulsed Light technology. Only hairs in the growing phase are susceptible to treatment with light.

2 Degradation phase (catagen phase)

Hair growth stops and the root shrinks before the hair sheds.


3 Resting phase (telogen phase)

The old hair is shed. The follicle remains at rest until the biological clock tells it to become active again and grow a new hair.

Working principle

- The appliance works by heating up the hair and root beneath the skin. The melanin in the hair and hair root absorbs the energy from the light pulses. The darker the colour of the hair, the more light will be absorbed. This process stimulates the hair to go into the resting phase.
- In between treatments, hair sheds naturally and regrowth is inhibited. After the treatment, it can take one or two weeks for the hair to shed.

- Optimal hair reduction cannot be obtained with a single session, since only hairs in the growing phase are susceptible to treatment with light. During a treatment, hairs on the treated area are usually in different stages of the hair growth cycle. To effectively treat all hairs you have to repeat the treatment once every two weeks for the first two months.

An integrated optical filter ensures that no UV light reaches the skin and ensures that the appliance is safe for the skin when used according to the guidelines.

General description (Fig. 5)

- 1 Intensity lights (1- 5)
- 2 On/off button
- 3 'Ready to flash' light
- 4 Air vents
- 5 Appliance socket
- 6 Flash button
- 7 Openings for electronic contacts
- 8 Light exit window with integrated UV filter

- 9 Body attachment (suitable for all body areas below the neck)
 - a Electronic contacts
 - b Metallic frame inside the body attachment
 - c Integrated safety system
 - d Skin tone sensor
- 10 Precision attachment (suitable for all body areas, including upper lip, chin and cheeks)
 - a Filter glass
 - b Integrated safety system
 - c Skin tone sensor

11 Adapter

12 Small plug

Also included:

- pouch
- cleaning cloth

Note: *The main filter is integrated in the appliance. It ensures safe use on body areas below the neck. Therefore the body attachment does not need an extra filter. However, the precision attachment has an additional filter for safe and pleasant treatment of sensitive skin, such as the face.*

Important

Read this user manual carefully before you use the appliance and save it for future reference.

Danger

- Water and electricity are a dangerous combination. Do not use this appliance in wet surroundings (e.g. near a filled bath, a running shower or a filled swimming pool).
- Keep the appliance and the adapter dry.
- If the appliance is broken, do not touch any inner part to avoid electric shock.
- Never insert any objects into the appliance.

Warning

- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the appliance.
- The appliance is not intended for children under the age of 15 years. Teenagers aged between 15 and 18 years can use the appliance with the consent and/or assistance of their parents or the persons who have parental authority over them. Adults from 19 years and older can use the appliance freely.
- Only use the adapter supplied.
- Do not use the appliance or the adapter if it is damaged.

- The adapter contains a transformer. Do not cut off the adapter to replace it with another plug, as this causes a hazardous situation.
- If the adapter is damaged, always have it replaced with one of the original type in order to avoid a hazard.
- Do not use the appliance if the glass of the light exit window is broken.
- Do not use the precision attachment if its reddish glass is broken.

To prevent damage

- Make sure nothing obstructs the air flow through the ventilation slots.
- Never subject the appliance to heavy shocks and do not shake or drop it.
- If you take the appliance from a very cold environment to a very warm environment or vice versa, wait approximately 3 hours before you use it.
- Let the appliance cool down after use and then store it in the pouch supplied, so that it does not become dusty.
- Never leave the appliance unattended when it is switched on. Always switch off the appliance after use.
- Do not expose the appliance to direct sunlight or UV light for extended periods, to prevent damage.
- Do not flash against any other surface than the skin. This can cause severe damage to the attachments and/or the light exit window. Only flash when the appliance is in contact with skin.

Caution

- Never use the appliance around the eyes or near the eyebrows to avoid the risk of eye damage.
- Use this appliance only at settings suitable for your skin type (see the table in section 'Recommended light intensities' in chapter 'Preparing for use'). Use at higher settings than those recommended can increase the risk of skin reactions and side effects.
- This appliance is only intended for removing unwanted human body hair from areas below the cheekbones with the two attachments supplied. Do not use it for any other purpose.
- Make sure the skin in the area to be treated is clean, dry and free from creams, lotions and other cosmetics. If the skin is not clean and dry, the integrated safety system may not work as intended and the appliance may not flash or may flash unintentionally.
- Do not use the appliance if any of the situations mentioned in section 'Contraindications' apply to you.

Keep the following things in mind to avoid compromising the lifetime of your Philips Lumea Essential:

- Do not expose the appliance to ambient temperatures lower than 15°C or higher than 35°C during use.
- The appliance is equipped with a built-in overheat protection. If the appliance overheats, it does not flash. Do not switch off the appliance, but let it cool down for approx. 15 minutes before you continue with your treatment. After 5 minutes the appliance switches off automatically. Then let it cool down for another 10 minutes.
- Always return the appliance to a service centre authorised by Philips for examination or repair. Repair by unqualified people could cause an extremely hazardous situation for the user.

Electromagnetic fields (EMF)

This Philips appliance complies with all applicable standards and regulations regarding exposure to electromagnetic fields.

Contraindications

For whom is Philips Lumea Essential NOT suitable?

Philips Lumea Essential is not suitable for everyone. If any of the following is true for you, then this appliance is **not** suitable for you to use!

Do not use the appliance if you are pregnant or breast feeding.

Philips Lumea Essential has not been tested on pregnant women, therefore we advise you not to use the appliance if you are pregnant or breast feeding. Hormonal changes could increase the sensitivity of your skin and the risk of injury to the skin.

Never use the appliance if you have an unsuitable skin and/or hair colour, i.e.:

- If you are of skin type V (you rarely get sunburnt, you tan very well and your skin tone is dark brown) or if you are of skin type VI (you rarely or never get sunburnt and you tan very darkly; your skin tone is brownish black or darker); in these cases you run a high risk of developing skin reactions such as discolourations, redness or burns, when you use Philips Lumea Essential. See also chapter 'Possible side effects and skin reactions' and consult the 'Recommended light intensities table' in the quick start guide.
- If your natural hair colour in the treatment areas is light blond, white, grey or red. The treatment will not be effective.

Do not use Philips Lumea Essential to treat any area more often than once every two weeks.

Never use the appliance if you suffer from any of the diseases listed below, i.e.:

- If you have a skin disease such as active skin cancer; you have a history of skin cancer or any other localised cancer in the areas to be treated.
- If you have pre-cancerous lesions or multiple atypical moles in the areas to be treated.
- If you have a history of collagen disorder; including a history of keloid scar formation or a history of poor wound healing.
- If you have a history of vascular disorder; such as the presence of varicose veins or vascular ectasia in the areas to be treated.
- If your skin is sensitive to light and easily develops a rash or an allergic reaction.
- If you have infections, eczema, burns, inflammation of hair follicles, open lacerations, abrasions, herpes simplex, wounds or lesions and haematomas in the areas to be treated.
- If you have had surgery in the areas to be treated.
- If you have epilepsy with flashlight sensitivity.
- If you have diabetes, lupus erythematoses, porphyria or congestive heart disease.
- If you have any bleeding disorder.
- If you have a history of immunosuppressive disease (including HIV infection or AIDS).

Never use the appliance if you take any of the medications listed below, i.e.:

- If your skin is currently being treated with or has recently been treated with Alpha-Hydroxy Acids (AHAs), Beta- Hydroxy Acids (BHAs), topical isotretinoin and azelaic acid.
- If you have taken any form of isotretinoin such as Accutane® or Roaccutane® in the last six months. This treatment can make skin more susceptible to tears, wounds and irritations.
- If you are on painkillers, which reduce the skin's sensitivity to heat.
- If you are taking photosensitising agents or medications, check the package insert of your medicine. Never use the appliance if the medication instructions include cautions regarding photo-allergic reactions, photo-toxic reactions or if you must avoid sun when taking this medicine.
- If you take anticoagulation medications, including heavy use of aspirin, in a manner that does not allow you to stop using them at least one week prior to each treatment with Philips Lumea Essential.
- If you take immunosuppressive medications.

Never use the appliance on the following areas:

- Around the eyes and near the eyebrows.
- On nipples, areolas, labia minora, vagina, anus and the inside of the nostrils and ears. See chapter 'Using the appliance', section 'Using the appliance in the bikini area'.
- Men must not use it on the face or on the scrotum. Men's facial hair differs from women's facial hair and the precision attachment is only designed for treatment on women's facial hair.
- Over or near anything artificial like silicone implants, pacemakers, subcutaneous injection ports (insulin dispenser) or piercings.
- On moles, freckles, large veins, darker pigmented areas, scars or other skin anomalies without consulting your doctor. This can result in a burn and a change in skin colour, which makes it potentially harder to identify skin-related diseases (e.g. skin cancer).
- On warts, tattoos or permanent make-up. This can result in a burn and a change in skin colour (hypopigmentation or hyperpigmentation). Usage on warts may even spread the warts, as warts are viral infections.
- On areas, where you use long-lasting deodorants. This can result in skin reactions. See chapter 'Possible side effects and skin reactions'.

Never use the appliance on sunburnt, recently tanned (last 2 weeks) or fake-tanned skin.

- All types of tanning can influence the safety of the treatment. This applies to exposure to natural sunlight as well as artificial tanning methods, such as tanning lotions, tanning appliances etc. If you are tanned or are planning to tan, see chapter 'Preparing for use', section 'Tanning advice'.

Note: This list is not exhaustive. If you are not sure whether you can use the appliance we advise you to consult a doctor.

Possible side effects and skin reactions

Common skin reactions

- Your skin may show slight redness and/or may prickle, tingle or feel warm. This reaction is absolutely harmless and disappears quickly.
- A skin reaction similar to sunburn may occur on the treated area for a short while. If this does not disappear within 3 days, we advise you to consult a doctor.
- Dry skin and itching may occur due to shaving or a combination of shaving and light treatment. This is harmless and disappears within a few days. You can cool the area with an ice pack or a wet facecloth. If dryness persists, you can apply a non-scented moisturiser on the treated area 24 hours after the treatment.

Rare side effects


- Burns, excessive redness and swelling: these reactions occur rarely and are usually the result of using a light intensity that is too high for your skin colour. If these reactions do not disappear within 3 days, we advise you to consult a doctor. Wait with the next treatment until the skin has healed completely and make sure you use a lower light intensity.
- Skin discolouration: this occurs very rarely. Skin discolouration manifests itself as either a darker or lighter patch than the surrounding area. This is the result of using a light intensity that is too high for your skin colour. If the discolouration does not disappear within 2 weeks, we advise you to consult a doctor. Do not treat discoloured areas until the discolouration has disappeared and your skin has regained its normal skin colour.
- Skin infection and inflammation: This occurs very rarely and may be caused by using the appliance on wounds or cuts that resulted from shaving, on already existing wounds or on ingrown hairs.
- Excessive pain: this can occur during or after treatment if you have used the appliance on unshaved skin, if you use the appliance at a light intensity that is too high for your skin colour, if you flash the same area more than once and if you use the appliance on open wounds, inflammations, infections, tattoos, burns, etc. See chapter 'Contraindications' for more information.

After use

- Do not apply any perfumed cosmetics to the areas treated immediately after treatment.
- Do not use deodorant immediately after treating the underarms. Wait until any skin redness has disappeared completely.
- See chapter 'Preparing for use', section 'Tanning advice' for information about sun exposure and artificial tanning.

Preparing for use

For optimal results, prepare the areas on which you intend to use the appliance in the manner described below.

Preparing treatment areas

- 1 Shave the areas you intend to treat as long as these areas still show hair regrowth.

Note: Do not use depilatory creams instead of shaving.

Note: If your last depilation method involved removing hairs by the roots (e.g. epilating), you can start the treatment right away. In this case, we advise you to shave in between treatments for optimal results.

Please note that treatment on unshaved or improperly shaved areas has a number of undesired effects:

- The light exit window and attachments may become contaminated with stray hairs and dirt. A contaminated light exit window and contaminated attachments can reduce the lifetime and efficacy of the appliance. Contamination of the light exit window can also cause the treatment to become painful or cause skin reactions like redness and discolouration.
- Hairs on your skin can burn, which can result in an unpleasant smell.

- 2 Clean your skin and make sure it is hair-free, entirely dry and free from oily substances (e.g. residues of deodorants, lotions, perfumes, tanning creams and sunscreens).

- 3 Check the appliance, especially the light exit window and attachments, for trapped hair, dust and fibres. If you see any contamination, clean the appliance according to the instructions in chapter 'Cleaning and maintenance'.

If shaving causes wounds or lesions, your skin may be more sensitive to IPL than that of others. This may cause the IPL treatment to feel uncomfortable or it may lead to other unwanted side effects. We advise you to wait until the irritation has disappeared.

- 4 Once the hairs stop growing back, which usually happens after 4-5 treatments, you no longer have to shave before you use the appliance.

Light intensity

Unlike professional light-based appliances, Philips Lumea Essential has been developed for consumer use and its light pulses are gentler. Therefore you can expect a **convenient** and **effective** treatment at a light intensity that you find comfortable.

Philips Lumea Essential provides five different light intensities. Adjust the light intensity to the colour of your skin, body hair and facial hair and to a level you find comfortable. The recommended light intensities table below helps you to establish the right light intensity.

Note: The appliance is automatically disabled when your skin is too dark, to prevent you from developing skin reactions. However, the integrated skin colour sensor does not tell you which light intensity setting to use, as preferences differ strongly from person to person. Philips Lumea Essential gives you the freedom to select the light intensity that you find most convenient, as long as the appliance is suitable for your skin colour. To establish the right light intensity for use on your skin, see the recommended light intensities table below and on the quick start guide.

- 1** Consult the table below to determine which light intensities are appropriate for your skin and the colour of your skin, body hair and facial hair (if there is no suitable light intensity for your skin/hair colour combination, this is indicated with x in the table).

Tip: You can put the table in the quick start guide on your skin and compare it with your skin colour.

Tip: We advise you to store the quick start guide in the pouch for future reference.

- 2** Reduce the light intensity if you experience any pain or discomfort. Persons with darker skin run a higher risk of skin irritation because their skin absorbs more light. Therefore lower intensities are recommended for persons with darker skin.

The light pulse can feel warm or hot on the skin but should never be painful. Different body areas may require different light intensities, see the table below. To prevent skin reactions, see chapter ‘Possible side effects and skin reactions’.

Recommended light intensities (1-5)

In the table below, the numbers indicate the light intensity to be selected and ‘x’ means that the appliance is not suitable for you.

Note: The colour of your skin may change, for example because your skin has tanned. It may then be darker than during the previous treatment. Therefore it is important that you check the setting to be used before every treatment and per area. Some areas may be darker than others.

Body hair colour:	Dark blond/light brown/ brown/dark brown/black	white/grey/red/light blond
<i>Skin colour and characteristics:</i>		
White (always sunburn, no tanning)	4/5	x
Beige (easy sunburn, minimal tanning)	4/5	x
Light brown (sometimes sunburn, good tanning)	4/5	x
Mid brown (rarely sunburn, fast and good tanning)	1/2/3	x
Dark brown (rarely sunburn, very good tanning)	x	x
Brownish black and darker (rarely to never sunburn, very dark tanning)	x	x

Important

To check whether Philips Lumea Essential is suitable for you, see chapter ‘Contraindications’.

Tanning advice

- **Avoid sunbathing when you use Philips Lumea Essential regularly (every two weeks).**

Sunbathing can cause skin to darken (tan) or can cause sunburn (reddening of the skin). Using Philips Lumea Essential on darkened (tanned) and/or sunburned skin increases the risk of adverse side effects and skin reactions. See chapter 'Possible side effects and skin reactions'.

- **Sunbathing before treatment: Wait 2 weeks and perform a skin test.**

Wait at least 2 weeks after sunbathing before you use the appliance. Perform a skin test (see chapter 'Using the appliance', section 'First use and skin test') on the area to be treated to determine the appropriate light intensity. Failure to do so increases the risk of developing skin reactions after treatment with Philips Lumea Essential.

- **Sunbathing after treatment: Wait 48 hours before sunbathing after treatment with Philips Lumea Essential.**

Wait at least 48 hours after treatment before you expose treated areas to the sun. Even after 48 hours, make sure that the treated skin shows no adverse reaction to the treatment before you expose it to the sun. Cover treated areas when you go out into the sun or use a sunblock (SPF 30+) in the two weeks after treatment. Exposing treated areas to sunlight without protection immediately after treatment increases the risk of adverse side effects and skin reactions. See chapter 'Possible side effects and skin reactions'.

- **Tanning with artificial light**

Follow the instructions in 'Sun exposure before treatment' and 'Sun exposure after treatment' above.

- **Tanning with creams**

If you have used an artificial tanning lotion, wait until the artificial tan has fully disappeared before you use Philips Lumea Essential.

Note: For more details on unwanted skin reactions, see chapter 'Possible side effects and skin reactions'.


Using the appliance

Check the recommended light intensities table in chapter 'Preparing for use', section 'Light intensity' to determine which light intensities are most suitable for your skin and body hair colour. Use the appliance in a well-lit room; this reduces the perceived brightness of the flash.

Shave the area you intend to treat before you use the appliance (see chapter 'Preparing for use', section 'Preparing treatment areas').

If the treatment becomes painful at any point in time, reduce the light intensity setting until the treatment is comfortable again. If you notice skin reactions stronger than slight redness, stop the treatment immediately and read section 'Possible side effects and skin reactions' in this chapter.


Choosing the right attachment

- Body attachment

The body attachment can be used to treat unwanted body hair on all body areas below the neck (for exceptions, see chapter 'Contraindications'). The body attachment can cover an area of 4cm² and has been specially designed for faster treatment of larger areas like legs.

Never use the body attachment on the face.


- Precision attachment

The precision attachment includes an additional integrated light filter and can be used for safe treatment of unwanted facial hair on the upper lip, chin and cheeks. It may also be used on other sensitive areas of the body, such as the underarms and bikini area or hard-to-reach areas.

Never use the appliance around the eyes or near the eyebrows.


First use and skin test


When you use the appliance for the first time or after recent tanning, we advise you to perform a skin test on **each** area to be treated to check your skin's reaction to the treatment and to get accustomed to light-based treatment.

- 1** Choose the attachment appropriate for the area you intend to treat.
- 2** Choose an area close to the area you intend to treat.
- 3** Insert the small plug into the appliance and put the adapter in the wall socket.

Note: Make sure the selected attachment is placed on the appliance.


- 4** Press the on/off button to switch on the appliance.

► Intensity light 1 goes on.

Note: The appliance always automatically selects the lowest light intensity setting. See section 'Light intensity' in chapter 'Preparing for use' to determine the appropriate light intensity.


- 5** To adjust the light intensity setting, press the on/off button one or more times until you have reached the required setting. The corresponding intensity light lights up.


- 6** Place the appliance at a 90° angle on the skin so that the attachment is in contact with the skin. Press the appliance firmly onto the skin to ensure skin contact.

Tip: For easier use, stretch the area of skin you want to treat.


- ▶ The appliance automatically performs a skin colour analysis.

- When the 'ready to flash' light on the back of the appliance lights up white, your skin is suitable for treatment with the appliance and you can proceed.


- When the 'ready to flash' light first flashes orange and then lights up orange continuously, your skin is not suitable for treatment with the appliance and the appliance is disabled. You might be able to use Philips Lumea Essential on less tanned or naturally lighter body parts. If this does not work, please contact Philips Consumer Care.

- 7** Press the flash button to release a flash.

*Note: The visible light produced by the appliance is the reflection of the flash on the skin and is harmless to your eyes. It is not necessary to wear goggles during use. Use the appliance in a **well-lit** room; this reduces the perceived brightness of the flash.*

- 8** Release the flash button once the appliance has flashed.

- ▶ The 'ready to flash' light starts blinking white to indicate that the appliance is getting ready to release a flash.
- ▶ When the appliance is ready to release a flash, the 'ready to flash' light lights up white continuously.

- 9** Slide the appliance over the skin to the next area to be treated. As long as it feels comfortable (not painful), increase the setting by one level within the recommended range for your skin type and apply one flash for each setting. **Do not apply more than one flash on the same spot.**

- 10** When you have finished the treatment, press and hold the on/off button for 2 seconds to switch off the appliance.
- 11** After the skin test, wait for 24 hours and check your skin for any reaction. If your skin shows reactions, choose the highest setting that did not result in any skin reaction for subsequent use.

Subsequent use

- 1** Before every treatment, clean the light exit window, the outside of the attachments, the metallic frames inside the attachments and the reddish filter glass inside the precision attachment. For more information see chapter 'Cleaning and maintenance'.
- 2** Start each treatment by selecting the light intensity very carefully, based on your skin's sensitivity and on how comfortable treatment feels. See section 'Recommended light intensities' to select the right setting for your skin type and hair colour combination.

Your skin may react differently on different days/occasions for a number of reasons. See chapter 'Possible side effects and skin reactions' in this chapter for more information.

- 3** Follow steps 2-7 of the previous section ('Before first use and skin test').

Note: It is normal that the appliance and the attachment become warm during use.

- 4** Slide the appliance over the skin to the next area to be treated and repeat steps 5 and 6 of section 'First use and skin test' in this chapter. After each flash, it takes up to 3.5 seconds until the appliance is ready to flash again. The 'ready to flash' light on the back of the appliance starts blinking white again to indicate that it is charging. Wait until this light lights up continuously before you continue.


Tip: For a more efficient treatment, keep the flash button pressed while you slide the appliance over your skin.

Do not treat exactly the same area multiple times during one session. This does not improve the effectiveness of the treatment, but increases the risk of skin reactions.

- 5** To avoid untreated areas, make sure the flashes are produced close to each other. Note that effective light only comes out of the light exit window. However, always make sure you do not flash the same area twice.
- 6** When you have finished the treatment, press and hold the on/off button for 2 seconds to switch off the appliance.
- 7** Check the light exit window and the attachments for hairs and dirt and clean the appliance after use (see chapter 'Cleaning and maintenance').

Treatment plan and effectiveness

A typical full hair growth cycle may take 18-24 months. During this time, multiple treatments with Philips Lumea Essential may be required in order to achieve permanent hair reduction. The effectiveness of hair reduction varies from person to person according to body area, hair colour and how the Philips Lumea Essential is used.

How to achieve optimal results

- It is normal that hairs are still visible after the first few treatments (see chapter 'Introduction', section 'Working principle' for more information).
- To remove all hairs successfully and prevent the hair follicle from becoming active again, the treatments have to be repeated once every two weeks for the first two months. Your skin should be smooth and hair-free by then. Please read section 'Maintenance phase' to find out how you can keep these results.

Maintenance phase

- To ensure that your skin stays smooth and hair-free, we advise you to repeat the treatment every four to six weeks. The time between treatments may vary based on your individual hair regrowth and also across different body areas. If too many hairs grow back between treatments, simply shorten the treatment interval again, but do not treat any area more often than once every two weeks. When you use the appliance more frequently than recommended, you do not increase the effectiveness but you do increase the risk of skin reactions.

Guidelines for treatment

Average treatment time

Area	Approx. treatment time	Recommended attachment
Underarms	1 minute per underarm	Body/Precision attachment
Bikini line	1 minute per side	Body/Precision attachment
Entire bikini area	4 minutes	Body/Precision attachment
One lower leg	8 minutes	Body attachment
One full leg	14 minutes	Body attachment
Upper lip/chin/cheeks	1 minute	Precision attachment

Note: The average treatment times above have been observed during research and tests. Individual treatment times may vary from person to person.

Intended treatment areas


The appliance is intended to treat unwanted hairs on the upper lip, chin, cheeks, underarms, arms, stomach, bikini area, back and legs.

Using the appliance on the legs

Check the recommended light intensities table in chapter 'Preparing for use', section 'Light intensity' to see which settings are suitable for your legs.

- You can move the appliance up and down or across the leg, whichever is more convenient to you. In order to maintain full contact with the skin, you can also place and hold the appliance on your skin sideways.

Note: Bony areas like shinbone and ankles are more sensitive than other areas of the body. When you press the attachment hard onto these areas, the skin can be more susceptible to skin reactions, such as skin irritation and redness. See chapter 'Possible side effects and skin reactions'. Press the attachment onto the skin just hard enough to make the 'ready to flash' light go on, but not harder.

Using the appliance in the bikini area

Check the recommended light intensities table in chapter 'Preparing for use', section 'Light intensity' to see which settings are suitable for your bikini area. In principle, Philips Lumea Essential is safe and gentle enough for treatment of the entire bikini area (e.g. for 'Brazilian' or 'Hollywood' style). Use a lower setting on darker pigmented, more sensitive skin in the bikini area. Do not use the appliance on inner labia, vagina and anus. Men who wish to treat their pubic area must not use the appliance on the scrotum.

Note: If you want to leave some hair in place, e.g. in case of a 'Brazilian', make sure that the skin colour sensor is not placed on the remaining hair during the treatment, as this can cause the appliance to block.

Note: For hygienic reasons, always clean the attachment before you start treating another body area (see chapter 'Cleaning and maintenance').

Using the appliance on the underarms

Check the recommended light intensities table in chapter 'Preparing for use', section 'Setting the light intensity' to see which settings are suitable for your underarms.

Make sure you remove any deodorant, perfume, cream or other skin care product from your underarms before treatment.

Tip: To achieve optimal results on your underarms, stand in front of a mirror to get a good view of the area to be treated and the 'ready to flash' light. Raise your arm and place your hand on the back of your neck to stretch the skin in the underarm area.

Note: The reflected light does not harm your eyes. If the scattered light is unpleasant, close your eyes when you release a flash.

Using the appliance on the face (upper lip, chin and cheeks)

Check the recommended light intensities table in chapter 'Preparing for use', section 'Setting the light intensity' to see which settings are suitable for your face (upper lip, chin and cheeks).

Always use only the **precision attachment** to treat the upper lip, chin and cheeks. Never use the appliance around the eyes or near the eyebrows. When you treat your upper lip, chin or cheeks, make sure you do not expose your face to direct sunlight in the 24 hours following the treatment. Apply a sunblock (SPF 30+) after the treatment for the next two weeks when you expose your face to direct sunlight.

Caution: When you treat the area above your upper lip, be careful that you do not flash on the lip.

Tip: Apply a thick layer of white eyeliner to your lip and lip margin. This white line ensures the energy of the flash is reflected and does not touch your lip if you flash the lip unintentionally. Only use a white eyeliner. Any other colour absorbs the light and may cause skin irritation.

Tip: To achieve optimal results on your upper lip, chin and cheeks, stand in front of a mirror to get a good view of the area to be treated and the 'ready to flash' light.

Tip: Because of the contours of the jaw, it may be difficult to bring the integrated safety system into full contact with the skin and to release a flash. Put your tongue between your teeth and upper lip or puff up your cheeks to make this easier.

Tip: The reflected light does not harm your eyes. If the scattered light is unpleasant, close your eyes when you release a flash.

Cleaning and maintenance

To ensure optimal results and a long life for the appliance, clean the appliance before and after each use and, if necessary, also during the treatment. The appliance loses its effectiveness if you do not clean it properly.

Never clean the appliance or any of its parts under the tap or in the dishwasher.

Never use scouring pads, abrasive cleaning agents or aggressive liquids such as petrol or acetone to clean the appliance.

Never scratch the light exit window or the metallic frame inside the body attachment.

Note: Stop using the appliance when it is no longer possible to clean the light exit window or the attachments.


- 1** Switch off the appliance, unplug it and let it cool down.


Note: The light exit window becomes hot during use. Make sure it has cooled down before you clean it.

- 2** To remove the attachment, grab the projections at the top and bottom of the attachment and pull gently.
- 3** If necessary, clean the outside of the appliance with the soft cleaning cloth supplied with the appliance.
- 4** Moisten the soft cloth supplied with the appliance with a few drops of water and use it to clean the following parts:
 - the light exit window
 - the outside surface of the attachments
 - the metallic frames inside the attachments
 - the reddish filter glass inside the precision attachment

Make sure that the electronic contacts do not get wet.

Note: If water does not clean effectively, use a few drops of high-percentage alcohol to clean the parts mentioned above.

Note: Make sure you clean the appliance before and after each use and, if necessary, also during the treatment according to the instructions in this section.

Storage

- 1** Switch off the appliance, unplug it and let it cool down.
- 2** Clean the appliance before you store it.
- 3** Store the appliance and the attachments in the pouch supplied at a temperature between -20°C and $+70^{\circ}\text{C}$.

Tip: We advise you to store the quick start guide in the pouch for future reference.

Ordering accessories

To buy accessories or spare parts, visit **www.shop.philips.com/service** or go to your Philips dealer. You can also contact the Philips Consumer Care Centre in your country (see the worldwide guarantee leaflet for contact details).

Spare parts

The following spare parts are available:

- Body attachment (type number 4222.100.1629)
- Precision attachment (type number 4222.100.1630)
- Adapter
- Cleaning cloth
- Pouch

Environment

- Do not throw away the appliance with the normal household waste at the end of its life, but hand it in at an official collection point for recycling. By doing this, you help to preserve the environment.

Guarantee and support

If you need information or support, please visit **www.philips.com/support** or read the separate worldwide guarantee leaflet.

Technical specifications

Model SC1996

Rated voltage	100V-240V
Rated frequency	50Hz-60Hz
Rated input	36W
Protection against electric shock	Class II 
Protection rating	IP 30 (EN 60529)
Operating conditions	Temperature: +15° to +35°C Relative humidity: 30% to 95%
Storage conditions	Temperature: -20°C to +70°C Relative humidity: 5% to 95%
Spectrum of body attachment	>500nm
Spectrum of precision attachment	>600nm

Adapter

Model: S036NU2400150

Rated	100-240V, 50/60Hz, 100mA
Protection against electric shock	Listed Class 2 power supply Non-suitable for use in the presence of flammable anesthetics with air, oxygen and nutritious oxide. Normal temperature and humidity for storage and transportation. Method of continuous operation with intermittent load.

- Do not expose to excessive moisture or rain.

- Protection against electric shock: Type BF Applied Part.

- Fragile device.

Troubleshooting

This chapter summarises the most common problems you could encounter with the appliance. If you are unable to solve the problem with the information below, visit www.philips.com/support for a list of frequently asked questions or contact the Consumer Care Centre in your country.

Problem	Possible cause	Solution
The adapter/ appliance becomes warm during charging.	This is normal.	No action required.
I have switched on the appliance but I cannot increase or decrease the light intensity.	The appliance has to be reset.	To reset the appliance, switch off the appliance and then switch it on again. If you still cannot adjust the light intensity, contact the Consumer Care Centre in your country, your Philips dealer or a Philips service centre.
The 'ready to flash' light does not work when I place the appliance on the skin.	You have not properly placed the appliance on the skin.	Place the appliance at a 90° angle onto the skin in such a way that all 4 sides of the integrated safety system touch the skin. Then check if the 'ready to flash' light lights up white continuously and press the flash button. If this does not work, place the appliance on a part of your body where it is easy to establish full skin contact, for example on your lower arm. Then check if the 'ready to flash' light works. If the 'ready to flash' light still does not start blinking when you place the appliance on the skin, contact the Consumer Care Centre in your country, your Philips dealer or a Philips service centre.

Problem	Possible cause	Solution
The 'ready to flash' light first blinks orange and then lights up orange continuously when I place the appliance on the skin and the appliance does not release a flash.	The skin you are trying to treat is too dark.	To prevent your skin from developing side effects, the appliance does not work on skin that is too dark. You might be able to use Philips Lumea Essential on lighter coloured body parts. If this does not work either, contact the Consumer Care Centre in your country, your Philips dealer or a Philips service centre.
The 'ready to flash' light does not work when I place the attachment on my skin, but the cooling fan works.	The overheat protection has been activated.	When the overheat protection has been activated, the fan still works. Do not switch off the appliance. Let it cool down for approx. 15 minutes before you continue to use it. After 5 minutes the appliance switches off automatically. Then let it cool down for 10 more minutes. If the 'ready to flash' light still does not work when you place the appliance on the skin, contact the Consumer Care Centre in your country, your Philips dealer or a Philips service centre.
The 'ready to flash' light lights up continuously, but the appliance does not produce a flash when I press the flash button.	The appliance needs to be reset.	To reset the appliance, switch off the appliance and switch it on again. If the 'ready to flash' light lights up continuously, but the appliance still does not produce a flash when you press the flash button, contact the Consumer Care Centre in your country, your Philips dealer or a Philips service centre.
Intensity light 1 continues to blink and I cannot change the setting.	The attachment is not placed on the appliance properly.	Place the attachment on the appliance properly.
All intensity lights go on alternately.	The appliance needs to be reset.	To reset the appliance, switch off the appliance and switch it on again. If the intensity lights continue to go on alternately, contact the Consumer Care Centre in your country, your Philips dealer or a Philips service centre.
All intensity lights are blinking simultaneously.	The appliance is defective.	Contact the Consumer Care Centre in your country, your Philips dealer or a Philips service centre.
The appliance produces a strange smell.	The attachment and/or the filter glass is dirty.	Clean the attachment and the light exit window carefully. If you are unable to clean the attachment properly, contact the Consumer Care Centre in your country, your Philips dealer or a Philips service centre to replace the attachment.
	You have not shaved the area to be treated properly.	If there are hairs on the area to be treated, these hairs may get burnt when you use the appliance. As a result, you notice a strange smell. Shave the area to be treated properly before you use the appliance. If shaving causes skin irritation, trim hairs as short as possible and use the Philips Lumea Essential at a comfortable setting.

Problem	Possible cause	Solution
	You stored the appliance in a dusty environment.	Clean the appliance and the attachments (see chapter 'Cleaning and maintenance').
The skin feels more sensitive than usual during treatment.	The light intensity you use is too high.	Check if you have selected the right light intensity. If necessary, select a lower light intensity.
	You did not shave the areas to be treated.	Shave the areas to be treated before you use the appliance. If shaving causes skin irritation, trim hairs as short as possible and use the Philips Lumea Essential at a comfortable setting.
	The light exit window is broken.	If the light exit window is broken, do not use the appliance anymore, to avoid the risk of skin reactions. Contact the Consumer Care Centre in your country, your Philips dealer or a Philips service centre.
I experience an unacceptable pain sensation when I use the appliance.	You did not shave the areas to be treated	Shave the areas to be treated before you use the appliance. If shaving causes skin irritation, trim hairs as short as possible and use the Philips Lumea Essential at a comfortable setting.
	You treated an area for which the appliance is not intended.	You must not use the appliance on your inner labia, vagina, anus, nipples and areolas. Men must not use it on the face and on their scrotum. Men's facial hair differs from women's facial hair and the precision attachment is only designed for treatment of women's facial hair.
	You have used a light intensity that is too high for you.	Reduce the light intensity to a level that is comfortable for you. See chapter 'Preparing for use', section 'Light intensity'.
	The light exit window of the body attachment or the filter glass of the precision attachment is broken.	If the light exit window of the body attachment or the filter glass of the precision attachment is broken, do not use the appliance with the attachments anymore. Go to your dealer or a Philips service centre or contact the Consumer Care Centre to have the attachment replaced.
The treated areas become red after the treatment.	Slight redness is harmless and normal and disappears quickly.	No action required.
The skin reaction after the treatment lasts longer than usual.	You have used a light intensity that is too high for you.	Select a lower light intensity the next time. See chapter 'Preparing for use', section 'Light intensity'. If the skin reaction lasts longer than 3 days, contact your doctor.
The hair-reduction results are not satisfactory.	You have used a light intensity that is too low for you.	Perform a skin test to determine if you can use a higher light intensity (see chapter 'Using the appliance', section 'First use and skin test').

Problem	Possible cause	Solution
	You did not make enough overlaps when you used the appliance.	For proper hair-reduction results, you have to make overlaps when you use the appliance. See chapter 'Using the appliance', section 'Subsequent use', step 4.
	You do not use the appliance as often as recommended.	To remove all hairs successfully and to prevent the hair root from becoming active again, the treatments have to be repeated once every two weeks for the first two months. To ensure that your skin stays smooth, we advise you to repeat the treatment every four to six weeks. The time between treatments may vary based on your individual hair regrowth and also across different body areas. If too many hairs grow back between treatments, simply shorten the treatment interval again, but do not treat any area more often than once every two weeks. When you use the appliance more frequently than recommended, you do not increase the effectiveness but you do increase the risk of skin reactions.
	You respond more slowly to IPL treatment than the average user.	Continue using the device for at least 6 months, as hair regrowth can still decrease over the course of this period. Also take into account the other solutions to this problem.
	The appliance is not suitable for your hair or skin colour.	Do not use the appliance if your body hair is light blond, red, grey or white. Do not use the appliance either if you have dark skin. Consult the body hair and skin colour table in the Quick Start Guide.
Hair starts to grow again in some spots on the treated areas.	You did not make enough overlaps when you used the appliance.	For proper hair-reduction results, you have to make overlaps when you use the appliance. See chapter 'Using the appliance', section 'Subsequent use', step 4.
	Hair starts to grow again as part of the natural cycle of hair growth.	This is normal. It is part of the natural cycle of hair growth that hair starts growing again. The hairs that come back, however, tend to be softer and thinner.
The appliance does not work.	The appliance is defective.	Contact the Consumer Care Centre in your country, your Philips dealer or a Philips service centre.

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