

PHILIPS

Lumea

SC2005

English	5
Czech	20
Polski	35
Romna	52
Slovenščina	68
繁體中文	83

Welcome

Welcome to the beauty world of Lumea! You are only a few weeks away from silky-smooth skin.

Philips Lumea uses Intense Pulsed Light (IPL) technology, known as one of the most effective methods to continuously prevent hair regrowth. In close cooperation with licensed dermatologists we adapted this light-based technology, originally used in professional beauty salons, for easy and effective use in the safety of your home. Philips Lumea is gentle and offers convenient and effective treatment at a light intensity that you find comfortable. Unwanted hairs are finally a thing of the past. Enjoy the feeling of being hair-free and look and feel amazing every day.

To fully benefit from the support that Philips offers, register your product at www.philips.com/welcome. For further information, please go to www.philips.com/lumea to find our experts' advice, tutorial videos and FAQs and make the most of your Lumea.

Device overview

- 1 Light exit window with integrated UV filter
- 2 Body attachment
- 3 Reflector inside the attachment
- 4 Safety ring
- 5 Flash button
- 6 Intensity lights (1-5)
- 7 Intensity increase button
- 8 Intensity decrease button
- 9 Charging light and battery low indication
- 10 On/off button
- 11 Device socket
- 12 'Ready to flash' light
- 13 Air vents
- 14 Adapter
- 15 Small plug

Who should not use Lumea? Contraindications

- Never use the device if you have skin type VI (You rarely to never sunburn, very dark tanning). In this case you run a high risk of developing skin reactions, such as hyperpigmentation and hypopigmentation, strong redness or burns.

- Never use the device if you are pregnant or breast feeding as the device was not tested on pregnant or breast feeding women.
- Never use the device if you suffer from any of the diseases listed below:
 - If you have a skin disease such as active skin cancer, you have a history of skin cancer or any other localized cancer in the areas to be treated.

- If you have pre-cancerous lesions or multiple atypical moles in the areas to be treated.
- If you have a history of collagen disorder, including a history of keloid scar formation or a history of poor wound healing.
- If you have a history of vascular disorder, such as the presence of varicose veins or vascular ectasia in the areas to be treated.
- If your skin is sensitive to light and easily develops a rash or an allergic reaction.
- If you have infections, eczema, burns, inflammation of hair follicles, open lacerations, abrasions, herpes simplex, wounds or lesions and haematomas in the areas to be treated.
- If you have had surgery in the areas to be treated in the last three weeks.
- If you have epilepsy with flashlight sensitivity.
- If you have diabetes, lupus erythematoses, porphyria or congestive heart disease.
- If you have any bleeding disorder.
- If you have a history of immunosuppressive disease (including HIV infection or AIDS)
- Never use the device if you take any of the medications listed below:
 - If your skin is currently being treated with or has recently been treated in the past week with Alpha-Hydroxy Acids (AHAs), Beta-Hydroxy Acids (BHAs), topical isotretinoin and azelaic acid.
 - If you have taken any form of isotretinoin Accutane or Roaccutane in the last six months. This treatment can make skin more susceptible to tears, wounds and irritations.
 - If you are on painkillers which reduce the skin's sensitivity to heat.
 - If you are taking photosensitising agents or medications, check the package insert of your medicine and never use the device if it is stated that it can cause photo-allergic reactions, photo-toxic reactions or if you have to avoid sun when taking this medicine.
 - If you take anticoagulation medications, including heavy use of aspirin, in a manner which does not allow for a minimum 1-week washout period prior to each treatment.
- If you take immunosuppressive medications.
- Never use the device on the following areas:
 - On the face.
 - On nipples, areolas, labia minora, vagina, anus and the inside of the nostrils and ears.
 - Men must not use it on the scrotum and face.
 - Over or near anything artificial like silicone implants, pacemakers, subcutaneous injection ports (insulin dispenser) or piercings.
 - On moles, freckles, large veins, darker pigmented areas, scars, skin anomalies without consulting your doctor. This can result in a burn and a change in skin color, which makes it potentially harder to identify skin-related diseases.
 - On warts, tattoos or permanent make-up. This can result in a burn and a change in skin color.

- On areas, where you use long-lasting deodorants. This can result in skin reactions.
- Never use the device on sunburned, recently tanned or fake-tanned skin.

Note: This list is not exhaustive. If you are not sure whether you can use the device, we advise you to consult a doctor.

Important

Danger

- Keep the device and the adapter dry.
- If the device is broken, do not touch any inner part to avoid electric shock.
- Never insert any objects into the device.

Warning

- This device is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, unless they have been given supervision or instruction concerning use of the device by a person responsible for their safety.
- Children should be supervised to ensure that they do not play with the device.
- The device is not intended for children under the age of 15 years. Teenagers aged between 15 and 18 years can use the device with the consent and/or assistance of their parents or the persons who have parental authority over them. Adults of 18 years and older can use the device freely.
- Only charge the device with the adapter supplied.
- Always check the device before you use it. Do not use the device or adapter if it is damaged. Always replace a damaged part with one of the original type.
- The adapter contains a transformer. Do not cut off the adapter to replace it with another plug, as this causes a hazardous situation.
- Do not use the device if the UV filter of the light exit window is broken.

To prevent damage

- Make sure nothing obstructs the airflow through the ventilation slots.
- Never subject the device to heavy shocks and do not shake or drop it.
- If you take the device from a very cold environment to a very warm environment and vice versa, wait approximately 3 hours before you use it.

- Do not flash against any other surface than the skin. This can cause severe damage to the attachment and/or the light exit window. Only flash when the device is in contact with skin.

Caution


- This device is not washable. Never immerse the device in water and do not rinse it under the tap.
- For hygienic reasons, the device should only be used by one person.
- Use the device only at settings suitable for your skin type. Use at higher settings than those recommended can increase the risk of skin reactions and side effects.
- Only use the device for its intended purpose as shown in the user manual.
- Never use compressed air, scouring pads, abrasive cleaning agents or aggressive liquids such as petrol or acetone to clean the device.
- Always return the device to a service centre authorised by Philips for examination or repair. Repair by unqualified people could cause an extremely hazardous situation for the user.


Electromagnetic fields (EMF)

This Philips device complies with all applicable standards and regulations regarding exposure to electromagnetic fields.


How IPL works

With IPL technology, gentle pulses of light are applied to the skin and absorbed by the hair root. The darker the hair, the better the pulses of light are absorbed.


The pulses of light stimulate the hair follicle to go into a resting phase. As a consequence, the hair sheds naturally and hair regrowth is prevented, leaving your skin continuously silky-smooth.

The cycle of hair growth consists of different phases. IPL technology is only effective when the hair is in its growing phase. Not all hairs are in the growing phase at the same time. This is why we recommend an 8-week treatment schedule followed by touch-ups to make sure all hairs are effectively treated in the growing phase.

Note: Treatment with Lumea is not effective if you have light blond, blond as light hairs do not absorb enough light. Below you can see the hair colors for which Lumea is suitable and effective.

Suitable body hair colors

What to expect

Immediately after the 1st treatment

After the first treatment, it can take 1 to 2 weeks for the hairs to fall out. In the first weeks following the initial treatments, you still see some hairs growing. These are likely to be hairs that were not in their growing phase during the first treatments.

After 2-3 treatments

After 2-3 treatments, you should see a noticeable reduction in hair growth. However, to effectively treat all hairs, it is important to keep on treating according to the recommended treatment schedule.

After 4-5 treatments

After 4-5 treatments, you should see a significant reduction of hair growth in the areas that you treated with Lumea. A reduction of hair density should be visible as well. Keep on treating with frequent touch-ups (every 4-8 weeks) to maintain the result.

Tanning advice

Tanning with natural or artificial sunlight

Intentionally exposing your skin to natural or artificial sunlight, with the aim of developing a tan influences the sensitivity and color of your skin.

Therefore the following is important:

- After each treatment, wait at least 24 hours before tanning. Even after 24 hours, make sure that the treated skin does not show any redness from the treatment anymore. Use a sunblock SPF 50+ on the treated areas in the 48 hours after treatment when exposing your skin to the sun. After this period, you can use your regular sunblock.
- After tanning, wait at least 48 hours before you use Lumea. Perform a skin test 48 hours after tanning to determine the appropriate light intensity setting.
- Do not use Lumea on sunburned body areas.


Note: Occasional and indirect sun exposure does not qualify as tanning.

Tanning with creams

If you have used an artificial tanning lotion, wait until the artificial tan has disappeared completely before you use the device.

Before you use your Lumea

Pretreating your skin

Before you use Lumea, you should pretreat your skin by removing hairs on the surface of your skin. This allows the light to be absorbed by the hair parts below the skin surface to ensure effective treatment. You can either shave, short-trim, epilate or wax. Do not use depilatory creams, as chemicals may cause skin reactions.

If you chose to wax, please wait 24 hours before using Lumea to let your skin rest. We recommend that you take a shower before the treatment to ensure that all possible residue of wax has been removed from your skin.

- 1 Pretreat the areas you intend to treat with Lumea.
- 2 Clean your skin and make sure it is hair-free, entirely dry and free from oily substances.

Note: Once the hairs stop growing back, which usually happens after 4-5 treatments, you no longer have to pretreat your skin before you use the device.


Skin test

When you use Lumea for the first time or after recent tanning, perform a skin test on each area to be treated. The skin test is necessary to check your skin's reaction to the treatment and to determine the correct light intensity setting for each body area.


- 1 Choose an area close to the area you intend to treat.
- 2 Turn on the device. Make sure you select setting 1.
- 3 Put the device against your skin and press the flash button to release a flash.
- 4 Slide the device over the skin to the next area to be treated.
- 5 Increase the setting by one level, apply a flash and slide the device to the next area. Repeat this for all levels within the recommended range for your skin type.
- 6 After the skin test, wait 24 hours and check your skin for any reaction. If your skin shows reactions, choose the highest setting that did not result in any skin reaction for subsequent use.

Using your Philips Lumea

Placing and removing

To place the attachment, simply snap it onto the light exit window.


To remove the attachment, pull it off the light exit window.

Selecting the right light intensity

Lumea provides 5 different light intensities. Depending on your skin type and light intensity level you find comfortable, you can select the right light intensity setting.

- 1 Consult the table below to select the right setting.

12 English


- 2 To adjust the light intensity setting, press the intensity increase button or the intensity decrease button until you have reached the required setting..
- 3 Lumea should never be painful. If you experience discomfort, reduce the light intensity setting.
- 4 After recent tanning, perform a skin test to determine the right light intensity setting.

Skin type	Skin tone	Light intensity setting
I	White; You always sunburn, never tan.	4/5
II	Beige; You easily sunburn, tan minimally.	4/5
III	Light brown; You sunburn easily, tan slowly to light brown.	4/5
IV	Mid brown; You rarely sunburn, tan easily.	3/4
V	Dark brown; You rarely sunburn, tan very easily.	1/2/3
VI	Brownish black or darker; You rarely or never sunburn, very dark tanning.	You cannot use the device


Note: Your skin may react differently on different days/occasions for a number of reasons.

Handling the device

- 1 Before use, charge the device, clean the attachments and the light exit window and select the right light intensity for your skin tone.
- 2 Place the device at a 90° angle on the skin so that the attachment and the safety ring are in contact with your skin.
The safety ring prevents unintentional flashing without skin contact.

- 3 Press the device firmly onto your skin to ensure proper skin contact.

The 'ready to flash' light on the back of the device lights up green to indicate that you can proceed with the treatment.

- 4 Press the flash button to release a flash. You should feel a warm sensation because of the flash.

Note: The light produced by the device is harmless to your eyes. It is not necessary to wear goggles during use. Use the device in a well-lit room so that the light is less glaring to your eyes.

- 5 Place it on the next area to be treated. After each flash, it takes up to 3.5 seconds until the device is ready to flash again. You can release a flash when the 'ready to flash' light lights up again.

- 6 To make sure that you have treated all areas, release the flashes close to each other. However, make sure you flash the same area only once. Flashing the same area twice doesn't improve the effectiveness of the treatment, but increases the risk of skin reactions

Two treatment modes: Stamp & Flash and Slide & Flash


Your Philips Lumea has two treatment modes for more convenient use on different body areas:

- The Stamp & Flash mode is ideal to treat small or curvy areas like knees and underarms. Simply press and release the flash button to release a single flash.
- The Slide & Flash mode offers convenient use on larger areas like legs. Keep the flash button pressed while you slide the device over your skin to release several flashes in a row.

Recommended treatment schedule

Initial phase

For the first 4 to 5 treatments, we advise you to use Philips Lumea once every two weeks to ensure that all hairs are treated.


Touch-up phase

After the initial treatment phase (4-5 treatments), we recommend touch-ups every 4 to 8 weeks, when you see hairs growing back. This is to maintain results and enjoy smooth skin for months. The time between treatments may vary based on your individual hair regrowth and also across different body areas.

Treatment time per area

After use

Common skin reactions


Your skin may show slight redness and/or may prickle, tingle or feel warm. This reaction is absolutely harmless and disappears quickly.

Dry skin and itching may occur due to shaving or a combination of shaving and light treatment. You can cool the area with an ice pack or a wet face cloth. If dryness persists, you can apply a non-scented moisturizer on the treated area 24 hours after the treatment.

Aftercare

After use, you can safely apply lotions, creams, deodorant, moisturizer or cosmetics to the treated areas. If you experience skin irritation or skin redness after treatment, wait until it disappears before applying any product to your skin. If you experience skin irritation after applying a product to your skin, wash it off with water.

Charging


Charge the batteries of the device before you use it for the first time. Charging takes approx. 1 hour and 40 minutes. Fully charge the batteries every 3 to 4 months, even if you do not use the device for a longer time. When the device is charging, the charging light flashes green.

Note: This device can only be used without cord.

- 1 Make sure the device is switched off.
- 2 Insert the small plug into the device and put the adapter in the wall socket.


Note: The device is equipped with battery-overheat protection and does not charge if the room temperature exceeds 40 °C.

Battery fully charged

When the batteries are fully charged, the charging light lights up solid green. Fully charged batteries provide at least 370 flashes at light intensity 5.


Battery low

When the battery is almost empty, the charging light starts to flash orange.

Cleaning & storage

- 1 After use, switch off the device, unplug it and let it cool down.

- 2 Remove the attachment.
- 3 Moisten the soft cloth supplied with the device with a few drops of water and use it to clean the following parts:
 - the light exit window
 - the outside surface of the attachment
 - the reflector inside the attachment
- 4 Let all parts air dry thoroughly.
- 5 Store the device in a dust-free and dry place at a temperature between 0 °C and 60 °C.

Recycling

- This symbol means that this product shall not be disposed of with normal household waste (2012/19/EU).

- This symbol means that this product contains a built-in rechargeable battery which shall not be disposed of with normal household waste (2006/66/EC). We strongly advise you to take your product to an official collection point or a Philips service centre to have a professional remove the rechargeable battery.
- Follow your country's rules for the separate collection of electrical and electronic products and rechargeable batteries. Correct disposal helps prevent negative consequences for the environment and human health.

Removing the rechargeable battery

Warning: Only remove the rechargeable battery when you discard the appliance. Make sure the battery is completely empty when you remove it.

- 1 Check if there are screws in the back or front of the appliance. If so, remove them.
- 2 Remove the back and/or front panel of the appliance with a screwdriver. If necessary, also remove additional screws and/or parts until you see the printed circuit board with the rechargeable battery.
- 3 Remove the rechargeable battery.

Guarantee and support

If you need information or support, please visit www.philips.com/support or read the separate worldwide guarantee leaflet.

Technical specifications

Rated voltage	100V-240V
Rated frequency	50Hz-60Hz
Rated input	7,5W
Protection against electric shock	Class II
Protections rating	IP 30 (EN 60529)
Operating conditions	Temperature: +15 °C to 35° C
	Relative humidity: 5% to 95%
Lithium-ion battery	2 x 3.7 Volt 1500 mAh

Troubleshooting

This chapter summarizes the most common problems you could encounter with the device. If you are unable to solve the problem with the information below, visit www.philips.com/support for a list of frequently asked questions or contact the Consumer Care Center in your country.

Problem	Possible cause	Solution
The device suddenly switches off.	The rechargeable batteries are empty.	Charge the device. See Chapter Charging.
The device/adaptor becomes warm during charging	This is normal.	No action required.
I cannot use the device from the mains.	Due to safety reasons it is not possible to use the device from the mains.	Philips Lumea can only be used without cord.
The 'ready to flash light' doesn't light up when I place the attachment on my skin, but the cooling fan works.	The overheat protection has been activated.	When the overheat protection has been activated, the fan still works. Do not switch off the device, but let it cool down for approx.15 minutes before you continue to use it.
The 'ready to flash light' lights up green but the device does not produce a flash when I press the flash button.	The device needs to be reset.	To reset the device, switch it off and switch it on again.
The device produces a strange smell.	The attachment or the filter glass is dirty.	Clean the attachment carefully.
	You have not removed the hairs on the area to be treated properly. These hairs may get burned and can cause the smell.	Pretreat your skin before you use Lumea.

Problem	Possible cause	Solution
The skin feels more sensitive than usual during treatment/ I experience an unacceptable pain sensation when I use the device.	The light intensity setting you use is too high.	Check if you have selected the right light intensity setting. If necessary, select a lower setting.
	You did not remove the hairs on the areas to be treated.	Pretreat your skin before you use Lumea.
	The UV filter of the light exit window is broken.	If the UV filter is broken, do not use the device anymore. Contact the Consumer Care Center in your country, your Philips dealer or a Philips service center.
	You treated an area for which the device is not intended.	Never use the device on the following areas; face, inner labia, anus, nipples, areolas, lips, inside the nostrils and ears, around the eyes and near the eyebrows. Men must not use if on face or on the scrotum.
The skin reaction after treatment lasts longer than usual.	The device is not suitable for your skin tone.	Do not use the device if you have skin type VI, a brownish black or darker skin tone.
	You have used a light intensity setting which is too high for you.	Select a lower intensity next time. See Chapter Using your Philips Lumea, section Selecting the right light intensity.
The results of the treatment are not satisfactory.	You have used a light intensity setting which is too low for you.	Select a higher setting next time.
	You did not flash an area adjacent to an area you treated before.	You have to release the flashes close to each other.
	The device is not effective on your body hair color.	If you have light blond, blond, or red hair, the treatment is not effective.
My body attachment does not have an integrated glass/filter/window.	This is normal.	The body attachment does not contain a glass/filter/window.

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