

PHILIPS

IPL – Hair removal device

Lumea Advanced

For use on body

15 minutes to treat lower legs

>250,000 lamp flashes

Extra-long cord



SC1993/00



Gentle, long-term hair reduction

Philips Lumea Advanced works wonders to prevent the reappearance of body hair. When applied regularly, the gentle pulses of light keep skin silky-smooth every day.



Long-lasting hair removal with IPL technology

- Enjoy smooth skin every day
- Effortlessly effective
- Adapted for safe and effective use in your own home
- Suitable for a wide variety of hair and skin types

Safe, fast and effective on your body

- For use on the body
- Large treatment window for fast application
- Gentle treatment

Maximum freedom, minimum hassle

- No replacement parts, no hidden costs
- Five adjustable light energy settings



asimpleswitch.com


Highlights

IPL technology



Philips Lumea uses an innovative light-based technology called IPL (Intense Pulsed Light), derived from the technology used in professional beauty salons. Philips has adapted this technology for safe and effective use in the comfort of your own home. Philips worked closely with leading dermatologists to develop our breakthrough hair removal system. For more than 10 years we carried out extensive consumer research with more than 2000 volunteers.

Enjoy smooth skin every day



Philips Lumea applies gentle pulses of light to the hair root. As a consequence, the amount of hair your body grows gradually decreases. Repeating the treatment leaves your skin beautifully hair free and touchably smooth.

Effortlessly effective

Clinical results after treatment: **example 1 (leg)**



Clinical results after treatment: **example 2 (bikini line)**



Our clinical studies have shown significant hair reduction in as little as four fortnightly treatments resulting in smooth hair free skin. To maintain these results, simply repeat the treatment when needed. Time between treatments may vary based on your individual hair regrowth. - To enlarge this picture, please click on the image in the image gallery on the top of this page

Suitable hair and skin types




Philips Lumea works effectively on (naturally) dark blonde, brown and black hairs. As with other light-based treatments, Philips Lumea is NOT effective on red, light-blond or white/grey hair. This product is also NOT suitable for dark skin.

Application areas



Use with confidence to remove hair on body areas below the neckline: legs, underarms, bikini area, arms and stomach.

Fast treatment on legs



Big treatment window for fast application on large areas like legs.

Gentle treatment




Gentle on your skin.

No replacement parts needed



Philips Lumea is a complete solution that comes ready to work straight out of the box. More importantly, it does not require any replacement lamps or gels.

Five light energy settings



Philips Lumea has five adjustable light energy settings to ensure gentle but effective treatment. Used correctly, Philips Lumea's IPL treatment is safe and gentle to use even on sensitive skin and sensitive body areas.



[asimpleswitch.com](https://www.asimpleswitch.com)

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Technical specifications

High performance lamp: Delivers > 250 000 flashes

Application areas

Body areas: Arms, Legs, Belly, Bikini, Underarms

Attachments

Body attachment (4 cm²): For use below the neckline

Items included

Instructions for use: User manual

Storage: Pouch

Adapter: 24 V/1500 mA

Accessory: Cleaning cloth

Power

Voltage: 100–240 V

Safety and adjustable settings

Integrated UV filter: Protects skin from UV light
5 light energy settings: Adjustable to your skin type

Integrated safety system: Prevent unintentional flashing

Service

Warranty: 2 years global warranty

Application time

Lower leg: 8 min

Underarm: 1 min

Bikini line: 1 min

Application mode

Corded / cordless use: Corded

Stamp and Flash: For treatment on curvy areas

