

# PHILIPS

IPL hair removal  
device

Lumea Prestige

With SenseIQ technology

Underarms, bikini, body, face

With SmartSkin sensor

Both cordless and corded use



BRI956/60



## Enjoy 6 months of hair-free smooth skin\* 4 intelligent attachments – for optimal results

Philips Lumea Prestige with SenseIQ technology is our most effective IPL. Designed for convenient treatment at home, Lumea's intelligent attachments perfectly fit every curve of your body and adapt programs tailored for each body area.

### Long-lasting smooth skin with Lumea IPL technology

- Expert IPL technology at-home, developed with dermatologists
- Proven gentle and effective treatment
- Suitable for a wide range of hair and skin types
- SenseIQ technology for personalized hair removal
- Use cordless or corded

### Professional IPL treatment for your body's areas

- Intelligent attachments adapt programs for each body area
- Underarm attachment with curved-out design
- Bikini area attachment with extra light filter
- Body attachment with curved-in design
- Precision face attachment with extra filter

### Maximum freedom, minimum hassle


- With SmartSkin sensor
- The Lumea IPL app. With you every step of the way



[asimpleswitch.com](https://asimpleswitch.com)


# Highlights

## Lumea prevents hair regrowth\*




IPL stands for Intense Pulsed Light. Philips Lumea applies gentle pulses of light to the hair root, putting the follicle into a resting phase. As a consequence, the amount of hair your body grows gradually decreases. Repeating the treatment leaves your skin beautifully hair free and touchably smooth. Treatment to prevent hair regrowth is safe and gentle, even on sensitive areas. Philips Lumea is clinically tested and developed with dermatologists for easy and effective treatment, from the comfort of your own home.

## Proven effectiveness



Objective studies show up to 92% hair reduction after three treatments\*\*. Carry out the first four treatments every two weeks, and the next eight treatments every four weeks. After 12 treatments you can enjoy six months of hair-free smooth skin\*.

## Suitable hair and skin types



Philips Lumea Prestige is effective on a wide range of hair and skin types. It works on naturally dark blonde, brown and black hair and on skin tones from very white to dark

brown. IPL needs contrast between the pigment in the hair color and the pigment in the skin tone, therefore (as with other IPL-based treatments) Lumea cannot be used to treat white, grey, light blonde or red hair. And is not suitable for very dark skin.

## SenseIQ technology



Senses. Lumea's SmartSkin sensor detects your skin tone and indicates the most comfortable light setting for you. Adapts. Intelligent attachments adapt programs tailored for each body area. Cares. Co-developed with dermatologists and inspired by you for gentle, long-lasting hair removal, even on sensitive areas\*\*\*.

## Treatment of body areas



Your body is unique. And its skin, curves and contours should be treated uniquely too. Lumea's total-body attachments are shaped to perfectly fit every curve, and will adapt programs tailored to each body area. Treatment windows are designed to maximize contact with your skin to ensure light doesn't escape, enabling effective and gentle treatments even on delicate or sensitive areas.

## Underarm attachment



Specially designed to treat hard-to-reach underarm areas, the intelligent underarm attachment is medium-sized with a 3 cm<sup>2</sup> window. When the attachment is clicked on to the device it adapts a tailored program for your underarm area. 86% of women satisfied with hair reduction on underarms\*\*\*\*.

## Bikini attachment



Click on the intelligent bikini attachment and it will adapt a tailored program specifically for that area. It has an optimal 3 cm<sup>2</sup> window with a transparent filter for extra comfort and the curved design compliments your body's contours for maximum contact, efficiency and comfort. 78% of women satisfied with hair reduction in bikini area\*\*\*\*.



[asimpleswitch.com](https://www.asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO<sub>2</sub> emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Application areas

**Body areas:** Arms, Legs, Stomach  
**Face areas:** Chin, Upper lip, Jawline  
**Sensitive areas:** Armpits, Bikini

## Technical specs. attachments

**Body treatment attachment:** Shape: Convex curved, Window size: 4.1 cm<sup>2</sup>, Tailored treatment for body  
**Face treatment attachment:** Shape: Flat, Window size: 2 cm<sup>2</sup>, extra filter, Tailored treatment for face  
**Armpits treatment attachment:** Shape: Concave curved, Window size: 3 cm<sup>2</sup>, Tailored treatment for armpits  
**Bikini treatment attachment:** Shape: Concave curved, Window size: 3 cm<sup>2</sup>, extra filter, Tailored treatment for bikini

## Application mode

**Slide and Flash:** For quick application  
**Stamp and Flash:** For treatment on small areas  
**Corded / cordless use:** Both cordless and corded use

## Safety and adjustable settings

**Integrated UV filter:** Protects skin from UV light  
**5 light energy settings:** Adjustable to your skin type  
**Integrated safety system:** Prevent unintentional flashing  
**Skin Tone sensor:** Detects your skin tone  
**SmartSkin sensor:** The right setting on demand

## Application time

**Lower legs:** 8,5 min  
**Armpit:** 2,5 min  
**Bikini line:** 2 min  
**Face areas:** 1,5 min

## Technical specifications

**High performance lamp:** Delivers >250.000 flashes\*\*\*\*\*

## Power

**Voltage:** 100-240 V

## Items included

**Instructions for use:** User manual  
**Storage:** Luxurious pouch  
**Adapter:** 19.5V / 4000mA

## Service

**Warranty:** 2 years global warranty

