PHILIPS

IPL - Hair removal device

Lumea Prestige

2 intelligent attachments Body and face With SmartSkin sensor Both cordless and corded use





BRI950/00

Enjoy 6 months* of hair-free smoothness Intelligent attachments designed for best results

Philips Lumea Prestige, our most powerful IPL yet, designed for your body's curves and the easiest at-home experience. Uniquely curved, intelligent attachments perfectly fit your curves and adapt the treatment programmes for each body area.

Continuously smooth skin with Lumea IPL technology

- Expert IPL technology at home, developed with dermatologists
- Proven safe and effective treatment
- · Suitable for a wide variety of hair and skin types
- Use either cordless or corded

Specialised treatment of different body areas

- Intelligent attachments adapt the treatment of each area
- Precise face attachment with extra filter
- Large body attachment with curved-in design

Maximum freedom, minimum hassle

- With SmartSkin sensor
- FREE Lumea App for a personal coach at your fingertips
- · No replacement parts, no hidden costs



Highlights

Lumea prevents hair regrowth*



IPL stands for Intense Pulsed Light. Philips Lumea applies gentle pulses of light to the hair root, putting the follicle into a resting phase. As a consequence, the amount of hair your body grows gradually decreases. Repeating the treatment leaves your skin beautifully hair free and touchably smooth. Treatment to prevent hair regrowth is safe and gentle, even on sensitive areas. Philips Lumea is clinically tested and developed with dermatologists for easy and effective treatment, from the comfort of your own home.

Proven effectiveness

Objective studies show up to 92% hair reduction in as little as 3 treatments^{**}. The first 4 treatments must be carried out every 2 weeks, after which you should already be able to see the desired results. To maintain your results, simply touch up every 4 weeks. After just 8 touch-up treatments you can be hair-free for 6 months* To enlarge this picture, please click on the image in the image gallery at the top of this page.

Suitable hair and skin types



Philips Lumea Prestige works effectively, easily and safely on a wide variety of hair and skin types. It works on naturally dark blonde, brown and black hair and on skin tones from very white to dark brown. As with other IPLbased treatments, Philips Lumea cannot be used to treat white/grey, light blonde or red hair and is not suitable for very dark skin. This is due to the high contrast required between the pigment in the hair colour and the pigment in the skin tone.

Treatment of body areas



Your hair and your curves vary across your body, so each area needs treating in a different way to ensure the best results. Philips Lumea Prestige has uniquely curved intelligent attachments for the body (legs, arms, stomach) and face (upper lip, chin). Each intelligent attachment differs by combinations of shape, window size and filter, and triggers a tailored treatment programme for specific body part when added to the device. This mix of considered design and bespoke treatment enables effective and safe treatment even on delicate or sensitive areas.

Face attachment



To gently but effectively treat hair on your face, the attachment has been carefully designed with a small, 2 cm2 window that's also flat for precise use on the upper lip, chin or jawline. It has an extra integrated filter and when added to the device the intelligent attachment will automatically adjust the light treatment for your face. 84% of women satisfied with hair reduction on face***.

Body attachment



The body attachment has our largest treatment window (4.1 cm2) so it's perfect for a fast treatment on large body areas like legs, arms and stomach. The curved-in design closely follows body contours for optimal skin contact. When added to the device this intelligent attachment automatically triggers a tailored treatment suitable for these body areas. 83% of women satisfied with hair reduction on legs***.

SmartSkin sensor

Philips Lumea Prestige has 5 energy settings and makes it easy for you to choose the right one thanks to the SmartSkin sensor. The SmartSkin sensor advises the best setting for your skin tone, any time you want.

Both cordless and corded use

Philips Lumea Prestige can be used both cordless and corded. For big areas like legs, we recommend using the device corded for a faster treatment without interruption. When treating smaller or hard to reach areas like the face, underarms or bikini area, use the device cordless so you can enjoy full freedom of movement.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Application areas

Body areas: Arms, Legs, Stomach Face areas: Chin, Upper lip, Jawline Sensitive areas: Underarms, Bikini

Technical specs. attachments

Body treatment attachment: Window size: 4.1 cm2, Shape: Convex curved, Tailored treatment for body

Face treatment attachment: Window size: 2 cm2, extra filter, Shape: Flat, Tailored treatment for face

Application mode

Slide and Flash: For quick application Stamp and Flash: For treatment on small areas Corded / cordless use: Both cordless and corded use

Safety and adjustable settings

5 light energy settings: Adjustable to your skin type

SmartSkin sensor: The right setting on demand Integrated safety system: Prevent unintentional flashing Integrated UV filter: Protects skin from UV light Skin Tone sensor: Detects your skin tone

Application time

Underarm: 2.5 min Face areas: 1.5 min Bikini line: 2 min Lower legs: 8.5 min High performance lamp: Delivers > 250 000 flashes

Power Voltage: 100-240 V

Items included

Adapter: 19.5 V / 4000 mA Instructions for use: User manual Storage: Luxurious pouch

Service

Warranty: 2 years global warranty

Specifications are subject to change

without notice. Trademarks are the

property of Koninklijke Philips N.V. or

their respective owners.

Issue date 2019-11-26 Version: 7.3.1

EAN: 08 71010 37790 56

www.philips.com



- objective study on 48 women in the Netherlands and Austria
- * * Up to 92% hair reduction after 3 treatments when following the treatment schedule, measured on legs; individual results may vary
- * ** Study conducted in the Netherlands and Austria, 46 women, after 3 treatments on legs, after 2 treatments on face

^{*} Median result of legs: 78%, after 12 treatments,