
PHILIPS

Lumea
IPL hair removal system

Quick Start Guide



ENGLISH 6

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General description (Fig. 3)

- A Light exit window with integrated UV filter
- B Body attachment (suitable for all body areas below the neck)
- C Metallic frame inside the attachment
- D Safety system (safety ring with contact switches)
- E Flash button
- F Intensity lights (1- 5)
- G ▲ Intensity increase button
- H ▼ Intensity decrease button
- I Charging light and battery low indication
- J On/off button
- K 'Ready to flash' light
- L Appliance socket
- M Air vents

Introduction

Congratulations on your purchase and welcome to Philips! Light-based hair removal is one of the most effective methods to treat unwanted body hair with longer-lasting results. You can now use this innovative light-based hair regrowth prevention technology comfortably in the privacy of your own home, whenever and wherever you want. Philips Lumea enables treatment of unwanted hair on the whole body below the neck using the tailored attachment.

Philips Lumea hair removal system applies gentle pulses of light to the hair and hair root beneath the skin. As a consequence, the hair sheds naturally and hair regrowth is inhibited. For more information about Philips Lumea IPL hair removal system and movies how to use Lumea, see www.philips.com/lumea. You can also register your product on www.philips.com/welcome.

Before using

1 Be aware!

This appliance is not suitable for everyone. This Quick Start Guide does not substitute the user manual. To ensure that Philips Lumea hair removal system is suitable for you and safe for you to use, ***read the user manual carefully before first use.***

Before you use Philips Lumea, perform a skin test: Choose a hair-free area close to the area to be treated. Apply one flash (see steps 3-7 in chapter 'Using') at the lowest recommended intensity for your skin colour. Increase the setting by one level within the recommended range of settings. After the skin test, wait for 24 hours and check your skin for any reactions. If you notice any skin reactions, choose the highest setting that did not result in any skin reactions for subsequent use.

Note: All types of tanning may influence the safety of the treatment. Wait at least 48 hours after sunbathing before you use the appliance. Do not use Philips Lumea on sunburnt body areas. See section "Tanning advice" in the user manual for more information.

2 **Recommended light intensities table**

Consult the table on the following page to determine which light intensities are most suitable for your skin and body hair colour and check if this method is suitable for you.

The method is more effective at higher intensities; however, you should reduce the light intensity if you experience any pain or discomfort. The light pulse can feel warm or hot on the skin but should never be painful.

3 **Check “contraindications” in the user manual; some medications and health conditions might cause adverse reactions when you use the appliance.**


Recommended light intensity table (1-5)

Skin tone \ Hair color	White	Beige	Light Brown	Medium Brown	Dark Brown	Brownish black and darker
White/Grey	⊗	⊗	⊗	⊗	⊗	⊗
Red	⊗	⊗	⊗	⊗	⊗	⊗
Light blond	⊗	⊗	⊗	⊗	⊗	⊗
Dark blond	4/5	4/5	4/5	3/4	1/2/3	⊗
Brown	4/5	4/5	4/5	3/4	1/2/3	⊗
Dark brown	4/5	4/5	4/5	3/4	1/2/3	⊗
Black	4/5	4/5	4/5	3/4	1/2/3	⊗


In the table above, 'x' means that the appliance is not suitable for you.

Preparing for use

- 1 Charge the appliance.



- 2 Shave the areas you intend to treat.


Using

- 1 Check the attachment and light exit window for any dirt, hair or fibres. If there is any contamination, clean it with the cloth supplied.
- 2 Put on the attachment



3 Press the on/off button to switch on the appliance.

4 Press the intensity increase/decrease button up or down to change the intensity (from 1 to 5).

5 Place the appliance at a 90° angle on the skin so that the attachment and the safety ring are in contact with the skin.

6 Press the safety ring fully onto the skin with slight pressure. When all contact switches of the appliance are in contact with the skin, the “ready to flash” light on the back of the appliance lights up green.

- 7 Press the flash button to release a flash. Now you have 2 options to release the next flash:
- A Keep the flash button pressed. The appliance releases one flash after another as long as the safety ring is in full contact with the skin. Immediately after the appliance has flashed, slide the appliance to the next area. Make sure the safety ring is fully pressed onto the skin.


This 'Slide & Flash mode' is especially convenient for treatment of larger areas like the legs.

- B Let go of the flash button after each flash. For the next flash, place the appliance on the area right next to the previously flashed area. Make sure the safety ring is fully pressed onto the skin.


This 'Step & Flash mode' guarantees a precise treatment of e.g. the upper lip, knees or ankles.

Note: After a few seconds, the appliance is ready to flash again, so make sure that you move it to the next area immediately after it has flashed.

Do not treat exactly the same area multiple times during one session. This does not improve the effectiveness of the treatment, but increases the risk of skin reactions.




- 8 To avoid untreated areas, always make sure there is some overlap with the previously treated area when you place the appliance on the skin.
- 9 Switch off the appliance after you have finished the treatment.
- 10 If there is any dirt on the attachment and/or the light exit window, clean it before you store the appliance in the case. If necessary, moisten the supplied cloth with water or high percentage alcohol for cleaning. See chapter 'Cleaning and maintenance' in the user manual for more information.

After treatment

Wait at least 24 hours or until all redness has disappeared before you expose treated areas to the sun. Cover treated areas when you go out into the sun or use a sunblock (SPF 30+) in the two weeks after treatment.

How to achieve optimal results

- To remove all hairs successfully and prevent the hair follicle from becoming active again, use Philips Lumea every two weeks for the first 4 to 5 treatments (more frequent treatments do not lead to faster or better results).
- In between treatments, the hair sheds naturally and hair regrowth is inhibited. Your skin should be smooth by then. To ensure that your skin stays smooth, simply repeat treatment when needed. The time between treatments may vary based on your individual hair regrowth and also may vary from body area to body area. To maintain the best results, we recommend repeating the treatment every 4 to 8 weeks.
- For an effective treatment, shave the area to be treated before you use the appliance as long as you still experience regrowth, which is usually the case for the first few treatments.

Average treatment time:

Area	Approx. treatment time
Underarms	1 minute per underarm
Bikini line	1 minute per side
Entire bikini area	4 minutes
One lower leg	8 minutes
One full leg	14 minutes

Monitor your steps to smooth skin every day!

Treatment no.	Date	Treated body area	Next treatment date	Used setting	Note

Treatment no.	Date	Treated body area	Next treatment date	Used setting	Note

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